



INSTRUCTION SHEET

ADJUST THE NEWTON GRADE AID WHEEL LOCKS

i This document explains how to properly **adjust** the **position** of the **Newton 3 positions grade Aid Wheel Lock**.

Wheelchair model(s)


- All wheelchairs

Tool(s) and materials required

- Hex keys (Allen keys): 4 mm, 5 mm
- Medium strength threadlocker adhesive (blue Loctite)
- Torque wrench


OPERATION

- **Position 1 - Free:** The rear wheel is free to turn in both directions
- **Position 2 - Grade Aid:** The wheel can only be turned forward. Backward movement is blocked by wheel lock
- **Position 3 - Locked:** The wheel is locked in both directions





ADJUSTING A GRADE AID WHEEL LOCK

- **Manually release** the **wheel lock** by pushing the lever backward and **place** it so that the **short pin** of the mechanism engages in the **adjustment notch**



- Loosen the **two screws of the wheel lock fixing clamps**:
 - **Two pieces** clamp (photo): **5 mm** hex key
 - **Ultralight** one piece clamp: **4 mm** hex key



- **Slide** the wheel lock **stem forward or backward** in the **clamp** to get **about a 1 mm gap** between the **tire** and the **wheel lock** (see photo)
 - **TIP**: You can use a **credit card** to set the **gap**
- **Tighten** the **screws of the clamp without forcing them too much**
- Apply medium strength **threadlocker adhesive** (blue Loctite) to the end of the screw threads:
 - Loosen one screw and keep the other one in place
 - Apply the threadlocker adhesive and then tighten the screw
 - Do the same for the other screw
- Apply the proper **tightening torque** according to the size of the key:
 - **4 mm** key: **7 Nm**
 - **5 mm** key: **12 Nm**
- **Repeat the operations** for the **other wheel lock**
- **Test the correct operation** of the grade aid wheel locks

