

ORDER FORM

User manual: English French

Account : _____ PO : _____ Ship to : _____
 Vendor : _____ Name : _____
 Date : _____ Address : _____
 Ordered by : _____ City: _____ Province : _____
 Marked for : _____ Postal Code : _____ Phone : _____

Description

1	Ultralight manual folding helio wheelchair (1)		2850\$	<input type="checkbox"/>
FRAME COLOR				
2	Onyx black <input type="checkbox"/> Monza silver <input type="checkbox"/> Sapphire blue <input type="checkbox"/> Ruby red <input type="checkbox"/>		NC	<input type="checkbox"/>
3	Custom color		300\$	<input type="checkbox"/>
4	Carbon fiber finish		450\$	<input type="checkbox"/>
SEAT WIDTH & DEPTH				
5	Seat width 14" <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/>		NC	<input type="checkbox"/>
6	Seat width 19" <input type="checkbox"/> 20" <input type="checkbox"/>		285\$	<input type="checkbox"/>
7	Seat depth 13" (2) <input type="checkbox"/> 14" (3) <input type="checkbox"/> 15" <input type="checkbox"/> 16" <input type="checkbox"/> 17" <input type="checkbox"/> 18" <input type="checkbox"/>		NC	<input type="checkbox"/>
8	Seat depth 19" <input type="checkbox"/> 20" <input type="checkbox"/>		215\$	<input type="checkbox"/>
FRAME SIZE (See drawing at the end of order form)				
9	Regular frame		NC	<input type="checkbox"/>
10	Extended frame + 1" (4) <input type="checkbox"/> Extended frame + 2" (5) <input type="checkbox"/>		NC	<input type="checkbox"/>

notes :

(1) : Maximum weight limit: 250lbs

(2) : Available with Extended frame +2" only

(3) : Available with Extended frame +1" or +2" only

(4) : Not available with 20" seat depth

(5) : Not available with 19" and 20" seat depth

STRAIGHT BACK, HEIGHT ADJUSTABLE				
11	9" - 12" <input type="checkbox"/> 12" - 15" <input type="checkbox"/> 15" - 18" <input type="checkbox"/> 18" - 21" <input type="checkbox"/>		NC	<input type="checkbox"/>
Set height to : _____				

8° BEND BACK, HEIGHT ADJUSTABLE				
12	13½" - 15" <input type="checkbox"/> 15" - 18" <input type="checkbox"/> 18" - 21" <input type="checkbox"/>		96\$	<input type="checkbox"/>
Set height to : _____				

ADJUSTABLE ANGLE BACK, HEIGHT ADJUSTABLE				
Height				
13	12" - 15" <input type="checkbox"/> 15" - 18" <input type="checkbox"/> 18" - 21" <input type="checkbox"/>		236\$	<input type="checkbox"/>
Angle				
14	85° <input type="checkbox"/> 90° <input type="checkbox"/> 95° <input type="checkbox"/> 100° <input type="checkbox"/> 105° <input type="checkbox"/> 110° <input type="checkbox"/>			
Set height to : _____				

PUSH HANDLES				
15	Push handles		NC	<input type="checkbox"/>
16	Omit push handle		NC	<input type="checkbox"/>

17	FOLDING STABILIZER BAR* (Only available with push handles)		175\$	<input type="checkbox"/>
----	---	--	-------	--------------------------

*Suggested with 19" or 20" seat width

MOTIONCOMPOSITES

519, J-Oswald Forest, suite 101, Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0
 T.450 588-6555 1 866 650-6555 F.450 588-0200
 www.motioncomposites.com



December 22nd, 2009

Description

SEATING UPHOLSTERY				
18	Nylon seat sling		NC	<input type="checkbox"/>
19	Omit seat sling		NC	<input type="checkbox"/>
SEATING CUSHION				
20	2" cushion		115\$	<input type="checkbox"/>
21	3" cushion		135\$	<input type="checkbox"/>
22	Omit cushion		NC	<input type="checkbox"/>
BACK UPHOLSTERY				
23	Slip on back upholstery		NC	<input type="checkbox"/>
24	Adjustable tension back upholstery		162\$	<input type="checkbox"/>
25	Omit back upholstery		-15\$	<input type="checkbox"/>

FLIP BACK ARMREST				
26	Desk length		348\$	<input type="checkbox"/>
27	Full length		348\$	<input type="checkbox"/>
Height				
28	6" – 8" <input type="checkbox"/> 8" – 10" <input type="checkbox"/> 9½" – 12" <input type="checkbox"/> 12" – 14" <input type="checkbox"/>			

"T" ARMREST				
29	Desk length		280\$	<input type="checkbox"/>
30	Full length		280\$	<input type="checkbox"/>
Height				
31	6" – 8" <input type="checkbox"/> 8" – 10" <input type="checkbox"/> 9½" – 12" <input type="checkbox"/> 12" – 14" <input type="checkbox"/>			

SWING AWAY ARMREST Height adjustable 8" - 11 ½"				
32	Tubular arm pad		250\$	<input type="checkbox"/>
33	Desk length arm pad		250\$	<input type="checkbox"/>
34	Rigid side guard		190\$	<input type="checkbox"/>

FOOTREST				
35	60° swing away footrest, 14" – 18" length adjustable		55\$	<input type="checkbox"/>
36	70° swing away footrest, 14" – 18" length adjustable		NC	<input type="checkbox"/>
37	Length modification, less than 14" or over 18"		65\$	<input type="checkbox"/>
38	Direct mount (<i>Footplate(s) mounted to the front frame</i>)		35\$	<input type="checkbox"/>
Set length to : _____				

39	ELEVATING LEGREST 14" – 18" length adjustable		329\$	<input type="checkbox"/>
40	Length modification, less than 14" or over 18"		65\$	<input type="checkbox"/>
41	Calf pad		NC	<input type="checkbox"/>
Set length to: _____				

FOOTPLATES				
42	Standard (<i>Composite</i>)		NC	<input type="checkbox"/>
43	Adjustable angle oversize (<i>Aluminum</i>)		138\$	<input type="checkbox"/>
44	One piece adjustable angle flip-up		190\$	<input type="checkbox"/>
45	Length adjustable heel loop		NC	<input type="checkbox"/>
46	Length adjustable heel loop with ankle strap		76\$	<input type="checkbox"/>
47	Calf strap		36\$	<input type="checkbox"/>
48	Padded calf strap		58\$	<input type="checkbox"/>

MOTIONCOMPOSITES

519, J-Oswald Forest, suite 101, Saint-Roch-de-l'Acigan, QC, Canada J0K 3H0
 T.450 588-6555 1 866 650-6555 F.450 588-0200
 www.motioncomposites.com



December 22nd, 2009

Description

CASTERS					
49	Composite wheel /w PU tire (1" width)	3" □ 4" □		50\$	□
50		5" □ 6" □ 7" □ 8" □		NC	□
51	Composite wheel Soft-Roll® by Frog Legs	4" x 1 1/2"		85\$	□
52		5" x 1 1/2"		95\$	□
53		6" x 1 1/2"		105\$	□
54	Pneumatic	6" x 1 1/4" □ 8" x 1 1/4" □		120\$	□

FRONT SEAT-TO-FLOOR HEIGHT MATRIX

Circle the desired configuration

Forks		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
3"	P1	13 1/2	14 1/2	15 1/2															
	P2	14	15	16	14 1/2	15 1/2	16 1/2												
4"	P1	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2									
	P2																		
	P3	14 1/2	15 1/2	16 1/2	15 1/2	16 1/2	17 1/2	16 1/2	17 1/2	18 1/2	16 1/2	17 1/2	18 1/2	17 1/2	18 1/2	19 1/2	17 1/2	18 1/2	19 1/2
5"	P1																		
	P2	14 1/2	15 1/2	16 1/2	15 1/2	16 1/2	17 1/2	16 1/2	17 1/2	18 1/2	16 1/2	17 1/2	18 1/2	17 1/2	18 1/2	19 1/2	17 1/2	18 1/2	19 1/2
	P3																		
	P4	14 1/2	15 1/2	16 1/2	15 1/2	16 1/2	17 1/2	16 1/2	17 1/2	18 1/2	16 1/2	17 1/2	18 1/2	17 1/2	18 1/2	19 1/2	17 1/2	18 1/2	19 1/2
7"	P1										17 1/2	18 1/2	19 1/2	17 1/2	18 1/2	19 1/2			
	P2							17 1/2	18 1/2	19 1/2	17 1/2	18 1/2	19 1/2	18 1/2	19 1/2	20 1/2			
	P3							17 1/2	18 1/2	19 1/2	18 1/2	19 1/2	20 1/2	18 1/2	19 1/2	20 1/2	19 1/2	20 1/2	21 1/2
	P4				17 1/2	18 1/2	19 1/2	18 1/2	19 1/2	20 1/2	18 1/2	19 1/2	20 1/2	19 1/2	20 1/2	21 1/2	19 1/2	20 1/2	21 1/2
	P5	17 1/2	18 1/2	19 1/2	18 1/2	19 1/2	20 1/2	18 1/2	19 1/2	20 1/2	19 1/2	20 1/2	21 1/2	19 1/2	20 1/2	21 1/2	20 1/2	21 1/2	22 1/2


REAR WHEELS

55	Mag (Available with soft urethane tires and vinyl handrim only)	22" □ 24" □		125\$	□
56	Spoke	22" □ 24" □		NC	□
57	Spinergy Wire®	20" □ 22" □ 24" □ 26" □		275\$	□
58	Spinergy Spox®	22" □ 24" □ 25" □ 26" □		800\$	□
59	Spinergy LX®	22" □ 24" □ 25" □ 26" □		900\$	□
TIRES					
60	Soft urethane	20" □ 22" □ 24" □		NC	□
61	Pneumatic 1 3/8"	20" □ 22" □ 24" □ 26" □		NC	□
62	Pyramid – hard urethane	22" □ 24" □		70\$	□
63	Kenda – Ironcap high pressure	20" □ 22" □ 24" □ 25" □ 26" □		80\$	□
64	Soft urethane – high performance	22" □ 24" □ 25" □ 26" □		120\$	□

REAR SEAT-TO-FLOOR HEIGHT MATRIX


Circle the desired configuration

	20" wheel	22" wheel	24" wheel	26" wheel
1	12 3/8	13 3/8	14 3/8	15 3/8
7	12 3/4	13 3/4	14 3/4	15 3/4
2	13 1/8	14 1/8	15 1/8	16 1/8
8	13 5/8	14 5/8	15 5/8	16 5/8
3	14	15	16	17
9	14 1/2	15 1/2	16 1/2	17 1/2
4	14 7/8	15 7/8	16 7/8	17 7/8
10	15 1/4	16 1/4	17 1/4	18 1/4
5	15 3/8	16 3/4	17 3/4	18 3/4
11	16 1/8	17 1/8	18 1/8	19 1/8
6	16 5/8	17 5/8	18 5/8	19 5/8
12	17	18	19	20




FRAME SIZE AVAILABLE


Seat depth



Regular Frame



Extended Frame +1"



Extended Frame +2"

Extended Frame recommendation

- Allows for future growth
- Increased stability
- Use with seat depth < 17" to permit forward center of gravity adjustment