

New Receptacle Support

for Swing Away Armrest



Old model



New model

The new receptacle support for the Swing Away Armrest has been greatly improved.

The R&D department of Motion Composites has completely rethought the design of this piece in order to optimize its performance, increase its strength and reduce weight.

Following a detailed engineering study, the team was able to significantly reduce the weight by eliminating all non-useful material.

In doing so, the new profile design is formed as a double gripper, which increases the clamping force on the frame without increasing the tightening torque of the bolts.

Furthermore, a plastic insert was added to the interior of the container to ensure optimal performance of the Swing Away Armrests.



Compatibilité :

Codes :

33538 Réceptacle en L