

Double hook latch system

for Flip-Up armrest

Codes:

Update kit (2 latches, screws and springs): **31552**

Compleat kit: **100613**

The new double latch system allows two options for use of your armrests:
Removable or **permanent**.

Removable:



By installing the latch in «Removable» position (diamond hole side), the armrest can be removed from the latch.

Permanent:



By installing the latch in «permanent» position (diamond screw side), the armrest will not detach from the latch.

Old version
Black Single Hook



New version Red Double Hook

Tip to remove armrest in permanent configuration.

