



Stroller Handle

New design

The R & D team has redesigned the stroller handle by changing the angle to maximize leverage making it easier to manipulate the wheelchair.


This new handle consists of a tubular aluminum frame and foam padding.

To order:
Note that the appropriate push handles are included in the order.

| | |
|-----|-------|
| 12" | 51525 |
| 13" | 51526 |
| 14" | 51527 |
| 15" | 51528 |
| 16" | 51529 |
| 17" | 51530 |
| 18" | 51531 |
| 19" | 51532 |
| 20" | 51533 |



Old model



Compatibilité :

- VELOCE ✓
- HELIO C2 ✓
- HELIO A ✓
- HELIO kids ✓
- MOVE ✓