



Tie Down Conversion Kits

For UNOCCUPIED Transit Tie Downs

Compatible with the Helio C2, Helio KIDS, Helio A7 and MOVE

Conversion kits are for adding tie downs on a chair that originally didn't have any tie downs.



Configuration	Qty	Part	Short Description
Flip Back Armrest "U"	2	90982	Longer Back Cane Screw
	2	90670	Shorter Back Cane Screw
	4	90659	Caster Housing Screw
	2	32066	Plastic Spacer
	4	22004	Tie Downs
Swing away Armrest "L"	4	90659	Caster Housing Screw
	4	22004	Tie Downs
	4	90670	Back Cane Screw
"T" Armrest/ No Armrest	4	90669	Back Cane Screw
	2	32066	Plastic Spacer
	4	22004	Tie Downs
	4	90659	Caster Housing Screw
Depth Adjustable back	4	90659	Caster Housing Screw
	4	90670	Back Cane Screw
	4	99562	Washer
	1	22076	Kids Tie Down Right
	1	22077	Kids Tie Down Left
	2	22004	Tie Downs

Depth Adjustable back + "T" armrest = Depth adjustable tie down kit

Depth Adjustable back + No Armrest = Depth adjustable tie down kit

- One tie down is installed on both front casters and at the rear of each side frame.
- The curved tie downs in figure 1.2 are used with the depth adjustable back as part of the depth adjustable tie down kit.
- The curved tie down are only used for the depth adjustable back.