


MOVE
ORDER FORM - CANADA

P.O. number:* _____ Date: _____

Company Name:* _____ Phone:* _____

Ordered by:* _____

Clinician/ATP/RTS: _____

Email: _____

Marked for: _____

Ship to:* _____

Name: _____

Address: _____

City: _____ Prov: _____

Postal Code: _____

Phone: _____

* Required MVMA050 MVMA051 Taxes not included

English Manual (STD) *Manuel Français*

MODEL									
1	MVWC01	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)				WM0001196	\$1 395		
2	MVWC09	Heavy duty kit (350 lb weight capacity)				WAMF	\$775		
SEAT WIDTH ³									
3	MVSW	Seat width:	14" <input type="checkbox"/>	16" <input type="checkbox"/>	18" ¹ <input type="checkbox"/>			NC	<input type="checkbox"/>
4	MVSW	Seat width:		20" ¹ <input type="checkbox"/>	22" ² <input type="checkbox"/>		WAMA	\$295	<input type="checkbox"/>
SEAT DEPTH									
5	MVSD	Seat depth:		16" <input type="checkbox"/>	18" <input type="checkbox"/>			NC	<input type="checkbox"/>
6	MVSD	Seat depth:		20" <input type="checkbox"/>			WAMB	\$225	<input type="checkbox"/>

1: Available with HD option 2: Only available with HD option [3. Width calculator](#)

Move FRONT Seat-to-floor Height - Please select your choice

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2	13 1/2	14 1/2												
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1															
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1							15 1/2	16 1/2	17 1/2	16	17	18			
	P2				15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3				16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available + 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
7	MVCA01	Composite wheel w/ PU tire (1" width)	4" <input type="checkbox"/>	5" <input type="checkbox"/>	6" <input type="checkbox"/>	7" <input type="checkbox"/>	8" <input type="checkbox"/>			NC	<input type="checkbox"/>	
8	MVCA10	Pneumatic 6" x 1 1/4"								\$160	<input type="checkbox"/>	
9	MVCA11	Pneumatic 8" x 1 1/4"								\$160	<input type="checkbox"/>	
10	MVCA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 1/2"								\$98	<input type="checkbox"/>
11	MVCA15		5" x 1 1/2"								\$98	<input type="checkbox"/>
12	MVCA16		6" x 1 1/2"								\$98	<input type="checkbox"/>
13	MVCA17		4" x 1 1/2"								\$175	<input type="checkbox"/>
14	MVCA18	Newton UltraCaster <i>Soft-Roll</i> Aluminum	5" x 1 1/2"								\$185	<input type="checkbox"/>
15	MVCA19		6" x 1 1/2"								\$195	<input type="checkbox"/>

Move REAR Seat-to-floor Height - Please circle your choice

20"	22"	24"	25"	26"
12 1/4"	13"	14"	14 1/4"	15"
13 1/4"	14"	15"	15 1/4"	16"
14 1/4"	15"	16"	16 1/4"	17"
15 1/4"	16"	17"	17 1/4"	18"
16 1/4"	17"	18"	18 1/4"	19"
17 1/4"	18"	19"	19 1/4"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

SWING-AWAY FOOTREST										
Set Length at:										
16	MVFR01	60° swing away footrest See line 18 (5 - 13.5") <input type="checkbox"/>	13.5 - 17.5" <input type="checkbox"/>	15.5 - 19.5" <input type="checkbox"/>	17 - 21" <input type="checkbox"/>	Set at:*			\$250	<input type="checkbox"/>
17	MVFR02	70° swing away footrest See line 18 (5 - 13.5") <input type="checkbox"/>	13.5 - 17.5" <input type="checkbox"/>	15.5 - 19.5" <input type="checkbox"/>	17 - 21" <input type="checkbox"/>	Set at:*			\$250	<input type="checkbox"/>
FOOTREST OPTIONS										
18	MVFR09	High mount (footplate attached to the hanger with clamps)	Range: 5 - 13.5"			Set at:*			\$90	<input type="checkbox"/>

* +/- 1/2"

ELEVATING LEGREST (unit)									
					The elevating legrest replace the footrest side selected				
19	MVFR10	Elevating legrest with calf pad, 14" - 19" Length adjustable	Left <input type="checkbox"/>	Right <input type="checkbox"/>	Set at:		WAMK	\$390	<input type="checkbox"/>

FOOTPLATES									
20	MVFP01	Standard - Composite						STD	<input type="checkbox"/>
21	MVFP02	Newton adjustable angle - Composite					WAMJ	\$145	<input type="checkbox"/>
22	MVFP03	Newton adjustable angle - Aluminum					WAMJ	\$165	<input type="checkbox"/>

FOOTPLATE ACCESSORIES									
23	MVFP08	Length adjustable heel loop						\$25	<input type="checkbox"/>
24	MVFP09	Calf strap					SEFND0010	\$52	<input type="checkbox"/>
25	MVFP30	Bodypoint Aeromesh Padded calf strap				New	SEFND0010	\$129	<input type="checkbox"/>

RESIDUAL LIMB SUPPORT (Height and angle adjustable) (unit)									
26	MVFR23	Swing away residual limb support - 7" x 9"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			SEFND0005	\$390	<input type="checkbox"/>
27	MVFR12	Swing away residual limb support - 8" x 11"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			SEFND0005	\$390	<input type="checkbox"/>

Back canes with integrated fixed height push handle

FIXED ANGLE BACK										
28	MVBC01	8° bend aluminum back cane	16" <input type="checkbox"/>	18" <input type="checkbox"/>	20" <input type="checkbox"/>	STD			NC	<input type="checkbox"/>
ADJUSTABLE ANGLE BACK										
29	MVBC02	8° bend aluminum back cane		16" <input type="checkbox"/>	20" <input type="checkbox"/>	STD			\$309	<input type="checkbox"/>
30		Select back cane angle:	90° <input type="checkbox"/>	95° <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>			

STROLLER HANDLE AND STABILIZER BAR *									
31	MVPH04	Stroller handle (not available with the HD option)					WAMW	\$210	<input type="checkbox"/>
32	MVPH03	Folding stabilizer bar						\$258	<input type="checkbox"/>

* Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY										
33	MVBH05	Slip-on Nylon back upholstery							NC	<input type="checkbox"/>
34	MVBH06	Adjustable tension Nylon back upholstery					WAMD/WAME	\$245	<input type="checkbox"/>	
35	MVBH99	Omit back upholstery						STD	<input type="checkbox"/>	
RIGID BACK										
36		nXt back series (See nXt back pricelist)				See pricelist				<input type="checkbox"/>

SEAT CUSHIONS										
37	MVSC01	2" Cushion							\$175	<input type="checkbox"/>
38	MVSC02	3" Cushion							\$195	<input type="checkbox"/>
39	MVSC99	Omit seat cushion							NC	<input type="checkbox"/>
40		nXt FIT cushion series (See nXt FIT pricelist)				See pricelist				<input type="checkbox"/>

REAR WHEELS										
41	MVRW01	Mag ¹	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>				NC	<input type="checkbox"/>
42	MVRW02	Newton One - spoke wheel	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>				NC	<input type="checkbox"/>
43	MVRW03	Newton Gravity - ultralight wheel	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>		\$360	<input type="checkbox"/>

1 - HP tires are not available on Mag Wheels - use only aluminum or plastic coated handrim only

TIRES										
							Available			
44	MVTI01	Soft Urethane 1 1/4 "			20", 22" 24"				STD	<input type="checkbox"/>
45	MVTI02	Soft Urethane knobby 1 3/8 "			22", 24", 25", 26"				\$130	<input type="checkbox"/>
46	MVTI04	Soft Urethane - Shox 1"			22", 24", 25", 26"				\$155	<input type="checkbox"/>
47	MVTI05	Pneumatic 1 3/8 "			20", 22", 24", 25", 26"				\$130	<input type="checkbox"/>
48	MVTI08	Pneumatic with Airless insert 1 3/8"			22", 24"				\$215	<input type="checkbox"/>
49	MVTI06	Pneumatic HP - 1" *			20", 22", 24", 25", 26"				\$150	<input type="checkbox"/>
50	MVTI07	Pneumatic HP - Schwalbe Marathon Plus 1" *			24", 25", 26"				\$230	<input type="checkbox"/>

* Not available with Mag Wheels

HANDRIM										
51	MVHR01	Aluminum anodized							STD	<input type="checkbox"/>
52	MVHR07	Newton Air Grip* Black <input type="checkbox"/>					WAMO	\$150	<input type="checkbox"/>	
53	MVHR08	Newton Air Grip* (08) Sunkissed Orange <input type="checkbox"/> (09) Burgundy <input type="checkbox"/> (10) Sapphire Blue <input type="checkbox"/> (11) Monster Green <input type="checkbox"/> (12) Charcoal <input type="checkbox"/> (13) Fuchsia <input type="checkbox"/> (14) Steel Blue <input type="checkbox"/> (15) Ferrari Red <input type="checkbox"/> (16) Acid Green <input type="checkbox"/> (17) White <input type="checkbox"/> (18) Sky Blue <input type="checkbox"/>					WAMO	\$330	<input type="checkbox"/>	
Assembly position for aluminum and Newton AirGrip handrims										
54	MVHR29	Standard Handrim position							STD	<input type="checkbox"/>
55	MVHR30	Narrow Handrim position - Short tab (not available in 22")							\$25	<input type="checkbox"/>
56	MVHR31	Super Narrow Handrim - Short tab (not available in 22")							\$25	<input type="checkbox"/>

HANDRIM										
57	MVHR02	Plastic coated ²					WAMO	\$160	<input type="checkbox"/>	
58	MVHR03	Q-Grip High friction Neoprene coated* ^{1,2}					WAMO	\$430	<input type="checkbox"/>	
59	MVHR28	Surge* ^{1,2} (1-1/2" Oval handrim w/ Gripton high friction strip)						\$470	<input type="checkbox"/>	
60	MVHR05	Surge LT* ^{1,2} (1-1/8" Oval handrim w/ Gripton high friction strip)						\$500	<input type="checkbox"/>	
61	MVHR04	Natural Fit* Standard Grip ^{1,2} (1-1/2" Oval handrim w/standard thumb piece)						\$500	<input type="checkbox"/>	
62	MVHR25	Natural Fit* Super Grip ^{1,2} (1-1/2" Oval handrim w/super grip thumb piece)						\$530	<input type="checkbox"/>	
63	MVHR99	Omit handrim							NC	<input type="checkbox"/>


* Not available with Mag wheels 1 - Not available in 20" 2- Only standard handrim position available

AXLES (0° camber)										
64	MVAX01	Quick release axle					WAMU	STD	<input type="checkbox"/>	
65	MVAX02	Permanent axle							NC	<input type="checkbox"/>
66	MVAX03	Quad release axle					WAMU	\$170	<input type="checkbox"/>	

AXLE PLATE									
67	MVAX04	Amputee adjustable axle plate *					WAMT	\$250	<input type="checkbox"/>

* May limit seat to floor height available

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable	Most responsive
AMPUTEE	3/4"
<input type="checkbox"/>	<input type="checkbox"/>
	2"
	<input type="checkbox"/>

WHEEL LOCKS

68	MVWL08	Newton Wheel lock - push to lock				STD	<input type="checkbox"/>
69	MVWL02	Pull to lock				NC	<input type="checkbox"/>
70	MVWL04	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>			WAMR	\$175	<input type="checkbox"/>
71	MVWL13 MVWL14	6" removable extension handle				\$67	<input type="checkbox"/>
72	MVWL06	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>				\$230	<input type="checkbox"/>
73	MVWL99	Omit wheel locks				NC	<input type="checkbox"/>

FLIP BACK ARMREST (Convertible to "T" armrest)*

		Height Adjustment		Removable function			
74	MVAR01	10" Desk length arm pad	8 - 12" <input type="checkbox"/> STD	10 - 14" <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/> STD	\$380 <input type="checkbox"/>
75	MVAR02	14" Full length arm pad	8 - 12" <input type="checkbox"/> STD	10 - 14" <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/> STD	\$380 <input type="checkbox"/>
		Height Adjustment					
"T" ARMREST*							
76	MVAR03	10" Desk length arm pad	8 - 12" <input type="checkbox"/> STD	10 - 14" <input type="checkbox"/>			\$360 <input type="checkbox"/>
77	MVAR04	14" Full length arm pad	8 - 12" <input type="checkbox"/> STD	10 - 14" <input type="checkbox"/>			\$360 <input type="checkbox"/>

* May have a 1/2 " adjustment restriction on both limits of the range

GEL ARM PAD - Option for "T" Armrest

78	MVAR12	14" Full Length Arm Pad - Gel Ovation			SEAND0030	\$150	<input type="checkbox"/>
79	MVAR99	Omit armrest				NC	<input type="checkbox"/>

POSITIONING BELTS

80	MVSB01	Velcro adjustable	Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>			\$45 <input type="checkbox"/>
81	MVSB02	Auto buckle	Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>			\$55 <input type="checkbox"/>
82	MVSB07	Hardware attachment for 2 point seatbelt (1 pair)					NC <input type="checkbox"/>
83	MVSB03	Additional hardware for 2 extra attachment points <i>(needed for the 4 point seat belt) (1 pair)</i>					\$13 <input type="checkbox"/>
84	MVSB04	BodyPoint Hip Belt - 2 point not padded	M <input type="checkbox"/>	L <input type="checkbox"/>			\$82 <input type="checkbox"/>
85	MVSB05	BodyPoint Hip Belt - 2 point padded	S <input type="checkbox"/>	M <input type="checkbox"/>	L <input type="checkbox"/>	SERND0030	\$170 <input type="checkbox"/>
86	MVSB06	BodyPoint Hip Belt - 4 point padded	S <input type="checkbox"/>	M <input type="checkbox"/>	L <input type="checkbox"/>	SERND0020	\$195 <input type="checkbox"/>
87	MVSB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	M <input type="checkbox"/>	L <input type="checkbox"/>		New	\$324 <input type="checkbox"/>
88	MVSB99	Omit positioning belt					NC <input type="checkbox"/>

HEADREST *

89	MVAC01	Profiled Headrest	Small (3.5"X6") <input type="checkbox"/>	Medium (4.5"X7") <input type="checkbox"/>	Large (5.5"X8") <input type="checkbox"/>		\$205 <input type="checkbox"/>
90	MVAC02	Headrest mounting hardware with stabiliser bar *					\$338 <input type="checkbox"/>


* Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES												
91	MVAC03	Anti tippers								\$115	<input type="checkbox"/>	
92	MVAC04	Transit Tie Down (<i>Unoccupied</i>)								NC	<input type="checkbox"/>	
93	MVAC06	Cane carrier		Left <input type="checkbox"/>		Right <input type="checkbox"/>				\$105	<input type="checkbox"/>	
94	MVAC07	Oxygen tank holder - D Type <i>Not compatible with swing away armrest or amputee axle plate</i>		Left					WAMX	\$375	<input type="checkbox"/>	
95	MVAC08	Spoke guards*	Available: 20", 22", 24", 26"							WAMP	\$145	<input type="checkbox"/>
96	MVAC09	Footrest neoprene impact guard (<i>pair</i>) - Black							SEFND0035	\$80	<input type="checkbox"/>	
97	MVAC10	Seat pouch								\$95	<input type="checkbox"/>	
98	MVAC26	BodyPoint Mobility Bag							New	\$135	<input type="checkbox"/>	
99	MVAC23	Newton Backpack								\$120	<input type="checkbox"/>	
100	MVAC12	Velcro Closing Strap								\$50	<input type="checkbox"/>	
101	MVAC13	Tool kit								\$50	<input type="checkbox"/>	
102	MVAC16	Motion Composites Polo shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>					\$32	<input type="checkbox"/>	
			M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>	M/XXXXL <input type="checkbox"/>			

* Not available with Mag wheels

FRAME COLOR (with black aluminum crossbrace)												
103	MVFC01	Black <input type="checkbox"/>								NC	<input type="checkbox"/>	
104	MVFC03	Sunkissed Orange <input type="checkbox"/>	Burgundy <input type="checkbox"/>	Sapphire Blue <input type="checkbox"/>	Monster Green <input type="checkbox"/>					\$270	<input type="checkbox"/>	
		Charcoal <input type="checkbox"/>	Fuchsia <input type="checkbox"/>	Steel Blue <input type="checkbox"/>	Ferrari Red <input type="checkbox"/>	Acid Green <input type="checkbox"/>						
105	MVFC99	Custom Color	Specify PPG Paintit color:								\$410	<input type="checkbox"/>

TOTAL										
--------------	--	--	--	--	--	--	--	--	--	--



Notes :