


MOVE
ORDER FORM - INTERNATIONAL

P.O. number:*	Date:	Ship to:*
Company Name:*	Phone:*	Name:
Ordered by:*		Address:
Clinician/ATP/RTS:		City: State/Prov:
Email:		Postal/Zip Code:
Marked for:	Phone:	* Required

MVMA053 <input type="checkbox"/>	MVMA054 <input type="checkbox"/>	MVMA055 <input type="checkbox"/>	MVMA056 <input type="checkbox"/>	MVMA057 <input type="checkbox"/>	MVMA058 <input type="checkbox"/>	
English Manual (STD)	German Manual	Italian Manual	Dutch Manual	Spanish Manual	French Manual	Taxes not included

MODEL					
1	MVWC01	Ultralight MOVE wheelchair aluminum. (120 kg/265 lb weight capacity)			\$1 350 <input type="checkbox"/>
2	MVWC09	Heavy duty kit (159 kg/350 lb weight capacity)			\$500 <input type="checkbox"/>
SEAT WIDTH ³					
3	MVSW	Seat width: 35.6 (14") <input type="checkbox"/> 40.6 (16") <input type="checkbox"/> 45.7 (18") ¹ <input type="checkbox"/>			NC <input type="checkbox"/>
4	MVSW	Seat width: 50.8 (20") ¹ <input type="checkbox"/> 55.9 (22") ² <input type="checkbox"/>			\$285 <input type="checkbox"/>
SEAT DEPTH					
5	MVSD	Seat depth: 40.6 (16") <input type="checkbox"/> 45.7 (18") <input type="checkbox"/>			NC <input type="checkbox"/>
6	MVSD	Seat depth: 50.8 (20") <input type="checkbox"/>			\$215 <input type="checkbox"/>

1: Available with HD option

2: Only available with HD option

[3. Width calculator](#)
Move FRONT Seat-to-floor Height - Please select your choice

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+ 2.54 cm	+ 5.08cm	STD	+ 2.54 cm	+ 5.08cm	STD	+ 2.54 cm	+ 5.08cm	STD	+ 2.54 cm	+ 5.08cm	STD	+ 2.54 cm	+ 5.08cm
4 in 10,2cm	P1	31.8cm 12.5"	34.3cm 13.5"	36.8cm 14.5"												
	P2	33cm 13"	35.6cm 14"	38.1cm 15"	34.3cm 13.5"	36.8cm 14.5"	39.4cm 15.5"									
5 in 12,7cm	P1															
	P2	33cm 13"	35.6cm 14"	38.1cm 15"	34.3cm 13.5"	36.8cm 14.5"	39.4cm 15.5"									
	P3	34.3cm 13.5"	36.8cm 14.5"	39.4cm 15.5"	35.6cm 14"	38.1cm 15"	40.6cm 16"	36.8cm 14.5"	39.4cm 15.5"	41.9cm 16.5"						
	P4	35.6cm 14"	38.1cm 15"	40.6cm 16"	36.8cm 14.5"	39.4cm 15.5"	41.9cm 16.5"	38.1cm 15"	40.6cm 16"	43.2cm 17"	39.4cm 15.5"	41.9cm 16.5"	44.5cm 17.5"			
7 in 17,8cm	P1							39.4cm 15.5"	41.9cm 16.5"	44.5cm 17.5"	40.6cm 16"	43.2cm 17"	45.7cm 18"			
	P2				39.4cm 15.5"	41.9cm 16.5"	44.5cm 17.5"	40.6cm 16"	43.2cm 17"	45.7cm 18"	41.9cm 16.5"	44.5cm 17.5"	47cm 18.5"			
	P3				40.6cm 16"	43.2cm 17"	45.7cm 18"	41.9cm 16.5"	44.5cm 17.5"	47cm 18.5"	43.2cm 17"	45.7cm 18"	48.3cm 19"	44.5cm 17.5"	47cm 18.5"	49.5cm 19.5"
	P4	40.6cm 16"	43.2cm 17"	45.7cm 18"	41.9cm 16.5"	44.5cm 17.5"	47cm 18.5"	43.2cm 17"	45.7cm 18"	48.3cm 19"	44.5cm 17.5"	47cm 18.5"	49.5cm 19.5"	45.7cm 18"	48.3cm 19"	50.8cm 20"
	P5	41.9cm 16.5"	44.5cm 17.5"	47cm 18.5"	43.2cm 17"	45.7cm 18"	48.3cm 19"	44.5cm 17.5"	47cm 18.5"	49.5cm 19.5"	45.7cm 18"	48.3cm 19"	50.8cm 20"	47cm 18.5"	49.5cm 19.5"	52.1cm 20.5"

Height available +/- 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

CASTERS												
7	MVCA01	Composite wheel w/ PU tire (1" width)	4" <input type="checkbox"/>	5" <input type="checkbox"/>	6" <input type="checkbox"/>	7" <input type="checkbox"/>	8" <input type="checkbox"/>			NC	<input type="checkbox"/>	
8	MVCA10	Pneumatic 6" x 1 1/4"								\$160	<input type="checkbox"/>	
9	MVCA11	Pneumatic 8" x 1 1/4"								\$160	<input type="checkbox"/>	
10	MVCA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 1/2"								\$80	<input type="checkbox"/>
11	MVCA15		5" x 1 1/2"								\$80	<input type="checkbox"/>
12	MVCA16		6" x 1 1/2"								\$80	<input type="checkbox"/>
13	MVCA17		4" x 1 1/2"								\$175	<input type="checkbox"/>
14	MVCA18	Newton UltraCaster <i>Soft-Roll</i> Aluminum	5" x 1 1/2"								\$185	<input type="checkbox"/>
15	MVCA19		6" x 1 1/2"								\$195	<input type="checkbox"/>

Move - REAR Seat-to-floor Height - Please select your choice

20"	22"	24"	25"	26"
31.1cm / 12.25in	33.7cm / 13.25in	35.6cm / 14in	36.8cm / 14.5in	38.1cm / 15in
33.7cm / 13.25in	36.2cm / 14.25in	38.1cm / 15in	39.4cm / 15.5in	40.6cm / 16in
36.2cm / 14.25in	38.7cm / 15.25in	40.6cm / 16in	41.9cm / 16.5in	43.2cm / 17in
38.7cm / 15.25in	41.3cm / 16.25in	43.2cm / 17in	44.5cm / 17.5in	45.7cm / 18in
41.3cm / 16.25in	43.8cm / 17.25in	45.7cm / 18in	47cm / 18.5in	48.3cm / 19in
43.8cm / 17.25in	46.4cm / 18.25in	48.3cm / 19in	49.5cm / 19.5in	50.8cm / 20in

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1.25 cm (1/2") of selected seat to floor height"

SWING-AWAY FOOTREST											
Set Length at:											
16	MVFR01	60° swing away footrest	See line 18 (12,7 - 34,3cm) <input type="checkbox"/>	34,3 - 44,5cm <input type="checkbox"/>	39,4-49,5 cm <input type="checkbox"/>	43,2 - 53,3cm <input type="checkbox"/>	Set at:*			NC	<input type="checkbox"/>
17	MVFR02	70° swing away footrest	See line 18 (12,7 - 34,3cm) <input type="checkbox"/>	34,3 - 44,5cm <input type="checkbox"/>	39,4-49,5 cm <input type="checkbox"/>	43,2 - 53,3cm <input type="checkbox"/>	Set at:*			STD	<input type="checkbox"/>
FOOTREST OPTIONS											
18	MVFR09	High mount (footplate attached to the hanger with clamps)				Range: 12,7 - 34,3cm	Set at:*			\$75	<input type="checkbox"/>
										* +/- 1/2"	

ELEVATING LEGREST (unit)											
19	MVFR10	Elevating legrest with calf pad, 14" - 19" Length adjustable			Left <input type="checkbox"/>	Right <input type="checkbox"/>	Set at:			\$240	<input type="checkbox"/>

FOOTPLATES											
20	MVFP01	Standard - Composite								STD	<input type="checkbox"/>
21	MVFP02	Newton adjustable angle - Composite								\$210	<input type="checkbox"/>
22	MVFP03	Newton adjustable angle - Aluminum								\$230	<input type="checkbox"/>

FOOTPLATE ACCESSORIES											
23	MVFP08	Length adjustable heel loop								\$40	<input type="checkbox"/>
24	MVFP09	Calf strap								\$35	<input type="checkbox"/>
25	MVFP30	Bodypoint Aeromesh Padded calf strap					New	SEFND0010	\$115	<input type="checkbox"/>	

RESIDUAL LIMB SUPPORT (Height and angle adjustable) (unit)										
26	MVFR23	Swing away residual limb support - 7" x 9"	17.8 x 22.8 cm	Left <input type="checkbox"/>	Right <input type="checkbox"/>				\$375	<input type="checkbox"/>
27	MVFR12	Swing away residual limb support - 8" x 11"	20.3 x 27.9 cm	Left <input type="checkbox"/>	Right <input type="checkbox"/>				\$375	<input type="checkbox"/>

Back canes with integrated fixed height push handle

FIXED ANGLE BACK										
28	MVBC01	8° bend aluminum back cane	16" <input type="checkbox"/>	18" <input type="checkbox"/>	20" <input type="checkbox"/>	STD			NC	<input type="checkbox"/>
ADJUSTABLE ANGLE BACK										
29	MVBC02	8° bend aluminum back cane	16" <input type="checkbox"/>	20" <input type="checkbox"/>	STD				\$300	<input type="checkbox"/>
30		Select back cane angle:	90° <input type="checkbox"/>	95° <input type="checkbox"/>	STD	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>		

STROLLER HANDLE AND STABILIZER BAR *						
31	MVPH04	Stroller handle <i>(not available with the HD option)</i>			\$160	<input type="checkbox"/>
32	MVPH03	Folding stabilizer bar			\$250	<input type="checkbox"/>

* Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY						
33	MVBH05	Slip-on Nylon back upholstery			NC	<input type="checkbox"/>
34	MVBH06	Adjustable tension Nylon back upholstery			\$165	<input type="checkbox"/>
35	MVBH99	Omit back upholstery			STD	<input type="checkbox"/>

SEAT CUSHIONS						
36	MVSC01	2" Cushion			\$165	<input type="checkbox"/>
37	MVSC02	3" Cushion			\$185	<input type="checkbox"/>
38	MVSC99	Omit seat cushion			NC	<input type="checkbox"/>

REAR WHEELS									
39	MVRW01	Mag ¹	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	\$50	<input type="checkbox"/>		
40	MVRW02	Newton One - spoke wheel	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	NC	<input type="checkbox"/>		
41	MVRW03	Newton Gravity - ultralight wheel	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>	\$300	<input type="checkbox"/>

1 - HP tires are not available on Mag Wheels - use only aluminum or plastic coated handrim only

TIRES		Available								
42	MVTI01	Soft Urethane 1 1/4 "	20", 22" 24"						STD	<input type="checkbox"/>
43	MVTI02	Soft Urethane knobby 1 3/8 "	22", 24", 25", 26"						\$95	<input type="checkbox"/>
44	MVTI04	Soft Urethane - Shox 1"	22", 24", 25", 26"						\$115	<input type="checkbox"/>
45	MVTI05	Pneumatic 1 3/8 "	20", 22", 24", 25", 26"						NC	<input type="checkbox"/>
46	MVTI08	Pneumatic with Airless insert 1 3/8"	22", 24"						\$80	<input type="checkbox"/>
47	MVTI06	Pneumatic HP - 1" *	20", 22", 24", 25", 26"						\$95	<input type="checkbox"/>
48	MVTI07	Pneumatic HP - Schwalbe Marathon Plus 1" *	24", 25", 26"						\$125	<input type="checkbox"/>

* Not available with Mag Wheels

HANDRIM									
49	MVHR01	Aluminum anodized				STD	<input type="checkbox"/>		
50	MVHR07	Newton Air Grip* Black <input type="checkbox"/>				\$95	<input type="checkbox"/>		
51	MVHR08	Newton Air Grip* (08) Sunkissed Orange <input type="checkbox"/>	(09) Burgundy <input type="checkbox"/>	(10) Sapphire Blue		\$240	<input type="checkbox"/>		
		(11) Monster Green <input type="checkbox"/>	(12) Charcoal <input type="checkbox"/>	(13) Fuchsia <input type="checkbox"/>	(14) Steel Blue				
		(15) Ferrari Red <input type="checkbox"/>	(16) Acid Green <input type="checkbox"/>	(17) White <input type="checkbox"/>	(18) Sky Blue <input type="checkbox"/>				
Assembly position for aluminum and Newton AirGrip handrims									
52	MVHR29	Standard Handrim position				STD	<input type="checkbox"/>		
53	MVHR30	Narrow Handrim position - Short tab <i>(not available in 22")</i>				\$25	<input type="checkbox"/>		
54	MVHR31	Super Narrow Handrim - Short tab <i>(not available in 22")</i>				\$25	<input type="checkbox"/>		

HANDRIM							
55	MVHR02	Plastic coated ²				\$95	<input type="checkbox"/>
56	MVHR03	Q-Grip High friction Neoprene coated* ^{1,2}				\$350	<input type="checkbox"/>
57	MVHR28	Surge * ^{1,2} <i>(1-1/2" Oval handrim w/ Gripton high friction strip)</i>				\$375	<input type="checkbox"/>
58	MVHR05	Surge LT* ^{1,2} <i>(1-1/8" Oval handrim w/ Gripton high friction strip)</i>				\$405	<input type="checkbox"/>
59	MVHR04	Natural Fit* Standard Grip ^{1,2} <i>(1-1/2" Oval handrim w/standard thumb piece)</i>				\$395	<input type="checkbox"/>
60	MVHR25	Natural Fit* Super Grip ^{1,2} <i>(1-1/2" Oval handrim w/super grip thumb piece)</i>				\$450	<input type="checkbox"/>
61	MVHR99	Omit handrim				NC	<input type="checkbox"/>

* Not available with Mag wheels


1 - Not available in 20"

2- Only standard handrim position available

AXLES (0° camber)						
62	MVAX01	Quick release axle			STD	<input type="checkbox"/>
63	MVAX02	Permanent axle			NC	<input type="checkbox"/>
64	MVAX03	Quad release axle			\$125	<input type="checkbox"/>
AXLE PLATE						
65	MVAX04	Amputee adjustable axle plate *			\$195	<input type="checkbox"/>

* May limit seat to floor height available

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.
Default setting is at 3/4"

Most stable	Most responsive
AMPUTEE	1,9 cm - 3/4"
<input type="checkbox"/>	<input type="checkbox"/>
	5,1 cm - 2"
	<input type="checkbox"/>

WHEEL LOCKS						
66	MVWL08	Newton Wheel lock - push to lock			STD	<input type="checkbox"/>
67	MVWL02	Pull to lock			NC	<input type="checkbox"/>
68	MVWL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)			\$150	<input type="checkbox"/>
69	MVWL13 MVWL14	6" removable extension handle			\$75	<input type="checkbox"/>
70	MVWL06	Attendant Lock (additional set of wheel lock mounted to the back)			\$200	<input type="checkbox"/>
71	MVWL99	Omit wheel locks			NC	<input type="checkbox"/>

FLIP BACK ARMREST (Convertible to "T" armrest)*		Height Adjustment		Removable function			
72	MVAR01	10" Desk length arm pad	20 - 30cm (8 - 12") <input type="checkbox"/> STD	25 - 36cm (10 - 14") <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/> STD	\$350 <input type="checkbox"/>
73	MVAR02	14" Full length arm pad	20 - 30cm (8 - 12") <input type="checkbox"/> STD	25 - 36cm (10 - 14") <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/> STD	\$350 <input type="checkbox"/>
"T" ARMREST*		Height Adjustment					
74	MVAR03	10" Desk length arm pad	20 - 30cm (8 - 12") <input type="checkbox"/> STD	25 - 36cm (10 - 14") <input type="checkbox"/>			\$325 <input type="checkbox"/>
75	MVAR04	14" Full length arm pad	20 - 30cm (8 - 12") <input type="checkbox"/> STD	25 - 36cm (10 - 14") <input type="checkbox"/>			\$325 <input type="checkbox"/>

* May have a 1/2" adjustment restriction on both limits of the range

GEL ARM PAD - Option for "T" Armrest						
76	MVAR12	36cm (14") Full Length Arm Pad - Gel Ovation			\$150	<input type="checkbox"/>

77	MVAR99	Omit armrest			NC	<input type="checkbox"/>
----	--------	--------------	--	--	----	--------------------------

POSITIONING BELTS						
78	MVSB01	Velcro adjustable	Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>		\$45 <input type="checkbox"/>
79	MVSB02	Auto buckle	Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>		\$85 <input type="checkbox"/>
80	MVSB07	Hardware attachment for 2 point seatbelt (1 pair)				NC <input type="checkbox"/>
81	MVSB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)				\$15 <input type="checkbox"/>
82	MVSB04	BodyPoint Hip Belt - 2 point not padded	M <input type="checkbox"/>	L <input type="checkbox"/>		\$80 <input type="checkbox"/>
83	MVSB05	BodyPoint Hip Belt - 2 point padded	S <input type="checkbox"/>	M <input type="checkbox"/>	L <input type="checkbox"/>	\$170 <input type="checkbox"/>
84	MVSB06	BodyPoint Hip Belt - 4 point padded	S <input type="checkbox"/>	M <input type="checkbox"/>	L <input type="checkbox"/>	\$190 <input type="checkbox"/>
85	MVSB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	M <input type="checkbox"/>	L <input type="checkbox"/>	New	\$285 <input type="checkbox"/>
86	MVSB99	Omit positioning belt				NC <input type="checkbox"/>

HEADREST						
87	MVAC01	Profiled Headrest	Small 8,8 x 15,2 cm <input type="checkbox"/>	Medium 11,4 x 17,8 cm <input type="checkbox"/>	Large 14 x 2,3 cm <input type="checkbox"/>	\$180 <input type="checkbox"/>
88	MVAC02	Headrest mounting hardware with stabiliser bar *				\$295 <input type="checkbox"/>


* Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES											
89	MVAC03	Anti tippers							\$140	<input type="checkbox"/>	
90	MVAC04	Transit Tie Down (<i>Unoccupied</i>)							NC	<input type="checkbox"/>	
91	MVAC06	Cane carrier	Left <input type="checkbox"/>	Right <input type="checkbox"/>					\$120	<input type="checkbox"/>	
92	MVAC07	Oxygen tank holder - D Type <i>Not compatible with swing away armrest or amputee axle plate</i>	Left						\$310	<input type="checkbox"/>	
93	MVAC08	Spoke guards*	Available: 20", 22", 24", 26"							\$115	<input type="checkbox"/>
94	MVAC09	Footrest neoprene impact guard (pair) - Black							\$90	<input type="checkbox"/>	
95	MVAC10	Seat pouch							\$75	<input type="checkbox"/>	
96	MVAC26	BodyPoint Mobility Bag						New	\$125	<input type="checkbox"/>	
97	MVAC23	Newton Backpack							\$120	<input type="checkbox"/>	
98	MVAC12	Velcro Closing Strap							\$50	<input type="checkbox"/>	
99	MVAC13	Tool kit							\$45	<input type="checkbox"/>	
100	A7AC16	Motion Composites Polo shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>					\$32	<input type="checkbox"/>
			M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>	M/XXXXL <input type="checkbox"/>		

* Not available with Mag wheels

FRAME COLOR (with black aluminum crossbrace)											
101	MVFC01	Black <input type="checkbox"/>							NC	<input type="checkbox"/>	
102	MVFC03	Sunkissed Orange <input type="checkbox"/>	Burgundy <input type="checkbox"/>	Sapphire Blue <input type="checkbox"/>	Monster Green <input type="checkbox"/>				\$270	<input type="checkbox"/>	
		Charcoal <input type="checkbox"/>	Fuchsia <input type="checkbox"/>	Steel Blue <input type="checkbox"/>	Ferrari Red <input type="checkbox"/>	Acid Green <input type="checkbox"/>					
103	MVFC99	Custom Color	Specify PPG Paintit color:							\$360	<input type="checkbox"/>

TOTAL										
--------------	--	--	--	--	--	--	--	--	--	--



Notes :