



160 Armand-Majeau Sud  
 Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0  
 T.+1-866-650-6555 F. +1-888-966-6555  
[orders@motioncomposites.com](mailto:orders@motioncomposites.com)  
[www.motioncomposites.com](http://www.motioncomposites.com)



<b>Move</b>	<b>MVqs</b>
<b>Quick Ship Order Form</b>	
<b>Canada</b>	<b>4.8</b>

P.O. number:\*

Company Name:\*

Ordered by:\*

Clinician: Phone:\*

Email:

Marked for:

Ship to:\*

Name:

Address:

City: Province

Postal Code:

Phone:

\* Required Taxes not included

MOVE		MVQSWC01		
Ultraight MOVE wheelchair aluminum	Basalt Black	WM0001196	\$1 395	<input type="checkbox"/>
Included: seating upholstery / Excluded back upholstery and cushion	265 lbs. weight capacity			

**Front Seat-to-Floor Height Configuration**

13"				14 1/2"				16 1/2"						
M6 16x16	NC	<input type="checkbox"/>	M1 18x16	NC	<input type="checkbox"/>	M11 16x18	NC	<input type="checkbox"/>	M7 16x18	NC	<input type="checkbox"/>	M12 18x18	NC	<input type="checkbox"/>
M9 18x18	NC	<input type="checkbox"/>	M2 18x18	NC	<input type="checkbox"/>	M13 18x20	WAMB \$225	<input type="checkbox"/>	M3 18x20	WAMB \$225	<input type="checkbox"/>	M14 20x18	WAMA \$295	<input type="checkbox"/>
M10 20x18	WAMA \$295	<input type="checkbox"/>	M4 20x20	WAMA WAMB \$520	<input type="checkbox"/>									
WHEELS Front, Rear & Handrim				WHEELS Front, Rear & Handrim				WHEELS Front, Rear & Handrim						
4"x1" Composites caster w/PU tire	NC	<input type="checkbox"/>	5"x1" Composites caster w/PU tire	NC	<input type="checkbox"/>	6"x1" Composites caster w/PU tire	NC	<input type="checkbox"/>	5" fork + Std stem bolt	NC	<input type="checkbox"/>	7" fork + Std stem bolt	NC	<input type="checkbox"/>
22" Newton One rear wheel	NC	<input type="checkbox"/>	22" Newton One rear wheel	NC	<input type="checkbox"/>	24" Newton One rear wheel	NC	<input type="checkbox"/>	Soft urethane tire 1 1/4"	NC	<input type="checkbox"/>	Soft urethane tire 1 1/4"	NC	<input type="checkbox"/>
Aluminum anodized handrim	NC	<input type="checkbox"/>	Aluminum anodized handrim	NC	<input type="checkbox"/>	Aluminum anodized handrim	NC	<input type="checkbox"/>						
SWING-AWAY FOOTREST				SWING-AWAY FOOTREST				SWING-AWAY FOOTREST						
60° Swing-away footrests	\$250	<input type="checkbox"/>	70° Swing-away footrests	\$250	<input type="checkbox"/>	70° Swing-away footrests	\$250	<input type="checkbox"/>	Newton Composite Adjustable angle footplate w/ heel loop	WAMJ \$170	<input type="checkbox"/>	Newton Composite Adjustable angle footplate w/ heel loop	WAMJ \$170	<input type="checkbox"/>
<b>This configuration:</b>				<b>This configuration:</b>				<b>This configuration:</b>						
Front STF Height	Rear STF Height		Front STF Height	Rear STF Height		Front STF Height	Rear STF Height		13" - 14"	13" - 18"		13 1/2" - 14 1/2"	13" - 18"	
				<b>Possible configurations:</b>										
Front STF Height				12 1/2" - 19 1/2"				Rear STF Height						
								12 1/4" - 20"						

PARTS & ACCESSORIES on all configurations			
<b>8° BEND BACK</b> (with fixed height push handles)			
Aluminum back cane	20" height	NC	<input type="checkbox"/>
<b>FLIP BACK ARMREST</b> (Convertible to "T" armrest)			
14" Full Length, 8 - 12" height adjustable		\$380	<input type="checkbox"/>
<b>AXLE</b> (with 0° camber)			
Quick release axle		WAMU NC	<input type="checkbox"/>
<b>WHEEL LOCK</b>			
Push to lock with 6" removable extension handle		\$67	<input type="checkbox"/>
<b>POSITIONING BELT</b>			
Additional attachment hardware for 4 point seatbelt		\$13	<input type="checkbox"/>
<b>ACCESSORIES</b>			
Anti tippers		\$115	<input type="checkbox"/>

# Quick Ship Parts and Accessories Program

Motion Composites' Quick Ship Parts and Accessories Program adds flexibility to the convenience of our Quick Ship wheelchairs

These specially discounted parts will be added to your Quick Ship wheelchair order for you to install and customize upon receipt

Invoice code					Invoice code				
QP04	Aluminum straight back cane 18-21"	PR	\$110	<input type="checkbox"/>	QP36	4" fork (2 holes)	PR	\$60	<input type="checkbox"/>
QP05	Aluminum back cane 8° bend 18-21"	PR	\$130	<input type="checkbox"/>	QP37	5" fork (4 holes)	PR	\$60	<input type="checkbox"/>
QP06	Angle adjustable - Aluminium straight back cane 18-21"	PR	\$230	<input type="checkbox"/>	QP38	7" fork (5 holes)	PR	\$60	<input type="checkbox"/>
QP98	Push handles	PR	\$100	<input type="checkbox"/>	QP114	Newton One wheels w/Urethane tire & Alu Pushrim - 22"	PR	\$420	<input type="checkbox"/>
QP07	Slip-on Active-Nylon back upholstery	EA	\$61	<input type="checkbox"/>	QP115	Newton One wheels w/Urethane tire & Alu Pushrim - 24"	PR	\$420	<input type="checkbox"/>
QP08	Adjustable tension Active-Mesh back upholstery	EA	\$245	<input type="checkbox"/>	QP39	Mag wheels w/Urethane tire & Alu Pushrim - 20"	PR	\$420	<input type="checkbox"/>
QP99	Swing away residual limb support - Left 7"x9"	EA	\$390	<input type="checkbox"/>	QP40	Mag wheels w/Urethane tire & Alu Pushrim - 22"	PR	\$420	<input type="checkbox"/>
QP100	Swing away residual limb support - Right 7"x9"	EA	\$390	<input type="checkbox"/>	QP41	Mag wheels w/Urethane tire & Alu Pushrim - 24"	PR	\$420	<input type="checkbox"/>
QP09	Swing away residual limb support - Left 8"x11"	EA	\$390	<input type="checkbox"/>	QP44	Wheel lock - pull to lock	PR	\$96	<input type="checkbox"/>
QP10	Swing away residual limb support - Right 8"x11"	EA	\$390	<input type="checkbox"/>	QP45	Wheel lock - pull to lock with removable extension handle	PR	\$163	<input type="checkbox"/>
QP121	Swing away elevating legrest/calf pad/ std footplate/heal loop, Left	EA	\$390	<input type="checkbox"/>	QP46	Wheel lock - Newton grade aid push to lock	PR	\$285	<input type="checkbox"/>
QP122	Swing away elevating legrest/calf pad/ std footplate/heal loop, Right	EA	\$390	<input type="checkbox"/>	QP48	Positioning belt - velcro adjustable	EA	\$45	<input type="checkbox"/>
QP13	Calf strap (wheelchair width: 15 to 17 in)	EA	\$52	<input type="checkbox"/>	QP49	Positioning belt - auto buckle	EA	\$55	<input type="checkbox"/>
QP14	Calf strap (wheelchair width: 18 to 20 in)	EA	\$52	<input type="checkbox"/>	QP50	Bodypoint 2 point Non Padded - medium (56")	EA	\$82	<input type="checkbox"/>
QP125	Bodypoint Aeromesh Padded calf strap (wheelchair Width: 12 to 16 in)	EA	\$129	<input type="checkbox"/>	QP51	Bodypoint 2 point Non Padded - large (61")	EA	\$82	<input type="checkbox"/>
QP126	Bodypoint Aeromesh Padded calf strap (wheelchair width: 17 to 20 in)	EA	\$129	<input type="checkbox"/>	QP52	Bodypoint 2 point Padded - small (51")	EA	\$155	<input type="checkbox"/>
QP127	Bodypoint Aeromesh Padded calf strap (wheelchair width: 21 to 24 in)	EA	\$129	<input type="checkbox"/>	QP53	Bodypoint 2 point Padded - medium (56")	EA	\$155	<input type="checkbox"/>
QP17	Composite caster 3" x 1"	PR	\$80	<input type="checkbox"/>	QP54	Bodypoint 2 point Padded - large (61")	EA	\$155	<input type="checkbox"/>
QP18	Composite caster 4" x 1"	PR	\$80	<input type="checkbox"/>	QP55	Bodypoint 4 point Padded - small (51")	EA	\$195	<input type="checkbox"/>
QP19	Composite caster 5" x 1"	PR	\$90	<input type="checkbox"/>	QP56	Bodypoint 4 point Padded - medium (56")	EA	\$195	<input type="checkbox"/>
QP20	Composite caster 6" x 1"	PR	\$90	<input type="checkbox"/>	QP57	Bodypoint 4 point Padded - large (61")	EA	\$195	<input type="checkbox"/>
QP21	Composite caster 7" x 1"	PR	\$100	<input type="checkbox"/>	QP58	Cane carrier - Left	EA	\$105	<input type="checkbox"/>
QP22	Composite caster 8" x 1"	PR	\$100	<input type="checkbox"/>	QP59	Cane carrier - Right	EA	\$105	<input type="checkbox"/>
QP77	Newton Ultracaster Composite caster 4" x 1½" Soft-Roll	PR	\$98	<input type="checkbox"/>	QP60	Spoke guards - 22"	PR	\$145	<input type="checkbox"/>
QP78	Newton Ultracaster Composite caster 5" x 1½" Soft-Roll	PR	\$98	<input type="checkbox"/>	QP61	Spoke guards - 24"	PR	\$145	<input type="checkbox"/>
QP79	Newton Ultracaster Composite caster 6" x 1½" Soft-Roll	PR	\$98	<input type="checkbox"/>	QP62	Footrest neoprene impact guard (pr) - Black	PR	\$80	<input type="checkbox"/>
QP80	Newton Ultracaster Aluminum caster 4" x 1½" Soft-Roll	PR	\$175	<input type="checkbox"/>	QP63	Seat pouch	EA	\$95	<input type="checkbox"/>
QP81	Newton Ultracaster Aluminum caster 5" x 1½" Soft-Roll	PR	\$185	<input type="checkbox"/>	QP103	Newton Backpack	EA	\$120	<input type="checkbox"/>
QP82	Newton Ultracaster Aluminum caster 6" x 1½" Soft-Roll	PR	\$195	<input type="checkbox"/>	QP106	Motion Composites Polo-Shirt Men small	EA	\$32	<input type="checkbox"/>
QP31	Pneumatic caster 6" x 1¼"	PR	\$160	<input type="checkbox"/>	QP107	Motion Composites Polo-Shirt Men medium	EA	\$32	<input type="checkbox"/>
QP32	Pneumatic caster 8" x 1¼"	PR	\$160	<input type="checkbox"/>	QP108	Motion Composites Polo-Shirt Men large	EA	\$32	<input type="checkbox"/>
QP33	Standard stem bolt	PR	\$18	<input type="checkbox"/>	QP109	Motion Composites Polo-Shirt Men X-large	EA	\$32	<input type="checkbox"/>
QP34	+1" stem bolt	PR	\$18	<input type="checkbox"/>	QP110	Motion Composites Polo-Shirt Men XX-large	EA	\$32	<input type="checkbox"/>
QP35	+2" stem bolt	PR	\$18	<input type="checkbox"/>	QP118	Motion Composites Polo-Shirt Men XXX-large	EA	\$32	<input type="checkbox"/>
QP76	Length adjustable heel loop	PR	\$40	<input type="checkbox"/>	QP119	Motion Composites Polo-Shirt Men XXXX-large	EA	\$32	<input type="checkbox"/>
QP96	Newton Arigrip 22" black handrim	PR	\$150	<input type="checkbox"/>	QP111	Motion Composites Polo-Shirt Women small	EA	\$32	<input type="checkbox"/>
QP97	Newton Arigrip 24" black handrim	PR	\$150	<input type="checkbox"/>	QP112	Motion Composites Polo-Shirt Women medium	EA	\$32	<input type="checkbox"/>
					QP113	Motion Composites Polo-Shirt Women large	EA	\$32	<input type="checkbox"/>
					QP71	Tool Kit	EA	\$50	<input type="checkbox"/>

### Move REAR Seat-to-floor Height

20"	22"	24"	25"	26"
12 1/4"	13"	14"	14 1/4"	15"
13 1/4"	14"	15"	15 1/4"	16"
14 1/4"	15"	16"	16 1/4"	17"
15 1/4"	16"	17"	17 1/4"	18"
16 1/4"	17"	18"	18 1/4"	19"
17 1/4"	18"	19"	19 1/4"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 3/8" of selected seat to floor height"

Table 1

SLIP ON BACK UPHOLSTERY		LENGTH
		inches (cm)
		18-21 (45.7-53.3)
Width Inches (cm)	16 (40.6)	9001832
	18 (45.7)	9001842
	20 (50.8)	9001852

### Move FRONT Seat-to-floor Height

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2	13 1/2	14 1/2												
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1															
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1							15 1/2	16 1/2	17 1/2	16	17	18			
	P2				15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3				16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +- 1/4"

Table 2

ADJUSTABLE TENSION BACK UPHOLSTERY		LENGTH
		inches (cm)
		18-21 (45.7-53.3)
Width Inches (cm)	16 (40.6)	9001777
	18 (45.7)	9001787
	20 (50.8)	9001797

Notes :