



MOVE

ORDER FORM - CANADA

P.O. number: * _____ Date: _____

Company Name: * _____ Phone: * _____

Ordered by: * _____

Clinician/ATP/RTS: _____

Email: _____

Marked for: _____

Ship to: * _____

Name: _____

Address: _____

City: _____ Prov: _____

Postal Code: _____

Phone: _____

* Required

MVMA050 English Manual (STD) MVMA051 Manuel Français Taxes not included

MODEL							
1	MVWC01	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)		WM0001196	\$1 395	<input type="checkbox"/>	
2	MVWC09	Heavy duty kit (350 lb weight capacity)		WAMF	\$775	<input type="checkbox"/>	
SEAT WIDTH ³							
3	MVSW	Seat width:	14" <input type="checkbox"/> 16" <input type="checkbox"/> 18" ¹ <input type="checkbox"/>			NC	<input type="checkbox"/>
4	MVSW	Seat width:	20" ¹ <input type="checkbox"/> 22" ² <input type="checkbox"/>		WAMA	\$295	<input type="checkbox"/>
SEAT DEPTH							
5	MVSD	Seat depth:	16" <input type="checkbox"/> 18" <input type="checkbox"/>			NC	<input type="checkbox"/>
6	MVSD	Seat depth:	20" <input type="checkbox"/>		WAMB	\$225	<input type="checkbox"/>

1 - Available with HD option

2 - Only available with HD option

3. [Width calculator](#)

Move FRONT Seat-to-floor Height - Please select your choice

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2	13 1/2	14 1/2												
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1															
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1							15 1/2	16 1/2	17 1/2	16	17	18			
	P2				15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3				16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
7	MVCA01	Composite wheel w/ PU tire (1" width)	4" <input type="checkbox"/>	5" <input type="checkbox"/>	6" <input type="checkbox"/>	7" <input type="checkbox"/>	8" <input type="checkbox"/>			NC	<input type="checkbox"/>	
8	MVCA10	Pneumatic 6" x 1 1/4"								\$160	<input type="checkbox"/>	
9	MVCA11	Pneumatic 8" x 1 1/4"								\$160	<input type="checkbox"/>	
10	MVCA14	Newton UltraCaster Soft-Roll Composite	4" x 1 1/2"								\$98	<input type="checkbox"/>
11	MVCA15		5" x 1 1/2"								\$98	<input type="checkbox"/>
12	MVCA16		6" x 1 1/2"								\$98	<input type="checkbox"/>
13	MVCA17		4" x 1 1/2"								\$175	<input type="checkbox"/>
14	MVCA18		Newton UltraCaster Soft-Roll Aluminum	5" x 1 1/2"								\$185
15	MVCA19	6" x 1 1/2"								\$195	<input type="checkbox"/>	

Move REAR Seat-to-floor Height - Please circle your choice

20"	22"	24"	25"	26"
12 1/4"	13"	14"	14 1/4"	15"
13 1/4"	14"	15"	15 1/4"	16"
14 1/4"	15"	16"	16 1/4"	17"
15 1/4"	16"	17"	17 1/4"	18"
16 1/4"	17"	18"	18 1/4"	19"
17 1/4"	18"	19"	19 1/4"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

SWING-AWAY FOOTREST										
16	MVFR01	60° swing away footrest	See line 18 (5 - 13.5") <input type="checkbox"/>	13.5 - 17.5" <input type="checkbox"/>	15.5 - 19.5" <input type="checkbox"/>	17 - 21" <input type="checkbox"/>	Set at:*		\$250	<input type="checkbox"/>
17	MVFR02	70° swing away footrest	See line 18 (5 - 13.5") <input type="checkbox"/>	13.5 - 17.5" <input type="checkbox"/>	15.5 - 19.5" <input type="checkbox"/>	17 - 21" <input type="checkbox"/>	Set at:*		\$250	<input type="checkbox"/>
FOOTREST OPTIONS										
18	MVFR09	High mount (footplate attached to the hanger with clamps)	Range: 5 - 13.5"			Set at:*		\$90	<input type="checkbox"/>	

* +/- 1/2"

ELEVATING LEGREST (unit)									
<i>The elevating legrest replace the footrest side selected</i>									
19	MVFR10	Elevating legrest with calf pad, 14" - 19" Length adjustable	Left <input type="checkbox"/>	Right <input type="checkbox"/>	Set at:	WAMK	\$390	<input type="checkbox"/>	

FOOTPLATES									
20	MVFP01	Standard - Composite						STD	<input type="checkbox"/>
21	MVFP02	Newton adjustable angle - Composite					WAMJ	\$145	<input type="checkbox"/>
22	MVFP03	Newton adjustable angle - Aluminum					WAMJ	\$165	<input type="checkbox"/>

FOOTPLATE ACCESSORIES									
23	MVFP08	Length adjustable heel loop						\$25	<input type="checkbox"/>
24	MVFP09	Calf strap					SEFND0010	\$52	<input type="checkbox"/>
25	MVFP30	Bodypoint Aeromesh Padded calf strap					SEFND0010	\$129	<input type="checkbox"/>

RESIDUAL LIMB SUPPORT (Height and angle adjustable) (unit)									
26	MVFR23	Swing away residual limb support - 7" x 9"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			SEFND0005	\$390	<input type="checkbox"/>
27	MVFR12	Swing away residual limb support - 8" x 11"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			SEFND0005	\$390	<input type="checkbox"/>

Back canes with integrated fixed height push handle

FIXED ANGLE BACK										
28	MVBC01	8° bend aluminum back cane	16" <input type="checkbox"/>	18" <input type="checkbox"/>	20" <input type="checkbox"/>	STD			NC	<input type="checkbox"/>
ADJUSTABLE ANGLE BACK										
29	MVBC02	8° bend aluminum back cane		16" <input type="checkbox"/>	20" <input type="checkbox"/>	STD			\$309	<input type="checkbox"/>
30		Select back cane angle:	90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>			

STROLLER HANDLE AND STABILIZER BAR ¹									
31	MVPH04	Stroller handle (not available with the HD option)					WAMW	\$210	<input type="checkbox"/>
32	MVPH03	Folding stabilizer bar						\$258	<input type="checkbox"/>

¹ - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY						
33	MVBH05	Slip-on Nylon back upholstery			NC	<input type="checkbox"/>
34	MVBH06	Adjustable tension Nylon back upholstery			WAMD/WAME	\$245 <input type="checkbox"/>
35	MVBH99	Omit back upholstery			STD	<input type="checkbox"/>
RIGID BACK						
36		nXt back series (See nXt back pricelist)			See pricelist	<input type="checkbox"/>

SEAT CUSHIONS						
37	MVSC01	2" Cushion			\$175	<input type="checkbox"/>
38	MVSC02	3" Cushion			\$195	<input type="checkbox"/>
39	MVSC99	Omit seat cushion			NC	<input type="checkbox"/>
40		nXt FIT cushion series (See nXt FIT pricelist)			See pricelist	<input type="checkbox"/>

REAR WHEELS									
41	MVRW01	Mag ¹	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	NC	<input type="checkbox"/>		
42	MVRW02	Newton One - spoke wheel ²	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	NC	<input type="checkbox"/>		
43	MVRW03	Newton Gravity - ultralight wheel	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>	\$360	<input type="checkbox"/>

1 - HP tires are not compatible with Mag Wheels -use only aluminum or plastic coated handrim 2- Not compatible with Shox tires

TIRES						Available	
44	MVTI01	Soft Urethane 1 1/4 "			20", 22" 24"	STD	<input type="checkbox"/>
45	MVTI02	Soft Urethane knobby 1 3/8 "			22", 24", 25", 26"	\$130	<input type="checkbox"/>
46	MVTI04	Soft Urethane - Shox 1"			22", 24", 25", 26"	\$155	<input type="checkbox"/>
47	MVTI05	Pneumatic 1 3/8 "			20", 22", 24", 25", 26"	\$130	<input type="checkbox"/>
48	MVTI08	Pneumatic with Airless insert 1 3/8"			22", 24"	\$215	<input type="checkbox"/>
49	MVTI06	Pneumatic HP - 1"			20", 22", 24", 25", 26"	\$150	<input type="checkbox"/>
50	MVTI07	Pneumatic HP - Schwalbe Marathon Plus 1"			24", 25", 26"	\$230	<input type="checkbox"/>

HANDRIM							
51	MVHR01				STD	<input type="checkbox"/>	
52	MVHR07	Newton Air Grip ¹			WAMO	\$150 <input type="checkbox"/>	
53	MVHR08	Newton Air Grip ¹	(8) Orange Sunkissed	(9) Burgundy <input type="checkbox"/>	(10) Sapphire blue <input type="checkbox"/>	WAMO	\$330 <input type="checkbox"/>
		(11) Monster green <input type="checkbox"/>	(12) Charcoal <input type="checkbox"/>	(13) Fuchsia <input type="checkbox"/>	(14) Steel blue <input type="checkbox"/>		
		(15) Ferrari red <input type="checkbox"/>	(16) Acid green <input type="checkbox"/>	(17) White <input type="checkbox"/>	(18) Sky blue <input type="checkbox"/>		
Assembly position for aluminum and Newton AirGrip handrims							
54	MVHR29	Standard Handrim position		Full tab	STD	<input type="checkbox"/>	
55	MVHR30	Narrow Handrim position (not available in 22")		Short tab 2 positions	\$25	<input type="checkbox"/>	
56	MVHR31	Super Narrow Handrim (not available in 22")		Short tab 2 positions Short tab 1 position <input type="checkbox"/>	\$25	<input type="checkbox"/>	

HANDRIM						
57	MVHR02	W	^{2, 3}	WAMO	\$160	<input type="checkbox"/>
58	MVHR03	Q-Grip	^{1,2,3}	WAMO	\$430	<input type="checkbox"/>
59	MVHR28	Surge ^{1, 2, 3}	(1-1/2")	\$470	<input type="checkbox"/>
60	MVHR05	Surge LT ^{1,2,3}	(1-1/8")	\$500	<input type="checkbox"/>
61	MVHR04	Natural Fit Standard Grip ^{1,2,3}	(1-1/2")	\$500	<input type="checkbox"/>
62	MVHR25	Natural Fit Super Grip ^{1,2,3}	(1-1/2")	\$530	<input type="checkbox"/>
63	MVHR99				NC	<input type="checkbox"/>

1- 2 - 3-


ONE ARM DRIVE								
64	MVHR22	One arm drive ⁴	24"	left <input type="checkbox"/>	right <input type="checkbox"/>	New	\$1 395	<input type="checkbox"/>

4 - Only available with Newton One wheels, aluminum or plastic handrims

AXLES (0° camber)							
65	MVAX01	Quick release axle			WAMU	STD	<input type="checkbox"/>
66	MVAX02	Permanent axle				NC	<input type="checkbox"/>
67	MVAX03	Quad release axle			WAMU	\$170	<input type="checkbox"/>
AXLE PLATE							
68	MVAX04	Amputee adjustable axle plate ¹			WAMT	\$250	<input type="checkbox"/>

1 - May limit seat to floor height available

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable	Most responsive
AMPUTEE	3/4" 2"
<input type="checkbox"/>	<input type="checkbox"/>

WHEEL LOCKS						
69	MVWL08	Newton Wheel lock - push to lock			STD	<input type="checkbox"/>
70	MVWL02	Pull to lock			NC	<input type="checkbox"/>
71	MVWL04	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>		WAMR	\$175	<input type="checkbox"/>
72	MVWL13 MVWL14	6" removable extension handle			\$67	<input type="checkbox"/>
73	MVWL06	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>			\$230	<input type="checkbox"/>
74	MVWL99	Omit wheel locks			NC	<input type="checkbox"/>

FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ¹		Height Adjustment		Removable function				
75	MVAR01	10" Desk length arm pad	8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/> STD	\$380	<input type="checkbox"/>
76	MVAR02	14" Full length arm pad	8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>	Locked <input type="checkbox"/>	Unlocked <input type="checkbox"/> STD	\$380	<input type="checkbox"/>
"T" ARMREST ¹		Height Adjustment						
77	MVAR03	10" Desk length arm pad	8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>			\$360	<input type="checkbox"/>
78	MVAR04	14" Full length arm pad	8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>			\$360	<input type="checkbox"/>

¹ - May have a 1/2" adjustment restriction on both limits of the range

GEL ARM PAD						
79	MVAR12	14" Full Length Arm Pad - Gel Ovation		SEAND0030	\$150	<input type="checkbox"/>
80	MVAR99	Omit armrest			NC	<input type="checkbox"/>

POSITIONING BELTS								
81	MVSB01	Velcro adjustable	Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>		\$45	<input type="checkbox"/>	
82	MVSB02	Auto buckle	Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>		\$55	<input type="checkbox"/>	
83	MVSB07	Hardware attachment for 2 point seatbelt (1 pair)				NC	<input type="checkbox"/>	
84	MVSB03	Additional hardware for 2 extra attachment points <i>(needed for the 4 point seat belt) (1 pair)</i>				\$13	<input type="checkbox"/>	
85	MVSB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>		\$82	<input type="checkbox"/>	
86	MVSB05	BodyPoint Hip Belt - 2 point padded	Small 51" <input type="checkbox"/>	Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>	SERND0030	\$170	<input type="checkbox"/>
87	MVSB06	BodyPoint Hip Belt - 4 point padded	Small 51" <input type="checkbox"/>	Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>	SERND0020	\$195	<input type="checkbox"/>
88	MVSB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" <input type="checkbox"/>	Long 24" <input type="checkbox"/>		\$324	<input type="checkbox"/>	
89	MVSB99	Omit positioning belt				NC	<input type="checkbox"/>	

HEADREST ¹								
90	MVAC01	Profiled Headrest	Small (3.5"X6") <input type="checkbox"/>	Medium (4.5"X7") <input type="checkbox"/>	Large (5.5"X8") <input type="checkbox"/>		\$205	<input type="checkbox"/>
91	MVAC02	Headrest mounting hardware with stabiliser bar					\$338	<input type="checkbox"/>

¹ - Not available with a stroller handle or a folding stabiliser bar


ACCESSORIES											
92	MVAC03	Anti tippers								\$115	<input type="checkbox"/>
93	MVAC27	Tip assist ¹							New	\$60	<input type="checkbox"/>
94	MVAC04	Transit Tie Down (Unoccupied)								NC	<input type="checkbox"/>
95	MVAC05	WC-19 Transit Tie Down w/Q'straint belt ²		Short belt 20-41" <input type="checkbox"/>		Long belt 20-46" <input type="checkbox"/>			New	\$345	<input type="checkbox"/>
96	MVAC06	Cane carrier		Left <input type="checkbox"/>		Right <input type="checkbox"/>				\$105	<input type="checkbox"/>
97	MVAC07	Oxygen tank holder - E Type <i>Not compatible with swing away armrest or amputee axle plate</i>		Left						WAMX	\$375 <input type="checkbox"/>
98	MVAC08	Spoke guards ³		Available: 20", 22", 24", 26"						WAMP	\$145 <input type="checkbox"/>
99	MVAC09	Footrest neoprene impact guard (pair) - Black								SEFND0035	\$80 <input type="checkbox"/>
100	MVAC10	Seat pouch									\$95 <input type="checkbox"/>
101	MVAC26	BodyPoint Mobility Bag									\$135 <input type="checkbox"/>
102	MVAC23	Newton Backpack									\$120 <input type="checkbox"/>
103	MVAC12	Velcro Closing Strap									\$50 <input type="checkbox"/>
104	MVAC13	Tool kit									\$50 <input type="checkbox"/>
105	MVAC16	Motion Composites Polo shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>						\$32 <input type="checkbox"/>
			M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>	M/XXXXL <input type="checkbox"/>		

1 - Cannot be combined with Anti Tippers 2 - Not available with the swing away armrest ; Not available with HD option 3 - Not available with Mag wheels

FRAME COLOR (with black aluminum crossbrace)												
106	MVFC01	Black <input type="checkbox"/>								NC	<input type="checkbox"/>	
107	MVFC03	Sunkissed Orange <input type="checkbox"/>	Burgundy <input type="checkbox"/>	Sapphire Blue <input type="checkbox"/>	Monster Green <input type="checkbox"/>					\$270	<input type="checkbox"/>	
		Charcoal <input type="checkbox"/>	Fuchsia <input type="checkbox"/>	Steel Blue <input type="checkbox"/>	Ferrari Red <input type="checkbox"/>	Acid Green <input type="checkbox"/>						
108	MVFC99	Custom Color	Specify PPG Paintit color:								\$410	<input type="checkbox"/>

TOTAL

Fork and Stem bolt combinations



Notes :