



# INSTALLATION INSTRUCTIONS


## RIGID SEAT PAN FOR RIGID WHEELCHAIR

COMPATIBLE WITH RIGID WHEELCHAIR

Material(s) required:

- Fixed alu seat pan
- Universal 1.25in alu clamp
- Button cap screw m6x1x30mm

1. Fix the anchors with a standard torque of 4.5Nm. Refer to the configuration of the wheelchair according to the available space for their height locations.



Note: the assembly must be delivered to the customer with their rigid seat pan velcro hook sticker kit. Make sure you have the kit in hand.