



1- Purchasing information

P.O. number*: _____ Date: _____
Company Name*: _____
Purchasing contact*: _____
Phone*: _____
Email*: _____
Contact preference*: Phone ☐ Email ☐
Marked for: _____

3- Shipping information

Ship to*: _____
Name: _____
Address: _____
City: _____ State: _____
Zip Code: _____
Phone: _____

2- d W I Clinician information

d W I Clinician*: _____
Phone*: _____
Email*: _____
Contact preference*: Phone ☐ Email ☐

Required		English Manual (STD)						Spanish Manual		Taxes not included
MODEL								USD		
WC03	☐	Ultralight folding Helio wheelchair (265lb weight capacity) A7 packageside frames, symmetrical hydroformed crossbrace and back canes made of 7000 series							K0005	\$2 295
WC08	☐	C2 carbon fiber upgrade side frames, symmetrical molded crossbrace and back canes made of ultralight carbon fiber								\$1 835
WC09	☐	HDHeavy Duty package (350 lb weight capacity)								\$525
WC10	☐	E-Motion adaptor Kit (2 adaptors)Only necessary with 6° camber								\$35
SEAT WIDTH ³										
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐			NC
	☐	20" ¹ ☐ 21" ² ☐ 22" ² ☐							E2201	\$495
SEAT DEPTH										
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐			NC
	☐	20" ☐							E2203	\$475

1 - Available with HD Kit (350 lb)

2 - Only with HD Kit (350 lb)

[3. Width calculator](#)

NEWTON ULTRALIGHT PACKAGE				Save \$225	\$490
PK01	☐	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels			

Helio FRONT Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2				14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						
	P4	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
7"	P1										16 1/2"	17 1/2"	18 1/2"	17"	18"	19"			
	P2							16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"			
	P3							17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
	P4	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"
	P5	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"	19 1/2"	20 1/2"	21 1/2"

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS						
CA01-03	☐	Composite wheel w/ PU tire	3" x 1"			NC
CA01-04	☐		4" x 1"			NC
CA01-05	☐		5" x 1"			NC
CA01-06	☐		6" x 1"			NC
CA01-07	☐		7" x 1"			NC
CA01-08	☐		8" x 1"			NC
CA10	☐	Pneumatic casters	6" x 1 ¼"		E2214/E2215	\$160
CA11	☐		8" x 1 ¼"		E2214/E2215	\$160
CA14	☐	Newton UltraCaster Composite	4" x 1 ½ "			\$85
CA15	☐		5" x 1 ½ "			\$85
CA16	☐		6" x 1 ½ "			\$85
CA17	☐	Newton UltraCaster Aluminum	4" x 1 ½"			\$185
CA18	☐		5" x 1 ½"			\$195
CA19	☐		6" x 1 ½"			\$205
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (available with +1" or +2" Stem bolt ONLY)			K0073	\$295
CA13	☐	Carbon fiber Frog Legs suspension forks (not available with caster pin locks)			E1015	\$475

1- Not available with WC-19; 2- Not available with 7" and 8" casters; limited height of 18.5" with 6" casters

Helio - REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height: _____
20" (451) wheel	☐	12 1/2	to 17 1/4	
22" (501) wheel	☐	13 1/2	to 18 1/4	
24" (540) wheel	☐	14 1/2	to 19 1/4	
25" (559) wheel	☐	15	to 19 3/4	
26" (590) wheel	☐	15 1/2	to 20 1/4	

E} \$ W d] E •] | U v \$ E } (P E A } \$ Ç • \$ \$ } v P v • \$ • o } % u Ç E • µ o \$ } v A E] \$ } v } (= | r * _ } (• o \$ • \$ \$ } (o)

Footrest					
SWING-IN/SWING-OUT FOOTREST					
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"		Set at: *	NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at: *	NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at: *	NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"		Set at: *	NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at: *	NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at: *	NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16"		Set at: *	NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at: *	NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at: *	NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension)	7.5" - 11"	Set at: *	NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension ^{1,2}	9.5"- 13"	Set at: *	NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension ^{1,2}	9.5" - 14.5"	Set at: *	NC
FOOTREST OPTIONS					
FR09	☐	High mount (footplate attached to the hanger with clamps) ³		Set at: *	\$80
FR10	☐	Elevating legrest with calf pad, Length adjustable Left ☐ Right ☐ 14" - 22"		Set at:	E0990 \$240
FR99	☐	Omit Footrest			NC

1- Available only with two pieces adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount; 4- The elevating legrest replace the footrest side selected; 1-1/2"

FOOTPLATES					
FP01	☐	Standard - Composite			STD
FP02	☐	Newton adjustable angle - Composite		K0040	\$210
FP03	☐	Newton adjustable angle - Aluminum		K0040	\$230
FP04	☐	Newton adjustable angle - Carbon fiber (with rubber inserts)			\$435
FP06	☐	One piece adjustable angle flip-up - Aluminum		K0040	\$190
FP07	☐	One piece adjustable angle flip-up - Carbon fiber (with rubber inserts)			\$455
FOOTPLATE ACCESSORIES					
FP08	☐	Length adjustable heel loop (not available with one piece footplate)		E0951	\$40
FP09	☐	Calf strap		K0038	\$35
FP30	☐	Bodypoint Aeromesh Padded calf strap		K0038	\$115

1- Not available with front mount

SWING AWAY RESIDUAL LIMB SUPPORT							
FR27	☐	Residual limb support - flat 6"x8"	Left ☐	Right ☐		E1020	\$405
FR28	☐	Residual limb support - flat 8"x10"	Left ☐	Right ☐		E1020	\$405
FR29	☐	Residual limb support - flat 10"x1	Left ☐	Right ☐		E1020	\$405
FR30	☐	Residual limb support - flat 10"x1	Left ☐	Right ☐		E1020	\$405
FR32	☐	Residual limb support - contour 6"x8"	Left ☐	Right ☐		E1020	\$405
FR33	☐	Residual limb support - contour 8"x10"	Left ☐	Right ☐		E1020	\$405

Backrest

FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE (not available with HD option)							
BC60	☐	Straight aluminum back cane 21" - 24"	increment of 1"		Set at ² :	New	\$ ó ñ
BC59	☐	8° bend aluminum back cane 21" - 24"	increment of 1"		Set at ² :	New	\$ ó ñ

FIXED ANGLE AND HEIGHT ADJUSTABLE BACK WITH PUSH HANDLE OPTION									
BC11	☐	Straight back cane*	12 - 15" ☐	15 - 18" ☐	18 - 21"STD ☐	increment of 1/2"	Set at²:		NC
BC06									
BC12	☐	8° Bend back cane*	13 ½ - 15" ☐	15 - 18" ☐	18 - 21"STD ☐	increment of 1/2"	Set at²:		\$55
BC08									

*** fixed height Carbon fiber back canes are unavailable until the end of 2022. If this option is selected, Aluminum back canes will be installed. Otherwise, carbon fiber angle and height adjustable are still available. See options below**

ANGLE AND HEIGHT ADJUSTABLE BACK WITH PUSH HANDLE OPTION								
BC13	☐	Straight back cane	15 - 18" ☐	18 - 21"STD ☐	increment of 1/2"	Set at²:		\$210
BC09								
BC14	☐	8° bend back cane	15 - 18" ☐	18 - 21"STD ☐	increment of 1/2"	Set at²:		\$265
BC10								
	Select back cane angle:		85°¹ ☐	90° ☐	95°STD ☐	100° ☐	105° ☐	110° ☐

PUSH HANDLE OPTIONS							
PH01	☐	Push handles					STD
PH05	☐	Fold down push handles ^{3,5}					\$126
PH06-920	☐	Clamp on height adjustable push handles 9" - 20"					\$315
PH06-923	☐	Clamp on height adjustable push handles 9" - 23"					\$315
PH02	☐	Omit push handles					NC

STROLLER HANDLE AND STABILIZER BAR							
PH04	☐	Stroller handle (only available with push handle) - (not available with the HD option)				K0108	\$160
PH03	☐	Newton folding stabilizer bar (not compatible with stroller handle)					\$265

¹- Not available with a 8° bend back cane; If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery
²- Not available with swing away armrests, 8° bend back canes and amputee axle plates; ³- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock
⁴- Not compatible with swing away armrests, 8° bend back canes and amputee axle plates; ⁵- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY							
BH07	☐	Slip-on Active-Mesh back upholstery					STD
BH08	☐	Adjustable tension Active-Mesh back upholstery				E2611	\$300
BH99	☐	Omit Back upholstery					NC
	☐	NXT back Support				See pricelist	

SEAT CUSHION							
SC01	☐	2" Cushion				E2601	\$165
SC02	☐	3" Cushion				E2601	\$185
SC99	☐	Omit seat cushion					STD
	☐	NXT cushion				See pricelist	

Rear wheels					
ONE ARM DRIVE			Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims		
HR22	☐	One arm drive ^{2 3}	24"	left ☐	right ☐
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)			NC
REAR WHEELS			Available		
RW01	☐	Mag ⁴	20", 22" 24"		\$55
RW02	☐	Newton One - spoke wheel	20", 22" 24"		NC
RW03	☐	Newton Gravity - ultralight wheel	20", 22", 24", 25", 26"		\$315
RW04	☐	Spinergy Spox Blk STD ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐	22", 24", 25", 26"		\$835
RW05	☐	Spinergy LX Blk STD ☐ Pur ⁶ ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐ Pink ⁶ ☐ Grn ⁶ ☐ Org ⁶ ☐	22", 24", 25", 26"		\$1 045
RW07	☐	Spinergy CLX Blk STD ☐ Pur ⁶ ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐ Pink ⁶ ☐ Grn ⁶ ☐ Org ⁶ ☐	24", 25"		\$1 555
TIRES			Available		
TI01	☐	Soft Urethane 1 3/8"	20", 22" 24"	E2218	STD
TI02	☐	Soft Urethane knobby1 3/8 "	22", 24", 25", 26"	E2218	\$95
TI04	☐	Soft Urethane Shox 1"	22", 24", 25", 26"	E2218	\$115
TI05	☐	Pneumatic 1 3/8"	20", 22", 24", 26"	E2211/E2212	NC
TI08	☐	Pneumatic with Airless insert - 1 3/8"	22", 24", 26"	E2211/E2213	\$80
TI06	☐	Pneumatic HP 1"	20", 22", 24", 25", 26"	E2211/E2212	\$95
TI07	☐	Pneumatic HPPuncture resistant, black 1"	24", 25", 26"	E2211/E2212	\$125

1- Offered with handrims on each side of the wheelchair 2- Only standard handrim position available 3 - Not available with HD option
4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 5 - Not compatible with Shox tires
6 - Color spokes require an additionnal 10 days lead-time 7- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims					
HANDRIM			Available		
HR01	☐	Aluminum anodized	20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized	20", 22", 24", 25"		\$70
HR02	☐	Plastic coated	20", 22", 24", 25", 26"		\$100
HR03	☐	Q-Grip High friction Neoprene coated	22", 24", 25", 26"		\$370
HR07	☐	Newton Air Grip ²	20", 22", 24", 25", 26"		\$100
HR	☐	Newton Air Grip ²	20", 22", 24", 25", 26"		\$250
HR28	☐	Surge ^{1,2} (1-1/2" Oval handrim w/ Gripton high friction strip)	**Limited availability	22", 24", 25", 26"	\$395
HR05	☐	Surge LT ^{1,2} (1-1/8" Oval handrim w/ Gripton high friction strip)	**Limited availability	22", 24", 25", 26"	\$425
HR04	☐	Natural Fit ^{1,2}	**Limited availability	22", 24", 25", 26"	\$325
HR99	☐	Omit handrim			NC

1- Only standard handrim position available 2- Not available with Mag wheels or OAD

Handrims option					
Assembly position for aluminum and Newton AirGrip handrims					
HR29	☐	Standard Handrim position	Full tab		STD
HR30	☐	Narrow Handrim position	Short tab 2 positions		\$27
HR31	☐	Super Narrow Handrim			\$27
Thumb piece selection for Natural Fit handrims					
HR23	☐	Standard grip	**Limited availability		\$90
HR27	☐	Super grip	**Limited availability		\$145
HR099	☐	Omit thumb piece			NC

3- Compatible with Surge, Surge LT handrims

AXLE					
AX01	☐	Quick release axle			STD
AX02	☐	Permanent axle			NC
AX03	☐	Quad release axle		K0108	\$125
AXLE PLATE					
AX04	☐	Amputee adjustable axle plate		E0959	\$195

1- Not required with one arm drive 2 - May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER					
CM00	☐	0° camber			STD
CM03	☐	3° camber			NC
CM06	☐	6° camber(not available with amputee axle plate)			\$160

CENTER OF GRAVITY

Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													Most responsive
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"	
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	

X = _____

If not indicated, default setting is at 1.5"

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

May interfere with some COG and other configurations. For more details, refer to the table at page 7 and the "How-to documents" page of our website: https://www.motioncomposites.com/en_us/support-and-education/how-to-documents

Armrests

HEIGHT ADJUSTABLE "T" ARMPAD ¹					
AR13-610	☐	10" Desk length armpad 6" - 10"	<i>New</i>	E0973	\$325
AR03-812	☐	10" Desk length armpad 8" - 12"		E0973	\$325
AR03-1014	☐	10" Desk length armpad 10" - 14"		E0973	\$325
AR14-610	☐	14" Full length armpad 6" - 10"	<i>New</i>	E0973	\$325
AR04-812	☐	14" Full length armpad 8" - 12"		E0973	\$325
AR04-1014	☐	14" Full length armpad 10" - 14"		E0973	\$325

ARMREST OPTIO					
AR12	☐	14" Full Length armrest - Gel Ovation			\$160
AR99	☐	Omit armrest			NC

1- May have 1/2" adjustment limit on both limits of the range; Could interfere with some rigid backrests; Not available with swing away and flip back armrest; Not available with attendant lock

5- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (Not available with Flip Back and "T" Armrest)					
SG01	☐	Plastic side guards (removable)		K0108	\$195
SG02	☐	Carbon fiber side guards (fixed)			\$265
SG03	☐	Carbon fiber side guards (removable)		K0108	\$290
SG06	☐	Plastic fender side guards (removable)			\$285
SG04	☐	Carbon fiber fender side guards (removable)		K0108	\$540
SG99	☐	Omit side guard			STD

SIDE GUARDS OPTION					
SG07	☐	Side guards clip (pair) ¹			NC

1- Not available with fixed carbon fiber side guard

POSITIONING BELT					
SB01-S	☐	Velcro adjustableStandard 35"		E0978	\$45
SB01-L	☐	Velcro adjustableLong 60"		E0978	\$45
SB02-S	☐	Auto buckle Standard 35"	Improved	E0978	\$85
SB02-L	☐	Auto buckle Long 60"	Improved	E0978	\$85
SB07	☐	Hardware attachment for 2 point seatbelt (1 pair)			NC
SB03	☐	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)			\$16
SB04-M	☐	BodyPoint Hip Belt - 2 point not paddedMedium 56"		E0978	\$80
SB04-L	☐	BodyPoint Hip Belt - 2 point not paddedLong 61"		E0978	\$80
SB05-S	☐	BodyPoint Hip Belt - 2 point paddedSmall 51"		E0978	\$170
SB05-M	☐	BodyPoint Hip Belt - 2 point paddedMedium 56"		E0978	\$170
SB05-L	☐	BodyPoint Hip Belt - 2 point paddedLong 61"		E0978	\$170
SB06-S	☐	BodyPoint Hip Belt - 4 point paddedSmall 51"		E0978	\$190
SB06-M	☐	BodyPoint Hip Belt - 4 point paddedMedium 56"		E0978	\$190
SB06-L	☐	BodyPoint Hip Belt - 4 point paddedLong 61"		E0978	\$190
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps-18"		K0108	\$300
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps-24"		K0108	\$300
SB99	☐	Omit positioning belt			STD

HEADREST					
AC01-S	☐	Profiled HeadrestSmall (3.5"X6")		E0955	\$180
AC01-M	☐	Profiled HeadrestMedium (4.5"X7")		E0955	\$180
AC01-L	☐	Profiled HeadrestLarge (5.5"X8")		E0955	\$180
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar		E1028	\$295

1 - Only available with push handles - Not available with a stroller handle or a folding stabilizer bar


ACCESSORIES					
AC03	☐	Anti-tippers		E0971	\$140
AC29	☐	Swing-away anti-tipper Left ☐ Right ☐			\$160
AC27	☐	Tip assist			\$65
AC04	☐	Transit Tie Down(Unoccupied)			NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt Short belt 20"-41"			\$115
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt Long belt 20"-46"			\$115
AC35	☐	WC-19 Transit Tie Down without belt	New		\$315
AC06	☐	Cane and crutch holder Left ☐ Right ☐		E2207	\$195
AC07	☐	Oxygen tank holder - D and E Type Left Not compatible with swing away armrest or amputee axle plate		E2208	\$310
AC08	☐	Spoke guards Available: 20", 22", 24", 26"		K0065	\$115
AC09	☐	Footrest neoprene impact guard (pair) - Black		K0108	\$90
AC10	☐	Seat pouch			\$80
AC30	☐	Newton multi-purpose pouch			\$60
AC26	☐	BodyPoint Mobility Bag			\$130
AC23	☐	Newton Backpack			\$125
AC12	☐	Velcro Closing Strap			\$55
AC13	☐	Tool kit			\$50
AC16	☐	Motion Composites Polo Shirt W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐			\$35
AC15	☐	Military Patches Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐ Combat Wounded Purple Heart ☐			\$95

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" a 2 - Cannot be combined with Anti Tippers

3 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD optic 4 - Not available with Mag wheels and Spinergy

FRAME COLOR					
FC01	☐	Black/Carbon Fiber Black	Black color is replaced by the Carbon fiber finish with the C2 upgrade)		NC
FC	☐	(03) Sunkissed orange ☐ (04) Burgundy ☐ (05) Sapphire Blue ☐ (06) Monster Green ☐ (07) Charcoal ☐ (08) Fuchsia ☐ (09) Steel Blue ☐ (17) Ferrari Red ☐ (18) Acid Green ☐			NC
FC99	☐	Custom Color Specify PPG Paintit color:	Additional 10 days lead-time		\$380


Fork and Stem bolt combinations



Unilateral wheel lock restrictions

Seat depth must be 16in and greater		
Rear wheel diameter 24 in and less		
"T" or flip-back armrest only available with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Wheel lock lever extension not available		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
22 in	16	1.0 to 2.75
	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :

BioFit™ CUSHION		WIDTH										
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
 NXCAF_____ (W" x D") HCPCS CODE: E2607-E2608	DEPTH	12"	☐ \$425	☐ \$425	☐ \$425							
		14"	☐ \$425	☐ \$425	☐ \$425							
		15"		☐ \$425	☐ \$425	☐ \$425	☐ \$425					
		16"		☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	☐ \$455	☐ \$455
		17"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	
		18"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	☐ \$455
		19"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	
		20"				☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	☐ \$455	☐ \$455
Optional Incontinent / Infection Control Outer Cover NXCAF____ V ☐ N/C												

THORACIC SHAPE

Optima™ THORACIC BACK SUPPORT



Optima™	Fits Width	Height	✓	\$
NOACB1613	16"-18"	13"	☐	\$660
NOACB1616	16"-18"	16"	☐	\$660
NOACB1619	16"-18"	19"	☐	\$660
NOACB1813	18"-20"	13"	☐	\$660
NOACB1816	18"-20"	16"	☐	\$660
NOACB1819	18"-20"	19"	☐	\$660
NOACB2016	20"-22"	16"	☐	\$660
NOACB2019	20"-22"	19"	☐	\$660
NOACB2213	22"-24"	13"	☐	\$710
NOACB2216	22"-24"	16"	☐	\$710
NOACB2219	22"-24"	19"	☐	\$710

HCPCS CODE:
E2620-E2621

Optima™ DEEP THORACIC BACK SUPPORT



Optima™ DEEP	Fits Width	Height	✓	\$
NOPCB1616	16"-18"	16"	☐	\$660
NOPCB1619	16"-18"	19"	☐	\$660
NOPCB1816	18"-20"	16"	☐	\$660
NOPCB1819	18"-20"	19"	☐	\$660
NOPCB2016	20"-22"	16"	☐	\$660
NOPCB2019	20"-22"	19"	☐	\$660
NOPCB2213	22"-24"	13"	☐	\$710
NOPCB2216	22"-24"	16"	☐	\$710
NOPCB2219	22"-24"	19"	☐	\$710

HCPCS CODE:
E2620-E2621

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT



Xtend™	Fits Width	Height	✓	\$
NHAACB1616/20	16"-18"	16"-20"	☐	\$660
NHAACB1816/20	18"-20"	16"-20"	☐	\$660
NHAACB2016/20	20"-22"	16"-20"	☐	\$660
NHAACB2216/20	22"-24"	16"-20"	☐	\$710

HCPCS CODE:
E2620-E2621

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT



Xtend™ DEEP	Fits Width	Height	✓	\$
NHAPCB1616/20	16"-18"	16"-20"	☐	\$660
NHAPCB1816/20	18"-20"	16"-20"	☐	\$660
NHAPCB2016/20	20"-22"	16"-20"	☐	\$660
NHAPCB2216/20	22"-24"	16"-20"	☐	\$710

HCPCS CODE:
E2620-E2621

Classic™ PELVIC BACK SUPPORT



Classic™	Fits Width	Height	✓	\$
NACB1616	16"-18"	16"	☐	\$660
NACB1618	16"-18"	18"	☐	\$660
NACB1816	18"-20"	16"	☐	\$660
NACB1818	18"-20"	18"	☐	\$660
NACB1820	18"-20"	20"	☐	\$660
NACB2020	20"-22"	20"	☐	\$660
NACB2216	22"-24"	16"	☐	\$710
NACB2218	22"-24"	18"	☐	\$710
NACB2220	22"-24"	20"	☐	\$710

HCPCS CODE:
E2620-E2621

Classic™ DEEP PELVIC BACK SUPPORT



Classic DEEP™	Fits Width	Height	✓	\$
NPCB1616	16"-18"	16"	☐	\$660
NPCB1618	16"-18"	18"	☐	\$660
NPCB1818	18"-20"	18"	☐	\$660
NPCB1820	18"-20"	20"	☐	\$660
NPCB2020	20"-22"	20"	☐	\$660
NPCB2216	22"-24"	16"	☐	\$710
NPCB2218	22"-24"	18"	☐	\$710
NPCB2220	22"-24"	20"	☐	\$710

HCPCS CODE:
E2620-E2621

PELVIC SHAPE