

ACCEPT NOTHING LESS THAN GREATER MOBILITY

The APEX Aluminum wheelchair by Motion Composites is the perfect combination of lightness, sturdiness, performance, adjustability and style. Its bold design and innovative approach have redefined perceptions of rigid wheelchairs. The triple-butted aluminum and strategic touch of carbon fiber strike a perfect balance, making this chair an example of modernity and lightness.

World-renowned, Motion Composites uses its knowhow and its passion for innovation to design unique, cutting-edge products inspired by you. As a visionary and passionate company, our ambition is driven by your experience.

We believe that high-performance wheelchairs enhance people's lives with components, design, features and quality. APEX Aluminum is our passion for performance, our commitment to quality and our vision of greater freedom brought to life.

GO BEYOND!







LIGHTWEIGHT AND INNOVATIVE

- 10.5 lb. A lightness that has redefined perceptions of rigid wheelchairs
- AluLite: a stronger and lighter aluminum alloy

Made with the industry's most advanced materials, the APEX Aluminum wheelchair is a true innovation. It combines the advantages of triple-butted aluminum with those of carbon fiber for unbeatable strength and lightness. The result is a wheelchair that is easy to handle and transport, with an ease of propulsion that makes it easier to get to where you want to go.

AGILE AND HIGH-PERFORMANCE

Unprecedented convenience Uncompromising stability

APEX Aluminum is an easily maneuverable rigid wheelchair that offers safe and high-performance handling. Thanks to its unparalleled design, it performs like its carbon fiber counterpart, APEX Carbon. The rigidizing bar unique to the APEX model maximizes stability while preserving reactivity. This rigidity reduces lateral movement and offers optimal propulsion for greater mobility.

ADJUSTABLE AND FLEXIBLE

- Infinitely adjustable
- Adapts to your lifestyle

The APEX Aluminum wheelchair offers a wide range of settings for a completely customized experience. Fully modular, the back angle, seat-to-floor height and rear wheel positions can all be adjusted for your comfort. Laser-etched markings make it easier to adjust key components and specify exact positioning. Simple. Quick. Adaptive.

STURDY AND DURABLE

- The industry's most perfected materials
- Exceptional reliability

AluLite aluminum alloy and triple-While technology and function butted aluminum make this chair are important, so are design and considerably lighter and aesthetics. Each detail of the APEX remarkably strong. Thanks to the Aluminum was carefully considered quality of the components and and thoughtfully designed to the touch of carbon fiber. create a head-turning chair the APEX Aluminum is highly that challenges all wheelchair durable to serve you for many stereotypes. years to come.

The Apex Aluminum made such a difference in my life. Everything is so easier now.



– Martha Whitefield, Georgetown, USA

DESIGN AND STYLE

- Bold, modern and refined design
- Attention to detail and aesthetics

FEATURES

CANTELIVER FRAME

Frame assembled by a precision fixing system offering exceptional rigidity and maximizing efficiency while allowing adjustments



BACK REST ADJUSTMENT Offers wide range of adjustments to fit any position



2 STANDARD CARBON FIBER REAR RIGIDIZER BAR

Allows for optimal propulsion and frame stability



3 STANDARD CARBON FIBER CAMBER TUBE

Adds carbon fiber strength and reduces overall weigh.



TRIPLE-BUTTED TUBE Reduces weight in non-critical sections and increases strength where needed



5 BUBBLE LEVELS Highly accurate and calibrated system for a quick and perfect adjustment



FRONT BUMPER Provides style and protection



7 REAR AXLE PLATE Provides a multitude of adjustments while maximizing stability and responsiveness

TECHNOLOGIES



ALULITE ALUMINUM ALLOY Stronger and lighter than 7005 aluminum.



TRIPLE-BUTTED ALUMINUM Reduces weight while maintaining strength. We put the strength where it counts!



CLAMPED RIGIDIZING SYSTEM Patented modular system that enhances versatility and rear frame stability.



BUBBLE LEVELS Built-in highly accurate and calibrated gauges make caster adjustments quick and precise.



STANDARD CARBON FIBER CAMBER TUBE

Provides a responsive ride without added weight.



LASER ETCHED MARKINGS Allow for guick and easy symmetrical adjustments.



NEWTON ACCESSORIES Parts and accessories designed to be lighter with improved functionality.





CENTER OF GRAVITY 0 in. to 6 in. | 0 cm to 15.2 cm



SEAT DEPTH 12 in. to 20 in. | 30.5 cm to 50.8 cm 12 in. to 20 in. | 30.5 cm to 50.8 cm



OVERALL WIDTH

18 ¾ in. min. to 31 ½ in. max.

47.6 cm min. to 80 cm max.

FRONT SEAT-TO-FLOOR HEIGHT 14 in. to 21 in. | 35.6 cm to 53.3 cm



ARMREST HEIGHT 6 in. to 14 in. | 15.2 cm to 35.6 cm



WHEEL CAMBER 0°, 2°, 4°, 6°, 8°



WEIGHT LIMIT 265 lb. | 120 kg

COLOR CHART COLORS ARE ON PROTECTORS. FRAME IS BLACK.



Looking for a carbon option? Check out our APEX Carbon > motioncomposites.com

70°, 75°, 80°, 85°, 90° 80° to 101°

REAR SEAT-TO-FLOOR HEIGHT BACK HEIGHT 14 in. to 20 in. | 35.6 cm to 50.8 cm 9 in. to 21 in. | 22.9 cm to 53.3 cm



TECHNICAL SPECIFICATIONS

Riaid

Triple-butted aluminum

265 lb. | 120 ka

10.6 lb. | 4.8 kg (16 x 16 without rear wheels, wheel locks, armrests, cushion and anti-tippers)

12 in. to 20 in. | 30.5 cm to 50.8 cm

12 in. to 20 in. | 30.5 cm to 50.8 cm

14 in. to 21 in. | 35.6 cm to 53.3 cm (with 3 in. to 6 in. caster)

14 in. to 20 in. | 35.6 cm to 50.8 cm (with 20 in. to 26 in. wheel)

18.0 lb. 8.2 kg, with wheel locks and wheels

FRAME

WIDTH

DEPTH

MATERIAL

WEIGHT LIMIT

TRANSPORT WEIGHT

LIGHTEST CONFIGURATION

FRONT SEAT-TO-FLOOR HEIGHT

REAR SEAT-TO-FLOOR HEIGHT

DIMENSIONS

SEAT WIDTH

FRONT FRAME ANGLE

BACK ANGLE

GO BEYOND!

Who we are: passionate innovators pushing the boundaries of what is possible. What we do: combine leading-edge engineering with innovative materials to create movement. What makes us proud: having our ultralight wheelchairs echo the love of life their users demonstrate each day. Through our unique approach and distinct attitude, our aim is to enhance lives, nourish a sense of freedom and revolutionize the way people see wheelchairs.

LEARN MORE To learn more about the APEX A and find more information regarding our products, visit our website at **motioncomposites.com**



T 1 (866) 650-6555 F 1 (888) 966-6555

motioncomposites.com

370689B-00-v4-APEX ALUMINUM BROCHURE - ENGLISH