



**1- Purchasing information**

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_

Company Name:\* \_\_\_\_\_

Purchasing contact:\* \_\_\_\_\_

Phone:\* \_\_\_\_\_

Email:\* \_\_\_\_\_

Contact preference:\* Phone ☐ Email ☐

Marked for: \_\_\_\_\_

**3- Shipping information**

Ship to:\* \_\_\_\_\_

Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ Province: \_\_\_\_\_

Postal Code: \_\_\_\_\_

Phone: \_\_\_\_\_

**2- Dealer/Representative information**

Dealer/Representative:\* \_\_\_\_\_

Phone:\* \_\_\_\_\_

Email:\* \_\_\_\_\_

Contact preference:\* Phone ☐ Email ☐

\* Required

MA040 ☐

English Manual (STD)

MA041 ☐

Manuel Français

Taxes not included

MODEL									
WC11	<input type="checkbox"/>	Ultralight <b>Aluminum APEX Rigid</b> Wheelchair (265 lb weight limit) Not available with painted frame - Frame will be black	Triple-buttet aluminum frame / Aluminum back canes and back rigidizer / Aluminum seat rigidizer / Ultralight carbon fiber camber tube and rear frame rigidizer						\$3,845
EVO-R	<input type="checkbox"/>	Evo Rigid - Chair resizing service*	To be used within 24 months of chair purchase						\$700

\*By selecting this option, I agree to the terms and conditions of the EVO-RIGID program.

[https://www.motioncomposites.com/en\\_us/terms-and-conditions](https://www.motioncomposites.com/en_us/terms-and-conditions)

Frame											
SEAT WIDTH											
SW	<input type="checkbox"/>	12" <input type="checkbox"/>	13" <input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <input type="checkbox"/>			NC
	<input type="checkbox"/>		19" <input type="checkbox"/>	20" <input type="checkbox"/>							\$360
SEAT DEPTH											
SD	<input type="checkbox"/>	12" <input type="checkbox"/>	13" <input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <input type="checkbox"/>			NC
	<input type="checkbox"/>		19" <input type="checkbox"/>	20" <input type="checkbox"/>							\$275
FRAME LENGTH											
SP01	<input type="checkbox"/>	Standard frame length configuration <sup>1</sup>									STD
SP02	<input type="checkbox"/>	Frame length 1" longer than standard configuration									NC
SP03	<input type="checkbox"/>	Frame length 2" longer than standard configuration									NC
SP04	<input type="checkbox"/>	Frame length 1" shorter than standard configuration									NC
SP05	<input type="checkbox"/>	Frame length 2" shorter than standard configuration									NC
TIE-DOWN											
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>2,3</sup>				Short belt 20"-41"			<i>Client cost</i>		\$365
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>2,3</sup>				Long belt 20"-46"			<i>Client cost</i>		\$365
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt <sup>3</sup>							<i>Client cost</i>		\$345

1- Distance of 4.37" from end of seat upholstery to bend

2- Not available with depth adjustable back; 3- May interfere with the swing away armrest and the clamp on height adjustable push handle; Add a 3" of frame length behind backrest


[4. Width calculator](#)

FRONT Seat-to-floor height		
Caster sizes		Height available +- 1/4"
<input type="checkbox"/>	3" Caster	14" to 21"
<input type="checkbox"/>	4" Caster	14.5" to 21"
<input type="checkbox"/>	5" Caster	15.5" to 21"
<input type="checkbox"/>	6" Caster	16.5" to 21" (N/A with 90° front frame)
		Specify desired height:  _____ (inches)



CASTERS				Available		
CA01-03 to CA01-06	<input type="checkbox"/>	Composite wheel w/ PU tire (1")		3", 4", 5", 6"		NC
CA10	<input type="checkbox"/>	Pneumatic casters (1 ¼")		6"		\$170
CA14 to CA16	<input type="checkbox"/>	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		\$105
CA17 to CA19	<input type="checkbox"/>	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		\$198
CASTER OPTIONS						
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks				\$525

REAR Seat-to-floor height			
Wheel sizes		Height available +/- 1/4"	Specify desired height:  <div>(inches)</div>
<input type="checkbox"/>	20" (451) wheel	14" to 17.5"	
<input type="checkbox"/>	22" (501) wheel	14" to 18.5"	
<input type="checkbox"/>	24" (540) wheel	15" to 19.5"	
<input type="checkbox"/>	25" (559) wheel	15.5" to 20"	
<input type="checkbox"/>	26" (590) wheel	16" to 20"	




FRONT FRAME ANGLE						
FA70	<input type="checkbox"/>	70° Front frame angle				NC
FA75	<input type="checkbox"/>	75° Front frame angle				NC
FA80	<input type="checkbox"/>	80° Front frame angle				NC
FA85	<input type="checkbox"/>	85° Front frame angle				NC
FA90	<input type="checkbox"/>	90° Front frame angle				NC

Seat upholstery to footplate length	
12" to 19" without high mount 7" to 13" with high mount <sup>1</sup>	Specify desired length:  (inches) <hr/> Required



1- Choose high mount in the footplates section

Footrest width / Taper selection				
Seat Width	Straight	1" Taper	2" Taper	
Outside Measurement	2 ½" Narrower than Seat Width	4 ½" Narrower than Seat Width	6 ½" Narrower than Seat Width	
12"	9.5"	7.5"		
13"	10.5"	8.5"		
14"	11.5"	9.5"		
15"	12.5"	10.5"		
16"	13.5"	11.5"	9.5"	
17"	14.5"	12.5"	10.5"	
18"	15.5"	13.5"	11.5"	
19"	16.5"	14.5"	12.5"	
20"	17.5"	15.5"	13.5"	

Note: higher than 20" Front STFH will reduce the footrest width up to 1" with a 2" taper

FOOTPLATES						
FP12	<input type="checkbox"/>	Tubular open loop <sup>1</sup>				STD
FP13	<input type="checkbox"/>	Tubular with ABS cover				NC
FP31	<input type="checkbox"/>	Newton adjustable angle - Aluminum ( <i>one piece</i> )				\$105
FP32	<input type="checkbox"/>	Newton adjustable angle - Carbon fiber ( <i>one piece with rubber inserts</i> )				\$210
FP14	<input type="checkbox"/>	Flip back, adjustable angle aluminum				\$240
FP15	<input type="checkbox"/>	Flip back, adjustable angle carbon fiber				\$355
FR25	<input type="checkbox"/>	High mount with adjustable angle footplate ( <i>Aluminum - one piece</i> )				\$200
FR26	<input type="checkbox"/>	High mount with adjustable angle footplate ( <i>Carbon fiber with rubber inserts - one piece</i> )				\$280
FOOTPLATES ACCESSORIES						
FP16	<input type="checkbox"/>	2" riser ( <i>only available with aluminum angle adjustable footplate</i> ) <sup>2</sup>				\$105
FP09	<input type="checkbox"/>	Calf strap				\$57
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap				\$140

1 - Available on: all straight frames; 1" taper with seat width of 14" and greater; 2" taper with seat width of 16" and greater

2- Not available with high mount

Backrest

FOLD DOWN BACKREST

BC20	<input type="checkbox"/>	Locking backrest			STD
BC21	<input type="checkbox"/>	Ultralight non locking backrest (not locking in upright and down position) <sup>2</sup>			\$105

FIXED HEIGHT, ADJUSTABLE ANGLE FOLD DOWN BACKREST (Not available with push handles)


BC	<input type="checkbox"/>	Fixed height without push handle	(30) 9" <input type="checkbox"/>	(31) 9.5" <input type="checkbox"/>	(32) 10" <input type="checkbox"/>	(33) 10.5" <input type="checkbox"/>	(34) 11" <input type="checkbox"/>		NC	
		(35) 11.5" <input type="checkbox"/>	(36) 12" <input type="checkbox"/>	(37) 12.5" <input type="checkbox"/>	(38) 13" <input type="checkbox"/>	(39) 13.5" <input type="checkbox"/>	(40) 14" <input type="checkbox"/>			(41) 14.5" <input type="checkbox"/>
		(42) 15" <input type="checkbox"/>	(43) 15.5" <input type="checkbox"/>	(44) 16" <input type="checkbox"/>	(45) 16.5" <input type="checkbox"/>	(46) 17" <input type="checkbox"/>	(47) 17.5" <input type="checkbox"/>			(48) 18" <input type="checkbox"/>
		(49) 18.5" <input type="checkbox"/>	(50) 19" <input type="checkbox"/>	(51) 19.5" <input type="checkbox"/>	(52) 20" <input type="checkbox"/>	(53) 20.5" <input type="checkbox"/>	(54) 21" <input type="checkbox"/>			

HEIGHT AND ANGLE ADJUSTABLE FOLD DOWN BACKREST WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)

BC	<input type="checkbox"/>	Straight back canes	(22) 9 - 12" <input type="checkbox"/>	(23) 12 - 15" <input type="checkbox"/>	(24) 15 - 18" <input type="checkbox"/>	(25) 18 - 21" <input type="checkbox"/>			NC
Set height to: *		( in increments of 1/2", measured from seat sling to top of back upholstery ) * If not indicated, default height setting is at the lowest position							

Select back cane angle: <sup>1</sup>

101° ☐ 99° ☐ 96° ☐ 93° ☐ 90° ☐ 87° ☐ 83° ☐ 80° ☐



STD

1- If not indicated default angle setting is at 93°; 2- Maximum back cane height of 19.5" with aluminum chair

PUSH HANDLES OPTIONS (ONE option at a time)


PH01	<input type="checkbox"/>	Push handles			STD
PH05	<input type="checkbox"/>	Fold down push handles			\$165
PH06-920	<input type="checkbox"/>	Clamp on height adjustable push handles 9 - 20" <sup>1</sup>			\$345
PH06-923	<input type="checkbox"/>	Clamp on height adjustable push handles 9 - 23" <sup>1</sup>			\$345
PH02	<input type="checkbox"/>	Omit push handles			NC

STROLLER HANDLE

PH04	<input type="checkbox"/>	Stroller handle (Push handles included)		WAMW	\$220
------	--------------------------	---	--	------	-------

1- Not compatible with swing away armrests;

BACKRESTS OPTIONS

BC61	<input type="checkbox"/>	Depth adjustable back - 1" depth adjustment <sup>1,4</sup>			\$260
BC62	<input type="checkbox"/>	Depth adjustable back - 2" depth adjustment <sup>2,4</sup>			\$260
BC56	<input type="checkbox"/>	Depth adjustable back - 3" depth adjustment <sup>3,4</sup>			\$260
BC57	<input type="checkbox"/>	Standard 3 inch deep rigidizer bar			STD
BC55	<input type="checkbox"/>	Four inch deep Aluminum rigidizer bar			NC
BC88	<input type="checkbox"/>	Custom height rigidizer bar set at <div>_____ (inches)</div> <div>Measured from the top of the seat sling to the center of the rigidizer - standard height is 5,7" Lowest position is 4,25"</div>	<div></div>		\$135

If used in combinaison with WC-19: 1- Maximum depth of 19"; 2- Maximum depth of 18"; 3- Maximum depth of 17"; 4- Only compatible with WC-19 transit tie down without Q'Straint belt

BACK UPHOLSTERY

BH03	<input type="checkbox"/>	Slip-on Syntech-Suede back upholstery			STD
BH04	<input type="checkbox"/>	Adjustable tension Syntech-Suede back upholstery	Seating benefit		\$260
BH99	<input type="checkbox"/>	Omit Back upholstery			NC

SEAT CUSHION

SC01	<input type="checkbox"/>	2" Cushion	Seating benefit		\$205
SC02	<input type="checkbox"/>	3" Cushion	Seating benefit		\$275
SC99	<input type="checkbox"/>	Omit seat cushion			STD

APEX Orderform - AADL 6.3

October 1st, 2024

3 of 6

Rear wheels														
REAR WHEELS							Available							
RW01	<input type="checkbox"/>	Mag <sup>1</sup>					20", 22" 24"			\$140				
RW02	<input type="checkbox"/>	Newton One - spoke wheel <sup>2</sup>					20", 22" 24"			STD				
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel					20", 22", 24", 25", 26"			\$390				
RW04	<input type="checkbox"/>	Spinergy Spox	Blk STD <input type="checkbox"/>	Wht <sup>3</sup> <input type="checkbox"/>	Blue <sup>3</sup> <input type="checkbox"/>	Red <sup>3</sup> <input type="checkbox"/>	Yel <sup>3</sup> <input type="checkbox"/>	22", 24", 25", 26"		\$875				
RW05	<input type="checkbox"/>	Spinergy LX	Blk STD <input type="checkbox"/>	Pur <sup>3</sup> <input type="checkbox"/>	Wht <sup>3</sup> <input type="checkbox"/>	Blue <sup>3</sup> <input type="checkbox"/>	Red <sup>3</sup> <input type="checkbox"/>	Yel <sup>3</sup> <input type="checkbox"/>	Pink <sup>3</sup> <input type="checkbox"/>	Grn <sup>3</sup> <input type="checkbox"/>	Org <sup>3</sup> <input type="checkbox"/>	22", 24", 25", 26"		\$1,095
RW07	<input type="checkbox"/>	Spinergy CLX <sup>4</sup>	Blk STD <input type="checkbox"/>	Pur <sup>3</sup> <input type="checkbox"/>	Wht <sup>3</sup> <input type="checkbox"/>	Blue <sup>3</sup> <input type="checkbox"/>	Red <sup>3</sup> <input type="checkbox"/>	Yel <sup>3</sup> <input type="checkbox"/>	Pink <sup>3</sup> <input type="checkbox"/>	Grn <sup>3</sup> <input type="checkbox"/>	Org <sup>3</sup> <input type="checkbox"/>	24", 25"		\$1,650
TIRES							Available							
AIRLESS														
TI01	<input type="checkbox"/>	Soft Urethane 1 3/8" - Low Tread, black						20", 22" 24"			NC			
TI02	<input type="checkbox"/>	Soft Urethane 1 3/8 " - Medium Tread, black <sup>5</sup>						22", 24", 25", 26"			\$120			
TI04	<input type="checkbox"/>	Soft Urethane 1" - No tread, black - <i>Shox</i>						22", 24", 25", 26"			\$140			
TI08	<input type="checkbox"/>	Pneumatic Airless insert 1 3/8" - Medium Tread <sup>6</sup>						22", 24", 26"			\$225			
PNEUMATIC														
TI05	<input type="checkbox"/>	Pneumatic 1 3/8" - Medium Tread						20", 22", 24", 25", 26"			STD			
TI06	<input type="checkbox"/>	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>						20", 22", 24", 25", 26"			\$120			
TI07	<input type="checkbox"/>	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>						24", 25", 26"			\$200			

1- Not compatible with HP tires; Only available with aluminum or plastic coated handrim in standard position; 2- Not compatible with Shox tires; 3- Colour spokes require an additional 10 days lead-time; 4- Only available with 1" tires; 5- The 22" is only available with Mags and Newton One rear wheels; 6- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels

Handrims					
HANDRIM				Available	
HR01	<input type="checkbox"/>	Aluminum anodized		20", 22", 24", 25", 26"	STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized <sup>2</sup>		20", 22", 24", 25"	\$85
HR02	<input type="checkbox"/>	Plastic coated <sup>1</sup>		20", 22", 24", 25", 26"	\$170
HR03	<input type="checkbox"/>	High friction coated <sup>1,2</sup>		22", 24", 25", 26"	\$460
HR07-HR18 and HR65-HR71	<input type="checkbox"/>	Newton Air Grip <sup>2</sup> <i>See page 6 for the color selection and price</i>		20", 22", 24", 25", 26"	
HR28	<input type="checkbox"/>	Surge <sup>1,2,3</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"	\$485
HR05	<input type="checkbox"/>	Surge LT <sup>1,2,3</sup> (1 1/8" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"	\$515
HR04	<input type="checkbox"/>	Natural Fit <sup>1,2,3</sup>		22", 24", 25", 26"	\$320
HR62	<input type="checkbox"/>	Simi H <sup>1,2</sup> (1 3/8" Oval handrim with high friction strip)		22", 24", 25", 26"	\$485
HR63	<input type="checkbox"/>	Nova H <sup>1,2</sup> (1 1/4" Oval handrim with high friction strip)		22", 24", 25", 26"	\$515
HR64	<input type="checkbox"/>	Optimum H <sup>1,2</sup> (1 1/8" heart-shaped handrim)		22", 24", 25", 26"	\$320
HR99	<input type="checkbox"/>	Omit handrim			NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels 3- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

Handrims option					
ASSEMBLY POSITION					
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab		STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	Full tab		NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut <sup>3</sup>	Short tab 2 positions		\$29
HR31	<input type="checkbox"/>	Super Narrow Handrim <sup>3</sup>	Short tab 2 positions STD <input type="checkbox"/> Short tab 1 position <input type="checkbox"/>		\$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	<input type="checkbox"/>	Standard grip			\$200
HR27	<input type="checkbox"/>	Super grip			\$230
HR099	<input type="checkbox"/>	Omit thumb piece			NC

3- Only available with Aluminum and Air Grip handrims;

AXLES					
AX01	<input type="checkbox"/>	Quick release axle			STD
AX03	<input type="checkbox"/>	Quad release axle			\$180
REAR WHEEL CAMBER					
CM00	<input type="checkbox"/>	0° camber			STD
CM02	<input type="checkbox"/>	2° camber			NC
CM04	<input type="checkbox"/>	4° camber			NC
CM06	<input type="checkbox"/>	6° camber			NC
CM08	<input type="checkbox"/>	8° camber			NC


CENTER OF GRAVITY

Specify the position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

0 - 6" range	Set At: _____ (inches)
--------------	------------------------

If not indicated, default setting is at 2.5" Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG



1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

HEIGHT ADJUSTABLE "T" ARMBREST <sup>1</sup>	
---	--

AR13-610	<input type="checkbox"/>	10" Desk length armpad <sup>3</sup>	6" - 10"		\$380
----------	--------------------------	-------------------------------------	----------	--	-------

HEIGHT ADJUSTABLE SWING AWAY ARMREST <sup>2</sup>	
---	--

AR05	<input type="checkbox"/>	Tubular armrest - short pad	8" - 11"			\$315
------	--------------------------	-----------------------------	----------	--	--	-------

[illegible]

AP12	<input type="checkbox"/>	14" Full Length armed Gel Quatern (pair) 4			\$160
------	--------------------------	--	--	--	-------

1- May have 1/2" adjustment restriction on both limits of the range; 2- not available with the Ultralight non-locking backrest; 3- May have some restriction with seat-to-floor height and rear wheel

4-Only available with "T" full length armrest; Not available with swing away and desk length armrests

**SIDE GUARDS (not available with "T" Armrest)**

SG01	<input type="checkbox"/>	Plastic side guards (removable)			\$235
------	--------------------------	---------------------------------	--	--	-------




























---

SP01-6		M.L. "L.A." #27"1		\$50
--------	--	-------------------	--	------

1- Not recommended for seat width of 18" and more

[illegible]

1. Not available with Mag wheels and Energizer.

Colors																						
Options <sup>1</sup>	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)		
																						
AirGrip	<input type="checkbox"/> \$160		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$350	
Camber Tube	<input type="checkbox"/> STD/NC		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$75	
																						
Highlight <sup>3</sup>	<input type="checkbox"/>		<input type="checkbox"/>					<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>			<input type="checkbox"/>					NC

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.;    Finish between Airgrip and frame might be different;

3- To see different highlight, visit our website: [https://www.motioncomposites.com/en\\_co/products/rigid-wheelchairs/apex-c](https://www.motioncomposites.com/en_co/products/rigid-wheelchairs/apex-c)

Notes :