



P.O. number:*	Date:	Ship to:*
Company Name:*	Phone:*	Name:
Ordered by:*		Address:
Sales rep/Clinician:		City: Prov:
Email:		Postal Code:
Marked for:		Phone:

* Required MA040 English Manual (STD) MA041 Manuel Français Taxes not included


MODEL					
1	WC11	Ultralight Aluminum APEX Rigid Wheelchair (265 lb weight limit)	Tripple-butted aluminum frame / Aluminum back canes and back rigidizer / Aluminum seat rigidizer / Ultralight carbon fiber camber tube and rear frame rigidizer	WM0001305	\$3 595
2	WC21	Apex Carbon fiber frame upgrade (High Gloss Carbon Fiber)	Ultralight carbon fiber frame, back canes, back rigidizer, seat rigidizer, camber tube and rear frame rigidizer.	WAMZ	\$900

SEAT WIDTH					
3	SW	Seat width:	12" ☐	13" ☐	14" ☐
4	SW	Seat width:	19" ☐	20" ☐	
SEAT DEPTH					
5	SD	Seat depth:	12" ☐	13" ☐	14" ☐
6	SD	Seat depth:	19" ☐	20" ☐	

FRAME LENGTH					
7	SP01	Standard frame length configuration ¹		STD	☐
8	SP02	Frame length 1" longer than standard configuration		NC	☐
9	SP03	Frame length 2" longer than standard configuration		NC	☐
10	SP04	Frame length 1" shorter than standard configuration		NC	☐
11	SP05	Frame length 2" shorter than standard configuration		NC	☐

1- Distance of 4.37" from seat end to bend

FRONT Seat-to-floor height		
Caster sizes	Height available +- 1/4"	Specify desired height: (inches)
3" Caster	14" to 21"	
4" Caster	14.5" to 21"	
5" Caster	15.5" to 21"	
6" Caster	16.5" to 21" (N/A with 90° front frame)	



CASTERS									
12	CA01	Composite wheel w/ PU tire (1" width)	3" ☐	4" ☐	5" ☐	6" ☐		NC	☐
13	CA10	Pneumatic casters 6" x 1 ¼"						\$160	☐
14	CA14	Newton UltraCaster Soft-Roll Composite	4" x 1 ½ "					\$98	☐
15	CA15		5" x 1 ½ "					\$98	☐
16	CA16		6" x 1 ½ "					\$98	☐
17	CA17	Newton UltraCaster Soft-Roll Aluminum	4" x 1 ½ "					\$175	☐
18	CA18		5" x 1 ½ "					\$185	☐
19	CA19		6" x 1 ½ "					\$195	☐
CASTER OPTIONS									
20	CA13	Carbon fiber Frog Legs suspension forks						\$489	☐

A line drawing of a wheelchair from a side profile. A red vertical double-headed arrow indicates the measurement from the floor to the center of the rear wheel. The arrow starts at a horizontal line representing the floor and extends upwards to a point on the wheel's axle.

A line drawing of a wheelchair from a side profile. A red double-headed vertical arrow is positioned between the floor and the top of the seat, indicating the measurement for the 'Seat to floor' dimension.

Required

Note: higher than 20" Front STFH will reduce the footrest width up to 1" with a 2" taper

FOOTPLATES					
26	FP12	Tubular open loop ¹		STD	□
27	FP13	Tubular with ABS cover		NC	□
28	FP31	Adjustable angle <i>(Aluminum - one piece)</i>	WAMI / SEFND0015	\$98	□
29	FP32	Newton adjustable angle <i>(Carbon fiber with rubber inserts - one piece)</i>	WAMI / SEFND0015	\$198	□
30	FP14	Flip back, adjustable angle <i>(Aluminum)</i>	WAMI / SEFND0015	\$230	□
31	FP15	Flip back, adjustable angle <i>(Carbon fiber)</i>	WAMI / SEFND0015	\$335	□
32	FR25	High mount with adjustable angle footplate <i>(Aluminum - one piece)</i>	New	190 \$	□
33	FR26	High mount with adjustable angle footplate <i>(Carbon fiber with rubber inserts - one piece)</i>	New	260 \$	□
FOOTPLATES ACCESSORIES					
34	FP16	2" riser <i>(only available with aluminum angle adjustable footplate)</i> ²	SEFND0030	\$95	□
35	FP09	Calf strap	SEFND0010	\$52	□
36	FP30	Bodypoint Aeromesh Padded calf strap	SEFND0010	\$129	□

Motion Composites

APEX Orderform CAN 5.1.xlsx

March 23rd, 2020 2 of 6

FOLD DOWN BACKREST					
37	BC20	Locking backrest		STD	☐
38	BC21	Ultralight non locking backrest (not locking in upright and down position)		\$98	☐

Back canes fixed height, angle adjustable, Without push handle

ADJUSTABLE ANGLE FOLD DOWN BACKREST										
39	BC	Fixed height without push handle		(30) 9" ☐	(31) 9.5" ☐	(32) 10" ☐	(33) 10.5" ☐	(34) 11" ☐	NC	☐
		(35) 11.5" ☐	(36) 12" ☐	(37) 12.5" ☐	(38) 13" ☐	(39) 13.5" ☐	(40) 14" ☐	(41) 14.5" ☐		
		(42) 15" ☐	(43) 15.5" ☐	(44) 16" ☐	(45) 16.5" ☐	(46) 17" ☐	(47) 17.5" ☐	(48) 18" ☐		
		(49) 18.5" ☐	(50) 19" ☐	(51) 19.5" ☐	(52) 20" ☐	(53) 20.5" ☐	(54) 21" ☐			

Back canes, height and angle adjustable, with option

ADJUSTABLE ANGLE FOLD DOWN BACKREST									
40	BC	Straight back canes						NC	☐
		(22) 9 - 12" ☐	(23) 12 - 15" ☐	(24) 15 - 18" ☐	(25) 18 - 21" ☐				


Set height to: *		(in increments of 1/2", measured from seat sling to top of back upholstery) * If not indicated, default height setting is at the lowest position											
41		Select back cane angle:	101° ☐	99° ☐	96° ☐	93° ¹ ☐	90° ☐	87° ☐	83° ☐	80° ☐		STD	☐

1 - If not indicated default angle setting is at 93°

PUSH HANDLE OPTION					
42	PH01	Push handles		STD	☐
43	PH05	Fold down push handles ¹		\$150	☐
44	PH06	Clamp on height adjustable push handle ^{1,2}	9 - 20" STD ☐	9-23" ☐	\$325 ☐
45	PH02	Without push handles		NC	☐

1- 265 pounds weight limit

2 - Not compatible with swing away armrests; Only available with 4" rigidizer bar

BACKRESTS OPTIONS					
46	BC57	Standard 3 inch deep rigidizer bar		STD	☐
47	BC55	Four inch deep Aluminum rigidizer bar ¹		NC	☐
48	BC56	Depth adjustable back (provides 3" of depth adjustment) ²		\$245	☐
49	BC88	Custom height rigidizer bar set at _____ (inches) Measured from the top of the seat sling to the center of the rigidizer - standard height is 5,7 " Lowest position is 4,25"			\$125 ☐

1- Available with seat depths of 14" and more

2- Not available with seat depths of 18" and more

BACK UPHOLSTERY					
50	BH03	Slip-on Syntech-Suede back upholstery		NC	☐
51	BH04	Adjustable tension Syntech-Suede back upholstery	WAMD/WAME	\$245	☐
52	BH99	Omit Back upholstery		STD	☐
RIGID BACK					
53		nXt back series (See nXt back pricelist)		See price list p.7	

SEAT CUSHION					
54	SC01	2" Cushion		\$175	☐
55	SC02	3" Cushion		\$195	☐
56	SC99	Omit seat cushion		STD	☐
57		nXt FIT cushion series (See nXt FIT pricelist)		See price list p.7	☐

Rear wheels									
REAR WHEELS									
58	RW01	Mag ¹	20" ☐	22" ☐	24" ☐			\$129	☐
59	RW02	Newton One - spoke wheel ²	20" ☐	22" ☐	24" ☐			STD	☐
60	RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐	\$360	☐
61	RW04	Spinergy Spox		22" ☐	24" ☐	25" ☐	26" ☐	\$800	☐
62		Spinergy Spox spokes colors ³	Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC	☐
63	RW05	Spinergy LX		22" ☐	24" ☐	25" ☐	26" ☐	\$995	☐
64		Spinergy LX spokes colors ³	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC
65	RW07	Spinergy CLX Carbon Blade			24" ☐	25" ☐		\$1 540	☐
66		Spinergy CLX spokes colors ³	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC
						Pink ☐	Grn ☐	Org ☐	☐
TIRES							Available		
67	TI01	Soft Urethane 1 1/4 "					20", 22" 24"	NC	☐
68	TI02	Soft Urethane knobby 1 3/8 "					22", 24", 25", 26"	\$110	☐
69	TI04	Soft Urethane - Shox 1"					22", 24", 25", 26"	\$130	☐
70	TI05	Pneumatic 1 3/8 "					20", 22", 24", 25", 26"	STD	☐
71	TI08	Pneumatic with Airless insert 1 3/8"					22", 24"	\$215	☐
72	TI06	Pneumatic HP - 1"					20", 22", 24", 25", 26"	\$110	☐
73	TI07	Pneumatic HP - Schwalbe Marathon Plus 1"					24", 25", 26"	\$190	☐

1 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position
2- Not compatible with Shox tires

3 - Colour spokes require an additionaln 10 days lead-time

Handrims									
HANDRIMS									
74	HR01	Aluminum anodized						STD	☐
75	HR02	Plastic coated ³					WAMO	\$160	☐
76	HR03	Q-Grip High friction Neoprene coated ^{1,2,3}					WAMO	\$430	☐
77	HR07	Newton Air Grip ¹	Black ☐				WAMO	\$150	☐
78	HR	Newton Air Grip ¹	(08) Sunkissed Orange ☐	(09) Burgundy ☐	(10) Sapphire Blue ☐		WAMO	\$330	☐
		(11) Monster Green ☐	(12) Charcoal ☐	(13) Fuchsia ☐	(14) Steel Blue ☐				
		(15) Ferrari Red ☐	(16) Acid Green ☐	(17) White ☐	(18) Sky Blue ☐				
79	HR28	Surge ^{1,2,3}	(1-1/2" Oval handrim w/ Gripton high friction strip)					\$470	☐
80	HR05	Surge LT ^{1,2,3}	(1-1/8" Oval handrim w/ Gripton high friction strip)					\$500	☐
81	HR04	Natural Fit ^{1,2,3}						\$310	☐
82	HR99	Omit handrim						NC	☐


1 - Not available with Mag wheels
2 - Not available in 20"
3- Only standard handrim position available

Handrims option							
Assembly position for aluminum and Newton AirGrip handrims							
83	HR29	Standard Handrim position	Full tab			STD	☐
84	HR30	Narrow Handrim position	Short tab 2 positions			\$25	☐
85	HR31	Super Narrow Handrim	Short tab 2 positions STD ☐		Short tab 1 position ☐		☐
Thumb piece ⁴ selection for Natural Fit handrims							
86	HR23	Standard grip				\$190	☐
87	HR27	Super grip				\$220	☐
88	HR099	Omit thumb piece				NC	☐

4- Compatible with Surge, Surge LT handrims

AXLES									
89	AX01	Quick release axle					WAMU	STD	☐
90	AX03	Quad release axle					WAMU	\$170	☐

REAR WHEEL CAMBER									
91	CM00	0° camber						STD	☐
92	CM02	2° camber						NC	☐
93	CM04	4° camber						NC	☐
94	CM06	6° camber						NC	☐
95	CM08	8° camber						NC	☐

CENTER OF GRAVITY									
Specify the position of rear wheel axle plate: Measure distance from front of back post to center of axle receiver. 0 - 6" range Set At: (inches) If not indicated, default setting is at 2,5" If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position. 									
WHEEL LOCK									
96	WL08	Newton wheel lock - push to lock						STD	☐
97	WL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)					WAMR	\$175	☐
98	WL13 WL14	6" removable extension handle						\$67	☐
99	WL05	Scissor Lock						\$175	☐
100	WL99	Omit wheel locks						NC	☐
Armrests									
"T" ARMREST ¹			Height Adjustment						
101	AR03	10" Desk Length	8 - 12" STD ☐	10 - 14" ☐			\$360	☐	
102	AR04	14" Full Length	8 - 12" STD ☐	10 - 14" ☐			\$360	☐	
SWING AWAY ARMREST ²			Height Adjustment						
103	AR05	Tubular arm pad	8 - 11"				\$295	☐	
104	AR06	10" Desk length arm pad	8 - 11"				\$295	☐	
105	AR99	Omit armrest					NC	☐	
1 - May have a 1/2 " adjustment restriction on both limits of the range 2 - not available with the Ultralight non-locking backrest									
SIDE GUARDS <i>(not available with "T" Armrest)</i>									
106	SG01	Plastic side guards <i>(removable)</i>					WAMV	\$225	☐
107	SG03	Carbon fiber side guards <i>(removable)</i>					WAMV	\$295	☐
108	SG06	Plastic fender side guards <i>(removable)</i> ¹					WAMV	\$320	☐
109	SG04	Carbon fiber fender side guards <i>(removable)</i>					WAMV	\$495	☐
110	SG99	Omit side guards						NC	☐
SIDE GUARDS OPTION									
111	SG07	Side guard clip <i>(pair)</i>					New	NC	☐
1 - Not compatible with 20" and 22" rear wheels									
POSITIONING BELT									
112	SB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐			\$45	☐	
113	SB02	Auto buckle	Standard 35" ☐	Long 60" ☐			\$55	☐	
114	SB07	Hardware attachment for 2 point seatbelt (1 pair)						NC	☐
115	SB03	Additional hardware for 2 extra attachment points <i>(needed for the 4 point seat belt)</i> (1 pair)						\$13	☐
116	SB04	BodyPoint Hip Belt - 2 point non padded		Medium 56" ☐	Long 61" ☐		\$82	☐	
117	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	SERND0030	\$170	☐	
118	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	SERND0020	\$195	☐	
119	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps		Medium 18" ☐	Long 24" ☐		\$324	☐	
120	SB99	Omit positioning belt						STD	☐
ACCESSORIES									
121	AC03	Anti-tippers						\$115	☐
122	AC05	WC-19 Transit Tie Down w/Q'straint belt ¹		Short belt 20-41" ☐	Long belt 20-46" ☐		\$345	☐	
123	AC08	Spoke guards ²		Available: 20", 22", 24", 26"		WAMP	\$145	☐	
124	AC10	Seat pouch						\$95	☐
125	AC26	BodyPoint Mobility Bag						\$135	☐
126	AC23	Newton Backpack						\$120	☐
127	AC13	Tool kit						\$50	☐
128	AC16	Motion Composites Polo Shirt	W/S ☐ M/S ☐	W/M ☐ M/M ☐	W/L ☐ M/L ☐	M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐	\$32	☐	
1 -May interfere with the swing away armrest and the Clamp on height adjustable push handle 2 - Not available with Mag wheels and Spinergy									

COLOR						
129	FC	Highlight color:	(10) Ferrari Red ☐ (11) White ☐ (12) Black ☐ (13) Sunkissed Orange ☐ (14) Sky Blue ☐ (15) Fuchsia ☐ (16) Acid Green ☐		NC	☐
				Total		

Notes :



BioFit™ CUSHION

	MODEL #	Fits Width	Depth	\$	✓	MODEL #	Fits Width	Depth	\$	✓
	NXCABF1616	15"-16"	16"	595		NXCABF2016	19" - 20"	16"	595	
	NXCABF1618	15"-16"	18"	595		NXCABF2018	19" - 20"	18"	595	
	NXCABF1620	15"-16"	20"	595		NXCABF2020	19" - 20"	20"	625	
	NXCABF1816	17"-18"	16"	595		NXCABF2216	21" - 22"	16"	625	
	NXCABF1818	17"-18"	18"	595		NXCABF2218	21" - 22"	18"	625	
	NXCABF1820	17"-18"	20"	595		NXCABF2220	21" - 22"	20"	625	

THORACIC SHAPE



Optima™
THORACIC BACK SUPPORT

CARBON FIBRE SHELL (WEIGHT SAVINGS – 1lb)

Model # Carbon Fibre	Model # Aluminum	Fits Width	Height	\$ Carbon	✓	\$ Aluminum	✓
	NOACB1307	13"-15"	7"			690	
	NOACB1309	13"-15"	9"			690	
	NOACB1313	13"-15"	13"			690	
	NOACB1316	13"-15"	16"			690	
	NOACB1407	14"-16"	7"			690	
NOCFACB1409	NOACB1409	14"-16"	9"	1095		690	
NOCFACB1413	NOACB1413	14"-16"	13"	1095		690	
NOCFACB1416	NOACB1416	14"-16"	16"	1095		690	
	NOACB1507	15"-17"	7"			690	
NOCFACB1509	NOACB1509	15"-17"	9"	1095		690	
NOCFACB1513	NOACB1513	15"-17"	13"	1095		690	
NOCFACB1516	NOACB1516	15"-17"	16"	1095		690	
	NOACB1607	16"-18"	7"			690	
NOCFACB1609	NOACB1609	16"-18"	9"	1095		690	
NOCFACB1613	NOACB1613	16"-18"	13"	1095		690	
NOCFACB1616	NOACB1616	16"-18"	16"	1095		690	
	NOACB1707	17"-19"	7"			690	
NOCFACB1709	NOACB1709	17"-19"	9"	1095		690	
NOCFACB1713	NOACB1713	17"-19"	13"	1095		690	
NOCFACB1716	NOACB1716	17"-19"	16"	1095		690	
	NOACB1807	18"-20"	7"			690	
NOCFACB1809	NOACB1809	18"-20"	9"	1095		690	
NOCFACB1813	NOACB1813	18"-20"	13"	1095		690	
NOCFACB1816	NOACB1816	18"-20"	16"	1095		690	
	NOACB2007	20"-22"	7"			690	
NOCFACB2009	NOACB2009	20"-22"	9"	1095		690	
NOCFACB2013	NOACB2013	20"-22"	13"	1095		690	
NOCFACB2016	NOACB2016	20"-22"	16"	1095		690	



Optima™ DEEP
THORACIC BACK SUPPORT

CARBON FIBRE SHELL (WEIGHT SAVINGS – 1lb)

Model # Carbon Fibre	Model # Aluminum	Fits Width	Height	\$ Carbon	✓	\$ Aluminum	✓
	NOPCB1307	13"-15"	7"			690	
	NOPCB1309	13"-15"	9"			690	
	NOPCB1313	13"-15"	13"			690	
	NOPCB1316	13"-15"	16"			690	
	NOPCB1407	14"-16"	7"			690	
NOCFPCB1409	NOPCB1409	14"-16"	9"	1195		690	
NOCFPCB1413	NOPCB1413	14"-16"	13"	1195		690	
NOCFPCB1416	NOPCB1416	14"-16"	16"	1195		690	
	NOPCB1507	15"-17"	7"			690	
NOCFPCB1509	NOPCB1509	15"-17"	9"	1195		690	
NOCFPCB1513	NOPCB1513	15"-17"	13"	1195		690	
NOCFPCB1516	NOPCB1516	15"-17"	16"	1195		690	
	NOPCB1607	16"-18"	7"			690	
NOCFPCB1609	NOPCB1609	16"-18"	9"	1195		690	
NOCFPCB1613	NOPCB1613	16"-18"	13"	1195		690	
NOCFPCB1616	NOPCB1616	16"-18"	16"	1195		690	
	NOPCB1707	17"-19"	7"			690	
NOCFPCB1709	NOPCB1709	17"-19"	9"	1195		690	
NOCFPCB1713	NOPCB1713	17"-19"	13"	1195		690	
NOCFPCB1716	NOPCB1716	17"-19"	16"	1195		690	
	NOPCB1807	18"-20"	7"			690	
NOCFPCB1809	NOPCB1809	18"-20"	9"	1195		690	
NOCFPCB1813	NOPCB1813	18"-20"	13"	1195		690	
NOCFPCB1816	NOPCB1816	18"-20"	16"	1195		690	
	NOPCB2007	20"-22"	7"			690	
NOCFPCB2009	NOPCB2009	20"-22"	9"	1195		690	
NOCFPCB2013	NOPCB2013	20"-22"	13"	1195		690	
NOCFPCB2016	NOPCB2016	20"-22"	16"	1195		690	

BACK SUPPORT HARDWARE OPTIONS

All backs come STANDARD with QuickRelease Mounting Hardware. The following are fixed hardware options available at no charge when ordered with a back support.

Model #	DESCRIPTION	WEIGHT OF HARDWARE ONLY	✓
	QuickFit Quick Release Hardware System	STANDARD (2 LBS)	STANDARD
NIB-1642	Micro Mount 2 point fixed hardware. Back width needs to match chair width	Weight savings of 1.4 lb over STANDARD	
NIB-1207	2 point fixed hardware. Width expandable to fit 3 sizes of chair	Weight savings of 1.0 lb over STANDARD	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com