



P.O. number:*	Date:
Company Name:*	Phone:*
Purchasing contact:*	
Sales rep/Clinician:	
Email:	
Marked for:	

Ship to:*	
Name:	
Address:	
City:	Prov:
Postal Code:	
Phone:	

MA040 ☐

English Ma

MA041 □
Manuel Français

Taxes not included

SEAT WIDTH												
3	SW	Seat width:	12" ☐	13" ☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐		NC	☐
4	SW	Seat width:	19" ☐	20" ☐						WAMA	\$350	☐

FRAME LENGTH					
7	SP01	Standard frame length configuration ¹		STD	☐
8	SP02	Frame length 1" longer than standard configuration		NC	☐
9	SP03	Frame length 2" longer than standard configuration		NC	☐
10	SP04	Frame length 1" shorter than standard configuration		NC	☐
11	SP05	Frame length 2" shorter than standard configuration		NC	☐

FRONT Seat-to-floor height		
Caster sizes	Height available +/- 1/4"	
3" Caster	14" to 21"	Specify desired height: (inches)
4" Caster	14.5" to 21"	
5" Caster	15.5" to 21"	
6" Caster	16.5" to 21" (N/A with 90° front frame)	



CASTERS


12	CA01	Composite wheel w/ PU tire (1" width)	3" □	4" □	5" □	6" □		NC	□
13	CA10	Pneumatic casters 6" x 1 ¼"						\$165	□
14	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "					\$98	□
15	CA15		5" x 1 ½ "					\$98	□
16	CA16		6" x 1 ½ "					\$98	□
17	CA17		4" x 1 ½ "					\$180	□
18	CA18	Newton UltraCaster <i>Soft-Roll</i> Aluminum	5" x 1 ½ "					\$190	□
19	CA19		6" x 1 ½ "					\$200	□

CASTER OPTIONS

20	CA13	Carbon fiber Frog Legs suspension forks	495 \$	☐
----	------	---	--------	--------------------------

REAR Seat-to-floor height

Wheel sizes	Height available +/- 1/4"	
20" (451) wheel ☐	14" to 17.5"	Specify desired height: (inches)
22" (501) wheel ☐	14" to 18.5"	
24" (540) wheel ☐	15" to 19.5"	
25" (559) wheel ☐	15.5" to 20"	
26" (590) wheel ☐	16" to 20"	




FRONT FRAME ANGLE

21	FA70	70° Front frame angle			NC	☐
22	FA75	75° Front frame angle			NC	☐
23	FA80	80° Front frame angle			NC	☐
24	FA85	85° Front frame angle			NC	☐
25	FA90	90° Front frame angle			NC	☐

Seat upholstery to footplate length

Range available	Specify desired length:
12" to 19" without High mount	
7" to 13" with High mount ¹	_____ (inches)
	<i>Required</i>



1- Select high mount in the footplates section

Footrest width / Taper selection

Seat Width	Straight	1" Taper	2" Taper
Outside Measurement	2 ½" Narrower than Seat Width	4 ½" Narrower than Seat Width	6 ½" Narrower than Seat Width
12"	9.5"	7.5"	
13"	10.5"	8.5"	
14"	11.5"	9.5"	
15"	12.5"	10.5"	
16"	13.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"

Note: higher than 20" Front STFH will reduce the footrest width up to 1" with a 2" taper

FOOTPLATES				
26	FP12	Tubular open loop ¹		STD
27	FP13	Tubular with ABS cover		NC
28	FP31	Adjustable angle <i>(Aluminum - one piece)</i>	WAMJ / SEFND0015	\$100
29	FP32	Newton adjustable angle <i>(Carbon fiber with rubber inserts - one piece)</i>	WAMJ / SEFND0015	\$205
30	FP14	Flip back, adjustable angle <i>(Aluminum)</i>	WAMJ / SEFND0015	\$235
31	FP15	Flip back, adjustable angle <i>(Carbon fiber)</i>	WAMJ / SEFND0015	\$345
32	FR25	High mount with adjustable angle footplate <i>(Aluminum - one piece)</i>		\$195
33	FR26	High mount with adjustable angle footplate <i>(Carbon fiber with rubber inserts - one piece)</i>		\$270
FOOTPLATES ACCESSORIES				
34	FP16	2" riser <i>(only available with aluminum angle adjustable footplate)</i> ²	SEFND0030	\$100
35	FP09	Calf strap	SEFND0010	\$55
36	FP30	Bodypoint Aeromesh Padded calf strap	SEFND0010	\$135

1 - Available on: all straight frames; 1" taper with seat width of 14" and greater; 2" taper with seat width of 16" and greater

2- Not available with high mount


Backrest												
FOLD DOWN BACKREST												
37	BC20	Locking backrest							STD	☐		
38	BC21	Ultralight non locking backrest (not locking in upright and down position)							\$100	☐		
FIXED HEIGHT, ADJUSTABLE ANGLE FOLD DOWN BACKREST (NO PUSH HANDLE)												
39	BC	Fixed height without push handle		(30) 9" ☐	(31) 9.5" ☐	(32) 10" ☐	(33) 10.5" ☐	(34) 11" ☐	NC	☐		
		(35) 11.5" ☐	(36) 12" ☐	(37) 12.5" ☐	(38) 13" ☐	(39) 13.5" ☐	(40) 14" ☐	(41) 14.5" ☐				
		(42) 15" ☐	(43) 15.5" ☐	(44) 16" ☐	(45) 16.5" ☐	(46) 17" ☐	(47) 17.5" ☐	(48) 18" ☐				
		(49) 18.5" ☐	(50) 19" ☐	(51) 19.5" ☐	(52) 20" ☐	(53) 20.5" ☐	(54) 21" ☐					
HEIGHT AND ANGLE ADJUSTABLE FOLD DOWN BACKREST WITH PUSH HANDLE OPTIONS												
40	BC	Straight back canes							NC	☐		
		(22) 9 - 12" ☐	(23) 12 - 15" ☐	(24) 15 - 18" ☐	(25) 18 - 21" ☐							
Set height to: *			(in increments of 1/2", measured from seat sling to top of back upholstery) * If not indicated, default height setting is at the lowest position									
41		Select back cane angle:						101° ☐ 99° ☐ 96° ☐ 93° ¹ ☐ 90° ☐ 87° ☐ 83° ☐ 80° ☐			STD	☐

1 - If not indicated default angle setting is at 93°

PUSH HANDLE OPTIONS				
42	PH01	Push handles		STD
43	PH05	Fold down push handles ¹		\$155
44	PH06	Clamp on height adjustable push handle ^{1,2}	9 - 20" STD ☐	9-23" ☐
45	PH02	Omit push handles		NC

1- 265 pounds weight limit

2 - Not compatible with swing away armrests;

BACKRESTS OPTIONS				
46	BC56	Depth adjustable back <i>(provides 3" of depth adjustment)</i> ^{2, 3}		\$250
47	BC57	Standard 3 inch deep rigidizer bar		STD
48	BC55	Four inch deep Aluminum rigidizer bar ¹		NC
49	BC88	Custom height rigidizer bar set at _____ (inches) Measured from the top of the seat sling to the center of the rigidizer - standard height is 5,7 " Lowest position is 4,25"		\$130

1- Available with seat width of 14" and more

2- Not available with seat depths of 18" and more

3- Not available with WC-19

BACK UPHOLSTERY				
50	BH03	Slip-on Syntech-Suede back upholstery		NC
51	BH04	Adjustable tension Syntech-Suede back upholstery	WAMD/WAME	\$250
52	BH99	Omit Back upholstery		STD
RIGID BACK				
53		nXt back series (See nXt back pricelist)		See price list

SEAT CUSHION				
54	SC01	2" Cushion		\$180
55	SC02	3" Cushion		\$200
56	SC99	Omit seat cushion		STD
57		nXt FIT cushion series (See nXt FIT pricelist)		See price list


Rear wheels									
REAR WHEELS									
58	RW01	Mag ¹	20" ☐	22" ☐	24" ☐			\$135	☐
59	RW02	Newton One - spoke wheel ²	20" ☐	22" ☐	24" ☐			STD	☐
60	RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐	\$370	☐
61	RW04	Spinergy Spox		22" ☐	24" ☐	25" ☐	26" ☐	\$825	☐
62		Spinergy Spox spokes colors ³	Bik STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC	☐
63	RW05	Spinergy LX		22" ☐	24" ☐	25" ☐	26" ☐	\$1 025	☐
64		Spinergy LX spokes colors ³	Bik STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐
				Grn ☐	Org ☐			NC	☐
65	RW07	Spinergy CLX Carbon Blade ²			24" ☐	25" ☐		1 585 \$	☐
66		Spinergy CLX spokes colors ³	Bik STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐
				Grn ☐	Org ☐			NC	☐
TIRES						Available			
67	TI01	Soft Urethane 1 1/4 "				20", 22" 24"		NC	☐
68	TI02	Soft Urethane knobby 1 3/8 "				22", 24", 25", 26"		\$115	☐
69	TI04	Soft Urethane - Shox 1"				22", 24", 25", 26"		\$135	☐
70	TI05	Pneumatic 1 3/8 "				20", 22", 24", 25", 26"		STD	☐
71	TI08	Pneumatic with Airless insert 1 3/8" ⁴				22", 24", 26"		\$220	☐
72	TI06	Pneumatic HP - 1"				20", 22", 24", 25", 26"		\$115	☐
73	TI07	Pneumatic HP - Schwalbe Marathon Plus 1"				24", 25", 26"		\$195	☐

1 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position
2- Not compatible with Shox tires

3 - Colour spokes require an additionnal 10 days lead-time
4- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims									
HANDRIMS									
74	HR01	Aluminum anodized						STD	☐
75	HR61	Black aluminum hard anodized ^{1,3}						\$77	☐
76	HR02	Plastic coated ²					WAMO	\$165	☐
77	HR03	Q-Grip High friction Neoprene coated ^{2,3,4}					WAMO	\$445	☐
78	HR07	Newton Air Grip ³	Black ☐				WAMO	\$155	☐
79	HR	Newton Air Grip ³	(08) Sunkissed Orange ☐	(09) Burgundy ☐	(10) Sapphire Blue ☐		WAMO	\$340	☐
			(11) Monster Green ☐	(12) Charcoal ☐	(13) Fuchsia ☐	(14) Steel Blue ☐			
			(15) Ferrari Red ☐	(16) Acid Green ☐	(17) White ☐	(18) Sky Blue ☐			
80	HR28	Surge ^{2,3,4}	(1-1/2" Oval handrim w/ Gripton high friction strip)					\$485	☐
81	HR05	Surge LT ^{2,3,4}	(1-1/8" Oval handrim w/ Gripton high friction strip)					\$515	☐
82	HR04	Natural Fit ^{2,3,4}						\$320	☐
83	HR99	Omit handrim						NC	☐

1 - Not available in 26"
2- Only standard handrim position available
3 - Not available with Mag wheels
4 - Not available in 20"

Handrims option					
Assembly position for aluminum and Newton AirGrip handrims					
84	HR29	Standard Handrim position	Full tab		STD ☐
85	HR30	Narrow Handrim position	Short tab 2 positions		\$26 ☐
86	HR31	Super Narrow Handrim	Short tab 2 positions STD ☐ Short tab 1 position ☐		\$26 ☐
Thumb piece ⁴ selection for Natural Fit handrims					
87	HR23	Standard grip			\$195 ☐
88	HR27	Super grip			\$225 ☐
89	HR099	Omit thumb piece			NC ☐
4- Compatible with Surge, Surge LT handrims					
AXLES					
90	AX01	Quick release axle		WAMU	STD ☐
91	AX03	Quad release axle		WAMU	\$175 ☐
REAR WHEEL CAMBER					
92	CM00	0° camber			STD ☐
93	CM02	2° camber			NC ☐
94	CM04	4° camber			NC ☐
95	CM06	6° camber			NC ☐
96	CM08	8° camber			NC ☐
CENTER OF GRAVITY					
Specify the position of rear wheel axle plate: Measure distance from front of back post to center of axle receiver. 0 - 6" range Set At: _____ (inches) <i>If not indicated, default setting is at 2,5"</i> <i>If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.</i>					
WHEEL LOCK					
97	WL08	Newton wheel lock - push to lock			STD ☐
98	WL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)		WAMR	\$180 ☐
99	WL13 WL14	6" removable extension handle			\$70 ☐
100	WL05	Scissor Lock			\$180 ☐
101	WL99	Omit wheel locks			NC ☐
Armrests					
"T" ARMREST ¹			Height Adjustment		
102	AR03	10" Desk Length	8 - 12" STD ☐	10 - 14" ☐	\$370 ☐
103	AR04	14" Full Length	8 - 12" STD ☐	10 - 14" ☐	\$370 ☐
SWING AWAY ARMREST ²			Height Adjustment		
104	AR05	Tubular arm pad	8 - 11"		\$305 ☐
105	AR06	10" Desk length arm pad	8 - 11"		\$305 ☐
106	AR99	Omit armrest			NC ☐
1 - May have a 1/2 " adjustment restriction on both limits of the range			2 -not available with the Ultralight non-locking backrest		

SIDE GUARDS <i>(not available with "T" Armrest)</i>						
107	SG01	Plastic side guards <i>(removable)</i>		WAMV	\$230	☐
108	SG03	Carbon fiber side guards <i>(removable)</i>		WAMV	\$305	☐
109	SG06	Plastic fender side guards <i>(removable)</i> ¹		WAMV	\$330	☐
110	SG04	Carbon fiber fender side guards <i>(removable)</i>		WAMV	\$495	☐
111	SG99	Omit side guards			NC	☐
		SIDE GUARDS OPTION				
112	SG07	Side guard clip <i>(pair)</i>			NC	☐

1 - Not compatible with 20" and 22" rear wheels

POSITIONING BELT						
113	SB01	Velcro adjustable	Standard 35" ☐ Long 60" ☐		\$47	☐
114	SB02	Auto buckle	Standard 35" ☐ Long 60 " ☐		\$57	☐
115	SB07	Hardware attachment for 2 point seatbelt (1 pair)			NC	☐
116	SB03	Additional hardware for 2 extra attachment points <i>(needed for the 4 point seat belt) (1 pair)</i>			\$14	☐
117	SB04	BodyPoint Hip Belt - 2 point non padded	Medium 56" ☐ Long 61" ☐		\$85	☐
118	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐ Medium 56" ☐ Long 61" ☐	SERND0030	\$175	☐
119	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐ Medium 56" ☐ Long 61" ☐	SERND0020	\$200	☐
120	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" ☐ Long 24" ☐		\$335	☐
121	SB99	Omit positioning belt			STD	☐

ACCESSORIES						
122	AC03	Anti-tippers			\$120	☐
123	AC05	WC-19 Transit Tie Down w/Q'straint belt ¹	Short belt 20-41" ☐ Long belt 20-46" ☐		\$355	☐
124	AC08	Spoke guards ²	Available: 20", 22", 24", 26"	WAMP	\$150	☐
125	AC10	Seat pouch			\$98	☐
126	AC30	Newton santizing multi-purpose pouch			\$72	☐
127	AC26	BodyPoint Mobility Bag			\$140	☐
128	AC23	Newton Backpack			\$125	☐
129	AC13	Tool kit			\$52	☐
130	AC16	Motion Composites Polo Shirt	W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐		\$35	☐


1 -May interfere with the swing away armrest and the Clamp on height adjustable push handle

2 - Not available with Mag wheels and Spinergy

COLOR						
131	FC	Highlight color:	(10) Ferrari Red ☐ (11) White ☐ (12) Black ☐ (13) Sunkissed Orange ☐ (14) Sky Blue ☐ (15) Fuchsia ☐ (16) Acid Green ☐		NC	☐

Notes :

BioFit™ CUSHION

	MODEL #	Fits Width	Depth	\$	✓	MODEL #	Fits Width	Depth	\$	✓
	NXCABF1616	15"-16"	16"	595		NXCABF2016	19" - 20"	16"	595	
	NXCABF1618	15"-16"	18"			NXCABF2018	19" - 20"	18"		
	NXCABF1620	15"-16"	20"			NXCABF2020	19" - 20"	20"	625	
	NXCABF1816	17"-18"	16"			NXCABF2216	21" - 22"	16"		
	NXCABF1818	17"-18"	18"			NXCABF2218	21" - 22"	18"		
	NXCABF1820	17"-18"	20"			NXCABF2220	21" - 22"	20"		

THORACIC SHAPE



Optima™
THORACIC BACK SUPPORT

CARBON FIBRE SHELL (WEIGHT SAVINGS – 1lb)

Model # Carbon Fibre	Model # Aluminum	Fits Width	Height	\$ Carbon	✓	\$ Aluminum	✓
	NOACB1307	13"-15"	7"				
	NOACB1309	13"-15"	9"				
	NOACB1313	13"-15"	13"				
	NOACB1316	13"-15"	16"				
	NOACB1407	14"-16"	7"				
NOCFACB1409	NOACB1409	14"-16"	9"	1095			
NOCFACB1413	NOACB1413	14"-16"	13"	1095			
NOCFACB1416	NOACB1416	14"-16"	16"	1095			
	NOACB1507	15"-17"	7"				
NOCFACB1509	NOACB1509	15"-17"	9"	1095			
NOCFACB1513	NOACB1513	15"-17"	13"	1095			
NOCFACB1516	NOACB1516	15"-17"	16"	1095			
	NOACB1607	16"-18"	7"				
NOCFACB1609	NOACB1609	16"-18"	9"	1095			
NOCFACB1613	NOACB1613	16"-18"	13"	1095			
NOCFACB1616	NOACB1616	16"-18"	16"	1095			
	NOACB1707	17"-19"	7"				
NOCFACB1709	NOACB1709	17"-19"	9"	1095			
NOCFACB1713	NOACB1713	17"-19"	13"	1095			
NOCFACB1716	NOACB1716	17"-19"	16"	1095			
	NOACB1807	18"-20"	7"				
NOCFACB1809	NOACB1809	18"-20"	9"	1095			
NOCFACB1813	NOACB1813	18"-20"	13"	1095			
NOCFACB1816	NOACB1816	18"-20"	16"	1095			
	NOACB2007	20"-22"	7"				
NOCFACB2009	NOACB2009	20"-22"	9"	1095			
NOCFACB2013	NOACB2013	20"-22"	13"	1095			
NOCFACB2016	NOACB2016	20"-22"	16"	1095			

690



Optima™ DEEP
THORACIC BACK SUPPORT

CARBON FIBRE SHELL (WEIGHT SAVINGS – 1lb)

Model # Carbon Fibre	Model # Aluminum	Fits Width	Height	\$ Carbon	✓	\$ Aluminum	✓
	NOPCB1307	13"-15"	7"				
	NOPCB1309	13"-15"	9"				
	NOPCB1313	13"-15"	13"				
	NOPCB1316	13"-15"	16"				
	NOPCB1407	14"-16"	7"				
NOCFPCB1409	NOPCB1409	14"-16"	9"	1195			
NOCFPCB1413	NOPCB1413	14"-16"	13"	1195			
NOCFPCB1416	NOPCB1416	14"-16"	16"	1195			
	NOPCB1507	15"-17"	7"				
NOCFPCB1509	NOPCB1509	15"-17"	9"	1195			
NOCFPCB1513	NOPCB1513	15"-17"	13"	1195			
NOCFPCB1516	NOPCB1516	15"-17"	16"	1195			
	NOPCB1607	16"-18"	7"				
NOCFPCB1609	NOPCB1609	16"-18"	9"	1195			
NOCFPCB1613	NOPCB1613	16"-18"	13"	1195			
NOCFPCB1616	NOPCB1616	16"-18"	16"	1195			
	NOPCB1707	17"-19"	7"				
NOCFPCB1709	NOPCB1709	17"-19"	9"	1195			
NOCFPCB1713	NOPCB1713	17"-19"	13"	1195			
NOCFPCB1716	NOPCB1716	17"-19"	16"	1195			
	NOPCB1807	18"-20"	7"				
NOCFPCB1809	NOPCB1809	18"-20"	9"	1195			
NOCFPCB1813	NOPCB1813	18"-20"	13"	1195			
NOCFPCB1816	NOPCB1816	18"-20"	16"	1195			
	NOPCB2007	20"-22"	7"				
NOCFPCB2009	NOPCB2009	20"-22"	9"	1195			
NOCFPCB2013	NOPCB2013	20"-22"	13"	1195			
NOCFPCB2016	NOPCB2016	20"-22"	16"	1195			

790

BACK SUPPORT HARDWARE OPTIONS

All backs come STANDARD with QuickRelease Mounting Hardware. The following are fixed hardware options available at no charge when ordered with a back support.

Model #	DESCRIPTION	WEIGHT OF HARDWARE ONLY	✓
	QuickFit Quick Release Hardware System	STANDARD (2 LBS)	STANDARD
NIB-1642	Micro Mount 2 point fixed hardware. Back width needs to match chair width	Weight savings of 1.4 lb over STANDARD	
NIB-1207	2 point fixed hardware. Width expandable to fit 3 sizes of chair	Weight savings of 1.0 lb over STANDARD	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

 For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

ver: jw-05-14-21