



1- Purchasing information

P.O. number: * _____ Date: _____

Company Name: * _____

Purchasing contact: * _____

Phone: * _____

Email: * _____

Contact preference: * Phone ☐ Email ☐

Marked for: _____

3- Shipping information

Ship to: * _____

Name: _____

Address: _____

City: _____ State: _____

Zip Code: _____

Phone: _____

2- Clinician information

Clinician: * _____

Phone: * _____

Email: * _____

Contact preference: * Phone ☐ Email ☐

		MA043 ☐	MA044 ☐	MA045 ☐	MA046 ☐	MA047 ☐	MA048 ☐	MA090 ☐		
		English Manual (STD)	German Manual	Italian Manual	Dutch Manual	Spanish Manual	French Manual	Portuguese Manual		
*Required										taxes not included
MODEL										USD
WC11	☐	Ultralight Aluminum APEX Rigid Wheelchair (120kg /265 lb weight limit)		Triple-buttet aluminum frame / Aluminum back canes and back rigidizer / Aluminum seat rigidizer / Ultralight carbon fiber camber tube and rear frame rigidizer						\$2,895
WC21	☐	Apex Carbon fiber frame upgrade (High Gloss Carbon Fiber)		Ultralight carbon fiber frame, back canes, back rigidizer, seat rigidizer, camber tube and rear frame rigidizer.						\$1,995
EVO-R	☐	Evo Rigid - Chair resizing service*		To be used within 24 months of chair purchase					Aluminum Carbon	\$750 \$1,000

*By selecting this option, I agree to the terms and conditions of the EVO-RIGID program.

https://www.motioncomposites.com/en_us/terms-and-conditions

Frame																		
SEAT WIDTH																		
SW	☐	30.48 (12")	☐	33.02 (13")	☐	35.56 (14")	☐	38.10(15")	☐	40.64(16")	☐	43.18 (17")	☐	45.72 (18")	☐		NC	
	☐	48.26 (19")	☐	50.80 (20")	☐											\$535		
SEAT DEPTH																		
SD	☐	30.48 (12")	☐	33.02 (13")	☐	35.56 (14")	☐	38.10(15")	☐	40.64(16")	☐	43.18 (17")	☐	45.72 (18")	☐		NC	
	☐	48.26 (19")	☐	50.80 (20")	☐											\$515		
FRAME LENGTH																		
SP01	☐	Standard frame length configuration ¹																STD
SP02	☐	Frame length 2.54cm (1") longer than standard configuration																NC
SP03	☐	Frame length 5.08cm (2") longer than standard configuration																NC
SP04	☐	Frame length 2.54cm (1") shorter than standard configuration																NC
SP05	☐	Frame length 5.08cm (2") shorter than standard configuration																NC
TIE-DOWN																		
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ^{2,3}				Short belt 50-104 (20-41")									\$365			
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ^{2,3}				Long belt 50-116 (20-46")									\$365			
AC35	☐	WC-19 Transit Tie Down without belt ³																\$335

1- Distance of 4.37" from end of seat upholstery to bend

2- Not available with depth adjustable back; 3- May interfere with the swing away armrest and the clamp on height adjustable push handle; Add a 3" of frame length behind backrest

4. Width calculator

FRONT Seat-to-floor height			
Caster sizes		Height available +- 1/4"	
☐	3" Caster	35.6cm to 53.3cm (14" to 21")	
☐	4" Caster	36.8cm to 53.3cm (14.5" to 21")	
☐	5" Caster	39.4cm to 53.3cm (15.5" to 21")	
☐	6" Caster	41.9cm to 53.3cm (16.5" to 21") (N/A with 90° front frame)	
		Specify desired height: (inches) _____	

Frog Legs FRONT Seat-to-floor height			
Caster Sizes		Height available +- 1/4"	
☐	3" Caster	38.7cm to 53.3cm (15.25" to 21")	
☐	4" Caster	40.0cm to 53.3cm (15.75" to 21")	
☐	5" Caster	42.54cm to 53.3cm (16.75" to 21")	
☐	6" Caster	45.08cm to 53.3cm (17.75" to 21") (N/A with 90° front frame)	
		Specify desired height: (inches) _____	



HANDRIMS

HANDRIM			Available		
HR01	☐	Aluminum anodized	20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²	20", 22", 24", 25"		\$80
HR02	☐	Plastic coated ¹	20", 22", 24", 25", 26"		\$145
HR03	☐	High friction coated ^{1,2}	22", 24", 25", 26"		\$380
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ² See page 6 for the color selection and price	20", 22", 24", 25", 26"		
HR28	☐	Surge ^{1,2,3} (1 1/2" Oval handrim w/ Gripton high friction strip)	22", 24", 25", 26"		\$425
HR05	☐	Surge LT ^{1,2,3} (1 1/8" Oval handrim w/ Gripton high friction strip)	22", 24", 25", 26"		\$465
HR04	☐	Natural Fit ^{1,2,3}	22", 24", 25", 26"		\$395
HR62	☐	Simi H ^{1,2} (1 3/8" Oval handrim with high friction strip)	22", 24", 25", 26"		\$395
HR63	☐	Nova H ^{1,2} (1 1/4" Oval handrim with high friction strip)	22", 24", 25", 26"		\$425
HR64	☐	Optimum H ^{1,2} (1 1/8" heart-shaped handrim)	22", 24", 25", 26"		\$325
HR99	☐	Omit handrim			NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels

3- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

HANDRIMS OPTIONS

ASSEMBLY POSITION					
HR29	☐	Standard Handrim position	Full tab		STD
HR030	☐	Narrow Handrim position no cut	Full tab		NC
HR30	☐	Narrow Handrim position with cut ³	Short tab 2 positions		\$32
HR31	☐	Super Narrow Handrim ³	Short tab 2 positions STD ☐ Short tab 1 position ☐		\$32
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	☐	Standard grip			\$95
HR27	☐	Super grip			\$150
HR099	☐	Omit thumb piece			NC

3- Only available with Aluminum and Air Grip handrims;

AXLES					
AX01	☐	Quick release axle			STD
AX03	☐	Quad release axle			\$135
REAR WHEEL CAMBER					
CM00	☐	0° camber			STD
CM02	☐	2° camber			NC
CM04	☐	4° camber			NC
CM06	☐	6° camber			NC
CM08	☐	8° camber			NC

CENTER OF GRAVITY

Specify the position of rear wheel axle plate:


Measure distance from front of back post to center of axle receiver.


0 - 6" range

Set At: (inches)

If not indicated, default setting is at 2.5"

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG





WHEEL LOCK					
WL08	☐	Newton wheel lock - push to lock			STD
WL16	☐	Aluminum push to lock with extensions			\$83
WL02	☐	Pull to lock			NC
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)			\$165
WL13	☐	6" removable extension handle ¹			\$83
WL14					
WL05	☐	Scissors lock			\$155
WL99	☐	Omit wheel locks			NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

FOOTPLATES					
FP12	☐	Tubular open loop ¹			STD
FP13	☐	Tubular with ABS cover			NC
FP31	☐	Newton adjustable angle - Aluminum (one piece)			\$105
FP32	☐	Newton adjustable angle - Carbon fiber (one piece with rubber inserts)			\$190
FP14	☐	Flip back, adjustable angle aluminum			\$210
FP15	☐	Flip back, adjustable angle carbon fiber			\$300
FR25	☐	High mount with adjustable angle footplate (Aluminum - one piece)			\$180
FR26	☐	High mount with adjustable angle footplate (Carbon fiber with rubber inserts - one piece)			\$250

FOOTPLATES ACCESSORIES					
FP16	☐	5.08cm (2") riser (only available with aluminum angle adjustable footplate) ²			\$78
FP09	☐	Calf strap			\$43
FP30	☐	Bodypoint Aeromesh Padded calf strap			\$125

1 - Available on: all straight frames; 2,5cm (1") taper with seat width of 14" and greater; 5cm (2") taper with seat width of 16" and greater; 2 - Not available with high mount

SEATING AND BACKREST					
BACK UPHOLSTERY					
BH03	☐	Slip-on Syntech-Suede back upholstery			STD
BH04	☐	Adjustable tension Syntech-Suede back upholstery			\$180
BH99	☐	Omit Back upholstery			NC
SEAT CUSHION					
SC01	☐	5.08cm (2") Cushion			\$194
SC02	☐	7.62cm (3") Cushion			\$200
SC99	☐	Omit seat cushion			STD
FOLD DOWN BACKREST					
BC20	☐	Locking backrest			STD
BC21	☐	Ultralight non locking backrest (not locking in upright and down position) ²			\$105
FIXED HEIGHT, ADJUSTABLE ANGLE FOLD DOWN BACKREST (Not available with push handles)					
BC	☐	Fixed height without push handle (30) 22.8cm/9" ☐ (31) 24.1cm/9.5" ☐ (32) 25.4cm/10" ☐ (33) 26.7cm/10.5" ☐ (34) 27.9 cm/11" ☐ (35) 29.2cm/11.5" ☐ (36) 30.5cm/12" ☐ (37) 31.7cm/12.5" ☐ (38) 33cm/13" ☐ (39) 34.3cm/13.5" ☐ (40) 35.5cm/14" ☐ (41) 35.8cm/14.5" ☐ (42) 38.1cm/15" ☐ (43) 39.4cm/15.5" ☐ (44) 40.6cm/16" ☐ (45) 41.9cm/16.5" ☐ (46) 43.2cm/17" ☐ (47) 44.5cm/17.5" ☐ (48) 45.7cm/18" ☐ (49) 47cm/18.5" ☐ (50) 48.3/19" ☐ (51) 49.5cm/19.5" ☐ (52) 50.8cm/20" ☐ (53) 52.7cm/20.5" ☐ (54) 53.3cm/21" ☐			NC
HEIGHT AND ANGLE ADJUSTABLE FOLD DOWN BACKREST WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)					
BC	☐	Straight back canes (22) 22.8 - 30.5cm ☐ (23) 30.5 - 38.1cm ☐ (24) 38.1 - 45.7cm ☐ (25) 45.7 - 53.3cm ☐			NC
Set height to: *		(in increments of 1/2", measured from seat sling to top of back upholstery) * If not indicated, default height setting is at the lowest position			
		Select back cane angle: ¹ 101° ☐ 99° ☐ 96° ☐ 93° ☐ 90° ☐ 87° ☐ 83° ☐ 80° ☐			STD

1 - If not indicated default angle setting is at 93°; 2- Maximum back cane height of 19.5" with aluminum chair

PUSH HANDLES OPTIONS (ONE option at a time)					
PH01	☐	Push handles			STD
PH05	☐	Fold down push handles			\$150
PH06-920	☐	Clamp on height adjustable push handles 22.9-50.80cm (9 - 20") ¹			\$330
PH06-923	☐	Clamp on height adjustable push handles 22.9-58.4cm (9 - 23") ¹			\$330
PH02	☐	Omit push handles			NC
STROLLER HANDLE					
PH04	☐	Stroller handle (Push handles included)			\$175

1- Not compatible with swing away armrests; Not available with 9" to 12" back canes

BACKRESTS OPTIONS					
BC61	☐	Depth adjustable back - 2.54cm (1") depth adjustment ^{1,4}			\$195
BC62	☐	Depth adjustable back - 5.08cm (2") depth adjustment ^{2,4}			\$195
BC56	☐	Depth adjustable back - 7.62cm (3") depth adjustment ^{3,4}			\$195
BC57	☐	Standard 7.62cm (3") deep rigidizer bar			STD
BC55	☐	Four inch deep Aluminum rigidizer bar			NC
BC88	☐	Custom height rigidizer bar set at _____ (inches) Measured from the top of the seat sling to the center of the rigidizer - standard height is 5,7" Lowest position is 4,25"			\$135

If used in combinaison with WC-19: 1- Maximum depth of 19"; 2- Maximum depth of 18"; 3- Maximum depth of 17"; 4- Only compatible with WC-19 transit tie down **without** Q'Straint belt

ARMREST					
HEIGHT ADJUSTABLE "T" ARMREST ¹					
AR13-610	☐	25.4cm (10") Desk Length armpad ³	15.25 - 25cm (6" - 10")		\$360
AR03-812	☐	25.4cm (10") Desk Length armpad	20.32 - 30.48cm (8" - 12")		\$360
AR03-1014	☐	25.4cm (10") Desk Length armpad	25.40 - 35.56cm (10" - 14")		\$360
AR14-610	☐	35.56cm (14") Full Length armpad ³	15.25 - 25cm (6" - 10")		\$360
AR04-812	☐	35.56cm (14") Full Length armpad	20.32 - 30.48cm (8" - 12")		\$360
AR04-1014	☐	35.56cm (14") Full Length armpad	25.40 - 35.56cm (10" - 14")		\$360
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²					
AR05	☐	Tubular armrest - Short pad	20.32 - 27.94cm (8 - 11")		\$250
AR23	☐	Tubular armrest - Long pad	20.32 - 27.94cm (8 - 11")		\$250
AR06	☐	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")		\$250
AR17	☐	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")		\$250
ARMREST OPTIONS					
AR12	☐	14" Full Length armpad - Gel Ovation (pair) ⁴			\$165
AR99	☐	Omit armrest			NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- not available with the Ultralight non-locking backrest; 3- May have some restriction with seat-to-floor height and rear wheel

4- Only available with "T" full length armrest; Not available with swing away and desk lenght armrests

ACCESSORIES							
SIDE GUARDS (not available with "T" Armrest)							
SG01	☐	Plastic side guards (removable)	Straight bracket -STD	Offset bracket ¹			\$210
SG03	☐	Carbon fiber side guards (removable)	Straight bracket -STD	Offset bracket ¹			\$315
SG06	☐	Plastic fender side guards (removable)	Straight bracket -STD	Offset bracket ¹			\$295
SG04	☐	Carbon fiber fender side guards (removable)	Straight bracket -STD	Offset bracket ¹			\$550
SG99	☐	Omit side guards					NC

1- only compatible with 0° and 2° camber; Increase width by 1/2" per side for a total of 1"

POSITIONING BELT							
SB01-S	☐	Velcro adjustable - Small 88.9cm (35") ¹					\$50
SB01-L	☐	Velcro adjustable - Long 152.40cm (60")					\$50
SB02-S	☐	Auto buckle - Small 88.9cm (35") ¹					\$95
SB02-M	☐	Auto buckle - Medium 121.92cm (48")					\$95
SB02-L	☐	Auto buckle - Long 152.40cm (60")					\$95
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")					\$95
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")					\$95
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")					\$185
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")					\$185
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")					\$185
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")					\$205
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")					\$205
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")					\$205
SB03	☐	Hardware attachment for 2 extra points (1 pair)					\$17
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66cm (18")					\$310
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")					\$310
SB99	☐	Omit positioning belt					STD

1- Not recommended for seat width of 45.72cm (18") and more

ACCESSORIES							
AC03	☐	Anti-tippers					\$155
AC06	☐	Cane and crutch holder	Left ☐	Right ☐			\$130
AC08	☐	Spoke guards ¹	Available: 20", 22", 24", 25" 26"				\$125
AC10	☐	Seat pouch					\$83
AC30	☐	Newton multi-purpose pouch					\$62
AC26	☐	BodyPoint Mobility Bag					\$135
AC23	☐	Newton Backpack					\$130
AC13	☐	Tool kit					\$50

1- Not available with Mag wheels and Spinergy; the 24" spoke guard will be installed on the 25" wheel

Colors																				
Options ¹	Black ⁴	High Gloss Carbon Fiber ⁴	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
																				
AirGrip	☐ \$125		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Camber Tube	☐ STD/NC			☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
																				
Highlight ³	☐ STD		☐					☐		☐	☐			☐			☐			
FC99 ²	☐	Custom Color	Specify PPG Paintit color: additional 10 days lead-time																	\$390

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- Only available with carbon fiber upgrade; 3- To see different highlight, visit our website: https://www.motioncomposites.com/en_intl/products/rigid-wheelchairs/apex-c

4- Black frame will be STD for the aluminum/High Gloss Carbon Fiber will be STD for the Carbon frame upgrade



CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the
EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The **EVO Custom Request Form** is the best place to start.

Notes :