

Two views of a black and red wheelchair. The left view shows the wheelchair from a side profile, highlighting the large rear wheels with red rims and the smaller front casters. The right view shows the wheelchair from a slightly different angle, emphasizing the red frame and the black seat and backrest.

1- Purchasing information

Marked for:

2- ATP/Clinician information

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Phone: _____

** Required*

MA089 ☐


Spanish Manual

Taxes not included

MODEL												USD
WC11	☐	Ultralight Aluminum APEX Rigid Wheelchair (265 lb weight limit)				Tripple-buttet aluminum frame / Aluminum back canes and back rigidizer / Aluminum seat rigidizer / Ultralight carbon fiber camber tube and rear frame rigidizer				K0005	\$2 595	
WC21	☐	Apex Carbon fiber frame upgrade (High Gloss Carbon Fiber)				Ultralight carbon fiber frame, back canes, back rigidizer, seat rigidizer, camber tube and rear frame rigidizer.					\$2 087	
SEAT WIDTH												
SW	☐	12" ☐	13" ☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐		NC	
	☐	20" ☐								E2201	\$495	
SEAT DEPTH												
SD	☐	12" ☐	13" ☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐		NC	
	☐	20" ☐								E2203	\$475	
FRAME LENGTH												
SP01	☐	Standard frame length configuration ¹									STD	
SP02	☐	Frame length 1" longer than standard configuration									NC	
SP03	☐	Frame length 2" longer than standard configuration									NC	
SP04	☐	Frame length 1" shorter than standard configuration									NC	
SP05	☐	Frame length 2" shorter than standard configuration									NC	


1- Distance of 4.37" from end of seat upholstery to bend

FRONT Seat-to-floor height		
Caster sizes	Height available +/- 1/4"	
3" Caster	14" to 21"	Specify desired height: (inches)
4" Caster	14.5" to 21"	
5" Caster	15.5" to 21"	
6" Caster	16.5" to 21" (N/A with 90° front frame)	




CASTERS						
CA01-03	☐	Composite wheel w/ PU tire	3" x 1"			NC
CA01-04	☐		4" x 1"			NC
CA01-05	☐		5" x 1"			NC
CA01-06	☐		6" x 1"			NC
CA10	☐	Pneumatic casters	6" x 1 ¼"			\$160
CA14	☐	Newton UltraCaster Composite	4" x 1 ½ "			\$85
CA15	☐		5" x 1 ½ "			\$85
CA16	☐		6" x 1 ½ "			\$85
CA17	☐	Newton UltraCaster Aluminum	4" x 1 ½ "			\$185
CA18	☐		5" x 1 ½ "			\$195
CA19	☐		6" x 1 ½ "			\$205
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks			E1015	\$475

REAR Seat-to-floor height		
Wheel sizes	Height available +/- 1/4"	Specify desired height: (inches)
20" (451) wheel ☐	14" to 17.5"	
22" (501) wheel ☐	14" to 18.5"	
24" (540) wheel ☐	15" to 19.5"	
25" (559) wheel ☐	15.5" to 20"	
26" (590) wheel ☐	16" to 20"	



FRONT FRAME ANGLE				
FA70	☐	70° Front frame angle		NC
FA75	☐	75° Front frame angle		NC
FA80	☐	80° Front frame angle		NC
FA85	☐	85° Front frame angle		NC
FA90	☐	90° Front frame angle		NC

Seat upholstery to footplate length	
12" to 19" without high mount	Specify desired length: (inches) <i>Required</i>
7" to 13" with high mount ¹	



1- Choose high mount in the footplates section


Footrest width / Taper selection				
Seat Width	Straight	1" Taper	2" Taper	
Outside Measurement	2 ½" Narrower than Seat Width	4 ½" Narrower than Seat Width	6 ½" Narrower than Seat Width	
12"	9.5"	7.5"		
13"	10.5"	8.5"		
14"	11.5"	9.5"		
15"	12.5"	10.5"		
16"	13.5"	11.5"	9.5"	
17"	14.5"	12.5"	10.5"	
18"	15.5"	13.5"	11.5"	
19"	16.5"	14.5"	12.5"	
20"	17.5"	15.5"	13.5"	

Note: higher than 20" Front STFH will reduce the footrest width up to 1" with a 2" taper

FOOTPLATES				
FP12	☐	Tubular open loop ¹		STD
FP13	☐	Tubular with ABS cover	K0040	NC
FP31	☐	Adjustable angle - Aluminum (<i>one piece</i>)	K0040	\$95
FP32	☐	Newton adjustable angle - Carbon fiber (<i>one piece with rubber inserts</i>)	K0040	\$185
FP14	☐	Flip back, adjustable angle aluminum	K0040	\$195
FP15	☐	Flip back, adjustable angle carbon fiber		\$290
FR25	☐	High mount with adjustable angle footplate (<i>Aluminum - one piece</i>)	K0040	\$165
FR26	☐	High mount with adjustable angle footplate (<i>Carbon fiber with rubber inserts - one piece</i>)		\$245
FOOTPLATES ACCESSORIES				
FP16	☐	2" riser (<i>only available with aluminum angle adjustable footplate</i>) ²		\$75
FP09	☐	Calf strap	K0038	\$35
FP30	☐	Bodypoint Aeromesh Padded calf strap	K0038	\$115

1 - Available on: all straight frames; 1" taper with seat width of 14" and greater; 2" taper with seat width of 16" and greater

2 - Not available with high mount


Backrest												
FOLD DOWN BACKREST												
BC20	☐	Locking backrest										STD
BC21	☐	Ultralight non locking backrest (not locking in upright and down position)										\$100
FIXED HEIGHT, ADJUSTABLE ANGLE FOLD DOWN BACKREST (NO PUSH HANDLE)												
BC	☐	Fixed height without push handle		(30) 9" ☐	(31) 9.5" ☐	(32) 10" ☐	(33) 10.5" ☐	(34) 11" ☐		NC		
		(35) 11.5" ☐	(36) 12" ☐	(37) 12.5" ☐	(38) 13" ☐	(39) 13.5" ☐	(40) 14" ☐	(41) 14.5" ☐				
		(42) 15" ☐	(43) 15.5" ☐	(44) 16" ☐	(45) 16.5" ☐	(46) 17" ☐	(47) 17.5" ☐	(48) 18" ☐				
		(49) 18.5" ☐	(50) 19" ☐	(51) 19.5" ☐	(52) 20" ☐	(53) 20.5" ☐	(54) 21" ☐					
HEIGHT AND ANGLE ADJUSTABLE FOLD DOWN BACKREST WITH PUSH HANDLE OPTIONS												
BC	☐	Straight back canes										NC
		(22) 9 - 12" ☐	(23) 12 - 15" ☐	(24) 15 - 18" ☐	(25) 18 - 21" ☐							
Set height to: *		<i>(in increments of 1/2", measured from seat sling to top of back upholstery)</i> <i>* If not indicated, default height setting is at the lowest position</i>										
Select back cane angle: ¹		101° ☐	99° ☐	96° ☐	93° ☐	90° ☐	87° ☐	83° ☐	80° ☐			STD

1 - If not indicated default angle setting is at 93°

PUSH HANDLE OPTIONS				
PH01	☐	Push handles		STD
PH05	☐	Fold down push handles ¹		\$126
PH06-920	☐	Clamp on height adjustable push handles 9 - 20" ^{1,2}		\$315
PH06-923	☐	Clamp on height adjustable push handles 9 - 23" ^{1,2}		\$315
PH02	☐	Omit push handles		NC

1- 265 pounds weight limit

2 - Not compatible with swing away armrests;

BACKRESTS OPTIONS				
BC61	☐	Depth adjustable back - 1" depth adjustment ^{1,4}		\$190
BC62	☐	Depth adjustable back - 2" depth adjustment ^{2,4}		\$190
BC56	☐	Depth adjustable back - 3" depth adjustment ^{3,4}		\$190
BC57	☐	Standard 3 inch deep rigidizer bar		STD
BC55	☐	Four inch deep Aluminum rigidizer bar		NC
BC88	☐	Custom height rigidizer bar set at _____ (inches) Measured from the top of the seat sling to the center of the rigidizer - standard height is 5,7" Lowest position is 4,25"		\$130

1- Maximum depth of 19"; 2- Maximum depth of 18"; 3- Maximum depth of 17"; 4- Not compatible with WC-19 transit tie down option that includes the Q'Straint belt;

BACK UPHOLSTERY						
BH03	☐	Slip-on Syntech-Suede back upholstery				STD
BH04	☐	Adjustable tension Syntech-Suede back upholstery			E2611	\$300
BH99	☐	Omit Back upholstery				NC
	☐	NXT back Support			See price list	
SEAT CUSHION						
SC01	☐	2" Cushion			E2601	\$165
SC02	☐	3" Cushion			E2601	\$185
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See price list	

Rear wheels						
REAR WHEELS				Available		
RW01	☐	Mag ¹		20", 22" 24"		\$55
RW02	☐	Newton One - spoke wheel ²		20", 22" 24"		STD
RW03	☐	Newton Gravity - ultralight wheel		20", 22", 24", 25", 26"		\$315
RW04	☐	Spinergy Spox Blk STD ☐ Wht ³ ☐ Blue ³ ☐ Red ³ ☐ Yel ³ ☐		22", 24", 25", 26"		\$835
RW05	☐	Spinergy LX Blk STD ☐ Pur ³ ☐ Wht ³ ☐ Blue ³ ☐ Red ³ ☐ Yel ³ ☐ Pink ³ ☐ Grn ³ ☐ Org ³ ☐		22", 24", 25", 26"		\$1 045
RW07	☐	Spinergy CLX ³ Blk STD ☐ Pur ³ ☐ Wht ³ ☐ Blue ³ ☐ Red ³ ☐ Yel ³ ☐ Pink ³ ☐ Grn ³ ☐ Org ³ ☐		24", 25"		\$1 555
TIRES				Available		
TI01	☐	Soft Urethane 1 3/8 "		20", 22" 24"	E2218	NC
TI02	☐	Soft Urethane knobby 1 3/8"		22", 24", 25", 26"	E2218	\$95
TI04	☐	Soft Urethane - <i>Shox 1"</i>		22", 24", 25", 26"	E2218	\$115
TI05	☐	Pneumatic 1 3/8 "		20", 22", 24", 26"	E2211/E2212	STD
TI08	☐	Pneumatic with Airless insert 1 3/8" ⁴		22", 24", 26"	E2211/E2213	\$80
TI06	☐	Pneumatic HP - 1"		20", 22", 24", 25", 26"	E2211/E2212	\$95
TI07	☐	Pneumatic HP - <i>Puncture resistant, black 1"</i>		24", 25", 26"	E2211/E2212	\$125

1- Not compatible with HP tires; Only available with aluminum or plastic coated handrim in standard position; 2- Not compatible with Shox tires
3- color spokes require an additional 10 days lead-time; 4- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims						
HANDRIM				Available		
HR01	☐	Aluminum anodized		20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²		20", 22", 24", 25"		\$70
HR02	☐	Plastic coated ¹		20", 22", 24", 25", 26"		\$100
HR03	☐	Q-Grip High friction Neoprene coated ^{1,2}		22", 24", 25", 26"		\$370
HR07	☐	Newton Air Grip ² <i>Black</i> ☐		20", 22", 24", 25", 26"		\$100
HR	☐	Newton Air Grip ² (08) <i>Sunkissed Orange</i> ☐ (09) <i>Burgundy</i> ☐ (10) <i>Sapphire Blue</i> ☐ (11) <i>Monster Green</i> ☐ (12) <i>Charcoal</i> ☐ (13) <i>Fuchsia</i> ☐ (14) <i>Steel Blue</i> ☐ (15) <i>Ferrari Red</i> ☐ (16) <i>Acid Green</i> ☐ (17) <i>White</i> ☐ (18) <i>Sky Blue</i> ☐		20", 22", 24", 25", 26"		\$250
HR28	☐	Surge ^{1,2} (1-1/2" Oval handrim w/ Gripton high friction strip) **Limited availability		22", 24", 25", 26"		\$395
HR05	☐	Surge LT ^{1,2} (1-1/8" Oval handrim w/ Gripton high friction strip) **Limited availability		22", 24", 25", 26"		\$425
HR04	☐	Natural Fit ^{1,2} **Limited availability		22", 24", 25", 26"		\$325
HR99	☐	Omit handrim				NC


1- Only standard handrim position available 2- Not available with Mag wheels

Handrims option						
Assembly position for aluminum and Newton AirGrip handrims						
HR29	☐	Standard Handrim position <i>Full tab</i>				STD
HR30	☐	Narrow Handrim position <i>Short tab 2 positions</i>				\$27
HR31	☐	Super Narrow Handrim <i>Short tab 2 positions STD</i> ☐ <i>Short tab 1 position</i> ☐				\$27
Thumb piece ⁴ selection for Natural Fit handrims						
HR23	☐	Standard grip **Limited availability				\$90
HR27	☐	Super grip **Limited availability				\$145
HR099	☐	Omit thumb piece				NC

4- Compatible with Surge, Surge LT handrims

AXLES					
AX01	☐	Quick release axle			STD
AX03	☐	Quad release axle		K0108	\$125
REAR WHEEL CAMBER					
CM00	☐	0° camber			STD
CM02	☐	2° camber			NC
CM04	☐	4° camber			NC
CM06	☐	6° camber			NC
CM08	☐	8° camber			NC

CENTER OF GRAVITY	
Specify the position of rear wheel axle plate: Measure distance from front of back post to center of axle receiver.	
0 - 6" range	Set At: _____ (inches)
<i>If not indicated, default setting is at 2.5"</i> <i>Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG</i>	



WHEEL LOCK					
WL08	☐	Newton wheel lock - push to lock			STD
WL04	☐	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)		E0974	\$150
WL13 WL14	☐	6" removable extension handle		E0961	\$75
WL05	☐	Scissors lock			\$150
WL99	☐	Omit wheel locks			NC

Armrests					
HEIGHT ADJUSTABLE "T" ARMREST <i>(improved)</i> ¹					
AR13-610	☐	10" Desk length armpad ³	6" - 10"	New	\$325
AR03-812	☐	10" Desk length armpad	8" - 12"		\$325
AR03-1014	☐	10" Desk length armpad	10" - 14"		\$325
AR14-610	☐	14" Full length armpad ³	6" - 10"	New	\$325
AR04-812	☐	14" Full length armpad	8" - 12"		\$325
AR04-1014	☐	14" Full length armpad	10" - 14"		\$325
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²					
AR05	☐	Tubular armrest without plastic side guards	8" - 11"		E0973 \$225
AR19	☐	Tubular armrest with plastic side guards	8" - 11"		E0973 \$385
AR06	☐	10" Desk length flat armpad without plastic side guards	8" - 11"		E0973 \$225
AR20	☐	10" Desk length flat armpad with plastic side guards	8" - 11"		E0973 \$385
AR17	☐	14" Full length flat armpad without plastic side guards	8" - 11"	New	E0973 \$225
AR21	☐	14" Full length flat armpad with plastic side guards	8" - 11"		E0973 \$385
AR99	☐	Omit armrest			NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- not available with the Ultralight non-locking backrest; 3- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS <i>(not available with "T" Armrest)</i>					
SG01	☐	Plastic side guards <i>(removable)</i>		K0108	\$195
SG03	☐	Carbon fiber side guards <i>(removable)</i>		K0108	\$290
SG06	☐	Plastic fender side guards <i>(removable)</i>			\$285
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>		K0108	\$540
SG99	☐	Omit side guards			NC
SIDE GUARDS OPTION					
SG07	☐	Side guard clip <i>(pair)</i>			NC

POSITIONING BELT					
SB01-S	☐	Velcro adjustable - Standard 35"		E0978	\$45
SB01-L	☐	Velcro adjustable - Long 60"		E0978	\$45
SB02-S	☐	Auto buckle - Standard 35"	Improved	E0978	\$85
SB02-L	☐	Auto buckle - Long 60"	Improved	E0978	\$85
SB07	☐	Hardware attachment for 2 point seatbelt (1 pair)			NC
SB03	☐	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)			\$16
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"		E0978	\$80
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"		E0978	\$80
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"		E0978	\$170
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"		E0978	\$170
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"		E0978	\$170
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"		E0978	\$190
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"		E0978	\$190
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"		E0978	\$190
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"		K0108	\$300
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"		K0108	\$300
SB99	☐	Omit positioning belt			STD


ACCESSORIES					
AC03	☐	Anti-tippers		E0971	\$140
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ^{1,2} Short belt 20"-41"			\$355
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ^{1,2} Long belt 20"-46"			\$355
AC35	☐	WC-19 Transit Tie Down without belt²	New		\$325
AC06	☐	Cane and crutch holder Left ☐ Right ☐	New	E2207	\$195
AC08	☐	Spoke guards ³ Available: 20", 22", 24", 26"		K0065	\$115
AC10	☐	Seat pouch			\$80
AC30	☐	Newton multi-purpose pouch			\$60
AC26	☐	BodyPoint Mobility Bag			\$130
AC23	☐	Newton Backpack			\$125
AC13	☐	Tool kit			\$50
AC16	☐	Motion Composites Polo Shirt W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐			\$35
AC15	☐	Military Patches Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐ Combat Wounded Purple Heart ☐			\$95

1- Not available with dept adjustable back; 2- May interfere with the swing away armrest and the clamp on height adjustable push handle; 3- Not available with Mag wheels and Spinergy

HIGHLIGHT COLOR					
FC	☐	(10) Ferrari Red ☐ (11) White ☐ (12) Black ☐ (13) Sunkissed Orange ☐ (14) Sky Blue ☐ (15) Fuchsia ☐ (16) Acid Green ☐			NC

Notes :

BioFit™ CUSHION		WIDTH									
DEPTH		12"	13"	14"	15"	16"	17"	18"	19"	20"	
	12"	☐ \$425	☐ \$425	☐ \$425							
	14"	☐ \$425	☐ \$425	☐ \$425							
	15"		☐ \$425	☐ \$425	☐ \$425						
	16"		☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	
	17"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	
	18"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	
	19"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	
	20"				☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	☐ \$455	
		Optional Incontinent / Infection Control Outer Cover NXCCAF__ __ __ V ☐ \$N/C									


NXCBF__ __ __
(W" x D")

HCPCS CODE: E2607-E2608

THORACIC SHAPE	Optima™ THORACIC BACK SUPPORT						Optima™ DEEP THORACIC BACK SUPPORT					
	CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)						CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)					
	HCPCS CODE: E2620-E2621						HCPCS CODE: E2620-E2621					
	Fits Width	Height	Model # CARBON FIBER	\$ CARBON	Model # Aluminum	\$ Aluminum	Fits Width	Height	Model # CARBON FIBER	\$ CARBON	Model # Aluminum	\$ Aluminum
	13"-15"	7"			NOACB1307	☐ \$660	13"-15"	7"			NOPCB1307	☐ \$660
	13"-15"	9"			NOACB1309	☐ \$660	13"-15"	9"			NOPCB1309	☐ \$660
	13"-15"	13"			NOACB1313	☐ \$660	13"-15"	13"			NOPCB1313	☐ \$660
	13"-15"	16"			NOACB1316	☐ \$660	13"-15"	16"			NOPCB1316	☐ \$660
	14"-16"	7"			NOACB1407	☐ \$660	14"-16"	7"			NOPCB1407	☐ \$660
	14"-16"	9"	NOCFACB1409	☐ \$870	NOACB1409	☐ \$660	14"-16"	9"	NOCFPCB1409	☐ \$965	NOPCB1409	☐ \$660
	14"-16"	13"	NOCFACB1413	☐ \$900	NOACB1413	☐ \$660	14"-16"	13"	NOCFPCB1413	☐ \$985	NOPCB1413	☐ \$660
	14"-16"	16"	NOCFACB1416	☐ \$930	NOACB1416	☐ \$660	14"-16"	16"	NOCFPCB1416	☐ \$1025	NOPCB1416	☐ \$660
	15"-17"	7"			NOACB1507	☐ \$660	15"-17"	7"			NOPCB1507	☐ \$660
	15"-17"	9"	NOCFACB1509	☐ \$870	NOACB1509	☐ \$660	15"-17"	9"	NOCFPCB1509	☐ \$965	NOPCB1509	☐ \$660
	15"-17"	13"	NOCFACB1513	☐ \$900	NOACB1513	☐ \$660	15"-17"	13"	NOCFPCB1513	☐ \$985	NOPCB1513	☐ \$660
	15"-17"	16"	NOCFACB1516	☐ \$930	NOACB1516	☐ \$660	15"-17"	16"	NOCFPCB1516	☐ \$1025	NOPCB1516	☐ \$660
	16"-18"	7"			NOACB1607	☐ \$660	16"-18"	7"			NOPCB1607	☐ \$660
	16"-18"	9"	NOCFACB1609	☐ \$870	NOACB1609	☐ \$660	16"-18"	9"	NOCFPCB1609	☐ \$965	NOPCB1609	☐ \$660
	16"-18"	13"	NOCFACB1613	☐ \$900	NOACB1613	☐ \$660	16"-18"	13"	NOCFPCB1613	☐ \$985	NOPCB1613	☐ \$660
	16"-18"	16"	NOCFACB1616	☐ \$930	NOACB1616	☐ \$660	16"-18"	16"	NOCFPCB1616	☐ \$1025	NOPCB1616	☐ \$660
	17"-19"	7"			NOACB1707	☐ \$660	17"-19"	7"			NOPCB1707	☐ \$660
	17"-19"	9"	NOCFACB1709	☐ \$870	NOACB1709	☐ \$660	17"-19"	9"	NOCFPCB1709	☐ \$965	NOPCB1709	☐ \$660
	17"-19"	13"	NOCFACB1713	☐ \$900	NOACB1713	☐ \$660	17"-19"	13"	NOCFPCB1713	☐ \$985	NOPCB1713	☐ \$660
	17"-19"	16"	NOCFACB1716	☐ \$930	NOACB1716	☐ \$660	17"-19"	16"	NOCFPCB1716	☐ \$1025	NOPCB1716	☐ \$660
	18"-20"	7"			NOACB1807	☐ \$660	18"-20"	7"			NOPCB1807	☐ \$660
	18"-20"	9"	NOCFACB1809	☐ \$870	NOACB1809	☐ \$660	18"-20"	9"	NOCFPCB1809	☐ \$965	NOPCB1809	☐ \$660
	18"-20"	13"	NOCFACB1813	☐ \$900	NOACB1813	☐ \$660	18"-20"	13"	NOCFPCB1813	☐ \$985	NOPCB1813	☐ \$660
	18"-20"	16"	NOCFACB1816	☐ \$930	NOACB1816	☐ \$660	18"-20"	16"	NOCFPCB1816	☐ \$1025	NOPCB1816	☐ \$660
	20"-22"	7"			NOACB2007	☐ \$660	20"-22"	7"			NOPCB2007	☐ \$660
	20"-22"	9"	NOCFACB2009	☐ \$870	NOACB2009	☐ \$660	20"-22"	9"	NOCFPCB2009	☐ \$965	NOPCB2009	☐ \$660
	20"-22"	13"	NOCFACB2013	☐ \$900	NOACB2013	☐ \$660	20"-22"	13"	NOCFPCB2013	☐ \$985	NOPCB2013	☐ \$660
	20"-22"	16"	NOCFACB2016	☐ \$930	NOACB2016	☐ \$660	20"-22"	16"	NOCFPCB2016	☐ \$1025	NOPCB2016	☐ \$660

BACK SUPPORT HARDWARE OPTIONS

All backs come STANDARD with QuickRelease Mounting Hardware. The following are fixed hardware options available at no charge when ordered with a back support.

Model #	DESCRIPTION	WEIGHT OF HARDWARE ONLY	✓
	QuickFit Quick Release Hardware System	STANDARD (2 LBS)	STANDARD
NIB-1642	Micro Mount 2 point fixed hardware. Back width needs to match chair width	Weight savings of 1.4 lb over STANDARD	☐
NIB-1207	2 point fixed hardware. Width expandable to fit 3 sizes of chair	Weight savings of 1.0 lb over STANDARD	☐

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com