



Two views of a black and red wheelchair. The left view shows the wheelchair from a side profile, highlighting the large rear wheels with red rims and the smaller front casters. The right view shows the wheelchair from a slightly different angle, emphasizing the red frame and the black seat and backrest.

1- Purchasing information

Marked for:

2- ATP/Clinician information

Contact preference:* Phone ☐ Email ☐

3- Shipping information


Phone: _____

Taxes not included

2- Not available with dept adjustable back; **3-** May interfere with the swing away armrest and the clamp on height adjustable push handle; Add a 3" of frame length behind backrest


CASTERS						
CA01-03	☐	Composite wheel w/ PU tire	3" x 1"			NC
CA01-04	☐		4" x 1"			NC
CA01-05	☐		5" x 1"			NC
CA01-06	☐		6" x 1"			NC
CA10	☐	Pneumatic casters	6" x 1 ¼"			\$165
CA14	☐	Newton UltraCaster Composite	4" x 1 ½ "			\$87
CA15	☐		5" x 1 ½ "			\$87
CA16	☐		6" x 1 ½ "			\$87
CA17	☐	Newton UltraCaster Aluminum	4" x 1 ½ "			\$190
CA18	☐		5" x 1 ½ "			\$200
CA19	☐		6" x 1 ½ "			\$210
CASTER OPTIONS						
CA13	☐	Carbon fiber Frog Legs suspension forks			E1015	\$485

REAR Seat-to-floor height		
Wheel sizes	Height available +/- 1/4"	Specify desired height: _____ (inches)
20" (451) wheel ☐	14" to 17.5"	
22" (501) wheel ☐	14" to 18.5"	
24" (540) wheel ☐	15" to 19.5"	
25" (559) wheel ☐	15.5" to 20"	
26" (590) wheel ☐	16" to 20"	



FRONT FRAME ANGLE				
FA70	☐	70° Front frame angle		NC
FA75	☐	75° Front frame angle		NC
FA80	☐	80° Front frame angle		NC
FA85	☐	85° Front frame angle		NC
FA90	☐	90° Front frame angle		NC

Seat upholstery to footplate length	
12" to 19" without high mount	Specify desired length: _____ (inches) <i>Required</i>
7" to 13" with high mount ¹	



1- Choose high mount in the footplates section


Footrest width / Taper selection				
Seat Width	Straight	1" Taper	2" Taper	
Outside Measurement	2 ½" Narrower than Seat Width	4 ½" Narrower than Seat Width	6 ½" Narrower than Seat Width	
12"	9.5"	7.5"		
13"	10.5"	8.5"		
14"	11.5"	9.5"		
15"	12.5"	10.5"		
16"	13.5"	11.5"	9.5"	
17"	14.5"	12.5"	10.5"	
18"	15.5"	13.5"	11.5"	
19"	16.5"	14.5"	12.5"	
20"	17.5"	15.5"	13.5"	

Note: higher than 20" Front STFH will reduce the footrest width up to 1" with a 2" taper

FOOTPLATES						
FP12	☐	Tubular open loop ¹				STD
FP13	☐	Tubular with ABS cover			K0040	NC
FP31	☐	Adjustable angle - Aluminum <i>(one piece)</i>			K0040	\$97
FP32	☐	Newton adjustable angle - Carbon fiber <i>(one piece with rubber inserts)</i>			K0040	\$190
FP14	☐	Flip back, adjustable angle aluminum			K0040	\$200
FP15	☐	Flip back, adjustable angle carbon fiber				\$295
FR25	☐	High mount with adjustable angle footplate <i>(Aluminum - one piece)</i>			K0040	\$170
FR26	☐	High mount with adjustable angle footplate <i>(Carbon fiber with rubber inserts - one piece)</i>				\$250
FOOTPLATES ACCESSORIES						
FP16	☐	2" riser <i>(only available with aluminum angle adjustable footplate)</i> ²				\$77
FP09	☐	Calf strap			K0038	\$36
FP30	☐	Bodypoint Aeromesh Padded calf strap			K0038	\$118

1 - Available on: all straight frames; 1" taper with seat width of 14" and greater; 2" taper with seat width of 16" and greater

2 - Not available with high mount


Backrest														
FOLD DOWN BACKREST														
BC20	☐	Locking backrest											STD	
BC21	☐	Ultralight non locking backrest (not locking in upright and down position)											\$102	
FIXED HEIGHT, ADJUSTABLE ANGLE FOLD DOWN BACKREST (Not available with push handles)														
BC	☐	Fixed height without push handle		(30) 9" ☐	(31) 9.5" ☐	(32) 10" ☐	(33) 10.5" ☐	(34) 11" ☐	NC					
		(35) 11.5" ☐	(36) 12" ☐	(37) 12.5" ☐	(38) 13" ☐	(39) 13.5" ☐	(40) 14" ☐	(41) 14.5" ☐						
		(42) 15" ☐	(43) 15.5" ☐	(44) 16" ☐	(45) 16.5" ☐	(46) 17" ☐	(47) 17.5" ☐	(48) 18" ☐						
		(49) 18.5" ☐	(50) 19" ☐	(51) 19.5" ☐	(52) 20" ☐	(53) 20.5" ☐	(54) 21" ☐							
HEIGHT AND ANGLE ADJUSTABLE FOLD DOWN BACKREST WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)														
BC	☐	Straight back canes		(22) 9 - 12" ☐ (23) 12 - 15" ☐ (24) 15 - 18" ☐ (25) 18 - 21" ☐								NC		
Set height to: *			(in increments of 1/2", measured from seat sling to top of back upholstery) * If not indicated, default height setting is at the lowest position											
Select back cane angle: ¹		101° ☐	99° ☐	96° ☐	93° ☐	90° ☐	87° ☐	83° ☐	80° ☐					STD

1 - If not indicated default angle setting is at 93°

PUSH HANDLES OPTIONS <i>(ONE option at a time)</i>						
PH01	☐	Push handles				STD
PH05	☐	Fold down push handles ¹				\$130
PH06-920	☐	Clamp on height adjustable push handles 9 - 20" ^{1,2}				\$320
PH06-923	☐	Clamp on height adjustable push handles 9 - 23" ^{1,2}				\$320
PH02	☐	Omit push handles				NC

1- 265 pounds weight limit

2 - Not compatible with swing away armrests;

BACKRESTS OPTIONS					
BC61	☐	Depth adjustable back - 1" depth adjustment ^{1,4}			\$195
BC62	☐	Depth adjustable back - 2" depth adjustment ^{2,4}			\$195
BC56	☐	Depth adjustable back - 3" depth adjustment ^{3,4}			\$195
BC57	☐	Standard 3 inch deep rigidizer bar			STD
BC55	☐	Four inch deep Aluminum rigidizer bar			NC
BC88	☐	Custom height rigidizer bar set at _____ (inches) Measured from the top of the seat sling to the center of the rigidizer - standard height is 5,7" Lowest position is 4,25"			\$135

If used in combination with WC-19: 1- Maximum depth of 19"; 2- Maximum depth of 18"; 3- Maximum depth of 17"; 4- Only compatible with WC-19 transit tie down **without** Q'Straint belt

BACK UPHOLSTERY						
BH03	☐	Slip-on Syntech-Suede back upholstery				STD
BH04	☐	Adjustable tension Syntech-Suede back upholstery			E2611	\$305
BH99	☐	Omit Back upholstery				NC
	☐	NXT back Support			See price list	
SEAT CUSHION						
SC01	☐	2" Cushion			E2601	\$170
SC02	☐	3" Cushion			E2601	\$190
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See price list	

Rear wheels						
REAR WHEELS				Available		
RW01	☐	Mag ¹		20", 22" 24"		\$57
RW02	☐	Newton One - spoke wheel ²		20", 22" 24"		STD
RW03	☐	Newton Gravity - ultralight wheel		20", 22", 24", 25", 26"		\$320
RW04	☐	Spinergy Spox Blk STD ☐ Wht ³ ☐ Blue ³ ☐ Red ³ ☐ Yel ³ ☐		22", 24", 25", 26"		\$850
RW05	☐	Spinergy LX Blk STD ☐ Pur ³ ☐ Wht ³ ☐ Blue ³ ☐ Red ³ ☐ Yel ³ ☐ Pink ³ ☐ Grn ³ ☐ Org ³ ☐		22", 24", 25", 26"		\$1 065
RW07	☐	Spinergy CLX ⁴ Blk STD ☐ Pur ³ ☐ Wht ³ ☐ Blue ³ ☐ Red ³ ☐ Yel ³ ☐ Pink ³ ☐ Grn ³ ☐ Org ³ ☐		24", 25"		\$1 585
TIRES				Available		
TI01	☐	Soft Urethane 1 3/8 "		20", 22" 24"	E2218	NC
TI02	☐	Soft Urethane knobby 1 3/8"		22", 24", 25", 26"	E2218	\$97
TI04	☐	Soft Urethane - Shox 1"		22", 24", 25", 26"	E2218	\$118
TI05	☐	Pneumatic 1 3/8 "		20", 22", 24", 25", 26"	E2211/E2212	STD
TI08	☐	Pneumatic with Airless insert 1 3/8" ⁵		22", 24", 26"	E2211/E2213	\$82
TI06	☐	Pneumatic HP - 1"		20", 22", 24", 25", 26"	E2211/E2212	\$97
TI07	☐	Pneumatic HP - Puncture resistant, black 1"		24", 25", 26"	E2211/E2212	\$130

1- Not compatible with HP tires; Only available with aluminum or plastic coated handrim in standard position; 2- Not compatible with Shox tires
3- color spokes require an additional 10 days lead-time; 4- Not available with Soft Urethane Knobby; 5- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims						
HANDRIM				Available		
HR01	☐	Aluminum anodized		20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²		20", 22", 24", 25"		\$72
HR02	☐	Plastic coated ¹		20", 22", 24", 25", 26"		\$102
HR03	☐	High friction coated ^{1,2}		22", 24", 25", 26"		\$375
HR07	☐	Newton Air Grip ² <i>Black</i> ☐		20", 22", 24", 25", 26"		\$102
HR	☐	Newton Air Grip ² (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐ (11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐ (15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐		20", 22", 24", 25", 26"		\$255
HR28	☐	Surge ^{1,2} (1 1/2" Oval handrim w/ Gripton high friction strip) **Limited availability		22", 24", 25", 26"		\$395
HR05	☐	Surge LT ^{1,2} (1 1/8" Oval handrim w/ Gripton high friction strip) **Limited availability		22", 24", 25", 26"		\$425
HR04	☐	Natural Fit ^{1,2} **Limited availability		22", 24", 25", 26"		\$325
HR62	☐	Tour H ^{1,2} (1 5/16" Oval handrim with high friction strip)		22", 24", 25", 26"	New	\$395
HR63	☐	Nova H ^{1,2} (1 1/4" Oval handrim with high friction strip)		22", 24", 25", 26"	New	\$425
HR64	☐	Optimum H ^{1,2} (1 1/8" heart-shaped handrim)		22", 24", 25", 26"	New	\$325
HR99	☐	Omit handrim				NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels

Handrims option						
ASSEMBLY POSITION						
HR29	☐	Standard Handrim position <i>Full tab</i>				STD
HR030	☐	Narrow Handrim position no cut <i>Full tab</i>				NC
HR30	☐	Narrow Handrim position with cut ³ <i>Short tab 2 positions</i>				\$30
HR31	☐	Super Narrow Handrim ³ <i>Short tab 2 positions STD</i> ☐ <i>Short tab 1 position</i> ☐				\$30
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip **Limited availability				\$92
HR27	☐	Super grip **Limited availability				\$150
HR099	☐	Omit thumb piece				NC

3- Only available with Aluminum and Air Grip handrims;

AXLES					
AX01	☐	Quick release axle			STD
AX03	☐	Quad release axle		K0108	\$130
REAR WHEEL CAMBER					
CM00	☐	0° camber			STD
CM02	☐	2° camber			NC
CM04	☐	4° camber			NC
CM06	☐	6° camber			NC
CM08	☐	8° camber			NC

CENTER OF GRAVITY


Specify the position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

0 - 6" range	Set At: _____ (inches)
--------------	------------------------

If not indicated, default setting is at 2.5"

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG



WHEEL LOCK					
WL08	☐	Newton wheel lock - push to lock			STD
WL04	☐	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)		E0974	\$155
WL13 WL14	☐	6" removable extension handle		E0961	\$77
WL05	☐	Scissors lock			\$155
WL99	☐	Omit wheel locks			NC

Armrests

HEIGHT ADJUSTABLE "T" ARMREST ¹					
AR13-610	☐	10" Desk length armpad ³	6" - 10"		\$330
AR03-812	☐	10" Desk length armpad	8" - 12"		\$330
AR03-1014	☐	10" Desk length armpad	10" - 14"		\$330
AR14-610	☐	14" Full length armpad ³	6" - 10"		\$330
AR04-812	☐	14" Full length armpad	8" - 12"		\$330
AR04-1014	☐	14" Full length armpad	10" - 14"		\$330
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²					
AR05	☐	Tubular armrest without plastic side guards - Short pad	8" - 11"	E0973	\$230
AR19	☐	Tubular armrest with plastic side guards - short pad	8" - 11"	E0973	\$395
AR23	☐	Tubular armrest without plastic side guards - long pad	8" - 11"	new E0973	\$230
AR24	☐	Tubular armrest with plastic side guards - long pad	8" - 11"	new E0973	\$395
AR06	☐	10" Desk length flat armpad without plastic side guards	8" - 11"	E0973	\$230
AR20	☐	10" Desk length flat armpad with plastic side guards	8" - 11"	E0973	\$395
AR17	☐	14" Full length flat armpad without plastic side guards	8" - 11"	E0973	\$230
AR21	☐	14" Full length flat armpad with plastic side guards	8" - 11"	E0973	\$395
AR99	☐	Omit armrest			NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- not available with the Ultralight non-locking backrest; 3- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with "T" Armrest)					
SG01	☐	Plastic side guards (removable)		K0108	\$200
SG03	☐	Carbon fiber side guards (removable)		K0108	\$295
SG06	☐	Plastic fender side guards - straight bracket (removable)			\$290
SG08	☐	Plastic fender side guards - offset bracket (removable) ¹			\$290
SG04	☐	Carbon fiber fender side guards (removable)		K0108	\$550
SG99	☐	Omit side guards			NC
SIDE GUARDS OPTION					
SG07	☐	Side guard clip (pair)			NC

1- only compatible with 0° and 2° camber; Increase distance between fender side guards by 1 1/8"

POSITIONING BELT					
SB01-S	☐	Velcro adjustable - Standard 35"		E0978	\$46
SB01-L	☐	Velcro adjustable - Long 60"		E0978	\$46
SB02-S	☐	Auto buckle - Standard 35"		E0978	\$87
SB02-L	☐	Auto buckle - Long 60"		E0978	\$87
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"		E0978	\$82
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"		E0978	\$82
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"		E0978	\$175
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"		E0978	\$175
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"		E0978	\$175
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"		E0978	\$195
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"		E0978	\$195
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"		E0978	\$195
SB03	☐	Hardware attachment for 2 extra points (1 pair)			\$17
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"		K0108	\$305
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"		K0108	\$305
SB99	☐	Omit positioning belt			STD

ACCESSORIES					
AC03	☐	Anti-tippers		E0971	\$145
AC06	☐	Cane and crutch holder	Left ☐ Right ☐	E2207	\$200
AC08	☐	Spoke guards ¹	Available: 20", 22", 24", 26"	K0065	\$118
AC10	☐	Seat pouch			\$82
AC30	☐	Newton multi-purpose pouch			\$61
AC26	☐	BodyPoint Mobility Bag			\$135
AC23	☐	Newton Backpack			\$130
AC13	☐	Tool kit			\$51
AC16	☐	Motion Composites Polo Shirt	W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐		\$36
AC15	☐	Military Patches	Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐ Combat Wounded Purple Heart ☐		\$97

1- Not available with Mag wheels and Spinergy

HIGHLIGHT COLOR					
FC	☐	(10) Ferrari Red ☐ (11) White ☐ (12) Black ☐ (13) Sunkissed Orange ☐ (14) Sky Blue ☐ (15) Fuchsia ☐ (16) Acid Green ☐			NC

Notes :

BioFit™ CUSHION		WIDTH									
DEPTH		12"	13"	14"	15"	16"	17"	18"	19"	20"	
	12"	☐ \$445	☐ \$445	☐ \$445							
	14"	☐ \$445	☐ \$445	☐ \$445							
	15"		☐ \$445	☐ \$445	☐ \$445						
	16"		☐ \$445	☐ \$445	☐ \$445	☐ \$445	☐ \$445	☐ \$445	☐ \$480	☐ \$480	
	17"				☐ \$445	☐ \$445	☐ \$445	☐ \$445	☐ \$480	☐ \$480	
	18"				☐ \$445	☐ \$445	☐ \$445	☐ \$445	☐ \$480	☐ \$480	
	19"				☐ \$445	☐ \$445	☐ \$445	☐ \$445	☐ \$480	☐ \$480	
	20"				☐ \$445	☐ \$445	☐ \$445	☐ \$480	☐ \$480	☐ \$480	
Option: Change cover for incontinent/infection control outer cover ☐ \$N/C											

<

SEATING AND CUSHION ACCESSORIES

	MODEL NO.	Description	Fits Back Width	\$	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6 & 43.2 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7 & 48.3 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	80	

THORACIC SHAPE	Optima™ THORACIC BACK SUPPORT						Optima™ DEEP THORACIC BACK SUPPORT					
	Fits Width	Height	Model # CARBON FIBER	\$ CARBON	Model # Aluminum	\$ Aluminum	Fits Width	Height	Model # CARBON FIBER	\$ CARBON	Model # Aluminum	\$ Aluminum
	13"-15"	7"			NOACB1307	☐ \$695	13"-15"	7"			NOPCB1307	☐ \$695
	13"-15"	9"			NOACB1309	☐ \$695	13"-15"	9"			NOPCB1309	☐ \$695
	13"-15"	13"			NOACB1313	☐ \$695	13"-15"	13"			NOPCB1313	☐ \$695
	13"-15"	16"			NOACB1316	☐ \$695	13"-15"	16"			NOPCB1316	☐ \$695
	14"-16"	7"			NOACB1407	☐ \$695	14"-16"	7"			NOPCB1407	☐ \$695
	14"-16"	9"	NOCFACB1409	☐ \$895	NOACB1409	☐ \$695	14"-16"	9"	NOCFPCB1409	☐ \$995	NOPCB1409	☐ \$695
	14"-16"	13"	NOCFACB1413	☐ \$925	NOACB1413	☐ \$695	14"-16"	13"	NOCFPCB1413	☐ \$1015	NOPCB1413	☐ \$695
	14"-16"	16"	NOCFACB1416	☐ \$960	NOACB1416	☐ \$695	14"-16"	16"	NOCFPCB1416	☐ \$1055	NOPCB1416	☐ \$695
	15"-17"	7"			NOACB1507	☐ \$695	15"-17"	7"			NOPCB1507	☐ \$695
	15"-17"	9"	NOCFACB1509	☐ \$895	NOACB1509	☐ \$695	15"-17"	9"	NOCFPCB1509	☐ \$995	NOPCB1509	☐ \$695
	15"-17"	13"	NOCFACB1513	☐ \$925	NOACB1513	☐ \$695	15"-17"	13"	NOCFPCB1513	☐ \$1015	NOPCB1513	☐ \$695
	15"-17"	16"	NOCFACB1516	☐ \$960	NOACB1516	☐ \$695	15"-17"	16"	NOCFPCB1516	☐ \$1055	NOPCB1516	☐ \$695
	16"-18"	7"			NOACB1607	☐ \$695	16"-18"	7"			NOPCB1607	☐ \$695
	16"-18"	9"	NOCFACB1609	☐ \$895	NOACB1609	☐ \$695	16"-18"	9"	NOCFPCB1609	☐ \$995	NOPCB1609	☐ \$695
	16"-18"	13"	NOCFACB1613	☐ \$925	NOACB1613	☐ \$695	16"-18"	13"	NOCFPCB1613	☐ \$1015	NOPCB1613	☐ \$695
	16"-18"	16"	NOCFACB1616	☐ \$960	NOACB1616	☐ \$695	16"-18"	16"	NOCFPCB1616	☐ \$1055	NOPCB1616	☐ \$695
	17"-19"	7"			NOACB1707	☐ \$695	17"-19"	7"			NOPCB1707	☐ \$695
	17"-19"	9"	NOCFACB1709	☐ \$895	NOACB1709	☐ \$695	17"-19"	9"	NOCFPCB1709	☐ \$995	NOPCB1709	☐ \$695
	17"-19"	13"	NOCFACB1713	☐ \$925	NOACB1713	☐ \$695	17"-19"	13"	NOCFPCB1713	☐ \$1015	NOPCB1713	☐ \$695
	17"-19"	16"	NOCFACB1716	☐ \$960	NOACB1716	☐ \$695	17"-19"	16"	NOCFPCB1716	☐ \$1055	NOPCB1716	☐ \$695
	18"-20"	7"			NOACB1807	☐ \$695	18"-20"	7"			NOPCB1807	☐ \$695
	18"-20"	9"	NOCFACB1809	☐ \$895	NOACB1809	☐ \$695	18"-20"	9"	NOCFPCB1809	☐ \$995	NOPCB1809	☐ \$695
	18"-20"	13"	NOCFACB1813	☐ \$925	NOACB1813	☐ \$695	18"-20"	13"	NOCFPCB1813	☐ \$1015	NOPCB1813	☐ \$695
	18"-20"	16"	NOCFACB1816	☐ \$960	NOACB1816	☐ \$695	18"-20"	16"	NOCFPCB1816	☐ \$1055	NOPCB1816	☐ \$695
	20"-22"	7"			NOACB2007	☐ \$695	20"-22"	7"			NOPCB2007	☐ \$695
	20"-22"	9"	NOCFACB2009	☐ \$895	NOACB2009	☐ \$695	20"-22"	9"	NOCFPCB2009	☐ \$995	NOPCB2009	☐ \$695
	20"-22"	13"	NOCFACB2013	☐ \$925	NOACB2013	☐ \$695	20"-22"	13"	NOCFPCB2013	☐ \$1015	NOPCB2013	☐ \$695
	20"-22"	16"	NOCFACB2016	☐ \$960	NOACB2016	☐ \$695	20"-22"	16"	NOCFPCB2016	☐ \$1055	NOPCB2016	☐ \$695

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

 For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

BACK SUPPORT HARDWARE OPTIONS

All backs come STANDARD with QuickRelease Mounting Hardware. The following are fixed hardware options available at no charge when ordered with a back support.

Model #	DESCRIPTION	WEIGHT OF HARDWARE ONLY	✓
	QuickFit Quick Release Hardware System	STANDARD (2 LBS)	STANDARD
NIB-1642	Micro Mount 2 point fixed hardware. Back width needs to match chair width	Weight savings of 1.4 lb over STANDARD	☐
NIB-1207	2 point fixed hardware. Width expandable to fit 3 sizes of chair	Weight savings of 1.0 lb over STANDARD	☐

BACK SUPPORT ACCESSORIES – PRIVACY FLAP

	Model #	Description	Fits Back Width	\$	QTY
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	70	