







MOTION
COMPOSITES

EVO PROGRAM



EVO
CUSTOM

EVO PROGRAM

The **EVO Program** gives wheelchair users peace of mind, knowing that their unique needs won't leave them with an ill-fitting chair, whether that means a frame replacement or customizations.



THE OPTIMAL SOLUTION TO ENSURE EVERY COMPONENT OF YOUR WHEELCHAIR IS TAILORED TO YOUR UNIQUE NEEDS.



A PRECISELY ADJUSTED REPLACEMENT BASE FRAME DESIGNED TO ACCOMMODATE YOUR CHANGING NEEDS.

CUSTOMIZATION STARTS HERE!

Our wheelchairs are designed for adjustability and meet a wide range of needs. **EVO Custom**, a part of the **EVO Program**, takes it even further! We know that the right customization solution can be life changing for wheelchair users. From subtle tweaks for comfort and mobility improvements to fully reimagined components tailored to your unique requirements, we're here to support you.

EXPRESS

OFF THE SHELF

- Recurrent customization requests that are kept in stock.
- Same quoting process as a regular order through customer service.
- Same lead times as a regular order.

TAILORED

MADE TO ORDER

- Recurrent customization requests that require manufacturing, but have already been designed.
- Same quoting process as a regular order through customer service.
- Minimal added lead times for manufacturing.

SIGNATURE

GO BEYOND

- Customization requests requiring a new design, a technical assessment, or further research.
- Quotes provided within a few extra days.
- Lead time vary depending on complexity.

HOW TO ORDER

GETTING STARTED WITH EVO CUSTOM IS SIMPLE.



SCAN HERE!

For **EXPRESS** and **TAILORED** customization levels, you can include your request or existing custom solution from our catalog in the notes section of your standard wheelchair order form. Prefer using a dedicated form? Our [EVO Custom Request Form](#) is also available for all customization levels and helps ensure nothing gets missed.

For **SIGNATURE** requests, or if you are not sure where your idea lands, the [EVO Custom Request Form](#) is strongly recommended to help us assess your needs and provide accurate timeline and pricing.

EVO CUSTOM



CONTENT

FRONT WHEELS	04
FOOT SUPPORTS	05
BACKRESTS	09
UPHOLSTERY & SEATING	12
REAR WHEELS	14
WHEEL LOCKS	16
ARMRESTS	18
SIDE GUARDS	20
ACCESSORIES	21
SIGNATURE CUSTOMS	23
PRICE LIST	30

FRONT WHEELS

CASTER LOCK PIN **EXPRESS**

0005330A-XX

APEX C | APEX A | VELOCE

Mechanism to prevent caster wheels from pivoting to help maintain direction. Can be activated or deactivated.

- Straight stem bolt
- 1" taper stem bolt
- 2" taper stem bolt

Approx. weight: 1.04 - 1.16 lbs (0.475 - 0.525 kg)

Available for standard length stem bolt only. Not compatible with Frog Legs.

- Facilitates unidirectional movement (straight trajectory).
- Prevents wheels from swiveling.
- Allows for safer transfers by preventing the front casters from swiveling.



FOOT SUPPORTS

CARBON PLATFORM FOOTPLATE

TAILORED

0005399A-XX

APEX C | APEX A | VELOCE

Custom sized carbon footplates that are not available on the order form.

Full-width plate only. Customer must specify desired width and depth.

- Accommodates lower extremity asymmetries.
- Accommodates range of motion limitations in the foot/ankle complex.
- Provides durable foot platform while being very lightweight.



ALUMINUM PLATFORM FOOTPLATE

TAILORED

0005400A-XX

APEX C | APEX A | VELOCE

Custom sized aluminum footplates that are not available on the order form.

Full-width plate only. Customer must specify desired width and depth.

- Accommodates lower extremity asymmetries.
- Accommodates range of motion limitations in the foot/ankle complex.



WINGED PLATFORM FOOTPLATE EXPRESS

0014753A-00

APEX C | APEX A

Custom 3D printed footplate with raised sides installed on an open loop structure.

Approx. weight: 1.04 lbs (0.47 kg)

Available only on the open loop.

- Facilitates optimal foot alignment.
- Helps prevent the lower extremities from falling laterally off the footplate.



ELEVATING LEG RESTS EXTENSION TUBE EXPRESS

0013514A-00

HELIO C2 | HELIO A6 | HELIO A7

Custom elongated tube for the elevated leg rests. Adds 2" to the longest extension available on order form.

Adds maximum 2" of length to elevating leg rests.

- Accommodates longer leg lengths.

SWING AWAY RESIDUAL LIMB SUPPORT EXPRESS

0006252A-XX

APEX C | APEX A | VELOCE

Sturdy and easily removable swing away residual limb support. Available for either left or right, in two shapes (flat or contoured), and in several different sizes.

- Swing out residual limb support - Flat 6x8"
- Swing out residual limb support - Flat 8x10"
- Swing out residual limb support - Flat 10x10"
- Swing out residual limb support - Flat 10x14"
- Swing out residual limb support - Contoured 6x8"
- Swing out residual limb support - Contoured 8x10"

Approx. weight: 2.79 lbs (1.26 kg)

Full-width plate only. Additional cushion sizes may be available as a Tailored or Signature.

- Provides support to the residual limb.
- Facilitates optimal alignment of the residual limb.
- Swing away feature helps facilitate transfers and ADLs.



HIGH MOUNT FLIP BACK FOOTPLATE TAILORED

0011380A-00

APEX C | APEX A

Full flip back footrest assembly installed as a high mount on APEX chairs. Easily adjustable to various heights.

Approx. weight: 1.04 lbs (0.47 kg)

Custom includes structure only. Footplate is delivered with the chair. Minimum chair width: 12".

- Accommodates a range of leg lengths at a desired height.
- Flip back feature facilitates transfers and ADLs.



APEX FOOTPLATE ASSEMBLY WITH LONG CLAMP

TAILORED

0012756A-00

APEX C | APEX A

Custom footrest clamp on APEX chairs that allows for height adjustability by splitting the tube in 2 and using a longer clamp to unite both ends. This allows freedom of movement going up and down even with a taper in the footrest.

Angle adjustable footrest only. Clamp and structure only (no footplate).

- To maximize footrest height adjustability



BACKRESTS

HEIGHT ADJUSTABLE SINGLE HANDLE **EXPRESS**

0009313C-00

APEX C | APEX A

Height adjustable and removable back cane installed on rear aluminum rigidizer. Single push handle mount on APEX rigidizer.

For APEX C, the carbon rigidizer needs to be replaced with an aluminum version.

- Facilitates caregiver assistance for dependent mobility.
- Handle can be quickly adjusted to the desired heights.



REINFORCED DEPTH ADJUSTABLE BACKREST STRUCTURE **TAILORED**

0013276A-XX

HELIO kids

Can help prevent breakage under demanding use.

- 9-12"
- 12-15"
- 15-18"
- 18-21"
- 8° - 9-12"
- 8° - 12-15"
- 8° - 15-18"
- 8° - 18-21"

Depth adjustable back only. Weight capacity 175 lb.

- Depth adjustability facilitates optimal postural alignment.
- The reinforced material helps accommodate abnormal tone and reduces risk of premature back cane failure.



FOLDING BACK CANES - FIXED ANGLE EXPRESS

MULTIPLE

VELDCE | HELIO C2 | HELIO A7 | HELIO A6
MOVE | HELIO kids

Allows back canes to fold down on folding chairs. Compatible with any fixed angle back cane configuration.

- Straight aluminum canes 9-12"
- Straight aluminum canes 12-15"
- Straight aluminum canes 15-18"
- Straight aluminum canes 18-21"
- 8° aluminum canes 13.5-15"
- 8° aluminum canes 15-18"
- 8° aluminum canes 18-21"
- Straight carbon canes 9-12"
- Straight carbon canes 12-15"
- Straight carbon canes 15-18"
- Straight carbon canes 18-21"
- 8° carbon canes 13.5-15"
- 8° carbon canes 15-18"
- 8° carbon canes 18-21"
- Integrated push handles aluminum canes 9-15"
- Integrated push handles aluminum canes 15-21"
- Integrated push handles aluminum canes 21-24"
- 8° integrated push handles aluminum canes 17-20"
- 8° integrated push handles aluminum canes 21-24"

Approx. weight: 1.10 - 1.70 lbs (0.50 - 0.80 kg)

Not compatible with ISO7176-19 or WC-19. Not compatible with flip back armrests.

- Facilitates ease of storage and transportation.



FOLDING BACKCANES - ADJUSTABLE ANGLE

EXPRESS

MULTIPLE

VELOCE | HELIO C2 HELIO A7 | HELIO A6
MOVE | HELIO kids

Allows back canes to fold down on folding chairs. Compatible with any adjustable angle back cane configuration.

- Aluminum canes 15-18"
- Aluminum canes 18-21"
- 8° aluminum canes 15-18"
- 8° aluminum canes 18-21"
- Carbon canes 14-15"
- Carbon canes 15-18"
- Carbon canes 18-21"
- 8° carbon canes 14-15"
- 8° carbon canes 15-18"
- 8° carbon canes 18-21"

Approx. weight: 1.10 - 1.70 lbs (0.50 - 0.80 kg)

Not compatible with ISO7176-19 or WC-19. Not compatible with flip back armrests.

- Facilitates ease of storage and transportation.
- Angle adjustable back canes facilitate optimal positional adjustments.



UPHOLSTERY & SEATING

ADDUCTOR TAILORED

0015196A-XX

**APEX C | APEX A | VELOCE | HELIO C2
HELIO A7 | HELIO A6 | MOVE | HELIO kids**

Custom external laterals that can be configured in any shape, length, height or thickness. Allows seating area to be reduced.

It is mandatory to specify the size and thickness of each lateral pad..

- Facilitates neutral alignment of the lower trunk and pelvis.
- Prevents excessive lateral movement of the lower trunk and pelvis.



KNEE ADDUCTOR TAILORED

0006009A-XX

**HELIO C2 | HELIO A7 | HELIO A6
MOVE | HELIO kids**

Knee adductor available left and right installed on chair to hold user's knees in place. Can be adjusted in height by sliding in the holder and is easily removable.

- Left side
- Right side

It is mandatory to specify the size and thickness of each pad.

- Facilitates neutral alignment of the knees.
- Prevents excessive abduction of the distal thighs.





CARBON FIBER SOLID SEAT PAN TAILORED

0010294A-XX

APEX C | APEX A

Carbon fiber solid seat pan to replace fabric seat slings, can fit any chair size, includes installation hardware.

- Provides durable and stable base of support for the seat cushion.
- Prevents hammocking and helps maintain optimal positioning.

NXT BACK SUPPORT VELCRO PATCH EXPRESS

MULTIPLE

AVAILABLE ON ALL NXT BACK SUPPORTS

USA Military service patches for any solid NXT backrest. Can be preinstalled, or sent with adhesive velcro for easy self installation.

- Air force patch with velcros
- Army patch with velcros
- Navy patch with velcros
- Marine corps patch with velcros
- Purple heart patch with velcros
- Coast guard patch with velcros
- Space force patch with velcros
- USA flag patch with velcros

Available for installation internally when ordering a NXT back support, or to retrofit on an existing back support.



REAR WHEELS

ONE ARM DRIVE KIT FOR 22" WHEEL

TAILORED

0005993A-XX

HELIO C2 | HELIO A7 | HELIO A6 | MOVE

One arm drive assembly for 22" rear wheels, comes complete with anodized or plastic handrims.

Newton One wheel only.

- Facilitates unilateral manual wheelchair propulsion.
- Facilitates independent mobility for individuals with hemiplegia or hemiparesis.



CUSTOM COLORED SPINERGY WHEELS

TAILORED

0014849A-XX

**APEX C | APEX A | VELOCE
HELIO C2 | HELIO A7 | HELIO A6 | HELIO kids**

Spinergy wheels in any hub and spokes color variations.

Customer must specify desired Spinergy configuration.

- Facilitates stability and rigidity while being lightweight.
- Offers vibration damping properties.
- Allows user to express their unique self.



VERTICAL PROJECTION HANDRIM EXPRESS

0011423A-00

APEX C | APEX A | VELOCE | HELIO C2
HELIO A7 | HELIO A6 | MOVE | HELIO kids

Vertical projections on handrims, includes protective covers on projections.

Approx. weight: 1.22 lbs - (0.55 kg).

Not compatible with mag wheels. Available for all Motion Composites and Newton handrims.

- Facilitates propulsion for users with limited hand function.
- The vertical direction of the projections does not impact the overall wheelchair width.



WHEEL LOCKS

UNILATERAL WHEEL LOCK WITH EXTENSION

EXPRESS

0004570A-XX

VELOC2 | HELIO C2 | HELIO A7 | HELIO A6
MOVE

Unilateral Pull-to-Lock wheel lock mechanism with extension installed on lever, to lock both wheels from one side only.

- 10-15" left + extension
- 10-15" right + extension
- 16-22" left + extension
- 16-22" right + extension

Approx. weight: 1.48 - 1.59 lbs (0.67-0.72 kg)

Includes unilateral wheel lock and extension installed when sold as parts.

- Facilitates unilateral locking of the manual wheelchair.
- Facilitates independent engagement of wheellocks for individuals with hemiplegia or hemiparesis.



EXTENSION FOR UNILATERAL WHEEL LOCK

EXPRESS

0014849A-XX

VELOC2 | HELIO C2 | HELIO A7 | HELIO A6
MOVE

Unilateral wheel lock extension sent as a add-on to the mechanism. Easy to install.

Extension only.

- Facilitates access to the unilateral wheel lock.
- Facilitates independent engagement of the unilateral wheel lock for clients with limited unilateral reach.



ATTENDANT WHEEL LOCK EXPRESS

0005422A-00

APEX C | APEX A | VELOCE

Attendant wheel lock installed on APEX or VELOCE chair frame.

Approx. weight: 1.11 lbs (0.50 kg)

Only one side included.

- Facilitates caregiver assisted engagement of the wheel locks for the dependent client.



APEX C | APEX A



VELOCE



ANTI-ROLLBACK SYSTEM EXPRESS

0006946A-01

VELOCE | HELIO C2 | HELIO A7
HELIO A6 | MOVE | HELIO kids

Anti-rollback mechanism that locks the rear wheels for added safety when user gets off the wheelchair seat, even when they fail to use standard wheel locks. Can be installed on chair at time of order, or sent to a dealer with proper instructions for easy installation.

- 22" PratiKo wheel lock system
- 24" PratiKo wheel lock system

0 deg camber, no amputee plate, 250 lb max., maximum width 20". Kit sold as parts for the customer to install (contains instructions).

- Prevents the wheelchair from rolling back during user transfers.
- Facilitates safety during user transfers or when standing.

ARMRESTS

HEIGHT ADJUSTABLE FLIP BACK ARMREST EXPRESS

0011495A-XX

APEX C | APEX A

Mechanism to allow armrest to pivot backwards and clear space on chair sides on rigid wheelchairs where it is normally not available.

- Desk length 8-12" height-adjustable flip back armrest on APEX
- Full length 8-12" height-adjustable flip back armrest on APEX

Additional heights may be available as a Tailored or Signature

- The flip back armrests can be flipped back to facilitate transfers, ADLs or proximity to desks and tables.



ARM TROUGH EXPRESS

0011719A-XX

APEX C | APEX A | VELOC2 | HELIO C2
HELIO A7 | HELIO A6 | MOVE | HELIO kids

Arm trough support installed on armrest structure, comes with arm and hand support.

- Rolko Arm trough assembly (left)
- Rolko Arm trough assembly (right)

- Facilitates optimal positioning of the upper extremity at the desired height.
- Facilitates optimal positioning of the upper extremity to abnormal tone or movement.



ARMREST FRONT EXTENSION EXPRESS

0012167A-00

APEX C | APEX A | VELOCE | HELIO C2
HELIO A7 | HELIO A6 | MOVE | HELIO kids

Three inch forward extension to single post or flip back armrests. Comes with a reinforced structure for added safety.

- Armrest front extension 8-12"
- Armrest front extension 10-14"

Approx. weight: 0.17 lbs (0.08 kg)

Full length armpad included.

- Facilitates transfers for users with limited reach.
- Facilitates safety and independence during transfers.



SIDE GUARDS

FOLD DOWN SIDE GUARD BRACKET EXPRESS

0005476A-00

APEX C | APEX A | VELOCE | HELIO C2
HELIO A7 | HELIO A6 | MOVE | HELIO kids

Bracket designed to allow side guards to fold inwards.

Folding wheelchair cannot fold with side guard folded. Bracket only.

- Efficiency for folding and storing the wheelchair.
- Eliminates the need of removing the sideguard to fold down back canes.



ACCESSORIES

APEX SWING AWAY ANTI-TIPPER AND TIP ASSIST

EXPRESS

0006589A-00

APEX C | APEX A

Anti-tipping device for APEX models easily flipped in or under when not in use. Can be activated by user or caregiver to facilitate transitions or obstacles. Equipped with tip assist for caregivers.

Approx. weight: 1.28 lbs (0.58 kg)

APEX only, only one side included.

- These anti-tippers can be flipped in or under when not in use.
- The tip assist facilitates caregiver assistance when traversing over uneven ground.



UNOCCUPIED TRANSIT TIE DOWN **EXPRESS**

0012717A-00

APEX C | APEX A | VELOCE

Complete tie down system for unoccupied transit. Hooks to attach and transport the wheelchair in a vehicle.

Approx. weight: 0.76 lbs (0.35 kg)

Contains all 4 tie downs and hardware.

- Facilitates safety during unoccupied wheelchair transportation.





IV POLE ASSEMBLY WITH CLAMP TAILORED

0017446A-00

**VELOCE | HELIO C2 | HELIO A7
HELIO A6 | MOVE | HELIO kids**

Height adjustable IV pole support, can be installed on either side of the wheelchair.

- Facilitates medication transportation and administration.

OXYGEN TANK HOLDER EXPRESS

0016719A-01

**VELOCE | HELIO C2 | HELIO A7
HELIO A6 | MOVE**

Oxygen tank holder installed on the right at the back of the chair. Pivots to clear space when necessary.

Approx. weight: 1.54 lbs (0.70 kg)

Compatible with D and E size tanks. Available for D or E size tanks.

- Facilitates transportation of oxygen for users with supplemental oxygen needs.





EVO
CUSTOM



SIGNATURE

AESTHETICS & PERSONALIZATION



PERFORMANCE &
MOBILITY ENHANCEMENTS



TITANIUM HANDRIMS –
SATIN FINISH



TC MAX STAINLESS STEEL
SPINERGY HANDRIM



X-CORE WHEELS



ALUMINUM CAMBER TUBE




SMALLER DIAMETER HANDRIM



TITANIUM QUICK RELEASE AXLE



QUAD RELEASE SIDE GUARDS



ALUMINUM SOLID SEAT PAN



BODYPOINT OPTIONS
& ACCESSORIES



THERAFIN POSITIONING ACCESSORIES

SEATING & POSITIONING



OFFSET ARMREST RECEIVERS



FOOTREST & LEG SUPPORT



GEL FOOTPLATE COVER



FOOTREST TRIGGER
EXTENSION



GRIPTAPE ON FOOTPLATE



CUSTOM SIZED ALUMINUM
FOOTPLATE (HELIO)



CUSTOM SIZED CARBON
FOOTPLATE (HELIO)



OFFSET FOOTREST
ANCHOR (20 DEG)




3" FOOTPLATE RISER


ACCESSORIES & ADD-ONS



REMOVABLE LUGGAGE CARRIER FOR HELIO



TRANSFER HANDLES FOR HELIO FRAME



CUP HOLDER

EVO
CUSTOM



PRICE LIST

PRICE LIST - USA

NAV ITEM	NAME	DESCRIPTION	LEVEL	UoM	PRICE
FRONT WHEELS					
0005330A-00	CASTER LOCK PIN 1" TAPER STD STEM BOLT	Mechanism to prevent caster wheels from pivoting to help maintain direction. Can be activated or deactivated.	Express	Pair	\$1,505
0005330A-01	CASTER LOCK PIN STRAIGHT STD STEM BOLT	Mechanism to prevent caster wheels from pivoting to help maintain direction. Can be activated or deactivated.	Express	Pair	\$1,505
0005330A-02	CASTER LOCK PIN 2" TAPER STD STEM BOLT	Mechanism to prevent caster wheels from pivoting to help maintain direction. Can be activated or deactivated.	Express	Pair	\$1,505
FOOTRESTS					
0005399A-XX	CUSTOM CARBON PLATFORM FOOTPLATE	Custom sized carbon footplates for any sizes not mentioned on order form.	Tailored	Each	\$230
0005400A-XX	CUSTOM ALUMINUM PLATFORM FOOTPLATE	Custom sized aluminum footplates for any sizes not mentioned on order form.	Tailored	Each	\$120
0014753A-00	WINGED PLATFORM FOOTPLATE	Custom 3D printed footplate with raised sides installed on an open loop structure.	Express	Each	\$150
0013514A-00	ELEVATING LEG RESTS EXTENSION TUBE	Custom elongated tube for the elevated leg rests. Adds 2" to the longest extension available on order form.	Express	Each	\$40
0006252A-00	SWING AWAY LEFT STUMP SUPPORT 6x8", FLAT (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-02	SWING AWAY LEFT STUMP SUPPORT 8x10", FLAT (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-01	SWING AWAY LEFT STUMP SUPPORT 10x10", FLAT (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-03	SWING AWAY LEFT STUMP SUPPORT 10x14", FLAT (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-04	SWING AWAY LEFT STUMP SUPPORT 6x8", CONTOURED (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-05	SWING AWAY LEFT STUMP SUPPORT 8x10", CONTOURED (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-06	SWING AWAY RIGHT STUMP SUPPORT 6x8", FLAT (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-07	SWING AWAY RIGHT STUMP SUPPORT 8x10", FLAT (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-08	SWING AWAY RIGHT STUMP SUPPORT 10x10", FLAT (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-09	SWING AWAY RIGHT STUMP SUPPORT 10x14", FLAT (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575

NAV ITEM	NAME	DESCRIPTION	LEVEL	UoM	PRICES
0006252A-10	SWING AWAY RIGHT STUMP SUPPORT 6x8", CONTOURED (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0006252A-11	SWING AWAY RIGHT STUMP SUPPORT 8x10", CONTOURED (APEX/VELOCE)	Sturdy and easily removable swing away residual limb support. Available in multiple cushion sizes.	Express	Each	\$575
0011380A-00	HIGH MOUNT FLIP BACK ASSEMBLY ON APEX	Full flip back footrest assembly installed as a high mount on APEX chairs. Easily adjustable to various heights.	Tailored	Each	\$315
0012756A-00	APEX FOOTPLATE ASSEMBLY WITH LONG CLAMP	Custom footrest clamp on apex chairs that allows for height adjustability on tapered chairs.	Tailored	Each	\$105
BACKRESTS					
0009313C-00	HEIGHT ADJUSTABLE SINGLE HANDLE FOR APEX MOUNTED ON THE BACK RIGIDIZER	Height adjustable and removable back cane installed on rear aluminum rigidizer. Single push handle mount on APEX rigidizer.	Express	Each	\$370
0013276A-00	REINFORCED DEPTH ADJUSTABLE BACK OPTION - 18-21"	Custom reinforced depth adjustable backrests to prevent breakage for users with spasms.	Tailored	Each	\$275
0013276A-01	REINFORCED DEPTH ADJUSTABLE BACK OPTION - 15-18"	Custom reinforced depth adjustable backrests to prevent breakage for users with spasms.	Tailored	Each	\$275
0013276A-02	REINFORCED DEPTH ADJUSTABLE BACK OPTION - 12-15"	Custom reinforced depth adjustable backrests to prevent breakage for users with spasms.	Tailored	Each	\$275
0013276A-03	REINFORCED DEPTH ADJUSTABLE BACK OPTION - 9-12"	Custom reinforced depth adjustable backrests to prevent breakage for users with spasms.	Tailored	Each	\$275
0013276A-04	REINFORCED DEPTH ADJUSTABLE BACK OPTION - 8° - 18-21"	Custom reinforced depth adjustable backrests to prevent breakage for users with spasms.	Tailored	Each	\$275
0013276A-05	REINFORCED DEPTH ADJUSTABLE BACK OPTION - 8° - 15-18"	Custom reinforced depth adjustable backrests to prevent breakage for users with spasms.	Tailored	Each	\$275
0013276A-06	REINFORCED DEPTH ADJUSTABLE BACK OPTION - 8° - 12-15"	Custom reinforced depth adjustable backrests to prevent breakage for users with spasms.	Tailored	Each	\$275
0012805A-00	FOLDING BACKREST - FIXED ANGLE - INTEGRATED PUSH HANDLES ALUMINUM CANES 9-15"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012805A-01	FOLDING BACKREST - FIXED ANGLE - INTEGRATED PUSH HANDLES ALUMINUM CANES 15-21"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012805A-02	FOLDING BACKREST - FIXED ANGLE - 8° INTEGRATED PUSH HANDLES ALUMINUM CANES 17-20"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012805A-03	FOLDING BACKREST - FIXED ANGLE - INTEGRATED PUSH HANDLES ALUMINUM CANES 21-24"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012805A-04	FOLDING BACKREST - FIXED ANGLE - 8° INTEGRATED PUSH HANDLES ALUMINUM CANES 21-24"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012541A-00	FOLDING BACKREST - FIXED ANGLE - STRAIGHT 9-12" ALUMINUM CANES	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012541A-01	FOLDING BACKREST - FIXED ANGLE - STRAIGHT 12-15" ALUMINUM CANES	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012541A-02	FOLDING BACKREST - FIXED ANGLE - STRAIGHT 15-18" ALUMINUM CANES	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575

NAV ITEM	NAME	DESCRIPTION	LEVEL	UoM	PRICES
0012541A-03	FOLDING BACKREST - FIXED ANGLE - STRAIGHT 18-21" ALUMINUM CANES	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012541A-04	FOLDING BACKREST - FIXED ANGLE - 8° ALUMINUM CANES 13.5-15"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012541A-05	FOLDING BACKREST - FIXED ANGLE - 8° ALUMINUM CANES 15-18"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012541A-06	FOLDING BACKREST - FIXED ANGLE - 8° ALUMINUM CANES 18-21"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012543A-00	FOLDING BACKREST - FIXED ANGLE - STRAIGHT 9-12" CARBON CANES	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012543A-01	FOLDING BACKREST - FIXED ANGLE - STRAIGHT 12-15" CARBON CANES	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012543A-02	FOLDING BACKREST - FIXED ANGLE - STRAIGHT 15-18" CARBON CANES	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012543A-03	FOLDING BACKREST - FIXED ANGLE - STRAIGHT 18-21" CARBON CANES	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012543A-04	FOLDING BACKREST - FIXED ANGLE - 8° CARBON CANES 13.5-15"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012543A-05	FOLDING BACKREST - FIXED ANGLE - 8° CARBON CANES 15-18"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012543A-06	FOLDING BACKREST - FIXED ANGLE - 8° CARBON CANES 18-21"	Allows backrest to fold down on folding chairs. Available with any back cane configuration. This version is not angle adjustable.	Express	Each	\$575
0012548A-00	FOLDING BACKREST - ADJUSTABLE ANGLE - ALUMINUM CANES 15-18"	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$575
0012548A-01	FOLDING BACKREST - ADJUSTABLE ANGLE - ALUMINUM CANES 18-21"	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$575
0012548A-02	FOLDING BACKREST - ADJUSTABLE ANGLE - ALUMINUM CANES 18-21"	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$575
0012548A-03	FOLDING BACKREST - ADJUSTABLE ANGLE - 8° ALUMINUM CANES 15-18"	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$575
0012573A-00	FOLDING BACKREST - ADJUSTABLE ANGLE - 8° ALUMINUM CANES 18-21"	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$590
0012573A-01	FOLDING BACKREST - ADJUSTABLE ANGLE - 14-15" CARBON CANES	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$590
0012573A-02	FOLDING BACKREST - ADJUSTABLE ANGLE - 15-18" CARBON CANES	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$590
0012573A-03	FOLDING BACKREST - ADJUSTABLE ANGLE - 21" CARBON CANES	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$590
0012573A-04	FOLDING BACKREST - ADJUSTABLE ANGLE - 14-15" 8° CARBON CANES	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$590
0012573A-05	FOLDING BACKREST - ADJUSTABLE ANGLE - 15-18" 8° CARBON CANES	Allows backrest to fold down on folding chairs. Comes with any back cane configuration. This version is angle adjustable.	Express	Each	\$590

NAV ITEM	NAME	DESCRIPTION	LEVEL	UoM	PRICES
UPHOLSTERY & SEATING					
0015196A-00	ADDUCTOR	Custom wheelchair side pads that can be configured in any shape, lengths, heights or thickness. Allows seating area to be reduced.	Tailored	Each	\$110
0006009A-00	LEFT SIDE KNEE ADDUCTOR	Knee adductor available left and right installed on chair to hold user's knees in place.	Tailored	Each	\$480
0006009A-01	RIGHT SIDE KNEE ADDUCTOR	Knee adductor available left and right installed on chair to hold user's knees in place.	Tailored	Each	\$480
0010294A-XX	CARBON FIBRE SOLID SEAT PAN	Carbon fiber solid seat pan to replace upholstery, can fit any chair size, comes with installation hardware	Tailored	Each	\$440
0018808A-00	AIR FORCE PATCH WITH VELCRO FOR NXT BACK SUPPORTS	Air force patch with velcro	Express	Each	\$52
0018808A-01	ARMY PATCH WITH VELCRO FOR NXT BACK SUPPORTS	Army patch with velcro	Express	Each	\$52
0018808A-02	NAVY PATCH WITH VELCRO FOR NXT BACK SUPPORTS	Navy patch with velcro	Express	Each	\$52
0018808A-03	MARINE CORPS PATCH WITH VELCRO FOR NXT BACK SUPPORTS	Marine corps patch with velcro	Express	Each	\$52
0018808A-04	PURPLE HEART PATCH WITH VELCRO FOR NXT BACK SUPPORTS	Purple heart patch with velcro	Express	Each	\$52
0018808A-05	COAST GUARD PATCH WITH VELCRO FOR NXT BACK SUPPORTS	Coast guard patch with velcro	Express	Each	\$52
0018808A-06	SPACE FORCE PATCH WITH VELCRO FOR NXT BACK SUPPORTS	Space force patch with velcro	Express	Each	\$52
0018808A-07	USA FLAG PATCH WITH VELCRO FOR NXT BACK SUPPORTS	Usa flag patch with velcro	Express	Each	\$52
REAR WHEELS					
0005993A-XX	ONE ARM DRIVE KIT WITH NEWTON ONE 22"	Custom wheelchair side pads that can be configured in any shape, lengths, heights or thickness. Allows seating area to be reduced.	Tailored	Each	\$1,660
0010899A-XX	CUSTOM COLOURED SPINERGY CLX WHEELS	Spinergy wheels in any hub and spokes color variations.	Tailored	Each	\$910
0013625A-XX	COLOURED SPINERGY LX WHEELS	Spinergy wheels in any hub and spokes color variations.	Tailored	Each	\$630
0014849A-XX	COLOURED SPINERGY SPOX SPORT WHEELS	Spinergy wheels in any hub and spokes color variations.	Tailored	Each	\$515
0011423A-00	VERTICAL PROJECTION HANDRIM	Vertical projections on handrims, includes protective covers on projections	Express	Each	\$415
WHEEL LOCKS					
0004570A-04	UNILATERAL WHEEL LOCK 16-22 LEFT + EXT (PARTS)	Unilateral wheel lock mechanism with extension installed on lever, to lock both wheels by users from one side only	Express	Each	\$354
0004570A-05	UNILATERAL WHEEL LOCK 16-22 RIGHT + EXT (PARTS)	Unilateral wheel lock mechanism with extension installed on lever, to lock both wheels by users from one side only	Express	Each	\$354
0004570A-06	UNILATERAL WHEEL LOCK 10-15 LEFT + EXT (PARTS)	Unilateral wheel lock mechanism with extension installed on lever, to lock both wheels by users from one side only	Express	Each	\$354
0004570A-07	UNILATERAL WHEEL LOCK 10-15 RIGHT + EXT (PARTS)	Unilateral wheel lock mechanism with extension installed on lever, to lock both wheels by users from one side only	Express	Each	\$354

NAV ITEM	NAME	DESCRIPTION	LEVEL	UoM	PRICES
0011150A-00	EXTENSION FOR UNILATERAL WHEEL LOCK	Unilateral wheel lock extension sent as a add-on to the mechanism. Easy to install.	Express	Each	\$35
0005422A-00	ATTENDANT WHEEL LOCK	Attendant wheel lock installed on apex chair frame. Comes in a pairs.	Express	Pair	\$180
0006946A-01	PRATIKO WHEEL LOCK SYSTEM 22"	Anti-rollback mechanism sold as parts, locks the rear wheels when user gets off the wheelchair seat. Comes with proper instructions for easy installation.	Express	Each	\$990
0006947A-01	PRATIKO WHEEL LOCK SYSTEM 24"	Anti-rollback mechanism sold as parts, locks the rear wheels when user gets off the wheelchair seat. Comes with proper instructions for easy installation.	Express	Each	\$990
ARMRESTS					
0011495A-00	SHORT 8-12 HEIGHT ADJUSTABLE FLIP BACK ARMREST ON APEX	Mechanism to allow armrest to pivot backwards and clear space on chair sides.	Express	Each	\$255
0011495A-01	LONG 8-12 HEIGHT ADJUSTABLE FLIP BACK ARMREST ON APEX	Mechanism to allow armrest to pivot backwards and clear space on chair sides.	Express	Each	\$255
0011495A-02	SHORT 10-14 HEIGHT ADJUSTABLE FLIP BACK ARMREST ON APEX	Mechanism to allow armrest to pivot backwards and clear space on chair sides.	Express	Each	\$255
0011495A-03	LONG 10-14 HEIGHT ADJUSTABLE FLIP BACK ARMREST ON APEX	Mechanism to allow armrest to pivot backwards and clear space on chair sides.	Express	Each	\$255
0011719A-00	LEFT SIDE ROLKO ARM TROUGH ASSEMBLY	Arm trough support installed on armrest structure, comes with arm and hand support.	Express	Each	\$420
0011719A-01	RIGHT SIDE ROLKO ARM TROUGH ASSEMBLY	Arm trough support installed on armrest structure, comes with arm and hand support.	Express	Each	\$420
0012167A-00	T-ARMREST FRONT EXTENSION 10-14	Extension in length to armrests, comes with a reinforcement structure for added safety.	Express	Each	\$60
0012167A-01	T-ARMREST FRONT EXTENSION 8-12	Extension in length to armrests, comes with a reinforcement structure for added safety.	Express	Each	\$60
0012977A-00	GEL PAD ARMREST DESK LENGTH	Custom made gel arm pad, comes in any thickness, width and length to add comfort for user.	Express	Each	\$220
0012977A-01	GEL PAD ARMREST FULL LENGTH	Custom made gel arm pad, comes in any thickness, width and length to add comfort for user.	Express	Each	\$220
SIDE GUARDS					
0005476A-00	FOLD DOWN SIDEGUARD BRACKET	Bracket designed to allow side guards to fold inwards.	Express	Each	\$300
ACCESSORIES					
0006589A-00	SWING AWAY ANTI-TIPPER (APEX)	Sturdy anti-tipping device easily flipped up when not in use. Activated by user to pass over bumps and equipped with tip assist for caregivers	Express	Each	\$270
0012717A-00	UNOCCUPIED TRANSIT TIE DOWN (APEX)	Complete tie-down system for chair transits without user. Hooks to attach and transport the wheelchair in a vehicule	Express	Each	\$95
0013816A-00	UNOCCUPIED TRANSIT TIE DOWN (VELOCE)	Complete tie-down system for chair transits without user.	Express	Each	\$95
0017446A-00	IV POLE ASSEMBLY WITH CLAMP	Adjustable iv pole support, can be installed on both sides of the chair and easily adujstable in height	Tailored	Each	\$410
0016719A-01	RIGHT MOUNTING OXYGEN TANK HOLDER	Custom oxygen tank holder for wheelchairs with certified transit tie down. Pivots easily to clear space when necessary.	Express	Each	\$695

GO BEYOND!

Who we are: passionate innovators pushing the boundaries of what is possible. What we do: combine leading-edge engineering with innovative materials to create movement. What makes us proud: having our ultralight wheelchairs echo the love of life their users demonstrate each day.

Through our unique approach and distinct attitude, our aim is to enhance lives, nourish a sense of freedom and revolutionize the way people see wheelchairs.

LEARN MORE To learn and find more information regarding our products, options and accessories, visit our website at motioncomposites.com



An updated version of this document is available online. Scan this QR code to access it.



T 1 866 650-6555

F 1 888 966-6555

motioncomposites.com

EVO CUSTOMS BROCHURE

