



## EZPZ-SP & EZPZ-AP

Self-Propel and Attendant Propel  
Mobile Shower Commode Chairs



EZPZ-AP  
Attendant Propel

EZPZ-SP  
Self Propel



**EZPZ Mobile Shower Commode Chairs** offer Raz Design quality, exceptional comfort and function at an affordable price. Standard features such as corrosion resistant aluminum frame, soft back support and locking casters make these commodes a great solution for residential and institutional use.

### EZPZ-SP & AP FEATURES

- Flip-up arm supports
- Contoured seat
- Tool-less height adjustment
- Ergonomic commode pan and lid

**EZPZ mobile shower commode chairs** are exceptionally durable and easy to clean... ideal for infection-control environments.

**FEATURES**



Corrosion-resistant aluminum frame with tool-less height adjustment



Height-adjustable without tools (EZPZ-AP) – for toilet bowl clearance and matching transfer heights



Ergonomically contoured seat for best-in-class comfort



Swing-away, removable foot supports – facilitates transfers



Soft and flexible back support – comfortable and easy to clean



Flip-up arm supports – aids in transfers



Dual-locking casters (swivel and roll) – stabilizes chair for toileting, bathing and transfers



Commode pan with lid – includes anti-splash baffle and high-gloss surface to ease cleaning

**OPTIONS**



Soft seat overlay – molded / contoured foam for added pressure distribution, comfort and positioning



Pelvic belt – side-release buckle (*standard webbing or infection control material*)



Chest belt – 1-piece with side-release buckle (*standard webbing or infection control material*)



Snap-on Infection Control Calf Strap – dries quickly



Arm support locks – provide more security during transfers



Swing-away laterals – provide lateral trunk support



Quick-release axles (EZPZ-SP) – makes transporting the chair easier

SPECIFICATIONS	EZPZ-AP	EZPZ-SP
WEIGHT CAPACITY	330 lbs	330 lbs
SEAT WIDTH	19.6"	19.6"
SEAT DEPTH (front of back post to front of seat)	16.9"	16.9"
BACK WIDTH	18.5"	18.5"
BACK HEIGHT (from top of seat to top of upholstery)	16.9"	16.9"
TOILET HEIGHT CLEARANCE	15.3" - 19.2"	15.3" - 19.2"
TOILET WIDTH CLEARANCE	16.9"	16.9"
ARM SUPPORT HEIGHT ABOVE SEAT	9"	9"
SEAT TO FLOOR HEIGHT – 1" (2.5 cm) increments	18.8" - 22.8"	18.8" - 22.8"
OVERALL WIDTH	22.8"	27.9"
PRODUCT WEIGHT	35.2 lbs	39.6 lbs
FRAME MATERIAL	ALUMINUM	ALUMINUM



Motion Composites is the exclusive Canadian Distributor of RAZ DESIGN INC.

motioncomposites.com  
 T 1 (866) 650-6555 | F 1 (888) 966-6555