



**Helio A6**

ORDER FORM - MANITOBA POSSIBLE

**1- Purchasing information**

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_  
 Company Name:\* \_\_\_\_\_  
 Purchasing contact:\* \_\_\_\_\_  
 Phone:\* \_\_\_\_\_  
 Email:\* \_\_\_\_\_  
 Contact preference:\* Phone  Email   
 Marked for: \_\_\_\_\_

**3- Shipping information**

Ship to:\* **Manitoba Wheelchair Program**  
 Name: \_\_\_\_\_  
 Address: **1857 Notre Dame Ave.**  
 City: **Wpg** Province: **MB**  
 Postal Code: **R3E 3E7**  
 Phone: **204-975-3250**

NC = NO CHARGE (option available through Manitoba Possible)  
 J = JUSTIFICATION REQUIRED (provide clinical rationale)  
 OTP = OPTION TO PURCHASE (client responsible for cost)  
 N/A = NOT AVAILABLE (option not available through Manitoba Possible)

MA080   
 English Manual (STD)

Taxes not included

MODEL								
WC19	<input type="checkbox"/>	Ultralight folding Helio A6 Wheelchair - 6000 series aluminum (265lb weight capacity)						NC
WC09	<input type="checkbox"/>	HD Heavy duty kit (350 lb weight capacity)						NC
<b>Frame</b>								
<b>SEAT WIDTH</b> <sup>3</sup>								
SW	<input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <sup>1</sup> <input type="checkbox"/>	NC	
	<input type="checkbox"/>	19" <sup>1</sup> <input type="checkbox"/>	20" <sup>1</sup> <input type="checkbox"/>	21" <sup>2</sup> <input type="checkbox"/>	22" <sup>2</sup> <input type="checkbox"/>		NC	
<b>SEAT DEPTH</b>								
SD	<input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <input type="checkbox"/>	NC	
	<input type="checkbox"/>	19" <input type="checkbox"/>	20" <input type="checkbox"/>				NC	
<b>TRANSIT TIE DOWN</b>								
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)						NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup> Short belt 20"-41"						J
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup> Long belt 20"-46"						J
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt <sup>4</sup>						J

1- Available with HD Kit (350 lb)      2- Only with HD Kit (350 lb)      3. [Width calculator](#)  
 4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

**A6 FRONT Seat-to-floor Height - Please select your choice**

Forks / Position	3" caster (OTP)			4" caster (OTP)			5" caster (J)			6" caster (NC)			7" caster (OTP)			8" caster (NC)			
	Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			
	STD	+1"	+2"																
4"	P1	12"	13"	14"	12 1/2"	13 1/2"	14 1/2"												
	P2	12 1/2"	13 1/2"	14 1/2"	13"	14"	15"		14 1/2"	15 1/2"									
5"	P1																		
	P2				13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
	P3	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"		15 1/2"	16 1/2"						
7"	P4	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"		16 1/2"	17 1/2"			
	P1										15 1/2"	16 1/2"	17 1/2"	16"	17"	18"			
	P2							15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
	P3							16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"
	P4	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"
P5	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form); Orange boxes are not available with a seat slope of more than 1.5in.

FRONT CASTERS

CASTERS			Availability		
CA01-03 to CA01-08	<input type="checkbox"/>	Composite wheel w/ PU tire (1")	3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	<input type="checkbox"/>	Pneumatic casters (1 ¼")	6", 8"		OTP
CA14 to CA16	<input type="checkbox"/>	Newton UltraCaster Composite (1 ½ ")	4", 5", 6"		J
CA17 to CA19	<input type="checkbox"/>	Newton UltraCaster Aluminum (1 ½ ")	4", 5", 6"		OTP
CASTER OPTIONS					
CA12	<input type="checkbox"/>	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) <sup>1</sup>			OTP
CA13	<input type="checkbox"/>	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)			OTP

1- Not available with WC-19

Helio A6 - REAR Seat-to-floor Height - Please select your choice

20" (451) (NC)	22" (501) (NC)	24" (540) (NC)	25" (559) (OTP)	26" (590) (OTP)
12"	13"	13 3/4"	14 1/4"	14 3/4"
12 1/2"	13 1/2"	14 1/4"	14 3/4"	15 1/4"
13"	14"	14 3/4"	15 1/4"	15 3/4"
13 1/2"	14 1/2"	15 1/4"	15 3/4"	16 1/4"
14"	15"	15 3/4"	16 1/4"	16 3/4"
14 1/2"	15 1/2"	16 1/4"	16 3/4"	17 1/4"
15"	16"	16 3/4"	17 1/4"	17 3/4"
15 1/2"	16 1/2"	17 1/4"	17 3/4"	18 1/4"
16"	17"	17 3/4"	18 1/4"	18 3/4"
16 1/2"	17 1/2"	18 1/4"	18 3/4"	19 1/4"
17"	18"	18 3/4"	19 1/4"	19 3/4"
17 1/2"	18 1/2"	19 1/4"	19 3/4"	20 1/4"

Configuration used: 1" seat slope with a center of gravity at 2"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

REAR WHEELS

ONE ARM DRIVE <sup>1</sup>			Only available with Newton One rear wheel, permanent axles and 0° camber			Available						
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim <sup>2,3</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			J				
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim <sup>2,3</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>			OTP				
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)						OTP				
REAR WHEELS						Available						
RW01	<input type="checkbox"/>	Mag <sup>4</sup>				20", 22" 24"		OTP				
RW02	<input type="checkbox"/>	Newton One - spoke wheel <sup>5</sup>				20", 22" 24"		NC				
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"		OTP				
RW04	<input type="checkbox"/>	Spinerigy Spox <b>Blk STD</b> <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>	22", 24", 25", 26"	OTP				
RW05	<input type="checkbox"/>	Spinerigy LX <b>Blk STD</b> <input type="checkbox"/>	Pur <sup>6</sup> <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>	Pink <sup>6</sup> <input type="checkbox"/>	Grn <sup>6</sup> <input type="checkbox"/>	Org <sup>6</sup> <input type="checkbox"/>	22", 24", 25", 26"	OTP
RW07	<input type="checkbox"/>	Spinerigy CLX <sup>7</sup> <b>Blk STD</b> <input type="checkbox"/>	Pur <sup>6</sup> <input type="checkbox"/>	Wht <sup>6</sup> <input type="checkbox"/>	Blue <sup>6</sup> <input type="checkbox"/>	Red <sup>6</sup> <input type="checkbox"/>	Yel <sup>6</sup> <input type="checkbox"/>	Pink <sup>6</sup> <input type="checkbox"/>	Grn <sup>6</sup> <input type="checkbox"/>	Org <sup>6</sup> <input type="checkbox"/>	24", 25"	OTP
TIRES			Solid			Available						
TI01	<input type="checkbox"/>	Soft Urethane 1 3/8"				20", 22" 24"		NC				
TI02	<input type="checkbox"/>	Soft Urethane knobby 1 3/8" <sup>8</sup>				22", 24", 25", 26"		J				
TI04	<input type="checkbox"/>	Soft Urethane 1" - Shox				22", 24", 25", 26"		OTP				
			Pneumatic									
TI05	<input type="checkbox"/>	Pneumatic 1 3/8"				20", 22", 24", 25", 26"		NC				
TI08	<input type="checkbox"/>	Pneumatic with Airless insert 1 3/8" <sup>9</sup>				22", 24", 26"		OTP				
TI06	<input type="checkbox"/>	Pneumatic HP, Low Tread, Puncture resistant, black 1" - SpeedLite				20", 22", 24", 25", 26"		OTP				
TI07	<input type="checkbox"/>	Pneumatic HP, Medium Tread, Puncture resistant, black 1" - TrailBlazer				24", 25", 26"		OTP				

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; 3- Not available with HD option

4- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 5- Not compatible with Shox tires; 6- color spokes require an additional 10 days lead-time;

7- Only available with 1" tires; 8- 22" is only available with Mags and Newton One rear wheels; 9- Not available with Spinerigy CLX, Spinerigy LX 25", 26" and Spox 25", 26" rear wheels

HANDRIMS

HANDRIM			Available					
HR01	<input type="checkbox"/>	Aluminum anodized				20", 22", 24", 25", 26"		NC
HR61	<input type="checkbox"/>	Aluminum black hard anodized <sup>2</sup>				20", 22", 24", 25"		OTP
HR02	<input type="checkbox"/>	Plastic coated <sup>1</sup>				20", 22", 24", 25", 26"		OTP
HR03	<input type="checkbox"/>	High friction coated <sup>1,2,3</sup>				22", 24", 25", 26"		OTP
HR07-HR18 and HR65-HR71	<input type="checkbox"/>	Newton Air Grip <sup>2</sup> <span style="color: red;">See page 6 for the color selection and price</span>				20", 22", 24", 25", 26"		J
HR28	<input type="checkbox"/>	Surge <sup>1,2,3</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"		OTP
HR05	<input type="checkbox"/>	Surge LT <sup>1,2,3</sup> (1 1/8" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"		OTP
HR04	<input type="checkbox"/>	Natural Fit <sup>1,2,3</sup>				22", 24", 25", 26"		OTP
HR62	<input type="checkbox"/>	Simi H <sup>1,2,3</sup> (1 1/4" Oval handrim with high friction strip)				22", 24", 25", 26"		OTP
HR63	<input type="checkbox"/>	Nova H <sup>1,2,3</sup> (1 1/4" Oval handrim with high friction strip)				22", 24", 25", 26"		OTP
HR64	<input type="checkbox"/>	Optimum H <sup>1,2,3</sup> (1 1/8" heart-shaped handrim)				22", 24", 25", 26"		OTP
HR99	<input type="checkbox"/>	Omit handrim						NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

## HANDRIMS OPTIONS

ASSEMBLY POSITION					
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab		NC
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	Full tab		NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut <sup>3</sup>	Short tab 2 positions		J
HR31	<input type="checkbox"/>	Super Narrow Handrim <sup>3</sup>	Short tab 2 positions STD <input type="checkbox"/> Short tab 1 position <input type="checkbox"/>		J
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	<input type="checkbox"/>	Standard grip			OTP
HR27	<input type="checkbox"/>	Super grip			OTP
HR099	<input type="checkbox"/>	Omit thumb piece			OTP


3- Only available with Aluminum and Air Grip handrims;

AXLE					
AX01	<input type="checkbox"/>	Quick release axle <sup>1</sup>			NC
AX02	<input type="checkbox"/>	Permanent axle			J
AX03	<input type="checkbox"/>	Quad release axle <sup>1</sup>			OTP
AXLE PLATE					
AX04	<input type="checkbox"/>	Amputee adjustable axle plate <sup>2</sup>			OTP

1- Not required with one arm drive; 2- May interfere with swing away armrest and limit seat to floor height

REAR WHEEL CAMBER					
CM00	<input type="checkbox"/>	0° camber			NC
CM03	<input type="checkbox"/>	3° camber			NC


### CENTER OF GRAVITY



#### Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 1 1/4"



Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1- Only available with amputee axle plate

WHEEL LOCK					
WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock			NC
WL14	<input type="checkbox"/>	6" removable extension handle for push to lock			J
WL16	<input type="checkbox"/>	Aluminum push to lock with extensions			NC
WL02	<input type="checkbox"/>	Pull to lock			NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock (works best with pneumatic or high pressure tire)			OTP
WL13	<input type="checkbox"/>	6" removable extension handle for pull to loc and grade aid			J
WL05	<input type="checkbox"/>	Scissors Lock			OTP
WL06	<input type="checkbox"/>	Attendant Lock (additional set of wheel lock mounted to the back)			OTP
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock <sup>1</sup>	Left <input type="checkbox"/> Right <input type="checkbox"/>		OTP
WL99	<input type="checkbox"/>	Omit wheel locks			NC

1- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

[https://www.motioncomposites.com/en\\_ca/support-and-education/how-to-documents](https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents)

## FOOTREST

SWING-IN/SWING-OUT FOOTREST					
FR03-01	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 13.5" - 17"		Set at: *	NC
FR03-02	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at: *	NC
FR03-03	<input type="checkbox"/>	60° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at: *	NC
FR04-01	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 13.5" - 17"		Set at: *	NC
FR04-02	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at: *	NC
FR04-03	<input type="checkbox"/>	70° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at: *	NC
FR15-01	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 13.5" - 16" <sup>5</sup>		Set at: *	NC
FR15-02	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 14.5" - 19.5" <sup>5</sup>		Set at: *	NC
FR15-03	<input type="checkbox"/>	80° Swing-in/swing-out footrest - 15.5" - 20.5" <sup>5</sup>		Set at: *	NC
FR06	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Standard (no extension) <sup>1</sup>	7.5" - 11"	Set at: *	NC
FR07	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Short extension <sup>2 1,2</sup>	9.5" - 13"	Set at: *	NC
FR08	<input type="checkbox"/>	90° swing-in/swing-out with front mount - Medium extension <sup>4 1,2</sup>	9.5" - 14.5"	Set at: *	NC
FOOTREST OPTIONS					
FR09	<input type="checkbox"/>	High mount (footplate attached to the hanger with clamps) <sup>3</sup>	5" - 13.5"	Set at: *	J
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, Length adjustable <sup>4</sup>	Left <input type="checkbox"/> Right <input type="checkbox"/> 14" - 22"	Set at:	J
FR99	<input type="checkbox"/>	Omit Footrest			N/A

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest; \* +/- 1/2"

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES						
FP01	<input type="checkbox"/>	Standard - Composite <sup>1</sup>				NC
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite				J
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum				OTP
FP06	<input type="checkbox"/>	One piece adjustable angle flip-up - Aluminum				OTP
FOOTPLATE ACCESSORIES						
FP08	<input type="checkbox"/>	Length adjustable heel loop <sup>2</sup>				NC
FP09	<input type="checkbox"/>	Calf strap				OTP
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap				OTP
<i>1- Not available with front mount; 2- not available with one piece footplate</i>						
SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	<input type="checkbox"/>	Residual limb support - flat 6"x8"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		OTP
FR28	<input type="checkbox"/>	Residual limb support - flat 8"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		OTP
FR29	<input type="checkbox"/>	Residual limb support - flat 10"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		OTP
FR30	<input type="checkbox"/>	Residual limb support - flat 10"x14"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		OTP
FR32	<input type="checkbox"/>	Residual limb support - contour 6"x8"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		OTP
FR33	<input type="checkbox"/>	Residual limb support - contour 8"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		OTP
SEATING AND BACKREST						
BACK UPHOLSTERY						
BH05	<input type="checkbox"/>	Slip-on Nylon back upholstery				NC
BH06	<input type="checkbox"/>	Adjustable tension Nylon back upholstery				OTP
BH99	<input type="checkbox"/>	Omit Back upholstery				NC
	<input type="checkbox"/>	NXT back Support			See pricelist	OTP
SEAT CUSHION						
SC01	<input type="checkbox"/>	2" Cushion				OTP
SC02	<input type="checkbox"/>	3" Cushion				OTP
SC99	<input type="checkbox"/>	Omit seat cushion				NC
	<input type="checkbox"/>	NXT cushion			See price list	OTP
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE (push handle included - no further selection needed)						
BC58-915	<input type="checkbox"/>	Straight aluminum back cane	9" - 15"	increment of 1/2"	Set at <sup>2</sup> :	NC
BC58-1521	<input type="checkbox"/>	Straight aluminum back cane	15" - 21"	increment of 1/2"	Set at <sup>2</sup> :	NC
BC60	<input type="checkbox"/>	Straight aluminum back cane	21" - 24"	increment of 1"	Set at <sup>2</sup> :	J
BC01-16	<input type="checkbox"/>	8° bend aluminum back cane	16"			NC
BC01-18	<input type="checkbox"/>	8° bend aluminum back cane	18"			NC
BC01-20	<input type="checkbox"/>	8° bend aluminum back cane	20"			NC
BC59	<input type="checkbox"/>	8° bend aluminum back cane	21" - 24"	increment of 1"	Set at <sup>2</sup> :	J
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)						
BC11-912	<input type="checkbox"/>	Straight aluminum back cane	9" - 12"	increment of 1/2"	Set at <sup>2</sup> :	J
BC11-1215	<input type="checkbox"/>	Straight aluminum back cane	12" - 15"	increment of 1/2"	Set at <sup>2</sup> :	J
BC11-1518	<input type="checkbox"/>	Straight aluminum back cane	15" - 18"	increment of 1/2"	Set at <sup>2</sup> :	J
BC11-1821	<input type="checkbox"/>	Straight aluminum back cane	18" - 21"	increment of 1/2"	Set at <sup>2</sup> :	J
BC12-1315	<input type="checkbox"/>	8° bend aluminum back cane	13.5" - 15"	increment of 1/2"	Set at <sup>2</sup> :	OTP
BC12-1518	<input type="checkbox"/>	8° bend aluminum back cane	15" - 18"	increment of 1/2"	Set at <sup>2</sup> :	OTP
BC12-1821	<input type="checkbox"/>	8° bend aluminum back cane	18" - 21"	increment of 1/2"	Set at <sup>2</sup> :	OTP
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)						
BC13-1518	<input type="checkbox"/>	Straight aluminum back cane	15" - 18"	increment of 1/2"	Set at <sup>2</sup> :	J
BC13-1821	<input type="checkbox"/>	Straight aluminum back cane	18" - 21"	increment of 1/2"	Set at <sup>2</sup> :	J
BC14-1518	<input type="checkbox"/>	8° bend aluminum back cane	15" - 18"	increment of 1/2"	Set at <sup>2</sup> :	OTP
BC14-1821	<input type="checkbox"/>	8° bend aluminum back cane	18" - 21"	increment of 1/2"	Set at <sup>2</sup> :	OTP
		Select back cane angle:	85° <sup>1</sup> <input type="checkbox"/>	90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>
					105° <input type="checkbox"/>	110° <input type="checkbox"/>
PUSH HANDLES OPTIONS (ONE option at a time)						
PH01	<input type="checkbox"/>	Push handles				NC
PH05	<input type="checkbox"/>	Fold down push handles <sup>3,5</sup>				OTP
PH06-920	<input type="checkbox"/>	Clamp on height adjustable push handles 9"-20" <sup>3,4,5</sup>				OTP
PH06-923	<input type="checkbox"/>	Clamp on height adjustable push handles 9"-23" <sup>3,4,5</sup>				OTP
PH02	<input type="checkbox"/>	Omit push handles				NC
STROLLER HANDLE AND STABILIZER BAR (only available with push handle)						
PH04	<input type="checkbox"/>	Stroller handle (not available with the HD option)				OTP
PH03	<input type="checkbox"/>	Newton folding stabilizer bar (not available with stroller handle)				OTP

1- Not compatible with 8 degree bent back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265lb weight limit

4- Not compatible with swing away armrests, 8° bend back canes, 21-24" back canes and amputee axle plates; 5- Not compatible with stroller handler, stabilizer bar, headrest and attendant lock;

ARMRESTS							
<b>HEIGHT ADJUSTABLE FLIP BACK ARMREST</b> (Convertible to "T" armrest) <sup>1,4</sup>							
AR01-812	<input type="checkbox"/>	10" Desk length armpad	8" - 12"	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>		NC
AR01-1014	<input type="checkbox"/>	10" Desk length armpad	10" - 14"	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>		NC
AR02-812	<input type="checkbox"/>	14" Full length armpad	8" - 12"	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>		NC
AR02-1014	<input type="checkbox"/>	14" Full length armpad	10" - 14"	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>		NC
<b>HEIGHT ADJUSTABLE "T" ARMREST</b> <sup>1</sup>							
AR13-610	<input type="checkbox"/>	10" Desk length armpad <sup>5</sup>	6" - 10"				NC
AR13-812	<input type="checkbox"/>	10" Desk length armpad	8" - 12"				NC
AR13-1014	<input type="checkbox"/>	10" Desk length armpad	10" - 14"				NC
AR14-610	<input type="checkbox"/>	14" Full length armpad <sup>5</sup>	6" - 10"				NC
AR14-812	<input type="checkbox"/>	14" Full length armpad	8" - 12"				NC
AR14-1014	<input type="checkbox"/>	14" Full length armpad	10" - 14"				NC
<b>SWING AWAY ARMREST</b> <sup>2</sup>							
AR05	<input type="checkbox"/>	Tubular armrest - Short pad	8" - 11"				J
AR23	<input type="checkbox"/>	Tubular armrest - long pad	8" - 11"				J
AR06	<input type="checkbox"/>	10" Desk length flat armpad	8" - 11"				J
AR17	<input type="checkbox"/>	14" Full length flat armpad	8" - 11"				J
<b>ARMREST OPTIONS</b>							
AR12	<input type="checkbox"/>	14" Full Length armpad - Gel Ovation (pair) <sup>3</sup>					OTP
AR99	<input type="checkbox"/>	Omit armrest					J

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;  
3- Only available with full length armrest; Not available with tubular (short and long) and short armrest; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

ACCESSORIES							
<b>SIDE GUARDS</b> (not available with Flip Back and "T" Armrest)							
SG01	<input type="checkbox"/>	Plastic side guards (removable)		Straight bracket - STD	Offset bracket <sup>1</sup>		OTP
SG02	<input type="checkbox"/>	Carbon fiber side guards (fixed)					OTP
SG03	<input type="checkbox"/>	Carbon fiber side guards (removable)		Straight bracket - STD	Offset bracket <sup>1</sup>		OTP
SG06	<input type="checkbox"/>	Plastic fender side guards (removable)		Straight bracket - STD	Offset bracket <sup>1</sup>		OTP
SG04	<input type="checkbox"/>	Carbon fiber fender side guards (removable)		Straight bracket - STD	Offset bracket <sup>1</sup>		OTP
SG99	<input type="checkbox"/>	Omit side guard					NC

1- only compatible with 0° camber; increase distance between fender side guards by 1 1/8";

POSITIONING BELT							
SB01-S	<input type="checkbox"/>	Velcro adjustable - Standard	35" <sup>1</sup>				OTP
SB01-L	<input type="checkbox"/>	Velcro adjustable - Long	60"				OTP
SB02-S	<input type="checkbox"/>	Auto buckle - Standard	35" <sup>1</sup>				NC
SB02-M	<input type="checkbox"/>	Auto buckle - Medium	48"				NC
SB02-L	<input type="checkbox"/>	Auto buckle - Long	60"				NC
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Medium	56"				OTP
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Long	61"				OTP
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Small	51"				OTP
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Medium	56"				OTP
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Long	61"				OTP
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Small	51"				OTP
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Medium	56"				OTP
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Long	61"				OTP
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)					OTP
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium	18"				OTP
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long	24"				OTP
SB99	<input type="checkbox"/>	Omit positioning belt					NC

1- Not available with wheelchair width of 18" and more

HEADREST							
AC01-S	<input type="checkbox"/>	Profiled Headrest - Small	(3.5"X6")				OTP
AC01-M	<input type="checkbox"/>	Profiled Headrest - Medium	(4.5"X7")				OTP
AC01-L	<input type="checkbox"/>	Profiled Headrest - Large	(5.5"X8")				OTP
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar	<sup>1</sup>				OTP

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar


ACCESSORIES						
AC03	<input checked="" type="checkbox"/>	Anti-tippers				NC
AC29	<input type="checkbox"/>	Swing-away anti-tipper <sup>1</sup>	Left <input type="checkbox"/> Right <input type="checkbox"/>			OTP
AC27	<input type="checkbox"/>	Tip assist <sup>2</sup>				OTP
AC06	<input type="checkbox"/>	Cane and crutch holder	Left <input type="checkbox"/> Right <input type="checkbox"/>			OTP
AC07	<input type="checkbox"/>	Oxygen tank holder <sup>4</sup> - D and E Type - <i>Installed to the left</i>				OTP
AC08	<input type="checkbox"/>	Spoke guards <sup>3</sup>	Available: 20", 22", 24", 26"			OTP
AC09	<input type="checkbox"/>	Footrest neoprene impact guard (pair) - Black				OTP
AC10	<input type="checkbox"/>	Seat pouch				OTP
AC30	<input type="checkbox"/>	Newton multi-purpose pouch				OTP
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag				OTP
AC23	<input type="checkbox"/>	Newton Backpack				OTP
AC12	<input type="checkbox"/>	Velcro Closing Strap				OTP
AC13	<input type="checkbox"/>	Tool kit				OTP

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate; 2- Cannot be combined with Anti Tipppers  
3- Not available with Mag wheels and Spinergy; 4- Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down

Colors																				
Options <sup>1</sup>	Black	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	<input type="checkbox"/>	OTP																		
Frame																				NC

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

### Fork and Stem bolt combinations



### Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.25 to 2.75
	16	0.75 to 2.25
22 in	17+	0.75 to 2.75
	16	0.75 to 1.75
24 in	17	0.75 to 2.25
	18+	0.75 to 2.75



## CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The EVO Custom Request Form is the best place to start.

Note :