



Helio A6

ORDER FORM - USA

P.O. number: * _____ Date: * _____
Company Name: * _____ Phone: * _____
Ordered by: * _____
Clinician/ATP/RTS: _____
Email: _____
Marked for: _____

Ship to: * _____
Name: _____
Address: _____
City: _____ State: _____
Zip Code: _____
Phone: _____

* Required								MA082	MA089	Taxes not included		
								English Manual USA (STD)	Spanish Manual USA			
MODEL								USD				
1	WC19	Ultralight folding Helio A6 Wheelchair - 6000 series aluminum (265lb weight capacity)							K0005	\$2 095	☐	
2	WC09	HD Heavy Duty package (350 lb weight capacity)								\$500	☐	
SEAT WIDTH ³												
3	SW	Seat width:	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐			NC	☐
4	SW	Seat width:		20" ¹ ☐	21" ² ☐	22" ² ☐				E2201	\$495	☐
SEAT DEPTH												
5	SD	Seat depth:	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐			NC	☐
6	SD	Seat depth:		20" ☐						E2203	\$475	☐

1 - Available with HD Kit (350 lb)

2 - Only with HD Kit (350 lb)

[3. Width calculator](#)

A6 FRONT Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
Fork/Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12	13	14	12 1/2	13 1/2	14 1/2												
	P2	12 1/2	13 1/2	14 1/2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1																		
	P2				13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13	14	15	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1										15 1/2	16 1/2	17 1/2	16	17	18			
	P2							15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3							16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
7	CA01	Composite wheel w/ PU tire (1" width)	3" ☐	4" ☐	5" ☐	6" ☐	7" ☐	8" ☐			NC	☐
8	CA10	Pneumatic casters 6" x 1 ¼"								E2214/E2215	\$160	☐
9	CA11	Pneumatic casters 8" x 1 ¼"								E2214/E2215	\$160	☐
10	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "								\$80	☐
11	CA15		5" x 1 ½ "								\$80	☐
12	CA16		6" x 1 ½ "								\$80	☐
13	CA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½ "								\$175	☐
14	CA18		5" x 1 ½ "								\$185	☐
15	CA19		6" x 1 ½ "								\$195	☐
CASTER OPTIONS												
16	CA12	Caster Pin Locks (<i>Available with +1" or +2" Stem bolt ONLY</i>)								K0073	\$295	☐
17	CA13	Carbon fiber Frog Legs suspension forks (<i>Not available with caster pin locks</i>)								E1015	\$475	☐

Helio A6 - REAR Seat-to-floor Height - Please select your choice

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12	13	13 3/4	14 1/4	14 3/4
12 1/2	13 1/2	14 1/4	14 3/4	15 1/4
13	14	14 3/4	15 1/4	15 3/4
13 1/2	14 1/2	15 1/4	15 3/4	16 1/4
14	15	15 3/4	16 1/4	16 3/4
14 1/2	15 1/2	16 1/4	16 3/4	17 1/4
15	16	16 3/4	17 1/4	17 3/4
15 1/2	16 1/2	17 1/4	17 3/4	18 1/4
16	17	17 3/4	18 1/4	18 3/4
16 1/2	17 1/2	18 1/4	18 3/4	19 1/4
17	18	18 3/4	19 1/4	19 3/4
17 1/2	18 1/2	19 1/4	19 3/4	20 1/4

Configuration used: 1" seat slope with a center of gravity at 2"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest										
SWING-IN/SWING-OUT FOOTREST										
18	FR03	60° Swing-in/swing-out footrest	13.5 - 17" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at: *			NC	□
19	FR04	70° Swing-in/swing-out footrest	13.5 - 17" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at: *			NC	□
20	FR15	80° Swing-in/swing-out footrest	13.5 - 16" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at: *			NC	□
90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT <i>(Footplate mounted to hanger tube)</i>										
21	FR06	90° swing-in/swing-out with front mount - Standard (no extension) ¹			7.5 - 11"	Set at: *			NC	□
22	FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}			9.5 - 13"	Set at: *			NC	□
23	FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}			9.5 - 14.5"	Set at: *			NC	□
FOOTREST OPTIONS ³										
24	FR09	High mount <i>(footplate attached to the hanger with clamps)</i>			5 - 13.5" □	Set at: *			\$75	□
25	FR99	Omit Footrest							NC	□

1 - Available only with two pieces adjustable footplates 2 - May interfere with front casters 3- Not available with front mount * +/- 1/2"

ELEVATING LEGREST <i>(unit)</i>			The elevating legrest replace the footrest side selected					
26	FR10	Elevating legrest with calf pad, 14" - 22" Length adjustable	Left ☐	Right ☐	Set at:	E0990	\$240	☐

FOOTPLATES											
27	FP01	Standard (<i>Composite</i>)								STD	☐
28	FP02	Newton adjustable angle (<i>Composite</i>)							K0040	\$210	☐
29	FP03	Newton adjustable angle (<i>Aluminum</i>)							K0040	\$230	☐
30	FP06	One piece adjustable angle flip-up (<i>Aluminum</i>)							K0040	\$190	☐
FOOTPLATE ACCESSORIES											
31	FP08	Length adjustable heel loop (<i>not available with one piece footplate</i>)							E0951	\$40	☐
32	FP09	Calf strap							K0038	\$35	☐
33	FP30	Bodypoint Aeromesh Padded calf strap							K0038	\$115	☐

RESIDUAL LIMB SUPPORT FLAT CUSHION (Swing away, depth, width, height and angle ajustable) - (unit)								
34	FR27	Residual limb support swing away- flat 6"x8"	Left ☐	Right ☐	New	E1020	\$405	☐
35	FR28	Residual limb support swing away - flat 8"x10"	Left ☐	Right ☐	New	E1020	\$405	☐
36	FR29	Residual limb support swing away - flat 10"x10"	Left ☐	Right ☐	New	E1020	\$405	☐
37	FR30	Residual limb support swing away - flat 10"x14"	Left ☐	Right ☐	New	E1020	\$405	☐
RESIDUAL LIMB SUPPORT CONTOUR CUSHION (Swing away, depth, width, height and angle ajustable) - (unit)								
38	FR32	Residual limb support swing away - contour 6"x8"	Left ☐	Right ☐	New	E1020	\$405	☐
39	FR33	Residual limb support swing away - contour 8"x10"	Left ☐	Right ☐	New	E1020	\$405	☐

Back canes with integrated push handle

FIXED ANGLE BACK								
40	BC58	Straight aluminum back cane	14" ☐	16" ☐	18" ☐	20" ☐	NC	☐
41	BC01	8° bend aluminum back cane	16" ☐	18" ☐	20" ☐		NC	☐

Back canes, height adjustable, with or without push handle

FIXED ANGLE BACK								
42	BC11	Straight aluminum back cane	12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐		NC	☐
43	BC12	8° bend aluminum back cane	13.5 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐		\$50	☐

ADJUSTABLE ANGLE BACK								
44	BC13	Adjustable angle, straight aluminum back cane	15 - 18" ☐	18 - 21" STD ☐			\$200	☐
45	BC14	Adjustable angle, 8° bend aluminum back cane	15 - 18" ☐	18 - 21" STD ☐			\$250	☐
46		Select back cane angle:	85° ¹ ☐	90° STD ☐	95° ☐	100° ☐	105° ☐	110° ☐

Set height to: ² (measured from seat sling to top of back upholstery)

PUSH HANDLE OPTION								
47	PH01	Push handles					STD	☐
48	PH05	Folding push handles ^{3,5}					\$120	☐
49	PH06	Clamp on height adjustable push handle ^{3,4,5}	9 - 20" STD ☐	9-23" ☐			\$300	☐
50	PH02	Without push handles					NC	☐

1 - Not compatible with 8 degree bent back 2 - If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit

4 - Not compatible with swing away armrests , 8° bend back canes and amputee axle plates.

5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

STROLLER HANDLE AND STABILIZER BAR ¹ (only available with push handle)								
51	PH04	Stroller handle (not available with the HD option)				K0108	\$160	☐
52	PH03	Newton folding stabilizer bar					\$250	☐

1 - Stroller handle and stabilizer bar can't be combined


BACK UPHOLSTERY								
53	BH05	Slip-on Nylon back upholstery					STD	☐
54	BH06	Adjustable tension Nylon back upholstery				E2611	\$300	☐
55	BH99	Omit Back upholstery					NC	☐
RIGID BACK								
56		NXT back series				See pricelist		

SEAT CUSHION								
57	SC01	2" Cushion					\$165	☐
58	SC02	3" Cushion					\$185	☐
59	SC99	Omit seat cushion					STD	☐
60		NXT FIT cushion series				See pricelist		

Rear wheels									
ONE ARM DRIVE ¹ <i>Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims</i>									
61	HR22	One arm drive ^{2,3} 24"	Left ☐	Right ☐				1 395 \$	☐
62	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)						NC	☐
REAR WHEELS									
63	RW01	Mag ⁴	20" ☐	22" ☐	24" ☐			\$50	☐
64	RW02	Newton One - spoke wheel ⁵	20" ☐	22" ☐	24" ☐			NC	☐
65	RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐	\$300	☐
66	RW04	Spinergy Spox		22" ☐	24" ☐	25" ☐	26" ☐	\$795	☐
67		Spinergy Spox spokes colors ⁶	Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC	☐
68	RW05	Spinergy LX		22" ☐	24" ☐	25" ☐	26" ☐	\$995	☐
69		Spinergy LX spokes colors ⁶	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC
70	RW07	Spinergy CLX Carbon Blade			24" ☐	25" ☐		\$1 480	☐
71		Spinergy CLX spokes colors ⁶	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC
						Pink ☐	Grn ☐	Org ☐	☐
TIRES					Available				
72	TI01	Soft Urethane 1 1/4"			20", 22" 24"		E2218	STD	☐
73	TI02	Soft Urethane knobby 1 3/8 "			22", 24", 25", 26"		E2218	\$95	☐
74	TI04	Soft Urethane - Shox 1"			22", 24", 25", 26"		E2218	\$115	☐
75	TI05	Pneumatic 1 3/8 "			20", 22", 24", 25", 26"		E2211/E2212	NC	☐
76	TI08	Pneumatic with Airless insert 1 3/8"			22", 24"		E2211/E2213	\$80	☐
77	TI06	Pneumatic HP - 1"			20", 22", 24", 25", 26"		E2211/E2212	\$95	☐
78	TI07	Pneumatic HP - Schwalbe Marathon Plus 1"			24", 25", 26"		E2211/E2212	\$125	☐
1 - Offered with handrims on each side of the wheelchair 2 - Only standard handrim position available 3 - Not available with HD option 4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 5 - Not compatible with Shox tires 6 - Colour spokes require an additionnal 10 days lead-time									
Handrims									
HANDRIM									
79	HR01	Aluminum anodized						STD	☐
80	HR02	Plastic coated ³					WAMO	\$95	☐
81	HR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}					WAMO	\$350	☐
82	HR07	Newton Air Grip ¹ Black ☐					WAMO	\$95	☐
83	HR08	Newton Air Grip ¹ (08) Sunkissed Orange ☐	(09) Burgundy ☐	(10) Sapphire Blue ☐			WAMO	\$240	☐
		(11) Monster Green ☐	(12) Charcoal ☐	(13) Fuchsia ☐	(14) Steel Blue ☐				
		(15) Ferrari Red ☐	(16) Acid Green ☐	(17) White ☐	(18) Sky Blue ☐				
84	HR28	Surge ^{1, 2, 3} (1-1/2" Oval handrim w/ Gripton high friction strip)						\$375	☐
85	HR05	Surge LT ^{1, 2, 3} (1-1/8" Oval handrim w/ Gripton high friction strip)						\$405	☐
86	HR04	Natural Fit ^{1, 2, 3}						\$310	☐
87	HR99	Omit handrim						NC	☐
1 - Not available with Mag wheels or OAD 2 - Not available in 20" 3- Only standard handrim position available									
Handrims option									
Assembly position for aluminum and Newton Air Grip handrim									
88	HR29	Standard Handrim position Full tab						STD	☐
89	HR30	Narrow Handrim position Short tab 2 positions						\$25	☐
90	HR31	Super Narrow Handrim Short tab 2 positions STD ☐	Short tab 1 position ☐					\$25	☐
Thumb piece ⁴ selection for Natural Fit handrims									
91	HR23	Standard grip						\$85	☐
92	HR27	Super grip						\$140	☐
93	HR099	Omit thumb piece						NC	☐
4- Compatible with Surge, Surge LT handrims									
AXLE									
94	AX01	Quick release axle						STD	☐
95	AX02	Permanent axle						NC	☐
96	AX03	Quad release axle					K0108	\$125	☐
AXLE PLATE									
97	AX04	Amputee adjustable axle plate (may interfere with swing away armrest)					E0959	\$195	☐

REAR WHEEL CAMBER					
98	CM00	0° camber			STD ☐
99	CM03	3° camber			NC ☐

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 1" 1/4"

Most stable						Most responsive	
AMPUTÉ	3/4"	1" 1/4	1" 3/4	2" 1/4	2" 3/4		
☐	☐	☐	☐	☐	☐		

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK					
100	WL08	Newton Wheel lock - push to lock			STD ☐
101	WL02	Pull to lock			NC ☐
102	WL04	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>		E0974	\$150 ☐
103	WL13 WL14	6" removable extension handle		E0961	\$75 ☐
104	WL05	Scissor Lock			\$140 ☐
105	WL06	Attendant Lock (additional set of wheel lock mounted to the back)			\$200 ☐
106	WL99	Omit wheel locks			NC ☐

Armrests

FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ¹						Height Adjustment		
107	AR01	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐	STD	E0973 \$365 ☐
108	AR02	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐	STD	E0973 \$365 ☐
"T" ARMREST ¹						Height Adjustment		
109	AR13	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐				E0973 \$325 ☐
110	AR14	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐				E0973 \$325 ☐
SWING AWAY ARMREST ²						Height Adjustment		
111	AR05	Tubular arm pad	8 - 11"					E0973 \$225 ☐
112	AR06	10" Desk length arm pad	8 - 11"					E0973 \$225 ☐
GEL ARM PAD ³								
113	AR12	14" Full Length Arm Pad - Gel Ovation <i>(pair)</i>						\$150 ☐
114	AR99	Omit armrest						NC ☐

1 - May have a 1/2" adjustment restriction on both limits of the range

2 - Could interfere with some rigid back rests

3 - Not available with the swing away armrest

SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>								
115	SG01	Plastic side guards <i>(removable)</i>						\$195 ☐
116	SG02	Carbon fiber side guards <i>(fixed)</i>						\$250 ☐
117	SG03	Carbon fiber side guards <i>(removable)</i>						\$290 ☐
118	SG06	Plastic fender side guards <i>(removable)</i> ¹						\$270 ☐
119	SG04	Carbon fiber fender side guards <i>(removable)</i>						\$540 ☐
120	SG99	Omit side guard						NC ☐
SIDE GUARDS OPTION								
121	SG07	Side guards clip <i>(pair)</i> ²				New		NC ☐

1- Not compatible with 20" and 22" rear wheels

2- Not available with fixed carbon fiber side guard

POSITIONING BELT									
122	SB01	Velcro adjustable		Standard 35" ☐	Long 60" ☐		E0978	\$45	☐
123	SB02	Auto buckle		Standard 35" ☐	Long 60" ☐		E0978	\$85	☐
124	SB07	Hardware attachment for 2 point seatbelt (1 pair)						NC	☐
125	SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)						\$15	☐
126	SB04	BodyPoint Hip Belt - 2 point not padded		Medium 56" ☐	Long 61" ☐		E0978	\$80	☐
127	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐		E0978	\$170	☐
128	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐		E0978	\$190	☐
129	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps		Medium 18" ☐	Long 24" ☐		E0978	\$285	☐
130	SB99	Omit positioning belt						STD	☐

HEADREST									
131	AC01	Profiled Headrest	Small (3.5"x6") ☐	Medium (4.5"x7") ☐	Large (5.5"x8") ☐		E0955	\$180	☐
132	AC28	Headrest mounting hardware with removable non folding stabilizer bar ¹					E1028	\$295	☐

¹ - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES									
133	AC03	Anti tippers					E0971	\$140	☐
134	AC27	Tip assist ¹						\$60	☐
135	AC04	Transit Tie Down (<i>Unoccupied</i>)						NC	☐
136	AC05	WC-19 Transit Tie Down w/Q'straint belt ²		Short belt 20-41" ☐	Long belt 20-46" ☐			\$320	☐
137	AC06	Cane Carrier	Left ☐	Right ☐		E2207	\$195	☐	
138	AC07	Oxygen tank holder - E Type	Left			E2208	\$310	☐	
139	AC08	Spoke guards ³		Available: 20", 22", 24", 26"		K0065	\$115	☐	
140	AC09	Footrest neoprene impact guard (pair) - Black				K0108	\$90	☐	
141	AC10	Seat pouch					\$75	☐	
142	AC26	BodyPoint Mobility Bag					\$125	☐	
143	AC23	Newton Backpack					\$120	☐	
144	AC12	Velcro Closing Strap					\$50	☐	
145	AC13	Tool kit					\$45	☐	
146	AC16	Motion Composites Polo Shirt	W/S ☐	W/M ☐	W/L ☐			\$32	☐
			M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	M/XXXL ☐	M/XXXXL ☐
147	AC15	Military Patches	Army ☐	Navy ☐	Air Force ☐	Marines ☐	Coast Guard ☐	\$90	☐
148	AC15	Combat Wounded Purple Heart Patch						\$90	☐


¹ - Cannot be combined with Anti Tippers ² - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

³ - Not available with Mag wheels and Spinergy

FRAME COLOR (with black hydroformed aluminum crossbrace)									
149	FC01	Black ☐						NC	☐
150	FC03	Sunkissed Orange ☐	Burgundy ☐	Sapphire Blue ☐	Monster Green ☐			NC	☐
		Charcoal ☐	Fuchsia ☐	Steel Blue ☐	Ferrari Red ☐	Acid Green ☐			
151	FC99	Custom Color	Specify PPG Paintit color:				additional 10 days lead-time	\$360	☐

TOTAL							
-------	--	--	--	--	--	--	--

Fork and Stem bolt combinations



Notes :



BioFit™ CUSHION

BioFit™	BioFit™	Fits Width	Depth	\$	✓
 HCPCS CODE: E2607-E2608	NXCABF1616	15"-16"	16"	425	
	NXCABF1618	15"-16"	18"	425	
	NXCABF1620	15"-16"	20"	425	
	NXCABF1816	17"-18"	16"	425	
	NXCABF1818	17"-18"	18"	425	
	NXCABF1820	17"-18"	20"	425	

BioFit™	Fits Width	Depth	\$	✓
NXCABF2016	19" - 20"	16"	425	
NXCABF2018	19" - 20"	18"	425	
NXCABF2020	19" - 20"	20"	425	
NXCABF2216	21" - 22"	16"	455	
NXCABF2218	21" - 22"	18"	455	
NXCABF2220	21" - 22"	20"	455	

Optima™ THORACIC BACK SUPPORT

Optima™	Optima™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NOACB1613	16"-18"	13"	660	
	NOACB1616	16"-18"	16"	660	
	NOACB1619	16"-18"	19"	660	
	NOACB1813	18"-20"	13"	660	
	NOACB1816	18"-20"	16"	660	
	NOACB1819	18"-20"	19"	660	
	NOACB2016	20"-22"	16"	660	
	NOACB2019	20"-22"	19"	660	
	NOACB2213	22"-24"	13"	710	
	NOACB2216	22"-24"	16"	710	
	NOACB2219	22"-24"	19"	710	

Optima™ DEEP THORACIC BACK SUPPORT

Optima™ DEEP	Optima™ DEEP	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NOPCB1616	16"-18"	16"	660	
	NOPCB1619	16"-18"	19"	660	
	NOPCB1816	18"-20"	16"	660	
	NOPCB1819	18"-20"	19"	660	
	NOPCB2016	20"-22"	16"	660	
	NOPCB2019	20"-22"	19"	660	
	NOPCB2213	22"-24"	13"	710	
	NOPCB2216	22"-24"	16"	710	
	NOPCB2219	22"-24"	19"	710	

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™	Xtend™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NHAACB1616/20	16"-18"	16"-20"	660	
	NHAACB1816/20	18"-20"	16"-20"	660	
	NHAACB2016/20	20"-22"	16"-20"	660	
	NHAACB2216/20	22"-24"	16"-20"	710	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™ DEEP	Xtend™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NHAACB1616/20	16"-18"	16"-20"	660	
	NHAACB1816/20	18"-20"	16"-20"	660	
	NHAACB2016/20	20"-22"	16"-20"	660	
	NHAPCB2216/20	22"-24"	16"-20"	710	

Classic™ PELVIC BACK SUPPORT

Classic™	Classic™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NACB1616	16"-18"	16"	660	
	NACB1618	16"-18"	18"	660	
	NACB1816	18"-20"	16"	660	
	NACB1818	18"-20"	18"	660	
	NACB1820	18"-20"	20"	660	
	NACB2020	20"-22"	20"	660	
	NACB2216	22"-24"	16"	710	
	NACB2218	22"-24"	18"	710	
	NACB2220	22"-24"	20"	710	

Classic™ DEEP PELVIC BACK SUPPORT

Classic™ DEEP	Classic DEEP™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NPCB1616	16"-18"	16"	660	
	NPCB1618	16"-18"	18"	660	
	NPCB1818	18"-20"	18"	660	
	NPCB1820	18"-20"	20"	660	
	NPCB2020	20"-22"	20"	660	
	NPCB2216	22"-24"	16"	710	
	NPCB2218	22"-24"	18"	710	
	NPCB2220	22"-24"	20"	710	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

 HCPCS CODE: E0955		
Model #NIB-0784	\$	✓
10" x 5" Large Rectangular	195	

 HCPCS CODE: E0955		
Model # NIB-0818	\$	✓
14" x 5" Extra-Large Rectangular	195	

MultiFit™ HEAD SUPPORT – FIXTURES

		
Model #NIB-0747	\$	✓
Multi-Axis Hardware	145	

		
Model # NIB-0823	\$	✓
Head Support Receiver Bracket	55	

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com