



2915 Ogletown RD # 2270
Newark, DE 19713
T.+1-866-650-6555 F. +1-888-966-6555
quotes@motioncomposites.com
orders@motioncomposites.com
www.motioncomposites.com



6.8

Helio A6

ORDER FORM - USA

1- Purchasing information

P.O. number:* _____ Date: _____
Company Name:* _____
Purchasing contact:* _____
Phone:* _____
Email:* _____
Contact preference:* Phone ☐ Email ☐
Marked for: _____

2- ATP/Clinician information

ATP/Clinician:* _____
Phone:* _____
Email:* _____
Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:* _____
Name: _____
Address: _____
City: _____ State: _____
Zip Code: _____
Phone: _____

* Required

MA082 ☐
English Manual (STD)

MA089 ☐
Spanish Manual

Taxes not included

MODEL		English Manual (SP1)		Spanish Manual		Prices Not Included					
						USD					
WC19	☐	Ultralight folding Helio A6 Wheelchair - 6000 series aluminum (265lb weight capacity)				K0005	\$2,195				
WC09	☐	HD Heavy duty kit (350 lb weight capacity) ⁵					\$550				
Frame											
SEAT WIDTH ³											
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐			NC	
	☐						20" ¹ ☐	21" ² ☐	22" ² ☐		E2201
SEAT DEPTH											
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐			NC	
	☐						20" ☐			E2203	\$575
FRAME LENGTH ⁶											
SP01	☐	Standard frame length configuration ⁷								STD	
SP02	☐	Frame length 1" longer than standard configuration								NC	
SP03	☐	Frame length 2" longer than standard configuration								NC	
SP04	☐	Frame length 1" shorter than standard configuration								NC	
TRANSIT TIE DOWN											
AC04	☐	Transit Tie Down (Unoccupied)								NC	
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 20"-41"								\$365	
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 20"-46"								\$365	
AC35	☐	WC-19 Transit Tie Down without belt ⁴								\$350	

1- Available with HD kit; 2- Only with HD Kit (350 lb);

4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

5- HD fee is required for clients weighing 300-350 lbs

6- See table on page 7; 7- Distance of 2 1/8" from end of seat upholstery to end of frame;

[3. Width calculator](#)



PACKAGES

VP01	☐	Vision Package : Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹	K0108	\$95
------	--------------------------	---	-------	------

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo *VP01; Recommended with black chair;
See picture of the product on page 7.

A6 FRONT Seat-to-floor Height - Please select your choice

Fork/Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12" ☐	13" ☐	14" ☐	12 1/2" ☐	13 1/2" ☐	14 1/2" ☐												
	P2	12 1/2" ☐	13 1/2" ☐	14 1/2" ☐	13" ☐	14" ☐	15" ☐		14 1/2" ☐	15 1/2" ☐									
5"	P1																		
	P2				13" ☐	14" ☐	15" ☐	13 1/2" ☐	14 1/2" ☐	15 1/2" ☐									
	P3	13" ☐	14" ☐	15" ☐	13 1/2" ☐	14 1/2" ☐	15 1/2" ☐	14" ☐	15" ☐	16" ☐		15 1/2" ☐	16 1/2" ☐						
	P4	13 1/2" ☐	14 1/2" ☐	15 1/2" ☐	14" ☐	15" ☐	16" ☐	14 1/2" ☐	15 1/2" ☐	16 1/2" ☐	15" ☐	16" ☐	17" ☐		16 1/2" ☐	17 1/2" ☐			
7"	P1										15 1/2" ☐	16 1/2" ☐	17 1/2" ☐	16" ☐	17" ☐	18" ☐			
	P2							15 1/2" ☐	16 1/2" ☐	17 1/2" ☐	16" ☐	17" ☐	18" ☐	16 1/2" ☐	17 1/2" ☐	18 1/2" ☐			
	P3							16" ☐	17" ☐	18" ☐	16 1/2" ☐	17 1/2" ☐	18 1/2" ☐	17" ☐	18" ☐	19" ☐	17 1/2" ☐	18 1/2" ☐	19 1/2" ☐
	P4	15 1/2" ☐	16 1/2" ☐	17 1/2" ☐	16" ☐	17" ☐	18" ☐	16 1/2" ☐	17 1/2" ☐	18 1/2" ☐	17" ☐	18" ☐	19" ☐	17 1/2" ☐	18 1/2" ☐	19 1/2" ☐	18" ☐	19" ☐	20" ☐
	P5	16" ☐	17" ☐	18" ☐	16 1/2" ☐	17 1/2" ☐	18 1/2" ☐	17" ☐	18" ☐	19" ☐	17 1/2" ☐	18 1/2" ☐	19 1/2" ☐	18" ☐	19" ☐	20" ☐	18 1/2" ☐	19 1/2" ☐	20 1/2" ☐
Froglegs	P1				13 1/2" ☐	14 1/2" ☐	15 1/2" ☐												
	P2	13 1/2" ☐	14 1/2" ☐	15 1/2" ☐	14" ☐	15" ☐	16" ☐	14 1/2" ☐	15 1/2" ☐	16 1/2" ☐									
	P3	14" ☐	15" ☐	16" ☐	14 1/2" ☐	15 1/2" ☐	16 1/2" ☐	15" ☐	16" ☐	17" ☐	15 1/2" ☐	16 1/2" ☐	17 1/2" ☐						

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form); Orange boxes are not available with a seat slope of more than 1.5in.

FRONT CASTERS					Availability		
CASTERS							
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")			3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 1/4")			6", 8"	E2214/E2215	\$165
CA14 to CA16	☐	Newton UltraCaster Composite (1 1/2")			4", 5", 6"		\$95
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 1/2")			4", 5", 6"		\$215
CASTER OPTIONS							
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹				K0073	\$300
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				E1015	\$490

1- Not available with WC-19

Helio A6 - REAR Seat-to-floor Height - Please select your choice

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12" ☐	13" ☐	13 3/4" ¹ ☐	14 1/4" ☐	14 3/4" ☐
12 1/2" ☐	13 1/2" ☐	14 1/4" ☐	14 3/4" ☐	15 1/4" ☐
13" ☐	14" ☐	14 3/4" ☐	15 1/4" ☐	15 3/4" ☐
13 1/2" ☐	14 1/2" ☐	15 1/4" ☐	15 3/4" ☐	16 1/4" ☐
14" ☐	15" ☐	15 3/4" ☐	16 1/4" ☐	16 3/4" ☐
14 1/2" ☐	15 1/2" ☐	16 1/4" ☐	16 3/4" ☐	17 1/4" ☐
15" ☐	16" ☐	16 3/4" ☐	17 1/4" ☐	17 3/4" ☐
15 1/2" ☐	16 1/2" ☐	17 1/4" ☐	17 3/4" ☐	18 1/4" ☐
16" ☐	17" ☐	17 3/4" ☐	18 1/4" ☐	18 3/4" ☐
16 1/2" ☐	17 1/2" ☐	18 1/4" ☐	18 3/4" ☐	19 1/4" ☐
17" ☐	18" ☐	18 3/4" ☐	19 1/4" ☐	19 3/4" ☐
17 1/2" ☐	18 1/2" ☐	19 1/4" ¹ ☐	19 3/4" ☐	20 1/4" ☐

Configuration used: 1" seat slope with a center of gravity at 2"

1- Not available with One Arm Drive

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

REAR WHEELS									
ONE ARM DRIVE ¹ <i>Only available with Newton One rear wheels permanent axles, 0° camber</i>									
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐		E0958	\$1,495		
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐		E0958	\$1,600		
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)					NC		
REAR WHEELS					Available				
RW01	☐	Mag ³			20", 22" 24"		\$57		
RW02	☐	Newton One - spoke wheel ⁴			20", 22" 24"		NC		
RW03	☐	Newton Gravity - ultralight wheel			20", 22", 24", 25", 26"		\$360		
RW04	☐	Spinergy Spox Blk STD ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐			22", 24", 25", 26"		\$915		
RW05	☐	Spinergy LX Blk STD ☐ Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐			22", 24", 25", 26"		\$1,145		
RW07	☐	Spinergy CLX ⁶ Blk STD ☐ Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐			24", 25"		\$1,625		
TIRES <i>(The brand and color may vary depending on parts availability)</i>					Available				
AIRLESS									
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black			20", 22" 24"	E2218	STD		
TI02	☐	Soft Urethane 1 3/8" - Medium Tread, black ⁷			22", 24", 25", 26"	E2218	\$97		
TI13	☐	Soft Urethane 1" - No tread, black - <i>Sentinel</i>			20", 22", 24", 25", 26"	E2218	\$118		
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸			22", 24", 26"	E2211/E2213	\$95		
PNEUMATICS									
TI05	☐	Pneumatic 1 3/8" - Medium Tread			20", 22", 24", 26"	E2211/E2212	NC		
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>			20", 22", 24", 25", 26"	E2211/E2212	\$105		
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>			24", 25", 26"	E2211/E2212	\$150		
TI03	☐	Pneumatic 2" - Large Tread, black ⁹			22", 24", 25"	E2211/E2212	\$115		

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide;
 2- Only standard handrim position available; 3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires;
 5- color spokes require an additional 10 days lead-time; 6- Only available with 1" tires; 7- 22" is only available with Mags, Newton One and Newton Gravity wheels
 8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels; 9- Not available with armrests when 3° camber is selected.

HANDRIMS									
HANDRIM					Available				
HR01	☐	Aluminum anodized			20", 22", 24", 25", 26"		STD		
HR61	☐	Aluminum black hard anodized ²			20", 22", 24", 25"		\$75		
HR02	☐	Plastic coated ¹			20", 22", 24", 25", 26"		\$110		
HR03	☐	High friction coated ^{1,2,3}			22", 24", 25", 26"		\$375		
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ² <i>See page 6 for the color selection and price</i>			20", 22", 24", 25", 26"				
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"		\$450		
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"		\$495		
HR04	☐	Natural Fit ^{1,2,3,4}			22", 24", 25", 26"		\$395		
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)			22", 24", 25", 26"		\$395		
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)			22", 24", 25", 26"		\$445		
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)			22", 24", 25", 26"		\$345		
HR99	☐	Omit handrim					NC		

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;
 4- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

HANDRIMS OPTION									
ASSEMBLY POSITION									
HR29	☐	Standard Handrim position	Full tab				STD		
HR030	☐	Narrow Handrim position no cut ⁶	Full tab				NC		
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions				\$30		
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐ Short tab 1 position ⁶ ☐				\$30		
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS									
HR23	☐	Standard grip					\$95		
HR27	☐	Super grip					\$150		
HR099	☐	Omit thumb piece					NC		


5- Only available with Aluminum and Air Grip handrims; 6- Not available with the combination of 24" Newton Gravity, Spinergy Spox, LX, wheels and a 1 3/8" Airless insert pneumatic tire with any handrim;

AXLE ¹									
AX01	☐	Quick release axle					STD		
AX02	☐	Permanent axle					NC		
AX03	☐	Quad release axle				K0108	\$130		
AXLE PLATE									
AX04	☐	Amputee adjustable axle plate ²				E0959	\$205		

1- Not required with one arm drive; 2- May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down and the swing away anti-tipper

REAR WHEEL CAMBER									
CM00	☐	0° camber					STD		
CM03	☐	3° camber					NC		

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 1" 1/4"

Most stable			Most responsive		
AMPUTEE ¹	3/4"	1" 1/4"	1" 3/4"	2" 1/4"	2" 3/4"
☐	☐	☐	☐	☐	☐

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1- Only available with amputee axle plate

WHEEL LOCK

WL08	☐	Newton Wheel lock - push to lock			STD
WL02	☐	Pull to lock			NC
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)		E0974	\$155
WL13	☐	6" removable extension handle ¹ *VP01		E0961	\$77
WL14	☐				
WL05	☐	Scissors Lock			\$155
WL06	☐	Attendant Lock (additional set of wheel lock mounted to the back)			\$215
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐	\$330
WL99	☐	Omit wheel locks			NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

FOOTREST

SWING-IN/SWING-OUT FOOTREST

FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"	Set at*:		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at*:		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at*:		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"	Set at*:		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at*:		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at*:		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" ⁵	Set at*:		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" ⁵	Set at*:		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" ⁵	Set at*:		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5" - 11"	Set at*:	NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5" - 13"	Set at*:	NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5" - 14.5"	Set at*:	NC

FOOTREST OPTIONS

FR09	☐	High mount (footplate attached to the hanger with clamps) ³	5" - 13.5"	Set at*:	\$85
FR10	☐	Elevating legrest with calf pad, Length adjustable ⁴	Left ☐ Right ☐ 14" - 22"	Set at:	E0990 \$245
FR99	☐	Omit Footrest			NC

1- Available only with two pieces adjustable footplate; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest; * +/- 1/2"

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES

FP01	☐	Standard - Composite ¹			STD
FP02	☐	Newton adjustable angle - Composite ^{*VP01}		K0040	\$215
FP03	☐	Newton adjustable angle - Aluminum		K0040	\$235
FP23	☐	Newton adjustable angle extra deep - Aluminum		K0040	\$235
FP06	☐	One piece adjustable angle flip-up - Aluminum		K0040	\$195

FOOTPLATE ACCESSORIES

FP08	☐	Length adjustable heel loop ²		E0951	\$41
FP09	☐	Calf strap		K0038	\$39
FP30	☐	Bodypoint Aeromesh Padded calf strap		K0038	\$118

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT

FR27	☐	Residual limb support - flat 6"x8"	Left ☐ Right ☐	E1020	\$420
FR28	☐	Residual limb support - flat 8"x10"	Left ☐ Right ☐	E1020	\$420
FR29	☐	Residual limb support - flat 10"x10"	Left ☐ Right ☐	E1020	\$420
FR30	☐	Residual limb support - flat 10"x14"	Left ☐ Right ☐	E1020	\$420
FR32	☐	Residual limb support - contour 6"x8"	Left ☐ Right ☐	E1020	\$420
FR33	☐	Residual limb support - contour 8"x10"	Left ☐ Right ☐	E1020	\$420

SEATING AND BACKREST									
BACK UPHOLSTERY									
BH05	☐	Slip-on Nylon back upholstery							STD
BH06	☐	Adjustable tension Nylon back upholstery					E2611/E2612		\$325
BH99	☐	Omit Back upholstery							NC
	☐	NXT back Support					See price list		
SEAT CUSHION									
SC01	☐	2" Cushion							\$170
SC02	☐	3" Cushion					E2601/E2602		\$190
SC99	☐	Omit seat cushion							STD
	☐	NXT cushion					See price list		
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE (push handle included - no further selection needed)									
BC58-915	☐	Straight aluminum back cane	9" - 15"		increment of 1/2"	Set at ² :			NC
BC58-1521	☐	Straight aluminum back cane	15" - 21"		increment of 1/2"	Set at ² :			NC
BC60	☐	Straight aluminum back cane	21" - 24"		increment of 1"	Set at ² :			\$77
BC01-16	☐	8° bend aluminum back cane	16"						NC
BC01-18	☐	8° bend aluminum back cane	18"						NC
BC01-20	☐	8° bend aluminum back cane	20"						NC
BC59	☐	8° bend aluminum back cane	21" - 24"		increment of 1"	Set at ² :			\$77
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)									
BC11-912	☐	Straight aluminum back cane	9" - 12"		increment of 1/2"	Set at ² :			NC
BC11-1215	☐	Straight aluminum back cane	12" - 15"		increment of 1/2"	Set at ² :			NC
BC11-1518	☐	Straight aluminum back cane	15" - 18"		increment of 1/2"	Set at ² :			NC
BC11-1821	☐	Straight aluminum back cane	18" - 21"		increment of 1/2"	Set at ² :			NC
BC12-1315	☐	8° bend aluminum back cane	13.5" - 15"		increment of 1/2"	Set at ² :			\$60
BC12-1518	☐	8° bend aluminum back cane	15" - 18"		increment of 1/2"	Set at ² :			\$60
BC12-1821	☐	8° bend aluminum back cane	18" - 21"		increment of 1/2"	Set at ² :			\$60
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)									
BC13-1518	☐	Straight aluminum back cane	15" - 18"		increment of 1/2"	Set at ² :			\$215
BC13-1821	☐	Straight aluminum back cane	18" - 21"		increment of 1/2"	Set at ² :			\$215
BC14-1518	☐	8° bend aluminum back cane	15" - 18"		increment of 1/2"	Set at ² :			\$270
BC14-1821	☐	8° bend aluminum back cane	18" - 21"		increment of 1/2"	Set at ² :			\$270
		Select back cane angle:	85° ☐	90° STD ☐	95° ☐	100° ☐	105° ☐	110° ☐	
PUSH HANDLES OPTIONS (ONE option at a time)									
PH01	☐	Push handles							STD
PH05	☐	Fold down push handles ^{3, 5}							\$145
PH06-920	☐	Clamp on height adjustable push handles 9"-20" ^{3, 4, 5}							\$325
PH06-923	☐	Clamp on height adjustable push handles 9"-23" ^{3, 4, 5}							\$325
PH02	☐	Omit push handles							NC
STROLLER HANDLE AND STABILIZER BAR (only available with push handle)									
PH04	☐	Stroller handle (not available with the HD option)					K0108		\$165
PH03	☐	Newton folding stabilizer bar (not available with stroller handle)							\$270

1 - Not compatible with 8 degree bend back cane 2 - If not indicated, default height setting is at the lowest position - measured from seat sling to top of back upholstery; 3 - 265 lb weight limit
4 - Not compatible with swing away armrests, 8° bend back canes, 9-12" and 21-24" back canes and amputee axle plates; 5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

ARMRESTS									
HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ^{1, 4} *VP01									
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		E0973		\$355
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		E0973		\$355
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		E0973		\$355
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		E0973		\$355
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01									
AR13-610	☐	10" Desk length armpad ⁵	6" - 10"				E0973		\$330
AR13-812	☐	10" Desk length armpad	8" - 12"				E0973		\$330
AR13-1014	☐	10" Desk length armpad	10" - 14"				E0973		\$330
AR14-610	☐	14" Full length armpad ⁵	6" - 10"				E0973		\$330
AR14-812	☐	14" Full length armpad	8" - 12"				E0973		\$330
AR14-1014	☐	14" Full length armpad	10" - 14"				E0973		\$330
SWING AWAY ARMREST ²									
AR05	☐	Tubular armrest without plastic side guards - Short pad	8" - 11"				E0973		\$230
AR19	☐	Tubular armrest with plastic side guards - Short pad	8" - 11"				E0973		\$395
AR23	☐	Tubular armrest without plastic side guards - Long pad	8" - 11"				E0973		\$230
AR24	☐	Tubular armrest with plastic side guards - Long pad	8" - 11"				E0973		\$395
AR06	☐	10" Desk length flat armpad without plastic side guards	8" - 11"				E0973		\$230
AR20	☐	10" Desk length flat armpad with plastic side guards	8" - 11"				E0973		\$395
AR17	☐	14" Full length flat armpad without plastic side guards	8" - 11"				E0973		\$230
AR21	☐	14" Full length flat armpad with plastic side guards	8" - 11"				E0973		\$395
ARMREST OPTION									
AR12	☐	14" Full Length armpad - Gel Ovation (pair) ³							\$165
AR99	☐	Omit armrest							NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;
3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

ACCESSORIES									
SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>									
SG01	☐	Plastic side guards <i>(removable)</i>	Straight bracket - STD ☐ Offset bracket ¹ ☐					K0108	\$200
SG02	☐	Carbon fiber side guards <i>(fixed)</i>							\$270
SG03	☐	Carbon fiber side guards <i>(removable)</i>	Straight bracket - STD ☐ Offset bracket ¹ ☐					K0108	\$305
SG06	☐	Plastic fender side guards <i>(removable)</i>	Straight bracket - STD ☐ Offset bracket ¹ ☐						\$290
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>	Straight bracket - STD ☐ Offset bracket ¹ ☐					K0108	\$550
SG99	☐	Omit side guard							NC

1- only compatible with 0" camber; increase width by 1/2" per side for a total of 1" ;

POSITIONING BELT									
SB01-S	☐	Velcro adjustable - Small 35" ¹						E0978	\$46
SB01-L	☐	Velcro adjustable - Long 60"						E0978	\$46
SB02-S	☐	Auto buckle - Small 35" ¹						E0978	\$87
SB02-M	☐	Auto buckle - Medium 48"						E0978	\$87
SB02-L	☐	Auto buckle - Long 60"						E0978	\$87
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"						E0978	\$85
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"						E0978	\$85
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"						E0978	\$180
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"						E0978	\$180
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"						E0978	\$180
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"						E0978	\$195
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"						E0978	\$195
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"						E0978	\$195
SB03	☐	Hardware attachment for 2 extra points (1 pair)							\$17
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"						K0108	\$305
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"						K0108	\$305
SB99	☐	Omit positioning belt							STD

1- Not recommended for seat width of 18" and more

HEADREST									
AC01-S	☐	Profiled Headrest - Small (3.5"X6")						E0955	\$185
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")						E0955	\$185
AC01-L	☐	Profiled Headrest - Large (5.5"X8")						E0955	\$185
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar ¹						E1028	\$300

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES									
AC03	☐	Anti-tippers						E0971	\$150
AC29	☐	Swing-away anti-tipper ¹	Left ☐	Right ☐					\$190
AC27	☐	Tip assist ²							\$67
AC06	☐	Cane and crutch holder	Left ☐	Right ☐				E2207	\$200
AC07	☐	Oxygen tank holder ⁴ - D and E Type - <i>Installed to the left</i>						E2208	\$315
AC08	☐	Spoke guards ³	Available: 20", 22", 24", 25", 26"					K0065	\$118
AC09	☐	Footrest neoprene impact guard <i>(pair)</i> - Black						K0108	\$92
AC10	☐	Seat pouch							\$82
AC30	☐	Newton multi-purpose pouch							\$61
AC26	☐	BodyPoint Mobility Bag							\$135
AC23	☐	Newton Backpack							\$130
AC12	☐	Velcro Closing Strap							\$57
AC13	☐	Tool kit							\$51
AC15	☐	Military Patches	Army ☐	Navy ☐	Air Force ☐	Marines ☐	Coast Guard ☐		
		Combat Wounded Purple Heart ☐	Space Force ☐	American Flag ☐					\$96

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate

2- Cannot be combined with Anti Tippers; 3- Not available with Mag wheels and Spinerger; the 24" spoke guard will be installed on the 25" wheel ;


4- Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down

Colors																			
Options ¹	Black	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
AirGrip	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	\$110																		
Frame ²	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
	STD																		
FC99	☐	Custom Color	Specify PPG Paint color:														<i>additional 10 days lead-time</i>		\$390

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- with black hydroformed aluminum crossbrace 3- The frame color will be confirmed with the quote conversion into an order (in absence of selection, the color will be black)

Fork and Stem bolt combinations



Frame length chart

	Seat Depth						
Frame length	14"	15"	16"	17"	18"	19"	20"
STD							
+1							N/A
+2						N/A	N/A
-1	N/A						

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater

Rear wheel: shall be 24 in and less

"T" or flip-back armrest only compatible with 0° camber

One arm drive: not compatible with "T"

One arm drive and flip-back armrest only available with 17 in depth frame

"T" or flip-back armrest may require to move the rear wheels outwards

Not compatible with wheel lock lever extension

Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.25 to 2.75
	16	0.75 to 2.25
22 in	17+	0.75 to 2.75
	16	0.75 to 1.75
24 in	17	0.75 to 2.25
	18+	0.75 to 2.75

CUSTOMIZATION STARTS HERE!



With EVO Custom, whether it's subtle tweaks for comfort and mobility improvements or custom colors that make your accessories uniquely yours, we're here to support you.



Explore the
EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.

Notes :



FOLDING ON CHAIR – USA ORDER FORM



Use wheelchair configuration
Or use this size:
NXCAF
(W" x D")
HPCPS CODE E2607-E2608

DEPTH	WIDTH											
	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
12"	\$468	\$468	\$468		N/A	N/A	N/A	N/A	N/A	N/A	N/A	
13"						N/A	N/A	N/A	N/A	N/A	N/A	
14"	\$468	\$468	\$468									
15"		\$468	\$468	\$468	\$468							
16"		\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	
17"	N/A			\$468	\$468	\$468	\$468	\$468	\$468			
18"	N/A			\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	
19"	N/A			\$468	\$468	\$468	\$468	\$468	\$468			
20"	N/A			\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	

Option: Change cover for incontinent/infection control outer cover ☐ ☐ N/C

Green boxes are stocked sizes; the red boxes require customization*; the N/A boxes are not available sizes.
*Customization (red boxes) is available with an upcharge of \$160.00 and has a 2 to 4 weeks additional lead time.

SEATING AND CUSHION ACCESSORIES

MODEL #	Description	Fits Back Width	Retail	QTY
RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	\$80	



Optima™
THORACIC BACK SUPPORT
HPCPS CODE: E2620-E2621

WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOACB1413	\$730	☐	NIB-2193V	\$130	☐
14"	16"	NOACB1416	\$730	☐	NIB-2194V	\$130	☐
15"	13"	NOACB1513	\$730	☐	NIB-2197V	\$130	☐
15"	16"	NOACB1516	\$730	☐	NIB-2198V	\$130	☐
16"	13"	NOACB1613	\$730	☐	NIB-2201V	\$130	☐
16"	16"	NOACB1616	\$730	☐	NIB-2202V	\$130	☐
16"	19"	NOACB1619	\$730	☐	NIB-2203V	\$130	☐
17"	13"	NOACB1713	\$730	☐	NIB-2206V	\$130	☐
17"	16"	NOACB1716	\$730	☐	NIB-2207V	\$130	☐
18"	13"	NOACB1813	\$730	☐	NIB-2208V	\$130	☐
18"	16"	NOACB1816	\$730	☐	NIB-2211V	\$130	☐
18"	19"	NOACB1819	\$730	☐	NIB-2212V	\$130	☐
19"	16"	NOACB1916	\$730	☐	NIB-2213V	\$130	☐
19"	19"	NOACB1919	\$730	☐	NIB-2215V	\$130	☐
20"	13"	NOACB2013	\$730	☐	NIB-2218V	\$130	☐
20"	16"	NOACB2016	\$730	☐	NIB-2219V	\$130	☐
20"	19"	NOACB2019	\$730	☐	NIB-2220V	\$130	☐
21"	16"	NOACB2116	\$730	☐	NIB-2221V	\$130	☐
21"	19"	NOACB2119	\$730	☐	NIB-2222V	\$130	☐
22"	13"	NOACB2213	\$783	☐	NIB-2223V	\$130	☐
22"	16"	NOACB2216	\$783	☐	NIB-2224V	\$130	☐
22"	19"	NOACB2219	\$783	☐	NIB-2225V	\$130	☐

PLEASE NOTE: Each back will fit widths up to 2" wider.



Optima Carbon™
THORACIC BACK SUPPORT
CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)
HPCPS CODE: E2620-E2621

WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOCFACB1413	\$972	☐	NIB-2193V	\$130	☐
14"	16"	NOCFACB1416	\$1,008	☐	NIB-2194V	\$130	☐
15"	13"	NOCFACB1513	\$972	☐	NIB-2197V	\$130	☐
15"	16"	NOCFACB1516	\$1,008	☐	NIB-2198V	\$130	☐
16"	13"	NOCFACB1613	\$972	☐	NIB-2201V	\$130	☐
16"	16"	NOCFACB1616	\$1,008	☐	NIB-2202V	\$130	☐
17"	13"	NOCFACB1713	\$972	☐	NIB-2206V	\$130	☐
17"	16"	NOCFACB1716	\$1,008	☐	NIB-2207V	\$130	☐
18"	13"	NOCFACB1813	\$972	☐	NIB-2211V	\$130	☐
18"	16"	NOCFACB1816	\$1,008	☐	NIB-2212V	\$130	☐
20"	13"	NOCFACB2013	\$972	☐	NIB-2218V	\$130	☐
20"	16"	NOCFACB2016	\$1,008	☐	NIB-2219V	\$130	☐

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

THORACIC SHAPE			Optima™ DEEP THORACIC BACK SUPPORT HCPCS CODE: E2620-E2621					
	WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
	14"	13"	NOPCB1413	\$730	☐	NIB-2246V	\$130	☐
	14"	16"	NOPCB1416	\$730	☐	NIB-2247V	\$130	☐
	15"	13"	NOPCB1513	\$730	☐	NIB-2250V	\$130	☐
	15"	16"	NOPCB1516	\$730	☐	NIB-2251V	\$130	☐
	16"	13"	NOPCB1613	\$730	☐	NIB-2254V	\$130	☐
	16"	16"	NOPCB1616	\$730	☐	NIB-2255V	\$130	☐
	16"	19"	NOPCB1619	\$730	☐	NIB-2256V	\$130	☐
	17"	13"	NOPCB1713	\$730	☐	NIB-2259V	\$130	☐
	17"	16"	NOPCB1716	\$730	☐	NIB-2260V	\$130	☐
	18"	13"	NOPCB1813	\$730	☐	NIB-2264V	\$130	☐
	18"	16"	NOPCB1816	\$730	☐	NIB-2265V	\$130	☐
	18"	19"	NOPCB1819	\$730	☐	NIB-2266V	\$130	☐
	19"	16"	NOPCB1916	\$730	☐	NIB-2267V	\$130	☐
	19"	19"	NOPCB1919	\$730	☐	NIB-2268V	\$130	☐
	20"	13"	NOPCB2013	\$730	☐	NIB-2271V	\$130	☐
	20"	16"	NOPCB2016	\$730	☐	NIB-2272V	\$130	☐
	20"	19"	NOPCB2019	\$730	☐	NIB-2273V	\$130	☐
	21"	16"	NOPCB2116	\$730	☐	NIB-2274V	\$130	☐
	21"	19"	NOPCB2119	\$730	☐	NIB-2275V	\$130	☐
	22"	13"	NOPCB2213	\$783	☐	NIB-2276V	\$130	☐
	22"	16"	NOPCB2216	\$783	☐	NIB-2277V	\$130	☐
	22"	19"	NOPCB2219	\$783	☐	NIB-2278V	\$130	☐

		Optima Carbon™ DEEP THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb) HCPCS CODE: E2620-E2621					
WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOCFPCB1413	\$1,066	☐	NIB-2246V	\$130	☐
14"	16"	NOCFPCB1416	\$1,108	☐	NIB-2247V	\$130	☐
15"	13"	NOCFPCB1513	\$1,066	☐	NIB-2250V	\$130	☐
15"	16"	NOCFPCB1516	\$1,108	☐	NIB-2251V	\$130	☐
16"	13"	NOCFPCB1613	\$1,066	☐	NIB-2254V	\$130	☐
16"	16"	NOCFPCB1616	\$1,108	☐	NIB-2255V	\$130	☐
17"	13"	NOCFPCB1713	\$1,066	☐	NIB-2259V	\$130	☐
17"	16"	NOCFPCB1716	\$1,108	☐	NIB-2260V	\$130	☐
18"	13"	NOCFPCB1813	\$1,066	☐	NIB-2264V	\$130	☐
18"	16"	NOCFPCB1816	\$1,108	☐	NIB-2265V	\$130	☐
20"	13"	NOCFPCB2013	\$1,066	☐	NIB-2271V	\$130	☐
20"	16"	NOCFPCB2016	\$1,108	☐	NIB-2272V	\$130	☐

THORACIC SHAPE			Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT HCPCS CODE E2620-E2621					
	xtend™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
	NHAACB(L)1213/17	12"	13-17"	\$730	☐	NIB-2148V	\$130	☐
	NHAACB(L)1413/17	14"	13-17"	\$730	☐	NIB-2149V	\$130	☐
	NHAACB(L)1613/17	16"	13-17"	\$730	☐	NIB-2150V	\$130	☐
	NHAACB(L)1813/17	18"	13-17"	\$730	☐	NIB-2151V	\$130	☐
	NHAACB (L)2013/17	20"	13-17"	\$730	☐	NIB-2152V	\$130	☐
	NHAACB (L)2213/17	22"	13-17"	\$783	☐	NIB-2153V	\$130	☐

		Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT HCPCS CODE E2620-E2621					
xtend™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NHAPCB1616/20	16"	16"-20"	\$730	☐	NIB-2143V	\$130	☐
NHAPCB1816/20	18"	16"-20"	\$730	☐	NIB-2144V	\$130	☐
NHAPCB2016/20	20"	16"-20"	\$730	☐	NIB-2145V	\$130	☐
NHAPCB2216/20	22"	16"-20"	\$783	☐	NIB-2146V	\$130	☐

PELVIC SHAPE			Classic™ PELVIC BACK SUPPORT HCPCS CODE E2620-E2621					
	Classic™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
	NACB1616	16"	16"	\$730	☐	NIB-2163V	\$130	☐
	NACB1618	16"	18"	\$730	☐	NIB-2164V	\$130	☐
	NACB1816	18"	16"	\$730	☐	NIB-2166V	\$130	☐
	NACB1818	18"	18"	\$730	☐	NIB-2167V	\$130	☐
	NACB1820	18"	20"	\$730	☐	NIB-2168V	\$130	☐
	NACB2020	20"	20"	\$730	☐	NIB-2171V	\$130	☐
	NACB2216	22"	16"	\$783	☐	NIB-2172V	\$130	☐
	NACB2218	22"	18"	\$783	☐	NIB-2173V	\$130	☐
	NACB2220	22"	20"	\$783	☐	NIB-2174V	\$130	☐

		Classic™ DEEP PELVIC BACK SUPPORT HCPCS CODE E2620-E2621					
Classic™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NPCB1616	16"	16"	\$730	☐	NIB-2175V	\$130	☐
NPCB1618	16"	18"	\$730	☐	NIB-2176V	\$130	☐
NPCB1818	18"	18"	\$730	☐	NIB-2179V	\$130	☐
NPCB1820	18"	20"	\$730	☐	NIB-2180V	\$130	☐
NPCB2020	18"	20"	\$730	☐	NIB-2183V	\$130	☐
NPCB2216	22"	16"	\$783	☐	NIB-2184V	\$130	☐
NPCB2218	22"	18"	\$783	☐	NIB-2185V	\$130	☐
NPCB2220	22"	20"	\$783	☐	NIB-2186V	\$130	☐

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

BACK SUPPORT ACCESSORIES

	MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	\$205	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	\$70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	\$70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	\$70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	\$70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	\$70	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

HCPCS CODE: E0955

	MODEL	Description	RETAIL	✓
	NIB-0784	10" x 5" Large Rectangular	\$205	
	NIB-0818	14" x 5" Extra-Large Rectangular	\$205	

MultiFit™ HEAD SUPPORT – FIXTURES

HCPCS CODE HARDWARE: E1033

	MODEL	Description	RETAIL	✓
	NIB-0747	Multi-Axis Hardware	\$150	
	NIB-0823	Head Support Receiver Bracket	\$60	

EXTRA COVERS

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	NXCCBF_____V	\$130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	NXCCBF_____V	\$150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	NXCCAF_____V	\$130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	NXCCAF_____V	\$150	

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com