



Helio A7 - FiveDay

ORDER FORM - CANADA

P.O. number: * _____ Date: _____
Company Name: * _____ Phone: * _____
Ordered by: * _____
Clinician/Sales rep: _____
Email: _____
Marked for: _____

Ship to: * _____
Name: _____
Address: _____
City: _____ Prov: _____
Postal Code: _____
Phone: _____

* Required MA080 ☐ English Manual (STD) MA081 ☐ Manuel Français Taxes not included

MODEL											
1	WC02	Ultralight Helio A7 wheelchair 7000 series aluminum (265 lb weight capacity)						WM0001020	\$2 645	☐	
2	WC09	HD Heavy duty kit (350 lb weight capacity)						WAMF	\$775	☐	
3	WC10	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber							\$30		
SEAT WIDTH ³											
4	SW	Seat width:	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐			NC	☐
5	SW	Seat width:	19" ¹ ☐	20" ¹ ☐	21" ² ☐	22" ² ☐		WAMA	\$295	☐	
SEAT DEPTH											
6	SD	Seat depth:	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐			NC	☐
7	SD	Seat depth:	19" ☐	20" ☐				WAMB	\$225	☐	

1 - Available with HD option

2 - Only available with HD option

[3. Width calculator](#)

A7 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13	14	15	13 1/2	14 1/2	15 1/2												
	P2	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2									
5"	P1																		
	P2				14	15	16	14 1/2	15 1/2	16 1/2									
	P3	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2						
	P4	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
7"	P1										16 1/2	17 1/2	18 1/2	17	18	19			
	P2							16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2			
	P3							17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2
	P4	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21
	P5	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21	19 1/2	20 1/2	21 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
8	CA01	Composite wheel w/ PU tire (1" width)	3" ☐	4" ☐	5" ☐	6" ☐	7" ☐	8" ☐			NC	☐
9	CA10	Pneumatic casters 6" x 1 ¼"									\$160	☐
10	CA11	Pneumatic casters 8" x 1 ¼"									\$160	☐
11	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "								\$98	☐
12	CA15		5" x 1 ½ "								\$98	☐
13	CA16		6" x 1 ½ "								\$98	☐
14	CA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½"								\$175	☐
15	CA18		5" x 1 ½"								\$185	☐
16	CA19		6" x 1 ½"								\$195	☐
CASTER OPTIONS												
17	CA12	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY)							WAMS	\$255	☐	
18	CA13	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)								\$489	☐	

A7 - REAR Seat-to-floor Height

WHEEL SIZES	Height available	Specify desired height:
20" wheel ☐	12 1/2 to 17 1/4	
22" wheel ☐	13 1/2 to 18 1/4	
24" wheel ☐	14 1/2 to 19 1/4	
25" wheel ☐	15 to 19 3/4	
26" wheel ☐	15 1/2 to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest											
SWING-IN/SWING-OUT FOOTREST											
19	FR03	60° Swing-in/swing-out footrest	13.5 - 17" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐	Set at:*			NC	☐	
20	FR04	70° Swing-in/swing-out footrest	13.5 - 17" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐	Set at:*			NC	☐	
21	FR15	80° Swing-in/swing-out footrest	13.5 - 16" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐	Set at:*			NC	☐	
90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT <i>(Footplate mounted to hanger tube)</i>											
22	FR06	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5 - 11"			Set at:*			NC	☐	
23	FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5 - 13"			Set at:*			NC	☐	
24	FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5 - 14.5"			Set at:*			NC	☐	
FOOTREST OPTIONS ³											
25	FR09	High mount <i>(footplate attached to the hanger with clamps)</i>	5 - 13.5" ☐			Set at:*			\$90	☐	
26	FR99	Omit Footrest							NC	☐	

1 - Available only with adjustable footplates

2 - May interfere with front casters

3- Not available with front mount

* +/- 1/2"

ELEVATING LEGREST (unit)			The elevating legrest replace the footrest side selected					
27	FR10	Elevating legrest with calf pad, 14" - 22" Length adjustable	Left ☐	Right ☐	Set at:	WAMK	\$290	☐

FOOTPLATES						
28	FP01	Standard - Composite			STD	☐
29	FP02	Newton adjustable angle - Composite		WAMJ	\$145	☐
30	FP03	Newton adjustable angle - Aluminum		WAMJ	\$165	☐
31	FP04	Newton adjustable angle - Carbon fiber <i>(with rubber inserts)</i>		WAMJ	\$380	☐
32	FP06	One piece adjustable angle flip-up - Aluminum	WAMJ	SEFND0015	\$206	☐
33	FP07	One piece adjustable angle flip-up - Carbon fiber <i>(with rubber inserts)</i>	WAMJ	SEFND0015	\$450	☐
FOOTPLATE ACCESSORIES						
34	FP08	Length adjustable heel loop <i>(not available with one piece footplate)</i>			NC	☐
35	FP09	Calf strap		SEFND0010	\$52	☐
36	FP30	Bodypoint Aeromesh Padded calf strap		SEFND0010	\$129	☐

RESIDUAL LIMB SUPPORT FLAT CUSHION <i>(Swing away, depth, width, height and angle adjustable) - (unit)</i>								
37	FR27	Residual limb support swing away - flat 6"x8"	Left ☐	Right ☐	New	SEFND0005	\$475	☐
38	FR28	Residual limb support swing away - flat 8"x10"	Left ☐	Right ☐	New	SEFND0005	\$475	☐
39	FR29	Residual limb support swing away - flat 10"x10"	Left ☐	Right ☐	New	SEFND0005	\$475	☐
40	FR30	Residual limb support swing away - flat 10"x14"	Left ☐	Right ☐	New	SEFND0005	\$475	☐
RESIDUAL LIMB SUPPORT CONTOUR CUSHION <i>(Swing away, depth, width, height and angle adjustable) - (unit)</i>								
41	FR32	Residual limb support swing away - contour 6"x8"	Left ☐	Right ☐	New	SEFND0005	\$475	☐
42	FR33	Residual limb support swing away - contour 8"x10"	Left ☐	Right ☐	New	SEFND0005	\$475	☐
Back canes with integrated fixed height push handle								
FIXED ANGLE BACK								
43	BC01	8° bend aluminum back cane	16" ☐	18" ☐	20" STD ☐		NC	☐
Back canes, height adjustable, with option								
FIXED ANGLE BACK								
44	BC11	Straight aluminum back cane	9-12" ☐	12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐	\$75	☐
45	BC12	8° bend aluminum back cane	13.5 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐		\$75	☐
ADJUSTABLE ANGLE BACK								
46	BC13	Adjustable angle, straight aluminum back cane	15 - 18" ☐	18 - 21" STD ☐			\$259	☐
47	BC14	Adjustable angle, 8° bend aluminum back cane	15 - 18" ☐	18 - 21" STD ☐			\$309	☐
48		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐
Set height to: ²		<i>(measured from seat sling to top of back upholstery)</i>						
PUSH HANDLE OPTION								
49	PH01	Push handles					STD	☐
50	PH05	Fold down push handles ^{3,5}					\$150	☐
51	PH06	Clamp on height adjustable push handle ^{3,4,5}	9 - 20" STD ☐	9-23" ☐			\$325	☐
52	PH02	Without push handles					NC	☐
1 - Not available with a 8 degree bend cane 2 - If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit								
4 - Not compatible with swing away armrests, 8° bend back canes and amputee axle plates.								
5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock								
STROLLER HANDLE AND STABILIZER BAR ¹ <i>(only available with push handle)</i>								
53	PH04	Stroller handle <i>(not available with the HD option)</i>				WAMW	\$210	☐
54	PH03	Newton folding stabilizer bar					\$258	☐
1 - Stroller handle and stabilizer bar can't be combined								
BACK UPHOLSTERY								
55	BH07	Slip-on Active-Mesh back upholstery					NC	☐
56	BH08	Adjustable tension Active-Mesh back upholstery				WAMD/WAME	\$245	☐
57	BH99	Omit Back upholstery					STD	☐
RIGID BACK								
58		NXT back series				See pricelist p. 7		☐
SEAT CUSHION								
59	SC01	2" Cushion					\$175	☐
60	SC02	3" Cushion					\$195	☐
61	SC99	Omit seat cushion					STD	☐
62		NXT FIT cushion series				See pricelist p. 7		☐

Rear wheels													
ONE ARM DRIVE ¹ <i>Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims</i>													
63	HR22	One arm drive ^{2, 3}	24"	left ☐	right ☐						\$1 395	☐	
64	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)									NC	☐	
REAR WHEELS													
65	RW01	Mag ⁴	20" ☐	22" ☐	24" ☐						\$129	☐	
66	RW02	Newton One - spoke wheel ⁵	20" ☐	22" ☐	24" ☐						NC	☐	
67	RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐				\$360	☐	
68	RW04	Spinergy Spox		22" ☐	24" ☐	25" ☐	26" ☐				\$800	☐	
69		Spinergy Spox spokes colors ⁶	Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐				NC	☐	
70	RW05	Spinergy LX		22" ☐	24" ☐	25" ☐	26" ☐				\$995	☐	
71		Spinergy LX spokes colors ⁶	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐	NC	☐
72	RW07	Spinergy CLX Carbon Blade			24" ☐	25" ☐					\$1 540	☐	
73		Spinergy CLX spokes colors ⁶	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐	NC	☐
TIRES							Available						
74	TI01	Soft Urethane 1 1/4 "				20", 22" 24"					STD	☐	
75	TI02	Soft Urethane knobby 1 3/8 "				22", 24", 25", 26"					\$110	☐	
76	TI04	Soft Urethane - Shox 1"				22", 24", 25", 26"					\$130	☐	
77	TI05	Pneumatic 1 3/8 "				20", 22", 24", 25", 26"					NC	☐	
78	TI08	Pneumatic with Airless insert 1 3/8"				22", 24"					\$215	☐	
79	TI06	Pneumatic HP - 1"				20", 22", 24", 25", 26"					\$110	☐	
80	TI07	Pneumatic HP - Schwalbe Marathon Plus 1"				24", 25", 26"					\$190	☐	

1- Offered with handrims on each side of the wheelchair

2- Only standard handrim position available

3 - Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position

5 - Not compatible with Shox tires

6 - Colour spokes require an additional 10 days lead-time

Handrims										
HANDRIMS										
81	HR01	Aluminum anodized						STD	☐	
82	HR02	Plastic coated ³					WAMO	\$160	☐	
83	HR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}					WAMO	\$430	☐	
84	HR07	Newton Air Grip ¹ Black ☐					WAMO	\$150	☐	
85	HR08	Newton Air Grip ¹ (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐ (11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐ (15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐				Additional 5 day lead-time	WAMO	\$330	☐	
86	HR28	Surge ^{1, 2, 3} (1-1/2" Oval handrim w/ Gripton high friction strip)							\$470	☐
87	HR05	Surge LT ^{1, 2, 3} (1-1/8" Oval handrim w/ Gripton high friction strip)							\$500	☐
88	HR04	Natural Fit Standard Grip ^{1, 2, 3} (1-1/2" Oval handrim w/standard thumb piece)						\$500	☐	
89	HR99	Omit handrim						NC	☐	

1 - Not available with Mag wheels or OAD

2 - Not available in 20"

3- Only standard handrim position available

Handrims option									
Assembly position for aluminum and Newton AirGrip handrims									
90	HR29	Standard Handrim position	Full tab					STD	☐
91	HR30	Narrow Handrim position	Short tab 2 positions					\$25	☐
92	HR31	Super Narrow Handrim	Short tab 2 positions STD ☐	Short tab 1 position ☐				\$25	☐
Thumb piece ⁴ selection Natural Fit handrims									
93	HR23	Standard grip						\$190	☐
94	HR27	Super grip						\$220	☐
95	HR099	Omit thumb piece						NC	☐

4- Compatible with Surge, Surge LT handrims

AXLE									
96	AX01	Quick release axle					WAMU	STD	☐
97	AX02	Permanent axle						NC	☐
98	AX03	Quad release axle					WAMU	\$170	☐
AXLE PLATE									
99	AX04	Amputee adjustable axle plate ¹					WAMT	\$250	☐

1 -May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER							
100	CM00	0° camber				STD	☐
101	CM03	3° camber				NC	☐
102	CM06	6° camber <i>(not available with amputee axle plate)</i>				\$80	☐

Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable								Most responsive			
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK							
103	WL08	Newton Wheel lock - push to lock				STD	☐
104	WL02	Pull to lock				NC	☐
105	WL04	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>			WAMR	\$175	☐
106	WL13 WL14	6" removable extension handle				\$67	☐
107	WL05	Scissor Lock				\$175	☐
108	WL06	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>				\$230	☐
109	WL99	Omit wheel locks				NC	☐

Armrests							
FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ^{1,4}			Height Adjustment		Removable function		
110	AR01	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐ STD	\$380 ☐
111	AR02	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐ STD	\$380 ☐
"T" ARMREST ¹			Height Adjustment				
112	AR03	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐			\$360 ☐
113	AR04	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐			\$360 ☐
SWING AWAY ARMREST ²			Height Adjustment				
114	AR05	Tubular arm pad	8 - 11"				\$295 ☐
115	AR06	10" Desk length arm pad	8 - 11"				\$295 ☐
GEL ARM PAD ³							
116	AR12	14" Full Length Arm Pad - Gel Ovation				SEAND0030	\$150 ☐
117	AR99	Omit armrest					NC ☐

1 - May have a 1/2" adjustment restriction on both limits of the range 2 - Could interfere with some rigid back rests 3 - Not available with swing away armrest 4 - Non available with attendant lock

SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>							
118	SG01	Plastic side guards <i>(removable)</i>			WAMV	\$225	☐
119	SG02	Carbon fiber side guards <i>(fixed)</i> ¹			WAMV	\$255	☐
120	SG03	Carbon fiber side guards <i>(removable)</i>			WAMV	\$295	☐
121	SG06	Plastic fender side guards <i>(removable)</i> ²			WAMV	\$320	☐
122	SG04	Carbon fiber fender side guards <i>(removable)</i>			WAMV	\$495	☐
123	SG99	Omit side guard				NC	☐
SIDE GUARDS OPTION							
124	SG07	Side guards clip <i>(pair)</i> ³			New		NC ☐

1- Not available with back upholstery 2- Not compatible with 20" and 22" rear wheels 3- Not available with fixed carbon fiber side guard

POSITIONING BELT							
125	SB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐		\$45	☐
126	SB02	Auto buckle	Standard 35" ☐	Long 60" ☐		\$55	☐
127	SB07	Hardware attachment for 2 point seatbelt (pair)				NC	☐
128	SB03	Additional hardware for 2 extra attachment points <i>(needed for the 4 point seat belt) (pair)</i>				\$13	☐
129	SB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" ☐	Long 61" ☐		\$82	☐
130	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	SERND0030	\$170 ☐
131	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	SERND0020	\$195 ☐
132	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" ☐	Long 24" ☐		\$324	☐
133	SB99	Omit positioning belt				STD	☐

HEADREST ¹							
134	AC01	Profiled Headrest	Small (3.5"X6") ☐	Medium (4.5"X7") ☐	Large (5.5"X8") ☐	\$205	☐
135	AC28	Headrest mounting hardware with removable non folding stabilizer bar				\$338	☐

1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES										
136	AC03	Anti tippers							\$115	☐
137	AC27	Tip assist ¹							\$60	☐
138	AC04	Transit Tie Down <i>(Unoccupied)</i>							NC	☐
139	AC05	WC-19 Transit Tie Down w/Q'straint belt ²		Short belt 20-41" ☐	Long belt 20-46" ☐				\$345	☐
140	AC06	Cane carrier		Left ☐	Right ☐				\$105	☐
141	AC07	Oxygen tank holder - E Type <i>Not compatible with swing away armrest or amputee axle plate</i>			Left			WAMX	\$375	☐
142	AC08	Spoke guards ³		Available: 20", 22", 24", 26"				WAMP	\$145	☐
143	AC09	Footrest neoprene impact guard (pair) - Black						SEFND0035	\$80	☐
144	AC10	Seat pouch							\$95	☐
145	AC26	BodyPoint Mobility Bag							\$135	☐
146	AC23	Newton Backpack							\$120	☐
147	AC12	Velcro Closing Strap							\$50	☐
148	AC13	Tool kit							\$50	☐
149	AC16	Motion Composites Polo Shirt	W/S ☐	W/M ☐	W/L ☐				\$32	☐
			M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	M/XXXL ☐	M/XXXXL ☐	

1 - Cannot be combined with Anti Tippers


2 - Not available with the swing away armrest ; Not available with HD option

3 - Not available with Mag wheels et Spinergy

FRAME COLOR							
150	FC01	Black ☐				NC	☐
151	FC03	Sunkissed Orange ☐	Burgundy ☐	Sapphire Blue ☐	Monster Green ☐	additional 5 day lead-time	\$270 ☐
		Charcoal ☐	Fuchsia ☐	Steel Blue ☐	Ferrari Red ☐	Acid Green ☐	
152	FC99	Custom Color	Specify PPG Paintit color:			additional 15 days lead-time	\$410 ☐

TOTAL							
-------	--	--	--	--	--	--	--

Fork and Stem bolt combinations



Notes :



BioFit™ CUSHION

	BioFit™	Fits Width	Depth	\$	✓
	NXCABF1616	15"-16"	16"	595	
	NXCABF1618	15"-16"	18"	595	
	NXCABF1620	15"-16"	20"	595	
	NXCABF1816	17"-18"	16"	595	
	NXCABF1818	17"-18"	18"	595	
	NXCABF1820	17"-18"	20"	595	

BioFit™	Fits Width	Depth	\$	✓
NXCABF2016	19" - 20"	16"	595	
NXCABF2018	19" - 20"	18"	595	
NXCABF2020	19" - 20"	20"	625	
NXCABF2216	21" - 22"	16"	625	
NXCABF2218	21" - 22"	18"	625	
NXCABF2220	21" - 22"	20"	625	

Optima™ THORACIC BACK SUPPORT

	Optima™	Fits Width	Height	\$	✓
	NOACB1613	16"-18"	13"	690	
	NOACB1616	16"-18"	16"	690	
	NOACB1619	16"-18"	19"	690	
	NOACB1813	18"-20"	13"	690	
	NOACB1816	18"-20"	16"	690	
	NOACB1819	18"-20"	19"	690	
	NOACB2016	20"-22"	16"	690	
	NOACB2019	20"-22"	19"	690	
	NOACB2213	22"-24"	13"	720	
	NOACB2216	22"-24"	16"	720	
	NOACB2219	22"-24"	19"	720	

Optima™ DEEP THORACIC BACK SUPPORT

	Optima™ DEEP	Fits Width	Height	\$	✓
	NOPCB1616	16"-18"	16"	790	
	NOPCB1619	16"-18"	19"	790	
	NOPCB1816	18"-20"	16"	790	
	NOPCB1819	18"-20"	19"	790	
	NOPCB2016	20"-22"	16"	790	
	NOPCB2019	20"-22"	19"	790	
	NOPCB2213	22"-24"	13"	820	
	NOPCB2216	22"-24"	16"	820	
	NOPCB2219	22"-24"	19"	820	

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	690	
	NHAACB1816/20	18"-20"	16"-20"	690	
	NHAACB2016/20	20"-22"	16"-20"	690	
	NHAACB2216/20	22"-24"	16"-20"	720	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™ DEEP	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	690	
	NHAACB1816/20	18"-20"	16"-20"	690	
	NHAACB2016/20	20"-22"	16"-20"	690	
	NHAPCB2216/20	22"-24"	16"-20"	720	

Classic™ PELVIC BACK SUPPORT

	Classic™	Fits Width	Height	\$	✓
	NACB1616	16"-18"	16"	690	
	NACB1618	16"-18"	18"	690	
	NACB1816	18"-20"	16"	690	
	NACB1818	18"-20"	18"	690	
	NACB1820	18"-20"	20"	690	
	NACB2020	20"-22"	20"	660	
	NACB2216	22"-24"	16"	720	
	NACB2218	22"-24"	18"	720	
	NACB2220	22"-24"	20"	720	

Classic™ DEEP PELVIC BACK SUPPORT

	Classic™ DEEP	Fits Width	Height	\$	✓
	NPCB1616	16"-18"	16"	690	
	NPCB1618	16"-18"	18"	690	
	NPCB1818	18"-20"	18"	690	
	NPCB1820	18"-20"	20"	690	
	NPCB2020	20"-22"	20"	690	
	NPCB2216	22"-24"	16"	720	
	NPCB2218	22"-24"	18"	720	
	NPCB2220	22"-24"	20"	720	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

		
Model #NIB-0784	\$	✓
10" x 5" Large Rectangular	195	

		
Model # NIB-0818	\$	✓
14" x 5" Extra-Large Rectangular	195	

MultiFit™ HEAD SUPPORT – FIXTURES

		
Model #NIB-0747	\$	✓
Multi-Axis Hardware	145	

		
Model # NIB-0823	\$	✓
Head Support Receiver Bracket	55	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com