



Helio A7 - FiveDay

ORDER FORM - CANADA

P.O. number: * _____ Date: _____

Company Name: * _____ Phone: * _____

Ordered by: * _____

Clinician/Sales rep: _____

Email: _____

Marked for: _____

Ship to: * _____

Name: _____

Address: _____

City: _____ Prov: _____

Postal Code: _____

Phone: _____

* Required MA080 English Manual (STD) MA081 Manuel Français Taxes not included

MODEL							
1	WC02	Ultralight Helio A7 wheelchair 7000 series aluminum (265 lb weight capacity)		WM0001020	\$2 645	<input type="checkbox"/>	
2	WC09	HD Heavy duty kit (350 lb weight capacity)		WAMF	\$775	<input type="checkbox"/>	
3	WC10	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber			\$30		
SEAT WIDTH ³							
4	SW	Seat width:	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" ¹ <input type="checkbox"/>
5	SW	Seat width:	19" ¹ <input type="checkbox"/>	20" ¹ <input type="checkbox"/>	21" ² <input type="checkbox"/>	22" ² <input type="checkbox"/>	WAMA \$295 <input type="checkbox"/>
SEAT DEPTH							
6	SD	Seat depth:	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <input type="checkbox"/>
7	SD	Seat depth:	19" <input type="checkbox"/>	20" <input type="checkbox"/>			WAMB \$225 <input type="checkbox"/>

1 - Available with HD option 2 - Only available with HD option [3. Width calculator](#)

A7 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"															
4"	P1	13	14	15	13 1/2	14 1/2	15 1/2												
	P2	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2									
5"	P1																		
	P2				14	15	16	14 1/2	15 1/2	16 1/2									
	P3	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2						
	P4	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
7"	P1									16 1/2	17 1/2	18 1/2	17	18	19				
	P2							16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2			
	P3							17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2
	P4	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21
	P5	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21	19 1/2	20 1/2	21 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS											
8	CA01	Composite wheel w/ PU tire (1" width)	3" <input type="checkbox"/>	4" <input type="checkbox"/>	5" <input type="checkbox"/>	6" <input type="checkbox"/>	7" <input type="checkbox"/>	8" <input type="checkbox"/>			
9	CA10	Pneumatic casters 6" x 1 1/4"								NC <input type="checkbox"/>	
10	CA11	Pneumatic casters 8" x 1 1/4"								\$160 <input type="checkbox"/>	
11	CA14	Newton UltraCaster Soft-Roll Composite	4" x 1 1/2"								\$98 <input type="checkbox"/>
12	CA15		5" x 1 1/2"								\$98 <input type="checkbox"/>
13	CA16		6" x 1 1/2"								\$98 <input type="checkbox"/>
14	CA17		4" x 1 1/2"								\$175 <input type="checkbox"/>
15	CA18	Newton UltraCaster Soft-Roll Aluminum	5" x 1 1/2"								\$185 <input type="checkbox"/>
16	CA19		6" x 1 1/2"								\$195 <input type="checkbox"/>
CASTER OPTIONS											
17	CA12	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY)							WAMS	\$255 <input type="checkbox"/>	
18	CA13	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)								\$489 <input type="checkbox"/>	

A7 - REAR Seat-to-floor Height

WHEEL SIZES	Height available	Specify desired height: _____
20" wheel <input type="checkbox"/>	12 1/2 to 17 1/4	
22" wheel <input type="checkbox"/>	13 1/2 to 18 1/4	
24" wheel <input type="checkbox"/>	14 1/2 to 19 1/4	
25" wheel <input type="checkbox"/>	15 to 19 3/4	
26" wheel <input type="checkbox"/>	15 1/2 to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

Footrest										
SWING-IN/SWING-OUT FOOTREST										
19	FR03	60° Swing-in/swing-out footrest	13.5 - 17" <input type="checkbox"/>	14.5 - 19.5" <input type="checkbox"/>	15.5 - 20.5" <input type="checkbox"/>				Set at:*	NC <input type="checkbox"/>
20	FR04	70° Swing-in/swing-out footrest	13.5 - 17" <input type="checkbox"/>	14.5 - 19.5" <input type="checkbox"/>	15.5 - 20.5" <input type="checkbox"/>				Set at:*	NC <input type="checkbox"/>
21	FR15	80° Swing-in/swing-out footrest	13.5 - 16" <input type="checkbox"/>	14.5 - 19.5" <input type="checkbox"/>	15.5 - 20.5" <input type="checkbox"/>				Set at:*	NC <input type="checkbox"/>
90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT (Footplate mounted to hanger tube)										
22	FR06	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5 - 11"						Set at:*	NC <input type="checkbox"/>
23	FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{1,2}	9.5 - 13"						Set at:*	NC <input type="checkbox"/>
24	FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{1,2}	9.5 - 14.5"						Set at:*	NC <input type="checkbox"/>
FOOTREST OPTIONS ³										
25	FR09	High mount (footplate attached to the hanger with clamps)	5 - 13.5" <input type="checkbox"/>						Set at:*	\$90 <input type="checkbox"/>
26	FR99	Omit Footrest								NC <input type="checkbox"/>

1 - Available only with adjustable footplates 2 - May interfere with front casters 3 - Not available with front mount * +/- 1/2"

ELEVATING LEGREST (unit) The elevating legrest replace the footrest side selected										
27	FR10	Elevating legrest with calf pad, 14" - 22" Length adjustable	Left <input type="checkbox"/>	Right <input type="checkbox"/>					Set at:	WAMK \$290 <input type="checkbox"/>

FOOTPLATES										
28	FP01	Standard - Composite								STD <input type="checkbox"/>
29	FP02	Newton adjustable angle - Composite						WAMJ		\$145 <input type="checkbox"/>
30	FP03	Newton adjustable angle - Aluminum						WAMJ		\$165 <input type="checkbox"/>
31	FP04	Newton adjustable angle - Carbon fiber (with rubber inserts)						WAMJ		\$380 <input type="checkbox"/>
32	FP06	One piece adjustable angle flip-up - Aluminum						WAMJ	SEFND0015	\$206 <input type="checkbox"/>
33	FP07	One piece adjustable angle flip-up - Carbon fiber (with rubber inserts)						WAMJ	SEFND0015	\$450 <input type="checkbox"/>

FOOTPLATE ACCESSORIES										
34	FP08	Length adjustable heel loop (not available with one piece footplate)								NC <input type="checkbox"/>
35	FP09	Calf strap							SEFND0010	\$52 <input type="checkbox"/>
36	FP30	Bodypoint Aeromesh Padded calf strap							SEFND0010	\$129 <input type="checkbox"/>

RESIDUAL LIMB SUPPORT FLAT CUSHION (Swing away, depth, width, height and angle adjustable) - (unit)								
37	FR27	Residual limb support swing away - flat 6"x8"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	New	SEFND0005	\$475	<input type="checkbox"/>
38	FR28	Residual limb support swing away - flat 8"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	New	SEFND0005	\$475	<input type="checkbox"/>
39	FR29	Residual limb support swing away - flat 10"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	New	SEFND0005	\$475	<input type="checkbox"/>
40	FR30	Residual limb support swing away - flat 10"x14"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	New	SEFND0005	\$475	<input type="checkbox"/>

RESIDUAL LIMB SUPPORT CONTOUR CUSHION (Swing away, depth, width, height and angle adjustable) - (unit)								
41	FR32	Residual limb support swing away - contour 6"x8"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	New	SEFND0005	\$475	<input type="checkbox"/>
42	FR33	Residual limb support swing away - contour 8"x10"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	New	SEFND0005	\$475	<input type="checkbox"/>

Back canes with integrated fixed height push handle

FIXED ANGLE BACK								
43	BC01	8° bend aluminum back cane	16" <input type="checkbox"/>	18" <input type="checkbox"/>	20" STD <input type="checkbox"/>			NC <input type="checkbox"/>

Back canes, height adjustable, with option

FIXED ANGLE BACK								
44	BC11	Straight aluminum back cane	9-12" <input type="checkbox"/>	12 - 15" <input type="checkbox"/>	15 - 18" <input type="checkbox"/>	18 - 21" STD <input type="checkbox"/>		\$75 <input type="checkbox"/>
45	BC12	8° bend aluminum back cane	13.5 - 15" <input type="checkbox"/>	15 - 18" <input type="checkbox"/>	18 - 21" STD <input type="checkbox"/>			\$75 <input type="checkbox"/>

ADJUSTABLE ANGLE BACK								
46	BC13	Adjustable angle, straight aluminum back cane		15 - 18" <input type="checkbox"/>	18 - 21" STD <input type="checkbox"/>			\$259 <input type="checkbox"/>
47	BC14	Adjustable angle, 8° bend aluminum back cane		15 - 18" <input type="checkbox"/>	18 - 21" STD <input type="checkbox"/>			\$309 <input type="checkbox"/>
48		Select back cane angle:	85° ¹ <input type="checkbox"/>	90° <input type="checkbox"/>	95° STD <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>

Set height to: ² (measured from seat sling to top of back upholstery)

PUSH HANDLE OPTION								
49	PH01	Push handles						STD <input type="checkbox"/>
50	PH05	Fold down push handles ^{3,5}						\$150 <input type="checkbox"/>
51	PH06	Clamp on height adjustable push handle ^{3,4,5}	9 - 20" STD <input type="checkbox"/>	9-23" <input type="checkbox"/>				\$325 <input type="checkbox"/>
52	PH02	Without push handles						NC <input type="checkbox"/>

1 - Not available with a 8 degree bend cane 2 - If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit

4 - Not compatible with swing away armrests, 8° bend back canes and amputee axle plates.

5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

STROLLER HANDLE AND STABILIZER BAR ¹ (only available with push handle)								
53	PH04	Stroller handle (not available with the HD option)					WAMW	\$210 <input type="checkbox"/>
54	PH03	Newton folding stabilizer bar						\$258 <input type="checkbox"/>

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY								
55	BH07	Slip-on Active-Mesh back upholstery						NC <input type="checkbox"/>
56	BH08	Adjustable tension Active-Mesh back upholstery					WAMD/WAME	\$245 <input type="checkbox"/>
57	BH99	Omit Back upholstery						STD <input type="checkbox"/>

RIGID BACK								
58		NXT back series					See pricelist p. 7	<input type="checkbox"/>

SEAT CUSHION								
59	SC01	2" Cushion						\$175 <input type="checkbox"/>
60	SC02	3" Cushion						\$195 <input type="checkbox"/>
61	SC99	Omit seat cushion						STD <input type="checkbox"/>
62		NXT FIT cushion series					See pricelist p. 7	<input type="checkbox"/>

Rear wheels											
ONE ARM DRIVE ¹ <i>Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims</i>											
63	HR22	One arm drive ^{2,3}	24"	left <input type="checkbox"/>	right <input type="checkbox"/>				\$1 395	<input type="checkbox"/>	
64	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)							NC	<input type="checkbox"/>	
REAR WHEELS											
65	RW01	Mag ⁴	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>				\$129	<input type="checkbox"/>	
66	RW02	Newton One - spoke wheel ⁵	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>				NC	<input type="checkbox"/>	
67	RW03	Newton Gravity - ultralight wheel	20" <input type="checkbox"/>	22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>		\$360	<input type="checkbox"/>	
68	RW04	Spinergy Spox		22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>		\$800	<input type="checkbox"/>	
69		Spinergy Spox spokes colors ⁶	Blk STD <input type="checkbox"/>	Wht <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Yel <input type="checkbox"/>		NC	<input type="checkbox"/>	
70	RW05	Spinergy LX		22" <input type="checkbox"/>	24" <input type="checkbox"/>	25" <input type="checkbox"/>	26" <input type="checkbox"/>		\$995	<input type="checkbox"/>	
71		Spinergy LX spokes colors ⁶	Blk STD <input type="checkbox"/>	Pur <input type="checkbox"/>	Wht <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Yel <input type="checkbox"/>	Pink <input type="checkbox"/>	Grn <input type="checkbox"/>	Org <input type="checkbox"/>
72	RW07	Spinergy CLX Carbon Blade			24" <input type="checkbox"/>	25" <input type="checkbox"/>			\$1 540	<input type="checkbox"/>	
73		Spinergy CLX spokes colors ⁶	Blk STD <input type="checkbox"/>	Pur <input type="checkbox"/>	Wht <input type="checkbox"/>	Blue <input type="checkbox"/>	Red <input type="checkbox"/>	Yel <input type="checkbox"/>	Pink <input type="checkbox"/>	Grn <input type="checkbox"/>	Org <input type="checkbox"/>
TIRES											
Available											
74	TI01	Soft Urethane 1 1/4 "				20", 22" 24"			STD	<input type="checkbox"/>	
75	TI02	Soft Urethane knobby 1 3/8 "				22", 24", 25", 26"			\$110	<input type="checkbox"/>	
76	TI04	Soft Urethane - Shox 1"				22", 24", 25", 26"			\$130	<input type="checkbox"/>	
77	TI05	Pneumatic 1 3/8 "				20", 22", 24", 25", 26"			NC	<input type="checkbox"/>	
78	TI08	Pneumatic with Airless insert 1 3/8"				22", 24"			\$215	<input type="checkbox"/>	
79	TI06	Pneumatic HP - 1"				20", 22", 24", 25", 26"			\$110	<input type="checkbox"/>	
80	TI07	Pneumatic HP - Schwalbe Marathon Plus 1"				24", 25", 26"			\$190	<input type="checkbox"/>	

1- Offered with handrims on each side of the wheelchair

2- Only standard handrim position available

3 - Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position

5 - Not compatible with Shox tires

6 - Colour spokes require an additional 10 days lead-time

Handrims																	
HANDRIMS																	
81	HR01	Aluminum anodized							STD	<input type="checkbox"/>							
82	HR02	Plastic coated ³						WAMO	\$160	<input type="checkbox"/>							
83	HR03	Q-Grip High friction Neoprene coated ^{1,2,3}						WAMO	\$430	<input type="checkbox"/>							
84	HR07	Newton Air Grip ¹	Black <input type="checkbox"/>					WAMO	\$150	<input type="checkbox"/>							
85	HR08	Newton Air Grip ¹	(08) Sunkissed Orange <input type="checkbox"/>	(09) Burgundy <input type="checkbox"/>	(10) Sapphire Blue <input type="checkbox"/>	(11) Monster Green <input type="checkbox"/>	(12) Charcoal <input type="checkbox"/>	(13) Fuchsia <input type="checkbox"/>	(14) Steel Blue <input type="checkbox"/>	(15) Ferrari Red <input type="checkbox"/>	(16) Acid Green <input type="checkbox"/>	(17) White <input type="checkbox"/>	(18) Sky Blue <input type="checkbox"/>	Additional 5 day lead-time	WAMO	\$330	<input type="checkbox"/>
86	HR28	Surge ^{1,2,3}	(1-1/2" Oval handrim w/ Gripton high friction strip)							\$470	<input type="checkbox"/>						
87	HR05	Surge LT ^{1,2,3}	(1-1/8" Oval handrim w/ Gripton high friction strip)							\$500	<input type="checkbox"/>						
88	HR04	Natural Fit Standard Grip ^{1,2,3}	(1-1/2" Oval handrim w/standard thumb piece)							\$500	<input type="checkbox"/>						
89	HR99	Omit handrim							NC	<input type="checkbox"/>							

1 - Not available with Mag wheels or OAD

2 - Not available in 20"

3- Only standard handrim position available

Handrims option											
Assembly position for aluminum and Newton AirGrip handrims											
90	HR29	Standard Handrim position	Full tab							STD	<input type="checkbox"/>
91	HR30	Narrow Handrim position	Short tab 2 positions							\$25	<input type="checkbox"/>
92	HR31	Super Narrow Handrim	Short tab 2 positions STD <input type="checkbox"/>		Short tab 1 position <input type="checkbox"/>				\$25	<input type="checkbox"/>	
Thumb piece⁴ selection Natural Fit handrims											
93	HR23	Standard grip							\$190	<input type="checkbox"/>	
94	HR27	Super grip							\$220	<input type="checkbox"/>	
95	HR099	Omit thumb piece							NC	<input type="checkbox"/>	


4- Compatible with Surge, Surge LT handrims

AXLE										
96	AX01	Quick release axle						WAMU	STD	<input type="checkbox"/>
97	AX02	Permanent axle							NC	<input type="checkbox"/>
98	AX03	Quad release axle						WAMU	\$170	<input type="checkbox"/>
AXLE PLATE										
99	AX04	Amputee adjustable axle plate ¹						WAMT	\$250	<input type="checkbox"/>

1 - May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER				
100	CM00	0° camber		STD <input type="checkbox"/>
101	CM03	3° camber		NC <input type="checkbox"/>
102	CM06	6° camber <i>(not available with amputee axle plate)</i>		\$80 <input type="checkbox"/>

CENTER OF GRAVITY



Choose position of rear wheel axle plate:
Measure distance X from front of back post to center of axle receiver.

X =	Most stable								Most responsive			
	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"
	<input type="checkbox"/>											

If not indicated, default setting is at 1.5"
Possible conflict between a seat depth under 16" and a center of gravity over 2.5".
Possible conflict between a 26" rear wheel and a center of gravity over 2.5".
If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK				
103	WL08	Newton Wheel lock - push to lock		STD <input type="checkbox"/>
104	WL02	Pull to lock		NC <input type="checkbox"/>
105	WL04	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>		WAMR \$175 <input type="checkbox"/>
106	WL13 WL14	6" removable extension handle		\$67 <input type="checkbox"/>
107	WL05	Scissor Lock		\$175 <input type="checkbox"/>
108	WL06	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>		\$230 <input type="checkbox"/>
109	WL99	Omit wheel locks		NC <input type="checkbox"/>

Armrests					
FLIP BACK ARMREST <i>(Convertible to "T" armrest)^{1,4}</i>		Height Adjustment		Removable function	
110	AR01	10" Desk length arm pad	8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>	Locked <input type="checkbox"/> Unlocked <input type="checkbox"/> STD \$380 <input type="checkbox"/>
111	AR02	14" Full length arm pad	8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>	Locked <input type="checkbox"/> Unlocked <input type="checkbox"/> STD \$380 <input type="checkbox"/>
"T" ARMREST ¹		Height Adjustment			
112	AR03	10" Desk length arm pad	8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>	\$360 <input type="checkbox"/>
113	AR04	14" Full length arm pad	8 - 12" STD <input type="checkbox"/>	10 - 14" <input type="checkbox"/>	\$360 <input type="checkbox"/>
SWING AWAY ARMREST ²		Height Adjustment			
114	AR05	Tubular arm pad	8 - 11"		\$295 <input type="checkbox"/>
115	AR06	10" Desk length arm pad	8 - 11"		\$295 <input type="checkbox"/>
GEL ARM PAD ³					
116	AR12	14" Full Length Arm Pad - Gel Ovation			SEAND0030 \$150 <input type="checkbox"/>
117	AR99	Omit armrest			NC <input type="checkbox"/>

1 - May have a 1/2" adjustment restriction on both limits of the range 2 - Could interfere with some rigid back rests 3 - Not available with swing away armrest 4 - Non available with attendant lock

SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>				
118	SG01	Plastic side guards <i>(removable)</i>		WAMV \$225 <input type="checkbox"/>
119	SG02	Carbon fiber side guards <i>(fixed)</i>		WAMV \$255 <input type="checkbox"/>
120	SG03	Carbon fiber side guards <i>(removable)</i>		WAMV \$295 <input type="checkbox"/>
121	SG06	Plastic fender side guards <i>(removable)¹</i>		WAMV \$320 <input type="checkbox"/>
122	SG04	Carbon fiber fender side guards <i>(removable)</i>		WAMV \$495 <input type="checkbox"/>
123	SG99	Omit side guard		NC <input type="checkbox"/>
SIDE GUARDS OPTION				
124	SG07	Side guards clip <i>(pair)²</i>		New NC <input type="checkbox"/>

1- Not compatible with 20" and 22" rear wheels 2- Not available with fixed carbon fiber side guard

POSITIONING BELT										
125	SB01	Velcro adjustable	Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>					\$45	<input type="checkbox"/>
126	SB02	Auto buckle	Standard 35" <input type="checkbox"/>	Long 60" <input type="checkbox"/>					\$55	<input type="checkbox"/>
127	SB07	Hardware attachment for 2 point seatbelt (pair)							NC	<input type="checkbox"/>
128	SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (pair)							\$13	<input type="checkbox"/>
129	SB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>					\$82	<input type="checkbox"/>
130	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" <input type="checkbox"/>	Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>			SERND0030	\$170	<input type="checkbox"/>
131	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" <input type="checkbox"/>	Medium 56" <input type="checkbox"/>	Long 61" <input type="checkbox"/>			SERND0020	\$195	<input type="checkbox"/>
132	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" <input type="checkbox"/>	Long 24" <input type="checkbox"/>					\$324	<input type="checkbox"/>
133	SB99	Omit positioning belt							STD	<input type="checkbox"/>

HEADREST ¹										
134	AC01	Profiled Headrest	Small (3.5"X6") <input type="checkbox"/>	Medium (4.5"X7") <input type="checkbox"/>	Large (5.5"X8") <input type="checkbox"/>				\$205	<input type="checkbox"/>
135	AC28	Headrest mounting hardware with removable non folding stabilizer bar							\$338	<input type="checkbox"/>

1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar


ACCESSORIES														
136	AC03	Anti tippers							\$115	<input type="checkbox"/>				
137	AC27	Tip assist ¹							\$60	<input type="checkbox"/>				
138	AC04	Transit Tie Down (Unoccupied)							NC	<input type="checkbox"/>				
139	AC05	WC-19 Transit Tie Down w/Q'straint belt ²	Short belt 20-41" <input type="checkbox"/>	Long belt 20-46" <input type="checkbox"/>					\$345	<input type="checkbox"/>				
140	AC06	Cane carrier	Left <input type="checkbox"/>	Right <input type="checkbox"/>					\$105	<input type="checkbox"/>				
141	AC07	Oxygen tank holder - E Type <i>Not compatible with swing away armrest or amputee axle plate</i>	Left					WAMX	\$375	<input type="checkbox"/>				
142	AC08	Spoke guards ³	Available: 20", 22", 24", 26"						WAMP	\$145	<input type="checkbox"/>			
143	AC09	Footrest neoprene impact guard (pair) - Black						SEFND0035	\$80	<input type="checkbox"/>				
144	AC10	Seat pouch							\$95	<input type="checkbox"/>				
145	AC26	BodyPoint Mobility Bag							\$135	<input type="checkbox"/>				
146	AC23	Newton Backpack							\$120	<input type="checkbox"/>				
147	AC12	Velcro Closing Strap							\$50	<input type="checkbox"/>				
148	AC13	Tool kit							\$50	<input type="checkbox"/>				
149	AC16	Motion Composites Polo Shirt	W/S <input type="checkbox"/>	W/M <input type="checkbox"/>	W/L <input type="checkbox"/>	M/S <input type="checkbox"/>	M/M <input type="checkbox"/>	M/L <input type="checkbox"/>	M/XL <input type="checkbox"/>	M/XXL <input type="checkbox"/>	M/XXXL <input type="checkbox"/>	M/XXXXL <input type="checkbox"/>	\$32	<input type="checkbox"/>

1 - Cannot be combined with Anti Tippers 2 - Not available with the swing away armrest ; Not available with HD option 3 - Not available with Mag wheels et Spinergy

FRAME COLOR													
150	FC01	Black <input type="checkbox"/>							NC	<input type="checkbox"/>			
151	FC03	Sunkissed Orange <input type="checkbox"/>	Burgundy <input type="checkbox"/>	Sapphire Blue <input type="checkbox"/>	Monster Green <input type="checkbox"/>	Charcoal <input type="checkbox"/>	Fuchsia <input type="checkbox"/>	Steel Blue <input type="checkbox"/>	Ferrari Red <input type="checkbox"/>	Acid Green <input type="checkbox"/>	additional 5 day lead-time	\$270	<input type="checkbox"/>
152	FC99	Custom Color	Specify PPG Paintit color:						additional 15 days lead-time	\$410	<input type="checkbox"/>		

TOTAL

Fork and Stem bolt combinations



Notes :



BioFit™ CUSHION

	BioFit™	Fits Width	Depth	\$	✓
	NXCABF1616	15"-16"	16"	595	
	NXCABF1618	15"-16"	18"	595	
	NXCABF1620	15"-16"	20"	595	
	NXCABF1816	17"-18"	16"	595	
	NXCABF1818	17"-18"	18"	595	
	NXCABF1820	17"-18"	20"	595	

BioFit™	Fits Width	Depth	\$	✓
NXCABF2016	19"-20"	16"	595	
NXCABF2018	19",-20"	18"	595	
NXCABF2020	19"-20"	20"	625	
NXCABF2216	21"-22"	16"	625	
NXCABF2218	21"-22"	18"	625	
NXCABF2220	21"-22"	20"	625	

Optima™ THORACIC BACK SUPPORT

	Optima™	Fits Width	Height	\$	✓
	NOACB1613	16"-18"	13"	690	
	NOACB1616	16"-18"	16"	690	
	NOACB1619	16"-18"	19"	690	
	NOACB1813	18"-20"	13"	690	
	NOACB1816	18"-20"	16"	690	
	NOACB1819	18"-20"	19"	690	
	NOACB2016	20"-22"	16"	690	
	NOACB2019	20"-22"	19"	690	
	NOACB2213	22"-24"	13"	720	
	NOACB2216	22"--24"	16"	720	
	NOACB2219	22"-24"	19"	720	

Optima™ DEEP THORACIC BACK SUPPORT

	Optima™ DEEP	Fits Width	Height	\$	✓
	NOPCB1616	16"-18"	16"	790	
	NOPCB1619	16"-18"	19"	790	
	NOPCB1816	18"-20"	16"	790	
	NOPCB1819	18"-20"	19"	790	
	NOPCB2016	20"-22"	16"	790	
	NOPCB2019	20"-22"	19"	790	
	NOPCB2213	22"-24"	13"	820	
	NOPCB2216	22"-24"	16"	820	
	NOPCB2219	22"-24"	19"	820	

THORACIC SHAPE

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	690	
	NHAACB1816/20	18"-20"	16"-20"	690	
	NHAACB2016/20	20"-22"	16"-20"	690	
	NHAACB2216/20	22"-24"	16"-20"	720	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™ DEEP	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	690	
	NHAACB1816/20	18"-20"	16"-20"	690	
	NHAACB2016/20	20"-22"	16"-20"	690	
	NHAPCB2216/20	22"-24"	16"-20"	720	

Classic™ PELVIC BACK SUPPORT

	Classic™	Fits Width	Height	\$	✓
	NACB1616	16"-18"	16"	690	
	NACB1618	16"-18"	18"	690	
	NACB1816	18"-20"	16"	690	
	NACB1818	18"-20"	18"	690	
	NACB1820	18"-20"	20"	690	
	NACB2020	20"-22"	20"	660	
	NACB2216	22"-24"	16"	720	
	NACB2218	22"-24"	18"	720	
	NACB2220	22"-24"	20"	720	

Classic™ DEEP PELVIC BACK SUPPORT

	Classic™ DEEP	Fits Width	Height	\$	✓
	NPCB1616	16"-18"	16"	690	
	NPCB1618	16"-18"	18"	690	
	NPCB1818	18"-20"	18"	690	
	NPCB1820	18"-20"	20"	690	
	NPCB2020	20"-22"	20"	690	
	NPCB2216	22"-24"	16"	720	
	NPCB2218	22"-24"	18"	720	
	NPCB2220	22"-24"	20"	720	

PELVIC SHAPE

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

	
Model #NIB-0784	Model # NIB-0818
\$ 195	\$ 195
10" x 5" Large Rectangular	14" x 5" Extra-Large Rectangular

MultiFit™ HEAD SUPPORT – FIXTURES

	
Model #NIB-0747	Model # NIB-0823
\$ 145	\$ 55
Multi-Axis Hardware	Head Support Receiver Bracket

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com