



Helio A7

ORDER FORM - INTERNATIONAL

1- Purchasing information

P.O. number:* _____ Date: _____
Company Name:* _____
Purchasing contact:* _____
Phone:* _____
Email:* _____
Contact preference:* Phone ☐ Email ☐
Marked for: _____

3- Shipping information

Ship to:* _____
Name: _____
Address: _____
City: _____ State: _____
Zip Code: _____
Phone: _____

2- ATP/Clinician information

ATP/Clinician:* _____
Phone:* _____
Email:* _____
Contact preference:* Phone ☐ Email ☐

MA083 ☐ English Manual (STD) MA084 ☐ German Manual MA085 ☐ Italian Manual MA086 ☐ Dutch Manual MA087 ☐ Spanish Manual MA088 ☐ French Manual MA089 ☐ Portuguese Manual

MODEL								USD (taxes not included)		
WC02	☐	Ultralight Helio A7 wheelchair 7000 series aluminum (120 kg/265 lb weight capacity)								\$2,625
WC09	☐	HD Heavy duty kit (159 kg/ 350 lb weight capacity)								\$540
Frame										
SEAT WIDTH ³ - cm (inches)										
SW	☐	35.56 (14") ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") ¹ ☐	48.26 (19") ¹ ☐			NC
	☐	50.80 (20") ¹ ☐ 53.34 (21") ² ☐ 55.88 (22") ² ☐								\$535
SEAT DEPTH - cm (inches)										
SD	☐	35.56 (14") ☐	38.10 (15") ☐	40.64 (16") ☐	43.18 (17") ☐	45.72 (18") ☐	48.26 (19") ☐			NC
	☐	50.80 (20") ☐								\$515
FRAME LENGTH ⁵										
SP01	☐	Standard frame length configuration ⁶								STD
SP02	☐	Frame length 1" longer than standard configuration								NC
SP03	☐	Frame length 2" longer than standard configuration								NC
SP04	☐	Frame length 1" shorter than standard configuration								NC
TRANSIT TIE DOWN										
AC04	☐	Transit Tie Down (Unoccupied)								NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 50-104 cm (20" - 41")								\$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 50-116 cm (20" - 46")								\$365
AC35	☐	WC-19 Transit Tie Down without belt ⁴								\$335

1- Available with HD option; 2- Only available with HD option

4- May interfere with swing away armrest and Clamp on height adjustable push handle;

5- See table on page 7; 6- Distance of 2 1/8" from end of seat upholstery to end of frame;

[3. Width calculator](#)


PACKAGES

VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹	\$95
------	--------------------------	--	------

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo *VP01; Recommended with black chair; See picture of the product on page 7.

A7 FRONT Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"	STD	+ 2,54 cm 1"	+ 5,08cm 2"
4 in 10,2 cm	P1	33.02cm 13"	35.56cm 14"	38.10cm 15"	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"												
	P2	34.29cm 13.5"	36.83cm 14.5"	39.37cm 15.5"	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
5 in 12,7 cm	P1																		
	P2				35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"									
	P3	35.56cm 14"	38.10cm 15"	40.64cm 16"	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"						
	P4	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"	40.64cm 16"	43.18cm 17"	45.72cm 18"	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"			
7 in 17,8 cm	P1										41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"			
	P2							41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"			
	P3							43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"
	P4	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"
	P5	43.18cm 17"	45.72cm 18"	48.26cm 19"	44.45cm 17.5"	46.99cm 18.5"	49.53cm 19.5"	45.72cm 18"	48.26cm 19"	50.80cm 20"	46.99cm 18.5"	49.53cm 19.5"	52.07cm 20.5"	48.26cm 19"	50.80cm 20"	53.3cm 21"	49.53cm 19.5"	52.07cm 20.5"	54.6cm 21.5"
Froglegs	P1				36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"												
	P2	36.83cm 14.5"	39.37cm 15.5"	41.91cm 16.5"	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"									
	P3	38.10cm 15"	40.64cm 16"	43.18cm 17"	39.37cm 15.5"	41.91cm 16.5"	44.45cm 17.5"	40.64cm 16"	43.18cm 17"	45.72cm 18"	41.91cm 16.5"	44.45cm 17.5"	46.99cm 18.5"						

Height available + 0,65cm (1/4") (see stem bolt and fork combination at the end of the order form)

FRONT CASTERS							Availability		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")					3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 1/4")					6", 8"		\$175
CA14 to CA16	☐	Newton UltraCaster Composite (1 1/2 ")					4", 5", 6"		\$105
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 1/2 ")					4", 5", 6"		\$210
CASTER OPTIONS									
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹							\$320
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)							\$515

1- Not Available with WC-19


A7 - REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height: _____
☐	20" (451) wheel	31.75cm/12.5 "	to 43.81cm/17.25"	
☐	22" (501) wheel	34.29cm/13.5 "	to 46.35cm/18.25"	
☐	24" (540) wheel	36.83cm/14.5 "	to 48.89cm/19.25"	
☐	25" (559) wheel	38.10cm/15 "	to 50.16cm/19.75"	
☐	26" (590) wheel	39.37cm/15.5 "	to 51.34cm/20.25"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1,25 cm (1/2") of selected Seat to Floor height"

REAR WHEELS						
ONE ARM DRIVE ¹				Only available with Newton One rear wheels permanent axles, 0° camber		
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐		\$1,510
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐		\$1,595
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)				NC
REAR WHEELS				Available		
RW01	☐	Mag ³			20", 22" 24"	\$57
RW02	☐	Newton One - spoke wheel ⁴			20", 22" 24"	NC
RW03	☐	Newton Gravity - ultralight wheel			20", 22", 24", 25", 26"	\$350
RW04	☐	Spinerger Spox	Blk STD ☐	Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐	22", 24", 25", 26"	\$895
RW05	☐	Spinerger LX	Blk STD ☐ Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐		22", 24", 25", 26"	\$1,125
RW07	☐	Spinerger CLX ⁶	Blk STD ☐ Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐		24", 25"	\$1,600
TIRES (The brand and color may vary depending on parts availability)				Available		
AIRLESS						
TI01	☐	Soft Urethane 3.49cm (1 3/8") - Low Tread, black ⁹			20", 22" 24"	STD
TI02	☐	Soft Urethane 3.49cm (1 3/8") - Medium Tread, black ^{7,9}			22", 24", 25", 26"	\$105
TI13	☐	Soft Urethane 2.54cm (1") - No tread, black - <i>Sentinel</i>			20", 22", 24", 25", 26"	\$125
TI08	☐	Pneumatic Airless insert 3.49cm (1 3/8") - Medium Tread ⁸			22", 24", 26"	\$90
PNEUMATIC						
TI05	☐	Pneumatic 3.49cm (1 3/8") Medium Tread ⁹			20", 22", 24", 26"	NC
TI06	☐	Pneumatic HP 2.54cm (1") - Low Tread, Puncture resistant, black - <i>SpeedLite</i>			20", 22", 24", 25", 26"	\$110
TI07	☐	Pneumatic HP 2.54cm (1") - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>			24", 25", 26"	\$145
TI03	☐	Pneumatic 5.08cm (2") - Large Tread, black ¹⁰			22", 24", 25"	\$115
1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available; 3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additional 10 days lead-time; 6- Only available with 1" tires; 7- 22" is only available with Mags, Newton One and Newton Gravity wheels 8- Not available with Spinerger CLX, Spinerger LX 25", 26" and Spox 25", 26" rear wheels; 9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;						
HANDRIMS						
HANDRIM				Available		
HR01	☐	Aluminum anodized			20", 22", 24", 25", 26"	STD
HR61	☐	Aluminum black hard anodized ²			20", 22", 24", 25"	\$80
HR02	☐	Plastic coated ¹			20", 22", 24", 25", 26"	\$145
HR03	☐	High friction coated ^{1,2,3}			22", 24", 25", 26"	\$380
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ²	See page 6 for the color selection and price		20", 22", 24", 25", 26"	
HR28	☐	Surge ^{1,2,3,4}	(1 1/2" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"	\$425
HR05	☐	Surge LT ^{1,2,3,4}	(1 1/8" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"	\$465
HR04	☐	Natural Fit ^{1,2,3,4}			22", 24", 25", 26"	\$395
HR62	☐	Simi H ^{1,2,3}	(1 3/8" Oval handrim with high friction strip)		22", 24", 25", 26"	\$395
HR63	☐	Nova H ^{1,2,3}	(1 1/4" Oval handrim with high friction strip)		22", 24", 25", 26"	\$425
HR64	☐	Optimum H ^{1,2,3}	(1 1/8" heart-shaped handrim)		22", 24", 25", 26"	\$325
HR99	☐	Omit handrim				NC
1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive; 4- 25" and 26" are not available in narrow position; 24" is not available with the Spinerger Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece						
HANDRIMS OPTION						
ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	Full tab			STD
HR030	☐	Narrow Handrim position no cut ⁶	Full tab			NC
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions			\$32
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐ Short tab 1 position ⁶ ☐			\$32
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip				\$95
HR27	☐	Super grip				\$150
HR099	☐	Omit thumb piece				NC
5- Only available with Aluminum and Air Grip handrims; 6- Not available with the combination of 24" Newton Gravity, Spinerger Spox, LX, wheels and a 1 3/8" Airless insert pneumatic tire with any handrim;						
AXLE						
AX01	☐	Quick release axle ¹				STD
AX02	☐	Permanent axle				NC
AX03	☐	Quad release axle ¹				\$135
AXLE PLATE						
AX04	☐	Amputee adjustable axle plate ²				\$210
1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down a nd the swing away anti-tipper						
REAR WHEEL CAMBER						
CM00	☐	0° camber				STD
CM03	☐	3° camber				NC
CM06	☐	6° camber (not available with amputee axle plate)				\$165

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

X =	Most stable								Most responsive			
	2,5 1"	3,2 1/4"	3,8 1/2"	4,4 3/4"	5,1 2"	5,7 1/4"	6,3 1/2"	7 3/4"	7,6 3"	8,3 1/4"	8,9 1/2"	9,5 3/4"
	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 3,8 cm (1.5") Possible conflict between a seat depth under 40.6 cm (16") and a center of gravity over 6.3 cm (2.5"). Possible conflict between a 26" rear wheel and a center of gravity over 6.3 cm (2.5"). Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

WHEEL LOCK

WL08	☐	Newton Wheel lock - push to lock			STD
WL02	☐	Pull to lock			NC
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)			\$165
WL13 WL14	☐	6" removable extension handle ¹ *VP01			\$83
WL05	☐	Scissors Lock			\$155
WL06	☐	Attendant Lock (additional set of wheel lock mounted to the back)			\$215
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐	\$335
WL99	☐	Omit wheel locks			NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

FOOTREST

SWING-IN/SWING-OUT FOOTREST

FR03-01	☐	60° Swing-in/swing-out footrest	34.29 - 43.18cm (13.5" - 17")	Set at*:	NC
FR03-02	☐	60° Swing-in/swing-out footrest	36.83 - 49.53cm (14.5" - 19.5")	Set at*:	NC
FR03-03	☐	60° Swing-in/swing-out footrest	39.29 - 52.07cm (15.5" - 20.5")	Set at*:	NC
FR04-01	☐	70° Swing-in/swing-out footrest	34.29 - 43.18cm (13.5" - 17")	Set at*:	NC
FR04-02	☐	70° Swing-in/swing-out footrest	36.83 - 49.53cm (14.5" - 19.5")	Set at*:	NC
FR04-03	☐	70° Swing-in/swing-out footrest	39.37 - 52.07cm (15.5" - 20.5")	Set at*:	NC
FR15-01	☐	80° Swing-in/swing-out footrest ⁵	34.29 - 40.64cm (13.5" - 16")	Set at*:	NC
FR15-02	☐	80° Swing-in/swing-out footrest ⁵	36.83 - 49.53cm (14.5" - 19.5")	Set at*:	NC
FR15-03	☐	80° Swing-in/swing-out footrest ⁵	39.37 - 52.07cm (15.5" - 20.5")	Set at*:	NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	19.05 - 27.94cm (7.5" - 11")	Set at*:	NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1,2}	24.13 - 33.02cm (9.5" - 13")	Set at*:	NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1,2}	24.13 - 36.83cm (9.5" - 14.5")	Set at*:	NC

FOOTREST OPTIONS

FR09	☐	High mount (footplate attached to the hanger with clamps) ³	12.70 - 34.29cm (5" - 13.5")	Set at*:	\$85
FR10	☐	Elevating legrest with calf pad, length adjustable ⁴	Left ☐ Right ☐	Set at:	\$248
FR99	☐	Omit Footrest			NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest;

* +/- 1.27 cm (1/2")

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES

FP01	☐	Standard - Composite ¹			NC
FP02	☐	Newton adjustable angle - Composite ^{**VP01}			\$230
FP03	☐	Newton adjustable angle - Aluminum			\$255
FP04	☐	Newton adjustable angle - Carbon fiber with rubber inserts			\$450
FP23	☐	Newton adjustable angle extra deep - Aluminum			\$235
FP06	☐	One piece adjustable angle flip-up - Aluminum			\$205
FP07	☐	One piece adjustable angle flip-up - Carbon fiber with rubber inserts			\$470

FOOTPLATE ACCESSORIES

FP08	☐	Length adjustable heel loop ²			\$51
FP09	☐	Calf strap			\$43
FP30	☐	Bodypoint Aeromesh Padded calf strap			\$125

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT

FR27	☐	Residual limb support - flat 15.24x20.32cm (6"x8")	Left ☐	Right ☐	\$440
FR28	☐	Residual limb support - flat 20.32x24.4cm (8"x10")	Left ☐	Right ☐	\$440
FR29	☐	Residual limb support - flat 25.4x25.4cm (10"x10")	Left ☐	Right ☐	\$440
FR30	☐	Residual limb support - flat 25.4x35.56cm (10"x14")	Left ☐	Right ☐	\$440
FR32	☐	Residual limb support - contour 15.24x20.32cm (6"x8")	Left ☐	Right ☐	\$440
FR33	☐	Residual limb support - contour 20.32x24.4cm (8"x10")	Left ☐	Right ☐	\$440

SEATING AND BACKREST								
BACK UPHOLSTERY								
BH07	☐	Slip-on Active-Mesh back upholstery					NC	
BH08	☐	Adjustable tension Active-Mesh back upholstery					\$190	
BH99	☐	Omit Back upholstery					STD	
	☐	NXT back Support					See price list	
SEAT CUSHION								
SC01	☐	5.08 cm (2") Cushion					\$194	
SC02	☐	7.62 cm (3") Cushion					\$200	
SC99	☐	Omit seat cushion					STD	
	☐	NXT cushion					See price list	
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE <i>(push handle included - no further selection needed)</i>								
BC58-915	☐	Straight aluminum back cane	22.86 - 38.10cm (9" - 15")	<i>increment of 1/2"</i>	Set at ² :		NC	
BC58-1521	☐	Straight aluminum back cane	38.10 - 53.34cm (15" - 21")	<i>increment of 1/2"</i>	Set at ² :		NC	
BC60	☐	Straight aluminum back cane	53.34 - 60.96cm (21" - 24")	<i>increment of 1"</i>	Set at ² :		\$78	
BC01-16	☐	8° bend aluminum back cane	40.64cm (16")				NC	
BC01-18	☐	8° bend aluminum back cane	45.72cm (18")				NC	
BC01-20	☐	8° bend aluminum back cane	50.08cm (20")				NC	
BC59	☐	8° bend aluminum back cane	53.34 - 60.96cm (21" - 24")	<i>increment of 1"</i>	Set at ² :		\$78	
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>								
BC11-912	☐	Straight aluminum back cane	22.86 - 30.48cm (9" - 12")	<i>increment of 1/2"</i>	Set at ² :		NC	
BC11-1215	☐	Straight aluminum back cane	30.48 - 38.10cm (12" - 15")	<i>increment of 1/2"</i>	Set at ² :		NC	
BC11-1518	☐	Straight aluminum back cane	38.10 - 45.72cm (15" - 18")	<i>increment of 1/2"</i>	Set at ² :		NC	
BC11-1821	☐	Straight aluminum back cane	45.72 - 53.34cm (18" - 21")	<i>increment of 1/2"</i>	Set at ² :		NC	
BC12-1315	☐	8° bend aluminum back cane	33.02 - 38.10cm (13.5" - 15")	<i>increment of 1/2"</i>	Set at ² :		NC	
BC12-1518	☐	8° bend aluminum back cane	38.10 - 45.72cm (15" - 18")	<i>increment of 1/2"</i>	Set at ² :		NC	
BC12-1821	☐	8° bend aluminum back cane	45.72 - 53.34cm (18" - 21")	<i>increment of 1/2"</i>	Set at ² :		NC	
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>								
BC13-1518	☐	Straight aluminum back cane	38.10 - 45.72cm (15" - 18")	<i>increment of 1/2"</i>	Set at ² :		\$275	
BC13-1821	☐	Straight aluminum back cane	45.72 - 53.34cm (18" - 21")	<i>increment of 1/2"</i>	Set at ² :		\$275	
BC14-1518	☐	8° bend aluminum back cane	38.10 - 45.72cm (15" - 18")	<i>increment of 1/2"</i>	Set at ² :		\$275	
BC14-1821	☐	8° bend aluminum back cane	45.72 - 53.34cm (18" - 21")	<i>increment of 1/2"</i>	Set at ² :		\$275	
		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐
PUSH HANDLES OPTIONS <i>(ONE option at a time)</i>								
PH01	☐	Push handles					STD	
PH05	☐	Fold down push handles ^{3, 5}					\$150	
PH06-920	☐	clamp on height adjustable push handles 22.86 - 50.8cm (9" - 20") ^{3, 4, 5}					\$330	
PH06-923	☐	clamp on height adjustable push handles 22.86 - 58.4cm (9" - 23") ^{3, 4, 5}					\$330	
PH02	☐	Omit push handles					NC	
STROLLER HANDLE AND STABILIZER BAR <i>(only available with push handle)</i>								
PH04	☐	Stroller handle <i>(not available with the HD option)</i>					\$175	
PH03	☐	Newton folding stabilizer bar <i>(not available with stroller handle)</i>					\$275	

1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265 lb weight limit; 4- Not compatible with swing away armrests, 8° bend back canes, 9-12" and 21-24" back canes and amputee axle plates; 5- Not compatible with stroller handle, the stabilizer bar, headrest and attendant lock

ARMRESTS							
HEIGHT ADJUSTABLE FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ^{1, 4} *VP01							
AR01-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐	Unlocked STD ☐		\$390
AR01-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐	Unlocked STD ☐		\$390
AR02-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")	Locked ☐	Unlocked STD ☐		\$390
AR02-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")	Locked ☐	Unlocked STD ☐		\$390
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01							
AR13-610	☐	25.4cm (10") Desk length armpad ⁵	15.25 - 25.40cm (6" - 10")				\$360
AR03-812	☐	25.4cm (10") Desk length armpad	20.32 - 30.48cm (8" - 12")				\$360
AR03-1014	☐	25.4cm (10") Desk length armpad	25.40 - 35.56cm (10" - 14")				\$360
AR14-610	☐	35.56cm (14") Full length armpad ⁵	15.25 - 25.40cm (6" - 10")				\$360
AR04-812	☐	35.56cm (14") Full length armpad	20.32 - 30.48cm (8" - 12")				\$360
AR04-1014	☐	35.56cm (14") Full length armpad	25.40 - 35.56cm (10" - 14")				\$360
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²							
AR05	☐	Tubular armrest - Short pad	20.32 - 27.94cm (8 - 11")				\$250
AR23	☐	Tubular armrest - Long pad	20.32 - 27.94cm (8 - 11")				\$250
AR06	☐	25.4cm (10") Desk Length flat armpad	20.32 - 27.94cm (8 - 11")				\$250
AR17	☐	35.56cm (14") Full Length flat armpad	20.32 - 27.94cm (8 - 11")				\$250
ARMREST OPTIONS							
AR12	☐	35.56cm (14") Full Length Armpad - Gel Ovation <i>(pair)</i> ³					\$165
AR99	☐	Omit armrest					NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards; 3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

ACCESSORIES							
SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>							
SG01	☐	Plastic side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹			\$210
SG02	☐	Carbon fiber side guards <i>(fixed)</i>					\$275
SG03	☐	Carbon fiber side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹			\$315
SG06	☐	Plastic fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹			\$295
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹			\$550
SG99	☐	Omit side guard					NC

1- only compatible with 0° camber; Increase width by 1/2" per side for a total of 1" ;

POSITIONING BELT							
SB01-S	☐	Velcro adjustable - Small 88.9cm (35") ¹					\$50
SB01-L	☐	Velcro adjustable - Long 152.40cm (60")					\$50
SB02-S	☐	Auto buckle - Small 88.9cm (35") ¹					\$95
SB02-M	☐	Auto buckle - Medium 121.92cm (48")					\$95
SB02-L	☐	Auto buckle - Long 152.40cm (60")					\$95
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 142.24cm (56")					\$95
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 154.94cm (61")					\$95
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 129.54cm (51")					\$185
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 142.24cm (56")					\$185
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 154.94cm (61")					\$185
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 129.54cm (51")					\$205
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 142.24cm (56")					\$205
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 154.94cm (61")					\$205
SB03	☐	Hardware attachment for 2 extra points (1 pair)					\$17
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 73.66 (18")					\$310
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 93.98cm (24")					\$310
SB99	☐	Omit positioning belt					STD

1- Not recommended for seat width of 45.72cm (18") and more

HEADREST ¹							
AC01-S	☐	Profiled Headrest - Small 8.8x15.2cm (3.5"x6")					\$195
AC01-M	☐	Profiled Headrest - Medium 11.4x17.8cm (4.5"x7")					\$195
AC01-L	☐	Profiled Headrest - Large 14x2.3cm(5.5"x8")					\$195
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar ¹					\$320

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES							
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber					\$35
AC03	☐	Anti-tippers					\$155
AC29	☐	Swing-away anti-tipper ¹ Left ☐ Right ☐					\$165
AC27	☐	Tip assist ²					\$67
AC06	☐	Cane and crutch holder Left ☐ Right ☐					\$130
AC07	☐	Oxygen tank holder ⁴ - D and E Type - Installed to the left					\$335
AC08	☐	Spoke guards ³ Available: 20", 22", 24", 25", 26"					\$125
AC09	☐	Footrest neoprene impact guard (pair) - Black					\$100
AC10	☐	Seat pouch					\$83
AC30	☐	Newton multi-purpose pouch					\$62
AC26	☐	BodyPoint Mobility Bag					\$135
AC23	☐	Newton Backpack					\$130
AC12	☐	Velcro Closing Strap					\$57
AC13	☐	Tool kit					\$50

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 34,3cm (13 1/2") and 50,2cm (19 3/4"); Not available with amputee axle plate

2- Cannot be combined with Anti Tippers; 3- Not available with Mag wheels and Spinergy; the 24" spoke guard will be installed on the 25" wheel;

4- Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down


Colors																			
Options ¹		Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
	Black																		
AirGrip	☐ \$125	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame ² (STD)	☐ STD	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
FC99	☐	Custom Color Specify PPG Paintit color: additional 10 days lead-time																	\$390

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen; Finish between Airgrip and frame might be different;

2- with black hydroformed aluminum crossbrace

3- The frame color will be confirmed with the quote conversion into an order (in absence of selection, the color will be black)

Fork and Stem bolt combinations



Frame length chart

Frame length	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
STD							
+1							N/A
+2						N/A	N/A
-1	N/A						

Unilateral wheel lock restrictions

Seat depth: shall be 40.64 cm and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° camber		
One arm drive: not compatible with "T"		
One arm drive and flip-back armrest only available with 17 in depth frame		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (cm)	COG Compatibility Range (cm)
20 in	40.64+	3.8 to 9.5
	40.64	2.5 to 7
22 in	43.18+	2.5 to 8.3
	45.72+	2.5 to 9.5
24 in	40.64	2.5 to 5.1
	43.18	2.5 to 7
	45.72+	2.5 to 9.5



CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the
EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.




Have a customization in mind or a specific request? The **EVO Custom Request Form** is the best place to start.

Notes :



FOLDING ON CHAIR – INTL ORDER FORM

BioFit™ CUSHION		WIDTH											
 Use wheelchair configuration Or use this size: NXCABF _____ (W" x D")	DEPTH	12"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
		12"	\$425	\$425	\$425		N/A	N/A	N/A	N/A	N/A	N/A	N/A
		13"						N/A	N/A	N/A	N/A	N/A	N/A
		14"	\$425	\$425	\$425								
		15"		\$425	\$425	\$425	\$425						
		16"		\$425	\$425	\$425	\$425	\$425	\$425	\$460	\$460	\$460	\$460
		17"	N/A			\$425	\$425	\$425	\$425	\$460	\$460		
		18"	N/A			\$425	\$425	\$425	\$425	\$460	\$460	\$460	\$460
		19"	N/A			\$425	\$425	\$425	\$425	\$460	\$460		
		20"	N/A			\$425	\$425	\$425	\$460	\$460	\$460	\$460	\$460
Option: Change cover for incontinent/infection control outer cover												N/C	

Green boxes are stocked sizes; the red boxes require customization*; the N/A boxes are not available sizes.
*Customization (red boxes) is available with an upcharge of \$160.00 and has a 2 to 4 weeks additional lead time.

SEATING AND CUSHION ACCESSORIES

	MODEL #	Description	Fits Back Width	Retail	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	\$77	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	\$77	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	\$77	



Optima™

THORACIC BACK SUPPORT

THORACIC SHAPE

WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOACB1413	\$695		NIB-2193V	\$125	
14"	16"	NOACB1416	\$695		NIB-2194V	\$125	
15"	13"	NOACB1513	\$695		NIB-2197V	\$125	
15"	16"	NOACB1516	\$695		NIB-2198V	\$125	
16"	13"	NOACB1613	\$695		NIB-2201V	\$125	
16"	16"	NOACB1616	\$695		NIB-2202V	\$125	
16"	19"	NOACB1619	\$695		NIB-2203V	\$125	
17"	13"	NOACB1713	\$695		NIB-2206V	\$125	
17"	16"	NOACB1716	\$695		NIB-2207V	\$125	
18"	13"	NOACB1813	\$695		NIB-2208V	\$125	
18"	16"	NOACB1816	\$695		NIB-2211V	\$125	
18"	19"	NOACB1819	\$695		NIB-2212V	\$125	
19"	16"	NOACB1916	\$695		NIB-2213V	\$125	
19"	19"	NOACB1919	\$695		NIB-2215V	\$125	
20"	13"	NOACB2013	\$695		NIB-2218V	\$125	
20"	16"	NOACB2016	\$695		NIB-2219V	\$125	
20"	19"	NOACB2019	\$695		NIB-2220V	\$125	
21"	16"	NOACB2116	\$695		NIB-2221V	\$125	
21"	19"	NOACB2119	\$695		NIB-2222V	\$125	
22"	13"	NOACB2213	\$745		NIB-2223V	\$125	
22"	16"	NOACB2216	\$745		NIB-2224V	\$125	
22"	19"	NOACB2219	\$745		NIB-2225V	\$125	

PLEASE NOTE: Each back will fit widths up to 2" wider.



Optima Carbon™

THORACIC BACK SUPPORT

CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)

WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover &	Retail	✓
14"	13"	NOCFACB1413	\$900		NIB-2193V	\$125	
14"	16"	NOCFACB1416	\$930		NIB-2194V	\$125	
15"	13"	NOCFACB1513	\$900		NIB-2197V	\$125	
15"	16"	NOCFACB1516	\$930		NIB-2198V	\$125	
16"	13"	NOCFACB1613	\$900		NIB-2201V	\$125	
16"	16"	NOCFACB1616	\$930		NIB-2202V	\$125	
17"	13"	NOCFACB1713	\$900		NIB-2206V	\$125	
17"	16"	NOCFACB1716	\$930		NIB-2207V	\$125	
18"	13"	NOCFACB1813	\$900		NIB-2211V	\$125	
18"	16"	NOCFACB1816	\$930		NIB-2212V	\$125	
20"	13"	NOCFACB2013	\$900		NIB-2218V	\$125	
20"	16"	NOCFACB2016	\$930		NIB-2219V	\$125	

ALL PRICES QUOTED IN USD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com



Optima™ DEEP

THORACIC BACK SUPPORT

WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOPCB1413	\$695		NIB-2246V	\$125	
14"	16"	NOPCB1416	\$695		NIB-2247V	\$125	
15"	13"	NOPCB1513	\$695		NIB-2250V	\$125	
15"	16"	NOPCB1516	\$695		NIB-2251V	\$125	
16"	13"	NOPCB1613	\$695		NIB-2254V	\$125	
16"	16"	NOPCB1616	\$695		NIB-2255V	\$125	
16"	19"	NOPCB1619	\$695		NIB-2256V	\$125	
17"	13"	NOPCB1713	\$695		NIB-2259V	\$125	
17"	16"	NOPCB1716	\$695		NIB-2260V	\$125	
18"	13"	NOPCB1813	\$695		NIB-2264V	\$125	
18"	16"	NOPCB1816	\$695		NIB-2265V	\$125	
18"	19"	NOPCB1819	\$695		NIB-2266V	\$125	
19"	16"	NOPCB1916	\$695		NIB-2267V	\$125	
19"	19"	NOPCB1919	\$695		NIB-2268V	\$125	
20"	13"	NOPCB2013	\$695		NIB-2271V	\$125	
20"	16"	NOPCB2016	\$695		NIB-2272V	\$125	
20"	19"	NOPCB2019	\$695		NIB-2273V	\$125	
21"	16"	NOPCB2116	\$695		NIB-2274V	\$125	
21"	19"	NOPCB2119	\$695		NIB-2275V	\$125	
22"	13"	NOPCB2213	\$745		NIB-2276V	\$125	
22"	16"	NOPCB2216	\$745		NIB-2277V	\$125	
22"	19"	NOPCB2219	\$745		NIB-2278V	\$125	

PLEASE NOTE: Each back will fit widths up to 2" wider.



Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

xtend™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NHAACB(L)1213/17	12"	13-17"	\$655		NIB-2148V	\$125	
NHAACB(L)1413/17	14"	13-17"	\$655		NIB-2149V	\$125	
NHAACB(L)1613/17	16"	13-17"	\$655		NIB-2150V	\$125	
NHAACB(L)1813/17	18"	13-17"	\$655		NIB-2151V	\$125	
NHAACB (L)2013/17	20"	13-17"	\$655		NIB-2152V	\$125	
NHAACB (L)2213/17	22"	13-17"	\$825		NIB-2153V	\$125	

PLEASE NOTE: Each back will fit widths up to 2" wider.



Classic™ PELVIC BACK SUPPORT

Classic™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NACB1616	16"	16"	\$655		NIB-2163V	\$125	
NACB1618	16"	18"	\$655		NIB-2164V	\$125	
NACB1816	18"	16"	\$655		NIB-2166V	\$125	
NACB1818	18"	18"	\$655		NIB-2167V	\$125	
NACB1820	18"	20"	\$655		NIB-2168V	\$125	
NACB2020	20"	20"	\$655		NIB-2171V	\$125	
NACB2216	22"	16"	\$825		NIB-2172V	\$125	
NACB2218	22"	18"	\$825		NIB-2173V	\$125	
NACB2220	22"	20"	\$825		NIB-2174V	\$125	

PLEASE NOTE: Each back will fit widths up to 2" wider.



Optima Carbon™ DEEP

THORACIC BACK SUPPORT

CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)

WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover &	Retail	✓
14"	13"	NOCFPCB1413	\$985		NIB-2246V	\$125	
14"	16"	NOCFPCB1416	\$1,025		NIB-2247V	\$125	
15"	13"	NOCFPCB1513	\$985		NIB-2250V	\$125	
15"	16"	NOCFPCB1516	\$1,025		NIB-2251V	\$125	
16"	13"	NOCFPCB1613	\$985		NIB-2254V	\$125	
16"	16"	NOCFPCB1616	\$1,025		NIB-2255V	\$125	
17"	13"	NOCFPCB1713	\$985		NIB-2259V	\$125	
17"	16"	NOCFPCB1716	\$1,025		NIB-2260V	\$125	
18"	13"	NOCFPCB1813	\$985		NIB-2264V	\$125	
18"	16"	NOCFPCB1816	\$1,025		NIB-2265V	\$125	
20"	13"	NOCFPCB2013	\$985		NIB-2271V	\$125	
20"	16"	NOCFPCB2016	\$1,025		NIB-2272V	\$125	



Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

xtend™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover &	Retail	✓
NHAPCB1616/20	16"	16"-20"	\$655		NIB-2143V	\$125	
NHAPCB1816/20	18"	16"-20"	\$655		NIB-2144V	\$125	
NHAPCB2016/20	20"	16"-20"	\$655		NIB-2145V	\$125	
NHAPCB2216/20	22"	16"-20"	\$825		NIB-2146V	\$125	



Classic™ DEEP PELVIC BACK SUPPORT

Classic™ Deep	WIDTH	HEIGHT	RETAIL	✓	repl. cover & feet	Retail	✓
NPCB1616	16"	16"	\$655		NIB-2175V	\$125	
NPCB1618	16"	18"	\$655		NIB-2176V	\$125	
NPCB1818	18"	18"	\$655		NIB-2179V	\$125	
NPCB1820	18"	20"	\$655		NIB-2180V	\$125	
NPCB2020	18"	20"	\$655		NIB-2183V	\$125	
NPCB2216	20"	16"	\$825		NIB-2184V	\$125	
NPCB2218	22"	18"	\$825		NIB-2185V	\$125	
NPCB2220	22"	20"	\$825		NIB-2186V	\$125	

ALL PRICES QUOTED IN USD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

BACK SUPPORT ACCESSORIES

	MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	\$160	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	\$65	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	\$65	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	\$65	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	\$65	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	\$65	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

	MODEL	Description	RETAIL	✓
	NIB-0784	10" x 5" Large Rectangular	\$195	
	NIB-0818	14" x 5" Extra-Large Rectangular	\$195	

MultiFit™ HEAD SUPPORT – FIXTURES

	MODEL	Description	RETAIL	✓
	NIB-0747	Multi-Axis Hardware	\$145	
	NIB-0823	Head Support Receiver Bracket	\$58	

EXTRA COVERS

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	NXCCBF _____ V	\$125	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	NXCCBF _____ V	\$145	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	NXCCAF _____ V	\$125	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	NXCCAF _____ V	\$145	

ALL PRICES QUOTED IN USD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com