



Helio C2

ORDER FORM - CANADA

1- Purchasing information

P.O. number: * _____ Date: _____
Company Name: * _____
Purchasing contact: * _____
Phone: * _____
Email: * _____
Contact preference: * Phone ☐ Email ☐
Marked for: _____

2- Dealer/Representative information

Dealer/Representative: * _____
Phone: * _____
Email: * _____
Contact preference: * Phone ☐ Email ☐

3- Shipping information

Ship to: * _____
Name: _____
Address: _____
City: _____ Province: _____
Postal Code: _____
Phone: _____

* Required

MA080 ☐
English Manual (STD)

MA081 ☐
Manuel Français

Taxes not included

MODEL							
WC03	☒	Helio ultralight folding heelchair (265lb weight capacity)				WM0000964	\$2,945
WC08	☒	C2 carbon fiber option				WAMZ	\$1,450
WC09	☐	HD Heavy Duty package (350 lb weight capacity)				WAMF	\$180
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6" camber					\$35
Frame							
FREE Growth Kit within the first 5 years of the Helio C2 wheelchair ownership. The growth kit includes carbon crossbraces, crossbrace links, seat sling, back upholstery and footplates.							
SEAT WIDTH ³							
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	NC
	☐	19" ¹ ☐	20" ¹ ☐	21" ² ☐	22" ² ☐	WAMA	\$360
SEAT DEPTH							
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	NC
	☐	19" ☐	20" ☐			WAMB	\$275
TRANSIT TIE DOWN							
AC04	☐	Transit Tie Down (Unoccupied)					NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 20"-41"					\$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 20"-46"					\$365
AC35	☐	WC-19 Transit Tie Down without belt ⁴					\$345

1- Available with HD option; 2- Only available with HD option

[3. Width calculator](#)

4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option



PACKAGES				
PK01	☐	Newton Ultralight Package (include: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels)	Save \$245	\$540
VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button caps, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹		\$95

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo **"VP01"**; Recommended with black chair;
See picture of the product on page 6.

C2 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2				14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						
	P4	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
7"	P1										16 1/2"	17 1/2"	18 1/2"	17"	18"	19"			
	P2							16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"			
	P3							17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
	P4	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"
	P5	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"	19 1/2"	20 1/2"	21 1/2"
Froglegs	P1				14 1/2"	15 1/2"	16 1/2"												
	P2	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"									
	P3	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"						

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS				Available		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 1/4")		6", 8"		170 \$
CA14 to CA16	☐	Newton UltraCaster Composite (1 1/2")		4", 5", 6"		105 \$
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 1/2")		4", 5", 6"		198 \$
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹			WAMS	\$275
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				\$525

1- Not available with WC-19

C2 - REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height:
☐	20" (451) wheel	12 1/2	to 17 1/4	
☐	22" (501) wheel	13 1/2	to 18 1/4	
☐	24" (540) wheel	14 1/2	to 19 1/4	
☐	25" (559) wheel	15	to 19 3/4	
☐	26" (590) wheel	15 1/2	to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/4" of selected seat to floor height"

Footrest						
SWING-IN/SWING-OUT FOOTREST						
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"		Set at*:		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*:		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*:		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"		Set at*:		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*:		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*:		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" ⁵		Set at*:		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" ⁵		Set at*:		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" ⁵		Set at*:		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5" - 11"	Set at*:		NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5" - 13"	Set at*:		NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5" - 14.5"	Set at*:		NC
FOOTREST OPTIONS ³						
FR09	☐	High mount (footplate attached to the hanger with clamps) ³	5" - 13.5"	Set at*:		\$105
FR10	☐	Elevating legrest with calf pad, Length adjustable ⁴	Left ☐ Right ☐ 14" - 22"	Set at:	WAMK	\$295
FR99	☐	Omit Footrest				NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest; * +/- 1/2"

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES						
FP01	☐	Standard - Composite ¹				STD
FP02	☐	Newton adjustable angle - Composite ^{**VP01}			WAMJ	\$155
FP03	☐	Newton adjustable angle - Aluminum			WAMJ	\$175
FP04	☐	Newton adjustable angle - Carbon fiber <i>(with rubber inserts)</i>			WAMJ	\$395
FP23	☐	Newton adjustable angle extra deep - Aluminum			K0040	\$175
FP06	☐	One piece adjustable angle flip-up - Aluminum		WAMJ	SEFND0015	\$215
FP07	☐	One piece adjustable angle flip-up - Carbon fiber <i>(with rubber inserts)</i>		WAMJ	SEFND0015	\$480
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop ²				NC
FP09	☐	Calf strap			SEFND0010	\$57
FP30	☐	Bodypoint Aeromesh Padded calf strap			SEFND0010	\$140

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	☐	Residual limb support - flat 6"x8"	Left ☐ Right ☐		SEFND0005	\$505
FR28	☐	Residual limb support - flat 8"x10"	Left ☐ Right ☐		SEFND0005	\$505
FR29	☐	Residual limb support - flat 10"x10"	Left ☐ Right ☐		SEFND0005	\$505
FR30	☐	Residual limb support - flat 10"x14"	Left ☐ Right ☐		SEFND0005	\$505
FR32	☐	Residual limb support - contour 6"x8"	Left ☐ Right ☐		SEFND0005	\$505
FR33	☐	Residual limb support - contour 8"x10"	Left ☐ Right ☐		SEFND0005	\$505

Backrest							
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>							
BC06-912	☐	Straight Carbon fiber back cane	9" - 12"	increment of 1/2"	Set at ²		NC
BC06-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ²		NC
BC06-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²		NC
BC06-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²		NC
BC08-1315	☐	8° Bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ²		\$170
BC08-1518	☐	8° Bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²		\$170
BC08-1821	☐	8° Bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²		\$170
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>							
BC09-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ²		\$330
BC09-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²		\$330
BC09-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²		\$330
BC10-1315	☐	8° bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ²		\$380
BC10-1518	☐	8° bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²		\$380
BC10-1821	☐	8° bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²		\$380
		Select back cane angle:	85° ¹ ☐ 90° ☐ 95° STD ☐ 100° ☐ 105° ☐ 110° ☐				
PUSH HANDLES OPTIONS <i>(ONE option at a time)</i>							
PH01	☐	Push handles					STD
PH05	☐	Fold down push handles ^{3,5}					\$165
PH06-920	☐	Clamp on height adjustable push handles 9"-20" ^{3,4,5}					\$345
PH06-923	☐	Clamp on height adjustable push handles 9"-23" ^{3,4,5}					\$345
PH02	☐	Omit push handles					NC
STROLLER HANDLE AND STABILIZER BAR <i>(only available with push handle)</i>							
PH04	☐	Stroller handle <i>(not available with the HD option)</i>			WAMW		\$220
PH03	☐	Newton folding stabilizer bar <i>(not available with stroller handle)</i>					\$275

1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265 lb weight limit;

4- Not compatible with swing away armrests, 8° bend back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				NC
BH08	☐	Adjustable tension Active-Mesh back upholstery			WAMD/WAME	\$260
BH99	☐	Omit Back upholstery				STD
	☐	NXT back Support			See price list	
SEAT CUSHION						
SC01	☐	2" Cushion				\$185
SC02	☐	3" Cushion				\$205
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See price list	

Rear wheels							
ONE ARM DRIVE ¹		Only available with Newton One rear wheels, permanent axles, 0° camber					
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐		WAML	\$1,480
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐		WAML	\$1,595
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)					NC
REAR WHEELS					Available		
RW01	☐	Mag ³			20", 22" 24"		\$140
RW02	☐	Newton One - spoke wheel ⁴			20", 22" 24"		NC
RW03	☐	Newton Gravity - ultralight wheel			20", 22", 24", 25", 26"		\$390
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐	22", 24", 25", 26"		\$875
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	22", 24", 25", 26"		\$1,095
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	24", 25"		\$1,650
TIRES					Available		
AIRLESS							
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black ⁹			20", 22" 24"		STD
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black ^{7,9}			22", 24", 25", 26"		\$120
TI04	☐	Soft Urethane 1" - No tread, black - <i>Shox</i>			22", 24", 25", 26"		\$140
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸			22", 24", 26"		\$225
PNEUMATIC							
TI05	☐	Pneumatic 1 3/8" - Medium Tread ⁹			20", 22", 24", 25" 26"		NC
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>			20", 22", 24", 25", 26"		\$120
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>			24", 25", 26"		\$200
TI03	☐	Pneumatic 2" - Large Tread, black ¹⁰			22", 24", 25"		\$130

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available;
3- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additional 10 days lead-time;
6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels
9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;

Handrims							
HANDRIM		Available					
HR01	☐	Aluminum anodized			20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²			20", 22", 24", 25"		\$85
HR02	☐	Plastic coated ¹			20", 22", 24", 25", 26"	WAMO	\$170
HR03	☐	High friction coated ^{1,2,3}			22", 24", 25", 26"	WAMO	\$460
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ²	See on page 6 for color choice and price		20", 22", 24", 25", 26"	WAMO	
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"		\$485
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"		\$515
HR04	☐	Natural Fit ^{1,2,3,4}			22", 24", 25", 26"		\$320
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)			22", 24", 25", 26"		\$485
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)			22", 24", 25", 26"		\$515
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)			22", 24", 25", 26"		\$320
HR99	☐	Omit handrim					NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;
4- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

Handrims option							
ASSEMBLY POSITION							
HR29	☐	Standard Handrim position	Full tab				STD
HR030	☐	Narrow Handrim position no cut	Full tab				NC
HR30	☐	Narrow Handrim position with cut ³	Short tab 2 positions				\$29
HR31	☐	Super Narrow Handrim ³	Short tab 2 positions STD ☐ Short tab 1 position ☐				\$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS							
HR23	☐	Standard grip					\$200
HR27	☐	Super grip					\$230
HR099	☐	Omit thumb piece					NC


3- Only available with Aluminum and Air Grip handrims;

AXLE ¹							
AX01	☐	Quick release axle				WAMU	STD
AX02	☐	Permanent axle					NC
AX03	☐	Quad release axle				WAMU	\$180
AXLE PLATE							
AX04	☐	Amputee adjustable axle plate ²				WAMT	\$275

1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

REAR WHEEL CAMBER							
CM00	☐	0° camber					STD
CM03	☐	3° camber					NC
CM06	☐	6° camber <i>(not available with amputee axle plate)</i>					\$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable														
X =	AMPUTEE	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

WHEEL LOCK							
WL08	☐	Newton Wheel lock - push to lock				PRODUCT IMPROVED	STD
WL16	☐	Aluminum push to lock with extensions					\$73
WL02	☐	Pull to lock					NC
WL04	☐	Newton grade aid push to lock <i>(works best with pneumatic or high pressure tire)</i>					WAMR \$185
WL13 WL14	☐	6" removable extension handle ¹ *VP01					\$73
WL05	☐	Scissors Lock					\$185
WL06	☐	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>					\$240
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐			WAMN \$335
WL99	☐	Omit wheel locks					NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests

HEIGHT ADJUSTABLE FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ^{1,4} *VP01							
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		\$395
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		\$395
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		\$395
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		\$395
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01							
AR13-610	☐	10" Desk length armpad ⁵	6" - 10"				\$380
AR03-812	☐	10" Desk length armpad	8" - 12"				\$380
AR03-1014	☐	10" Desk length armpad	10" - 14"				\$380
AR14-610	☐	14" Full length armpad ⁵	6" - 10"				\$380
AR04-812	☐	14" Full length armpad	8" - 12"				\$380
AR04-1014	☐	14" Full length armpad	10" - 14"				\$380
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²							
AR05	☐	Tubular armrest - short pad	8" - 11"				\$315
AR23	☐	Tubular armrest - long pad	8" - 11"				\$315
AR06	☐	10" Desk length flat armpad	8" - 11"				\$315
AR17	☐	14" Full length flat armpad	8" - 11"				\$315
ARMRESTS OPTIONS							
AR12	☐	14" Full Length armpad - Gel Ovation (pair) ³				SEAND0030	\$160
AR99	☐	Omit armrest					NC


1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>							
SG01	☐	Plastic side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		WAMV	\$235
SG02	☐	Carbon fiber side guards <i>(fixed)</i>				WAMV	\$285
SG03	☐	Carbon fiber side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		WAMV	\$320
SG06	☐	Plastic fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		WAMV	\$340
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		WAMV	\$510
SG99	☐	Omit side guard					NC

1- only compatible with 0° camber; Increase width by 1/2" per side for a total of 1" ;

Fork and Stem bolt combinations



Vision Package Accents



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
22 in	16	1.0 to 2.75
	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :

FOLDING ON CHAIR – CDN ORDER FORM

EXTRA COVERS

Extra Outer smartx^{3D} Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	\$	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	__ __ X __ __	130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	__ __ X __ __	150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	__ __ X __ __	130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	__ __ X __ __	150	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com