



Helio C2

ORDER FORM - CANADA

1- Purchasing information

P.O. number: * _____ Date: _____

Company Name: * _____

Purchasing contact: * _____

Phone: * _____

Email: * _____

Contact preference: * Phone ☐ Email ☐

Marked for: _____

2- Dealer/Representative information

Dealer/Representative: * _____

Phone: * _____

Email: * _____

Contact preference: * Phone ☐ Email ☐

3- Shipping information

Ship to: * _____

Name: _____

Address: _____

City: _____ Province: _____

Postal Code: _____

Phone: _____

* Required

MA080 ☐
English Manual (STD)

MA081 ☐
Manuel Français

Taxes not included

MODEL								
WC03	☒	Helio ultralight folding wheelchair (265lb weight capacity)						\$2,995
WC08	☒	C2 carbon fiber option						\$1,500
WC09	☐	HD Heavy Duty package (350 lb weight capacity)						\$205
FRAME								
FREE Growth Kit within the first 5 years of the Helio C2 wheelchair ownership. The growth kit includes carbon crossbraces, crossbrace links, seat sling, back upholstery and footplates.								
SEAT WIDTH ³								
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐		NC
	☐	19" ¹ ☐	20" ¹ ☐	21" ² ☐	22" ² ☐			\$365
SEAT DEPTH								
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐		NC
	☐	19" ☐	20" ☐					\$280
FRAME LENGTH ⁵								
SP01	☐	Standard frame length configuration						STD
SP02	☐	Frame length 1" longer than standard configuration						NC
SP03	☐	Frame length 2" longer than standard configuration						NC
SP04	☐	Frame length 1" shorter than standard configuration						NC
SP05	☐	Frame length 2" shorter than standard configuration						NC
TRANSIT TIE DOWN								
AC04	☐	Transit Tie Down (Unoccupied)						NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 20"-41"						\$370
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 20"-46"						\$370
AC35	☐	WC-19 Transit Tie Down without belt ⁴						\$345

1- Available with HD option; 2- Only available with HD option

4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

5- See chart on page 7

3. Width calculator



PACKAGES				
PK01	☐	Newton Ultralight Package (include: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels)		Save \$245 \$540
VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button caps, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹		\$95

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo *VP01; Recommended with black chair;
See picture of the product on page 6.

C2 FRONT Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2				14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						
	P4	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
7"	P1										16 1/2"	17 1/2"	18 1/2"	17"	18"	19"			
	P2							16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"			
	P3							17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
	P4	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"
	P5	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"	19 1/2"	20 1/2"	21 1/2"
Froglegs	P1				14 1/2"	15 1/2"	16 1/2"												
	P2	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"									
	P3	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"						

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

FRONT CASTERS

CASTERS				Available		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 ¼")		6", 8"		170 \$
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		120 \$
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		198 \$
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹				\$275
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				\$525

1- Not available with WC-19

C2 - REAR Seat-to-floor Height

WHEEL SIZES		Height available	Specify desired height: _____
☐	20" (451) wheel	12 1/2 to 17 1/4	
☐	22" (501) wheel	13 1/2 to 18 1/4	
☐	24" (540) wheel	14 1/2 to 19 1/4	
☐	25" (559) wheel	15 to 19 3/4	
☐	26" (590) wheel	15 1/2 to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

REAR WHEELS						
ONE ARM DRIVE ¹		Only available with Newton One rear wheels, permanent axles, 0° camber				
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐		\$1,480
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐		\$1,595
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)				NC
REAR WHEELS					Available	
RW01	☐	Mag ³			20", 22" 24"	\$150
RW02	☐	Newton One - spoke wheel ⁴			20", 22" 24"	NC
RW03	☐	Newton Gravity - ultralight wheel			20", 22", 24", 25", 26"	\$390
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐	22", 24", 25", 26"	\$915
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	22", 24", 25", 26"	\$1,135
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐ Wht ⁵ ☐ Blue ⁵ ☐ Red ⁵ ☐ Yel ⁵ ☐ Pink ⁵ ☐ Grn ⁵ ☐ Org ⁵ ☐	24", 25"	\$1,690
TIRES					Available	
AIRLESS						
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black ⁹			20", 22" 24"	STD
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black ^{7,9}			22", 24", 25", 26"	\$120
TI04	☐	Soft Urethane 1" - No tread, black - <i>Shox</i>			22", 24", 25", 26"	\$150
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸			22", 24", 26"	\$235
PNEUMATIC						
TI05	☐	Pneumatic 1 3/8" - Medium Tread ⁹			20", 22", 24", 25" 26"	NC
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>			20", 22", 24", 25", 26"	\$130
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>			24", 25", 26"	\$210
TI03	☐	Pneumatic 2" - Large Tread, black ¹⁰			22", 24", 25"	\$140

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available;

3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additional 10 days lead-time;

6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 25" and 26" are not available with Newton Gravity wheels

8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels 9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;

HANDRIMS						
HANDRIM		Available				
HR01	☐	Aluminum anodized			20", 22", 24", 25", 26"	STD
HR61	☐	Aluminum black hard anodized ²			20", 22", 24", 25"	\$95
HR02	☐	Plastic coated ¹			20", 22", 24", 25", 26"	\$180
HR03	☐	High friction coated ^{1,2,3}			22", 24", 25", 26"	\$460
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ²	See on page 6 for color choice and price			
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"	\$485
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"	\$515
HR04	☐	Natural Fit ^{1,2,3,4}			22", 24", 25", 26"	\$320
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)			22", 24", 25", 26"	\$485
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)			22", 24", 25", 26"	\$515
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)			22", 24", 25", 26"	\$320
HR99	☐	Omit handrim				NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

4 - 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

HANDRIMS OPTION						
ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	Full tab			STD
HR030	☐	Narrow Handrim position no cut	Full tab			NC
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions			\$29
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐ Short tab 1 position ☐			\$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip				\$200
HR27	☐	Super grip				\$230
HR099	☐	Omit thumb piece				NC


5- Only available with Aluminum and Air Grip handrims;

AXLE ¹						
AX01	☐	Quick release axle				STD
AX02	☐	Permanent axle				NC
AX03	☐	Quad release axle				\$180
AXLE PLATE						
AX04	☐	Amputee adjustable axle plate ²				\$280

1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

REAR WHEEL CAMBER							
CM00	☐	0° camber					STD
CM03	☐	3° camber					NC
CM06	☐	6° camber (not available with amputee axle plate)					\$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													
AMPUTEE	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

X =

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

WHEEL LOCK							
WL08	☐	Newton Wheel lock - push to lock					STD
WL16	☐	Aluminum push to lock with extensions					\$85
WL02	☐	Pull to lock					NC
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)					\$200
WL13	☐	6" removable extension handle ¹ *VP01					\$75
WL14	☐						
WL05	☐	Scissors Lock					\$200
WL06	☐	Attendant Lock (additional set of wheel lock mounted to the back)					\$240
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐			\$335
WL99	☐	Omit wheel locks					NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

FOOTREST

SWING-IN/SWING-OUT FOOTREST							
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"			Set at*:		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"			Set at*:		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"			Set at*:		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"			Set at*:		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"			Set at*:		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"			Set at*:		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" ⁵			Set at*:		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" ⁵			Set at*:		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" ⁵			Set at*:		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹		7.5" - 11"	Set at*:		NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}		9.5" - 13"	Set at*:		NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}		9.5" - 14.5"	Set at*:		NC
FOOTREST OPTIONS ³							
FR09	☐	High mount (footplate attached to the hanger with clamps) ³		5" - 13.5"	Set at*:		\$105
FR10	☐	Elevating legrest with calf pad, Length adjustable ⁴	Left ☐	Right ☐	Set at:		\$295
FR99	☐	Omit Footrest					NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest; * +/- 1/2"

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES							
FP01	☐	Standard - Composite ¹					STD
FP02	☐	Newton adjustable angle - Composite *VP01					\$155
FP03	☐	Newton adjustable angle - Aluminum					\$175
FP04	☐	Newton adjustable angle - Carbon fiber (with rubber inserts)					\$395
FP23	☐	Newton adjustable angle extra deep - Aluminum					\$175
FP06	☐	One piece adjustable angle flip-up - Aluminum					\$215
FP07	☐	One piece adjustable angle flip-up - Carbon fiber (with rubber inserts)					\$480
FOOTPLATE ACCESSORIES							
FP08	☐	Length adjustable heel loop ²					NC
FP09	☐	Calf strap					\$62
FP30	☐	Bodypoint Aeromesh Padded calf strap					\$145

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT								
FR27	☐	Residual limb support - flat 6"x8"	Left ☐	Right ☐				\$505
FR28	☐	Residual limb support - flat 8"x10"	Left ☐	Right ☐				\$505
FR29	☐	Residual limb support - flat 10"x10"	Left ☐	Right ☐				\$505
FR30	☐	Residual limb support - flat 10"x14"	Left ☐	Right ☐				\$505
FR32	☐	Residual limb support - contour 6"x8"	Left ☐	Right ☐				\$505
FR33	☐	Residual limb support - contour 8"x10"	Left ☐	Right ☐				\$505

SEATING AND BACKREST

BACK UPHOLSTERY								
BH07	☐	Slip-on Active-Mesh back upholstery						NC
BH08	☐	Adjustable tension Active-Mesh back upholstery						\$265
BH99	☐	Omit Back upholstery						STD
	☐	NXT back Support				See price list		
SEAT CUSHION								
SC01	☐	2" Cushion						\$185
SC02	☐	3" Cushion						\$205
SC99	☐	Omit seat cushion						STD
	☐	NXT cushion				See price list		
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)								
BC06-912	☐	Straight Carbon fiber back cane	9" - 12"	increment of 1/2"	Set at ²			NC
BC06-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ²			NC
BC06-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²			NC
BC06-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²			NC
BC08-1315	☐	8° Bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ²			\$170
BC08-1518	☐	8° Bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²			\$170
BC08-1821	☐	8° Bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²			\$170
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)								
BC09-1215	☐	Straight Carbon fiber back cane	12" - 15"	increment of 1/2"	Set at ²			\$330
BC09-1518	☐	Straight Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²			\$330
BC09-1821	☐	Straight Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²			\$330
BC10-1315	☐	8° bend Carbon fiber back cane	13.5" - 15"	increment of 1/2"	Set at ²			\$380
BC10-1518	☐	8° bend Carbon fiber back cane	15" - 18"	increment of 1/2"	Set at ²			\$380
BC10-1821	☐	8° bend Carbon fiber back cane	18" - 21"	increment of 1/2"	Set at ²			\$380
		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐
PUSH HANDLES OPTIONS (ONE option at a time)								
PH01	☐	Push handles						STD
PH05	☐	Fold down push handles ^{3,5}						\$165
PH06-920	☐	Clamp on height adjustable push handles 9"-20" ^{3,4,5}						\$345
PH06-923	☐	Clamp on height adjustable push handles 9"-23" ^{3,4,5}						\$345
PH02	☐	Omit push handles						NC
STROLLER HANDLE AND STABILIZER BAR (only available with push handle)								
PH04	☐	Stroller handle (not available with the HD option)						\$220
PH03	☐	Newton folding stabilizer bar (not available with stroller handle)						\$275

1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265 lb weight limit;

4- Not compatible with swing away armrests, 8° bend back canes, 9-12" back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

ARMRESTS

HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ^{1,4} *VP01								
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐			\$405
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐			\$405
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐			\$405
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐			\$405
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01								
AR13-610	☐	10" Desk length armpad ⁵	6" - 10"					\$385
AR03-812	☐	10" Desk length armpad	8" - 12"					\$385
AR03-1014	☐	10" Desk length armpad	10" - 14"					\$385
AR14-610	☐	14" Full length armpad ⁵	6" - 10"					\$385
AR04-812	☐	14" Full length armpad	8" - 12"					\$385
AR04-1014	☐	14" Full length armpad	10" - 14"					\$385
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²								
AR05	☐	Tubular armrest - short pad	8" - 11"					\$325
AR23	☐	Tubular armrest - long pad	8" - 11"					\$325
AR06	☐	10" Desk length flat armpad	8" - 11"					\$325
AR17	☐	14" Full length flat armpad	8" - 11"					\$325
ARMRESTS OPTIONS								
AR12	☐	14" Full Length armpad - Gel Ovation (pair) ³						\$160
AR99	☐	Omit armrest						NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

ACCESSORIES						
SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>						
SG01	☐	Plastic side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$235
SG02	☐	Carbon fiber side guards <i>(fixed)</i>				\$285
SG03	☐	Carbon fiber side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$320
SG06	☐	Plastic fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$340
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$535
SG99	☐	Omit side guard				NC

POSITIONING BELT

SB01-S	☐	Velcro adjustable - Small 35" ¹			\$50
SB01-L	☐	Velcro adjustable - Long 60"			\$50
SB02-S	☐	Auto buckle - Small 35" ¹			\$80
SB02-M	☐	Auto buckle - Medium 48"			\$80
SB02-L	☐	Auto buckle - Long 60"			\$80
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"			\$105
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"			\$105
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"			\$195
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"			\$195
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"			\$195
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"			\$220
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"			\$220
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"			\$220
SB03	☐	Hardware attachment for 2 extra points (1 pair)			\$15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"			\$365
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"			\$365
SB99	☐	Omit positioning belt			STD

HEADREST¹

AC01-S	☐	Profiled Headrest - Small (3.5"X6")			\$215
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")			\$215
AC01-L	☐	Profiled Headrest - Large (5.5"X8")			\$215
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar ¹			\$365

ACCESSORIES

WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6" camber			\$35
AC03	☐	Anti-tippers			\$125
AC29	☐	Swing-away anti-tipper ¹ Left ☐ Right ☐			\$190
AC27	☐	Tip assist ²			\$70
AC06	☐	Cane and crutch holder Left ☐ Right ☐			\$115
AC07	☐	Oxygen tank holder ⁴ - D and E Type - installed on the left			\$395
AC08	☐	Spoke guards ³ Available: 20", 22", 24", 25", 26"			\$165
AC09	☐	Footrest neoprene impact guard (<i>pair</i>) - Black			\$90
AC10	☐	Seat pouch			\$105
AC30	☐	Newton multi-purpose pouch			\$80
AC26	☐	BodyPoint Mobility Bag			\$150
AC23	☐	Newton Backpack			\$135
AC12	☐	Velcro Closing Strap			\$60
AC13	☐	Tool kit			\$60


3- Not available with Mag wheels and Spinergy; the 24" spoke guard will be installed on the 25" wheel; **4-** Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down

Colors

Options ¹		Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
																						
AirGrip	☐ \$160		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	\$350
Frame ²		STD	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	NC
FC99	☐	Custom Color	Specify PPG Paintit color: additional 10 days lead-time																		\$422	

2- With carbon fiber crossbrace:

Fork and Stem bolt combinations



Frame length chart

Frame length	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
STD							
+1							N/A
+2						N/A	N/A
-1	N/A						
-2	N/A	N/A					

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75



CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the
EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The **EVO Custom Request Form** is the best place to start.

Notes :



Green boxes are stocked sizes; the red boxes require customization*; the N/A boxes are not available sizes.
*Customization (red boxes) is available with an upcharge of \$160.00 and has a 4 to 6 weeks additional lead time.

	MODEL #	Description	Fits Back Width	Retail	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	\$80	

THORACIC SHAPE	Optima™						Optima™ DEEP					
	 THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)						 THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)					
	WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM	WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM
	14"	13"	NOCFACB1413	\$1,150	NOACB1413	\$725	14"	13"	NOCFPCB1413	\$1,250	NOPCB1413	\$825
	14"	16"	NOCFACB1416	\$1,150	NOACB1416	\$725	14"	16"	NOCFPCB1416	\$1,250	NOPCB1416	\$825
	15"	13"	NOCFACB1513	\$1,150	NOACB1513	\$725	15"	13"	NOCFPCB1513	\$1,250	NOPCB1513	\$825
	15"	16"	NOCFACB1516	\$1,150	NOACB1516	\$725	15"	16"	NOCFPCB1516	\$1,250	NOPCB1516	\$825
	16"	13"	NOCFACB1613	\$1,150	NOACB1613	\$725	16"	13"	NOCFPCB1613	\$1,250	NOPCB1613	\$825
	16"	16"	NOCFACB1616	\$1,150	NOACB1616	\$725	16"	16"	NOCFPCB1616	\$1,250	NOPCB1616	\$825
	16"	19"			NOACB1619	\$725	16"	19"			NOPCB1619	\$825
	17"	13"	NOCFACB1713	\$1,150	NOACB1713	\$725	17"	13"	NOCFPCB1713	\$1,250	NOPCB1713	\$825
	17"	16"	NOCFACB1716	\$1,150	NOACB1716	\$725	17"	16"	NOCFPCB1716	\$1,250	NOPCB1716	\$825
	18"	13"	NOCFACB1813	\$1,150	NOACB1813	\$725	18"	13"	NOCFPCB1813	\$1,250	NOPCB1813	\$825
	18"	16"	NOCFACB1816	\$1,150	NOACB1816	\$725	18"	16"	NOCFPCB1816	\$1,250	NOPCB1816	\$825
	18"	19"			NOACB1819	\$725	18"	19"			NOPCB1819	\$825
	19"	16"			NOACB1916	\$725	19"	16"			NOPCB1916	\$825
	19"	19"			NOACB1919	\$725	19"	19"			NOPCB1919	\$825
	20"	13"	NOCFACB2013	\$1,150	NOACB2013	\$725	20"	13"	NOCFPCB2013	\$1,250	NOPCB2013	\$825
	20"	16"	NOCFACB2016	\$1,150	NOACB2016	\$725	20"	16"	NOCFPCB2016	\$1,250	NOPCB2016	\$825
	20"	19"			NOACB2019	\$725	20"	19"			NOPCB2019	\$825
	21"	16"			NOACB2116	\$755	21"	16"			NOPCB2116	\$860
	21"	19"			NOACB2119	\$755	21"	19"			NOPCB2119	\$860
	22"	13"			NOACB2213	\$755	22"	13"			NOPCB2213	\$860
	22"	16"			NOACB2216	\$755	22"	16"			NOPCB2216	\$860
	22"	19"			NOACB2219	\$755	22"	19"			NOPCB2219	\$860
	ITEM		WIDTH X DEPTH		RETAIL	QTY	ITEM		WIDTH X DEPTH		RETAIL	QTY
	Replacement Cover & Foam		NOACC		\$130		Replacement Cover & Foam		NOPCC		\$130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

THORACIC SHAPE	xtend™	WIDTH	HEIGHT	RETAIL	✓
	NHAACB (L)1213/17	12"	13"-17"	\$725	
	NHAACB (L)1413/17	14"	13"-17"	\$725	
	NHAACB (L)1613/17	16"	13"-17"	\$725	
	NHAACB (L)1813/17	18"	13"-17"	\$755	
	NHAACB (L)2013/17	20"	13"-17"	\$755	
	NHAACB (L)2213/17	22"	13"-17"	\$755	
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		N3LHAACC		\$130	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	xtend™ Deep	WIDTH	HEIGHT	RETAIL	✓
	NHAPCB1616/20	16"	16"-20"	\$725	
	NHAPCB1816/20	18"	16"-20"	\$725	
	NHAPCB2016/20	20"	16"-20"	\$725	
	NHAPCB2216/20	22"	16"-20"	\$755	
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		N3HAPCC		\$130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

Classic™ PELVIC BACK SUPPORT

PELVIC SHAPE	Classic™	WIDTH	HEIGHT	RETAIL	✓
	NACB1616	16"	16"	\$725	
	NACB1618	16"	18"	\$725	
	NACB1816	18"	16"	\$725	
	NACB1818	18"	18"	\$725	
	NACB1820	18"	20"	\$725	
	NACB2020	20"	20"	\$725	
	NACB2216	22"	16"	\$755	
	NACB2218	22"	18"	\$755	
	NACB2220	22"	20"	\$755	
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		NACC		\$130	

Classic™ DEEP PELVIC BACK SUPPORT

	Classic™ Deep	WIDTH	HEIGHT	RETAIL	✓
	NPCB1616	16"	16"	\$830	
	NPCB1618	16"	18"	\$830	
	NPCB1818	18"	18"	\$830	
	NPCB1820	18"	20"	\$830	
	NPCB2020	18"	20"	\$830	
	NPCB2216	20"	16"	\$860	
	NPCB2218	22"	18"	\$860	
	NPCB2220	22"	20"	\$860	
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		NPCC		\$130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

BACK SUPPORT ACCESSORIES

	MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	\$205	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	\$70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	\$70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	\$70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	\$70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	\$70	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

MODEL #	Description	RETAIL	✓
NIB-0784	10" x 5" Large Rectangular	\$205	
NIB-0818	14" x 5" Extra-Large Rectangular	\$205	

MultiFit™ HEAD SUPPORT – FIXTURES

MODEL	Description	RETAIL	✓
NIB-0747	Multi-Axis Hardware	\$150	
NIB-0823	Head Support Receiver Bracket	\$60	

EXTRA COVERS

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	NXCCBF _____ v	\$130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	NXCCBF _____ v	\$150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	NXCCAF _____ v	\$130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	NXCCAF _____ v	\$150	

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com