



2915 Ogletown RD # 2270
Newark, DE 19713
T. +1-866-650-6555 F. +1-888-966-6555
orders@motioncomposites.com
www.motioncomposites.com



**AWARD
HARDING
INNOVATIVE
PRODUCT**
2010
2013
2014
2016



5.1

Helio C2

ORDER FORM - USA

P.O. number: * _____ Date: _____

Company Name: * _____ Phone: * _____

Ordered by: * _____

Clinician/ATP/RTS: _____

Email: _____

Marked for: _____

Ship to: * _____

Name: _____

Address: _____

City: _____ State: _____

Zip Code: _____

Phone: _____

* Required MA082 English Manual (STD) MA089 Spanish Manual Taxes not included

MODEL										USD		
1	WC05	Ultralight Folding Carbon Fiber Helio C2 Wheelchair (265 lb weight capacity)								K0005	\$3 790	☐
2	WC09	HD Heavy Duty package (350 lb weight capacity)									\$500	☐
3	WC10	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6" camber									\$30	☐
4	SW	Seat width:	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐			NC	☐
5	SW	Seat width:	20" ¹ ☐	21" ² ☐	22" ² ☐					E2201	\$495	☐
SEAT DEPTH												
6	SD	Seat depth:	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐			NC	☐
7	SD	Seat depth:			20" ☐					E2203	\$475	☐

1 - Available with HD option

2 - Only available with HD option

[3. Width calculator](#)

NEWTON ULTRALIGHT PACKAGE												
8	PK01	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels								Save \$235	\$480	☐

C2 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13	14	15	13 1/2	14 1/2	15 1/2												
	P2	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2									
5"	P1																		
	P2				14	15	16	14 1/2	15 1/2	16 1/2									
	P3	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2						
	P4	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
7"	P1										16 1/2	17 1/2	18 1/2	17	18	19			
	P2							16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2			
	P3							17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2
	P4	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21
	P5	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21	19 1/2	20 1/2	21 1/2

Height available +- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
9	CA01	Composite wheel w/ PU tire (1" width)	3" ☐	4" ☐	5" ☐	6" ☐	7" ☐	8" ☐			NC	☐
10	CA10	Pneumatic casters 6" x 1 ¼"								E2214/E2215	\$160	☐
11	CA11	Pneumatic casters 8" x 1 ¼"								E2214/E2215	\$160	☐
12	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "								\$80	☐
13	CA15		5" x 1 ½ "								\$80	☐
14	CA16		6" x 1 ½ "								\$80	☐
15	CA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½ "								\$175	☐
16	CA18		5" x 1 ½ "								\$185	☐
17	CA19		6" x 1 ½ "								\$195	☐
CASTER OPTIONS												
18	CA12	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY)								K0073	\$295	☐
19	CA13	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)								E1015	\$475	☐

C2 - REAR Seat-to-floor Height

WHEEL SIZES		Height available			
20" (451) wheel ☐		12 1/2	to	17 1/4	Specify desired height:
22" (501) wheel ☐		13 1/2	to	18 1/4	
24" (540) wheel ☐		14 1/2	to	19 1/4	
25" (559) wheel ☐		15	to	19 3/4	
26" (590) wheel ☐		15 1/2	to	20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest											
SWING-IN/SWING-OUT FOOTREST											
20	FR03	60° Swing-in/swing-out footrest	13.5 - 17" □	14.5 - 19.5" □	15.5 - 20.5" □		Set at: *			NC	□
21	FR04	70° Swing-in/swing-out footrest	13.5 - 17" □	14.5 - 19.5" □	15.5 - 20.5" □		Set at: *			NC	□
22	FR15	80° Swing-in/swing-out footrest	13.5 - 16" □	14.5 - 19.5" □	15.5 - 20.5" □		Set at: *			NC	□
90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT <i>(Footplate mounted to hanger tube)</i>											
23	FR06	90° swing-in/swing-out with front mount - Standard <i>(no extension)</i> ¹	7.5 - 11"				Set at: *			NC	□
24	FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5 - 13"				Set at: *			NC	□
25	FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5 - 14.5"				Set at: *			NC	□
FOOTREST OPTIONS ³											
26	FR09	High mount <i>(footplate attached to the hanger with clamps)</i>	5 - 13.5" □				Set at: *			\$75	□
27	FR99	Omit Footrest								NC	□

1 - Available only with two pieces adjustable footplates 2 - May interfere with front casters 3- Not available with front mount * + or - 1/2"

ELEVATING LEGREST (unit)				The elevating legrest replace the footrest side selected							
28	FR10	Elevating legrest with calf pad, 14" - 22" Length adjustable	Left ☐	Right ☐	Set at:		E0990	\$240		☐	

FOOTPLATES											
29	FP01	Standard - Composite								STD	☐
30	FP02	Newton adjustable angle - Composite							K0040	\$210	☐
31	FP03	Newton adjustable angle - Aluminum							K0040	\$230	☐
32	FP04	Newton adjustable angle - Carbon fiber (with rubber inserts)							K0040	\$415	☐
33	FP06	One piece adjustable angle flip-up - Aluminum							K0040	\$190	☐
34	FP07	One piece adjustable angle flip-up - Carbon fiber (with rubber inserts)							K0040	\$435	☐

FOOTPLATE ACCESSORIES											
35	FP08	Length adjustable heel loop (not available with one piece footplate)							E0951	\$40	☐
36	FP09	Calf strap							K0038	\$35	☐
37	FP30	Bodypoint Aeromesh Padded calf strap							K0038	\$115	☐

RESIDUAL LIMB SUPPORT FLAT CUSHION (Swing away, depth, width, height and angle ajustable) - (unit)											
38	FR27	Residual limb support swing away - flat 6"x8"	Left ☐	Right ☐	New	E1020	\$405	☐			
39	FR28	Residual limb support swing away - flat 8"x10"	Left ☐	Right ☐	New	E1020	\$405	☐			
40	FR29	Residual limb support swing away - flat 10"x10"	Left ☐	Right ☐	New	E1020	\$405	☐			
41	FR30	Residual limb support swing away - flat 10"x14"	Left ☐	Right ☐	New	E1020	\$405	☐			
RESIDUAL LIMB SUPPORT CONTOUR CUSHION (Swing away, depth, width, height and angle ajustable) - (unit)											
42	FR32	Residual limb support swing away - contour 6"x8"	Left ☐	Right ☐	New	E1020	\$405	☐			
43	FR33	Residual limb support swing away - contour 8"x10"	Left ☐	Right ☐	New	E1020	\$405	☐			

Back canes, height adjustable, with option										
FIXED ANGLE BACK										
44	BC06	Straight Carbon fiber back cane	9 - 12" ☐	12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐			NC	☐
45	BC08	8° Bend Carbon fiber back cane		13 ½ - 15" ☐	15 - 18" ☐	18 - 21" STD ☐			\$50	☐
ADJUSTABLE ANGLE BACK										
46	BC09	Adjustable angle, straight Carbon fiber back cane		12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐			\$200	☐
47	BC10	Adjustable angle, 8° bend Carbon fiber back cane		13 ½ - 15" ☐	15 - 18" ☐	18 - 21" STD ☐			\$250	☐
48		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐		
Set height to: ²		(measured from seat sling to top of back upholstery)								
		PUSH HANDLE OPTION								
49	PH01	Push handles							STD	☐
50	PH05	Folding push handles ^{3, 5}							\$120	☐
51	PH06	clamp on height adjustable push handle ^{3, 4, 5}		9 - 20" STD ☐	9-23" ☐				\$300	☐
52	PH02	Without push handles							NC	☐

1- Not available with a 8 degree bend cane 2 - If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit

4 - Not compatible with swing away armrests, 8° bend back canes a nd ampute axle plates

5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

STROLLER HANDLE AND STABILIZER BAR ¹ (only available with push handle)									
53	PH04	Stroller handle (not available with the HD option)					K0108	\$160	☐
54	PH03	Newton folding stabilizer bar						\$250	☐

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY									
55	BH07	Slip-on Active-Mesh back upholstery						STD	☐
56	BH08	Adjustable tension Active-Mesh back upholstery					E2611	\$300	☐
57	BH99	Omit Back upholstery						NC	☐
58		NXT back series				see pricelist p.7			

SEAT CUSHION									
59	SC01	2" Cushion						\$165	☐
60	SC02	3" Cushion						\$185	☐
61	SC99	Omit seat cushion						STD	☐
62		NXT FIT cushion series				see pricelist p.7			

Rear wheels											
ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims											
63	HR22	One arm drive ^{2,3}			24"	left ☐	right ☐			\$1 395	☐
64	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)								NC	☐
REAR WHEELS											
65	RW01	Mag ⁴	20" ☐	22" ☐	24" ☐					\$50	☐
66	RW02	Newton One - spoke wheel ⁵	20" ☐	22" ☐	24" ☐					NC	☐
67	RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐			\$300	☐
68	RW04	Spinergy Spox			22" ☐	24" ☐	25" ☐	26" ☐		\$795	☐
69		Spinergy Spox spokes colors ⁶	Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐			NC	☐
70	RW05	Spinergy LX			22" ☐	24" ☐	25" ☐	26" ☐		\$995	☐
71		Spinergy LX spokes colors ⁶	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐
72	RW07	Spinergy CLX Carbon Blade				24" ☐	25" ☐			\$1 480	☐
73		Spinergy CLX spokes colors ⁶	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐
TIRES							Available				
74	TI01	Soft Urethane 1 1/4 "					20", 22" 24"		E2218	STD	☐
75	TI02	Soft Urethane knobby 1 3/8 "					22", 24", 25", 26"		E2218	\$95	☐
76	TI04	Soft Urethane - Shox 1"					22", 24", 25", 26"		E2218	\$115	☐
77	TI05	Pneumatic 1 3/8 "					20", 22", 24", 25", 26"		E2211/E2212	NC	☐
78	TI08	Pneumatic with Airless insert 1 3/8"					22", 24"		E2211/E2213	\$80	☐
79	TI06	Pneumatic HP - 1"					20", 22", 24", 25", 26"		E2211/E2212	\$95	☐
80	TI07	Pneumatic HP - Schwalbe Marathon Plus 1"					24", 25", 26"		E2211/E2212	\$125	☐

1- Offered with handrims on each side of the wheelchair

2- Only standard handrim position available

3 - Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position

5 - Not compatible with Shox tires

6 - Colour spokes require an additionnal 10 days lead-time

Motion Composites

T.450 588-6555 1 866 650-6555 F.1 888 966-6555

HELIO C2 Orderform USA 5.1.xlsx

March 23rd, 2020 3 of 6

Handrims						
HANDRIM						
81	HR01	Aluminum anodized			STD	☐
82	HR02	Plastic coated ³			\$95	☐
83	HR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}			\$350	☐
84	HR07	Newton Air Grip ¹ Black ☐			\$95	☐
85	HR08	Newton Air Grip ¹ (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐			\$240	☐
		(11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐				
		(15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐				
86	HR28	Surge ^{1, 2, 3} (1-1/2" Oval handrim w/ Gripton high friction strip)			\$375	☐
87	HR05	Surge LT ^{1, 2, 3} (1-1/8" Oval handrim w/ Gripton high friction strip)			\$405	☐
88	HR04	Natural Fit ^{1, 2, 3}			\$310	☐
89	HR99	Omit handrim			NC	☐

1 - Not available with Mag wheels or OAD
2 - Not available in 20"
3 - Only standard handrim position available

Handrims option						
Assembly position for aluminum and Newton AirGrip handrims						
90	HR29	Standard Handrim position Full tab			STD	☐
91	HR30	Narrow Handrim position Short tab 2 positions			\$25	☐
92	HR31	Super Narrow Handrim Short tab 2 positions STD ☐ Short tab 1 position ☐			\$25	☐
Thumb piece ⁴ selection for Natural Fit handrims						
93	HR23	Standard grip			\$85	☐
94	HR27	Super grip			\$140	☐
95	HR099	Omit thumb piece			NC	☐


4- Compatible with Surge, Surge LT handrims

AXLE						
96	AX01	Quick release axle			STD	☐
97	AX02	Permanent axle			NC	☐
98	AX03	Quad release axle		K0108	\$125	☐
AXLE PLATE						
99	AX04	Amputee adjustable axle plate ¹		E0959	\$195	☐

1 - May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER						
100	CM00	0° camber			NC	☐
101	CM03	3° camber			NC	☐
102	CM06	6° camber <i>(not available with amputee axle plate)</i>			\$150	☐

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable								Most responsive				
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK						
103	WL08	Newton Wheel lock - push to lock			STD	☐
104	WL02	Pull to lock			NC	☐
105	WL04	Newton Grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>		E0974	\$150	☐
106	WL13 WL14	6" removable extension handle		E0961	\$75	☐
107	WL05	Scissor Lock			\$140	☐
108	WL06	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>			\$200	☐
109	WL99	Omit wheel locks			NC	☐

Armrest									
FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ^{1,4}				Height Adjustment			Removable function		
110	AR01	10" Desk length arm pad		8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐ <i>STD</i>	E0973	\$350 ☐
111	AR02	14" Full length arm pad		8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked ☐ <i>STD</i>	E0973	\$350 ☐
"T" ARMREST ¹				Height Adjustment					
112	AR03	10" Desk length arm pad		8 - 12" STD ☐	10 - 14" ☐			E0973	\$325 ☐
113	AR04	14" Full length arm pad		8 - 12" STD ☐	10 - 14" ☐			E0973	\$325 ☐
SWING AWAY ARMREST ²				Height Adjustment					
114	AR05	Tubular arm pad		8 - 11"				E0973	\$225 ☐
115	AR06	10" Desk length arm pad		8 - 11"				E0973	\$225 ☐
GEL ARM PAD ³									
116	AR12	14" Full Length Arm Pad - <i>Gel Ovation</i>						SEAND0030	\$150 ☐
117	AR99	Omit armrest							NC ☐
1- May have a 1/2 " adjustment restriction on both limits of the range 2- Could interfere with some rigid back rests 3- Not available with swing away armrest; 4- Not available with attendant lock									
SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>									
118	SG01	Plastic side guards <i>(removable)</i>						\$195	☐
119	SG02	Carbon fiber side guards <i>(fixed)</i> ¹						\$250	☐
120	SG03	Carbon fiber side guards <i>(removable)</i>						\$290	☐
121	SG06	Plastic fender side guards <i>(removable)</i> ²						\$270	☐
122	SG04	Carbon fiber fender side guards <i>(removable)</i>						\$540	☐
123	SG99	Omit side guard						NC	☐
				SIDE GUARDS OPTION					
124	SG07	Side guards clip <i>(pair)</i> ³					<i>New</i>		NC ☐
1- Not available with back upholstery 2- Not compatible with 20" and 22" rear wheels 3- Not available with fixed carbon fiber side guard									
POSITIONING BELT									
125	SB01	Velcro adjustable		Standard 35" ☐	Long 60" ☐			E0978	\$45 ☐
126	SB02	Auto buckle		Standard 35" ☐	Long 60 " ☐			E0978	\$85 ☐
127	SB07	Hardware attachment for 2 point seatbelt (1 pair)						NC	☐
128	SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)						\$15	☐
129	SB04	BodyPoint Hip Belt - 2 point not padded		Medium 56" ☐	Long 61" ☐			E0978	\$80 ☐
130	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐			E0978	\$170 ☐
131	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐			E0978	\$190 ☐
132	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps		Medium 18" ☐	Long 24" ☐			E0978	\$285 ☐
133	SB99	Omit positioning belt						STD	☐
HEADREST ¹									
134	AC01	Profiled Headrest	Small (3.5"x6") ☐	Medium (4.5"x7") ☐	Large (5.5"x8") ☐			E0955	\$180 ☐
135	AC28	Headrest mounting hardware with removable non folding stabilizer bar						E1028	\$295 ☐
1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar									

ACCESSORIES					
136	AC03	Anti tippers		E0971	\$140 ☐
137	AC27	Tip assist ¹			\$60 ☐
138	AC04	Transit Tie Down (<i>Unoccupied</i>)			NC ☐
139	AC05	WC-19 Transit Tie Down w/Q'straint belt ²	Short belt 20-41" ☐ Long belt 20-46" ☐		\$320 ☐
140	AC06	Cane carrier	Left ☐ Right ☐	E2207	\$195 ☐
141	AC07	Oxygen tank holder - E Type <i>Not compatible with swing away armrest or amputee axle plate</i>	Left	E2208	\$310 ☐
142	AC08	Spoke guards ³	Available: 20", 22", 24", 26"	K0065	\$115 ☐
143	AC09	Footrest neoprene impact guard (pair) - Black		K0108	\$90 ☐
144	AC10	Seat pouch			\$75 ☐
145	AC26	BodyPoint Mobility Bag			\$125 ☐
146	AC23	Newton Backpack			\$120 ☐
147	AC12	Velcro Closing Strap			\$50 ☐
148	AC13	Tool kit			\$45 ☐
149	AC16	Motion Composites Polo Shirt	W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐		\$32 ☐
150	AC15	Military Patches	Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐		\$90 ☐
151	AC15	Combat Wounded Purple Heart Patch			\$90 ☐

1 - Cannot be combined with Anti Tippers 2 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option
3 - Not available with Mag wheels and Spinergy

FRAME COLOR					
152	FC01	Carbon Fiber Black ☐			NC ☐
153	FC03	Sunkissed Orange ☐ Burgundy ☐ Sapphire Blue ☐ Monster Green ☐ Charcoal ☐ Fuchsia ☐ Steel Blue ☐ Ferrari Red ☐ Acid Green ☐			NC ☐
154	FC99	Custom Color	Specify PPG Paintit color:	Additional 10 days lead-time	\$360 ☐

TOTAL

Fork and Stem bolt combinations

STD

+1"

+2"

STD

+1"

+2"

STD

+1"

+2"

4" 5" 7"

Notes :



BioFit™ CUSHION

BioFit™	BioFit™	Fits Width	Depth	\$	✓
 HCPCS CODE: E2607-E2608	NXCABF1616	15"-16"	16"	425	
	NXCABF1618	15"-16"	18"	425	
	NXCABF1620	15"-16"	20"	425	
	NXCABF1816	17"-18"	16"	425	
	NXCABF1818	17"-18"	18"	425	
	NXCABF1820	17"-18"	20"	425	

BioFit™	Fits Width	Depth	\$	✓
NXCABF2016	19" - 20"	16"	425	
NXCABF2018	19" - 20"	18"	425	
NXCABF2020	19" - 20"	20"	425	
NXCABF2216	21" - 22"	16"	455	
NXCABF2218	21" - 22"	18"	455	
NXCABF2220	21" - 22"	20"	455	

Optima™ THORACIC BACK SUPPORT

Optima™	Optima™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NOACB1613	16"-18"	13"	660	
	NOACB1616	16"-18"	16"	660	
	NOACB1619	16"-18"	19"	660	
	NOACB1813	18"-20"	13"	660	
	NOACB1816	18"-20"	16"	660	
	NOACB1819	18"-20"	19"	660	
	NOACB2016	20"-22"	16"	660	
	NOACB2019	20"-22"	19"	660	
	NOACB2213	22"-24"	13"	710	
	NOACB2216	22"-24"	16"	710	
	NOACB2219	22"-24"	19"	710	

Optima™ DEEP THORACIC BACK SUPPORT

Optima™ DEEP	Optima™ DEEP	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NOPCB1616	16"-18"	16"	660	
	NOPCB1619	16"-18"	19"	660	
	NOPCB1816	18"-20"	16"	660	
	NOPCB1819	18"-20"	19"	660	
	NOPCB2016	20"-22"	16"	660	
	NOPCB2019	20"-22"	19"	660	
	NOPCB2213	22"-24"	13"	710	
	NOPCB2216	22"-24"	16"	710	
	NOPCB2219	22"-24"	19"	710	

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™	Xtend™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NHAACB1616/20	16"-18"	16"-20"	660	
	NHAACB1816/20	18"-20"	16"-20"	660	
	NHAACB2016/20	20"-22"	16"-20"	660	
	NHAACB2216/20	22"-24"	16"-20"	710	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™ DEEP	Xtend™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NHAACB1616/20	16"-18"	16"-20"	660	
	NHAACB1816/20	18"-20"	16"-20"	660	
	NHAACB2016/20	20"-22"	16"-20"	660	
	NHAPCB2216/20	22"-24"	16"-20"	710	

Classic™ PELVIC BACK SUPPORT

Classic™	Classic™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NACB1616	16"-18"	16"	660	
	NACB1618	16"-18"	18"	660	
	NACB1816	18"-20"	16"	660	
	NACB1818	18"-20"	18"	660	
	NACB1820	18"-20"	20"	660	
	NACB2020	20"-22"	20"	660	
	NACB2216	22"-24"	16"	710	
	NACB2218	22"-24"	18"	710	
	NACB2220	22"-24"	20"	710	

Classic™ DEEP PELVIC BACK SUPPORT

Classic™ DEEP	Classic DEEP™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NPCB1616	16"-18"	16"	660	
	NPCB1618	16"-18"	18"	660	
	NPCB1818	18"-20"	18"	660	
	NPCB1820	18"-20"	20"	660	
	NPCB2020	20"-22"	20"	660	
	NPCB2216	22"-24"	16"	710	
	NPCB2218	22"-24"	18"	710	
	NPCB2220	22"-24"	20"	710	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

 HCPCS CODE: E0955		
Model #NIB-0784	\$	✓
10" x 5" Large Rectangular	195	

 HCPCS CODE: E0955		
Model # NIB-0818	\$	✓
14" x 5" Extra-Large Rectangular	195	

MultiFit™ HEAD SUPPORT – FIXTURES

		
Model #NIB-0747	\$	✓
Multi-Axis Hardware	145	

		
Model # NIB-0823	\$	✓
Head Support Receiver Bracket	55	

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com