



HELIO kids

ENHANCED AUTONOMY -

GROWS WITH YOUR KIDS



HELIO kids C



HELIO kids A

A WHEELCHAIR SO NICE, WE MADE IT TWICE!

Introducing the HELIO Kids folding pediatric wheelchair, available in carbon fiber and aluminum, shaped by years of experience working with children and families.

Whether it's playing with friends, going to school, having fun and getting around easily... kids should be free to be kids!

Easier propulsion and handling → Reduced strain and improved everyday mobility thanks to our renowned performance-driven frame design, featuring a one-piece side frame, ultra-rigid folding system, and symmetrical cross brace.

Adapts as the child grows → Long-term usability, enhanced by a depth-adjustable back, growth kit and other features for even greater fit and flexibility over time.

Ultralight and durable → Lightweight for greater freedom of movement. Engineered to withstand kids being kids.

Safety-oriented by design → A choice of accessories, options and configurations, including WC19 compliant transit anchors, support secure transportation and peace of mind.

Easy personalization through colors and accessories → Kids can express themselves and feel a sense of ownership over their wheelchair through independent color selections for the frame and foot hangers, optional handrim colors, or additional customization solutions.



MEET THE STARS



HELIO kids C



HELIO kids A

MATERIAL	High-Modulus Carbon <ul style="list-style-type: none"> Made of carbon fiber, a material renowned for its ultralight properties and excellent vibration damping. 	Ultralight Aluminum 7005 T6 <ul style="list-style-type: none"> Made of ultralight aluminum 7005 T6, an alloy known for its high strength-to-weight ratio.
FEATURES	<ul style="list-style-type: none"> Symmetrical molded cross brace and sturdy unibody frame maximize propulsion efficiency and energy conservation, helping kids go further with less effort while reducing weight and maintenance. Includes access to a growth kit for the expected lifetime of the HELIO Kids wheelchair to help ensure the chair fits as the child grows. Personalized solutions through EVO Custom, because the right customization can be life changing for wheelchair users. Augmented mobility, thanks to a wide variety of NEWTON parts and accessories. Compatible with highly functional NXT seat cushions and back supports, built with premium materials for enhanced comfort, support, and durability. 	
TRANSPORT WEIGHT ^{1,3}	11.3 lbs (5.12 kg)	13.1 lbs (5.93 kg)
LIGHTEST FUNCTIONAL CONFIGURATION ^{2,3}	21.1 lbs (9.56 kg)	22.9 lbs (10.37 kg)
WEIGHT CAPACITY	165 lbs (74.8 kg)	165 lbs (74.8 kg)
SEAT WIDTH	12" (30.5 cm) to 16" (40.6 cm)	12" (30.5 cm) to 16" (40.6 cm)
SEAT DEPTH	12" (30.5 cm) to 19" (48.3 cm)	12" (30.5 cm) to 19" (48.3 cm)
FRONT SEAT HEIGHT	13" (33 cm) to 21.5" (54.6 cm)	13" (33 cm) to 21.5" (54.6 cm)
REAR SEAT HEIGHT	12.5" (30.5 cm) to 19.25" (48.9 cm)	12.5" (30.5 cm) to 19.25" (48.9 cm)
COG	COG range: 1" - 4" COG increments: 0.25"	COG range: 1" - 3.75" COG increments: 0.25"

¹ Without easily removable parts such as : Rear wheels, Armrests / Sideguards, Seating products (backrest and cushion), Anti-tippers, Depth adjustable back and Footrests.

² With: Rear wheels, Wheel locks, Footrests, Back upholstery and Seat cushion.

³ Based on a 12"x12" wheelchair, without depth adjustable back. Please note that weights may slightly vary.

A KID FRIENDLY, CAREGIVER TRUSTED PLATFORM.

Motion Composites continues to lead the industry in delivering thoughtful, safety-driven solutions that now support children through many stages of growth. We believe that high-performance wheelchairs enhance people's lives with components, design, features and quality. The HELIO Kids folding pediatric wheelchairs represent our passion for performance, our commitment to quality, and our vision of greater freedom brought to life.

LEARN MORE ABOUT
THE HELIO KIDS [HERE](#).

