



Helio KIDS

ORDER FORM - CANADA

P.O. number:* _____ Date: _____

Company Name:* _____ Phone:* _____

Ordered by:* _____

Clinician/Sales rep: _____

Email: _____

Marked for: _____

Ship to:* _____

Name: _____

Address: _____

City: _____ Prov: _____

Postal Code: _____

Phone: _____

* Required


MODEL		English Manual (STD)	Manuel Français	Taxes not included
1	WC06	Ultralight Folding Carbon Fiber Helio KIDS Wheelchair Including angle & depth adjustable back (175 lb weight capacity)		WM0000563/WAMZ \$3 545 ☐
2	WC10	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6" camber		\$30 ☐

FREE Growth Kit within the first 5 years of the Helio KIDS wheelchair ownership
The growth kit includes carbon crossbraces, crossbrace links, seat sling, back upholstery, front riggings and footplates.

SEAT WIDTH¹

3	SW12	Seat width:	12" ☐	13" ☐	14" ☐	15" ☐	16" ☐	NC	☐
---	------	-------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	----	--------------------------

SEAT DEPTH (Seat sling and frame have a fixed depth for improved rigidity and reduced weight)

4	SD14	12" - 15" (14" seat sling/frame)	Set at:	Seat depth adjustment achieved by changing back cane position 	NC	☐
5	SD15	13" - 16" (15" seat sling/frame)	Set at:		NC	☐
6	SD16	14" - 17" (16" seat sling/frame)	Set at:		NC	☐
7	SD17	15" - 18" (17" seat sling/frame)	Set at:		NC	☐
8	SD18	16" - 19" (18" seat sling/frame)	Set at:		NC	☐

[1. Width Calculator](#)

Helio Kids FRONT Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13	14	15	13 1/2	14 1/2	15 1/2												
	P2	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2									
5"	P1																		
	P2				14	15	16	14 1/2	15 1/2	16 1/2									
	P3	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2						
	P4	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
7"	P1										16 1/2	17 1/2	18 1/2	17	18	19			
	P2							16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2			
	P3							17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2
	P4	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21
	P5	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21	19 1/2	20 1/2	21 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
9	CA01	Composite wheel w/ PU tire (1" width)	3" ☐	4" ☐	5" ☐	6" ☐	7" ☐	8" ☐			NC	☐
10	CA10	Pneumatic 6" x 1 ¼"									\$160	☐
11	CA11	Pneumatic 8" x 1 ¼"									\$160	☐
12	CA14	Newton UltraCaster Soft-Roll Composite	4" x 1 ½ "								\$98	☐
13	CA15		5" x 1 ½ "								\$98	☐
14	CA16		6" x 1 ½ "								\$98	☐
15	CA17	Newton UltraCaster Soft-Roll Aluminum	4" x 1 ½ "								\$175	☐
16	CA18		5" x 1 ½ "								\$185	☐
17	CA19		6" x 1 ½ "								\$195	☐
CASTER OPTIONS												
18	CA12	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY)								WAMS	\$255	☐
19	CA13	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)									\$489	☐

KIDS REAR Seat-to-floor Height

WHEEL SIZES		Height available			
20" wheel ☐		12 1/2	to	17 1/4	Specify desired height:
22" wheel ☐		13 1/2	to	18 1/4	
24" wheel ☐		14 1/2	to	19 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ¼" of selected seat to floor height"

Footrest											
SWING-IN/SWING-OUT FOOTREST											
20	FR03	60° Swing-in/swing-out footrest	13.5 - 17" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at:*			NC	□	
21	FR04	70° Swing-in/swing-out footrest	13.5 - 17" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at:*			NC	□	
22	FR15	80° Swing-in/swing-out footrest	13.5 - 16" □	14.5 - 19.5" □	15.5 - 20.5" □	Set at:*			NC	□	
90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT <i>(Footplate mounted to hanger tube)</i>											
23	FR06	90° swing-in/swing-out with front mount - Standard <i>(no extension)</i> ¹	7.5 - 11"			Set at:*			NC	□	
24	FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5 - 13"			Set at:*			\$130	□	
25	FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5 - 14.5"			Set at:*			\$130	□	
FOOTREST OPTIONS ³											
26	FR09	High mount <i>(footplate attached to the hanger with clamps)</i>	5 - 13.5" □			Set at:*			\$90	□	
27	FR99	Omit Footrest							NC	□	

1 - Available only with two pieces adjustable footplates 2 - May interfere with front casters 3- Not available with front mount * + or - 1/2"

ELEVATING LEGREST (unit)				The elevating legrest replace the footrest side selected							
28	FR10	Elevating legrest with calf pad, 14" - 22" Length adjustable			Left ☐	Right ☐	Set at:		WAMK	\$290	☐

FOOTPLATES						
29	FP01	Standard - Composite			STD	☐
30	FP02	Newton adjustable angle - Composite		WAMJ	\$145	☐
31	FP03	Newton adjustable angle - Aluminum ¹		WAMJ	\$165	☐
32	FP04	Newton adjustable angle - Carbon fiber <i>(with rubber inserts)</i>		WAMJ	\$380	☐
33	FP05	Newton adjustable angle extra deep - Aluminum		WAMJ	\$165	☐
34	FP06	One piece adjustable angle flip-up - <i>Aluminum (Not available with 90° Front Mount)</i>	WAMJ	SEFND0015	\$206	☐
35	FP07	One piece adjustable angle flip-up - Carbon fiber <i>(Not available with 90° Front Mount)</i>	WAMJ	SEFND0015	\$450	☐
FOOTPLATE ACCESSORIES						
36	FP08	Length adjustable heel loop <i>(not available with one piece footplate)</i>			NC	☐
37	FP09	Calf strap		SEFND0010	\$52	☐
38	FP30	Bodypoint Aeromesh Padded calf strap		SEFND0010	\$129	☐

1 - The Newton adjustable angle (Aluminum) is the only available option with a 12 and 13" wide frame

RESIDUAL LIMB SUPPORT FLAT CUSHION (Swing away, depth, width, height and angle ajustable) - (unit)											
39	FR27	Residual limb support swing away - flat 6"x8"	Left ☐	Right ☐	New	SEFND0005	\$475	☐			
40	FR28	Residual limb support swing away - flat 8"x10"	Left ☐	Right ☐	New	SEFND0005	\$475	☐			
41	FR29	Residual limb support swing away - flat 10"x10"	Left ☐	Right ☐	New	SEFND0005	\$475	☐			
42	FR30	Residual limb support swing away - flat 10"x14"	Left ☐	Right ☐	New	SEFND0005	\$475	☐			
RESIDUAL LIMB SUPPORT CONTOUR CUSHION (Swing away depth, width, height and angle ajustable) - (unit)											
43	FR32	Residual limb support swing away - contour 6"x8"	Left ☐	Right ☐	New	SEFND0005	\$475	☐			
44	FR33	Residual limb support swing away - contour 8"x10"	Left ☐	Right ☐	New	SEFND0005	\$475	☐			

Back canes, depth ajustable and angle ajustable, with option									
BACK CANES - <i>Depth ajustable and angle ajustable</i>									
45	BC15	Straight aluminum back cane	9 - 12" ☐	12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐			NC ☐
46	BC16	8° Bend Aluminum back cane		12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐			\$105 ☐
47		Select back cane angle:	90° ☐	95° STD ☐	100° ☐	105° ☐			
Set height to: ¹		<i>(measured from seat sling to top of back upholstery)</i>							
		PUSH HANDLE OPTION							
48	PH01	Push handles						STD	☐
49	PH05	Fold down push handles ²						\$150	☐
50	PH02	Without push handles						NC	☐

1 - If not indicated, default height setting is at the lowest position 2 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

STROLLER HANDLE AND STABILIZER BAR ¹ (only available with push handle)								
51	PH04	Stroller handle				WAMW	\$210	☐
52	PH03	Newton folding stabilizer bar					\$258	☐

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY						
53	BH07	Slip-on Active-Mesh back upholstery			NC	☐
54	BH08	Adjustable tension Active-Mesh back upholstery		WAMD/WAME	\$245	☐
55	BH99	Omit Back upholstery			STD	☐
RIGID BACK						
56		NXT back series	See pricelist p.7			☐

SEAT CUSHION							
57	SC01	2" Cushion				\$175	☐
58	SC02	3" Cushion				\$195	☐
59	SC99	Omit Cushion				STD	☐
60		NXT FIT cushion series			See pricelist p.7		☐

Rear wheels													
ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims													
61	HR22	One arm drive ^{2,3}	24"	left ☐	right ☐						\$1 395	☐	
62	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)									NC	☐	
REAR WHEELS													
63	RW01	Mag ⁴	20" ☐	22" ☐	24" ☐						\$129	☐	
64	RW02	Newton One - spoke wheel ⁵	20" ☐	22" ☐	24" ☐						NC	☐	
65	RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐						\$360	☐	
66	RW04	Spinergy Spox		22" ☐	24" ☐						\$800	☐	
67		Spinergy Spox spokes colors ⁶	Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐				NC	☐	
68	RW05	Spinergy LX		22" ☐	24" ☐						\$995	☐	
69		Spinergy LX spokes colors ⁶	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐	NC	☐
70	RW07	Spinergy CLX Carbon Blade			24" ☐						\$1 540	☐	
71		Spinergy CLX spokes colors ⁶	Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐	NC	☐
TIRES						Available							
72	TI01	Soft Urethane 1 1/4 "				20", 22" 24"					STD	☐	
73	TI02	Soft Urethane knobby 1 3/8 "				22", 24"					\$110	☐	
74	TI04	Soft Urethane - Shox 1"				22", 24"					\$130	☐	
75	TI05	Pneumatic 1 3/8 "				20", 22" 24"					NC	☐	
76	TI08	Pneumatic with Airless insert 1 3/8"				22", 24"					\$215	☐	
77	TI06	Pneumatic HP - 1"				20", 22" 24"					\$110	☐	
78	TI07	Pneumatic HP - Schwalbe Marathon Plus 1"				24"					\$190	☐	

1- Offered with handrims on each side of the wheelchair

2- Only standard handrim position available

3- Not available with HD option

4- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position

5 - Not compatible with Shox tires

6 - Colour spokes require an additionnal 10 days lead-time

Handrims					
HANDRIM					
79	HR01	Aluminum anodized			STD ☐
80	HR02	Plastic coated ³		WAMO	\$160 ☐
81	HR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}		WAMO	\$430 ☐
82	HR07	Newton Air Grip ¹ <i>Black</i> ☐		WAMO	\$150 ☐
83	HR08	Newton Air Grip ¹ (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐ (11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐ (15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐		WAMO	\$330 ☐
84	HR28	Surge ^{1, 2, 3} (1-1/2" Oval handrim w/ Gripton high friction strip)			\$470 ☐
85	HR05	Surge LT ^{1, 2, 3} (1-1/8" Oval handrim w/ Gripton high friction strip)			\$500 ☐
86	HR04	Natural Fit ^{1, 2, 3}			\$310 ☐
87	HR99	Omit handrim			NC ☐

1 - Not available with Mag wheels or OAD 2 - Not available in 20" 3- Only standard handrim position available

Handrims option					
Assembly position for aluminum and Newton AirGrip handrims					
88	HR29	Standard Handrim position <i>Full tab</i>			STD ☐
89	HR30	Narrow Handrim position <i>Short tab 2 positions</i>			\$25 ☐
90	HR31	Super Narrow Handrim <i>Short tab 2 positions STD</i> ☐ <i>Short tab 1 position</i> ☐			\$25 ☐
Thumb piece ⁴ selection Natural Fit handrims					
91	HR23	Standard grip			\$190 ☐
92	HR27	Super grip			\$220 ☐
93	HR099	Omit thumb piece			NC ☐


4- Compatible with Surge, Surge LT handrims

AXLE					
94	AX01	Quick release axle		WAMU	STD ☐
95	AX02	Permanent axle			NC ☐
96	AX03	Quad release axle		WAMU	\$170 ☐
AXLE PLATE					
97	AX04	Amputee adjustable axle plate ¹		WAMT	\$250 ☐

1 - May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER					
98	CM00	0° camber			STD ☐
99	CM03	3° camber			NC ☐
100	CM06	6° camber (not available with amputee axle plate)			\$80 ☐

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable												Most responsive			
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"			
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"


Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

WHEEL LOCK					
101	WL08	Newton Wheel lock - push to lock			STD ☐
102	WL02	Pull to lock			NC ☐
103	WL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)		WAMR	\$175 ☐
104	WL13 WL14	6" removable extension handle			\$67 ☐
105	WL05	Scissor Lock			\$175 ☐
106	WL06	Attendant Lock (additional set of wheel lock mounted to the back)			\$230 ☐
107	WL99	Omit wheel locks			NC ☐

Armrests										
"T" ARMREST ¹			Height Adjustment							
108	AR03	10" Desk length arm pad	8 - 12" STD ☐		10 - 14" ☐				\$360	☐
109	AR04	14" Full length arm pad	8 - 12" STD ☐		10 - 14" ☐				\$360	☐
LOCKING FLIP-UP ARMREST ²										
110	AR07	Locking Flip-Up Armrest	Set height to:						\$330	☐
GEL ARM PAD ³										
111	AR12	14" Full Length Arm Pad - <i>Gel Ovation</i>						SEAND0030	\$150	☐
112	AR99	Omit armrest							NC	☐
<i>1 - May have a 1/2 " adjustment restriction on both limits of the range</i>										
<i>2 - May interfere with rigid back hardware 2 - Only available with 18" - 21" straight aluminum back canes 2 - Not available with a back upholstery</i>										
<i>3- Not available with locking flip-up armrest</i>										
SIDE GUARDS <i>(not available with "T" Armrest)</i>										
113	SG01	Plastic side guards <i>(removable)</i>						WAMV	\$225	☐
114	SG03	Carbon fiber side guards <i>(removable)</i>						WAMV	\$295	☐
115	SG06	Plastic fender side guards <i>(removable)</i> ¹						WAMV	\$320	☐
116	SG04	Carbon fiber fender side guards <i>(removable)</i>						WAMV	\$495	☐
117	SG99	Omit side guard							NC	☐
		SIDE GUARDS OPTION								
118	SG07	Side guards clip <i>(pair)</i>					New		NC	☐
<i>1 - Not compatible with 20" and 22" rear wheels</i>										
POSITIONING BELT										
119	SB01	Velcro adjustable	Standard 35" ☐		Long 60" ☐				\$45	☐
120	SB02	Auto buckle	Standard 35" ☐		Long 60" ☐				\$55	☐
121	SB07	Hardware attachment for 2 point seatbelt <i>(1 pair)</i>							NC	☐
122	SB03	Additional hardware for 2 extra attachment points <i>(needed for the 4 point seat belt) (1 pair)</i>							\$13	☐
123	SB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" ☐		Long 61" ☐				\$82	☐
124	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐		SERN0030	\$170	☐	
125	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐		SERN0020	\$195	☐	
126	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" ☐		Long 24" ☐			\$324	☐	
127	SB99	Omit positioning belt							STD	☐
HEADREST ¹										
128	AC01	Profiled Headrest	Small (3.5"X6") ☐		Medium (4.5"X7") ☐	Large (5.5"X8") ☐			\$205	☐
129	AC28	Headrest mounting hardware with removable non folding stabilizer bar							\$338	☐
<i>1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar</i>										
ACCESSORIES										
130	AC03	Anti tippers							\$115	☐
131	AC27	Tip assist ¹							\$60	☐
132	AC04	Transit Tie Down <i>(Unoccupied)</i>							NC	☐
133	AC06	Cane carrier	Left ☐		Right ☐				\$105	☐
134	AC07	Oxygen tank holder - E Type <i>Not compatible with swing away armrest or amputee axle plate</i>	Left				WAMX	\$375	☐	
135	AC08	Spoke guards ²	Available: 20", 22", 24"					WAMP	\$145	☐
136	AC09	Footrest neoprene impact guard (pair) - Black					SEFND0035	\$80	☐	
137	AC10	Seat pouch							\$95	☐
138	AC26	BodyPoint Mobility Bag							\$135	☐
139	AC23	Newton Backpack							\$120	☐
140	AC12	Velcro Closing Strap							\$50	☐
141	AC13	Tool kit							\$50	☐
142	AC16	Motion Composites Polo Shirt	W/S ☐	W/M ☐	W/L ☐				\$32	☐
			M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	M/XXXL ☐	M/XXXXL ☐	
<i>1 - Cannot be combined with Anti Tippers 2 - Not available with Mag wheels and Spinergy</i>										
FRAME COLOR <i>(with carbon finish crossbrace)</i>										
143	FC01	Carbon Fiber Black ☐							NC	☐
144	FC03	Sunkissed Orange ☐	Burgundy ☐	Sapphire Blue ☐	Monster Green ☐				NC	☐
		Charcoal ☐	Fuchsia ☐	Steel Blue ☐	Ferrari Red ☐	Acid Green ☐				
145	FC99	Custom Color	Specify PPG Paintit color:					10 days additional lead-time	\$410	☐
TOTAL										

Fork and Stem bolt combinations



Notes :



BioFit™ CUSHION

	BioFit™	Fits Width	Depth	\$	✓
	NXCABF1616	15"-16"	16"	595	
	NXCABF1618	15"-16"	18"	595	
	NXCABF1620	15"-16"	20"	595	
	NXCABF1816	17"-18"	16"	595	
	NXCABF1818	17"-18"	18"	595	
	NXCABF1820	17"-18"	20"	595	

BioFit™	Fits Width	Depth	\$	✓
NXCABF2016	19" - 20"	16"	595	
NXCABF2018	19" - 20"	18"	595	
NXCABF2020	19" - 20"	20"	625	
NXCABF2216	21" - 22"	16"	625	
NXCABF2218	21" - 22"	18"	625	
NXCABF2220	21" - 22"	20"	625	

Optima™ THORACIC BACK SUPPORT

	Optima™	Fits Width	Height	\$	✓
	NOACB1613	16"-18"	13"	690	
	NOACB1616	16"-18"	16"	690	
	NOACB1619	16"-18"	19"	690	
	NOACB1813	18"-20"	13"	690	
	NOACB1816	18"-20"	16"	690	
	NOACB1819	18"-20"	19"	690	
	NOACB2016	20"-22"	16"	690	
	NOACB2019	20"-22"	19"	690	
	NOACB2213	22"-24"	13"	720	
	NOACB2216	22"-24"	16"	720	
	NOACB2219	22"-24"	19"	720	

Optima™ DEEP THORACIC BACK SUPPORT

	Optima™ DEEP	Fits Width	Height	\$	✓
	NOPCB1616	16"-18"	16"	790	
	NOPCB1619	16"-18"	19"	790	
	NOPCB1816	18"-20"	16"	790	
	NOPCB1819	18"-20"	19"	790	
	NOPCB2016	20"-22"	16"	790	
	NOPCB2019	20"-22"	19"	790	
	NOPCB2213	22"-24"	13"	820	
	NOPCB2216	22"-24"	16"	820	
	NOPCB2219	22"-24"	19"	820	

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	690	
	NHAACB1816/20	18"-20"	16"-20"	690	
	NHAACB2016/20	20"-22"	16"-20"	690	
	NHAACB2216/20	22"-24"	16"-20"	720	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™ DEEP	Fits Width	Height	\$	✓
	NHAACB1616/20	16"-18"	16"-20"	690	
	NHAACB1816/20	18"-20"	16"-20"	690	
	NHAACB2016/20	20"-22"	16"-20"	690	
	NHAPCB2216/20	22"-24"	16"-20"	720	

Classic™ PELVIC BACK SUPPORT

	Classic™	Fits Width	Height	\$	✓
	NACB1616	16"-18"	16"	690	
	NACB1618	16"-18"	18"	690	
	NACB1816	18"-20"	16"	690	
	NACB1818	18"-20"	18"	690	
	NACB1820	18"-20"	20"	690	
	NACB2020	20"-22"	20"	660	
	NACB2216	22"-24"	16"	720	
	NACB2218	22"-24"	18"	720	
	NACB2220	22"-24"	20"	720	

Classic™ DEEP PELVIC BACK SUPPORT

	Classic™ DEEP	Fits Width	Height	\$	✓
	NPCB1616	16"-18"	16"	690	
	NPCB1618	16"-18"	18"	690	
	NPCB1818	18"-20"	18"	690	
	NPCB1820	18"-20"	20"	690	
	NPCB2020	20"-22"	20"	690	
	NPCB2216	22"-24"	16"	720	
	NPCB2218	22"-24"	18"	720	
	NPCB2220	22"-24"	20"	720	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

		
Model #NIB-0784	\$	✓
10" x 5" Large Rectangular	195	

		
Model # NIB-0818	\$	✓
14" x 5" Extra-Large Rectangular	195	

MultiFit™ HEAD SUPPORT – FIXTURES

		
Model #NIB-0747	\$	✓
Multi-Axis Hardware	145	

		
Model # NIB-0823	\$	✓
Head Support Receiver Bracket	55	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com