



Helio KIDS

ORDER FORM - CANADA

1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

Marked for: _____

2- Dealer/Representative information

Dealer/Representative:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:* _____

Name: _____

Address: _____

City: _____ Province: _____

Postal Code: _____

Phone: _____

* Required

MA080 ☐
English Manual (STD)

MA081 ☐
Manuel Français

Taxes not included

MODEL							
WC06	☒	Ultralight Folding Helio KIDS Wheelchair Including fixed seat depth OR depth adjustable back (175 lb weight capacity)	WM0000563/WAMZ			\$2,945	
WC08	☒	Carbon fiber option	WAMZ			\$1,450	
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber				\$35	
FREE Growth Kit within the first 5 years of the Helio KIDS wheelchair ownership The growth kit includes carbon crossbraces, crossbrace links, seat sling, back upholstery, front riggings and footplates.							
SEAT WIDTH ¹							
SW	☐	12" ☐ 13" ☐ 14" ☐ 15" ☐ 16" ☐					NC


[1. Width Calculator](#)



↓ Caution: complete only one of the following two sections ↓

SEAT DEPTH (Available with WC-19)			
SD12	☐	12" Seat sling - frame of 14"	NC
SD13	☐	13" Seat sling - frame of 14"	NC
SD14	☐	14" Seat sling - frame of 14"	NC
SD15	☐	15" Seat sling - frame of 15"	NC
SD16	☐	16" Seat sling - frame of 16"	NC
FIXED BACK CANE (MUST select ONE push handle options)			
BC11-1215	☐	Aluminium straight back cane (90°)	12 - 15" NC
BC11-1518	☐	Aluminium straight back cane (90°)	15 - 18" NC
BC11-1821	☐	Aluminium straight back cane (90°) - STD	18 - 21" NC
BC12-1315	☐	Aluminum back cane 8° bend	13.5 - 15" NC
BC12-1518	☐	Aluminum back cane 8° bend	15 - 18" NC
BC12-1821	☐	Aluminum back cane 8° bend -STD	18 - 21" NC
Specify height: ⁵ _____ measured from seat sling to top of back upholstery			
ANGLE ADJUSTABLE BACK CANE (MUST select ONE push handle options)			
BC13-1518	☐	Aluminium straight back cane (90°)	15 - 18" NC
BC13-1821	☐	Aluminium straight back cane (90°) - STD	18 - 21" NC
BC14-1518	☐	Aluminum back cane 8° bend	15 - 18" NC
BC14-1821	☐	Aluminum back cane 8° bend -STD	18 - 21" NC
Select angle for the back cane:		85° ☐ 90° ☐ 95° STD ☐ 100° ☐ 105° ☐ 110° ☐	
Specify height: ⁵ _____			

5- If not indicated, default height setting is at the lowest position

SEAT DEPTH- Depth and angle adustable (Not available with WC-19) ³			
SD14	☐	12" to 15" (seat sling/frame 14")	NC
SD15	☐	13" to 16" (seat sling/frame 15")	NC
SD16	☐	14" to 17" (seat sling/frame 16")	NC
SD17	☐	15" to 18" (seat sling/frame 17")	NC
SD18	☐	16" to 19" (seat sling/frame 18")	NC
Specify depth: ⁴ _____			
3 - Maximum weight of 175 lb with a depth adjustable backrest			
4 - If not indicated, the depth will be equal to the seat sling			
			
DEPTH AND ANGLE ADJUSTABLE BACK CANE (MUST select ONE push handle options)			
BC15-912	☐	Aluminium straight back cane 9 - 12"	NC
BC15-1215	☐	Aluminium straight back cane 12 - 15"	NC
BC15-1518	☐	Aluminium straight back cane 15 - 18"	NC
BC15-1821	☐	Aluminium straight back cane 18 - 21" STD	NC
BC16-1215	☐	Aluminum back cane 8° bend 12 - 15"	NC
BC16-1518	☐	Aluminum back cane 8° bend 15 - 18"	NC
BC16-1821	☐	Aluminum back cane 8° bend STD 18 - 21"	NC
Select angle for the back cane:		90° ☐ 95° STD ☐ 100° ☐ 105° ☐	
Specify height: ⁵ _____			

TRANSIT TIE DOWN					
AC04	☐	Transit Tie Down (<i>Unoccupied</i>)			NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶ Short belt 20"-41"			\$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁶ Long belt 20"-46"			\$365
AC35	☐	WC-19 Transit Tie Down without belt ⁶			\$345

6- Not available with depth adjustable back

PACKAGES					
VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹			\$95

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo ***VP01** ; Recommended with black chair;
See picture of the product on page 6.

Helio Kids FRONT
 Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2				14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						
	P4	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
7"	P1										16 1/2"	17 1/2"	18 1/2"	17"	18"	19"			
	P2							16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"			
	P3							17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
	P4	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"
	P5	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"	19 1/2"	20 1/2"	21 1/2"
Froglegs	P1				14 1/2"	15 1/2"	16 1/2"												
	P2	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"									
	P3	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"						

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS				Availability		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 ¼")		6", 8"		\$170
CA14 to CA16	☐	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		\$105
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"		\$198
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹			WAMS	\$275
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				\$525

1- Not available with WC-19

KIDS REAR
 Seat-to-floor Height

WHEEL SIZES		Height available		
☐	20" (451) wheel	12 1/2	to	17 1/4
☐	22" (501) wheel	13 1/2	to	18 1/4
☐	24" (540) wheel	14 1/2	to	19 1/4
Specify desired height: _____				

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest						
SWING-IN/SWING-OUT FOOTREST						
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"		Set at*:		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*:		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*:		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"		Set at*:		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*:		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*:		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" ⁵		Set at*:		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" ⁵		Set at*:		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" ⁵		Set at*:		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5" - 11"	Set at*:		NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5" - 13"	Set at*:		NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5" - 14.5"	Set at*:		NC
FOOTREST OPTIONS						
FR09	☐	High mount (footplate attached to the hanger with clamps) ³	5" - 13.5"	Set at*:		\$105
FR10	☐	Elevating legrest with calf pad, Length adjustable ⁴ Left ☐ Right ☐	14" - 22"	Set at:	WAMK	\$295
FR99	☐	Omit Footrest				NC
¹ - Available only with two piece adjustable footplates; ² - May interfere with front casters; ³ - Not available with front mount and the elevating legrest; * + or - 1/2" ⁴ - The elevating legrest replace the footrest side selected; Not compatible with seat width of 12" and 13" and one piece footplate; ⁵ - Not available with 8" caster						
FOOTPLATES						
FP01	☐	Standard - Composite ¹				STD
FP02	☐	Newton adjustable angle - Composite ^{*VP01}			WAMJ	\$155
FP03	☐	Newton adjustable angle - Aluminum ²			WAMJ	\$175
FP04	☐	Newton adjustable angle - Carbon fiber (with rubber inserts)			WAMJ	\$395
FP23	☐	Newton adjustable angle extra deep - Aluminum			WAMJ	\$175
FP06	☐	One piece adjustable angle flip-up - Aluminum		WAMJ	SEFND0015	\$215
FP07	☐	One piece adjustable angle flip-up - Carbon fiber (with rubber inserts)		WAMJ	SEFND0015	\$480
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop (not available with one piece footplate)				NC
FP09	☐	Calf strap			SEFND0010	\$57
FP30	☐	Bodypoint Aeromesh Padded calf strap			SEFND0010	\$140
¹ - Not available with front mount; ² - The Newton adjustable angle (Aluminum) is the only available option with a 12 and 13" wide frame						
SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	☐	Residual limb support - flat 6"x8" Left ☐ Right ☐			SEFND0005	\$505
FR28	☐	Residual limb support - flat 8"x10" Left ☐ Right ☐			SEFND0005	\$505
FR29	☐	Residual limb support - flat 10"x10" Left ☐ Right ☐			SEFND0005	\$505
FR30	☐	Residual limb support - flat 10"x14" Left ☐ Right ☐			SEFND0005	\$505
FR32	☐	Residual limb support - contour 6"x8" Left ☐ Right ☐			SEFND0005	\$505
FR33	☐	Residual limb support - contour 8"x10" Left ☐ Right ☐			SEFND0005	\$505
PUSH HANDLES OPTIONS (ONE option at a time)						
PH01	☐	Push handles				STD
PH05	☐	Fold down push handles ¹				\$165
PH06-920	☐	Clamp on height adjustable push handles 9" - 20" ^{1,2}				\$345
PH06-923	☐	Clamp on height adjustable push handles 9" - 23" ^{1,2}				\$345
PH02	☐	Omit push handles				NC
STROLLER HANDLE AND STABILIZER BAR (only available with push handle)						
PH04	☐	Stroller handle			WAMW	\$220
PH03	☐	Newton folding stabilizer bar (not available with stroller handle)				\$275
¹ - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock; ² - Not compatible with depth adjustable back, 8" bend back canes and amputee axle plates						
BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				NC
BH08	☐	Adjustable tension Active-Mesh back upholstery			WAMD/WAME	\$260
BH99	☐	Omit Back upholstery				STD
	☐	NXT back Support			See price list	
SEAT CUSHION						
SC01	☐	2" Cushion				\$185
SC02	☐	3" Cushion				\$205
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See price list	

Rear wheels													
ONE ARM DRIVE ¹						Only available with Newton One rear wheels, permanent axles, 0° camber							
HR22	☐	One arm drive - aluminum inner handrim ²				Left ☐	Right ☐			WAML	\$1,480		
HR60	☐	One arm drive - plastic coated inner handrim ²				Left ☐	Right ☐			WAML	\$1,595		
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)									NC		
REAR WHEELS									Available				
RW01	☐	Mag ³							20", 22" 24"		\$140		
RW02	☐	Newton One - spoke wheel ⁴							20", 22" 24"		NC		
RW03	☐	Newton Gravity - ultralight wheel							20", 22", 24"		\$390		
RW04	☐	Spineray Spox	Blk STD ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	22", 24"		\$875			
RW05	☐	Spineray LX	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐	22", 24"	\$1,095
RW07	☐	Spineray CLX ⁶	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐	24"	\$1,650
TIRES									Available				
AIRLESS													
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black							20", 22" 24"		STD		
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black ⁷							22", 24"		\$120		
TI04	☐	Soft Urethane 1" - No tread, black - <i>Shox</i>							22", 24"		\$140		
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸							22", 24"		\$225		
PNEUMATIC													
TI05	☐	Pneumatic 1 3/8" - Medium Tread							20", 22", 24"		NC		
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>							20", 22", 24"		\$120		
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>							24"		\$200		
TI03	☐	Pneumatic 2" - Large Tread, black ⁹							22", 24"		\$130		


Handrims					
HANDRIM			Available		
HR01	☐	Aluminum anodized	20", 22", 24"		STD
HR61	☐	Aluminum black hard anodized ²	20", 22", 24"		\$85
HR02	☐	Plastic coated ¹	20", 22", 24"	WAMO	\$170
HR03	☐	High friction coated ^{1,2,3}	22", 24"	WAMO	\$460
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ² <i>See on page 6 for color choice and price</i>	20", 22", 24"	WAMO	
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)	22", 24"		\$485
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)	22", 24"		\$515
HR04	☐	Natural Fit ^{1,2,3,4}	22", 24"		\$320
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)	22", 24"		\$485
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)	22", 24"		\$515
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)	22", 24"		\$320
HR99	☐	Omit handrim			NC

Handrims option					
ASSEMBLY POSITION					
HR29	☐	Standard Handrim position	Full tab		STD
HR030	☐	Narrow Handrim position no cut	Full tab		NC
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions		\$29
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐ Short tab 1 position ☐		\$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	☐	Standard grip			\$200
HR27	☐	Super grip			\$230
HR099	☐	Omit thumb piece			NC

AXLE ¹						
AX01	☐	Quick release axle		WAMU	STD	
AX02	☐	Permanent axle			NC	
AX03	☐	Quad release axle		WAMU	\$180	
AXLE PLATE						
AX04	☐	Amputee adjustable axle plate ²		WAMT	\$275	

1- Not required with one arm drive, 2- May interfere with swing away armrest and limit seat to floor height, Not available with WC-19 transit tie down				
REAR WHEEL CAMBER				
CM00	☐	0° camber		STD
CM03	☐	3° camber		NC
CM06	☐	6° camber (not available with amputee axle plate)		\$85

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													
AMPUTEE ²	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
X =	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

WHEEL LOCK								
WL08	☐	Newton Wheel lock - push to lock	PRODUCT IMPROVED					STD
WL16	☐	Aluminum push to lock with extensions						\$73
WL02	☐	Pull to lock						NC
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)					WAMR	\$185
WL13	☐	6" removable extension handle ¹ *VP01						\$73
WL14	☐	Scissors Lock						\$185
WL05	☐	Attendant Lock (additional set of wheel lock mounted to the back)						\$240
WL06	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐			WAMN	\$335
WL15	☐	Omit wheel locks						NC
WL99	☐							

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests						
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01						
AR13-610	☐	10" Desk length armpad ⁴	6" - 10"			\$380
AR03-812	☐	10" Desk length armpad	8" - 12"			\$380
AR03-1014	☐	10" Desk length armpad	10" - 14"			\$380
AR14-610	☐	14" Full length armpad ⁴	6" - 10"			\$380
AR04-812	☐	14" Full length armpad	8" - 12"			\$380
AR04-1014	☐	14" Full length armpad	10" - 14"			\$380
LOCKING FLIP-UP ARMREST ²						
AR07	☐	Locking Flip-Up Armrest				\$350
ARMREST OPTIONS						
AR12	☐	14" Full Length armpad - Gel Ovation (pair) ³			SEAND0030	\$160
AR99	☐	Omit armrest				NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- May interfere with rigid back hardware; Only available with 18"-21" straight aluminum back canes; Not available with back upholstery;

3- Not available with locking flip-up armrest and short "T" armrest; 4- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with "T" Armrest)								
SG01	☐	Plastic side guards (removable)	Straight bracket - STD	Offset bracket ¹			WAMV	\$235
SG03	☐	Carbon fiber side guards (removable)	Straight bracket - STD	Offset bracket ¹			WAMV	\$320
SG06	☐	Plastic fender side guards (removable)	Straight bracket - STD	Offset bracket ¹			WAMV	\$340
SG04	☐	Carbon fiber fender side guards (removable)	Straight bracket - STD	Offset bracket ¹			WAMV	\$510
SG99	☐	Omit side guard						NC



1- only compatible with 0° camber; Increase width by 1/2" per side for a total of 1"

POSITIONING BELT								
SB01-S	☐	Velcro adjustable - Small 35"						\$50
SB01-L	☐	Velcro adjustable - Long 60"						\$50
SB02-S	☐	Auto buckle - Small 35"						\$60
SB02-M	☐	Auto buckle - Medium 48"						\$60
SB02-L	☐	Auto buckle - Long 60"						\$60
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"						\$95
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"						\$95
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"				SERND0030		\$185
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"				SERND0030		\$185
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"				SERND0030		\$185
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"				SERND0020		\$210
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"				SERND0020		\$210
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"				SERND0020		\$210
SB03	☐	Hardware attachment for 2 extra points (1 pair)						\$15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"						\$355
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"						\$355
SB99	☐	Omit positioning belt						STD

HEADREST ¹						
AC01-S	☐	Profiled Headrest - Small (3.5"X6")				\$215
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")				\$215
AC01-L	☐	Profiled Headrest - Large (5.5"X8")				\$215
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar ¹				\$360
1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar						
ACCESSORIES						
AC03	☐	Anti-tippers				\$120
AC29	☐	Swing-away anti-tipper ¹	Left ☐ Right ☐			\$185
AC27	☐	Tip assist ²				\$65
AC06	☐	Cane and crutch holder	Left ☐ Right ☐			\$115
AC07	☐	Oxygen tank holder ⁴ - D and E Type - installed on the left				\$395
AC08	☐	Spoke guards ³	Available: 20", 22", 24", 26"			\$155
AC09	☐	Footrest neoprene impact guard (<i>pair</i>) - Black				\$85
AC10	☐	Seat pouch				\$100
AC30	☐	Newton multi-purpose pouch				\$75
AC26	☐	BodyPoint Mobility Bag				\$145
AC23	☐	Newton Backpack				\$130
AC12	☐	Velcro Closing Strap				\$55
AC13	☐	Tool kit				\$55
AC16	☐	Motion Composites Polo Shirt	W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐			\$36

1- Must select the pair (left and right) with HD optionè Only available with rear seat to floor height between 13 1/2" and 19 3/4"
2- Cannot be combined with Anti Tippers


3- Not available with Mag wheels and Spinerigy; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down

Colors																					
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
																					
AirGrip	☐ \$160		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	\$350
Frame ²		☐ STD	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	NC
FC99	☐	Custom Color Specify PPG Paintit color: additional 10 days lead-time																			\$422

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- with black hydroformed aluminum crossbrace

Fork and Stem bolt combinations



Vision Package Accents



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :

BioFit™ CUSHION		WIDTH					
		12"		13"		14"	
DEPTH	12"	\$625	\$625	\$625			
	14"	\$625	\$625	\$625			
	15"		\$625	\$625	\$625	\$625	
	16"		\$625	\$625	\$625	\$625	
	17"				\$625	\$625	
	18"				\$625	\$625	
	19"				\$625	\$625	
Option: Change cover for incontinent / infection control outer cover							N/C

Green boxes are stocked sizes and the red boxes require customization.

SEATING AND CUSHION ACCESSORIES

	Model #	Description	Fits Back Width	\$	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6 & 43.2 cm) and 14"-18" deep (35.6-45.7 cm)	80	

	Optima™ THORACIC BACK SUPPORT
	CARBON FIBER SHELL (WEIGHT SAVINGS - 1lb)
	Optima™ DEEP THORACIC BACK SUPPORT
	CARBON FIBER SHELL (WEIGHT SAVINGS - 1lb)

THORACIC SHAPE

WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM
13"	07"	NOCFACB1307	\$1150	NOACB1307	\$720
13"	09"	NOCFACB1309	\$1150	NOACB1309	\$720
13"	13"	NOCFACB1313	\$1150	NOACB1313	\$720
13"	16"	NOCFACB1316	\$1150	NOACB1316	\$720
14"	13"	NOCFACB1413	\$1150	NOACB1413	\$720
14"	16"	NOCFACB1416	\$1150	NOACB1416	\$720
15"	13"	NOCFACB1513	\$1150	NOACB1513	\$720
15"	16"	NOCFACB1516	\$1150	NOACB1516	\$720
16"	13"	NOCFACB1613	\$1150	NOACB1613	\$720
16"	16"	NOCFACB1616	\$1150	NOACB1616	\$720
16"	19"	NOCFACB1619	\$1150	NOACB1619	\$720

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	NOACC _ _ _ _	130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	xtend™	WIDTH	HEIGHT	✓	\$
	**NHAACB (L)1213/17	12"	13"-17"		\$725
	NHAACB (L)1413/17	14"	13"-17"		\$725
	NHAACB (L)1613/17	16"	13"-17"		\$725

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	N3LHAACC _ _ _ _	130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

HELIO FOLDING ON CHAIR – CDN ORDER FORM

BACK SUPPORT ACCESSORIES

	MODEL #	DESCRIPTION	FITS BACK WIDTH	\$	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	205	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	70	


MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

											
MODEL # NIB-0784	✓	\$	MODEL # NIB-0818	✓	\$	MODEL # NIB-0747	✓	\$	MODEL # NIB-0823	✓	\$
10" x 5" Large Rectangular		\$205	14" x 5" Extra-Large Rectangular		\$205	Multi-Axis Hardware		\$150	Head Support Receiver Bracket		\$60

MultiFit™ HEAD SUPPORT – FIXTURES

EXTRA COVERS

			
Extra Outer smartx ^{3D} Cover		Extra Outer Incontinent Cover	
ITEM	WIDTH X DEPTH	\$	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	___ x ___	130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	___ X ___	150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	___ X ___	130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	___ X ___	150	



Solid Seat Insert



Pelvic Obliquity Wedges

SEATING AND CUSHION ACCESSORIES

MODEL NO.	DESCRIPTION	FITS BACK WIDTH	\$	QTY
NIB-0711	Pair of 0.5" (1.2 cm) pelvic build up wedges	All Sizes	80	
RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6 & 43.2 cm) and 14"-18" deep (35.6-45.7 cm)	80	
RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7 & 48.3 cm) and 14"-18" deep (35.6-45.7 cm)	80	
RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	80	

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com