


1- Purchasing information

P.O. number: * Date:
Company Name: *
Purchasing contact: *
Phone: *
Email: *
Contact preference: * ☐ Phone ☐ Email ☐
Marked for:
2- Dealer/Representative information

Dealer/Representative: *
Phone: *
Email: *
Contact preference: * ☐ Phone ☐ Email ☐
3- Shipping information

Ship to: *
Name:
Address:
City: Province:
Postal Code:
Phone:

* Required

MA080 ☐
English Manual (STD)

MA081 ☐
Manuel Français

Taxes not included

MODEL								
WC03	☐	Ultralight Folding Carbon Fiber Helio C2 Wheelchair (265 lb weight capacity)						\$3,545
WC09	☐	HD Heavy Duty package (350 lb weight capacity)						\$170
FRAME								
SEAT WIDTH ³								
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐		NC
	☐	19" ¹ ☐	20" ¹ ☐	21" ² ☐	22" ² ☐			\$338
SEAT DEPTH								
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐		NC
	☐	19" ☐	20" ☐					\$255
FRAME LENGTH ⁵								
SP01	☐	Standard frame length configuration						STD
SP02	☐	Frame length 1" longer than standard configuration						NC
SP03	☐	Frame length 2" longer than standard configuration						NC
SP04	☐	Frame length 1" shorter than standard configuration						NC
SP05	☐	Frame length 2" shorter than standard configuration						NC
TRANSIT TIE DOWN								
AC04	☐	Transit Tie Down (Unoccupied)						NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴		Short belt 20"-41"				\$345
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴		Long belt 20"-46"				\$345
AC35	☐	WC-19 Transit Tie Down without belt ⁴						\$325

1- Available with HD option; 2- Only available with HD option

4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

5- See chart on page 7

[3. Width calculator](#)


PACKAGES							
PK01	☐	Newton Ultralight Package (include: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels)				Save \$245	\$540

C2 FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2				14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						
	P4	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
7"	P1										16 1/2"	17 1/2"	18 1/2"	17"	18"	19"			
	P2							16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"			
	P3							17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
	P4	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"
	P5	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"	19 1/2"	20 1/2"	21 1/2"
Froglegs	P1				14 1/2"	15 1/2"	16 1/2"												
	P2	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"									
	P3	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"						

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

FRONT CASTERS						
CASTERS				Available		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 ¼")		6", 8"		160 \$
CA14	☐	Newton UltraCaster Composite	4" x 1 ½ "			120 \$
CA15	☐		5" x 1 ½ "			130 \$
CA16	☐		6" x 1 ½ "			140 \$
CA17	☐	Newton UltraCaster Aluminum	4" x 1 ½ "			175 \$
CA18	☐		5" x 1 ½ "			185 \$
CA19	☐		6" x 1 ½ "			195 \$
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹				\$255
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)				\$489

1- Not available with WC-19

C2 - REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height: _____
☐	20" (451) wheel	12 1/2	to 17 1/4	
☐	22" (501) wheel	13 1/2	to 18 1/4	
☐	24" (540) wheel	14 1/2	to 19 1/4	
☐	25" (559) wheel	15	to 19 3/4	
☐	26" (590) wheel	15 1/2	to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

REAR WHEELS							Available						
RW01	☐	Mag ¹					20", 22" 24"		\$129				
RW02	☐	Newton One - spoke wheel ²					20", 22" 24"		NC				
RW03	☐	Newton Gravity - ultralight wheel					20", 22", 24", 25", 26"		\$360				
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ³ ☐	Blue ³ ☐	Red ³ ☐	Yel ³ ☐	22", 24", 25", 26"	\$800				
RW05	☐	Spinergy LX	Blk STD ☐	Pur ³ ☐	Wht ³ ☐	Blue ³ ☐	Red ³ ☐	Yel ³ ☐	Pink ³ ☐	Grn ³ ☐	Org ³ ☐	22", 24", 25", 26"	\$995
TIRES							Available						
AIRLESS													
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black ⁶					20", 22" 24"		STD				
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black ^{4,6}					22", 24", 25", 26"		\$110				
TI04	☐	Soft Urethane 1" - No tread, black - <i>Shox</i>					22", 24", 25", 26"		\$130				
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁵					22", 24", 26"		\$215				
PNEUMATIC													
TI05	☐	Pneumatic 1 3/8" - Medium Tread ⁶					20", 22", 24", 25" 26"		NC				
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>					20", 22", 24", 25", 26"		\$110				
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>					24", 25", 26"		\$190				

1- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 2- Not compatible with Shox tires; 3- color spokes require an additional 10 days lead-time;

4- 22" is only available with Mags and Newton One rear wheels; 25" and 26" are not available with Newton Gravity wheels

5- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels 6- Not available with flip-back armrest when the 6" camber is selected;

HANDRIMS					
HANDRIM				Available	
HR01	☐	Aluminum anodized		20", 22", 24", 25", 26"	STD
HR02	☐	Plastic coated ¹		20", 22", 24", 25", 26"	\$160
HR03	☐	High friction coated ^{1,2,3}		22", 24", 25", 26"	\$430
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ²	See on page 6 for color choice and price	20", 22", 24", 25", 26"	
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"	\$470
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)		22", 24", 25", 26"	\$500
HR04	☐	Natural Fit ^{1,2,3,4}		22", 24", 25", 26"	\$310
HR99	☐	Omit handrim			NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;

4 - 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

HANDRIMS OPTION					
ASSEMBLY POSITION					
HR29	☐	Standard Handrim position	Full tab		STD
HR030	☐	Narrow Handrim position no cut	Full tab		NC
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions		\$25
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐ Short tab 1 position ☐		\$25
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS					
HR23	☐	Standard grip			\$190
HR27	☐	Super grip			\$220
HR099	☐	Omit thumb piece			NC


5- Only available with Aluminum and Air Grip handrims;

AXLE ¹					
AX01	☐	Quick release axle			STD
AX02	☐	Permanent axle			NC
AX03	☐	Quad release axle			\$170
AXLE PLATE					
AX04	☐	Amputee adjustable axle plate ²			\$250

1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

REAR WHEEL CAMBER					
CM00	☐	0° camber			STD
CM03	☐	3° camber			NC
CM06	☐	6° camber (not available with amputee axle plate)			\$80

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

	Most stable														
X =	AMPUTEE ¹	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"	
	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2.5".

Possible conflict between a 26" rear wheel and a center of gravity over 2.5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

WHEEL LOCK					
WL08	☐	Newton Wheel lock - push to lock			STD
WL02	☐	Pull to lock			NC
WL04	☐	Newton grade aid push to lock (works best with pneumatic or high pressure tire)			\$175
WL13 WL14	☐	6" removable extension handle			\$67
WL05	☐	Scissors Lock			\$300
WL06	☐	Attendant Lock (additional set of wheel lock mounted to the back)			\$230
WL99	☐	Omit wheel locks			NC

FOOTREST						
SWING-IN/SWING-OUT FOOTREST						
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"		Set at*:		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*:		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*:		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"		Set at*:		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"		Set at*:		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"		Set at*:		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" ⁵		Set at*:		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" ⁵		Set at*:		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" ⁵		Set at*:		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5" - 11"	Set at*:		NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5"- 13"	Set at*:		NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5" - 14.5"	Set at*:		NC
FOOTREST OPTIONS ³						
FR09	☐	High mount <i>(footplate attached to the hanger with clamps)</i> ³	5" - 13.5"	Set at*:		\$90
FR10	☐	Elevating legrest with calf pad, Length adjustable ⁴	Left ☐ Right ☐ 14" - 22"	Set at:		\$290
FR99	☐	Omit Footrest				NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating leagrest; * +/- 1/2"
4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES						
FP01	☐	Standard - Composite ¹				STD
FP02	☐	Newton adjustable angle - Composite				\$145
FP03	☐	Newton adjustable angle - Aluminum				\$165
FP04	☐	Newton adjustable angle - Carbon fiber <i>(with rubber inserts)</i>				\$380
FP06	☐	One piece adjustable angle flip-up - Aluminum				\$206
FP07	☐	One piece adjustable angle flip-up - Carbon fiber <i>(with rubber inserts)</i>				\$450
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop ²				NC
FP09	☐	Calf strap				\$52
FP10	☐	Standard Padded Calf strap				\$80
FP30	☐	Bodypoint Aeromesh Padded calf strap				\$120

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	☐	Residual limb support - flat 6"x8"	Left ☐ Right ☐			\$390
FR28	☐	Residual limb support - flat 8"x10"	Left ☐ Right ☐			\$390
FR29	☐	Residual limb support - flat 10"x10"	Left ☐ Right ☐			\$390
FR30	☐	Residual limb support - flat 10"x14"	Left ☐ Right ☐			\$390
FR32	☐	Residual limb support - contour 6"x8"	Left ☐ Right ☐			\$390
FR33	☐	Residual limb support - contour 8"x10"	Left ☐ Right ☐			\$390

SEATING AND BACKREST						
BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				NC
BH08	☐	Adjustable tension Active-Mesh back upholstery				\$245
BH99	☐	Omit Back upholstery				STD
	☐	NXT back Support			See price list	
SEAT CUSHION						
SC01	☐	2" Cushion				\$175
SC02	☐	3" Cushion				\$195
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See price list	
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>						
BC06-912	☐	Straight Carbon fiber back cane 9" - 12"	<i>increment of 1/2"</i>	Set at ²		NC
BC06-1215	☐	Straight Carbon fiber back cane 12" - 15"	<i>increment of 1/2"</i>	Set at ²		NC
BC06-1518	☐	Straight Carbon fiber back cane 15" - 18"	<i>increment of 1/2"</i>	Set at ²		NC
BC06-1821	☐	Straight Carbon fiber back cane 18" - 21"	<i>increment of 1/2"</i>	Set at ²		NC
BC08-1315	☐	8° Bend Carbon fiber back cane 13.5" - 15"	<i>increment of 1/2"</i>	Set at ²		\$100
BC08-1518	☐	8° Bend Carbon fiber back cane 15" - 18"	<i>increment of 1/2"</i>	Set at ²		\$100
BC08-1821	☐	8° Bend Carbon fiber back cane 18" - 21"	<i>increment of 1/2"</i>	Set at ²		\$100
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>						
BC09-1215	☐	Straight Carbon fiber back cane 12" - 15"	<i>increment of 1/2"</i>	Set at ²		\$250
BC09-1518	☐	Straight Carbon fiber back cane 15" - 18"	<i>increment of 1/2"</i>	Set at ²		\$250
BC09-1821	☐	Straight Carbon fiber back cane 18" - 21"	<i>increment of 1/2"</i>	Set at ²		\$250
BC10-1315	☐	8° bend Carbon fiber back cane 13.5" - 15"	<i>increment of 1/2"</i>	Set at ²		\$325
BC10-1518	☐	8° bend Carbon fiber back cane 15" - 18"	<i>increment of 1/2"</i>	Set at ²		\$325
BC10-1821	☐	8° bend Carbon fiber back cane 18" - 21"	<i>increment of 1/2"</i>	Set at ²		\$325
		Select back cane angle:	85° ¹ ☐ 90° ☐ 95° STD ☐ 100° ☐ 105° ☐ 110° ☐			

PUSH HANDLES OPTIONS <i>(ONE option at a time)</i>						
PH01	☐	Push handles				STD
PH05	☐	Fold down push handles ^{3, 5}				\$150
PH06-920	☐	Clamp on height adjustable push handles 9"-20" ^{3,4, 5}				\$325
PH06-923	☐	Clamp on height adjustable push handles 9"-23" ^{3,4, 5}				\$325
PH02	☐	Omit push handles				NC
STROLLER HANDLE AND STABILIZER BAR <i>(only available with push handle)</i>						
PH04	☐	Stroller handle <i>(not available with the HD option)</i>				\$210
PH03	☐	Newton folding stabilizer bar <i>(not available with stroller handle)</i>				\$258

1- Not available with 8 degree bend cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265 lb weight limit;
4- Not compatible with swing away armrests, 8" bend back canes, 9-12" back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

ARMRESTS						
HEIGHT ADJUSTABLE FLIP BACK ARMREST <i>(Convertible to "T" armrest) ^{1, 4}</i>						
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐	\$380
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐	\$380
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐	\$380
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐	\$380
HEIGHT ADJUSTABLE "T" ARMREST ¹						
AR13-610	☐	10" Desk length armpad ⁵	6" - 10"			\$360
AR03-812	☐	10" Desk length armpad	8" - 12"			\$360
AR03-1014	☐	10" Desk length armpad	10" - 14"			\$360
AR14-610	☐	14" Full length armpad ⁵	6" - 10"			\$360
AR04-812	☐	14" Full length armpad	8" - 12"			\$360
AR04-1014	☐	14" Full length armpad	10" - 14"			\$360
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²						
AR05	☐	Tubular armrest - short pad	8" - 11"			\$295
AR23	☐	Tubular armrest - long pad	8" - 11"			\$295
AR06	☐	10" Desk length flat armpad	8" - 11"			\$295
AR17	☐	14" Full length flat armpad	8" - 11"			\$295
ARMRESTS OPTIONS						
AR12	☐	14" Full Length armpad - Gel Ovation (pair) ³				\$150
AR99	☐	Omit armrest				NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;
3- Only available with full length armrest; Not available with swing away and desk lenght armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

ACCESSORIES						
SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>						
SG01	☐	Plastic side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$225
SG02	☐	Carbon fiber side guards <i>(fixed)</i>				\$255
SG03	☐	Carbon fiber side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$295
SG06	☐	Plastic fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$320
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$495
SG99	☐	Omit side guard				NC

1- only compatible with 0° camber;Increase width by 1/2" per side for a total of 1" ;

POSITIONING BELT						
SB01-S	☐	Velcro adjustable - Small 35" ¹				\$45
SB01-L	☐	Velcro adjustable - Long 60"				\$45
SB02-S	☐	Auto buckle - Small 35" ¹				\$55
SB02-M	☐	Auto buckle - Medium 48"				\$55
SB02-L	☐	Auto buckle - Long 60"				\$55
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"				\$82
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"				\$82
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"				\$155
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"				\$155
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"				\$155
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"				\$195
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"				\$195
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"				\$195
SB03	☐	Hardware attachment for 2 extra points (1 pair)				\$13
SB99	☐	Omit positioning belt				STD

1- Not recommended for seat width of 18" and more

HEADREST ¹						
AC01-S	☐	Profiled Headrest - Small (3.5"X6")				\$205
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")				\$205
AC01-L	☐	Profiled Headrest - Large (5.5"X8")				\$205
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar ¹				\$338

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES						
WC10	☐	E-Motion adaptor Kit (2 adaptors) - Only necessary with 6° camber				\$75
AC03	☐	Anti-tippers				\$115
AC06	☐	Cane and crutch holder	Left ☐	Right ☐		\$105
AC07	☐	Oxygen tank holder ¹ - D and E Type - installed on the left				\$375
AC08	☐	Spoke guards ² Available: 20", 22", 24", 25", 26"				\$145
AC09	☐	Footrest neoprene impact guard (pair) - Black				\$80
AC10	☐	Seat pouch				\$95
AC23	☐	Newton Backpack				\$120
AC12	☐	Velcro Closing Strap				\$50
AC13	☐	Tool kit				\$50


1- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down 2- Not available with Mag wheels and Spinergy; the 24" spoke guard will be installed on the 25" wheel ;

Colors																					
Options ¹																					
	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
																					
AirGrip	☐		☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	\$330
Frame ²		STD	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	NC
FC99	☐	Custom Color Specify PPG Paintit color: additional 10 days lead-time																			\$422

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- With carbon fiber crossbrace;

Fork and Stem bolt combinations



Frame length	Seat Depth						
	14"	15"	16"	17"	18"	19"	20"
STD							
+1							N/A
+2						N/A	N/A
-1	N/A						
-2	N/A	N/A					

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :