



2915 Ogletown RD # 2270
Newark, DE 19713
T.+1-866-650-6555 F. +1-888-966-6555
orders@motioncomposites.com
www.motioncomposites.com



5.1



AWARD 2010
HARDING 2013
INNOVATIVE 2014
PRODUCT 2016



Helio
ORDER FORM - USA

P.O. number: *	Date:	Ship to: *
Company Name: *	Phone: *	Name:
Ordered by: *		Address:
Clinician/ATP/RTS:		City: State:
Email:		Zip Code:
Marked for:		Phone:

* Required MA082 ☐ English Manual (STD) MA089 ☐ Spanish Manual Taxes not included

MODEL								USD			
1	WC03	Ultralight folding Helio wheelchair (265lb weight capacity) A7 package: side frames, symmetrical hydroformed crossbrace and back canes made of 7000 series aluminum						K0005	\$2 295	☐	
2	WC08	C2 carbon fiber upgrade: side frames, symmetrical molded crossbrace and back canes made of ultralight carbon fiber							\$1 495	☐	
3	WC09	HD Heavy Duty package (350 lb weight capacity)							\$500	☐	
4	WC10	E-Motion adaptor Kit (2 adaptors) - <i>Only necessary with 6° camber</i>							\$30	☐	
SEAT WIDTH ³											
5	SW	Seat width:	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐		NC	☐
6	SW	Seat width:		20" ¹ ☐	21" ² ☐	22" ² ☐			E2201	\$495	☐
SEAT DEPTH											
7	SD	Seat depth:	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐		NC	☐
8	SD	Seat depth:		20" ☐					E2203	\$475	☐

1 - Available with HD Kit (350 lb)

2 - Only with HD Kit (350 lb)

[3. Width calculator](#)

NEWTON ULTRALIGHT PACKAGE													
9	PK01	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels								Save \$235	\$480	☐	

Helio FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13	14	15	13 1/2	14 1/2	15 1/2												
	P2	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2									
5"	P1																		
	P2				14	15	16	14 1/2	15 1/2	16 1/2									
	P3	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2						
	P4	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
7"	P1										16 1/2	17 1/2	18 1/2	17	18	19			
	P2							16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2			
	P3							17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2
	P4	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21
	P5	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2	19	20	21	19 1/2	20 1/2	21 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS												
10	CA01	Composite wheel w/ PU tire (1" width)	3" ☐	4" ☐	5" ☐	6" ☐	7" ☐	8" ☐			NC	☐
11	CA10	Pneumatic casters 6" x 1 ¼"								E2214/E2215	\$160	☐
12	CA11	Pneumatic casters 8" x 1 ¼"								E2214/E2215	\$160	☐
13	CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "								\$80	☐
14	CA15		5" x 1 ½ "								\$80	☐
15	CA16		6" x 1 ½ "								\$80	☐
16	CA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½ "								\$175	☐
17	CA18		5" x 1 ½ "								\$185	☐
18	CA19		6" x 1 ½ "								\$195	☐
CASTER OPTIONS												
19	CA12	Caster Pin Locks (<i>Available with +1" or +2" Stem bolt ONLY</i>)							K0073	\$295	☐	
20	CA13	Carbon fiber Frog Legs suspension forks (<i>Not available with caster pin locks</i>)							E1015	\$475	☐	

Helio - REAR Seat-to-floor Height

WHEEL SIZES	Height available	Specify desired height:
20" (451) wheel ☐	12 1/2 to 17 1/4	
22" (501) wheel ☐	13 1/2 to 18 1/4	
24" (540) wheel ☐	14 1/2 to 19 1/4	
25" (559) wheel ☐	15 to 19 3/4	
26" (590) wheel ☐	15 1/2 to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest												
SWING-IN/SWING-OUT FOOTREST												
21	FR03	60° Swing-in/swing-out footrest	13.5 - 17" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐			Set at: *			NC	☐
22	FR04	70° Swing-in/swing-out footrest	13.5 - 17" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐			Set at: *			NC	☐
23	FR15	80° Swing-in/swing-out footrest	13.5 - 16" ☐	14.5 - 19.5" ☐	15.5 - 20.5" ☐			Set at: *			NC	☐
90° SWING-IN/SWING-OUT FOOTREST WITH FRONT MOUNT (Footplate mounted to hanger tube)												
24	FR06	90° swing-in/swing-out with front mount - Standard (no extension) ¹			7.5 - 11"			Set at: *			NC	☐
25	FR07	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}			9.5 - 13"			Set at: *			NC	☐
26	FR08	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}			9.5 - 14.5"			Set at: *			NC	☐
FOOTREST OPTIONS ³												
27	FR09	High mount (footplate attached to the hanger with clamps)		5 - 13.5" ☐				Set at: *			\$75	☐
28	FR99	Omit Footrest									NC	☐

1 - Available only with two pieces adjustable footplates 2 - May interfere with front casters 3- Not available with front mount * + / - 1/2"

ELEVATING LEGREST (unit)								The elevating legrest replace the footrest side selected				
29	FR10	Elevating legrest with calf pad, 14" - 22" Length adjustable		Left ☐		Right ☐		Set at:		E0990	\$240	☐

FOOTPLATES												
30	FP01	Standard - Composite									STD	☐
31	FP02	Newton adjustable angle - Composite								K0040	\$210	☐
32	FP03	Newton adjustable angle - Aluminum								K0040	\$230	☐
33	FP04	Newton adjustable angle - Carbon fiber (with rubber inserts)									\$415	☐
34	FP06	One piece adjustable angle flip-up - Aluminum								K0040	\$190	☐
35	FP07	One piece adjustable angle flip-up - Carbon fiber (with rubber inserts)									\$435	☐

FOOTPLATE ACCESSORIES												
36	FP08	Length adjustable heel loop (not available with one piece footplate)								E0951	\$40	☐
37	FP09	Calf strap								K0038	\$35	☐
38	FP30	Bodypoint Aeromesh Padded calf strap								SEFND0010	\$115	☐

RESIDUAL LIMB SUPPORT FLAT CUSHION (Swing away, depth, width, height and angle ajustable) - (unit)												
39	FR27	Residual limb support swing away - flat 6"x8"		Left ☐		Right ☐		New		E1020	\$405	☐
40	FR28	Residual limb support swing away - flat 8"x10"		Left ☐		Right ☐		New		E1020	\$405	☐
41	FR29	Residual limb support swing away - flat 10"x10"		Left ☐		Right ☐		New		E1020	\$405	☐
42	FR30	Residual limb support swing away - flat 10"x14"		Left ☐		Right ☐		New		E1020	\$405	☐

RESIDUAL LIMB SUPPORT CONTOUR CUSHION (Swing away, depth, width, height and angle ajustable) - (unit)												
43	FR32	Residual limb support swing away - contour 6"x8"		Left ☐		Right ☐		New		E1020	\$405	☐
44	FR33	Residual limb support swing away - contour 8"x10"		Left ☐		Right ☐		New		E1020	\$405	☐

Back canes, height adjustable, with option									
FIXED ANGLE BACK									
45	BC11	Straight height adjustable back cane	12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐			NC	☐
	BC06								
46	BC12	8° Bend height adjustable back cane	13 ½ - 15" ☐	15 - 18" ☐	18 - 21" STD ☐			\$50	☐
	BC08								
ADJUSTABLE ANGLE BACK									
47	BC13	Adjustable angle, straight back cane		15 - 18" ☐	18 - 21" STD ☐			\$200	☐
	BC09								
48	BC14	Adjustable angle, 8° bend back cane		15 - 18" ☐	18 - 21" STD ☐			\$250	☐
	BC10								
49		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐	
Set height to: ²		(measured from seat sling to top of back upholstery)							
		PUSH HANDLE OPTION							
50	PH01	Push handles						STD	☐
51	PH05	Folding push handles ^{3, 5}						\$120	☐
52	PH06	clamp on height adjustable push handle ^{3, 4, 5}	9 - 20" STD ☐	9-23" ☐				\$300	☐
53	PH02	Without push handles						NC	☐

1- Not available with a 8 degree bend cane 2- If not indicated, default height setting is at the lowest position 3 - 265 pounds weight limit

4 - Not compatible with swing away armrests, 8° bend back canes and amputee axle plates 5 - Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

STROLLER HANDLE AND STABILIZER BAR ¹									
54	PH04	Stroller handle (only available with push handle) - (not available with the HD option)					K0108	\$160	☐
55	PH03	Newton folding stabilizer bar						\$250	☐

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY									
56	BH07	Slip-on Active-Mesh back upholstery						STD	☐
57	BH08	Adjustable tension Active-Mesh back upholstery					E2611	\$300	☐
58	BH99	Omit Back upholstery						NC	☐
59		NXT back series				see pricelist p.7			

SEAT CUSHION									
60	SC01	2" Cushion						\$165	☐
61	SC02	3" Cushion						\$185	☐
62	SC99	Omit seat cushion						STD	☐
63		NXT FIT cushion series				see pricelist p.7			

Rear wheels									
ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims									
64	HR22	One arm drive ^{2, 3}	24" ☐	left ☐	right ☐			\$1 395	☐
65	HR98	Omit handrim on opposite side of the One Arm Drive (OAD)						NC	☐
REAR WHEELS									
66	RW01	Mag ⁴	20" ☐	22" ☐	24" ☐			\$50	☐
67	RW02	Newton One - spoke wheel ⁵	20" ☐	22" ☐	24" ☐			NC	☐
68	RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐	\$300	☐
69	RW04	Spinergy Spox	22" ☐	24" ☐	25" ☐	26" ☐		\$795	☐
70		Spinergy Spox spokes colors ⁶ Blk STD ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐		NC	☐
71	RW05	Spinergy LX	22" ☐	24" ☐	25" ☐	26" ☐		\$995	☐
72		Spinergy LX spokes colors ⁶ Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC	☐
73	RW07	Spinergy CLX Carbon Blade		24" ☐	25" ☐			\$1 480	☐
74		Spinergy CLX spokes colors ⁶ Blk STD ☐	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC	☐
TIRES					Available				
75	TI01	Soft Urethane 1 1/4 "			20", 22" 24"		E2218	STD	☐
76	TI02	Soft Urethane knobby 1 3/8 "			22", 24", 25", 26"		E2218	\$95	☐
77	TI04	Soft Urethane - Shox 1"			22", 24", 25", 26"		E2218	\$115	☐
78	TI05	Pneumatic 1 3/8 "			20", 22", 24", 25", 26"		E2211/E2212	NC	☐
79	TI08	Pneumatic with Airless insert 1 3/8 "			22", 24"		E2211/E2213	\$80	☐
80	TI06	Pneumatic HP - 1"			20", 22", 24", 25", 26"		E2211/E2212	\$95	☐
81	TI07	Pneumatic HP - Schwalbe Marathon Plus 1"			24", 25", 26"		E2211/E2212	\$125	☐

1 - Offered with handrims on each side of the wheelchair 2 - Only standard handrim position available 3 - Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 5 - Not compatible with Shox tires 6 - Colour spokes require an additionnal 10 days lead-time

Handrims						
HANDRIM						
82	HR01	Aluminum anodized			STD	☐
83	HR02	Plastic coated ³			\$95	☐
84	HR03	Q-Grip High friction Neoprene coated ^{1, 2, 3}			\$350	☐
85	HR07	Newton Air Grip ¹ <i>Black</i> ☐			\$95	☐
86	HR08	Newton Air Grip ¹ (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐			\$240	☐
		(11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐				
		(15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐				
87	HR28	Surge ^{1, 2, 3} (1-1/2" Oval handrim w/ Gripton high friction strip)			\$375	☐
88	HR05	Surge LT ^{1, 2, 3} (1-1/8" Oval handrim w/ Gripton high friction strip)			\$405	☐
89	HR04	Natural Fit ^{1, 2, 3}			\$310	☐
90	HR99	Omit handrim			NC	☐

1 - Not available with Mag wheels or OAD 2 - Not available in 20" 3- Only standard handrim position available

Handrims option						
Assembly position for aluminum and Newton AirGrip handrims						
91	HR29	Standard Handrim position	Full tab		STD	☐
92	HR30	Narrow Handrim position	Short tab 2 positions		\$25	☐
93	HR31	Super Narrow Handrim	Short tab 2 positions STD ☐ Short tab 1 position ☐		\$25	☐
Thumb piece ⁴ selection for Natural Fit handrims						
94	HR23	Standard grip			\$85	☐
95	HR27	Super grip			\$140	☐
96	HR099	Omit thumb piece			NC	☐


4- Compatible with Surge, Surge LT handrims

AXLE						
97	AX01	Quick release axle			STD	☐
98	AX02	Permanent axle			NC	☐
99	AX03	Quad release axle		K0108	\$125	☐
AXLE PLATE						
100	AX04	Amputee adjustable axle plate ¹		E0959	\$195	☐

1 - May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER						
101	CM00	0° camber			STD	☐
102	CM03	3° camber			NC	☐
103	CM06	6° camber (not available with amputee axle plate)			\$150	☐

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable												Most responsive	
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"	
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2,5".

Possible conflict between a 26" rear wheel and a center of gravity over 2,5".

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.

Only on carbon frame

WHEEL LOCK						
104	WL08	Newton Wheel lock - push to lock			STD	☐
105	WL02	Pull to lock			NC	☐
106	WL04	Newton grade aid push to lock (<i>suggested with pneumatic or high pressure tire</i>)		E0974	\$150	☐
107	WL13 WL14	6" removable extension handle		E0961	\$75	☐
108	WL05	Scissor Lock			\$140	☐
109	WL06	Attendant Lock (<i>additional set of wheel lock mounted to the back</i>)			\$200	☐
110	WL99	Omit wheel locks			NC	☐

Armrests						
FLIP BACK ARMREST (<i>Convertible to "T" armrest</i>) ^{1,4}			Height Adjustment		Removable function	
111	AR01	10" Desk length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked STD ☐
112	AR02	14" Full length arm pad	8 - 12" STD ☐	10 - 14" ☐	Locked ☐	Unlocked STD ☐
"T" ARMREST ¹			Height Adjustment			
113	AR03	10" Desk length arm pad	8 - 12" ☐ STD	10 - 14" ☐		E0973
114	AR04	14" Full length arm pad	8 - 12" ☐ STD	10 - 14" ☐		E0973
SWING AWAY ARMREST ²			Height Adjustment			
115	AR05	Tubular arm pad	8 - 11"			E0973
116	AR06	10" Desk length arm pad	8 - 11"			E0973
GEL ARM PAD ³						
117	AR12	14" Full Length Arm Pad - Gel Ovation				\$150
118	AR99	Omit armrest				NC

1 - May have a 1/2 " adjustment restriction on both limits of the range 2 - Could interfere with some rigid back rests 3 - Not available with swing away armrest 4 - Not available with attendant lock

SIDE GUARDS (<i>not available with Flip Back and "T" Armrest</i>)						
119	SG01	Plastic side guards (<i>removable</i>)				\$195
120	SG02	Carbon fiber side guards (<i>fixed</i>)				\$250
121	SG03	Carbon fiber side guards (<i>removable</i>)				\$290
122	SG06	Plastic fender side guards (<i>removable</i>) ¹				\$270
123	SG04	Carbon fiber fender side guards (<i>removable</i>)				\$540
124	SG99	Omit side guard				STD
SIDE GUARDS OPTION						
125	SG07	Side guards clip (<i>pair</i>) ²		New		NC

1 - Not compatible with 20" and 22" rear wheels 2 - Not available with fixed carbon fiber side guard

POSITIONING BELT						
126	SB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐		E0978
127	SB02	Auto buckle	Standard 35" ☐	Long 60" ☐		E0978
128	SB07	Hardware attachment for 2 point seatbelt (1 pair)				NC
129	SB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)				\$15
130	SB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" ☐	Long 61" ☐		E0978
131	SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	E0978
132	SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	E0978
133	SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" ☐	Long 24" ☐		E0978
134	SB99	Omit positioning belt				STD

HEADREST ¹						
135	AC01	Profiled Headrest	Small (3.5"X6") ☐	Medium (4.5"X7") ☐	Large (5.5"X8") ☐	E0955
136	AC28	Headrest mounting hardware with removable non folding stabilizer bar				E1028

1 - Only available with push handles - Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES					
137	AC03	Anti tippers		E0971	\$140 ☐
138	AC27	Tip assist ¹			\$60 ☐
139	AC04	Transit Tie Down (Unoccupied)			NC ☐
140	AC05	WC-19 Transit Tie Down w/Q'straint belt ²	Short belt 20-41" ☐ Long belt 20-46" ☐		\$320 ☐
141	AC06	Cane carrier	Left ☐ Right ☐	E2207	\$195 ☐
142	AC07	Oxygen tank holder - E Type <i>Not compatible with swing away armrest or amputee axle plate</i>	Left	E2208	\$310 ☐
143	AC08	Spoke guards ³	Available: 20", 22", 24", 26"	K0065	\$115 ☐
144	AC09	Footrest neoprene impact guard (pair) - Black		K0108	\$90 ☐
145	AC10	Seat pouch			\$75 ☐
146	AC26	BodyPoint Mobility Bag			\$125 ☐
147	AC23	Newton Backpack			\$120 ☐
148	AC12	Velcro Closing Strap			\$50 ☐
149	AC13	Tool kit			\$45 ☐
150	AC16	Motion Composites Polo Shirt	W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐		\$32 ☐
151	AC15	Military Patches	Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐		\$90 ☐
152	AC15	Combat Wounded Purple Heart Patch			\$90 ☐


1 - Cannot be combined with Anti Tippers 2 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

3 - Not available with Mag wheels and Spinerger

FRAME COLOR					
153	FC01	Black/Carbon Fiber Black ☐ (Black color is replaced by the Carbon fiber finish with the C2 upgrade)			NC ☐
154	FC03	Sunkissed Orange ☐ Burgundy ☐ Sapphire Blue ☐ Monster Green ☐ Charcoal ☐ Fuchsia ☐ Steel Blue ☐ Ferrari Red ☐ Acid Green ☐			NC ☐
155	FC99	Custom Color Specify PPG Paintit color:		Additional 10 days lead-time	\$360 ☐

TOTAL					
-------	--	--	--	--	--

Fork and Stem bolt combinations



Notes :



BioFit™ CUSHION

BioFit™	BioFit™	Fits Width	Depth	\$	✓
 HCPCS CODE: E2607-E2608	NXCABF1616	15"-16"	16"	425	
	NXCABF1618	15"-16"	18"	425	
	NXCABF1620	15"-16"	20"	425	
	NXCABF1816	17"-18"	16"	425	
	NXCABF1818	17"-18"	18"	425	
	NXCABF1820	17"-18"	20"	425	

BioFit™	Fits Width	Depth	\$	✓
NXCABF2016	19" - 20"	16"	425	
NXCABF2018	19" - 20"	18"	425	
NXCABF2020	19" - 20"	20"	425	
NXCABF2216	21" - 22"	16"	455	
NXCABF2218	21" - 22"	18"	455	
NXCABF2220	21" - 22"	20"	455	

Optima™ THORACIC BACK SUPPORT

Optima™	Optima™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NOACB1613	16"-18"	13"	660	
	NOACB1616	16"-18"	16"	660	
	NOACB1619	16"-18"	19"	660	
	NOACB1813	18"-20"	13"	660	
	NOACB1816	18"-20"	16"	660	
	NOACB1819	18"-20"	19"	660	
	NOACB2016	20"-22"	16"	660	
	NOACB2019	20"-22"	19"	660	
	NOACB2213	22"-24"	13"	710	
	NOACB2216	22"-24"	16"	710	
	NOACB2219	22"-24"	19"	710	

Optima™ DEEP THORACIC BACK SUPPORT

Optima™ DEEP	Optima™ DEEP	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NOPCB1616	16"-18"	16"	660	
	NOPCB1619	16"-18"	19"	660	
	NOPCB1816	18"-20"	16"	660	
	NOPCB1819	18"-20"	19"	660	
	NOPCB2016	20"-22"	16"	660	
	NOPCB2019	20"-22"	19"	660	
	NOPCB2213	22"-24"	13"	710	
	NOPCB2216	22"-24"	16"	710	
	NOPCB2219	22"-24"	19"	710	

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™	Xtend™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NHAACB1616/20	16"-18"	16"-20"	660	
	NHAACB1816/20	18"-20"	16"-20"	660	
	NHAACB2016/20	20"-22"	16"-20"	660	
	NHAACB2216/20	22"-24"	16"-20"	710	

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™ DEEP	Xtend™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NHAACB1616/20	16"-18"	16"-20"	660	
	NHAACB1816/20	18"-20"	16"-20"	660	
	NHAACB2016/20	20"-22"	16"-20"	660	
	NHAPCB2216/20	22"-24"	16"-20"	710	

Classic™ PELVIC BACK SUPPORT

Classic™	Classic™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NACB1616	16"-18"	16"	660	
	NACB1618	16"-18"	18"	660	
	NACB1816	18"-20"	16"	660	
	NACB1818	18"-20"	18"	660	
	NACB1820	18"-20"	20"	660	
	NACB2020	20"-22"	20"	660	
	NACB2216	22"-24"	16"	710	
	NACB2218	22"-24"	18"	710	
	NACB2220	22"-24"	20"	710	

Classic™ DEEP PELVIC BACK SUPPORT

Classic™ DEEP	Classic DEEP™	Fits Width	Height	\$	✓
 HCPCS CODE: E2620-E2621	NPCB1616	16"-18"	16"	660	
	NPCB1618	16"-18"	18"	660	
	NPCB1818	18"-20"	18"	660	
	NPCB1820	18"-20"	20"	660	
	NPCB2020	20"-22"	20"	660	
	NPCB2216	22"-24"	16"	710	
	NPCB2218	22"-24"	18"	710	
	NPCB2220	22"-24"	20"	710	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

 HCPCS CODE: E0955		
Model #NIB-0784	\$	✓
10" x 5" Large Rectangular	195	

 HCPCS CODE: E0955		
Model # NIB-0818	\$	✓
14" x 5" Extra-Large Rectangular	195	

MultiFit™ HEAD SUPPORT – FIXTURES

		
Model #NIB-0747	\$	✓
Multi-Axis Hardware	145	

		
Model # NIB-0823	\$	✓
Head Support Receiver Bracket	55	

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com