



1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

Marked for: _____

2- ATP/Clinician information

ATP/Clinician:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:* _____

Name: _____

Address: _____

City: _____ State: _____

Zip Code: _____

Phone: _____

* Required							English Manual (STD)	Spanish Manual	Taxes not included		
MODEL							USD				
WC03	☐	Ultralight folding Helio wheelchair (265lb weight capacity) A7 package: side frames, symmetrical hydroformed crossbrace and back canes made of 7000 series						K0005	\$2 295		
WC08	☐	C2 carbon fiber upgrade: side frames, symmetrical molded crossbrace and back canes made of ultralight carbon fiber							\$1 835		
WC09	☐	HD Heavy Duty package (350 lb weight capacity)							\$525		
WC10	☐	E-Motion adaptor Kit (2 adaptors) - <i>Only necessary with 6° camber</i>							\$35		
SEAT WIDTH ³											
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐			NC	
	☐	20" ¹ ☐		21" ² ☐	22" ² ☐				E2201	\$495	
SEAT DEPTH											
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐			NC	
	☐	20" ☐							E2203	\$475	

1 - Available with HD Kit (350 lb)

2 - Only with HD Kit (350 lb)

[3. Width calculator](#)

NEWTON ULTRALIGHT PACKAGE										
PK01	☐	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels					Save \$225		\$490	

Helio FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
4"	P1	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2				14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						
	P4	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
7"	P1										16 1/2"	17 1/2"	18 1/2"	17"	18"	19"			
	P2							16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"			
	P3							17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
	P4	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"
	P5	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"	19 1/2"	20 1/2"	21 1/2"

Height available + 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS						
CA01-03	☐	Composite wheel w/ PU tire	3" x 1"			NC
CA01-04	☐		4" x 1"			NC
CA01-05	☐		5" x 1"			NC
CA01-06	☐		6" x 1"			NC
CA01-07	☐		7" x 1"			NC
CA01-08	☐		8" x 1"			NC
CA10	☐	Pneumatic casters	6" x 1 ¼"		E2214/E2215	\$160
CA11	☐		8" x 1 ¼"		E2214/E2215	\$160
CA14	☐	Newton UltraCaster Composite	4" x 1 ½ "			\$85
CA15	☐		5" x 1 ½ "			\$85
CA16	☐		6" x 1 ½ "			\$85
CA17	☐	Newton UltraCaster Aluminum	4" x 1 ½"			\$185
CA18	☐		5" x 1 ½"			\$195
CA19	☐		6" x 1 ½"			\$205
CASTER OPTIONS						
CA12	☐	Caster Pin Locks <i>(Available with +1" or +2" Stem bolt ONLY)</i> ¹			K0073	\$295
CA13	☐	Carbon fiber Frog Legs suspension forks <i>(Not available with caster pin locks)</i> ²			E1015	\$475

1- Not available with WC-19; 2- Not available with 7" and 8" casters; limited height of 18.5" with 6" casters

Helio - REAR Seat-to-floor Height

WHEEL SIZES	Height available	
20" (451) wheel ☐	12 1/2 to 17 1/4	Specify desired height: _____
22" (501) wheel ☐	13 1/2 to 18 1/4	
24" (540) wheel ☐	14 1/2 to 19 1/4	
25" (559) wheel ☐	15 to 19 3/4	
26" (590) wheel ☐	15 1/2 to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest					
SWING-IN/SWING-OUT FOOTREST					
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"	Set at: *		NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at: *		NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at: *		NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"	Set at: *		NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at: *		NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at: *		NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16"	Set at: *		NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at: *		NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at: *		NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5" - 11"	Set at: *	NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5"- 13"	Set at: *	NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5" - 14.5"	Set at: *	NC
FOOTREST OPTIONS					
FR09	☐	High mount <i>(footplate attached to the hanger with clamps)</i> ³	5" - 13.5"	Set at:*	\$80
FR10	☐	Elevating legrest with calf pad, Length adjustable ⁴	Left ☐ Right ☐ 14" - 22"	Set at:	E0990 \$240
FR99	☐	Omit Footrest			NC

1- Available only with two pieces adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount; 4- The elevating legrest replace the footrest side selected; * +/-1/2"

FOOTPLATES					
FP01	☐	Standard - Composite ¹			STD
FP02	☐	Newton adjustable angle - Composite		K0040	\$210
FP03	☐	Newton adjustable angle - Aluminum		K0040	\$230
FP04	☐	Newton adjustable angle - Carbon fiber <i>(with rubber inserts)</i>			\$435
FP06	☐	One piece adjustable angle flip-up - Aluminum		K0040	\$190
FP07	☐	One piece adjustable angle flip-up - Carbon fiber <i>(with rubber inserts)</i>			\$455
FOOTPLATE ACCESSORIES					
FP08	☐	Length adjustable heel loop <i>(not available with one piece footplate)</i>		E0951	\$40
FP09	☐	Calf strap		K0038	\$35
FP30	☐	Bodypoint Aeromesh Padded calf strap		K0038	\$115

1- Not available with front mount

SWING AWAY RESIDUAL LIMB SUPPORT							
FR27	☐	Residual limb support - flat 6"x8"	Left ☐	Right ☐		E1020	\$405
FR28	☐	Residual limb support - flat 8"x10"	Left ☐	Right ☐		E1020	\$405
FR29	☐	Residual limb support - flat 10"x10"	Left ☐	Right ☐		E1020	\$405
FR30	☐	Residual limb support - flat 10"x14"	Left ☐	Right ☐		E1020	\$405
FR32	☐	Residual limb support - contour 6"x8"	Left ☐	Right ☐		E1020	\$405
FR33	☐	Residual limb support - contour 8"x10"	Left ☐	Right ☐		E1020	\$405

Backrest							
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE <i>(for aluminum chair)</i>							

BC60	☐	Straight aluminum back cane	21" - 24"	increment of 1"	Set at ² :	New	\$75
BC59	☐	8° bend aluminum back cane	21" - 24"	increment of 1"	Set at ² :	New	\$75

FIXED ANGLE AND HEIGHT ADJUSTABLE BACK WITH PUSH HANDLE OPTION									
BC11	☐	Straight back cane*	12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐	increment of 1/2"	Set at ² :		NC
BC06									
BC12	☐	8° Bend back cane*	13 ½ - 15" ☐	15 - 18" ☐	18 - 21" STD ☐	increment of 1/2"	Set at ² :		\$55
BC08									

*** fixed height Carbon fiber back canes are unavailable until the end of 2022. If this option is selected, Aluminum back canes will be installed. Otherwise, carbon fiber angle and height adjustable are still available. See options below**

ANGLE AND HEIGHT ADJUSTABLE BACK WITH PUSH HANDLE OPTION								
BC13	☐	Straight back cane	15 - 18" ☐	18 - 21" STD ☐	increment of 1/2"	Set at ² :		\$210
BC09								
BC14	☐	8° bend back cane	15 - 18" ☐	18 - 21" STD ☐	increment of 1/2"	Set at ² :		\$265
BC10								
	Select back cane angle:		85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐

PUSH HANDLE OPTIONS							
PH01	☐	Push handles					STD
PH05	☐	Fold down push handles ^{3, 5}					\$126
PH06-920	☐	Clamp on height adjustable push handles 9" - 20" ^{3, 4, 5}					\$315
PH06-923	☐	Clamp on height adjustable push handles 9" - 23" ^{3, 4, 5}					\$315
PH02	☐	Omit push handles					NC

STROLLER HANDLE AND STABILIZER BAR							
PH04	☐	Stroller handle <i>(only available with push handle) - (not available with the HD option)</i>				K0108	\$160
PH03	☐	Newton folding stabilizer bar <i>(not compatible with stroller handle)</i>					\$265

1- Not available with a 8° bend back cane; **2-** If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; **3-** Not available with HD option
4- Not compatible with swing away armrests, 8° bend back canes and amputee axle plates; **5-** Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY							
BH07	☐	Slip-on Active-Mesh back upholstery					STD
BH08	☐	Adjustable tension Active-Mesh back upholstery				E2611	\$300
BH99	☐	Omit Back upholstery					NC
	☐	NXT back Support				See pricelist	

SEAT CUSHION							
SC01	☐	2" Cushion				E2601	\$165
SC02	☐	3" Cushion				E2601	\$185
SC99	☐	Omit seat cushion					STD
	☐	NXT cushion				See pricelist	

Rear wheels						
ONE ARM DRIVE ¹		Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims				
HR22	☐	One arm drive ^{2, 3}	24"	left ☐	right ☐	\$1 465
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)				NC
REAR WHEELS			Available			
RW01	☐	Mag ⁴	20", 22" 24"			\$55
RW02	☐	Newton One - spoke wheel ⁵	20", 22" 24"			NC
RW03	☐	Newton Gravity - ultralight wheel	20", 22", 24", 25", 26"			\$315
RW04	☐	Spinergy Spox Blk STD ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐	22", 24", 25", 26"			\$835
RW05	☐	Spinergy LX Blk STD ☐ Pur ⁶ ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐ Pink ⁶ ☐ Grn ⁶ ☐ Org ⁶ ☐	22", 24", 25", 26"			\$1 045
RW07	☐	Spinergy CLX ⁵ Blk STD ☐ Pur ⁶ ☐ Wht ⁶ ☐ Blue ⁶ ☐ Red ⁶ ☐ Yel ⁶ ☐ Pink ⁶ ☐ Grn ⁶ ☐ Org ⁶ ☐	24", 25"			\$1 555
TIRES			Available			
TI01	☐	Soft Urethane - 1 3/8"	20", 22" 24"	E2218	STD	
TI02	☐	Soft Urethane - knobby 1 3/8 "	22", 24", 25", 26"	E2218	\$95	
TI04	☐	Soft Urethane - Shox 1"	22", 24", 25", 26"	E2218	\$115	
TI05	☐	Pneumatic - 1 3/8"	20", 22", 24", 26"	E2211/E2212	NC	
TI08	☐	Pneumatic with Airless insert - 1 3/8" ⁷	22", 24", 26"	E2211/E2213	\$80	
TI06	☐	Pneumatic HP - 1"	20", 22", 24", 25", 26"	E2211/E2212	\$95	
TI07	☐	Pneumatic HP - Puncture resistant, black 1"	24", 25", 26"	E2211/E2212	\$125	

1- Offered with handrims on each side of the wheelchair
2- Only standard handrim position available
3 - Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position
5 - Not compatible with Shox tires

6 - Color spokes require an additionnal 10 days lead-time
7- Only available with Mag, Newton One and Newton Gravity rear wheels

Handrims						
HANDRIM			Available			
HR01	☐	Aluminum anodized	20", 22", 24", 25", 26"		STD	
HR61	☐	Aluminum black hard anodized ²	20", 22", 24", 25"		\$70	
HR02	☐	Plastic coated ¹	20", 22", 24", 25", 26"		\$100	
HR03	☐	Q-Grip High friction Neoprene coated ^{1,2}	22", 24", 25", 26"		\$370	
HR07	☐	Newton Air Grip ² Black ☐	20", 22", 24", 25", 26"		\$100	
HR	☐	Newton Air Grip ² (08) Sunkissed Orange ☐ (09) Burgundy ☐ (10) Sapphire Blue ☐ (11) Monster Green (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐ (15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐	20", 22", 24", 25", 26"		\$250	
HR28	☐	Surge ^{1,2} (1-1/2" Oval handrim w/ Gripton high friction strip)	**Limited availability	22", 24", 25", 26"		\$395
HR05	☐	Surge LT ^{1,2} (1-1/8" Oval handrim w/ Gripton high friction strip)	**Limited availability	22", 24", 25", 26"		\$425
HR04	☐	Natural Fit ^{1,2}	**Limited availability	22", 24", 25", 26"		\$325
HR99	☐	Omit handrim				NC

1- Only standard handrim position available
2- Not available with Mag wheels or OAD

Handrims option						
Assembly position for aluminum and Newton AirGrip handrims						
HR29	☐	Standard Handrim position Full tab				STD
HR30	☐	Narrow Handrim position Short tab 2 positions				\$27
HR31	☐	Super Narrow Handrim Short tab 2 positions STD ☐ Short tab 1 position ☐				\$27
Thumb piece ³ selection for Natural Fit handrims						
HR23	☐	Standard grip	**Limited availability			\$90
HR27	☐	Super grip	**Limited availability			\$145
HR099	☐	Omit thumb piece				NC


3- Compatible with Surge, Surge LT handrims

AXLE ¹						
AX01	☐	Quick release axle				STD
AX02	☐	Permanent axle				NC
AX03	☐	Quad release axle		K0108		\$125
AXLE PLATE						
AX04	☐	Amputee adjustable axle plate ²		E0959		\$195

1- Not required with one arm drive
2 - May interfere with swing away armrest and limit seat to floor height available

REAR WHEEL CAMBER						
CM00	☐	0° camber				STD
CM03	☐	3° camber				NC
CM06	☐	6° camber (not available with amputee axle plate)				\$160

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable												Most responsive			
1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"			
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐			

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2,5".

Possible conflict between a 26" rear wheel and a center of gravity over 2,5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

Only on carbon frame

WHEEL LOCK						
WL08	☐	Newton Wheel lock - push to lock				STD
WL02	☐	Pull to lock				NC
WL04	☐	Newton grade aid push to lock <i>(suggested with pneumatic or high pressure tire)</i>			E0974	\$150
WL13 WL14	☐	6" removable extension handle			E0961	\$75
WL05	☐	Scissors Lock				\$150
WL06	☐	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>				\$210
WL15	☐	Unilateral wheel lock - pull to lock ¹		Left ☐ Right ☐	New	\$325
WL99	☐	Omit wheel locks				NC

1- Only available with: seat depths of ≥ 16"; rear wheels of 20", 22" and 24"; Not available with the following combinations: T or U armrests with One Arm Drive or >0° camber;

May interfere with some COG and other configurations. For more details, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_us/support-and-education/how-to-documents

<http://www.motioncomposites.com/en-us/support-and-education/how-to-documents>

Armrests								
HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) (improved) ^{1,4}								
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		E0973	\$350
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		E0973	\$350
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		E0973	\$350
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		E0973	\$350
HEIGHT ADJUSTABLE "T" ARMREST (improved) ¹								
AR13-610	☐	10" Desk length armpad ⁵	6" - 10"			New	E0973	\$325
AR03-812	☐	10" Desk length armpad	8" - 12"				E0973	\$325
AR03-1014	☐	10" Desk length armpad	10" - 14"				E0973	\$325
AR14-610	☐	14" Full length armpad ⁵	6" - 10"			New	E0973	\$325
AR04-812	☐	14" Full length armpad	8" - 12"				E0973	\$325
AR04-1014	☐	14" Full length armpad	10" - 14"				E0973	\$325
HEIGHT ADJUSTABLE SWING AWAY ARMREST ^{2,4}								
AR05	☐	Tubular armrest without plastic side guards		8" - 11"			E0973	\$225
AR19	☐	Tubular armrest with plastic side guards		8" - 11"			E0973	\$385
AR06	☐	10" Desk length flat armpad without plastic side guards		8" - 11"			E0973	\$225
AR20	☐	10" Desk length flat armpad with plastic side guards		8" - 11"			E0973	\$385
AR17	☐	14" Full length flat armpad without plastic side guards		8" - 11"		New	E0973	\$225
AR21	☐	14" Full length flat armpad with plastic side guards		8" - 11"			E0973	\$385
ARMREST OPTIONS								
AR12	☐	14" Full Length armpad ³ - Gel Ovation						\$160
AR99	☐	Omit armrest						NC

1- May have 1/2" adjustment limit on both limits of the range; 2- Could interfere with some rigid backrests; 3- Not available with swing away and flip back armrest; 4- Not available with attendant lock

5- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS (not available with Flip Back and "T" Armrest)						
SG01	☐	Plastic side guards (removable)			K0108	\$195
SG02	☐	Carbon fiber side guards (fixed)				\$265
SG03	☐	Carbon fiber side guards (removable)			K0108	\$290
SG06	☐	Plastic fender side guards (removable)				\$285
SG04	☐	Carbon fiber fender side guards (removable)			K0108	\$540
SG99	☐	Omit side guard				STD
SIDE GUARDS OPTION						
SG07	☐	Side guards clip (pair) ¹				NC

1- Not available with fixed carbon fiber side guard

POSITIONING BELT						
SB01-S	☐	Velcro adjustable - Standard 35"		E0978	\$45	
SB01-L	☐	Velcro adjustable - Long 60"		E0978	\$45	
SB02-S	☐	Auto buckle - Standard 35"	<i>Improved</i>	E0978	\$85	
SB02-L	☐	Auto buckle - Long 60"	<i>Improved</i>	E0978	\$85	
SB07	☐	Hardware attachment for 2 point seatbelt (1 pair)			NC	
SB03	☐	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)			\$16	
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"		E0978	\$80	
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"		E0978	\$80	
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"		E0978	\$170	
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"		E0978	\$170	
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"		E0978	\$170	
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"		E0978	\$190	
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"		E0978	\$190	
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"		E0978	\$190	
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"		K0108	\$300	
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"		K0108	\$300	
SB99	☐	Omit positioning belt			STD	

HEADREST ¹						
AC01-S	☐	Profiled Headrest - Small (3.5"X6")		E0955	\$180	
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")		E0955	\$180	
AC01-L	☐	Profiled Headrest - Large (5.5"X8")		E0955	\$180	
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar		E1028	\$295	

1 - Only available with push handles - Not available with a stroller handle or a folding stabilizer bar


ACCESSORIES						
AC03	☐	Anti-tippers		E0971	\$140	
AC29	☐	Swing-away anti-tipper ¹	Left ☐ Right ☐		\$160	
AC27	☐	Tip assist ²			\$65	
AC04	☐	Transit Tie Down (Unoccupied)			NC	
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ³	Short belt 20"-41"		\$355	
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ³	Long belt 20"-46"		\$355	
AC35	☐	WC-19 Transit Tie Down without belt ³		<i>New</i>	\$325	
AC06	☐	Cane and crutch holder	Left ☐ Right ☐	E2207	\$195	
AC07	☐	Oxygen tank holder - D and E Type <i>Not compatible with swing away armrest or amputee axle plate</i>	Left	E2208	\$310	
AC08	☐	Spoke guards ⁴	Available: 20", 22", 24", 26"	K0065	\$115	
AC09	☐	Footrest neoprene impact guard (pair) - Black		K0108	\$90	
AC10	☐	Seat pouch			\$80	
AC30	☐	Newton multi-purpose pouch			\$60	
AC26	☐	BodyPoint Mobility Bag			\$130	
AC23	☐	Newton Backpack			\$125	
AC12	☐	Velcro Closing Strap			\$55	
AC13	☐	Tool kit			\$50	
AC16	☐	Motion Composites Polo Shirt	W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐		\$35	
AC15	☐	Military Patches	Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐ Combat Wounded Purple Heart ☐		\$95	

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4" 2 - Cannot be combined with Anti Tippers

3 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option 4 - Not available with Mag wheels and Spinergy

FRAME COLOR						
FC01	☐	Black/Carbon Fiber Black (Black color is replaced by the Carbon fiber finish with the C2 upgrade)				NC
FC	☐	(03) Sunkissed orange ☐ (04) Burgundy ☐ (05) Sapphire Blue ☐ (06) Monster Green ☐ (07) Charcoal ☐ (08) Fuchsia ☐ (09) Steel Blue ☐ (17) Ferrari Red ☐ (18) Acid Green ☐				NC
FC99	☐	Custom Color	Specify PPG Paintit color:	<i>Additional 10 days lead-time</i>		\$380


Fork and Stem bolt combinations



Unilateral wheel lock restrictions

Seat depth must be 16in and greater		
Rear wheel diameter 24 in and less		
"T" or flip-back armrest only available with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Wheel lock lever extension not available		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
22 in	16	1.0 to 2.75
	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :

BioFit™ CUSHION		WIDTH											
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
 NXCAF _____ (W" x D") HCPCS CODE: E2607-E2608	DEPTH	12"	☐ \$425	☐ \$425	☐ \$425								
		14"	☐ \$425	☐ \$425	☐ \$425								
		15"		☐ \$425	☐ \$425	☐ \$425	☐ \$425						
		16"		☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	☐ \$455	☐ \$455	
		17"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455		
		18"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	☐ \$455	☐ \$455
		19"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455		
		20"				☐ \$425	☐ \$425	☐ \$425	☐ \$425	☐ \$455	☐ \$455	☐ \$455	☐ \$455
Optional Incontinent / Infection Control Outer Cover NXCAF _____ V ☐ N/C													

Optima™ THORACIC BACK SUPPORT

THORACIC SHAPE	Optima™	Fits Width	Height	✓	\$
	NOACB1613	16"-18"	13"	☐	\$660
	NOACB1616	16"-18"	16"	☐	\$660
	NOACB1619	16"-18"	19"	☐	\$660
	NOACB1813	18"-20"	13"	☐	\$660
	NOACB1816	18"-20"	16"	☐	\$660
	NOACB1819	18"-20"	19"	☐	\$660
	NOACB2016	20"-22"	16"	☐	\$660
	NOACB2019	20"-22"	19"	☐	\$660
	NOACB2213	22"-24"	13"	☐	\$710
	NOACB2216	22"-24"	16"	☐	\$710
	NOACB2219	22"-24"	19"	☐	\$710

HCPCS CODE: E2620-E2621

Optima™ DEEP THORACIC BACK SUPPORT

Optima™ DEEP	Fits Width	Height	✓	\$
NOPCB1616	16"-18"	16"	☐	\$660
NOPCB1619	16"-18"	19"	☐	\$660
NOPCB1816	18"-20"	16"	☐	\$660
NOPCB1819	18"-20"	19"	☐	\$660
NOPCB2016	20"-22"	16"	☐	\$660
NOPCB2019	20"-22"	19"	☐	\$660
NOPCB2213	22"-24"	13"	☐	\$710
NOPCB2216	22"-24"	16"	☐	\$710
NOPCB2219	22"-24"	19"	☐	\$710

HCPCS CODE: E2620-E2621

Xtend™ HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™	Fits Width	Height	✓	\$
NHAACB1616/20	16"-18"	16"-20"	☐	\$660
NHAACB1816/20	18"-20"	16"-20"	☐	\$660
NHAACB2016/20	20"-22"	16"-20"	☐	\$660
NHAACB2216/20	22"-24"	16"-20"	☐	\$710

HCPCS CODE: E2620-E2621

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™ DEEP	Fits Width	Height	✓	\$
NHAPCB1616/20	16"-18"	16"-20"	☐	\$660
NHAPCB1816/20	18"-20"	16"-20"	☐	\$660
NHAPCB2016/20	20"-22"	16"-20"	☐	\$660
NHAPCB2216/20	22"-24"	16"-20"	☐	\$710

HCPCS CODE: E2620-E2621

Classic™ PELVIC BACK SUPPORT

Classic™	Fits Width	Height	✓	\$
NACB1616	16"-18"	16"	☐	\$660
NACB1618	16"-18"	18"	☐	\$660
NACB1816	18"-20"	16"	☐	\$660
NACB1818	18"-20"	18"	☐	\$660
NACB1820	18"-20"	20"	☐	\$660
NACB2020	20"-22"	20"	☐	\$660
NACB2216	22"-24"	16"	☐	\$710
NACB2218	22"-24"	18"	☐	\$710
NACB2220	22"-24"	20"	☐	\$710

HCPCS CODE: E2620-E2621

Classic™ DEEP PELVIC BACK SUPPORT

Classic DEEP™	Fits Width	Height	✓	\$
NPCB1616	16"-18"	16"	☐	\$660
NPCB1618	16"-18"	18"	☐	\$660
NPCB1818	18"-20"	18"	☐	\$660
NPCB1820	18"-20"	20"	☐	\$660
NPCB2020	20"-22"	20"	☐	\$660
NPCB2216	22"-24"	16"	☐	\$710
NPCB2218	22"-24"	18"	☐	\$710
NPCB2220	22"-24"	20"	☐	\$710

HCPCS CODE: E2620-E2621

MultiFit™ HEAD SUPPORT - MULTI-ADJUSTABLE PADS

Model #	✓	\$
10" x 5" Large Rectangular	☐	\$195
14" x 5" Extra-Large Rectangular	☐	\$195

HCPCS CODE: E0955

MultiFit™ HEAD SUPPORT - FIXTURES

Model #	✓	\$
Multi-Axis Hardware	☐	\$145
Head Support Receiver Bracket	☐	\$55

HCPCS CODE: E0955

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com