



1- Purchasing information

P.O. number:* Date:

Company Name:*

Purchasing contact:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

Marked for:

2- ATP/Clinician information

ATP/Clinician:*

Phone:*

Email:*

Contact preference:* Phone ☐ Email ☐

3- Shipping information

Ship to:*

Name:

Address:

City: State:

Zip Code:

Phone:

* Required		English Manual (STD)				Spanish Manual		Taxes not included	
MODEL						USD			
WC03	☐	Ultralight folding Helio wheelchair (265lb weight capacity) A7 package: side frames, symmetrical hydroformed crossbrace and back canes made of 7000 series aluminum					K0005	\$2,395	
WC08	☐	C2 carbon fiber upgrade: side frames, symmetrical molded crossbrace and back canes made of ultralight carbon fiber						\$2,245	
WC09	☐	HD Heavy Duty package (350 lb weight capacity) ⁵						\$535	
WC10	☐	E-Motion adaptor Kit (2 adaptors) - <i>Only necessary with 6° camber</i>						\$36	
Frame									
SEAT WIDTH ³									
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐		NC
	☐		20" ¹ ☐	21" ² ☐	22" ² ☐			E2201	\$505
SEAT DEPTH									
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐		NC
	☐		20" ☐					E2203	\$485
TRANSIT TIE DOWN									
AC04	☐	Transit Tie Down (Unoccupied)					E1022	NC	
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 20"-41"					E1023	\$365	
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 20"-46"					E1023	\$365	
AC35	☐	WC-19 Transit Tie Down without belt ⁴					E1023	\$335	

1- Available with HD Kit (350 lb); 2- Only with HD Kit (350 lb)

[3. Width calculator](#)

4- May interfere with swing away armrest and Clamp on height adjustable push handle; not available with HD option when carbon upgrade is selected.

5- HD fee is required for clients weighing 300-350 lbs



NEWTON ULTRALIGHT PACKAGE			
PK01	☐	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels	Save \$275 \$530
VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹	K0108 \$95

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo *VP01; Recommended with black chair; See picture of the product on page 6.

Helio FRONT Seat-to-floor Height - Please select your choice

Forks / Position		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"												
	P2	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
5"	P1																		
	P2				14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						
	P4	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
7"	P1										16 1/2"	17 1/2"	18 1/2"	17"	18"	19"			
	P2							16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"			
	P3							17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
	P4	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"
	P5	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"	19"	20"	21"	19 1/2"	20 1/2"	21 1/2"
Froglegs	P1				14 1/2"	15 1/2"	16 1/2"												
	P2	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"									
	P3	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"						

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS				Available		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 1/4")		6", 8"	E2214/E2215	\$165
CA14 to CA16	☐	Newton UltraCaster Composite (1 1/4")		4", 5", 6"		\$95
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 1/2")		4", 5", 6"		\$215
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹			K0073	\$300
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)			E1015	\$490

1- Not available with WC-19

Helio - REAR Seat-to-floor Height

WHEEL SIZES		Height available		Specify desired height: _____
☐	20" (451) wheel	12 1/2"	to 17 1/4"	
☐	22" (501) wheel	13 1/2"	to 18 1/4"	
☐	24" (540) wheel	14 1/2"	to 19 1/4"	
☐	25" (559) wheel	15"	to 19 3/4"	
☐	26" (590) wheel	15 1/2"	to 20 1/4"	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/4" of selected seat to floor height"

Footrest				
SWING-IN/SWING-OUT FOOTREST				
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"	Set at*:	NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at*:	NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at*:	NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"	Set at*:	NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at*:	NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at*:	NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" ⁵	Set at*:	NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" ⁵	Set at*:	NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" ⁵	Set at*:	NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹	7.5" - 11"	Set at*:
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2}	9.5" - 13"	Set at*:
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2}	9.5" - 14.5"	Set at*:
FOOTREST OPTIONS				
FR09	☐	High mount (footplate attached to the hanger with clamps) ³	5" - 13.5"	Set at*:
FR10	☐	Elevating legrest with calf pad, Length adjustable ⁴ Left ☐ Right ☐	14" - 22"	Set at:
FR99	☐	Omit Footrest		NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest; * +/- 1/2"

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES						
FP01	☐	Standard - Composite ¹				STD
FP02	☐	Newton adjustable angle - Composite ^{*VP01}			K0040	\$215
FP03	☐	Newton adjustable angle - Aluminum			K0040	\$235
FP04	☐	Newton adjustable angle - Carbon fiber <i>(with rubber inserts)</i>				\$445
FP23	☐	Newton adjustable angle extra deep - Aluminum			K0040	\$235
FP06	☐	One piece adjustable angle flip-up - Aluminum			K0040	\$195
FP07	☐	One piece adjustable angle flip-up - Carbon fiber <i>(with rubber inserts)</i>				\$465
FOOTPLATE ACCESSORIES						
FP08	☐	Length adjustable heel loop ²			E0951	\$41
FP09	☐	Calf strap			K0038	\$36
FP30	☐	Bodypoint Aeromesh Padded calf strap			K0038	\$118

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	☐	Residual limb support - flat 6"x8"	Left ☐ Right ☐		E1020	\$420
FR28	☐	Residual limb support - flat 8"x10"	Left ☐ Right ☐		E1020	\$420
FR29	☐	Residual limb support - flat 10"x10"	Left ☐ Right ☐		E1020	\$420
FR30	☐	Residual limb support - flat 10"x14"	Left ☐ Right ☐		E1020	\$420
FR32	☐	Residual limb support - contour 6"x8"	Left ☐ Right ☐		E1020	\$420
FR33	☐	Residual limb support - contour 8"x10"	Left ☐ Right ☐		E1020	\$420

Backrest								
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLES - For Aluminum chair <i>(Not available with push handles options)</i>								
BC58	☐	Straight aluminum back cane	9 - 15" ☐	15 - 21" ☐		NC		
BC60	☐	Straight aluminum back cane	21" - 24" <i>increment of 1"</i>		Set at ² :	\$77		
BC01	☐	8° bend aluminum back cane	16" ☐	18" ☐	20" ☐	NC		
BC59	☐	8° bend aluminum back cane	21" - 24" <i>increment of 1"</i>		Set at ² :	\$77		
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>								
BC11	☐	Straight back cane	12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐	<i>increment of 1/2"</i>	Set at ² :	NC
BC06								
BC12	☐	8° Bend back cane	13 ½ - 15" ☐	15 - 18" ☐	18 - 21" STD ☐	<i>increment of 1/2"</i>	Set at ² :	\$60
BC08								
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS <i>(MUST select ONE push handle options)</i>								
BC13	☐	Straight back cane	15 - 18" ☐	18 - 21" STD ☐	<i>increment of 1/2"</i>		Set at ² :	\$215
BC09								
BC14	☐	8° bend back cane	15 - 18" ☐	18 - 21" STD ☐	<i>increment of 1/2"</i>		Set at ² :	\$270
BC10								
		Select back cane angle:	85° ¹ ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐
PUSH HANDLES OPTIONS <i>(ONE option at a time)</i>								
PH01	☐	Push handles						STD
PH05	☐	Fold down push handles ^{3,5}						\$145
PH06-920	☐	Clamp on height adjustable push handles 9" - 20" ^{3,4,5}						\$325
PH06-923	☐	Clamp on height adjustable push handles 9" - 23" ^{3,4,5}						\$325
PH02	☐	Omit push handles						NC
STROLLER HANDLE AND STABILIZER BAR								
PH04	☐	Stroller handle <i>(only available with push handle) - (not available with the HD option)</i>					K0108	\$165
PH03	☐	Newton folding stabilizer bar <i>(not compatible with stroller handle)</i>						\$270

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265lb weight limit

4- Not compatible with swing away armrests, 8° bend back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

BACK UPHOLSTERY						
BH07	☐	Slip-on Active-Mesh back upholstery				STD
BH08	☐	Adjustable tension Active-Mesh back upholstery			E2611/E2612	\$305
BH99	☐	Omit Back upholstery				NC
	☐	NXT back Support			See pricelist	
SEAT CUSHION						
SC01	☐	2" Cushion				\$170
SC02	☐	3" Cushion			E2601/E2602	\$190
SC99	☐	Omit seat cushion				STD
	☐	NXT cushion			See pricelist	

Rear wheels												
ONE ARM DRIVE ¹		Only available with Newton One rear wheels, permanent axles, 0° camber										
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐			E0958	\$1,495				
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐			E0958	\$1,600				
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)						NC				
REAR WHEELS					Available							
RW01	☐	Mag ³				20", 22" 24"		\$57				
RW02	☐	Newton One - spoke wheel ⁴				20", 22" 24"		NC				
RW03	☐	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"		\$360				
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	\$915				
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐	\$1,145
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	Pink ⁵ ☐	Grn ⁵ ☐	Org ⁵ ☐	\$1,625
TIRES					Available							
AIRLESS												
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black ⁹				20", 22" 24"	E2218	STD				
TI02	☐	Soft Urethane 1 3/8 " - Medium Tread, black ^{7,9}				22", 24", 25", 26"	E2218	\$97				
TI04	☐	Soft Urethane 1" - No tread, black - <i>Shox</i>				22", 24", 25", 26"	E2218	\$118				
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸				22", 24", 26"	E2211/E2213	\$82				
PNEUMATIC												
TI05	☐	Pneumatic 1 3/8" - Medium Tread ⁹				20", 22", 24", 25", 26"	E2211/E2212	NC				
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>				20", 22", 24", 25", 26"	E2211/E2212	\$105				
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>				24", 25", 26"	E2211/E2212	\$150				
TI03	☐	Pneumatic 2" - Large Tread, black ¹⁰				22", 24", 25"	E2211/E2212	\$115				

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available;
3- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additionnal 10 days lead-time;
6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels
9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;

Handrims							
HANDRIM		Available					
HR01	☐	Aluminum anodized			20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²			20", 22", 24", 25"		\$75
HR02	☐	Plastic coated ¹			20", 22", 24", 25", 26"		\$110
HR03	☐	High friction coated ^{1,2,3}			22", 24", 25", 26"		\$375
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ² See on page 6 for color choice and price			20", 22", 24", 25", 26"		
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"		\$450
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)			22", 24", 25", 26"		\$495
HR04	☐	Natural Fit ^{1,2,3,4}			22", 24", 25", 26"		\$395
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)			22", 24", 25", 26"		\$395
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)			22", 24", 25", 26"		\$445
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)			22", 24", 25", 26"		\$345
HR99	☐	Omit handrim					NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;
4- 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

Handrims option							
ASSEMBLY POSITION							
HR29	☐	Standard Handrim position	Full tab				STD
HR030	☐	Narrow Handrim position no cut	Full tab				NC
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions				\$30
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐	Short tab 1 position ☐			\$30
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS							
HR23	☐	Standard grip					\$95
HR27	☐	Super grip					\$150
HR099	☐	Omit thumb piece					NC

5- Only available with Aluminum and Air Grip handrims;

AXLE ¹							
AX01	☐	Quick release axle					STD
AX02	☐	Permanent axle					NC
AX03	☐	Quad release axle			K0108		\$130
AXLE PLATE							
AX04	☐	Amputee adjustable axle plate ²			E0959		\$205

1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

REAR WHEEL CAMBER						
CM00	☐	0° camber				STD
CM03	☐	3° camber				NC
CM06	☐	6° camber <i>(not available with amputee axle plate)</i>				\$165

Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

Most stable													
AMPUTEE	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2,5".

Possible conflict between a 26" rear wheel and a center of gravity over 2,5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

Only on carbon frame

WHEEL LOCK						
WL08	☐	Newton Wheel lock - push to lock		PRODUCT IMPROVED		STD
WL16	☐	Aluminum push to lock with extensions			E0961	\$77
WL02	☐	Pull to lock				NC
WL04	☐	Newton grade aid push to lock <i>(works best with pneumatic or high pressure tire)</i>			E0974	\$155
WL13 WL14	☐	6" removable extension handle ¹ *VP01			E0961	\$77
WL05	☐	Scissors Lock				\$155
WL06	☐	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>				\$215
WL15	☐	Unilateral wheel lock - pull to lock ²	Left ☐	Right ☐		\$330
WL99	☐	Omit wheel locks				NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:
https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests						
HEIGHT ADJUSTABLE FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i> ^{1,4} *VP01						
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐	E0973 \$355
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐	E0973 \$355
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐	E0973 \$355
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐	E0973 \$355
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01						
AR13-610	☐	10" Desk length armpad ⁵	6" - 10"			E0973 \$330
AR03-812	☐	10" Desk length armpad	8" - 12"			E0973 \$330
AR03-1014	☐	10" Desk length armpad	10" - 14"			E0973 \$330
AR14-610	☐	14" Full length armpad ⁵	6" - 10"			E0973 \$330
AR04-812	☐	14" Full length armpad	8" - 12"			E0973 \$330
AR04-1014	☐	14" Full length armpad	10" - 14"			E0973 \$330
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²						
AR05	☐	Tubular armrest without plastic side guards - Short pad	8" - 11"			E0973 \$230
AR19	☐	Tubular armrest with plastic side guards - Short pad	8" - 11"			E0973 \$395
AR23	☐	Tubular armrest without plastic side guards - Long pad	8" - 11"			E0973 \$230
AR24	☐	Tubular armrest with plastic side guards - Long pad	8" - 11"			E0973 \$395
AR06	☐	10" Desk length flat armpad without plastic side guards	8" - 11"			E0973 \$230
AR20	☐	10" Desk length flat armpad with plastic side guards	8" - 11"			E0973 \$395
AR17	☐	14" Full length flat armpad without plastic side guards	8" - 11"			E0973 \$230
AR21	☐	14" Full length flat armpad with plastic side guards	8" - 11"			E0973 \$395
ARMREST OPTIONS						
AR12	☐	14" Full Length armpad ³ - Gel Ovation				\$165
AR99	☐	Omit armrest				NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;

3- Only available with full length armrest; Not available with swing away and desk lenght armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

SIDE GUARDS *(not available with Flip Back and "T" Armrest)*

SG01	☐	Plastic side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹	K0108	\$200
SG02	☐	Carbon fiber side guards <i>(fixed)</i>				\$270
SG03	☐	Carbon fiber side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹	K0108	\$305
SG06	☐	Plastic fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹		\$290
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>	Straight bracket - STD	Offset bracket ¹	K0108	\$550
SG99	☐	Omit side guard				STD

POSITIONING BELT					
SB01-S	☐	Velcro adjustable - Small 35" ¹		E0978	\$46
SB01-L	☐	Velcro adjustable - Long 60"		E0978	\$46
SB02-S	☐	Auto buckle - Small 35" ¹		E0978	\$87
SB02-M	☐	Auto buckle - Medium 48"		E0978	\$87
SB02-L	☐	Auto buckle - Long 60"		E0978	\$87
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"		E0978	\$85
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"		E0978	\$85
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"		E0978	\$180
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"		E0978	\$180
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"		E0978	\$180
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"		E0978	\$195
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"		E0978	\$195
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"		E0978	\$195
SB03	☐	Hardware attachment for 2 extra points (1 pair)			\$17
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"		K0108	\$305
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"		K0108	\$305
SB99	☐	Omit positioning belt			STD

HEADREST ¹					
AC01-S	☐	Profiled Headrest - Small (3.5"X6")		E0955	\$185
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")		E0955	\$185
AC01-L	☐	Profiled Headrest - Large (5.5"X8")		E0955	\$185
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar		E1028	\$300

ACCESSORIES									
AC03	☐	Anti-tippers						E0971	\$145
AC29	☐	Swing-away anti-tipper ¹	Left ☐	Right ☐					\$190
AC27	☐	Tip assist ²							\$67
AC06	☐	Cane and crutch holder	Left ☐	Right ☐				E2207	\$200
AC07	☐	Oxygen tank holder ⁴ - D and E Type - installed on the left						E2208	\$315
AC08	☐	Spoke guards ³	Available: 20", 22", 24", 26"					K0065	\$118
AC09	☐	Footrest neoprene impact guard (pair) - Black						K0108	\$92
AC10	☐	Seat pouch							\$82
AC30	☐	Newton multi-purpose pouch							\$61
AC26	☐	BodyPoint Mobility Bag							\$135
AC23	☐	Newton Backpack							\$130
AC12	☐	Velcro Closing Strap							\$57
AC13	☐	Tool kit							\$51
AC16	☐	Motion Composites Polo Shirt	W/S ☐	W/M ☐	W/L ☐				\$36
			M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	M/XXXL ☐	
AC15	☐	Military Patches	Army ☐	Navy ☐	Air Force ☐	Marines ☐	Coast Guard ☐		\$96
			Combat Wounded Purple Heart ☐	Space Force ☐	American Flag ☐				


2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spineray; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down

Colors																					
Options ¹	Black	High Gloss Carbon Fiber	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
																					
AirGrip	☐ \$110	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	\$265	
Frame ²	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	NC	
FC99	☐	Custom Color	Specify PPG PaintIt color:																	Additional 10 days lead-time	\$390

2- Black color is replaced by the Carbon fiber finish with the C2 upgrade:

Fork and Stem bolt combinations



Vision Package Accents



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

Notes :

BioFit™ CUSHION		WIDTH										
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
DEPTH	12"	\$445	\$445	\$445								
	14"	\$445	\$445	\$445								
	15"		\$445	\$445	\$445	\$445						
	16"		\$445	\$445	\$445	\$445	\$445	\$445	\$480	\$480	\$480	\$480
	17"				\$445	\$445	\$445	\$445	\$480	\$480		
	18"				\$445	\$445	\$445	\$445	\$480	\$480	\$480	\$480
	19"				\$445	\$445	\$445	\$445	\$480	\$480		
	20"				\$445	\$445	\$445	\$480	\$480	\$480	\$480	\$480
Option: Change cover for incontinent/infection control outer cover												N/C

Green boxes are stocked sizes and the red boxes require customization.

SEATING AND CUSHION ACCESSORIES

	MODEL NO.	Description	Fits Back Width	\$	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6 & 43.2 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7 & 48.3 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	80	

THORACIC SHAPE	Optima™ THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb) HCPCS CODE: E2620-E2621						Optima™ DEEP THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb) HCPCS CODE: E2620-E2621					
	WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM	WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM
	14"	13"	NOCFACB1413	\$925	NOACB1413	\$695	14"	13"	NOCFPCB1413	\$1015	NOPCB1413	\$695
	14"	16"	NOCFACB1416	\$960	NOACB1416	\$695	14"	16"	NOCFPCB1416	\$1055	NOPCB1416	\$695
	15"	13"	NOCFACB1513	\$925	NOACB1513	\$695	15"	13"	NOCFPCB1513	\$1015	NOPCB1513	\$695
	15"	16"	NOCFACB1516	\$960	NOACB1516	\$695	15"	16"	NOCFPCB1516	\$1055	NOPCB1516	\$695
	16"	13"	NOCFACB1613	\$925	NOACB1613	\$695	16"	13"	NOCFPCB1613	\$1015	NOPCB1613	\$695
	16"	16"	NOCFACB1616	\$960	NOACB1616	\$695	16"	16"	NOCFPCB1616	\$1055	NOPCB1616	\$695
	16"	19"			NOACB1619	\$695	16"	19"			NOPCB1619	\$695
	17"	13"	NOCFACB1713	\$925	NOACB1713	\$695	17"	13"	NOCFPCB1713	\$1015	NOPCB1713	\$695
	17"	16"	NOCFACB1716	\$960	NOACB1716	\$695	17"	16"	NOCFPCB1716	\$1055	NOPCB1716	\$695
	18"	13"	NOCFACB1813	\$925	NOACB1813	\$695	18"	13"	NOCFPCB1813	\$1015	NOPCB1813	\$695
	18"	16"	NOCFACB1816	\$960	NOACB1816	\$695	18"	16"	NOCFPCB1816	\$1055	NOPCB1816	\$695
	18"	19"			NOACB1819	\$695	18"	19"			NOPCB1819	\$695
	19"	16"			NOACB1916	\$695	19"	16"			NOPCB1916	\$695
	19"	19"			NOACB1919	\$695	19"	19"			NOPCB1919	\$695
	20"	13"	NOCFACB2013	\$925	NOACB2013	\$695	20"	13"	NOCFPCB2013	\$1015	NOPCB2013	\$695
	20"	16"	NOCFACB2016	\$960	NOACB2016	\$695	20"	16"	NOCFPCB2016	\$1055	NOPCB2016	\$695
	20"	19"			NOACB2019	\$695	20"	19"			NOPCB2019	\$695
	21"	16"			NOACB2116	\$745	21"	16"			NOPCB2116	\$745
	21"	19"			NOACB2119	\$745	21"	19"			NOPCB2119	\$745
	22"	13"			NOACB2213	\$745	22"	13"			NOPCB2213	\$745
	22"	16"			NOACB2216	\$745	22"	16"			NOPCB2216	\$745
	22"	19"			NOACB2219	\$745	22"	19"			NOPCB2219	\$745

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	NOACC _ _ _ _	130	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	NOPCC _ _ _ _	130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

FOLDING ON CHAIR - USA ORDER FORM

Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

THORACIC SHAPE

	xtend™	WIDTH	HEIGHT	✓	\$
	**NHAACB (L)1213/17	12"	13"-17"		\$695
	NHAACB (L)1413/17	14"	13"-17"		\$695
	NHAACB (L)1613/17	16"	13"-17"		\$695
	NHAACB (L)1813/17	18"	13"-17"		\$745
	NHAACB (L)2013/17	20"	13"-17"		\$745
	NHAACB (L)2213/17	22"	13"-17"		\$745

HCPCS CODE:
E2620-E2621

	xtend™ DEEP	WIDTH	HEIGHT	✓	\$
	NHAPCB1616/20	16"	16"-20"		\$695
	NHAPCB1816/20	18"	16"-20"		\$695
	NHAPCB2016/20	20"	16"-20"		\$695
	NHAPCB2216/20	22"	16"-20"		\$745

HCPCS CODE:
E2620-E2621

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	N3LHAACC _ _ _ _	130	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	N3HAPCC _ _ _ _	130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

Classic™ PELVIC BACK SUPPORT

Classic™ DEEP PELVIC BACK SUPPORT

PELVIC SHAPE

	Classic™	WIDTH	HEIGHT	✓	\$
	NACB1616	16"-18"	16"		\$695
	NACB1618	16"-18"	18"		\$695
	NACB1816	18"-20"	16"		\$695
	NACB1818	18"-20"	18"		\$695
	NACB1820	18"-20"	20"		\$695
	NACB2020	20"-22"	20"		\$695
	NACB2216	22"-24"	16"		\$745
	NACB2218	22"-24"	18"		\$745
	NACB2220	22"-24"	20"		\$745

HCPCS CODE:
E2620-E2621

	Classic DEEP™	WIDTH	HEIGHT	✓	\$
	NPCB1616	16"-18"	16"		\$695
	NPCB1618	16"-18"	18"		\$695
	NPCB1818	18"-20"	18"		\$695
	NPCB1820	18"-20"	20"		\$695
	NPCB2020	20"-22"	20"		\$695
	NPCB2216	22"-24"	16"		\$745
	NPCB2218	22"-24"	18"		\$745
	NPCB2220	22"-24"	20"		\$745

HCPCS CODE:
E2620-E2621

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	NACC _ _ _ _	130	

ITEM	WIDTH X DEPTH	\$	QTY
Replacement Cover & Foam	NPCC _ _ _ _	130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

BACK SUPPORT ACCESSORIES

	MODEL #	DESCRIPTION	FITS BACK WIDTH	\$	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	205	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	70	

MultiFit™ HEAD SUPPORT - MULTI-ADJUSTABLE PADS

MultiFit™ HEAD SUPPORT - FIXTURES

		
HCPCS CODE: E0955		
MODEL # NIB-0784	✓	\$
10" x 5" Large Rectangular		\$205

		
HCPCS CODE: E0955		
MODEL # NIB-0818	✓	\$
14" x 5" Extra-Large Rectangular		\$205

		
HCPCS CODE: E1033		
MODEL # NIB-0747	✓	\$
Multi-Axis Hardware		\$150

		
HCPCS CODE: E1033		
MODEL # NIB-0823	✓	\$
Head Support Receiver Bracket		\$60

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

EXTRA COVERS

Extra Outer
smartx^{3D} Cover



Extra Outer
Incontinent Cover



ITEM	WIDTH X DEPTH	\$	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	__ __ x __ __	130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	__ __ X __ __	150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	__ __ X __ __	130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	__ __ X __ __	150	

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com