



**1- Purchasing information**

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_

Company Name:\* \_\_\_\_\_

Purchasing contact:\* \_\_\_\_\_

Phone:\* \_\_\_\_\_

Email:\* \_\_\_\_\_

Contact preference:\* Phone ☐ Email ☐

Marked for: \_\_\_\_\_

**3- Shipping information**

Ship to:\* \_\_\_\_\_

Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

Zip Code: \_\_\_\_\_

Phone: \_\_\_\_\_

**2- ATP/Clinician information**

ATP/Clinician:\* \_\_\_\_\_

Phone:\* \_\_\_\_\_

Email:\* \_\_\_\_\_

Contact preference:\* Phone ☐ Email ☐

|                           |                          |                                                                                                                                                                               |                                           |                                           |                                           |                                           |                                           |                                |       |                    |
|---------------------------|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------|-------------------------------------------|-------------------------------------------|-------------------------------------------|-------------------------------------------|--------------------------------|-------|--------------------|
| * Required                |                          | MA082 <input type="checkbox"/>                                                                                                                                                |                                           |                                           |                                           |                                           |                                           | MA089 <input type="checkbox"/> |       | Taxes not included |
|                           |                          | English Manual (STD)                                                                                                                                                          |                                           |                                           |                                           |                                           |                                           | Spanish Manual                 |       |                    |
| MODEL                     |                          |                                                                                                                                                                               |                                           |                                           |                                           |                                           |                                           | USD                            |       |                    |
| WC03                      | <input type="checkbox"/> | Ultralight folding Helio wheelchair (265lb weight capacity)<br><b>A7 package:</b> side frames, symmetrical hydroformed crossbrace and back canes made of 7000 series aluminum |                                           |                                           |                                           |                                           |                                           |                                | K0005 | \$2,395            |
| WC08                      | <input type="checkbox"/> | <b>C2 carbon fiber upgrade:</b> side frames, symmetrical molded crossbrace and back canes made of ultralight carbon fiber                                                     |                                           |                                           |                                           |                                           |                                           |                                |       | \$2,295            |
| WC09                      | <input type="checkbox"/> | <b>HD Heavy Duty package</b> (350 lb weight capacity )                                                                                                                        |                                           |                                           |                                           |                                           |                                           |                                |       | \$550              |
| Frame                     |                          |                                                                                                                                                                               |                                           |                                           |                                           |                                           |                                           |                                |       |                    |
| SEAT WIDTH <sup>3</sup>   |                          |                                                                                                                                                                               |                                           |                                           |                                           |                                           |                                           |                                |       |                    |
| SW                        | <input type="checkbox"/> | 14" <input type="checkbox"/>                                                                                                                                                  | 15" <input type="checkbox"/>              | 16" <input type="checkbox"/>              | 17" <input type="checkbox"/>              | 18" <sup>1</sup> <input type="checkbox"/> | 19" <sup>1</sup> <input type="checkbox"/> |                                |       | NC                 |
|                           | <input type="checkbox"/> |                                                                                                                                                                               | 20" <sup>1</sup> <input type="checkbox"/> | 21" <sup>2</sup> <input type="checkbox"/> | 22" <sup>2</sup> <input type="checkbox"/> |                                           |                                           |                                | E2201 | \$525              |
| SEAT DEPTH                |                          |                                                                                                                                                                               |                                           |                                           |                                           |                                           |                                           |                                |       |                    |
| SD                        | <input type="checkbox"/> | 14" <input type="checkbox"/>                                                                                                                                                  | 15" <input type="checkbox"/>              | 16" <input type="checkbox"/>              | 17" <input type="checkbox"/>              | 18" <input type="checkbox"/>              | 19" <input type="checkbox"/>              |                                |       | NC                 |
|                           | <input type="checkbox"/> |                                                                                                                                                                               | 20" <input type="checkbox"/>              |                                           |                                           |                                           |                                           |                                | E2203 | \$575              |
| FRAME LENGTH <sup>5</sup> |                          |                                                                                                                                                                               |                                           |                                           |                                           |                                           |                                           |                                |       |                    |
| SP01                      | <input type="checkbox"/> | Standard frame length configuration                                                                                                                                           |                                           |                                           |                                           |                                           |                                           |                                |       | STD                |
| SP02                      | <input type="checkbox"/> | Frame length 1" longer than standard configuration                                                                                                                            |                                           |                                           |                                           |                                           |                                           |                                |       | NC                 |
| SP03                      | <input type="checkbox"/> | Frame length 2" longer than standard configuration                                                                                                                            |                                           |                                           |                                           |                                           |                                           |                                |       | NC                 |
| SP04                      | <input type="checkbox"/> | Frame length 1" shorter than standard configuration                                                                                                                           |                                           |                                           |                                           |                                           |                                           |                                |       | NC                 |
| SP05                      | <input type="checkbox"/> | Frame length 2" shorter than standard configuration                                                                                                                           |                                           |                                           |                                           |                                           |                                           |                                |       | NC                 |
| TRANSIT TIE DOWN          |                          |                                                                                                                                                                               |                                           |                                           |                                           |                                           |                                           |                                |       |                    |
| AC04                      | <input type="checkbox"/> | Transit Tie Down (Unoccupied)                                                                                                                                                 |                                           |                                           |                                           |                                           |                                           |                                |       | NC                 |
| AC05-S                    | <input type="checkbox"/> | WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup> Short belt 20"-41"                                                                                                       |                                           |                                           |                                           |                                           |                                           |                                |       | \$365              |
| AC05-L                    | <input type="checkbox"/> | WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup> Long belt 20"-46"                                                                                                        |                                           |                                           |                                           |                                           |                                           |                                |       | \$365              |
| AC35                      | <input type="checkbox"/> | WC-19 Transit Tie Down without belt <sup>4</sup>                                                                                                                              |                                           |                                           |                                           |                                           |                                           |                                |       | \$350              |

1- Available with HD Kit (350 lb); 2- Only with HD Kit (350 lb)

4- May interfere with swing away armrest and Clamp on height adjustable push handle; not available with HD option when carbon upgrade is selected.

5- See table on page 8

[3. Width calculator](#)



**NEWTON ULTRALIGHT PACKAGE**

|      |                          |                                                                                                                                                                                                                                                  |  |  |  |            |  |       |
|------|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|------------|--|-------|
| PK01 | <input type="checkbox"/> | Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels                                                                                                                                                       |  |  |  | Save \$275 |  | \$530 |
| VP01 | <input type="checkbox"/> | Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates <sup>1</sup> |  |  |  | K0108      |  | \$95  |

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo **\*VP01**; Recommended with black chair; See picture of the product on page 6.

## Helio FRONT Seat-to-floor Height - Please select your choice

| Forks / Position |    | 3" caster        |         |         | 4" caster        |         |         | 5" caster        |         |         | 6" caster        |         |         | 7" caster        |         |         | 8" caster        |         |         |
|------------------|----|------------------|---------|---------|------------------|---------|---------|------------------|---------|---------|------------------|---------|---------|------------------|---------|---------|------------------|---------|---------|
|                  |    | Stem bolt Length |         |         | Stem bolt Length |         |         | Stem bolt Length |         |         | Stem bolt Length |         |         | Stem bolt Length |         |         | Stem bolt Length |         |         |
|                  |    | STD              | +1"     | +2"     | STD              | +1"     | +2"     | STD              | +1"     | +2"     | STD              | +1"     | +2"     | STD              | +1"     | +2"     | STD              | +1"     | +2"     |
| 4"               | P1 | 13"              | 14"     | 15"     | 13 1/2"          | 14 1/2" | 15 1/2" |                  |         |         |                  |         |         |                  |         |         |                  |         |         |
|                  | P2 | 13 1/2"          | 14 1/2" | 15 1/2" | 14"              | 15"     | 16"     | 14 1/2"          | 15 1/2" | 16 1/2" |                  |         |         |                  |         |         |                  |         |         |
| 5"               | P1 |                  |         |         |                  |         |         |                  |         |         |                  |         |         |                  |         |         |                  |         |         |
|                  | P2 |                  |         |         | 14"              | 15"     | 16"     | 14 1/2"          | 15 1/2" | 16 1/2" |                  |         |         |                  |         |         |                  |         |         |
|                  | P3 | 14"              | 15"     | 16"     | 14 1/2"          | 15 1/2" | 16 1/2" | 15"              | 16"     | 17"     | 15 1/2"          | 16 1/2" | 17 1/2" |                  |         |         |                  |         |         |
|                  | P4 | 14 1/2"          | 15 1/2" | 16 1/2" | 15"              | 16"     | 17"     | 15 1/2"          | 16 1/2" | 17 1/2" | 16"              | 17"     | 18"     | 16 1/2"          | 17 1/2" | 18 1/2" |                  |         |         |
| 7"               | P1 |                  |         |         |                  |         |         |                  |         |         | 16 1/2"          | 17 1/2" | 18 1/2" | 17"              | 18"     | 19"     |                  |         |         |
|                  | P2 |                  |         |         |                  |         |         | 16 1/2"          | 17 1/2" | 18 1/2" | 17"              | 18"     | 19"     | 17 1/2"          | 18 1/2" | 19 1/2" |                  |         |         |
|                  | P3 |                  |         |         |                  |         |         | 17"              | 18"     | 19"     | 17 1/2"          | 18 1/2" | 19 1/2" | 18"              | 19"     | 20"     | 18 1/2"          | 19 1/2" | 20 1/2" |
|                  | P4 | 16 1/2"          | 17 1/2" | 18 1/2" | 17"              | 18"     | 19"     | 17 1/2"          | 18 1/2" | 19 1/2" | 18"              | 19"     | 20"     | 18 1/2"          | 19 1/2" | 20 1/2" | 19"              | 20"     | 21"     |
|                  | P5 | 17"              | 18"     | 19"     | 17 1/2"          | 18 1/2" | 19 1/2" | 18"              | 19"     | 20"     | 18 1/2"          | 19 1/2" | 20 1/2" | 19"              | 20"     | 21"     | 19 1/2"          | 20 1/2" | 21 1/2" |
| Froglegs         | P1 |                  |         |         | 14 1/2"          | 15 1/2" | 16 1/2" |                  |         |         |                  |         |         |                  |         |         |                  |         |         |
|                  | P2 | 14 1/2"          | 15 1/2" | 16 1/2" | 15"              | 16"     | 17"     | 15 1/2"          | 16 1/2" | 17 1/2" |                  |         |         |                  |         |         |                  |         |         |
|                  | P3 | 15"              | 16"     | 17"     | 15 1/2"          | 16 1/2" | 17 1/2" | 16"              | 17"     | 18"     | 16 1/2"          | 17 1/2" | 18 1/2" |                  |         |         |                  |         |         |

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

### FRONT CASTERS

| CASTERS            |                          |                                                                               |  | Available              |             |       |
|--------------------|--------------------------|-------------------------------------------------------------------------------|--|------------------------|-------------|-------|
| CA01-03 to CA01-08 | <input type="checkbox"/> | Composite wheel w/ PU tire (1")                                               |  | 3", 4", 5", 6", 7", 8" |             | NC    |
| CA10 to CA11       | <input type="checkbox"/> | Pneumatic casters (1 1/4")                                                    |  | 6", 8"                 | E2214/E2215 | \$165 |
| CA14 to CA16       | <input type="checkbox"/> | Newton UltraCaster Composite (1 1/2")                                         |  | 4", 5", 6"             |             | \$95  |
| CA17 to CA19       | <input type="checkbox"/> | Newton UltraCaster Aluminum (1 1/2")                                          |  | 4", 5", 6"             |             | \$215 |
| CASTER OPTIONS     |                          |                                                                               |  |                        |             |       |
| CA12               | <input type="checkbox"/> | Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) <sup>1</sup>      |  |                        | K0073       | \$300 |
| CA13               | <input type="checkbox"/> | Carbon fiber Frog Legs suspension forks (Not available with caster pin locks) |  |                        | E1015       | \$490 |

1- Not available with WC-19

## Helio - REAR Seat-to-floor Height

| WHEEL SIZES              |                 | Height available |           | Specify desired height:<br><br>_____ |
|--------------------------|-----------------|------------------|-----------|--------------------------------------|
| <input type="checkbox"/> | 20" (451) wheel | 12 1/2           | to 17 1/4 |                                      |
| <input type="checkbox"/> | 22" (501) wheel | 13 1/2           | to 18 1/4 |                                      |
| <input type="checkbox"/> | 24" (540) wheel | 14 1/2           | to 19 1/4 |                                      |
| <input type="checkbox"/> | 25" (559) wheel | 15               | to 19 3/4 |                                      |
| <input type="checkbox"/> | 26" (590) wheel | 15 1/2           | to 20 1/4 |                                      |

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

| REAR WHEELS                |                          |                                                                                |                                                                        |                                                                                                                                                                                                                                                                                                                                                   |                         |  |             |         |  |
|----------------------------|--------------------------|--------------------------------------------------------------------------------|------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|--|-------------|---------|--|
| ONE ARM DRIVE <sup>1</sup> |                          |                                                                                | Only available with Newton One rear wheels, permanent axles, 0° camber |                                                                                                                                                                                                                                                                                                                                                   |                         |  |             |         |  |
| HR22                       | <input type="checkbox"/> | One arm drive - aluminum inner handrim <sup>2</sup>                            | Left <input type="checkbox"/>                                          | Right <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                    |                         |  | E0958       | \$1,495 |  |
| HR60                       | <input type="checkbox"/> | One arm drive - plastic coated inner handrim <sup>2</sup>                      | Left <input type="checkbox"/>                                          | Right <input type="checkbox"/>                                                                                                                                                                                                                                                                                                                    |                         |  | E0958       | \$1,600 |  |
| HR98                       | <input type="checkbox"/> | Omit handrim on opposite side of the One Arm Drive (OAD)                       |                                                                        |                                                                                                                                                                                                                                                                                                                                                   |                         |  |             | NC      |  |
| REAR WHEELS                |                          |                                                                                |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | Available               |  |             |         |  |
| RW01                       | <input type="checkbox"/> | Mag <sup>3</sup>                                                               |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 20", 22" 24"            |  |             | \$57    |  |
| RW02                       | <input type="checkbox"/> | Newton One - spoke wheel <sup>4</sup>                                          |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 20", 22" 24"            |  |             | NC      |  |
| RW03                       | <input type="checkbox"/> | Newton Gravity - ultralight wheel                                              |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 20", 22", 24", 25", 26" |  |             | \$360   |  |
| RW04                       | <input type="checkbox"/> | Spinergy Spox                                                                  | Blk STD <input type="checkbox"/>                                       | Wht <sup>5</sup> <input type="checkbox"/> Blue <sup>5</sup> <input type="checkbox"/> Red <sup>5</sup> <input type="checkbox"/> Yel <sup>5</sup> <input type="checkbox"/>                                                                                                                                                                          | 22", 24", 25", 26"      |  |             | \$915   |  |
| RW05                       | <input type="checkbox"/> | Spinergy LX                                                                    | Blk STD <input type="checkbox"/>                                       | Pur <sup>5</sup> <input type="checkbox"/> Wht <sup>5</sup> <input type="checkbox"/> Blue <sup>5</sup> <input type="checkbox"/> Red <sup>5</sup> <input type="checkbox"/> Yel <sup>5</sup> <input type="checkbox"/> Pink <sup>5</sup> <input type="checkbox"/> Grn <sup>5</sup> <input type="checkbox"/> Org <sup>5</sup> <input type="checkbox"/> | 22", 24", 25", 26"      |  |             | \$1,145 |  |
| RW07                       | <input type="checkbox"/> | Spinergy CLX <sup>6</sup>                                                      | Blk STD <input type="checkbox"/>                                       | Pur <sup>5</sup> <input type="checkbox"/> Wht <sup>5</sup> <input type="checkbox"/> Blue <sup>5</sup> <input type="checkbox"/> Red <sup>5</sup> <input type="checkbox"/> Yel <sup>5</sup> <input type="checkbox"/> Pink <sup>5</sup> <input type="checkbox"/> Grn <sup>5</sup> <input type="checkbox"/> Org <sup>5</sup> <input type="checkbox"/> | 24", 25"                |  |             | \$1,625 |  |
| TIRES                      |                          |                                                                                |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | Available               |  |             |         |  |
| AIRLESS                    |                          |                                                                                |                                                                        |                                                                                                                                                                                                                                                                                                                                                   |                         |  |             |         |  |
| TI01                       | <input type="checkbox"/> | Soft Urethane 1 3/8" - Low Tread, black <sup>9</sup>                           |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 20", 22" 24"            |  | E2218       | STD     |  |
| TI02                       | <input type="checkbox"/> | Soft Urethane 1 3/8 " - Medium Tread, black <sup>7,9</sup>                     |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 22", 24", 25", 26"      |  | E2218       | \$97    |  |
| TI04                       | <input type="checkbox"/> | Soft Urethane 1" - No tread, black - <i>Shox</i>                               |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 22", 24", 25", 26"      |  | E2218       | \$118   |  |
| TI08                       | <input type="checkbox"/> | Pneumatic Airless insert 1 3/8" - Medium Tread <sup>8</sup>                    |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 22", 24", 26"           |  | E2211/E2213 | \$95    |  |
| PNEUMATIC                  |                          |                                                                                |                                                                        |                                                                                                                                                                                                                                                                                                                                                   |                         |  |             |         |  |
| TI05                       | <input type="checkbox"/> | Pneumatic 1 3/8" - Medium Tread <sup>9</sup>                                   |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 20", 22", 24", 25", 26" |  | E2211/E2212 | NC      |  |
| TI06                       | <input type="checkbox"/> | Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>      |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 20", 22", 24", 25", 26" |  | E2211/E2212 | \$105   |  |
| TI07                       | <input type="checkbox"/> | Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i> |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 24", 25", 26"           |  | E2211/E2212 | \$150   |  |
| TI03                       | <input type="checkbox"/> | Pneumatic 2" - Large Tread, black <sup>10</sup>                                |                                                                        |                                                                                                                                                                                                                                                                                                                                                   | 22", 24", 25"           |  | E2211/E2212 | \$115   |  |

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available; 3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additional 10 days lead-time; 6- Only available with 1" tires; 7- 22" is only available with Mags and Newton One rear wheels; 25" and 26" are not available with Newton Gravity wheels 8 - Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels 9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;

| HANDRIMS                |                          |                                                                                  |  |                         |  |       |
|-------------------------|--------------------------|----------------------------------------------------------------------------------|--|-------------------------|--|-------|
| HANDRIM                 |                          |                                                                                  |  | Available               |  |       |
| HR01                    | <input type="checkbox"/> | Aluminum anodized                                                                |  | 20", 22", 24", 25", 26" |  | STD   |
| HR61                    | <input type="checkbox"/> | Aluminum black hard anodized <sup>2</sup>                                        |  | 20", 22", 24", 25"      |  | \$75  |
| HR02                    | <input type="checkbox"/> | Plastic coated <sup>1</sup>                                                      |  | 20", 22", 24", 25", 26" |  | \$110 |
| HR03                    | <input type="checkbox"/> | High friction coated <sup>1,2,3</sup>                                            |  | 22", 24", 25", 26"      |  | \$375 |
| HR07-HR18 and HR65-HR71 | <input type="checkbox"/> | Newton Air Grip <sup>2</sup> <i>See on page 6 for color choice and price</i>     |  | 20", 22", 24", 25", 26" |  |       |
| HR28                    | <input type="checkbox"/> | Surge <sup>1,2,3,4</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)    |  | 22", 24", 25", 26"      |  | \$450 |
| HR05                    | <input type="checkbox"/> | Surge LT <sup>1,2,3,4</sup> (1 1/8" Oval handrim w/ Gripton high friction strip) |  | 22", 24", 25", 26"      |  | \$495 |
| HR04                    | <input type="checkbox"/> | Natural Fit <sup>1,2,3,4</sup>                                                   |  | 22", 24", 25", 26"      |  | \$395 |
| HR62                    | <input type="checkbox"/> | Simi H <sup>1,2,3</sup> (1 3/8" Oval handrim with high friction strip)           |  | 22", 24", 25", 26"      |  | \$395 |
| HR63                    | <input type="checkbox"/> | Nova H <sup>1,2,3</sup> (1 1/4" Oval handrim with high friction strip)           |  | 22", 24", 25", 26"      |  | \$445 |
| HR64                    | <input type="checkbox"/> | Optimum H <sup>1,2,3</sup> (1 1/8" heart-shaped handrim)                         |  | 22", 24", 25", 26"      |  | \$345 |
| HR99                    | <input type="checkbox"/> | Omit handrim                                                                     |  |                         |  | NC    |

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive; 4 - 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

| HANDRIMS OPTION                                                 |                          |                                                                                                                                                  |  |           |  |       |
|-----------------------------------------------------------------|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|--|-----------|--|-------|
| ASSEMBLY POSITION                                               |                          |                                                                                                                                                  |  | Available |  |       |
| HR29                                                            | <input type="checkbox"/> | Standard Handrim position <i>Full tab</i>                                                                                                        |  |           |  | STD   |
| HR030                                                           | <input type="checkbox"/> | Narrow Handrim position no cut <i>Full tab</i>                                                                                                   |  |           |  | NC    |
| HR30                                                            | <input type="checkbox"/> | Narrow Handrim position with cut <sup>5</sup> <i>Short tab 2 positions</i>                                                                       |  |           |  | \$30  |
| HR31                                                            | <input type="checkbox"/> | Super Narrow Handrim <sup>5</sup> <i>Short tab 2 positions</i> STD <input type="checkbox"/> <i>Short tab 1 position</i> <input type="checkbox"/> |  |           |  | \$30  |
| THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS |                          |                                                                                                                                                  |  | Available |  |       |
| HR23                                                            | <input type="checkbox"/> | Standard grip                                                                                                                                    |  |           |  | \$95  |
| HR27                                                            | <input type="checkbox"/> | Super grip                                                                                                                                       |  |           |  | \$150 |
| HR099                                                           | <input type="checkbox"/> | Omit thumb piece                                                                                                                                 |  |           |  | NC    |


5- Only available with Aluminum and Air Grip handrims;

| AXLE <sup>1</sup> |                          |                                            |  | Available |       |       |
|-------------------|--------------------------|--------------------------------------------|--|-----------|-------|-------|
| AX01              | <input type="checkbox"/> | Quick release axle                         |  |           |       | STD   |
| AX02              | <input type="checkbox"/> | Permanent axle                             |  |           |       | NC    |
| AX03              | <input type="checkbox"/> | Quad release axle                          |  |           | K0108 | \$130 |
| AXLE PLATE        |                          |                                            |  | Available |       |       |
| AX04              | <input type="checkbox"/> | Amputee adjustable axle plate <sup>2</sup> |  |           | E0959 | \$205 |

1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

| REAR WHEEL CAMBER |                          |                                                          |       |
|-------------------|--------------------------|----------------------------------------------------------|-------|
| CM00              | <input type="checkbox"/> | 0° camber                                                | STD   |
| CM03              | <input type="checkbox"/> | 3° camber                                                | NC    |
| CM06              | <input type="checkbox"/> | 6° camber <i>(not available with amputee axle plate)</i> | \$165 |

#### CENTER OF GRAVITY



**Choose position of rear wheel axle plate:**

Measure distance X from front of back post to center of axle receiver.

| Most stable                                                                               |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |
|-------------------------------------------------------------------------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| X =                                                                                       | AMPUTEE                  | 1"                       | 1/4"                     | 1/2"                     | 3/4"                     | 2"                       | 1/4"                     | 1/2"                     | 3/4"                     | 3"                       | 1/4"                     | 1/2"                     | 3/4"                     | 4"                       |
|                                                                                           | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| <i>If not indicated, default setting is at 1.5"</i>                                       |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |
| <i>Possible conflict between a seat depth under 16" and a center of gravity over 2,5"</i> |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |
| Only on carbon frame                                                                      |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |                          |

**If not indicated, default setting is at 1.5"**

Possible conflict between a seat depth under 16" and a center of gravity over 2,5".

Possible conflict between a 26" rear wheel and a center of gravity over 2,5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

Only on carbon frame

#### WHEEL LOCK

|      |                          |                                                                                                                |  |       |       |
|------|--------------------------|----------------------------------------------------------------------------------------------------------------|--|-------|-------|
| WL08 | <input type="checkbox"/> | Newton Wheel lock - push to lock                                                                               |  |       | STD   |
| WL16 | <input type="checkbox"/> | Aluminum push to lock <b>with</b> extensions                                                                   |  | E0961 | \$77  |
| WL02 | <input type="checkbox"/> | Pull to lock                                                                                                   |  |       | NC    |
| WL04 | <input type="checkbox"/> | Newton grade aid push to lock <i>(works best with pneumatic or high pressure tire)</i>                         |  | E0974 | \$155 |
| WL13 | <input type="checkbox"/> | 6" removable extension handle <sup>1</sup> *VP01                                                               |  | E0961 | \$77  |
| WL05 | <input type="checkbox"/> | Scissors Lock                                                                                                  |  |       | \$155 |
| WL06 | <input type="checkbox"/> | Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>                                       |  |       | \$215 |
| WL15 | <input type="checkbox"/> | Unilateral wheel lock - pull to lock <sup>2</sup> Left <input type="checkbox"/> Right <input type="checkbox"/> |  |       | \$330 |
| WL99 | <input type="checkbox"/> | Omit wheel locks                                                                                               |  |       | NC    |

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

[https://www.motioncomposites.com/en\\_ca/support-and-education/how-to-documents](https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents)

#### FOOTREST

| SWING-IN/SWING-OUT FOOTREST |                          |                                                                                                                              |              |  |            |          |                |
|-----------------------------|--------------------------|------------------------------------------------------------------------------------------------------------------------------|--------------|--|------------|----------|----------------|
| FR03-01                     | <input type="checkbox"/> | 60° Swing-in/swing-out footrest - 13.5" - 17"                                                                                |              |  | Set at*:   |          | NC             |
| FR03-02                     | <input type="checkbox"/> | 60° Swing-in/swing-out footrest - 14.5" - 19.5"                                                                              |              |  | Set at*:   |          | NC             |
| FR03-03                     | <input type="checkbox"/> | 60° Swing-in/swing-out footrest - 15.5" - 20.5"                                                                              |              |  | Set at*:   |          | NC             |
| FR04-01                     | <input type="checkbox"/> | 70° Swing-in/swing-out footrest - 13.5" - 17"                                                                                |              |  | Set at*:   |          | NC             |
| FR04-02                     | <input type="checkbox"/> | 70° Swing-in/swing-out footrest - 14.5" - 19.5"                                                                              |              |  | Set at*:   |          | NC             |
| FR04-03                     | <input type="checkbox"/> | 70° Swing-in/swing-out footrest - 15.5" - 20.5"                                                                              |              |  | Set at*:   |          | NC             |
| FR15-01                     | <input type="checkbox"/> | 80° Swing-in/swing-out footrest - 13.5" - 16" <sup>5</sup>                                                                   |              |  | Set at*:   |          | NC             |
| FR15-02                     | <input type="checkbox"/> | 80° Swing-in/swing-out footrest - 14.5" - 19.5" <sup>5</sup>                                                                 |              |  | Set at*:   |          | NC             |
| FR15-03                     | <input type="checkbox"/> | 80° Swing-in/swing-out footrest - 15.5" - 20.5" <sup>5</sup>                                                                 |              |  | Set at*:   |          | NC             |
| FR06                        | <input type="checkbox"/> | 90° swing-in/swing-out with front mount - Standard (no extension) <sup>1</sup>                                               | 7.5" - 11"   |  | Set at*:   |          | NC             |
| FR07                        | <input type="checkbox"/> | 90° swing-in/swing-out with front mount - Short extension 2" <sup>1, 2</sup>                                                 | 9.5" - 13"   |  | Set at*:   |          | NC             |
| FR08                        | <input type="checkbox"/> | 90° swing-in/swing-out with front mount - Medium extension 4" <sup>1, 2</sup>                                                | 9.5" - 14.5" |  | Set at*:   |          | NC             |
| FOOTREST OPTIONS            |                          |                                                                                                                              |              |  |            |          |                |
| FR09                        | <input type="checkbox"/> | High mount <i>(footplate attached to the hanger with clamps)</i> <sup>3</sup>                                                |              |  | 5" - 13.5" | Set at*: | \$85           |
| FR10                        | <input type="checkbox"/> | Elevating legrest with calf pad, Length adjustable <sup>4</sup> Left <input type="checkbox"/> Right <input type="checkbox"/> |              |  | 14" - 22"  | Set at:  | E0990    \$245 |
| FR99                        | <input type="checkbox"/> | Omit Footrest                                                                                                                |              |  |            |          | NC             |

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest; \* +/-1/2"

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

#### FOOTPLATES

|      |                          |                                                                                |  |       |       |
|------|--------------------------|--------------------------------------------------------------------------------|--|-------|-------|
| FP01 | <input type="checkbox"/> | Standard - Composite <sup>1</sup>                                              |  |       | STD   |
| FP02 | <input type="checkbox"/> | Newton adjustable angle - Composite *VP01                                      |  | K0040 | \$215 |
| FP03 | <input type="checkbox"/> | Newton adjustable angle - Aluminum                                             |  | K0040 | \$235 |
| FP04 | <input type="checkbox"/> | Newton adjustable angle - Carbon fiber <i>(with rubber inserts)</i>            |  |       | \$445 |
| FP23 | <input type="checkbox"/> | Newton adjustable angle extra deep - Aluminum                                  |  | K0040 | \$235 |
| FP06 | <input type="checkbox"/> | One piece adjustable angle flip-up - Aluminum                                  |  | K0040 | \$195 |
| FP07 | <input type="checkbox"/> | One piece adjustable angle flip-up - Carbon fiber <i>(with rubber inserts)</i> |  |       | \$465 |

#### FOOTPLATE ACCESSORIES

|      |                          |                                          |  |       |       |
|------|--------------------------|------------------------------------------|--|-------|-------|
| FP08 | <input type="checkbox"/> | Length adjustable heel loop <sup>2</sup> |  | E0951 | \$41  |
| FP09 | <input type="checkbox"/> | Calf strap                               |  | K0038 | \$39  |
| FP30 | <input type="checkbox"/> | Bodypoint Aeromesh Padded calf strap     |  | K0038 | \$118 |

1- Not available with front mount; 2- not available with one piece footplate

#### SWING AWAY RESIDUAL LIMB SUPPORT

|      |                          |                                        |                                                              |  |       |       |
|------|--------------------------|----------------------------------------|--------------------------------------------------------------|--|-------|-------|
| FR27 | <input type="checkbox"/> | Residual limb support - flat 6"x8"     | Left <input type="checkbox"/> Right <input type="checkbox"/> |  | E1020 | \$420 |
| FR28 | <input type="checkbox"/> | Residual limb support - flat 8"x10"    | Left <input type="checkbox"/> Right <input type="checkbox"/> |  | E1020 | \$420 |
| FR29 | <input type="checkbox"/> | Residual limb support - flat 10"x10"   | Left <input type="checkbox"/> Right <input type="checkbox"/> |  | E1020 | \$420 |
| FR30 | <input type="checkbox"/> | Residual limb support - flat 10"x14"   | Left <input type="checkbox"/> Right <input type="checkbox"/> |  | E1020 | \$420 |
| FR32 | <input type="checkbox"/> | Residual limb support - contour 6"x8"  | Left <input type="checkbox"/> Right <input type="checkbox"/> |  | E1020 | \$420 |
| FR33 | <input type="checkbox"/> | Residual limb support - contour 8"x10" | Left <input type="checkbox"/> Right <input type="checkbox"/> |  | E1020 | \$420 |

## SEATING AND BACKREST

| BACK UPHOLSTERY                                                                                                    |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
|--------------------------------------------------------------------------------------------------------------------|--------------------------|----------------------------------------------------------------------------------------|--|------------------------------------------------------------|--------------------------------------------------------------------------------------------|-------------------------------|-------------------------------------------------------------|-------|
| BH07                                                                                                               | <input type="checkbox"/> | Slip-on Active-Mesh back upholstery                                                    |  |                                                            |                                                                                            |                               |                                                             | STD   |
| BH08                                                                                                               | <input type="checkbox"/> | Adjustable tension Active-Mesh back upholstery                                         |  |                                                            |                                                                                            |                               | E2611/E2612                                                 | \$305 |
| BH99                                                                                                               | <input type="checkbox"/> | Omit Back upholstery                                                                   |  |                                                            |                                                                                            |                               |                                                             | NC    |
|                                                                                                                    | <input type="checkbox"/> | NXT back Support                                                                       |  |                                                            |                                                                                            |                               | See pricelist                                               |       |
| SEAT CUSHION                                                                                                       |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| SC01                                                                                                               | <input type="checkbox"/> | 2" Cushion                                                                             |  |                                                            |                                                                                            |                               |                                                             | \$170 |
| SC02                                                                                                               | <input type="checkbox"/> | 3" Cushion                                                                             |  |                                                            |                                                                                            |                               | E2601/E2602                                                 | \$190 |
| SC99                                                                                                               | <input type="checkbox"/> | Omit seat cushion                                                                      |  |                                                            |                                                                                            |                               |                                                             | STD   |
|                                                                                                                    | <input type="checkbox"/> | NXT cushion                                                                            |  |                                                            |                                                                                            |                               | See pricelist                                               |       |
| FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLES - For Aluminum chair (Not available with push handles options) |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| BC58                                                                                                               | <input type="checkbox"/> | Straight aluminum back cane                                                            |  | 9 - 15" <input type="checkbox"/>                           | 15 - 21" <input type="checkbox"/>                                                          |                               |                                                             | NC    |
| BC60                                                                                                               | <input type="checkbox"/> | Straight aluminum back cane                                                            |  | 21" - 24" <span style="float:right">increment of 1"</span> |                                                                                            | Set at <sup>2</sup> :         |                                                             | \$77  |
| BC01                                                                                                               | <input type="checkbox"/> | 8° bend aluminum back cane                                                             |  | 16" <input type="checkbox"/>                               | 18" <input type="checkbox"/> <span style="float:right">20" <input type="checkbox"/></span> |                               |                                                             | NC    |
| BC59                                                                                                               | <input type="checkbox"/> | 8° bend aluminum back cane                                                             |  | 21" - 24" <span style="float:right">increment of 1"</span> |                                                                                            | Set at <sup>2</sup> :         |                                                             | \$77  |
| FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)           |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| BC11                                                                                                               | <input type="checkbox"/> | Straight back cane                                                                     |  | 12 - 15" <input type="checkbox"/>                          | 15 - 18" <input type="checkbox"/> 18 - 21" STD <input type="checkbox"/>                    | increment of 1/2"             | Set at <sup>2</sup> :                                       | NC    |
| BC06                                                                                                               |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| BC12                                                                                                               | <input type="checkbox"/> | 8° Bend back cane                                                                      |  | 13 ½ - 15" <input type="checkbox"/>                        | 15 - 18" <input type="checkbox"/> 18 - 21" STD <input type="checkbox"/>                    | increment of 1/2"             | Set at <sup>2</sup> :                                       | \$60  |
| BC08                                                                                                               |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)               |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| BC13                                                                                                               | <input type="checkbox"/> | Straight back cane                                                                     |  | 15 - 18" <input type="checkbox"/>                          | 18 - 21" STD <input type="checkbox"/>                                                      | increment of 1/2"             | Set at <sup>2</sup> :                                       | \$215 |
| BC09                                                                                                               |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| BC14                                                                                                               | <input type="checkbox"/> | 8° bend back cane                                                                      |  | 15 - 18" <input type="checkbox"/>                          | 18 - 21" STD <input type="checkbox"/>                                                      | increment of 1/2"             | Set at <sup>2</sup> :                                       | \$270 |
| BC10                                                                                                               |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
|                                                                                                                    | Select back cane angle:  |                                                                                        |  | 85° <sup>1</sup> <input type="checkbox"/>                  | 90° <input type="checkbox"/> 95° STD <input type="checkbox"/>                              | 100° <input type="checkbox"/> | 105° <input type="checkbox"/> 110° <input type="checkbox"/> |       |
| PUSH HANDLES OPTIONS (ONE option at a time)                                                                        |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| PH01                                                                                                               | <input type="checkbox"/> | Push handles                                                                           |  |                                                            |                                                                                            |                               |                                                             | STD   |
| PH05                                                                                                               | <input type="checkbox"/> | Fold down push handles <sup>3,5</sup>                                                  |  |                                                            |                                                                                            |                               |                                                             | \$145 |
| PH06-920                                                                                                           | <input type="checkbox"/> | Clamp on height adjustable push handles 9" - 20" <sup>3,4,5</sup>                      |  |                                                            |                                                                                            |                               |                                                             | \$325 |
| PH06-923                                                                                                           | <input type="checkbox"/> | Clamp on height adjustable push handles 9" - 23" <sup>3,4,5</sup>                      |  |                                                            |                                                                                            |                               |                                                             | \$325 |
| PH02                                                                                                               | <input type="checkbox"/> | Omit push handles                                                                      |  |                                                            |                                                                                            |                               |                                                             | NC    |
| STROLLER HANDLE AND STABILIZER BAR                                                                                 |                          |                                                                                        |  |                                                            |                                                                                            |                               |                                                             |       |
| PH04                                                                                                               | <input type="checkbox"/> | Stroller handle (only available with push handle) - (not available with the HD option) |  |                                                            |                                                                                            |                               | K0108                                                       | \$165 |
| PH03                                                                                                               | <input type="checkbox"/> | Newton folding stabilizer bar (not compatible with stroller handle)                    |  |                                                            |                                                                                            |                               |                                                             | \$270 |

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265lb weight limit

4- Not compatible with swing away armrests, 8° bend back canes, 9-12" back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

## ARMRESTS

| HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) <sup>1, 4</sup> *VP01 |                          |                                                                |           |                                 |                                       |       |       |
|----------------------------------------------------------------------------------------|--------------------------|----------------------------------------------------------------|-----------|---------------------------------|---------------------------------------|-------|-------|
| AR01-812                                                                               | <input type="checkbox"/> | 10" Desk length armpad                                         | 8" - 12"  | Locked <input type="checkbox"/> | Unlocked STD <input type="checkbox"/> | E0973 | \$355 |
| AR01-1014                                                                              | <input type="checkbox"/> | 10" Desk length armpad                                         | 10" - 14" | Locked <input type="checkbox"/> | Unlocked STD <input type="checkbox"/> | E0973 | \$355 |
| AR02-812                                                                               | <input type="checkbox"/> | 14" Full length armpad                                         | 8" - 12"  | Locked <input type="checkbox"/> | Unlocked STD <input type="checkbox"/> | E0973 | \$355 |
| AR02-1014                                                                              | <input type="checkbox"/> | 14" Full length armpad                                         | 10" - 14" | Locked <input type="checkbox"/> | Unlocked STD <input type="checkbox"/> | E0973 | \$355 |
| HEIGHT ADJUSTABLE "T" ARMREST <sup>1</sup> *VP01                                       |                          |                                                                |           |                                 |                                       |       |       |
| AR13-610                                                                               | <input type="checkbox"/> | 10" Desk length armpad <sup>5</sup>                            | 6" - 10"  |                                 |                                       | E0973 | \$330 |
| AR03-812                                                                               | <input type="checkbox"/> | 10" Desk length armpad                                         | 8" - 12"  |                                 |                                       | E0973 | \$330 |
| AR03-1014                                                                              | <input type="checkbox"/> | 10" Desk length armpad                                         | 10" - 14" |                                 |                                       | E0973 | \$330 |
| AR14-610                                                                               | <input type="checkbox"/> | 14" Full length armpad <sup>5</sup>                            | 6" - 10"  |                                 |                                       | E0973 | \$330 |
| AR04-812                                                                               | <input type="checkbox"/> | 14" Full length armpad                                         | 8" - 12"  |                                 |                                       | E0973 | \$330 |
| AR04-1014                                                                              | <input type="checkbox"/> | 14" Full length armpad                                         | 10" - 14" |                                 |                                       | E0973 | \$330 |
| HEIGHT ADJUSTABLE SWING AWAY ARMREST <sup>2</sup>                                      |                          |                                                                |           |                                 |                                       |       |       |
| AR05                                                                                   | <input type="checkbox"/> | Tubular armrest <b>without</b> plastic side guards - Short pad | 8" - 11"  |                                 |                                       | E0973 | \$230 |
| AR19                                                                                   | <input type="checkbox"/> | Tubular armrest <b>with</b> plastic side guards - Short pad    | 8" - 11"  |                                 |                                       | E0973 | \$395 |
| AR23                                                                                   | <input type="checkbox"/> | Tubular armrest <b>without</b> plastic side guards - Long pad  | 8" - 11"  |                                 |                                       | E0973 | \$230 |
| AR24                                                                                   | <input type="checkbox"/> | Tubular armrest <b>with</b> plastic side guards - Long pad     | 8" - 11"  |                                 |                                       | E0973 | \$395 |
| AR06                                                                                   | <input type="checkbox"/> | 10" Desk length flat armpad <b>without</b> plastic side guards | 8" - 11"  |                                 |                                       | E0973 | \$230 |
| AR20                                                                                   | <input type="checkbox"/> | 10" Desk length flat armpad <b>with</b> plastic side guards    | 8" - 11"  |                                 |                                       | E0973 | \$395 |
| AR17                                                                                   | <input type="checkbox"/> | 14" Full length flat armpad <b>without</b> plastic side guards | 8" - 11"  |                                 |                                       | E0973 | \$230 |
| AR21                                                                                   | <input type="checkbox"/> | 14" Full length flat armpad <b>with</b> plastic side guards    | 8" - 11"  |                                 |                                       | E0973 | \$395 |
| ARMREST OPTIONS                                                                        |                          |                                                                |           |                                 |                                       |       |       |
| AR12                                                                                   | <input type="checkbox"/> | 14" Full Length armpad <sup>3</sup> - Gel Ovation              |           |                                 |                                       |       | \$165 |
| AR99                                                                                   | <input type="checkbox"/> | Omit armrest                                                   |           |                                 |                                       |       | NC    |

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

| ACCESSORIES                                                       |                          |                                                    |                        |                             |  |       |       |
|-------------------------------------------------------------------|--------------------------|----------------------------------------------------|------------------------|-----------------------------|--|-------|-------|
| SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i> |                          |                                                    |                        |                             |  |       |       |
| SG01                                                              | <input type="checkbox"/> | Plastic side guards <i>(removable)</i>             | Straight bracket - STD | Offset bracket <sup>1</sup> |  | K0108 | \$200 |
| SG02                                                              | <input type="checkbox"/> | Carbon fiber side guards <i>(fixed)</i>            |                        |                             |  |       | \$270 |
| SG03                                                              | <input type="checkbox"/> | Carbon fiber side guards <i>(removable)</i>        | Straight bracket - STD | Offset bracket <sup>1</sup> |  | K0108 | \$305 |
| SG06                                                              | <input type="checkbox"/> | Plastic fender side guards <i>(removable)</i>      | Straight bracket - STD | Offset bracket <sup>1</sup> |  |       | \$290 |
| SG04                                                              | <input type="checkbox"/> | Carbon fiber fender side guards <i>(removable)</i> | Straight bracket - STD | Offset bracket <sup>1</sup> |  | K0108 | \$550 |
| SG99                                                              | <input type="checkbox"/> | Omit side guard                                    |                        |                             |  |       | STD   |

1- only compatible with 0° camber; Increase width by 1/2" per side for a total of 1" ;

| POSITIONING BELT |                          |                                                                                     |  |  |  |       |       |
|------------------|--------------------------|-------------------------------------------------------------------------------------|--|--|--|-------|-------|
| SB01-S           | <input type="checkbox"/> | Velcro adjustable - <b>Small 35"</b> <sup>1</sup>                                   |  |  |  | E0978 | \$46  |
| SB01-L           | <input type="checkbox"/> | Velcro adjustable - <b>Long 60"</b>                                                 |  |  |  | E0978 | \$46  |
| SB02-S           | <input type="checkbox"/> | Auto buckle - <b>Small 35"</b> <sup>1</sup>                                         |  |  |  | E0978 | \$87  |
| SB02-M           | <input type="checkbox"/> | Auto buckle - <b>Medium 48"</b>                                                     |  |  |  | E0978 | \$87  |
| SB02-L           | <input type="checkbox"/> | Auto buckle - <b>Long 60"</b>                                                       |  |  |  | E0978 | \$87  |
| SB04-M           | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point not padded - <b>Medium 56"</b>                         |  |  |  | E0978 | \$85  |
| SB04-L           | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point not padded - <b>Long 61"</b>                           |  |  |  | E0978 | \$85  |
| SB05-S           | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point padded - <b>Small 51"</b>                              |  |  |  | E0978 | \$180 |
| SB05-M           | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point padded - <b>Medium 56"</b>                             |  |  |  | E0978 | \$180 |
| SB05-L           | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point padded - <b>Long 61"</b>                               |  |  |  | E0978 | \$180 |
| SB06-S           | <input type="checkbox"/> | BodyPoint Hip Belt - 4 point padded - <b>Small 51"</b>                              |  |  |  | E0978 | \$195 |
| SB06-M           | <input type="checkbox"/> | BodyPoint Hip Belt - 4 point padded - <b>Medium 56"</b>                             |  |  |  | E0978 | \$195 |
| SB06-L           | <input type="checkbox"/> | BodyPoint Hip Belt - 4 point padded - <b>Long 61"</b>                               |  |  |  | E0978 | \$195 |
| SB03             | <input type="checkbox"/> | Hardware attachment for 2 extra points (1 pair)                                     |  |  |  |       | \$17  |
| SB08-M           | <input type="checkbox"/> | BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Medium 18"</b> |  |  |  | K0108 | \$305 |
| SB08-L           | <input type="checkbox"/> | BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Long 24"</b>   |  |  |  | K0108 | \$305 |
| SB99             | <input type="checkbox"/> | Omit positioning belt                                                               |  |  |  |       | STD   |

1- Not recommended for seat width of 18" and more


| HEADREST <sup>1</sup> |                          |                                                                      |  |  |  |       |       |
|-----------------------|--------------------------|----------------------------------------------------------------------|--|--|--|-------|-------|
| AC01-S                | <input type="checkbox"/> | Profiled Headrest - <b>Small (3.5"x6")</b>                           |  |  |  | E0955 | \$185 |
| AC01-M                | <input type="checkbox"/> | Profiled Headrest - <b>Medium (4.5"x7")</b>                          |  |  |  | E0955 | \$185 |
| AC01-L                | <input type="checkbox"/> | Profiled Headrest - <b>Large (5.5"x8")</b>                           |  |  |  | E0955 | \$185 |
| AC28                  | <input type="checkbox"/> | Headrest mounting hardware with removable non folding stabilizer bar |  |  |  | E1028 | \$300 |

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

| ACCESSORIES |                          |                                                                                                                                                                                                                                                                                                                                |  |  |  |       |              |
|-------------|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|--|--|-------|--------------|
| WC10        | <input type="checkbox"/> | E-Motion adaptor Kit (2 adaptors) - <i>Only necessary with 6° camber</i>                                                                                                                                                                                                                                                       |  |  |  |       | \$36         |
| AC03        | <input type="checkbox"/> | Anti-tippers                                                                                                                                                                                                                                                                                                                   |  |  |  | E0971 | <b>\$150</b> |
| AC29        | <input type="checkbox"/> | Swing-away anti-tipper <sup>1</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>                                                                                                                                                                                                                               |  |  |  |       | \$190        |
| AC27        | <input type="checkbox"/> | Tip assist <sup>2</sup>                                                                                                                                                                                                                                                                                                        |  |  |  |       | \$67         |
| AC06        | <input type="checkbox"/> | Cane and crutch holder      Left <input type="checkbox"/> Right <input type="checkbox"/>                                                                                                                                                                                                                                       |  |  |  | E2207 | \$200        |
| AC07        | <input type="checkbox"/> | Oxygen tank holder <sup>4</sup> - D and E Type - installed on the left                                                                                                                                                                                                                                                         |  |  |  | E2208 | \$315        |
| AC08        | <input type="checkbox"/> | Spoke guards <sup>3</sup> Available: 20", 22", 24", 25", 26"                                                                                                                                                                                                                                                                   |  |  |  | K0065 | \$118        |
| AC09        | <input type="checkbox"/> | Footrest neoprene impact guard (pair) - Black                                                                                                                                                                                                                                                                                  |  |  |  | K0108 | \$92         |
| AC10        | <input type="checkbox"/> | Seat pouch                                                                                                                                                                                                                                                                                                                     |  |  |  |       | \$82         |
| AC30        | <input type="checkbox"/> | Newton multi-purpose pouch                                                                                                                                                                                                                                                                                                     |  |  |  |       | \$61         |
| AC26        | <input type="checkbox"/> | BodyPoint Mobility Bag                                                                                                                                                                                                                                                                                                         |  |  |  |       | \$135        |
| AC23        | <input type="checkbox"/> | Newton Backpack                                                                                                                                                                                                                                                                                                                |  |  |  |       | \$130        |
| AC12        | <input type="checkbox"/> | Velcro Closing Strap                                                                                                                                                                                                                                                                                                           |  |  |  |       | \$57         |
| AC13        | <input type="checkbox"/> | Tool kit                                                                                                                                                                                                                                                                                                                       |  |  |  |       | \$51         |
| AC15        | <input type="checkbox"/> | Military Patches      Army <input type="checkbox"/> Navy <input type="checkbox"/> Air Force <input type="checkbox"/> Marines <input type="checkbox"/> Coast Guard <input type="checkbox"/><br>Combat Wounded Purple Heart <input type="checkbox"/> Space Force <input type="checkbox"/> American Flag <input type="checkbox"/> |  |  |  |       | \$96         |

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate; 2- Cannot be combined with Anti Tippers


3- Not available with Mag wheels and Spinerigy; the 24" spoke guard will be installed on the 25" wheel; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down

| Colors               |                                                                                     |                                                                                     |                                                                                          |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                      |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |
|----------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Options <sup>1</sup> | Black                                                                               | High Gloss Carbon Fiber                                                             | Sunkissed Orange                                                                         | Burgundy                                                                            | Sapphire Blue                                                                       | Monster Green                                                                       | Charcoal                                                                            | Fuchsia                                                                             | Steel Blue                                                                          | Ferrari Red                                                                         | Acid Green                                                                          | Electric Purple                                                                      | Caribbean Blue                                                                        | White                                                                                 | Pink Pearl                                                                            | Plum Purple                                                                           | Sky Blue                                                                              | Black Cherry                                                                          | Khaki Green (matte)                                                                   | Desert Sand (matte)                                                                   |
|                      |  |  |       |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| AirGrip              | <input type="checkbox"/> \$110                                                      | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                                 | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                             | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              |
| Frame <sup>2</sup>   | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                                 | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                            | <input type="checkbox"/>                                                             | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              | <input type="checkbox"/>                                                              |
| FC99                 | <input type="checkbox"/>                                                            | Custom Color                                                                        | Specify PPG Paintit color: <span style="color: red;">Additional 10 days lead-time</span> |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                     |                                                                                      |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       |                                                                                       | \$390                                                                                 |

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- Black color is replaced by the Carbon fiber finish with the C2 upgrade;

## Fork and Stem bolt combinations



## Frame length chart

| Frame length | Seat Depth |     |     |     |     |     |     |
|--------------|------------|-----|-----|-----|-----|-----|-----|
|              | 14"        | 15" | 16" | 17" | 18" | 19" | 20" |
| STD          |            |     |     |     |     |     |     |
| +1           |            |     |     |     |     |     | N/A |
| +2           |            |     |     |     |     | N/A | N/A |
| -1           | N/A        |     |     |     |     |     |     |
| -2           | N/A        | N/A |     |     |     |     |     |

## Unilateral wheel lock restrictions

| Seat depth: shall be 16 in and greater                                |                 |                              |
|-----------------------------------------------------------------------|-----------------|------------------------------|
| Rear wheel: shall be 24 in and less                                   |                 |                              |
| "T" or flip-back armrest only compatible with 0° Camber               |                 |                              |
| One arm drive: Not compatible with "T" or flip-back armrest           |                 |                              |
| "T" or flip-back armrest may require to move the rear wheels outwards |                 |                              |
| Not compatible with wheel lock lever extension                        |                 |                              |
| Rear wheel diameter                                                   | Seat depth (in) | COG Compatibility Range (in) |
| 20 in                                                                 | 16+             | 1.5 to 3.75                  |
| 22 in                                                                 | 16              | 1.0 to 2.75                  |
|                                                                       | 17              | 1.0 to 3.25                  |
|                                                                       | 18+             | 1.0 to 3.75                  |
| 24 in                                                                 | 16              | 1.0 to 2.0                   |
|                                                                       | 17              | 1.0 to 2.75                  |
|                                                                       | 18+             | 1.0 to 3.75                  |



## CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the  
EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The **EVO Custom Request Form** is the best place to start.

Notes :




## FOLDING ON CHAIR – USA ORDER FORM

|                                                                                   |     |       |                                                                    |       |       |       |       |       |       |       |       |       |       |  |
|-----------------------------------------------------------------------------------|-----|-------|--------------------------------------------------------------------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|--|
| BioFit™ CUSHION                                                                   |     | DEPTH | WIDTH                                                              |       |       |       |       |       |       |       |       |       |       |  |
|  |     |       | 12"                                                                | 13"   | 14"   | 15"   | 16"   | 17"   | 18"   | 19"   | 20"   | 21"   | 22"   |  |
|                                                                                   | 12" |       | \$445                                                              | \$445 | \$445 |       | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   |  |
|                                                                                   | 13" |       |                                                                    |       |       |       |       | N/A   | N/A   | N/A   | N/A   | N/A   | N/A   |  |
|                                                                                   | 14" |       | \$445                                                              | \$445 | \$445 |       |       |       |       |       |       |       |       |  |
|                                                                                   | 15" |       |                                                                    | \$445 | \$445 | \$445 | \$445 |       |       |       |       |       |       |  |
|                                                                                   | 16" |       |                                                                    | \$445 | \$445 | \$445 | \$445 | \$445 | \$445 | \$480 | \$480 | \$480 | \$480 |  |
|                                                                                   | 17" |       | N/A                                                                |       |       | \$445 | \$445 | \$445 | \$445 | \$480 | \$480 |       |       |  |
|                                                                                   | 18" |       | N/A                                                                |       |       | \$445 | \$445 | \$445 | \$445 | \$480 | \$480 | \$480 | \$480 |  |
|                                                                                   | 19" |       | N/A                                                                |       |       | \$445 | \$445 | \$445 | \$445 | \$480 | \$480 |       |       |  |
|                                                                                   | 20" |       | N/A                                                                |       |       | \$445 | \$445 | \$445 | \$480 | \$480 | \$480 | \$480 | \$480 |  |
| Use wheelchair configuration<br>Or use this size:<br>NXCABF _____<br>(W" x D")    |     |       | Option: Change cover for incontinent/infection control outer cover |       |       |       |       |       |       |       |       |       |       |  |
| HCPCS CODE E2607-E2608                                                            |     |       |                                                                    |       |       |       |       |       |       |       |       |       |       |  |

Green boxes are stocked sizes; the red boxes require customization\*; the N/A boxes are not available sizes.

\*Customization (red boxes) is available with an upcharge of \$160.00 and has a 4 to 6 weeks additional lead time.

### SEATING AND CUSHION ACCESSORIES

|  | MODEL #   | Description       | Fits Back Width                                                             | Retail | QTY |
|-----------------------------------------------------------------------------------|-----------|-------------------|-----------------------------------------------------------------------------|--------|-----|
|                                                                                   | RSSI 1600 | Solid Seat Insert | Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm) | \$80   |     |
|                                                                                   | RSSI 1800 | Solid Seat Insert | Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm) | \$80   |     |
|                                                                                   | RSSI 2000 | Solid Seat Insert | Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)   | \$80   |     |

| <b>THORACIC SHAPE</b> |  |        | <b>Optima™</b>                            |           |                  |             |
|-----------------------|------------------------------------------------------------------------------------|--------|-------------------------------------------|-----------|------------------|-------------|
|                       |                                                                                    |        | THORACIC BACK SUPPORT                     |           |                  |             |
|                       |                                                                                    |        | CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb) |           |                  |             |
|                       |                                                                                    |        | HCPCS CODE: E2620-E2621                   |           |                  |             |
|                       | WIDTH                                                                              | HEIGHT | MODEL # CARBON FIBER                      | \$ CARBON | MODEL # ALUMINUM | \$ ALUMINUM |
|                       | 14"                                                                                | 13"    | NOCFACB1413                               | \$925     | NOACB1413        | \$695       |
|                       | 14"                                                                                | 16"    | NOCFACB1416                               | \$960     | NOACB1416        | \$695       |
|                       | 15"                                                                                | 13"    | NOCFACB1513                               | \$925     | NOACB1513        | \$695       |
|                       | 15"                                                                                | 16"    | NOCFACB1516                               | \$960     | NOACB1516        | \$695       |
|                       | 16"                                                                                | 13"    | NOCFACB1613                               | \$925     | NOACB1613        | \$695       |
|                       | 16"                                                                                | 16"    | NOCFACB1616                               | \$960     | NOACB1616        | \$695       |
|                       | 16"                                                                                | 19"    |                                           |           | NOACB1619        | \$695       |
|                       | 17"                                                                                | 13"    | NOCFACB1713                               | \$925     | NOACB1713        | \$695       |
|                       | 17"                                                                                | 16"    | NOCFACB1716                               | \$960     | NOACB1716        | \$695       |
|                       | 18"                                                                                | 13"    | NOCFACB1813                               | \$925     | NOACB1813        | \$695       |
|                       | 18"                                                                                | 16"    | NOCFACB1816                               | \$960     | NOACB1816        | \$695       |
|                       | 18"                                                                                | 19"    |                                           |           | NOACB1819        | \$695       |
|                       | 19"                                                                                | 16"    |                                           |           | NOACB1916        | \$695       |
|                       | 19"                                                                                | 19"    |                                           |           | NOACB1919        | \$695       |
|                       | 20"                                                                                | 13"    | NOCFACB2013                               | \$925     | NOACB2013        | \$695       |
|                       | 20"                                                                                | 16"    | NOCFACB2016                               | \$960     | NOACB2016        | \$695       |
|                       | 20"                                                                                | 19"    |                                           |           | NOACB2019        | \$695       |
|                       | 21"                                                                                | 16"    |                                           |           | NOACB2116        | \$745       |
|                       | 21"                                                                                | 19"    |                                           |           | NOACB2119        | \$745       |
|                       | 22"                                                                                | 13"    |                                           |           | NOACB2213        | \$745       |
|                       | 22"                                                                                | 16"    |                                           |           | NOACB2216        | \$745       |
|                       | 22"                                                                                | 19"    |                                           |           | NOACB2219        | \$745       |
|                       | ITEM                                                                               |        | WIDTH X DEPTH                             |           | RETAIL           | QTY         |
|                       | Replacement Cover & Foam                                                           |        | NOACC                                     |           | \$130            |             |

|  |        | <b>Optima™ DEEP</b>                       |           |                  |             |
|------------------------------------------------------------------------------------|--------|-------------------------------------------|-----------|------------------|-------------|
|                                                                                    |        | THORACIC BACK SUPPORT                     |           |                  |             |
|                                                                                    |        | CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb) |           |                  |             |
|                                                                                    |        | HCPCS CODE: E2620-E2621                   |           |                  |             |
| WIDTH                                                                              | HEIGHT | MODEL # CARBON FIBER                      | \$ CARBON | MODEL # ALUMINUM | \$ ALUMINUM |
| 14"                                                                                | 13"    | NOCFPCB1413                               | \$1,015   | NOPCB1413        | \$695       |
| 14"                                                                                | 16"    | NOCFPCB1416                               | \$1,055   | NOPCB1416        | \$695       |
| 15"                                                                                | 13"    | NOCFPCB1513                               | \$1,015   | NOPCB1513        | \$695       |
| 15"                                                                                | 16"    | NOCFPCB1516                               | \$1,055   | NOPCB1516        | \$695       |
| 16"                                                                                | 13"    | NOCFPCB1613                               | \$1,015   | NOPCB1613        | \$695       |
| 16"                                                                                | 16"    | NOCFPCB1616                               | \$1,055   | NOPCB1616        | \$695       |
| 16"                                                                                | 19"    |                                           |           | NOPCB1619        | \$695       |
| 17"                                                                                | 13"    | NOCFPCB1713                               | \$1,015   | NOPCB1713        | \$695       |
| 17"                                                                                | 16"    | NOCFPCB1716                               | \$1,055   | NOPCB1716        | \$695       |
| 18"                                                                                | 13"    | NOCFPCB1813                               | \$1,015   | NOPCB1813        | \$695       |
| 18"                                                                                | 16"    | NOCFPCB1816                               | \$1,055   | NOPCB1816        | \$695       |
| 18"                                                                                | 19"    |                                           |           | NOPCB1819        | \$695       |
| 19"                                                                                | 16"    |                                           |           | NOPCB1916        | \$695       |
| 19"                                                                                | 19"    |                                           |           | NOPCB1919        | \$695       |
| 20"                                                                                | 13"    | NOCFPCB2013                               | \$1,015   | NOPCB2013        | \$695       |
| 20"                                                                                | 16"    | NOCFPCB2016                               | \$1,055   | NOPCB2016        | \$695       |
| 20"                                                                                | 19"    |                                           |           | NOPCB2019        | \$695       |
| 21"                                                                                | 16"    |                                           |           | NOPCB2116        | \$745       |
| 21"                                                                                | 19"    |                                           |           | NOPCB2119        | \$745       |
| 22"                                                                                | 13"    |                                           |           | NOPCB2213        | \$745       |
| 22"                                                                                | 16"    |                                           |           | NOPCB2216        | \$745       |
| 22"                                                                                | 19"    |                                           |           | NOPCB2219        | \$745       |
| ITEM                                                                               |        | WIDTH X DEPTH                             |           | RETAIL           | QTY         |
| Replacement Cover & Foam                                                           |        | NOPCC                                     |           | \$130            |             |

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at [www.motioncomposites.com](http://www.motioncomposites.com)



**Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT**

| THORACIC SHAPE           |  | xtend™            | WIDTH | HEIGHT        | RETAIL | ✓      |
|--------------------------|-----------------------------------------------------------------------------------|-------------------|-------|---------------|--------|--------|
|                          |                                                                                   | NHAACB (L)1213/17 | 12"   | 13"-17"       | \$695  |        |
|                          |                                                                                   | NHAACB (L)1413/17 | 14"   | 13"-17"       | \$695  |        |
|                          |                                                                                   | NHAACB (L)1613/17 | 16"   | 13"-17"       | \$695  |        |
|                          |                                                                                   | NHAACB (L)1813/17 | 18"   | 13"-17"       | \$745  |        |
|                          |                                                                                   | NHAACB (L)2013/17 | 20"   | 13"-17"       | \$745  |        |
|                          |                                                                                   | NHAACB (L)2213/17 | 22"   | 13"-17"       | \$745  |        |
|                          | HCPCS<br>E2620-<br>E2601                                                          | ITEM              |       | WIDTH X DEPTH |        | RETAIL |
| Replacement Cover & Foam |                                                                                   | N3LHAACC          |       | \$130         |        |        |

**Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT**

| xtend™ Deep              | WIDTH         | HEIGHT  | RETAIL | ✓   |
|--------------------------|---------------|---------|--------|-----|
| NHAPCB1616/20            | 16"           | 16"-20" | \$695  |     |
| NHAPCB1816/20            | 18"           | 16"-20" | \$695  |     |
| NHAPCB2016/20            | 20"           | 16"-20" | \$695  |     |
| NHAPCB2216/20            | 22"           | 16"-20" | \$745  |     |
| HCPCS CODE E2620-E2621   |               |         |        |     |
| ITEM                     | WIDTH X DEPTH |         | RETAIL | QTY |
| Replacement Cover & Foam | N3HAPCC       |         | \$130  |     |

PLEASE NOTE: Each back will fit widths up to 2" wider.

**Classic™ PELVIC BACK SUPPORT**

| PELVIC SHAPE             | <br>HCPCS<br>CODE<br>E2620-<br>E2621 | Classic™      | WIDTH | HEIGHT | RETAIL | ✓ |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------|---------------|-------|--------|--------|---|
|                          |                                                                                                                       | NACB1616      | 16"   | 16"    | \$695  |   |
|                          |                                                                                                                       | NACB1618      | 16"   | 18"    | \$695  |   |
|                          |                                                                                                                       | NACB1816      | 18"   | 16"    | \$695  |   |
|                          |                                                                                                                       | NACB1818      | 18"   | 18"    | \$695  |   |
|                          |                                                                                                                       | NACB1820      | 18"   | 20"    | \$695  |   |
|                          |                                                                                                                       | NACB2020      | 20"   | 20"    | \$695  |   |
|                          |                                                                                                                       | NACB2216      | 22"   | 16"    | \$745  |   |
|                          |                                                                                                                       | NACB2218      | 22"   | 18"    | \$745  |   |
|                          |                                                                                                                       | NACB2220      | 22"   | 20"    | \$745  |   |
| ITEM                     |                                                                                                                       | WIDTH X DEPTH |       | RETAIL | QTY    |   |
| Replacement Cover & Foam |                                                                                                                       | NACC          |       | \$130  |        |   |

**Classic™ DEEP PELVIC BACK SUPPORT**

| Classic™ Deep            | WIDTH         | HEIGHT | RETAIL | ✓   |
|--------------------------|---------------|--------|--------|-----|
| NPCB1616                 | 16"           | 16"    | \$695  |     |
| NPCB1618                 | 16"           | 18"    | \$695  |     |
| NPCB1818                 | 18"           | 18"    | \$695  |     |
| NPCB1820                 | 18"           | 20"    | \$695  |     |
| NPCB2020                 | 18"           | 20"    | \$695  |     |
| NPCB2216                 | 20"           | 16"    | \$745  |     |
| NPCB2218                 | 22"           | 18"    | \$745  |     |
| NPCB2220                 | 22"           | 20"    | \$745  |     |
| HCPCS CODE E2620-E2621   |               |        |        |     |
| ITEM                     | WIDTH X DEPTH |        | RETAIL | QTY |
| Replacement Cover & Foam | NPCC          |        | \$130  |     |

PLEASE NOTE: Each back will fit widths up to 2" wider.

**BACK SUPPORT ACCESSORIES**

|  | MODEL #        | DESCRIPTION                                               | FITS BACK WIDTH                    | RETAIL | QTY |
|-------------------------------------------------------------------------------------|----------------|-----------------------------------------------------------|------------------------------------|--------|-----|
|                                                                                     | Custom - Cover | Optional Infection Control Wipe Down Back Support Cover   | Specify back model # when ordering | \$205  |     |
|                                                                                     | NIB-1938       | Extra Small Privacy Flap for 11" - 13" wide back supports | 11" - 13" (27.9 - 33.0)            | \$70   |     |
|                                                                                     | NIB-1078       | Small Privacy Flap for 14"-16" wide back supports         | 14" - 16" (35.6-40.6 cm)           | \$70   |     |
|                                                                                     | NIB-1079       | Medium Flap for 17"-19" wide back supports                | 17" - 19" (43.2-48.3 cm)           | \$70   |     |
|                                                                                     | NIB-1080       | Large Privacy Flap for 20"-22" wide back supports         | 20" - 22" (50.8-55.9 cm)           | \$70   |     |
|                                                                                     | NIB-1081       | Extra Large Privacy Flap for 23" - 26" wide back supports | 23" - 26" (58.4 - 66.0 cm)         | \$70   |     |

**MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS**

HCPCS CODE: E0955

| MODEL #  | Description                      | RETAIL | ✓ |
|----------|----------------------------------|--------|---|
| NIB-0784 | 10" x 5" Large Rectangular       | \$205  |   |
| NIB-0818 | 14" x 5" Extra-Large Rectangular | \$205  |   |

**MultiFit™ HEAD SUPPORT – FIXTURES**

HCPCS CODE HARDWARE: E1033

| MODEL #  | Description                   | RETAIL | ✓ |
|----------|-------------------------------|--------|---|
| NIB-0747 | Multi-Axis Hardware           | \$150  |   |
| NIB-0823 | Head Support Receiver Bracket | \$60   |   |

**EXTRA COVERS**

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



| ITEM                                                                     | WIDTH X DEPTH | RETAIL | QTY |
|--------------------------------------------------------------------------|---------------|--------|-----|
| Extra Outer smartx <sup>3D</sup> Cover – For a cover up to 18" x 18"     | NXCCBF V      | \$130  |     |
| Extra Outer smartx <sup>3D</sup> Cover – For a cover 18" x 20" and above | NXCCBF V      | \$150  |     |
| Extra Outer Incontinent Cover – For a cover up to 18" x 18"              | NXCCAF V      | \$130  |     |
| Extra Outer Incontinent Cover – For a cover 18" x 20" and above          | NXCCAF V      | \$150  |     |

ALL PRICES QUOTED IN USD. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at [www.motioncomposites.com](http://www.motioncomposites.com)