



1- Purchasing information

P.O. number: * _____ Date: _____
Company Name: * _____
Purchasing contact: * _____
Phone: * _____
Email: * _____
Contact preference: * Phone ☐ Email ☐
Marked for: _____

2- ATP/Clinician information

ATP/Clinician: * _____
Phone: * _____
Email: * _____
Contact preference: * Phone ☐ Email ☐

3- Shipping information

Ship to: * _____
Name: _____
Address: _____
City: _____ State: _____
Zip Code: _____
Phone: _____

* Required

MODEL		English Manual (STD)		Spanish Manual		Taxes not included	
						USD	
WC03	☐	Ultralight folding Helio wheelchair (265lb weight capacity) A7 package: side frames, symmetrical hydroformed crossbrace and back canes made of 7000 series aluminum				K0005	\$2,395
WC08	☐	C2 carbon fiber upgrade: side frames, symmetrical molded crossbrace and back canes made of ultralight carbon fiber					\$2,295
WC09	☐	HD Heavy Duty package (350 lb weight capacity)					\$550
Frame							
SEAT WIDTH ³							
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐	19" ¹ ☐
	☐		20" ¹ ☐	21" ² ☐	22" ² ☐		
SEAT DEPTH							
SD	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐
	☐		20" ☐				
FRAME LENGTH ⁵							
SP01	☐	Standard frame length configuration ⁶					STD
SP02	☐	Frame length 1" longer than standard configuration					NC
SP03	☐	Frame length 2" longer than standard configuration					NC
SP04	☐	Frame length 1" shorter than standard configuration					NC
TRANSIT TIE DOWN							
AC04	☐	Transit Tie Down (Unoccupied)					NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 20"-41"					\$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 20"-46"					\$365
AC35	☐	WC-19 Transit Tie Down without belt ⁴					\$350

1- Available with HD Kit (350 lb); 2- Only with HD Kit (350 lb)

4- May interfere with swing away armrest and Clamp on height adjustable push handle; not available with HD option when carbon upgrade is selected.

5- See table on page 7; 6- Distance of 2 1/8" from end of seat upholstery to end of frame;

[3. Width calculator](#)



NEWTON ULTRALIGHT PACKAGE			
PK01	☐	Includes: Carbon fiber angle adjustable footplates + Newton Gravity Ultralight rear wheels	Save \$275 \$530
VP01	☐	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹	K0108 \$95

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo ***VP01**; Recommended with black chair; See picture of the product on page 6.

Helio FRONT Seat-to-floor Height - Please select your choice

		3" caster			4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length			Stem bolt Length		
Forks / Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	☐ 13"	☐ 14"	☐ 15"	☐ 13 1/2"	☐ 14 1/2"	☐ 15 1/2"												
	P2	☐ 13 1/2"	☐ 14 1/2"	☐ 15 1/2"	☐ 14"	☐ 15"	☐ 16"	☐ 14 1/2"	☐ 15 1/2"	☐ 16 1/2"									
5"	P1																		
	P2				☐ 14"	☐ 15"	☐ 16"	☐ 14 1/2"	☐ 15 1/2"	☐ 16 1/2"									
	P3	☐ 14"	☐ 15"	☐ 16"	☐ 14 1/2"	☐ 15 1/2"	☐ 16 1/2"	☐ 15"	☐ 16"	☐ 17"	☐ 15 1/2"	☐ 16 1/2"	☐ 17 1/2"						
	P4	☐ 14 1/2"	☐ 15 1/2"	☐ 16 1/2"	☐ 15"	☐ 16"	☐ 17"	☐ 15 1/2"	☐ 16 1/2"	☐ 17 1/2"	☐ 16"	☐ 17"	☐ 18"	☐ 16 1/2"	☐ 17 1/2"	☐ 18 1/2"			
7"	P1										☐ 16 1/2"	☐ 17 1/2"	☐ 18 1/2"	☐ 17"	☐ 18"	☐ 19"			
	P2							☐ 16 1/2"	☐ 17 1/2"	☐ 18 1/2"	☐ 17"	☐ 18"	☐ 19"	☐ 17 1/2"	☐ 18 1/2"	☐ 19 1/2"			
	P3							☐ 17"	☐ 18"	☐ 19"	☐ 17 1/2"	☐ 18 1/2"	☐ 19 1/2"	☐ 18"	☐ 19"	☐ 20"	☐ 18 1/2"	☐ 19 1/2"	☐ 20 1/2"
	P4	☐ 16 1/2"	☐ 17 1/2"	☐ 18 1/2"	☐ 17"	☐ 18"	☐ 19"	☐ 17 1/2"	☐ 18 1/2"	☐ 19 1/2"	☐ 18"	☐ 19"	☐ 20"	☐ 18 1/2"	☐ 19 1/2"	☐ 20 1/2"	☐ 19"	☐ 20"	☐ 21"
	P5	☐ 17"	☐ 18"	☐ 19"	☐ 17 1/2"	☐ 18 1/2"	☐ 19 1/2"	☐ 18"	☐ 19"	☐ 20"	☐ 18 1/2"	☐ 19 1/2"	☐ 20 1/2"	☐ 19"	☐ 20"	☐ 21"	☐ 19 1/2"	☐ 20 1/2"	☐ 21 1/2"
Froglegs	P1				☐ 14 1/2"	☐ 15 1/2"	☐ 16 1/2"												
	P2	☐ 14 1/2"	☐ 15 1/2"	☐ 16 1/2"	☐ 15"	☐ 16"	☐ 17"	☐ 15 1/2"	☐ 16 1/2"	☐ 17 1/2"									
	P3	☐ 15"	☐ 16"	☐ 17"	☐ 15 1/2"	☐ 16 1/2"	☐ 17 1/2"	☐ 16"	☐ 17"	☐ 18"	☐ 16 1/2"	☐ 17 1/2"	☐ 18 1/2"						

Height available + - 1/4" (see stem bolt and fork combination at the end of the order form)

FRONT CASTERS

CASTERS				Available		
CA01-03 to CA01-08	☐	Composite wheel w/ PU tire (1")		3", 4", 5", 6", 7", 8"		NC
CA10 to CA11	☐	Pneumatic casters (1 1/4")		6", 8"	E2214/E2215	\$165
CA14 to CA16	☐	Newton UltraCaster Composite (1 1/2 ")		4", 5", 6"		\$95
CA17 to CA19	☐	Newton UltraCaster Aluminum (1 1/2 ")		4", 5", 6"		\$215
CASTER OPTIONS						
CA12	☐	Caster Pin Locks (Available with +1" or +2" Stem bolt ONLY) ¹			K0073	\$300
CA13	☐	Carbon fiber Frog Legs suspension forks (Not available with caster pin locks)			E1015	\$490

1- Not available with WC-19

Helio - REAR Seat-to-floor Height


WHEEL SIZES		Height available		Specify desired height: _____
☐	20" (451) wheel	12 1/2	to 17 1/4	
☐	22" (501) wheel	13 1/2	to 18 1/4	
☐	24" (540) wheel	14 1/2	to 19 1/4	
☐	25" (559) wheel	15	to 19 3/4	
☐	26" (590) wheel	15 1/2	to 20 1/4	

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

REAR WHEELS									
ONE ARM DRIVE ¹				Only available with Newton One rear wheels, permanent axles, 0° camber					
HR22	☐	One arm drive - aluminum inner handrim ²	Left ☐	Right ☐			E0958	\$1,495	
HR60	☐	One arm drive - plastic coated inner handrim ²	Left ☐	Right ☐			E0958	\$1,600	
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)						NC	
REAR WHEELS				Available					
RW01	☐	Mag ³				20", 22" 24"		\$57	
RW02	☐	Newton One - spoke wheel ⁴				20", 22" 24"		NC	
RW03	☐	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"		\$360	
RW04	☐	Spinergy Spox	Blk STD ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	Yel ⁵ ☐	\$915	
RW05	☐	Spinergy LX	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	\$1,145	
RW07	☐	Spinergy CLX ⁶	Blk STD ☐	Pur ⁵ ☐	Wht ⁵ ☐	Blue ⁵ ☐	Red ⁵ ☐	\$1,625	
TIRES (The brand and color may vary depending on parts availability)				Available					
AIRLESS									
TI01	☐	Soft Urethane 1 3/8" - Low Tread, black ⁹				20", 22" 24"	E2218	STD	
TI02	☐	Soft Urethane 1 3/8" - Medium Tread, black ^{7,9}				22", 24", 25", 26"	E2218	\$97	
TI13	☐	Soft Urethane 1" - No tread, black - Sentinel				20", 22", 24", 25", 26"	E2218	\$118	
TI08	☐	Pneumatic Airless insert 1 3/8" - Medium Tread ⁸				22", 24", 26"	E2211/E2213	\$95	
PNEUMATIC									
TI05	☐	Pneumatic 1 3/8" - Medium Tread ⁹				20", 22", 24", 26"	E2211/E2212	NC	
TI06	☐	Pneumatic HP 1" - Low Tread, Puncture resistant, black - SpeedLite				20", 22", 24", 25", 26"	E2211/E2212	\$105	
TI07	☐	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - TrailBlazer				24", 25", 26"	E2211/E2212	\$150	
TI03	☐	Pneumatic 2" - Large Tread, black ¹⁰				22", 24", 25"	E2211/E2212	\$115	
1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide; 2- Only standard handrim position available; 3- Not compatible with HP tires; Only available with aluminum or plastic coated handrim, standard position; 4- Not compatible with Shox tires; 5- color spokes require an additional 10 days lead-time; 6- Only available with 1" tires; 7- 22" is only available with Mags, Newton One rear wheels and Newton Gravity wheels; 8- Not available with Spinergy CLX, Spinergy LX 25", 26" and Spox 25", 26" rear wheels 9- Not available with flip-back armrest when the 6° camber is selected; 10- Not available with armrests when the 3° or 6° camber are selected;									
HANDRIMS									
HANDRIM				Available					
HR01	☐	Aluminum anodized				20", 22", 24", 25", 26"		STD	
HR61	☐	Aluminum black hard anodized ²				20", 22", 24", 25"		\$75	
HR02	☐	Plastic coated ¹				20", 22", 24", 25", 26"		\$110	
HR03	☐	High friction coated ^{1,2,3}				22", 24", 25", 26"		\$375	
HR07-HR18 and HR65-HR71	☐	Newton Air Grip ² See on page 6 for color choice and price				20", 22", 24", 25", 26"			
HR28	☐	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"		\$450	
HR05	☐	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"		\$495	
HR04	☐	Natural Fit ^{1,2,3,4}				22", 24", 25", 26"		\$395	
HR62	☐	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)				22", 24", 25", 26"		\$395	
HR63	☐	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)				22", 24", 25", 26"		\$445	
HR64	☐	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)				22", 24", 25", 26"		\$345	
HR99	☐	Omit handrim						NC	
1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive; 4 - 25" and 26" are not available in narrow position; 24" is not available with the Spinergy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece									
HANDRIMS OPTION									
ASSEMBLY POSITION									
HR29	☐	Standard Handrim position	Full tab					STD	
HR030	☐	Narrow Handrim position no cut ⁶	Full tab					NC	
HR30	☐	Narrow Handrim position with cut ⁵	Short tab 2 positions					\$30	
HR31	☐	Super Narrow Handrim ⁵	Short tab 2 positions STD ☐	Short tab 1 position ⁶ ☐				\$30	
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS									
HR23	☐	Standard grip						\$95	
HR27	☐	Super grip						\$150	
HR099	☐	Omit thumb piece						NC	
5- Only available with Aluminum and Air Grip handrims; 6- Not available with the combination of 24" Newton Gravity, Spinergy Spox, LX, wheels and a 1 3/8" Airless insert pneumatic tire with any handrim;									
AXLE ¹									
AX01	☐	Quick release axle						STD	
AX02	☐	Permanent axle						NC	
AX03	☐	Quad release axle					K0108	\$130	
AXLE PLATE									
AX04	☐	Amputee adjustable axle plate ²					E0959	\$205	
1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down and the swing away anti-tipper									

REAR WHEEL CAMBER			
CM00	☐	0° camber	STD
CM03	☐	3° camber	NC
CM06	☐	6° camber <i>(not available with amputee axle plate)</i>	\$165

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance X from front of back post to center of axle receiver.

X =

Most stable	1"	1/4"	1/2"	3/4"	2"	1/4"	1/2"	3/4"	3"	1/4"	1/2"	3/4"	4"
AMPUTEE	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

If not indicated, default setting is at 1.5"

Possible conflict between a seat depth under 16" and a center of gravity over 2,5".

Possible conflict between a 26" rear wheel and a center of gravity over 2,5".

Motion Composites will select the closest possible center of gravity (COG) if an issue arises with the COG

1-Only available with the amputee axle plate

Only on carbon frame

WHEEL LOCK			
WL08	☐	Newton Wheel lock - push to lock	STD
WL02	☐	Pull to lock	NC
WL04	☐	Newton grade aid push to lock <i>(works best with pneumatic or high pressure tire)</i>	E0974 \$155
WL13	☐	6" removable extension handle ¹ *VP01	E0961 \$77
WL14	☐		
WL05	☐	Scissors Lock	\$155
WL06	☐	Attendant Lock <i>(additional set of wheel lock mounted to the back)</i>	\$215
WL15	☐	Unilateral wheel lock - pull to lock ² Left ☐ Right ☐	\$330
WL99	☐	Omit wheel locks	NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

FOOTREST

SWING-IN/SWING-OUT FOOTREST			
FR03-01	☐	60° Swing-in/swing-out footrest - 13.5" - 17"	Set at*: NC
FR03-02	☐	60° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at*: NC
FR03-03	☐	60° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at*: NC
FR04-01	☐	70° Swing-in/swing-out footrest - 13.5" - 17"	Set at*: NC
FR04-02	☐	70° Swing-in/swing-out footrest - 14.5" - 19.5"	Set at*: NC
FR04-03	☐	70° Swing-in/swing-out footrest - 15.5" - 20.5"	Set at*: NC
FR15-01	☐	80° Swing-in/swing-out footrest - 13.5" - 16" ⁵	Set at*: NC
FR15-02	☐	80° Swing-in/swing-out footrest - 14.5" - 19.5" ⁵	Set at*: NC
FR15-03	☐	80° Swing-in/swing-out footrest - 15.5" - 20.5" ⁵	Set at*: NC
FR06	☐	90° swing-in/swing-out with front mount - Standard (no extension) ¹ 7.5" - 11"	Set at*: NC
FR07	☐	90° swing-in/swing-out with front mount - Short extension 2" ^{1, 2} 9.5" - 13"	Set at*: NC
FR08	☐	90° swing-in/swing-out with front mount - Medium extension 4" ^{1, 2} 9.5" - 14.5"	Set at*: NC
FOOTREST OPTIONS			
FR09	☐	High mount <i>(footplate attached to the hanger with clamps)</i> ³ 5" - 13.5"	Set at*: \$85
FR10	☐	Elevating legrest with calf pad, Length adjustable ⁴ Left ☐ Right ☐ 14" - 22"	Set at: E0990 \$245
FR99	☐	Omit Footrest	NC

1- Available only with two piece adjustable footplates; 2- May interfere with front casters; 3- Not available with front mount and the elevating legrest; * +/-1/2"

4- The elevating legrest replace the footrest side selected; Not available with one piece footplate; 5- Not available with 8" caster

FOOTPLATES			
FP01	☐	Standard - Composite ¹	STD
FP02	☐	Newton adjustable angle - Composite *VP01	K0040 \$215
FP03	☐	Newton adjustable angle - Aluminum	K0040 \$235
FP04	☐	Newton adjustable angle - Carbon fiber <i>(with rubber inserts)</i>	\$445
FP23	☐	Newton adjustable angle extra deep - Aluminum	K0040 \$235
FP06	☐	One piece adjustable angle flip-up - Aluminum	K0040 \$195
FP07	☐	One piece adjustable angle flip-up - Carbon fiber <i>(with rubber inserts)</i>	\$465

FOOTPLATE ACCESSORIES			
FP08	☐	Length adjustable heel loop ²	E0951 \$41
FP09	☐	Calf strap	K0038 \$39
FP30	☐	Bodypoint Aeromesh Padded calf strap	K0038 \$118

1- Not available with front mount; 2- not available with one piece footplate

SWING AWAY RESIDUAL LIMB SUPPORT			
FR27	☐	Residual limb support - flat 6"x8" Left ☐ Right ☐	E1020 \$420
FR28	☐	Residual limb support - flat 8"x10" Left ☐ Right ☐	E1020 \$420
FR29	☐	Residual limb support - flat 10"x10" Left ☐ Right ☐	E1020 \$420
FR30	☐	Residual limb support - flat 10"x14" Left ☐ Right ☐	E1020 \$420
FR32	☐	Residual limb support - contour 6"x8" Left ☐ Right ☐	E1020 \$420
FR33	☐	Residual limb support - contour 8"x10" Left ☐ Right ☐	E1020 \$420

SEATING AND BACKREST

BACK UPHOLSTERY									
BH07	☐	Slip-on Active-Mesh back upholstery							STD
BH08	☐	Adjustable tension Active-Mesh back upholstery					E2611/E2612		\$305
BH99	☐	Omit Back upholstery							NC
	☐	NXT back Support					See pricelist		
SEAT CUSHION									
SC01	☐	2" Cushion							\$170
SC02	☐	3" Cushion					E2601/E2602		\$190
SC99	☐	Omit seat cushion							STD
	☐	NXT cushion					See pricelist		
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLES - For Aluminum chair (Not available with push handles options)									
BC58	☐	Straight aluminum back cane	9 - 15" ☐	15 - 21" ☐					NC
BC60	☐	Straight aluminum back cane	21" - 24"		increment of 1"	Set at ² :			\$77
BC01	☐	8° bend aluminum back cane	16" ☐	18" ☐	20" ☐				NC
BC59	☐	8° bend aluminum back cane	21" - 24"		increment of 1"	Set at ² :			\$77
FIXED ANGLE, HEIGHT ADJUSTABLE BACK CANES WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)									
BC11	☐	Straight back cane	12 - 15" ☐	15 - 18" ☐	18 - 21" STD ☐	increment of 1/2"	Set at ² :		NC
BC06									
BC12	☐	8° Bend back cane	13 ½ - 15" ☐	15 - 18" ☐	18 - 21" STD ☐	increment of 1/2"	Set at ² :		\$60
BC08									
ANGLE AND HEIGHT ADJUSTABLE BACK CANE WITH PUSH HANDLE OPTIONS (MUST select ONE push handle options)									
BC13	☐	Straight back cane	15 - 18" ☐	18 - 21" STD ☐		increment of 1/2"	Set at ² :		\$215
BC09									
BC14	☐	8° bend back cane	15 - 18" ☐	18 - 21" STD ☐		increment of 1/2"	Set at ² :		\$270
BC10									
		Select back cane angle:	85° ☐	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐	
PUSH HANDLES OPTIONS (ONE option at a time)									
PH01	☐	Push handles							STD
PH05	☐	Fold down push handles ^{3, 5}							\$145
PH06-920	☐	Clamp on height adjustable push handles 9" - 20" ^{3, 4, 5}							\$325
PH06-923	☐	Clamp on height adjustable push handles 9" - 23" ^{3, 4, 5}							\$325
PH02	☐	Omit push handles							NC
STROLLER HANDLE AND STABILIZER BAR									
PH04	☐	Stroller handle (only available with push handle) - (not available with the HD option)					K0108		\$165
PH03	☐	Newton folding stabilizer bar (not compatible with stroller handle)							\$270

1- Not available with a 8° bend back cane; 2- If not indicated, default height setting is at the lowest position; measured from seat sling to top of back upholstery; 3- 265lb weight limit

4- Not compatible with swing away armrests, 8° bend back canes, 9-12" back canes and amputee axle plates; 5- Not compatible with stroller handle, stabilizer bar, headrest and attendant lock

ARMRESTS

HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ^{1, 4} *VP01									
AR01-812	☐	10" Desk length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		E0973		\$355
AR01-1014	☐	10" Desk length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		E0973		\$355
AR02-812	☐	14" Full length armpad	8" - 12"	Locked ☐	Unlocked STD ☐		E0973		\$355
AR02-1014	☐	14" Full length armpad	10" - 14"	Locked ☐	Unlocked STD ☐		E0973		\$355
HEIGHT ADJUSTABLE "T" ARMREST ¹ *VP01									
AR13-610	☐	10" Desk length armpad ⁵	6" - 10"				E0973		\$330
AR03-812	☐	10" Desk length armpad	8" - 12"				E0973		\$330
AR03-1014	☐	10" Desk length armpad	10" - 14"				E0973		\$330
AR14-610	☐	14" Full length armpad ⁵	6" - 10"				E0973		\$330
AR04-812	☐	14" Full length armpad	8" - 12"				E0973		\$330
AR04-1014	☐	14" Full length armpad	10" - 14"				E0973		\$330
HEIGHT ADJUSTABLE SWING AWAY ARMREST ²									
AR05	☐	Tubular armrest without plastic side guards - Short pad	8" - 11"				E0973		\$230
AR19	☐	Tubular armrest with plastic side guards - Short pad	8" - 11"				E0973		\$395
AR23	☐	Tubular armrest without plastic side guards - Long pad	8" - 11"				E0973		\$230
AR24	☐	Tubular armrest with plastic side guards - Long pad	8" - 11"				E0973		\$395
AR06	☐	10" Desk length flat armpad without plastic side guards	8" - 11"				E0973		\$230
AR20	☐	10" Desk length flat armpad with plastic side guards	8" - 11"				E0973		\$395
AR17	☐	14" Full length flat armpad without plastic side guards	8" - 11"				E0973		\$230
AR21	☐	14" Full length flat armpad with plastic side guards	8" - 11"				E0973		\$395
ARMREST OPTIONS									
AR12	☐	14" Full Length armpad ³ - Gel Ovation							\$165
AR99	☐	Omit armrest							NC

1- May have 1/2" adjustment restriction on both limits of the range; 2- Could interfere with some rigid backrests; Not available with attendant lock and fixed carbon fiber side guards;

3- Only available with full length armrest; Not available with swing away and desk length armrests; 4- Not available with attendant lock; 5- May have some restriction with seat-to-floor height and rear wheel

ACCESSORIES									
SIDE GUARDS <i>(not available with Flip Back and "T" Armrest)</i>									
SG01	☐	Plastic side guards <i>(removable)</i>	Straight bracket - STD ☐		Offset bracket ¹ ☐			K0108	\$200
SG02	☐	Carbon fiber side guards <i>(fixed)</i>							\$270
SG03	☐	Carbon fiber side guards <i>(removable)</i>	Straight bracket - STD ☐		Offset bracket ¹ ☐			K0108	\$305
SG06	☐	Plastic fender side guards <i>(removable)</i>	Straight bracket - STD ☐		Offset bracket ¹ ☐				\$290
SG04	☐	Carbon fiber fender side guards <i>(removable)</i>	Straight bracket - STD ☐		Offset bracket ¹ ☐			K0108	\$550
SG99	☐	Omit side guard							STD

1- only compatible with 0° camber; Increase width by 1/2" per side for a total of 1" ;

POSITIONING BELT									
SB01-S	☐	Velcro adjustable - Small 35" ¹						E0978	\$46
SB01-L	☐	Velcro adjustable - Long 60"						E0978	\$46
SB02-S	☐	Auto buckle - Small 35" ¹						E0978	\$87
SB02-M	☐	Auto buckle - Medium 48"						E0978	\$87
SB02-L	☐	Auto buckle - Long 60"						E0978	\$87
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"						E0978	\$85
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"						E0978	\$85
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"						E0978	\$180
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"						E0978	\$180
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"						E0978	\$180
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"						E0978	\$195
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"						E0978	\$195
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"						E0978	\$195
SB03	☐	Hardware attachment for 2 extra points (1 pair)							\$17
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"						K0108	\$305
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"						K0108	\$305
SB99	☐	Omit positioning belt							STD

1- Not recommended for seat width of 18" and more

HEADREST ¹									
AC01-S	☐	Profiled Headrest - Small (3.5"X6")						E0955	\$185
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")						E0955	\$185
AC01-L	☐	Profiled Headrest - Large (5.5"X8")						E0955	\$185
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar						E1028	\$300

1- Only available with push handles - Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES									
WC10	☐	E-Motion adaptor Kit (2 adaptors) - <i>Only necessary with 6° camber</i>							\$36
AC03	☐	Anti-tippers						E0971	\$150
AC29	☐	Swing-away anti-tipper ¹	Left ☐	Right ☐					\$190
AC27	☐	Tip assist ²							\$67
AC06	☐	Cane and crutch holder	Left ☐	Right ☐				E2207	\$200
AC07	☐	Oxygen tank holder ⁴ - D and E Type - installed on the left						E2208	\$315
AC08	☐	Spoke guards ³	Available: 20", 22", 24", 25", 26"					K0065	\$118
AC09	☐	Footrest neoprene impact guard (pair) - Black						K0108	\$92
AC10	☐	Seat pouch							\$82
AC30	☐	Newton multi-purpose pouch							\$61
AC26	☐	BodyPoint Mobility Bag							\$135
AC23	☐	Newton Backpack							\$130
AC12	☐	Velcro Closing Strap							\$57
AC13	☐	Tool kit							\$51
AC15	☐	Military Patches	Army ☐	Navy ☐	Air Force ☐	Marines ☐	Coast Guard ☐		\$96
		Combat Wounded Purple Heart ☐	Space Force ☐	American Flag ☐					

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate; 2- Cannot be combined with Anti Tipppers


3- Not available with Mag wheels and Spinerigy; the 24" spoke guard will be installed on the 25" wheel; 4- Not compatible with swing away armrest or amputee axle plate and WC-19 transit tie-down

Colors																				
Options ¹	Black ²	High Gloss Carbon Fiber ²	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
																				
AirGrip	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame (NC)	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
FC99	☐	Custom Color	Specify PPG Paintit color:																	Additional 10 days lead-time
																				\$390

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- Black frame will be STD for the aluminum/High Gloss Carbon Fiber will be STD for the Carbon frame upgrade; 3- The frame color will be confirmed with the quote conversion into an order (in absence of selection, the color will be black)

Fork and Stem bolt combinations



Frame length chart

	Seat Depth						
Frame length	14"	15"	16"	17"	18"	19"	20"
STD							
+1							N/A
+2						N/A	N/A
-1	N/A						

Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater

Rear wheel: shall be 24 in and less

"T" or flip-back armrest only compatible with 0° camber

One arm drive: not compatible with "T"

One arm drive and flip-back armrest only available with 17 in depth frame

"T" or flip-back armrest may require to move the rear wheels outwards

Not compatible with wheel lock lever extension

Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	1.5 to 3.75
	16	1.0 to 2.75
22 in	17+	1.0 to 3.25
	18+	1.0 to 3.75
24 in	16	1.0 to 2.0
	17	1.0 to 2.75
	18+	1.0 to 3.75

CUSTOMIZATION STARTS HERE!



With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the
EVO Custom Catalog.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The **EVO Custom Request Form** is the best place to start.

Notes :



FOLDING ON CHAIR – USA ORDER FORM



Use wheelchair configuration
Or use this size:
NXCAF
(W" x D")
HCPCS CODE E2607-E2608

DEPTH	WIDTH											
	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	
12"	\$468	\$468	\$468		N/A	N/A	N/A	N/A	N/A	N/A	N/A	
13"						N/A	N/A	N/A	N/A	N/A	N/A	
14"	\$468	\$468	\$468									
15"		\$468	\$468	\$468	\$468							
16"		\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	
17"	N/A			\$468	\$468	\$468	\$468	\$468	\$468			
18"	N/A			\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	
19"	N/A			\$468	\$468	\$468	\$468	\$468	\$468			
20"	N/A			\$468	\$468	\$468	\$468	\$468	\$468	\$468	\$504	

Option: Change cover for incontinent/infection control outer cover ☐ ☐ N/C

Green boxes are stocked sizes; the red boxes require customization*; the N/A boxes are not available sizes.
*Customization (red boxes) is available with an upcharge of \$160.00 and has a 2 to 4 weeks additional lead time.

SEATING AND CUSHION ACCESSORIES

	MODEL #	Description	Fits Back Width	Retail	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	\$80	



Optima™

THORACIC BACK SUPPORT

HCPCS CODE: E2620-E2621

THORACIC SHAPE

WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOACB1413	\$730	☐	NIB-2193V	\$130	☐
14"	16"	NOACB1416	\$730	☐	NIB-2194V	\$130	☐
15"	13"	NOACB1513	\$730	☐	NIB-2197V	\$130	☐
15"	16"	NOACB1516	\$730	☐	NIB-2198V	\$130	☐
16"	13"	NOACB1613	\$730	☐	NIB-2201V	\$130	☐
16"	16"	NOACB1616	\$730	☐	NIB-2202V	\$130	☐
16"	19"	NOACB1619	\$730	☐	NIB-2203V	\$130	☐
17"	13"	NOACB1713	\$730	☐	NIB-2206V	\$130	☐
17"	16"	NOACB1716	\$730	☐	NIB-2207V	\$130	☐
18"	13"	NOACB1813	\$730	☐	NIB-2208V	\$130	☐
18"	16"	NOACB1816	\$730	☐	NIB-2211V	\$130	☐
18"	19"	NOACB1819	\$730	☐	NIB-2212V	\$130	☐
19"	16"	NOACB1916	\$730	☐	NIB-2213V	\$130	☐
19"	19"	NOACB1919	\$730	☐	NIB-2215V	\$130	☐
20"	13"	NOACB2013	\$730	☐	NIB-2218V	\$130	☐
20"	16"	NOACB2016	\$730	☐	NIB-2219V	\$130	☐
20"	19"	NOACB2019	\$730	☐	NIB-2220V	\$130	☐
21"	16"	NOACB2116	\$730	☐	NIB-2221V	\$130	☐
21"	19"	NOACB2119	\$730	☐	NIB-2222V	\$130	☐
22"	13"	NOACB2213	\$783	☐	NIB-2223V	\$130	☐
22"	16"	NOACB2216	\$783	☐	NIB-2224V	\$130	☐
22"	19"	NOACB2219	\$783	☐	NIB-2225V	\$130	☐

PLEASE NOTE: Each back will fit widths up to 2" wider.



Optima Carbon™

THORACIC BACK SUPPORT
CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)

HCPCS CODE: E2620-E2621

WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOCFACB1413	\$972	☐	NIB-2193V	\$130	☐
14"	16"	NOCFACB1416	\$1,008	☐	NIB-2194V	\$130	☐
15"	13"	NOCFACB1513	\$972	☐	NIB-2197V	\$130	☐
15"	16"	NOCFACB1516	\$1,008	☐	NIB-2198V	\$130	☐
16"	13"	NOCFACB1613	\$972	☐	NIB-2201V	\$130	☐
16"	16"	NOCFACB1616	\$1,008	☐	NIB-2202V	\$130	☐
17"	13"	NOCFACB1713	\$972	☐	NIB-2206V	\$130	☐
17"	16"	NOCFACB1716	\$1,008	☐	NIB-2207V	\$130	☐
18"	13"	NOCFACB1813	\$972	☐	NIB-2211V	\$130	☐
18"	16"	NOCFACB1816	\$1,008	☐	NIB-2212V	\$130	☐
20"	13"	NOCFACB2013	\$972	☐	NIB-2218V	\$130	☐
20"	16"	NOCFACB2016	\$1,008	☐	NIB-2219V	\$130	☐

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

THORACIC SHAPE			Optima™ DEEP THORACIC BACK SUPPORT HCPCS CODE: E2620-E2621					
	WIDTH	HEIGHT	MODEL # ALUMINUM	RETAIL	✓	Repl. cover & foam	Retail	✓
	14"	13"	NOPCB1413	\$730	☐	NIB-2246V	\$130	☐
	14"	16"	NOPCB1416	\$730	☐	NIB-2247V	\$130	☐
	15"	13"	NOPCB1513	\$730	☐	NIB-2250V	\$130	☐
	15"	16"	NOPCB1516	\$730	☐	NIB-2251V	\$130	☐
	16"	13"	NOPCB1613	\$730	☐	NIB-2254V	\$130	☐
	16"	16"	NOPCB1616	\$730	☐	NIB-2255V	\$130	☐
	16"	19"	NOPCB1619	\$730	☐	NIB-2256V	\$130	☐
	17"	13"	NOPCB1713	\$730	☐	NIB-2259V	\$130	☐
	17"	16"	NOPCB1716	\$730	☐	NIB-2260V	\$130	☐
	18"	13"	NOPCB1813	\$730	☐	NIB-2264V	\$130	☐
	18"	16"	NOPCB1816	\$730	☐	NIB-2265V	\$130	☐
	18"	19"	NOPCB1819	\$730	☐	NIB-2266V	\$130	☐
	19"	16"	NOPCB1916	\$730	☐	NIB-2267V	\$130	☐
	19"	19"	NOPCB1919	\$730	☐	NIB-2268V	\$130	☐
	20"	13"	NOPCB2013	\$730	☐	NIB-2271V	\$130	☐
	20"	16"	NOPCB2016	\$730	☐	NIB-2272V	\$130	☐
	20"	19"	NOPCB2019	\$730	☐	NIB-2273V	\$130	☐
	21"	16"	NOPCB2116	\$730	☐	NIB-2274V	\$130	☐
	21"	19"	NOPCB2119	\$730	☐	NIB-2275V	\$130	☐
	22"	13"	NOPCB2213	\$783	☐	NIB-2276V	\$130	☐
	22"	16"	NOPCB2216	\$783	☐	NIB-2277V	\$130	☐
	22"	19"	NOPCB2219	\$783	☐	NIB-2278V	\$130	☐

		Optima Carbon™ DEEP THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb) HCPCS CODE: E2620-E2621					
WIDTH	HEIGHT	MODEL # CARBON	RETAIL	✓	Repl. cover & foam	Retail	✓
14"	13"	NOCFPCB1413	\$1,066	☐	NIB-2246V	\$130	☐
14"	16"	NOCFPCB1416	\$1,108	☐	NIB-2247V	\$130	☐
15"	13"	NOCFPCB1513	\$1,066	☐	NIB-2250V	\$130	☐
15"	16"	NOCFPCB1516	\$1,108	☐	NIB-2251V	\$130	☐
16"	13"	NOCFPCB1613	\$1,066	☐	NIB-2254V	\$130	☐
16"	16"	NOCFPCB1616	\$1,108	☐	NIB-2255V	\$130	☐
17"	13"	NOCFPCB1713	\$1,066	☐	NIB-2259V	\$130	☐
17"	16"	NOCFPCB1716	\$1,108	☐	NIB-2260V	\$130	☐
18"	13"	NOCFPCB1813	\$1,066	☐	NIB-2264V	\$130	☐
18"	16"	NOCFPCB1816	\$1,108	☐	NIB-2265V	\$130	☐
20"	13"	NOCFPCB2013	\$1,066	☐	NIB-2271V	\$130	☐
20"	16"	NOCFPCB2016	\$1,108	☐	NIB-2272V	\$130	☐

THORACIC SHAPE			Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT HCPCS CODE E2620-E2621					
	xtend™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
	NHAACB(L)1213/17	12"	13-17"	\$730	☐	NIB-2148V	\$130	☐
	NHAACB(L)1413/17	14"	13-17"	\$730	☐	NIB-2149V	\$130	☐
	NHAACB(L)1613/17	16"	13-17"	\$730	☐	NIB-2150V	\$130	☐
	NHAACB(L)1813/17	18"	13-17"	\$730	☐	NIB-2151V	\$130	☐
	NHAACB (L)2013/17	20"	13-17"	\$730	☐	NIB-2152V	\$130	☐
	NHAACB (L)2213/17	22"	13-17"	\$783	☐	NIB-2153V	\$130	☐

		Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT HCPCS CODE E2620-E2621					
xtend™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NHAPCB1616/20	16"	16"-20"	\$730	☐	NIB-2143V	\$130	☐
NHAPCB1816/20	18"	16"-20"	\$730	☐	NIB-2144V	\$130	☐
NHAPCB2016/20	20"	16"-20"	\$730	☐	NIB-2145V	\$130	☐
NHAPCB2216/20	22"	16"-20"	\$783	☐	NIB-2146V	\$130	☐

PELVIC SHAPE			Classic™ PELVIC BACK SUPPORT HCPCS CODE E2620-E2621					
	Classic™	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
	NACB1616	16"	16"	\$730	☐	NIB-2163V	\$130	☐
	NACB1618	16"	18"	\$730	☐	NIB-2164V	\$130	☐
	NACB1816	18"	16"	\$730	☐	NIB-2166V	\$130	☐
	NACB1818	18"	18"	\$730	☐	NIB-2167V	\$130	☐
	NACB1820	18"	20"	\$730	☐	NIB-2168V	\$130	☐
	NACB2020	20"	20"	\$730	☐	NIB-2171V	\$130	☐
	NACB2216	22"	16"	\$783	☐	NIB-2172V	\$130	☐
	NACB2218	22"	18"	\$783	☐	NIB-2173V	\$130	☐
	NACB2220	22"	20"	\$783	☐	NIB-2174V	\$130	☐

		Classic™ DEEP PELVIC BACK SUPPORT HCPCS CODE E2620-E2621					
Classic™ Deep	WIDTH	HEIGHT	RETAIL	✓	Repl. cover & foam	Retail	✓
NPCB1616	16"	16"	\$730	☐	NIB-2175V	\$130	☐
NPCB1618	16"	18"	\$730	☐	NIB-2176V	\$130	☐
NPCB1818	18"	18"	\$730	☐	NIB-2179V	\$130	☐
NPCB1820	18"	20"	\$730	☐	NIB-2180V	\$130	☐
NPCB2020	18"	20"	\$730	☐	NIB-2183V	\$130	☐
NPCB2216	22"	16"	\$783	☐	NIB-2184V	\$130	☐
NPCB2218	22"	18"	\$783	☐	NIB-2185V	\$130	☐
NPCB2220	22"	20"	\$783	☐	NIB-2186V	\$130	☐

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

BACK SUPPORT ACCESSORIES

	MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	\$205	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	\$70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	\$70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	\$70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	\$70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	\$70	

MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS

HCPCS CODE: E0955

	MODEL	Description	RETAIL	✓
	NIB-0784	10" x 5" Large Rectangular	\$205	
	NIB-0818	14" x 5" Extra-Large Rectangular	\$205	

MultiFit™ HEAD SUPPORT – FIXTURES

HCPCS CODE HARDWARE: E1033

	MODEL	Description	RETAIL	✓
	NIB-0747	Multi-Axis Hardware	\$150	
	NIB-0823	Head Support Receiver Bracket	\$60	

EXTRA COVERS

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx ^{3D} Cover – For a cover up to 18" x 18"	NXCCBF_____V	\$130	
Extra Outer smartx ^{3D} Cover – For a cover 18" x 20" and above	NXCCBF_____V	\$150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	NXCCAF_____V	\$130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	NXCCAF_____V	\$150	

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com