



MOVE

ORDER FORM - AADL #W114, W114HD, W114P

1- Purchasing information

P.O. number: * _____ Date: _____

Company Name: * _____

Purchasing contact: * _____

Phone: * _____

Email: * _____

Contact preference: * Phone ☐ Email ☐

Marked for: _____

3- Shipping information

Ship to: * _____

Name: _____

Address: _____

City: _____ Province: _____

Postal Code: _____

Phone: _____

2- Dealer/Representative information

Dealer/Representative: * _____

Phone: * _____

Email: * _____

Contact preference: * Phone ☐ Email ☐

* Required

MA080 ☐

English Manual (STD)

MA081 ☐

Manuel Français

Taxes not included

MODEL								
WC01	<input type="checkbox"/>	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)						\$1,425
WC09	<input type="checkbox"/>	Heavy duty kit (350 lb weight capacity) ⁵						\$1,070
Frame								
SEAT WIDTH ^{3,5}								
SW	<input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" ¹ <input type="checkbox"/>		NC
	<input type="checkbox"/>	19" ¹ <input type="checkbox"/>	20" ¹ <input type="checkbox"/>	21" ²	22" ² <input type="checkbox"/>			\$315
SEAT DEPTH ⁶								
SD	<input type="checkbox"/>	16" <input type="checkbox"/>		17" <input type="checkbox"/>	18" <input type="checkbox"/>			NC
	<input type="checkbox"/>	20" <input type="checkbox"/>						\$235
TRANSIT TIE DOWN								
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)						NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 20"-41"				Client cost		\$365
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 20"-46"				Client cost		\$365
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt ⁴				Client cost		\$345

1- Available with HD option; 2- Only available with HD option

3- Width calculator

4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

5- Seat width of 20" to 22" are no charge with HD option; 6- Seat depth of 20" is no charge with HD option

PACKAGES			
VP01	<input type="checkbox"/>	Vision Package Kit: Orange accent added to point of contact accessories including: footrest latches, "T" and "U" armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates ¹	\$95

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; Recommended with black chair; See picture of the product on page 5.

Move FRONT Seat-to-floor Height - Please select your choice

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
4"	P1	12 1/2"	13 1/2"	14 1/2"												
	P2	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
5"	P1															
	P2	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
	P3	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"						
	P4	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"			
7"	P1							15 1/2"	16 1/2"	17 1/2"	16"	17"	18"			
	P2				15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
	P3				16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"
	P4	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"
	P5	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS				Available		
CA01-04 to CA10-08	<input type="checkbox"/>	Composite wheel w/ PU tire (1")		4", 5", 6", 7", 8"		NC
CA10 to CA11	<input type="checkbox"/>	Pneumatic casters (1 ¼")		6", 8"		\$170
CA14 to CA16	<input type="checkbox"/>	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		\$105
CA17 to CA19	<input type="checkbox"/>	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"	A1*	\$198

**A1 is clients' responsibility for costs/A2 requires clinical justification*

Move REAR Seat-to-floor Height - Please circle your choice

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12"	13"	14"	14 1/2"	15"
13"	14"	15"	15 1/2"	16"
14"	15"	16"	16 1/2"	17"
15"	16"	17"	17 1/2"	18"
16"	17"	18"	18 1/2"	19"
17"	18"	19"	19 1/2"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest						
SWING-AWAY FOOTREST						
FR01-01	<input type="checkbox"/>	60° swing away footrest	13.5 - 17.5"	Set at:*		\$270
FR01-02	<input type="checkbox"/>	60° swing away footrest	15.5 - 19.5"	Set at:*		\$270
FR01-03	<input type="checkbox"/>	60° swing away footrest	17 - 21"	Set at:*		\$270
FR02-01	<input type="checkbox"/>	70° swing away footrest	13.5 - 17.5"	Set at:*		\$270
FR02-02	<input type="checkbox"/>	70° swing away footrest	15.5 - 19.5"	Set at:*		\$270
FR02-03	<input type="checkbox"/>	70° swing away footrest	17 - 21"	Set at:*		\$270
FOOTREST OPTIONS						
FR09	<input type="checkbox"/>	High mount <i>(footplate attached to the hanger with clamps)</i>	5" - 13.5"	Set at:*		\$105
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, Length adjustable ¹	Left <input type="checkbox"/> Right <input type="checkbox"/> 14" - 22"	Set at:		\$395
FR99	<input type="checkbox"/>	Omit Footrest				NC

1- The elevating legrest replace the footrest side selected; Not available with high mount

** +/- 1/2"*

FOOTPLATES						
FP01	<input type="checkbox"/>	Standard - Composite				STD
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite			W107	\$155
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum			W107	\$175
FOOTPLATE ACCESSORIES						
FP08	<input type="checkbox"/>	Length adjustable heel loop				\$26
FP09	<input type="checkbox"/>	Calf strap				\$57
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap			<i>basic seating benefit</i>	\$140

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	<input type="checkbox"/>	Residual limb support - flat 6"x8"	Left <input type="checkbox"/> Right <input type="checkbox"/>			\$505
FR28	<input type="checkbox"/>	Residual limb support - flat 8"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>			\$505
FR29	<input type="checkbox"/>	Residual limb support - flat 10"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>			\$505
FR30	<input type="checkbox"/>	Residual limb support - flat 10"x14"	Left <input type="checkbox"/> Right <input type="checkbox"/>			\$505
FR32	<input type="checkbox"/>	Residual limb support - contour 6"x8"	Left <input type="checkbox"/> Right <input type="checkbox"/>			\$505
FR33	<input type="checkbox"/>	Residual limb support - contour 8"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>			\$505

Backrest						
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE						
BC01-16	<input type="checkbox"/>	8° bend aluminum back cane	16"			NC
BC01-18	<input type="checkbox"/>	8° bend aluminum back cane	18"			NC
BC01-20	<input type="checkbox"/>	8° bend aluminum back cane	20"			NC
BC59	<input type="checkbox"/>	8° bend aluminum back cane	21" - 24"	<i>increment of 1"</i>	Set at ² :	\$100
ADJUSTABLE ANGLE BACK WITH INTEGRATED PUSH HANDLE						
BC02-16	<input type="checkbox"/>	8° bend aluminum back cane 16"				\$209
BC02-18	<input type="checkbox"/>	8° bend aluminum back cane 18"				\$209
BC02-20	<input type="checkbox"/>	8° bend aluminum back cane 20"				\$209
	Select back cane angle:		90° <input type="checkbox"/> 95° STD <input type="checkbox"/> 100° <input type="checkbox"/> 105° <input type="checkbox"/> 110° <input type="checkbox"/>			

STROLLER HANDLE AND STABILIZER BAR ¹						
PH04	<input type="checkbox"/>	Stroller handle <i>(not available with the HD option)</i>			A1*	\$220
PH03	<input type="checkbox"/>	Newton folding stabilizer bar <i>(not available with stroller handle)</i>				\$275

1- Stroller handle and stabilizer bar can't be combined

2- If not indicated, will be at the lowest position

BACK UPHOLSTERY						
BH05	<input type="checkbox"/>	Slip-on Nylon back upholstery				NC
BH06	<input type="checkbox"/>	Adjustable tension Nylon back upholstery			<i>Basic seating benefit</i>	\$260
BH99	<input type="checkbox"/>	Omit back upholstery				STD
SEAT CUSHIONS						
SC01	<input type="checkbox"/>	2" Cushion			<i>Basic seating benefit</i>	\$185
SC02	<input type="checkbox"/>	3" Cushion			<i>Basic seating benefit</i>	\$205
SC99	<input type="checkbox"/>	Omit seat cushion				NC

**A1 is clients' responsibility for costs/A2 requires clinical justification*

Rear wheels									
ONE ARM DRIVE ¹						Only available with Newton One rear wheels, permanent axles, 0° camber			
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim ²	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		W104		\$1,480
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim ²	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>		W104		\$1,595
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)							NC
REAR WHEELS						Available			
RW01	<input type="checkbox"/>	Mag ³				20", 22" 24"			NC
RW02	<input type="checkbox"/>	Newton One - spoke wheel ⁴				20", 22" 24"			NC
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"	A1*		\$390
TIRES						Available			
AIRLESS									
TI01	<input type="checkbox"/>	Soft Urethane 1 3/8" - Low Tread, black				20", 22" 24"			STD
TI02	<input type="checkbox"/>	Soft Urethane 1 3/8 " - Medium Tread, black ⁵				22", 24", 25", 26"	A1*		\$140
TI04	<input type="checkbox"/>	Soft Urethane 1" - No tread, black - <i>Shox</i>				22", 24", 25", 26"	A1*		\$165
TI08	<input type="checkbox"/>	Pneumatic Airless insert 1 3/8" - Medium Tread				22", 24", 26"	A1*		\$225
PNEUMATIC									
TI05	<input type="checkbox"/>	Pneumatic 1 3/8" - Medium Tread				20", 22", 24", 25", 26"			\$140
TI06	<input type="checkbox"/>	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>				20", 22", 24", 25", 26"			\$160
TI07	<input type="checkbox"/>	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>				24", 25", 26"			\$240

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available;
3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position; 4- Not compatible with Shox tires; 5- The 22" is only available with Mags and Newton One rear wheels


Handrims							
HANDRIM						Available	
HR01	<input type="checkbox"/>	Aluminum anodized				20", 22", 24", 25", 26"	STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized ²				20", 22", 24", 25"	\$85
HR02	<input type="checkbox"/>	Plastic coated ¹				20", 22", 24", 25", 26"	\$170
HR03	<input type="checkbox"/>	High friction coated ^{1,2,3}				22", 24", 25", 26"	A1* \$460
HR07-HR18 and HR65-HR71	<input type="checkbox"/>	Newton Air Grip ² <i>See page 5 for the color selection and price</i>				20", 22", 24", 25", 26"	
HR28	<input type="checkbox"/>	Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"	A1* \$485
HR05	<input type="checkbox"/>	Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"	A1* \$515
HR04	<input type="checkbox"/>	Natural Fit ^{1,2,3,4}				22", 24", 25", 26"	A1* \$320
HR62	<input type="checkbox"/>	Simi H ^{1,2,3} (1 3/8" Oval handrim with high friction strip)				22", 24", 25", 26"	A1* \$485
HR63	<input type="checkbox"/>	Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip)				22", 24", 25", 26"	A1* \$515
HR64	<input type="checkbox"/>	Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim)				22", 24", 25", 26"	A1* \$320
HR99	<input type="checkbox"/>	Omit handrim					NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;
4- 25" and 26" are not available in narrow position; 24" is not available with the Spinygro Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

Handrims option							
ASSEMBLY POSITION							
HR29	<input type="checkbox"/>	Standard Handrim position	Full tab				STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	Full tab				NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut ⁵	Short tab 2 positions				\$29
HR31	<input type="checkbox"/>	Super Narrow Handrim ⁵	Short tab 2 positions STD <input type="checkbox"/> Short tab 1 position <input type="checkbox"/>				\$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS							
HR23	<input type="checkbox"/>	Standard grip				A1*	\$200
HR27	<input type="checkbox"/>	Super grip				A1*	\$230
HR099	<input type="checkbox"/>	Omit thumb piece					NC

AXLES (0° camber) ¹							
AX01	<input type="checkbox"/>	Quick release axle					STD
AX02	<input type="checkbox"/>	Permanent axle					NC
AX03	<input type="checkbox"/>	Quad release axle				A2*	\$180
AXLE PLATE							
AX04	<input type="checkbox"/>	Amputee adjustable axle plate ²				A2*	\$275

1- Not required with one arm drive; 2- May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down
*A1 is clients' responsibility for costs/A2 requires clinical justification



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable

3/4"

Most responsive

AMPUTEE ¹

2"

☐
☐
☐

1- Only available with amputee axle plate

WHEEL LOCKS					
WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock	PRODUCT IMPROVED		STD
WL16	<input type="checkbox"/>	Aluminum push to lock with extensions			\$73
WL02	<input type="checkbox"/>	Pull to lock			NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock (<i>works best with pneumatic or high pressure tire</i>)		A2*	\$185
WL13	<input type="checkbox"/>	6" removable extension handle ¹			\$73
WL14					
WL06	<input type="checkbox"/>	Attendant Lock (<i>additional set of wheel lock mounted to the back</i>)		A1*	\$240
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock ² Left <input type="checkbox"/> Right <input type="checkbox"/>		A1*	\$335
WL99	<input type="checkbox"/>	Omit wheel locks			NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:
https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests					
HEIGHT ADJUSTABLE FLIP BACK ARMREST (<i>Convertible to "T" armrest</i>) ^{1,2}					
AR01-812	<input type="checkbox"/>	10" Desk length armpad 8" - 12"	Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>		\$395
AR01-1014	<input type="checkbox"/>	10" Desk length armpad 10" - 14"	Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>		\$395
AR02-812	<input type="checkbox"/>	14" Full length armpad 8" - 12"	Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>		\$395
AR02-1014	<input type="checkbox"/>	14" Full length armpad 10" - 14"	Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/>		\$395
HEIGHT ADJUSTABLE "T" ARMREST ¹					
AR13-610	<input type="checkbox"/>	10" Desk length armpad ⁴ 6" - 10"			\$380
AR03-812	<input type="checkbox"/>	10" Desk length armpad 8" - 12"			\$380
AR03-1014	<input type="checkbox"/>	10" Desk length armpad 10" - 14"			\$380
AR14-610	<input type="checkbox"/>	14" Full length armpad ⁴ 6" - 10"			\$380
AR04-812	<input type="checkbox"/>	14" Full length armpad 8" - 12"			\$380
AR04-1014	<input type="checkbox"/>	14" Full length armpad 10" - 14"			\$380
ARMREST OPTIONS					
AR12	<input type="checkbox"/>	14" Full Length armpad - Gel Ovation (pair) ³			\$160
AR99	<input type="checkbox"/>	Omit armrest			NC

1- May have 1/2" adjust. restriction on both limits of the range; 2- Not available with attendant lock; 3- Only available with full length armrest; Not available with swing away and desk length armrests

4- May have restriction with seat-to-floor height and rear wheel

POSITIONING BELTS					
SB01-S	<input type="checkbox"/>	Velcro adjustable - Small 35" ¹	Basic seating benefit		\$50
SB01-L	<input type="checkbox"/>	Velcro adjustable - Long 60"	Basic seating benefit		\$50
SB02-S	<input type="checkbox"/>	Auto buckle - Small 35" ¹			\$22
SB02-M		Auto buckle - Medium 48"			\$22
SB02-L	<input type="checkbox"/>	Auto buckle - Long 60"			\$22
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Medium 56"	Basic seating benefit		\$95
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - Long 61"	Basic seating benefit		\$95
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Small 51"	Basic seating benefit		\$185
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Medium 56"	Basic seating benefit		\$185
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - Long 61"	Basic seating benefit		\$185
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Small 51"	Basic seating benefit		\$210
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Medium 56"	Basic seating benefit		\$210
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - Long 61"	Basic seating benefit		\$210
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)	Basic seating benefit		\$15
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"	Basic seating benefit		\$355
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"	Basic seating benefit		\$355
SB99	<input type="checkbox"/>	Omit positioning belt			STD

1- Not recommended for seat width of 18" and more

AC01-S	<input type="checkbox"/>	Profiled Headrest - Small (3.5"X6")	basic seating benefit		\$215
AC01-M	<input type="checkbox"/>	Profiled Headrest - Medium (4.5"X7")	Basic seating benefit		\$215
AC01-L	<input type="checkbox"/>	Profiled Headrest - Large (5.5"X8")	Basic seating benefit		\$215
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar ¹	Basic seating benefit		\$360

1- Not available with a stroller handle or a folding stabilizer bar

AC03	<input type="checkbox"/>	Anti-tippers			\$120
AC29	<input type="checkbox"/>	Swing-away anti-tipper ¹ Left <input type="checkbox"/> Right <input type="checkbox"/>			\$185
AC27	<input type="checkbox"/>	Tip assist ²			\$65
AC06	<input type="checkbox"/>	Cane and crutch holder Left <input type="checkbox"/> Right <input type="checkbox"/>			\$115
AC07	<input type="checkbox"/>	Oxygen tank holder ⁴ - D and E Type - Installed to the left			\$398
AC08	<input type="checkbox"/>	Spoke guards ³ Available: 20", 22", 24", 26"			\$155
AC09	<input type="checkbox"/>	Footrest neoprene impact guard (<i>pair</i>) - Black			\$85
AC10	<input type="checkbox"/>	Seat pouch			\$100
AC30	<input type="checkbox"/>	Newton multi-purpose pouch			\$75
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag			\$145
AC23	<input type="checkbox"/>	Newton Backpack			\$130
AC12	<input type="checkbox"/>	Velcro Closing Strap			\$55
AC13	<input type="checkbox"/>	Tool kit			\$55
AC16	<input type="checkbox"/>	Motion Composites Polo Shirt W/S <input type="checkbox"/> W/M <input type="checkbox"/> W/L <input type="checkbox"/> M/S <input type="checkbox"/> M/M <input type="checkbox"/> M/L <input type="checkbox"/> M/XL <input type="checkbox"/> M/XXL <input type="checkbox"/> M/XXXL <input type="checkbox"/>			\$36

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"


2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinerig; 4- Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-

Colors																			
Options ¹	Black*	Sunkissed Orange	Burgundy	Sapphire Blue*	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
AirGrip	<input type="checkbox"/> \$160	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Frame ²	<input type="checkbox"/> STD																		\$50
FC99	<input type="checkbox"/> Custom Color	Specify PPG Paintit color: additional 10 days lead-time																	\$422

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;
2- with black crossbrace;

Fork and Stem bolt combinations



Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	2
22 in	16	0.75 to 2
	17+	0.75 to 2
24 in	16	0.75
	17	0.75 to 2
	18+	0.75 to 2

Notes :