



**MOVE**

ORDER FORM - AADL #W114, W114HD, W114P

**1- Purchasing information**

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_  
 Company Name:\* \_\_\_\_\_  
 Purchasing contact:\* \_\_\_\_\_  
 Phone:\* \_\_\_\_\_  
 Email:\* \_\_\_\_\_  
 Contact preference:\* Phone  Email   
 Marked for: \_\_\_\_\_

**3- Shipping information**

Ship to:\* \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ Province: \_\_\_\_\_  
 Postal Code: \_\_\_\_\_  
 Phone: \_\_\_\_\_

**2- Dealer/Representative information**

Dealer/Representative:\* \_\_\_\_\_  
 Phone:\* \_\_\_\_\_  
 Email:\* \_\_\_\_\_  
 Contact preference:\* Phone  Email

MA080  English Manual (STD)      MA081  Manuel Français      Taxes not included

MODEL							
WC01	<input type="checkbox"/>	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)					\$1,425
WC09	<input type="checkbox"/>	Heavy duty kit (350 lb weight capacity) <sup>5</sup>					\$1,070
<b>Frame</b>							
SEAT WIDTH <sup>3,5</sup>							
SW	<input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <sup>1</sup> <input type="checkbox"/>	NC
	<input type="checkbox"/>	19" <sup>1</sup> <input type="checkbox"/>	20" <sup>1</sup> <input type="checkbox"/>	21" <sup>2</sup> <input type="checkbox"/>	22" <sup>2</sup> <input type="checkbox"/>		\$315
SEAT DEPTH <sup>6</sup>							
SD	<input type="checkbox"/>	16" <input type="checkbox"/>		17" <input type="checkbox"/>	18" <input type="checkbox"/>		NC
	<input type="checkbox"/>	20" <input type="checkbox"/>					\$235
TRANSIT TIE DOWN							
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)					NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup> Short belt 20"-41"			Client cost		\$365
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup> Long belt 20"-46"			Client cost		\$365
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt <sup>4</sup>			Client cost		\$345

1- Available with HD option; 2- Only available with HD option      3- Width calculator      4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option  
 5- Seat width of 20" to 22" are no charge with HD option; 6- Seat depth of 20" is no charge with HD option

PACKAGES							
VP01	<input type="checkbox"/>	Vision Package Kit: Orange accent added to point of contact accessories including: footrest latches, "T" and "U" armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates <sup>1</sup>					\$95

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; Recommended with black chair; See picture of the product on page 5.

**Move FRONT Seat-to-floor Height - Please select your choice**

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"												
4"	P1	12 1/2"	13 1/2"	14 1/2"												
	P2	13"	14"	15"		14 1/2"	15 1/2"									
5"	P1															
	P2	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
	P3	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"		15 1/2"	16 1/2"						
	P4	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"		16 1/2"	17 1/2"			
7"	P1							15 1/2"	16 1/2"	17 1/2"	16"	17"	18"			
	P2				15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
	P3				16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"
	P4	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"
	P5	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form); Orange boxes are not available with a seat slope of more than 1.5in.

**FRONT CASTERS**

CASTERS				Available		
CA01-04 to CA01-08	<input type="checkbox"/>	Composite wheel w/ PU tire (1")		4", 5", 6", 7", 8"		NC
CA10 to CA11	<input type="checkbox"/>	Pneumatic casters (1 ¼")		6", 8"		\$170
CA14 to CA16	<input type="checkbox"/>	Newton UltraCaster Composite (1 ½ ")		4", 5", 6"		\$105
CA17 to CA19	<input type="checkbox"/>	Newton UltraCaster Aluminum (1 ½ ")		4", 5", 6"	<b>A1*</b>	\$198

*\*A1 is clients' responsibility for costs/A2 requires clinical justification*

**Move REAR Seat-to-floor Height - Please circle your choice**

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12"	13"	14"	14 1/2"	15"
13"	14"	15"	15 1/2"	16"
14"	15"	16"	16 1/2"	17"
15"	16"	17"	17 1/2"	18"
16"	17"	18"	18 1/2"	19"
17"	18"	19"	19 1/2"	20"

*Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"*

**REAR WHEELS**

ONE ARM DRIVE <sup>1</sup>				Available			
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim <sup>2</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	<b>W104</b>	\$1,480
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim <sup>2</sup>	24"	Left <input type="checkbox"/>	Right <input type="checkbox"/>	<b>W104</b>	\$1,595
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)					NC
REAR WHEELS				Available			
RW01	<input type="checkbox"/>	Mag <sup>3</sup>				20", 22" 24"	NC
RW02	<input type="checkbox"/>	Newton One - spoke wheel <sup>4</sup>				20", 22" 24"	NC
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"	<b>A1*</b> \$390
TIRES				Available			
AIRLESS							
TI01	<input type="checkbox"/>	Soft Urethane 1 3/8" - Low Tread, black				20", 22" 24"	STD
TI02	<input type="checkbox"/>	Soft Urethane 1 3/8" - Medium Tread, black <sup>5</sup>				22", 24", 25", 26"	<b>A1*</b> \$140
TI04	<input type="checkbox"/>	Soft Urethane 1" - No tread, black - <i>Shox</i>				22", 24", 25", 26"	<b>A1*</b> \$165
TI08	<input type="checkbox"/>	Pneumatic Airless insert 1 3/8" - Medium Tread				22", 24", 26"	<b>A1*</b> \$225
PNEUMATIC							
TI05	<input type="checkbox"/>	Pneumatic 1 3/8" - Medium Tread				20", 22", 24", 25", 26"	\$140
TI06	<input type="checkbox"/>	Pneumatic HP 1" - Low Tread, Puncture resistant, black - <i>SpeedLite</i>				20", 22", 24", 25", 26"	\$160
TI07	<input type="checkbox"/>	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - <i>TrailBlazer</i>				24", 25", 26"	\$240

*1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; Not available with 21" and 22" wide*

*3- Not compatible with HP tires; Only available with aluminum or plastic coated handrim in standard position; 4- Not compatible with Shox tires; 5- The 22" is only available with Mags and Newton One rear wheels*

**HANDRIMS**

HANDRIM				Available			
HR01	<input type="checkbox"/>	Aluminum anodized				20", 22", 24", 25", 26"	STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized <sup>2</sup>				20", 22", 24", 25"	\$85
HR02	<input type="checkbox"/>	Plastic coated <sup>1</sup>				20", 22", 24", 25", 26"	\$170
HR03	<input type="checkbox"/>	High friction coated <sup>1,2,3</sup>				22", 24", 25", 26"	<b>A1*</b> \$460
HR07-HR18 and HR65-HR71	<input type="checkbox"/>	Newton Air Grip <sup>2</sup>	<i>See page 5 for the color selection and price</i>			20", 22", 24", 25", 26"	
HR28	<input type="checkbox"/>	Surge <sup>1,2,3</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"	<b>A1*</b> \$485
HR05	<input type="checkbox"/>	Surge LT <sup>1,2,3</sup> (1 1/8" Oval handrim w/ Gripton high friction strip)				22", 24", 25", 26"	<b>A1*</b> \$515
HR04	<input type="checkbox"/>	Natural Fit <sup>1,2,3</sup>				22", 24", 25", 26"	<b>A1*</b> \$320
HR62	<input type="checkbox"/>	Simi H <sup>1,2,3</sup> (1 3/8" Oval handrim with high friction strip)				22", 24", 25", 26"	<b>A1*</b> \$485
HR63	<input type="checkbox"/>	Nova H <sup>1,2,3</sup> (1 1/4" Oval handrim with high friction strip)				22", 24", 25", 26"	<b>A1*</b> \$515
HR64	<input type="checkbox"/>	Optimum H <sup>1,2,3</sup> (1 1/8" heart-shaped handrim)				22", 24", 25", 26"	<b>A1*</b> \$320
HR99	<input type="checkbox"/>	Omit handrim					NC

*1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;*

**HANDRIMS OPTION**

ASSEMBLY POSITION				Available			
HR29	<input type="checkbox"/>	Standard Handrim position	<i>Full tab</i>				STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut	<i>Full tab</i>				NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut <sup>3</sup>	<i>Short tab 2 positions</i>				\$29
HR31	<input type="checkbox"/>	Super Narrow Handrim <sup>3</sup>	<i>Short tab 2 positions STD</i> <input type="checkbox"/>	<i>Short tab 1 position</i>	<input type="checkbox"/>		\$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS							
HR23	<input type="checkbox"/>	Standard grip				<b>A1*</b>	\$200
HR27	<input type="checkbox"/>	Super grip				<b>A1*</b>	\$230
HR099	<input type="checkbox"/>	Omit thumb piece					NC


*3- Only available with Aluminum and Air Grip handrims;*

AXLES (0° camber) <sup>1</sup>						
AX01	<input type="checkbox"/>	Quick release axle				STD
AX02	<input type="checkbox"/>	Permanent axle				NC
AX03	<input type="checkbox"/>	Quad release axle			A2*	\$180
AXLE PLATE						
AX04	<input type="checkbox"/>	Amputee adjustable axle plate <sup>2</sup>			A2*	\$275

1- Not required with one arm drive; 2- May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

\*A1 is clients' responsibility for costs/A2 requires clinical justification

#### CENTER OF GRAVITY



#### Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable		Most responsive
AMPUTEE <sup>1</sup>	3/4"	2"
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

1- Only available with amputee axle plate

#### WHEEL LOCKS

WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock				STD
WL16	<input type="checkbox"/>	Aluminum push to lock <b>with</b> extensions				\$73
WL02	<input type="checkbox"/>	Pull to lock				NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock ( <i>works best with pneumatic or high pressure tire</i> )			A2*	\$185
WL13	<input type="checkbox"/>	6" removable extension handle <sup>1</sup>				\$73
WL14	<input type="checkbox"/>					
WL06	<input type="checkbox"/>	Attendant Lock ( <i>additional set of wheel lock mounted to the back</i> )			A1*	\$240
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock <sup>2</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>			A1*	\$335
WL99	<input type="checkbox"/>	Omit wheel locks				NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

[https://www.motioncomposites.com/en\\_ca/support-and-education/how-to-documents](https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents)

#### FOOTREST

SWING-AWAY FOOTREST						
FR01-01	<input type="checkbox"/>	60° swing away footrest	13.5 - 17.5"		Set at:*	\$270
FR01-02	<input type="checkbox"/>	60° swing away footrest	15.5 - 19.5"		Set at:*	\$270
FR01-03	<input type="checkbox"/>	60° swing away footrest	17 - 21"		Set at:*	\$270
FR02-01	<input type="checkbox"/>	70° swing away footrest	13.5 - 17.5"		Set at:*	\$270
FR02-02	<input type="checkbox"/>	70° swing away footrest	15.5 - 19.5"		Set at:*	\$270
FR02-03	<input type="checkbox"/>	70° swing away footrest	17 - 21"		Set at:*	\$270
FOOTREST OPTIONS						
FR09	<input type="checkbox"/>	High mount ( <i>footplate attached to the hanger with clamps</i> )		5" - 13.5"	Set at:*	\$105
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, Length adjustable <sup>1</sup>		Left <input type="checkbox"/> Right <input type="checkbox"/>	Set at:	\$395
FR99	<input type="checkbox"/>	Omit Footrest				NC

1- The elevating legrest replace the footrest side selected; Not available with high mount

\* +/- 1/2"

FOOTPLATES						
FP01	<input type="checkbox"/>	Standard - Composite				STD
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite			W107	\$155
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum			W107	\$175

FOOTPLATE ACCESSORIES						
FP08	<input type="checkbox"/>	Length adjustable heel loop				\$26
FP09	<input type="checkbox"/>	Calf strap				\$57
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap			basic seating benefit	\$140

SWING AWAY RESIDUAL LIMB SUPPORT						
FR27	<input type="checkbox"/>	Residual limb support - flat 6"x8"		Left <input type="checkbox"/> Right <input type="checkbox"/>		\$505
FR28	<input type="checkbox"/>	Residual limb support - flat 8"x10"		Left <input type="checkbox"/> Right <input type="checkbox"/>		\$505
FR29	<input type="checkbox"/>	Residual limb support - flat 10"x10"		Left <input type="checkbox"/> Right <input type="checkbox"/>		\$505
FR30	<input type="checkbox"/>	Residual limb support - flat 10"x14"		Left <input type="checkbox"/> Right <input type="checkbox"/>		\$505
FR32	<input type="checkbox"/>	Residual limb support - contour 6"x8"		Left <input type="checkbox"/> Right <input type="checkbox"/>		\$505
FR33	<input type="checkbox"/>	Residual limb support - contour 8"x10"		Left <input type="checkbox"/> Right <input type="checkbox"/>		\$505

**SEATING AND BACKREST**

BACK UPHOLSTERY							
BH05	<input type="checkbox"/>	Slip-on Nylon back upholstery					NC
BH06	<input type="checkbox"/>	Adjustable tension Nylon back upholstery				<i>Basic seating benefit</i>	\$260
BH99	<input type="checkbox"/>	Omit back upholstery					STD
SEAT CUSHIONS							
SC01	<input type="checkbox"/>	2" Cushion				<i>Basic seating benefit</i>	\$185
SC02	<input type="checkbox"/>	3" Cushion				<i>Basic seating benefit</i>	\$205
SC99	<input type="checkbox"/>	Omit seat cushion					NC
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE							
BC01-16	<input type="checkbox"/>	8° bend aluminum back cane	16"			NC	
BC01-18	<input type="checkbox"/>	8° bend aluminum back cane	18"			NC	
BC01-20	<input type="checkbox"/>	8° bend aluminum back cane	20"			NC	
BC59	<input type="checkbox"/>	8° bend aluminum back cane	21" - 24"	<i>increment of 1"</i>	Set at <sup>2</sup> :	\$100	
ADJUSTABLE ANGLE BACK WITH INTEGRATED PUSH HANDLE							
BC02-16	<input type="checkbox"/>	8° bend aluminum back cane 16"					\$209
BC02-18	<input type="checkbox"/>	8° bend aluminum back cane 18"					\$209
BC02-20	<input type="checkbox"/>	8° bend aluminum back cane 20"					\$209
Select back cane angle:		90° <input type="checkbox"/>	95° <b>STD</b> <input type="checkbox"/>	100° <input type="checkbox"/>	105° <input type="checkbox"/>	110° <input type="checkbox"/>	
STROLLER HANDLE AND STABILIZER BAR <sup>1</sup>							
PH04	<input type="checkbox"/>	Stroller handle (not available with the HD option)					\$220
PH03	<input type="checkbox"/>	Newton folding stabilizer bar (not available with stroller handle)					\$275

1- Stroller handle and stabilizer bar can't be combined

2- If not indicated, will be at the lowest position

*\*A1 is clients' responsibility for costs/A2 requires clinical justification*

**ARMRESTS**

HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) <sup>1,2</sup>							
AR01-812	<input type="checkbox"/>	10" Desk length armpad	8" - 12"	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>	\$395	
AR01-1014	<input type="checkbox"/>	10" Desk length armpad	10" - 14"	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>	\$395	
AR02-812	<input type="checkbox"/>	14" Full length armpad	8" - 12"	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>	\$395	
AR02-1014	<input type="checkbox"/>	14" Full length armpad	10" - 14"	Locked <input type="checkbox"/>	Unlocked STD <input type="checkbox"/>	\$395	
HEIGHT ADJUSTABLE "T" ARMREST <sup>1</sup>							
AR13-610	<input type="checkbox"/>	10" Desk length armpad <sup>4</sup>	6" - 10"			\$380	
AR03-812	<input type="checkbox"/>	10" Desk length armpad	8" - 12"			\$380	
AR03-1014	<input type="checkbox"/>	10" Desk length armpad	10" - 14"			\$380	
AR14-610	<input type="checkbox"/>	14" Full length armpad <sup>4</sup>	6" - 10"			\$380	
AR04-812	<input type="checkbox"/>	14" Full length armpad	8" - 12"			\$380	
AR04-1014	<input type="checkbox"/>	14" Full length armpad	10" - 14"			\$380	
ARMREST OPTIONS							
AR12	<input type="checkbox"/>	14" Full Length armpad - Gel Ovation (pair) <sup>3</sup>					\$160
AR99	<input type="checkbox"/>	Omit armrest					NC

1- May have 1/2" adjust. restriction on both limits of the range; 2- Not available with attendant lock; 3- Only available with full length armrest; Not available with swing away and desk length armrests

4- May have restriction with seat-to-floor height and rear wheel

**ACCESSORIES**

POSITIONING BELTS							
SB01-S	<input type="checkbox"/>	Velcro adjustable - <b>Small 35"</b> <sup>1</sup>				<i>Basic seating benefit</i>	\$50
SB01-L	<input type="checkbox"/>	Velcro adjustable - <b>Long 60"</b>				<i>Basic seating benefit</i>	\$50
SB02-S	<input type="checkbox"/>	Auto buckle - <b>Small 35"</b> <sup>1</sup>				<i>Basic seating benefit</i>	\$22
SB02-M	<input type="checkbox"/>	Auto buckle - <b>Medium 48"</b>				<i>Basic seating benefit</i>	\$22
SB02-L	<input type="checkbox"/>	Auto buckle - <b>Long 60"</b>				<i>Basic seating benefit</i>	\$22
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - <b>Medium 56"</b>				<i>Basic seating benefit</i>	\$95
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - <b>Long 61"</b>				<i>Basic seating benefit</i>	\$95
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Small 51"</b>				<i>Basic seating benefit</i>	\$185
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Medium 56"</b>				<i>Basic seating benefit</i>	\$185
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Long 61"</b>				<i>Basic seating benefit</i>	\$185
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Small 51"</b>				<i>Basic seating benefit</i>	\$210
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Medium 56"</b>				<i>Basic seating benefit</i>	\$210
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Long 61"</b>				<i>Basic seating benefit</i>	\$210
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)				<i>Basic seating benefit</i>	\$15
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Medium 18"</b>				<i>Basic seating benefit</i>	\$355
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Long 24"</b>				<i>Basic seating benefit</i>	\$355
SB99	<input type="checkbox"/>	Omit positioning belt					STD

1- Not recommended for seat width of 18" and more

HEADREST <sup>1</sup>			
AC01-S	<input type="checkbox"/>	Profiled Headrest - <b>Small (3.5"X6")</b>	<i>basic seating benefit</i> \$215
AC01-M	<input type="checkbox"/>	Profiled Headrest - <b>Medium (4.5"X7")</b>	<i>Basic seating benefit</i> \$215
AC01-L	<input type="checkbox"/>	Profiled Headrest - <b>Large (5.5"X8")</b>	<i>Basic seating benefit</i> \$215
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar <sup>1</sup>	<i>Basic seating benefit</i> \$360

1- Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES			
AC03	<input type="checkbox"/>	Anti-tippers	\$120
AC29	<input type="checkbox"/>	Swing-away anti-tipper <sup>1</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>	\$185
AC27	<input type="checkbox"/>	Tip assist <sup>2</sup>	\$65
AC06	<input type="checkbox"/>	Cane and crutch holder      Left <input type="checkbox"/> Right <input type="checkbox"/>	\$115
AC07	<input type="checkbox"/>	Oxygen tank holder <sup>4</sup> - D and E Type - Installed to the left	\$398
AC08	<input type="checkbox"/>	Spoke guards <sup>3</sup> Available: 20", 22", 24", 26"	\$155
AC09	<input type="checkbox"/>	Footrest neoprene impact guard ( <i>pair</i> ) - Black	\$85
AC10	<input type="checkbox"/>	Seat pouch	\$100
AC30	<input type="checkbox"/>	Newton multi-purpose pouch	\$75
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag	\$145
AC23	<input type="checkbox"/>	Newton Backpack	\$130
AC12	<input type="checkbox"/>	Velcro Closing Strap	\$55
AC13	<input type="checkbox"/>	Tool kit	\$55

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate; 2- Cannot be combined with Anti Tippers


3- Not available with Mag wheels and Spinergy; 4- Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down

Colors																				
Options <sup>1</sup>	Black*	Sunkissed Orange	Burgundy	Sapphire Blue*	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)	
AirGrip	<input type="checkbox"/> \$160	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$350
Frame <sup>2</sup>	<input type="checkbox"/> STD																		\$50	
FC99	<input type="checkbox"/> Custom Color	<a href="#">Specify PPG Paintit color:</a>														<i>additional 10 days lead-time</i>			\$422	

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- with black crossbrace;

### Fork and Stem bolt combinations



### Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	2
	17+	0.75 to 2
22 in	16	0.75 to 2
	17	0.75 to 2
	18+	0.75 to 2
24 in	16	0.75
	17	0.75 to 2
	18+	0.75 to 2

Notes :