



MOVE ORDER FORM - CANADA

P.O. number:*	Date:	Ship to:*
Company Name: *	Phone:*	Name:
Purchasing contact:*		Address:
Clinician/Sales rep:		City: Prov:
Email:		Postal Code:
Marked for:		Phone:

* Required MA080 ☐ English Manual (STD) MA081 ☐ Manuel Français Taxes not included

MODEL					
WC01	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)		WM0001196	\$1 425	☐
WC09	Heavy duty kit (350 lb weight capacity)		WAMF	\$800	☐
SEAT WIDTH ³					
SW	Seat width: 14" ☐ 15" ☐ 16" ☐ 17" ☐ 18" ¹ ☐			NC	☐
SW	Seat width: 19" ¹ ☐ 20" ¹ ☐ 22" ² ☐		WAMA	\$305	☐
SEAT DEPTH					
SD	Seat depth: 16" ☐ 18" ☐			NC	☐
SD	Seat depth: 20" ☐		WAMB	\$230	☐

1 - Available with HD option

2 - Only available with HD option

[3. Width calculator](#)

Move FRONT Seat-to-floor Height - Please select your choice

		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
Fork/Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2	13 1/2	14 1/2												
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1															
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1							15 1/2	16 1/2	17 1/2	16	17	18			
	P2				15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3				16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS										
CA01	Composite wheel w/ PU tire (1" width)	4" ☐	5" ☐	6" ☐	7" ☐	8" ☐			NC	☐
CA10	Pneumatic 6" x 1 ¼"								\$165	☐
CA11	Pneumatic 8" x 1 ¼"								\$165	☐
CA14	Newton UltraCaster <i>Soft-Roll</i> Composite	4" x 1 ½ "							\$98	☐
CA15		5" x 1 ½ "							\$98	☐
CA16		6" x 1 ½ "							\$98	☐
CA17	Newton UltraCaster <i>Soft-Roll</i> Aluminum	4" x 1 ½ "							\$180	☐
CA18		5" x 1 ½ "							\$190	☐
CA19		6" x 1 ½ "							\$200	☐

Move REAR Seat-to-floor Height - Please circle your choice

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12"	13"	14"	14 1/2"	15"
13"	14"	15"	15 1/2"	16"
14"	15"	16"	16 1/2"	17"
15"	16"	17"	17 1/2"	18"
16"	17"	18"	18 1/2"	19"
17"	18"	19"	19 1/2"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest									
SWING-AWAY FOOTREST									
FR01	60° swing away footrest	13.5 - 17.5" ☐	15.5 - 19.5" ☐	17 - 21" ☐	Set at:*			\$260	☐
FR02	70° swing away footrest	13.5 - 17.5" ☐	15.5 - 19.5" ☐	17 - 21" ☐	Set at:*			\$260	☐
FOOTREST OPTIONS									
FR09	High mount (footplate attached to the hanger with clamps)	5 - 13.5" ☐			Set at:*			\$95	☐
FR10	Elevating legrest with calf pad, 14" - 22" Length adjustable ¹	Left ☐	Right ☐		Set at:	WAMK		\$390	☐

1- The elevating legrest replace the footrest side selected

* +/- 1/2"

FOOTPLATES									
FP01	Standard - Composite							STD	☐
FP02	Newton adjustable angle - Composite					WAMJ		\$150	☐
FP03	Newton adjustable angle - Aluminum					WAMJ		\$170	☐
FOOTPLATE ACCESSORIES									
FP08	Length adjustable heel loop							\$25	☐
FP09	Calf strap					SEFND0010		\$55	☐
FP30	Bodypoint Aeromesh Padded calf strap					SEFND0010		\$135	☐

SWING AWAY RESIDUAL LIMB SUPPORT									
FR27	Residual limb support - flat 6"x8"	Left ☐	Right ☐			SEFND0005		\$490	☐
FR28	Residual limb support - flat 8"x10"	Left ☐	Right ☐			SEFND0005		\$490	☐
FR29	Residual limb support - flat 10"x10"	Left ☐	Right ☐			SEFND0005		\$490	☐
FR30	Residual limb support - flat 10"x14"	Left ☐	Right ☐			SEFND0005		\$490	☐
FR32	Residual limb support - contour 6"x8"	Left ☐	Right ☐			SEFND0005		\$490	☐
FR33	Residual limb support - contour 8"x10"	Left ☐	Right ☐			SEFND0005		\$490	☐

Backrest									
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE									
BC01	8° bend aluminum back cane	16-20"	increment of 2"	Set at:				NC	☐
ADJUSTABLE ANGLE BACK									
BC02	8° bend aluminum back cane	16" ☐	20" STD ☐					\$320	☐
	Select back cane angle:	90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐			

STROLLER HANDLE AND STABILIZER BAR ¹									
PH04	Stroller handle <i>(not available with the HD option)</i>						WAMW	\$215	☐
PH03	Folding stabilizer bar							\$265	☐

1 - Stroller handle and stabilizer bar can't be combined

BACK UPHOLSTERY					
BH05	Slip-on Nylon back upholstery			NC	☐
BH06	Adjustable tension Nylon back upholstery		WAMD/WAME	\$250	☐
BH99	Omit back upholstery			STD	☐
	NXT back series	See pricelist			☐
SEAT CUSHIONS					
SC01	2" Cushion			\$180	☐
SC02	3" Cushion			\$200	☐
SC99	Omit seat cushion			NC	☐
	NXT FIT cushion series	See pricelist			☐

Rear wheels										
ONE ARM DRIVE ¹ Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims										
HR22	One arm drive ^{2,3}	24"	left ☐	right ☐				\$1 435	☐	
HR98	Omit handrim on opposite side of the One Arm Drive (OAD)							NC	☐	
REAR WHEELS										
RW01	Mag ⁴	20" ☐	22" ☐	24" ☐				NC	☐	
RW02	Newton One - spoke wheel ⁵	20" ☐	22" ☐	24" ☐				NC	☐	
RW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐		\$370	☐	
TIRES						Available				
TI01	Soft Urethane 1 1/4 "				20", 22" 24"				STD	☐
TI02	Soft Urethane knobby 1 3/8 "				22", 24", 25", 26"				\$135	☐
TI04	Soft Urethane - Shox 1"				22", 24", 25", 26"				\$160	☐
TI05	Pneumatic 1 3/8 "				20", 22", 24", 25", 26"				\$135	☐
TI08	Pneumatic with Airless insert 1 3/8"				22", 24", 26"				\$220	☐
TI06	Pneumatic HP - 1"				20", 22", 24", 25", 26"				\$155	☐
TI07	Pneumatic HP - Schwalbe Marathon Plus 1"				24", 25", 26"				\$235	☐

1- Offered with handrims on each side of the wheelchair 2- Only standard handrim position available 3- Not available with HD option

4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 5- Not compatible with Shox tires

Handrims									
HANDRIM									
HR01	Aluminum anodized							STD	☐
HR61	Aluminum black hard anodized ^{1,3}							\$77	☐
HR02	Plastic coated ²						WAMO	\$165	☐
HR03	Q-Grip High friction Neoprene coated ^{2,3,4}						WAMO	\$445	☐
HR07	Newton Air Grip ³ <i>Black</i> ☐						WAMO	\$155	☐
HR	Newton Air Grip ³	(08) <i>Sunkissed Orange</i> ☐	(09) <i>Burgundy</i> ☐	(10) <i>Sapphire Blue</i> ☐			WAMO	\$340	☐
	(11) <i>Monster Green</i> ☐	(12) <i>Charcoal</i> ☐	(13) <i>Fuchsia</i> ☐	(14) <i>Steel Blue</i> ☐					
	(15) <i>Ferrari Red</i> ☐	(16) <i>Acid Green</i> ☐	(17) <i>White</i> ☐	(18) <i>Sky Blue</i> ☐					
HR28	Surge ^{2,3,4} <i>(1-1/2" Oval handrim w/ Gripton high friction strip)</i>							\$485	☐
HR05	Surge LT ^{2,3,4} <i>(1-1/8" Oval handrim w/ Gripton high friction strip)</i>							\$515	☐
HR04	Natural Fit ^{2,3,4}							\$320	☐
HR99	Omit handrim							NC	☐

1 - Not available in 26" 2- Only standard handrim position available 3 - Not available with Mag wheels or OAD 4 - Not available in 20"


Handrims option					
Assembly position for aluminum and Newton AirGrip handrims					
HR29	Standard Handrim position	Full tab		STD	☐
HR30	Narrow Handrim position	Short tab 2 positions		\$26	☐
HR31	Super Narrow Handrim	Short tab 2 positions STD ☐ Short tab 1 position ☐		\$26	☐
Thumb piece ⁴ selection Natural Fit handrims					
HR23	Standard grip			\$195	☐
HR27	Super grip			\$225	☐
HR099	Omit thumb piece			NC	☐

4- Compatible with Surge, Surge LT handrims

AXLES (0° camber) ¹					
AX01	Quick release axle		WAMU	STD	☐
AX02	Permanent axle			NC	☐
AX03	Quad release axle		WAMU	\$175	☐
AXLE PLATE					
AX04	Amputee adjustable axle plate ²		WAMT	\$260	☐

1- Not required with one arm drive

2- May limit seat to floor height available



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.
Default setting is at 3/4"

Most stable	Most responsive
AMPUTEE	3/4" 2"
☐	☐

WHEEL LOCKS					
WL08	Newton Wheel lock - push to lock			STD	☐
WL02	Pull to lock			NC	☐
WL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)		WAMR	\$180	☐
WL13	6" removable extension handle			\$70	☐
WL14					
WL06	Attendant Lock (additional set of wheel lock mounted to the back)			\$235	☐
WL99	Omit wheel locks			NC	☐

Armrests					
FLIP BACK ARMREST (Convertible to "T" armrest) ^{1,2}			Removable function		
Height Adjustment					
AR01	10" Desk length arm pad	8 - 12" STD ☐ 10 - 14" ☐	Locked ☐ Unlocked ☐	STD	\$390 ☐
AR02	14" Full length arm pad	8 - 12" STD ☐ 10 - 14" ☐	Locked ☐ Unlocked ☐	STD	\$390 ☐
"T" ARMREST ¹					
Height Adjustment					
AR03	10" Desk length arm pad	8 - 12" STD ☐ 10 - 14" ☐			\$370 ☐
AR04	14" Full length arm pad	8 - 12" STD ☐ 10 - 14" ☐			\$370 ☐
GEL ARM PAD ³					
AR12	14" Full Length Arm Pad - Gel Ovation		SEAND0030	\$155	☐
AR99	Omit armrest			NC	☐

1- May have a 1/2" adjustment restriction on both limits of the range 2- Not available with attendant lock 3- Not available with flip back armrest

POSITIONING BELTS							
SB01	Velcro adjustable	Standard 35" ☐	Long 60" ☐			\$47	☐
SB02	Auto buckle	Standard 35" ☐	Long 60" ☐			\$57	☐
SB07	Hardware attachment for 2 point seatbelt (1 pair)					NC	☐
SB03	Additional hardware for 2 extra attachment points <i>(needed for the 4 point seat belt) (1 pair)</i>					\$14	☐
SB04	BodyPoint Hip Belt - 2 point not padded	Medium 56" ☐	Long 61" ☐			\$85	☐
SB05	BodyPoint Hip Belt - 2 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	SERND0030	\$175	☐
SB06	BodyPoint Hip Belt - 4 point padded	Small 51" ☐	Medium 56" ☐	Long 61" ☐	SERND0020	\$200	☐
SB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	Medium 18" ☐	Long 24" ☐			\$335	☐
SB99	Omit positioning belt					NC	☐

HEADREST ¹							
AC01	Profiled Headrest	Small (3.5"X6") ☐	Medium (4.5"X7") ☐	Large (5.5"X8") ☐		\$210	☐
AC28	Headrest mounting hardware with stabiliser bar					\$350	☐

1 - Not available with a stroller handle or a folding stabiliser bar

ACCESSORIES													
AC03	Anti tippers							\$120	☐				
AC29	Swing-away anti-tippers ¹	Left ☐	Right ☐				\$175	☐					
AC27	Tip assist ²						\$62	☐					
AC04	Transit Tie Down <i>(Unoccupied)</i>						NC	☐					
AC05	WC-19 Transit Tie Down w/Q'straint belt ³	Short belt 20-41" ☐	Long belt 20-46" ☐				\$355	☐					
AC06	Cane carrier	Left ☐	Right ☐				\$110	☐					
AC07	Oxygen tank holder - D and E Type	Left				WAMX	\$385	☐					
AC08	Spoke guards ⁴	Available: 20", 22", 24", 26"				WAMP	\$150	☐					
AC09	Footrest neoprene impact guard <i>(pair)</i> - Black					SEFND0035	\$82	☐					
AC10	Seat pouch						\$98	☐					
AC30	Newton sanitizing multi-purpose pouch						\$72	☐					
AC26	BodyPoint Mobility Bag						\$140	☐					
AC23	Newton Backpack						\$125	☐					
AC12	Velcro Closing Strap						\$52	☐					
AC13	Tool kit						\$52	☐					
AC16	Motion Composites Polo shirt	W/S ☐	W/M ☐	W/L ☐	M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	M/XXXL ☐	M/XXXXL ☐	\$35	☐

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height of 13 1/2" and 19 3/4"


2 - Cannot be combined with Anti Tippers

3 - May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

4 - Not available with Mag wheels and Spinergy

FRAME COLOR <i>(with black aluminum crossbrace)</i>							
FC01	Black ☐					NC	☐
FC	(03) Sunkissed Orange ☐	(04) Burgundy ☐	(05) Sapphire Blue ☐	(06) Monster Green ☐		\$270	☐
	(07) Charcoal ☐	(08) Fuchsia ☐	(09) Steel Blue ☐	(17) Ferrari Red ☐	(18) Acid Green ☐		
FC99	Custom Color Specify PPG Paintit color:			Additional 10 days lead-time		\$410	☐

Fork and Stem bolt combinations



Notes :

