



1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

Marked for: _____

3- Shipping information

Ship to:* _____

Name: _____

Address: _____

City: _____ Province: _____

Postal Code: _____

Phone: _____

2- Dealer/Representative information

Dealer/Representative:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

* Required

MA080 ☐
English Manual (STD)

MA081 ☐
Manuel Français

Taxes not included

MODEL								
WC01	☐	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)					WM0001196	\$1,425
WC09	☐	Heavy duty kit (350 lb weight capacity)					WAMF	\$825
Frame								
SEAT WIDTH ³								
SW	☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ¹ ☐		NC
	☐	19" ¹ ☐	20" ¹ ☐	22" ² ☐			WAMA	\$315
SEAT DEPTH								
SD	☐	16" ☐			18" ☐			NC
	☐	20" ☐					WAMB	\$235
TRANSIT TIE DOWN								
AC04	☐	Transit Tie Down (Unoccupied)						NC
AC05-S	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Short belt 20"-41"						\$365
AC05-L	☐	WC-19 Transit Tie Down w/Q'straint belt ⁴ Long belt 20"-46"						\$365
AC35	☐	WC-19 Transit Tie Down without belt ⁴						\$335

1- Available with HD option; 2- Only available with HD option

[3. Width calculator](#)

4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

Move FRONT Seat-to-floor Height - Please select your choice

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
4"	P1	12 1/2"	13 1/2"	14 1/2"												
	P2	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
5"	P1															
	P2	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
	P3	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"						
	P4	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"			
7"	P1							15 1/2"	16 1/2"	17 1/2"	16"	17"	18"			
	P2				15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
	P3				16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"
	P4	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"
	P5	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form)

CASTERS						
CA01-04	☐	Composite wheel w/ PU tire	4" x 1"			NC
CA01-05	☐		5" x 1"			NC
CA01-06	☐		6" x 1"			NC
CA01-07	☐		7" x 1"			NC
CA01-08	☐		8" x 1"			NC
CA10	☐	Pneumatic casters	6" x 1 ¼"			\$170
CA11	☐		8" x 1 ¼"			\$170
CA14	☐	Newton UltraCaster Composite	4" x 1 ½ "			\$100
CA15	☐		5" x 1 ½ "			\$100
CA16	☐		6" x 1 ½ "			\$100
CA17	☐	Newton UltraCaster Aluminum	4" x 1 ½ "			\$185
CA18	☐		5" x 1 ½ "			\$195
CA19	☐		6" x 1 ½ "			\$205

Move REAR Seat-to-floor Height - Please circle your choice

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12"	13"	14"	14 1/2"	15"
13"	14"	15"	15 1/2"	16"
14"	15"	16"	16 1/2"	17"
15"	16"	17"	17 1/2"	18"
16"	17"	18"	18 1/2"	19"
17"	18"	19"	19 1/2"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Footrest						
SWING-AWAY FOOTREST						
FR01-01	☐	60° swing away footrest	13.5 - 17.5"	Set at:*		\$270
FR01-02	☐	60° swing away footrest	15.5 - 19.5"	Set at:*		\$270
FR01-03	☐	60° swing away footrest	17 - 21"	Set at:*		\$270
FR02-01	☐	70° swing away footrest	13.5 - 17.5"	Set at:*		\$270
FR02-02	☐	70° swing away footrest	15.5 - 19.5"	Set at:*		\$270
FR02-03	☐	70° swing away footrest	17 - 21"	Set at:*		\$270
FOOTREST OPTIONS						
FR09	☐	High mount <i>(footplate attached to the hanger with clamps)</i>		5" - 13.5"	Set at:*	\$100
FR10	☐	Elevating legrest with calf pad, Length adjustable ¹	Left ☐ Right ☐	14" - 22"	Set at:	WAMK \$395
FR99	☐	Omit Footrest				NC

1- The elevating legrest replace the footrest side selected

* +/- 1/2"

FOOTPLATES								
FP01	☐	Standard - Composite						STD
FP02	☐	Newton adjustable angle - Composite					WAMJ	\$155
FP03	☐	Newton adjustable angle - Aluminum					WAMJ	\$175
FOOTPLATE ACCESSORIES								
FP08	☐	Length adjustable heel loop						\$26
FP09	☐	Calf strap					SEFND0010	\$57
FP30	☐	Bodypoint Aeromesh Padded calf strap					SEFND0010	\$140
SWING AWAY RESIDUAL LIMB SUPPORT								
FR27	☐	Residual limb support - flat 6"x8"	Left ☐	Right ☐		SEFND0005	\$505	
FR28	☐	Residual limb support - flat 8"x10"	Left ☐	Right ☐		SEFND0005	\$505	
FR29	☐	Residual limb support - flat 10"x10"	Left ☐	Right ☐		SEFND0005	\$505	
FR30	☐	Residual limb support - flat 10"x14"	Left ☐	Right ☐		SEFND0005	\$505	
FR32	☐	Residual limb support - contour 6"x8"	Left ☐	Right ☐		SEFND0005	\$505	
FR33	☐	Residual limb support - contour 8"x10"	Left ☐	Right ☐		SEFND0005	\$505	

Backrest									
FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE									
BC01-16	☐	8° bend aluminum back cane	16"				NC		
BC01-18	☐	8° bend aluminum back cane	18"				NC		
BC01-20	☐	8° bend aluminum back cane	20"				NC		
BC59	☐	8° bend aluminum back cane	21" - 24"	increment of 1"	Set at ² :		\$100		
ADJUSTABLE ANGLE BACK WITH INTEGRATED PUSH HANDLE									
BC02-16	☐	8° bend aluminum back cane	16"				\$330		
BC02-20	☐	8° bend aluminum back cane	20"				\$330		
	Select back cane angle:			90° ☐	95° STD ☐	100° ☐	105° ☐	110° ☐	
STROLLER HANDLE AND STABILIZER BAR ¹									
PH04	☐	Stroller handle (not available with the HD option)				WAMW	\$220		
PH03	☐	Newton folding stabilizer bar (not available with stroller handle)					\$275		
1 - Stroller handle and stabilizer bar can't be combined 2- If not indicated, will be at the lowest position									
BACK UPHOLSTERY									
BH05	☐	Slip-on Nylon back upholstery					NC		
BH06	☐	Adjustable tension Nylon back upholstery				WAMD/WAME	\$260		
BH99	☐	Omit back upholstery					STD		
	☐	NXT back Support			See pricelist				
SEAT CUSHIONS									
SC01	☐	2" Cushion					\$185		
SC02	☐	3" Cushion					\$205		
SC99	☐	Omit seat cushion					NC		
	☐	NXT cushion				See price list			

Rear wheels								
ONE ARM DRIVE ¹ <i>Only available with Newton One rear wheels with Aluminum, Plastic coated Handrims, permanent axles, 0° camber</i>								
HR22	☐	One arm drive - aluminum inner handrim ^{2, 3}	24"	Left ☐	Right ☐		WAML	\$1,480
HR60	☐	One arm drive - plastic coated inner handrim ^{2, 3}	24"	Left ☐	Right ☐		WAML	\$1,480
HR98	☐	Omit handrim on opposite side of the One Arm Drive (OAD)						NC
REAR WHEELS						Available		
RW01	☐	Mag ⁴				20", 22" 24"		NC
RW02	☐	Newton One - spoke wheel ⁵				20", 22" 24"		NC
RW03	☐	Newton Gravity - ultralight wheel				20", 22", 24", 25", 26"		\$380
TIRES						Available		
TI01	☐	Soft Urethane - 1 3/8"				20", 22" 24"		STD
TI02	☐	Soft Urethane - knobby 1 3/8 "				22", 24", 25", 26"		\$140
TI04	☐	Soft Urethane - Shox 1"				22", 24", 25", 26"		\$165
TI05	☐	Pneumatic - 1 3/8"				20", 22", 24", 25", 26"		\$140
TI08	☐	Pneumatic with Airless insert - 1 3/8"				22", 24", 26"		\$225
TI06	☐	Pneumatic HP - 1"				20", 22", 24", 25", 26"		\$160
TI07	☐	Pneumatic HP - Puncture resistant, black 1"				24", 25", 26"		\$240

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; 3- Not available with HD option
4 - Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position 5- Not compatible with Shox tires

Handrims						
HANDRIM				Available		
HR01	☐	Aluminum anodized		20", 22", 24", 25", 26"		STD
HR61	☐	Aluminum black hard anodized ²		20", 22", 24", 25"		\$80
HR02	☐	Plastic coated ¹		20", 22", 24", 25", 26"	WAMO	\$170
HR03	☐	High friction coated ^{1,2}		22", 24", 25", 26"	WAMO	\$460
HR07 to HR18	☐	Newton Air Grip ²	See page 5 for the color selection and price	20", 22", 24", 25", 26"	WAMO	
HR28	☐	Surge ^{1,2}	(1 1/2" Oval handrim w/ Gripton high friction strip)	**Limited availability	22", 24", 25", 26"	\$485
HR05	☐	Surge LT ^{1,2}	(1 1/8" Oval handrim w/ Gripton high friction strip)	**Limited availability	22", 24", 25", 26"	\$515
HR04	☐	Natural Fit ^{1,2}		**Limited availability	22", 24", 25", 26"	\$320
HR63	☐	Nova H ^{1,2}	(1 1/4" Oval handrim with high friction strip)		22", 24", 25", 26"	\$515
HR64	☐	Optimum H ^{1,2}	(1 1/8" heart-shaped handrim)		22", 24", 25", 26"	\$320
HR99	☐	Omit handrim				NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels

Handrims option						
ASSEMBLY POSITION						
HR29	☐	Standard Handrim position	Full tab			STD
HR030	☐	Narrow Handrim position no cut	Full tab			NC
HR30	☐	Narrow Handrim position with cut ³	Short tab 2 positions			\$27
HR31	☐	Super Narrow Handrim ³	Short tab 2 positions STD ☐ Short tab 1 position ☐			\$27
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS						
HR23	☐	Standard grip	**Limited availability			\$200
HR27	☐	Super grip	**Limited availability			\$230
HR099	☐	Omit thumb piece				NC

3- Only available with Aluminum and Air Grip handrims;

AXLES (0° camber) ¹						
AX01	☐	Quick release axle			WAMU	STD
AX02	☐	Permanent axle				NC
AX03	☐	Quad release axle			WAMU	\$180
AXLE PLATE						
AX04	☐	Amputee adjustable axle plate ²			WAMT	\$270

1- Not required with one arm drive

2 - May limit seat to floor height available

CENTER OF GRAVITY

Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable

AMPUTEE ¹

3/4"

2"

Most responsive

☐
☐
☐

1- Only available with amputee axle plate

WHEEL LOCKS						
WL08	☐	Newton Wheel lock - push to lock				STD
WL16		Aluminum push to lock with extension				\$73
WL02	☐	Pull to lock				NC
WL04	☐	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)			WAMR	\$185
WL13 WL14	☐	6" removable extension handle				\$73
WL06	☐	Attendant Lock (additional set of wheel lock mounted to the back)				\$240
WL15	☐	Unilateral wheel lock - pull to lock ¹ Left ☐ Right ☐			WAMN	\$335
WL99	☐	Omit wheel locks				NC

1- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website: https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests						
HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ^{1,2}						
AR01-812	☐	10" Desk length armpad 8" - 12"	Locked ☐ Unlocked STD ☐			\$395
AR01-1014	☐	10" Desk length armpad 10" - 14"	Locked ☐ Unlocked STD ☐			\$395
AR02-812	☐	14" Full length armpad 8" - 12"	Locked ☐ Unlocked STD ☐			\$395
AR02-1014	☐	14" Full length armpad 10" - 14"	Locked ☐ Unlocked STD ☐			\$395
HEIGHT ADJUSTABLE "T" ARMREST ¹						
AR13-610	☐	10" Desk length armpad ⁵ 6" - 10"				\$380
AR03-812	☐	10" Desk length armpad 8" - 12"				\$380
AR03-1014	☐	10" Desk length armpad 10" - 14"				\$380
AR14-610	☐	14" Full length armpad ⁵ 6" - 10"				\$380
AR04-812	☐	14" Full length armpad 8" - 12"				\$380
AR04-1014	☐	14" Full length armpad 10" - 14"				\$380
ARMREST OPTIONS						
AR12	☐	14" Full Length armpad - Gel Ovation (pair) ³			SEAND0030	\$160
AR99	☐	Omit armrest				NC

1- May have 1/2" adjust. restriction on both limits of the range; 2- Not available with attendant lock; 3- Not available with flip back armrest; 4- May have restriction with seat-to-floor height and rear wheel

POSITIONING BELTS							
SB01-S	☐	Velcro adjustable - Standard 35"					\$50
SB01-L	☐	Velcro adjustable - Long 60"					\$50
SB02-S	☐	Auto buckle - Standard 35"					\$60
SB02-L	☐	Auto buckle - Long 60"					\$60
SB04-M	☐	BodyPoint Hip Belt - 2 point not padded - Medium 56"					\$90
SB04-L	☐	BodyPoint Hip Belt - 2 point not padded - Long 61"					\$90
SB05-S	☐	BodyPoint Hip Belt - 2 point padded - Small 51"			SERND0030		\$180
SB05-M	☐	BodyPoint Hip Belt - 2 point padded - Medium 56"			SERND0030		\$180
SB05-L	☐	BodyPoint Hip Belt - 2 point padded - Long 61"			SERND0030		\$180
SB06-S	☐	BodyPoint Hip Belt - 4 point padded - Small 51"			SERND0020		\$205
SB06-M	☐	BodyPoint Hip Belt - 4 point padded - Medium 56"			SERND0020		\$205
SB06-L	☐	BodyPoint Hip Belt - 4 point padded - Long 61"			SERND0020		\$205
SB03	☐	Hardware attachment for 2 extra points (1 pair)					\$15
SB08-M	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18"					\$345
SB08-L	☐	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24"					\$345
SB99	☐	Omit positioning belt					STD

HEADREST ¹							
AC01-S	☐	Profiled Headrest - Small (3.5"X6")					\$215
AC01-M	☐	Profiled Headrest - Medium (4.5"X7")					\$215
AC01-L	☐	Profiled Headrest - Large (5.5"X8")					\$215
AC28	☐	Headrest mounting hardware with removable non folding stabilizer bar ¹					\$360

1- Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES								
AC03	☐	Anti-tippers						\$120
AC29	☐	Swing-away anti-tipper ¹	Left ☐	Right ☐				\$180
AC27	☐	Tip assist ²						\$65
AC06	☐	Cane and crutch holder	Left ☐	Right ☐				\$115
AC07	☐	Oxygen tank holder ⁴ - D and E Type - Installed to the left						\$398
AC08	☐	Spoke guards ³	Available: 20", 22", 24", 26"					\$155
AC09	☐	Footrest neoprene impact guard (<i>pair</i>) - Black						\$85
AC10	☐	Seat pouch						\$100
AC30	☐	Newton multi-purpose pouch						\$75
AC26	☐	BodyPoint Mobility Bag						\$145
AC23	☐	Newton Backpack						\$130
AC12	☐	Velcro Closing Strap						\$55
AC13	☐	Tool kit						\$55
AC16	☐	Motion Composites Polo Shirt	W/S ☐	W/M ☐	W/L ☐			\$36
			M/S ☐	M/M ☐	M/L ☐	M/XL ☐	M/XXL ☐	

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"

2 - Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinerigy; 4- Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down


Colors																			
Options ¹	Black	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)
																			
AirGrip	☐ \$160	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
Frame ²	NC	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
FC99	☐	Custom Color Specify PPG Paintit color: <i>additional 10 days lead-time</i>																	\$422

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.;

Finish between Airgrip and frame might be different;

2- with black crossbrace

Fork and Stem bolt combinations



Unilateral wheel lock restrictions


Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	2
22 in	16	0.75 to 2
	17+	0.75 to 2
24 in	16	0.75
	17	0.75 to 2
	18+	0.75 to 2

Notes :

FOLDING ON CHAIR - CDN ORDER FORM

BioFit™ CUSHION		WIDTH										
		12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"
		☐ \$625	☐ \$625	☐ \$625								
DEPTH	12"	☐ \$625	☐ \$625	☐ \$625								
	14"	☐ \$625	☐ \$625	☐ \$625								
	15"		☐ \$625	☐ \$625	☐ \$625	☐ \$625						
	16"		☐ \$625	☐ \$625	☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655	☐ \$655	☐ \$655	
	17"				☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655			
	18"				☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655	☐ \$655	☐ \$655	
	19"				☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655			
	20"				☐ \$625	☐ \$625	☐ \$625	☐ \$655	☐ \$655	☐ \$655	☐ \$655	☐ \$655
Option: Change cover for incontinent / infection control outer cover												☐ N/C

SEATING AND CUSHION ACCESSORIES

	Model #	Description	Fits Back Width	\$	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6 & 43.2 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7 & 48.3 cm) and 14"-18" deep (35.6-45.7 cm)	80	
	RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	80	

Optima™ THORACIC BACK SUPPORT

	Optima™	Fits Width	Height	✓	\$
	NOACB1613	16"-18"	13"	☐	\$725
	NOACB1616	16"-18"	16"	☐	\$725
	NOACB1619	16"-18"	19"	☐	\$725
	NOACB1813	18"-20"	13"	☐	\$725
	NOACB1816	18"-20"	16"	☐	\$725
	NOACB1819	18"-20"	19"	☐	\$725
	NOACB2016	20"-22"	16"	☐	\$725
	NOACB2019	20"-22"	19"	☐	\$725
	NOACB2213	22"-24"	13"	☐	\$755
	NOACB2216	22"-24"	16"	☐	\$755
	NOACB2219	22"-24"	19"	☐	\$755

Optima™ DEEP THORACIC BACK SUPPORT

	Optima™ DEEP	Fits Width	Height	✓	\$
	NOPCB1616	16"-18"	16"	☐	\$825
	NOPCB1619	16"-18"	19"	☐	\$825
	NOPCB1816	18"-20"	16"	☐	\$825
	NOPCB1819	18"-20"	19"	☐	\$825
	NOPCB2016	20"-22"	16"	☐	\$825
	NOPCB2019	20"-22"	19"	☐	\$825
	NOPCB2213	22"-24"	13"	☐	\$860
	NOPCB2216	22"-24"	16"	☐	\$860
	NOPCB2219	22"-24"	19"	☐	\$860

Optima™ CARBON THORACIC BACK SUPPORT

	Optima™ Carbon	Fits Width	Height	✓	\$
	NOCFACB1413	14"-16"	13"	☐	\$1150
	NOCFACB1416	14"-16"	16"	☐	\$1150
	NOCFACB1513	15"-17"	13"	☐	\$1150
	NOCFACB1516	15"-17"	16"	☐	\$1150
	NOCFACB1613	16"-18"	13"	☐	\$1150
	NOCFACB1616	16"-18"	16"	☐	\$1150
	NOCFACB1713	17"-19"	13"	☐	\$1150
	NOCFACB1716	17"-19"	16"	☐	\$1150
	NOCFACB1813	18"-20"	13"	☐	\$1150
	NOCFACB1816	18"-20"	16"	☐	\$1150
	NOCFACB2013	20"-22"	13"	☐	\$1150
	NOCFACB2016	20"-22"	16"	☐	\$1150

Optima™ CARBON DEEP THORACIC BACK SUPPORT

	Optima™ Carbon DEEP	Fits Width	Height	✓	\$
	NOCFPCB1413	14"-16"	13"	☐	\$1250
	NOCFPCB1416	14"-16"	16"	☐	\$1250
	NOCFPCB1513	15"-17"	13"	☐	\$1250
	NOCFPCB1516	15"-17"	16"	☐	\$1250
	NOCFPCB1613	16"-18"	13"	☐	\$1250
	NOCFPCB1616	16"-18"	16"	☐	\$1250
	NOCFPCB1713	17"-19"	13"	☐	\$1250
	NOCFPCB1716	17"-19"	16"	☐	\$1250
	NOCFPCB1813	18"-20"	13"	☐	\$1250
	NOCFPCB1816	18"-20"	16"	☐	\$1250
	NOCFPCB2013	20"-22"	13"	☐	\$1250
	NOCFPCB2016	20"-22"	16"	☐	\$1250

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com

FOLDING ON CHAIR - CDN ORDER FORM

Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

THORACIC SHAPE		Xtend™	Fits Width	Height	✓	\$
		**NHAACB (L)1213/17	12"-14"	13"-17"	☐	\$725
		NHAACB (L)1413/17	14"-16"	13"-17"	☐	\$725
		NHAACB (L)1613/17	16"-18"	13"-17"	☐	\$725
		NHAACB (L)1813/17	18"-20"	13"-17"	☐	\$755
		NHAACB (L)2013/17	20"-22"	13"-17"	☐	\$755
		NHAACB (L)2213/17	22"-24"	13"-17"	☐	\$755

Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT

	Xtend™ DEEP	Fits Width	Height	✓	\$
	NHAPCB1616/20	16"-18"	16"-20"	☐	\$725
	NHAPCB1816/20	18"-20"	16"-20"	☐	\$725
	NHAPCB2016/20	20"-22"	16"-20"	☐	\$725
	NHAPCB2216/20	22"-24"	16"-20"	☐	\$755

Classic™ PELVIC BACK SUPPORT

PELVIC SHAPE		Classic™	Fits Width	Height	✓	\$
		NACB1616	16"-18"	16"	☐	\$720
		NACB1618	16"-18"	18"	☐	\$720
		NACB1816	18"-20"	16"	☐	\$720
		NACB1818	18"-20"	18"	☐	\$720
		NACB1820	18"-20"	20"	☐	\$720
		NACB2020	20"-22"	20"	☐	\$720
		NACB2216	22"-24"	16"	☐	\$755
		NACB2218	22"-24"	18"	☐	\$755
		NACB2220	22"-24"	20"	☐	\$755

Classic™ DEEP PELVIC BACK SUPPORT

	Classic DEEP™	Fits Width	Height	✓	\$
	NPCB1616	16"-18"	16"	☐	\$825
	NPCB1618	16"-18"	18"	☐	\$825
	NPCB1818	18"-20"	18"	☐	\$825
	NPCB1820	18"-20"	20"	☐	\$825
	NPCB2020	20"-22"	20"	☐	\$825
	NPCB2216	22"-24"	16"	☐	\$860
	NPCB2218	22"-24"	18"	☐	\$860
	NPCB2220	22"-24"	20"	☐	\$860

BACK SUPPORT ACCESSORIES - PRIVACY FLAP

	Model #	Description	Fits Back Width	\$	QTY
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	70	

MultiFit™ HEAD SUPPORT - MULTI-ADJUSTABLE PADS

			
Model # NIB-0784	✓	\$	
10" x 5" Large Rectangular	☐	\$205	

			
Model # NIB-0818	✓	\$	
14" x 5" Extra-Large Rectangular	☐	\$205	

MultiFit™ HEAD SUPPORT - FIXTURES

		
Model # NIB-0747	✓	\$
Multi-Axis Hardware	☐	\$150

		
Model # NIB-0823	✓	\$
Head Support Receiver Bracket	☐	\$60

ALL PRICES QUOTED IN CDN. PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at www.motioncomposites.com