



**1- Purchasing information**

P.O. number:\*  Date:

Company Name:\*

Purchasing contact:\*

Phone:\*

Email:\*

Contact preference:\* Phone  Email

Marked for:

**3- Shipping information**

Ship to:\*

Name:

Address:

City:  Province:

Postal Code:

Phone:

**2- Dealer/Representative information**

Dealer/Representative:\*

Phone:\*

Email:\*

Contact preference:\* Phone  Email

MA080  English Manual (STD)      MA081  Manuel Français      Taxes not included

MODEL		Frame							
WC01	<input type="checkbox"/>	Ultralight MOVE wheelchair aluminum. (265 lb weight capacity)							\$1,425
WC09	<input type="checkbox"/>	Heavy duty kit (350 lb weight capacity)							\$850
SEAT WIDTH <sup>3</sup>									
SW	<input type="checkbox"/>	14" <input type="checkbox"/>	15" <input type="checkbox"/>	16" <input type="checkbox"/>	17" <input type="checkbox"/>	18" <sup>1</sup> <input type="checkbox"/>			NC
	<input type="checkbox"/>	19" <sup>1</sup> <input type="checkbox"/>	20" <sup>1</sup> <input type="checkbox"/>	21" <sup>2</sup> <input type="checkbox"/>	22" <sup>2</sup> <input type="checkbox"/>				\$320
SEAT DEPTH									
SD	<input type="checkbox"/>	16" <input type="checkbox"/>		18" <input type="checkbox"/>					NC
	<input type="checkbox"/>	20" <input type="checkbox"/>							\$280
FRAME LENGTH <sup>5</sup>									
SP01	<input type="checkbox"/>	Standard frame length configuration							STD
SP02	<input type="checkbox"/>	Frame length 1" longer than standard configuration							NC
SP03	<input type="checkbox"/>	Frame length 2" longer than standard configuration							NC
SP04	<input type="checkbox"/>	Frame length 1" shorter than standard configuration							NC
SP05	<input type="checkbox"/>	Frame length 2" shorter than standard configuration							NC
TRANSIT TIE DOWN									
AC04	<input type="checkbox"/>	Transit Tie Down (Unoccupied)							NC
AC05-S	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup>		Short belt 20"-41"					\$370
AC05-L	<input type="checkbox"/>	WC-19 Transit Tie Down w/Q'straint belt <sup>4</sup>		Long belt 20"-46"					\$370
AC35	<input type="checkbox"/>	WC-19 Transit Tie Down without belt <sup>4</sup>							\$345

1- Available with HD option; 2- Only available with HD option  
 4- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option  
 5- See table on page 6

[3. Width calculator](#)



**PACKAGES**

VP01	<input type="checkbox"/>	Vision Package: Orange accent added to point of contact accessories including: footrest latches, "T" and flip back armrest button cap, 6" Push to lock extension and rubber inserts on Newton composite angle adjustable footplates <sup>1</sup>	\$95
------	--------------------------	--	------

1- Please select footplates, footrest, armrest, push to lock wheel lock and the extensions in the appropriate sections; To know the compatible parts, look for this logo \*VP01; Recommended with black chair; See picture of the product on page 5.

**Move FRONT Seat-to-floor Height - Please select your choice**

Fork/Position		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
		STD	+1"	+2"												
4"	P1	12 1/2"	13 1/2"	14 1/2"												
	P2	13"	14"	15"		14 1/2"	15 1/2"									
5"	P1															
	P2	13"	14"	15"	13 1/2"	14 1/2"	15 1/2"									
	P3	13 1/2"	14 1/2"	15 1/2"	14"	15"	16"		15 1/2"	16 1/2"						
	P4	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"		16 1/2"	17 1/2"			
7"	P1							15 1/2"	16 1/2"	17 1/2"	16"	17"	18"			
	P2				15 1/2"	16 1/2"	17 1/2"	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"			
	P3				16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"
	P4	16"	17"	18"	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"
	P5	16 1/2"	17 1/2"	18 1/2"	17"	18"	19"	17 1/2"	18 1/2"	19 1/2"	18"	19"	20"	18 1/2"	19 1/2"	20 1/2"
Froglegs	P1	13 1/2"	14 1/2"	15 1/2"												
	P2	14"	15"	16"	14 1/2"	15 1/2"	16 1/2"									
	P3	14 1/2"	15 1/2"	16 1/2"	15"	16"	17"	15 1/2"	16 1/2"	17 1/2"						

Height available +/- 1/4" (see stem bolt and fork combination at the end of the order form); Orange boxes are not available with a seat slope of more than 1.5in.

**FRONT CASTERS**

CASTERS				Available	
CA01-04 to CA01-08	<input type="checkbox"/>	Composite wheel w/ PU tire (1")		4", 5", 6", 7", 8"	NC
CA10 to CA11	<input type="checkbox"/>	Pneumatic casters (1 1/4")		6", 8"	\$170
CA14 to CA16	<input type="checkbox"/>	Newton UltraCaster Composite (1 1/2")		4", 5", 6"	\$120
CA17 to CA19	<input type="checkbox"/>	Newton UltraCaster Aluminum (1 1/2")		4", 5", 6"	\$198

**Move REAR Seat-to-floor Height - Please circle your choice**

20" (451)	22" (501)	24" (540)	25" (559)	26" (590)
12"	13"	14"	14 1/2"	15"
13"	14"	15"	15 1/2"	16"
14"	15"	16"	16 1/2"	17"
15"	16"	17"	17 1/2"	18"
16"	17"	18"	18 1/2"	19"
17"	18"	19"	19 1/2"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

**REAR WHEELS**

ONE ARM DRIVE <sup>1</sup> Only available with Newton One rear wheels, permanent axles, 0° camber				Available	
HR22	<input type="checkbox"/>	One arm drive - aluminum inner handrim <sup>2</sup>	Left <input type="checkbox"/> Right <input type="checkbox"/>		\$1,480
HR60	<input type="checkbox"/>	One arm drive - plastic coated inner handrim <sup>2</sup>	Left <input type="checkbox"/> Right <input type="checkbox"/>		\$1,595
HR98	<input type="checkbox"/>	Omit handrim on opposite side of the One Arm Drive (OAD)			NC
REAR WHEELS				Available	
RW01	<input type="checkbox"/>	Mag <sup>3</sup>		20", 22" 24"	NC
RW02	<input type="checkbox"/>	Newton One - spoke wheel <sup>4</sup>		20", 22" 24"	NC
RW03	<input type="checkbox"/>	Newton Gravity - ultralight wheel		20", 22", 24", 25", 26"	\$390
TIRES				Available	
AIRLESS					
TI01	<input type="checkbox"/>	Soft Urethane 1 3/8" - Low Tread, black		20", 22" 24"	STD
TI02	<input type="checkbox"/>	Soft Urethane 1 3/8" - Medium Tread, black <sup>5</sup>		22", 24", 25", 26"	\$140
TI04	<input type="checkbox"/>	Soft Urethane 1" - No tread, black - Shox		22", 24", 25", 26"	\$175
TI08	<input type="checkbox"/>	Pneumatic Airless insert 1 3/8" - Medium Tread		22", 24", 26"	\$235
PNEUMATIC					
TI05	<input type="checkbox"/>	Pneumatic 1 3/8" - Medium Tread		20", 22", 24", 25", 26"	\$150
TI06	<input type="checkbox"/>	Pneumatic HP 1" - Low Tread, Puncture resistant, black - SpeedLite		20", 22", 24", 25", 26"	\$170
TI07	<input type="checkbox"/>	Pneumatic HP 1" - Medium Tread, Puncture resistant, black - TrailBlazer		24", 25", 26"	\$250
TI03	<input type="checkbox"/>	Pneumatic 2" - Large Tread, black		22", 24", 25"	\$140

1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; Only available with 24" wheel size; Not available with 21" and 22" wide 2- Only standard handrim position available;

3- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position; 4- Not compatible with Shox tires;

5- The 22" is only available with Mags and Newton One rear wheels; 25" and 26" are not available with Newton Gravity wheels

## HANDRIMS

HANDRIM		Available		
HR01	<input type="checkbox"/>	Aluminum anodized	20", 22", 24", 25", 26"	STD
HR61	<input type="checkbox"/>	Aluminum black hard anodized <sup>2</sup>	20", 22", 24", 25"	\$95
HR02	<input type="checkbox"/>	Plastic coated <sup>1</sup>	20", 22", 24", 25", 26"	\$180
HR03	<input type="checkbox"/>	High friction coated <sup>1,2,3</sup>	22", 24", 25", 26"	\$460
HR07-HR18 and HR65-HR71	<input type="checkbox"/>	Newton Air Grip <sup>2</sup> <span style="color: red;">See page 5 for the color selection and price</span>	20", 22", 24", 25", 26"	
HR28	<input type="checkbox"/>	Surge <sup>1,2,3,4</sup> (1 1/2" Oval handrim w/ Gripton high friction strip)	22", 24", 25", 26"	\$485
HR05	<input type="checkbox"/>	Surge LT <sup>1,2,3,4</sup> (1 1/8" Oval handrim w/ Gripton high friction strip)	22", 24", 25", 26"	\$515
HR04	<input type="checkbox"/>	Natural Fit <sup>1,2,3,4</sup>	22", 24", 25", 26"	\$320
HR62	<input type="checkbox"/>	Simi H <sup>1,2,3</sup> (1 3/8" Oval handrim with high friction strip)	22", 24", 25", 26"	\$485
HR63	<input type="checkbox"/>	Nova H <sup>1,2,3</sup> (1 1/4" Oval handrim with high friction strip)	22", 24", 25", 26"	\$515
HR64	<input type="checkbox"/>	Optimum H <sup>1,2,3</sup> (1 1/8" heart-shaped handrim)	22", 24", 25", 26"	\$320
HR99	<input type="checkbox"/>	Omit handrim		NC

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive; 4 - 25" and 26" are not available in narrow position; 24" is not available with the Spinegy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

## HANDRIMS OPTION

ASSEMBLY POSITION				
HR29	<input type="checkbox"/>	Standard Handrim position <span style="float: right;">Full tab</span>		STD
HR030	<input type="checkbox"/>	Narrow Handrim position no cut <span style="float: right;">Full tab</span>		NC
HR30	<input type="checkbox"/>	Narrow Handrim position with cut <sup>5</sup> <span style="float: right;">Short tab 2 positions</span>		\$29
HR31	<input type="checkbox"/>	Super Narrow Handrim <sup>5</sup> <span style="float: right;">Short tab 2 positions STD <input type="checkbox"/> Short tab 1 position <input type="checkbox"/></span>		\$29
THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS				
HR23	<input type="checkbox"/>	Standard grip		\$200
HR27	<input type="checkbox"/>	Super grip		\$230
HR099	<input type="checkbox"/>	Omit thumb piece		NC

5- Only available with Aluminum and Air Grip handrims;

### AXLES (0° camber) <sup>1</sup>


AX01	<input type="checkbox"/>	Quick release axle		STD
AX02	<input type="checkbox"/>	Permanent axle		NC
AX03	<input type="checkbox"/>	Quad release axle		\$180

### AXLE PLATE

AX04	<input type="checkbox"/>	Amputee adjustable axle plate <sup>2</sup>		\$280
------	--------------------------	--	--	-------

1- Not required with one arm drive; 2 - May interfere with swing away armrest and limit seat to floor height; Not available with WC-19 transit tie down

## CENTER OF GRAVITY



### Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable	3/4"	Most responsive
AMPUTEES <sup>1</sup>	3/4"	2"
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

1- Only available with amputee axle plate

## WHEEL LOCKS

WL08	<input type="checkbox"/>	Newton Wheel lock - push to lock		STD
WL16	<input type="checkbox"/>	Aluminum push to lock <b>with</b> extensions		\$85
WL02	<input type="checkbox"/>	Pull to lock		NC
WL04	<input type="checkbox"/>	Newton grade aid push to lock (works best with pneumatic or high pressure tire)		\$200
WL13 WL14	<input type="checkbox"/>	6" removable extension handle <sup>1</sup> *VP01		\$75
WL06	<input type="checkbox"/>	Attendant Lock (additional set of wheel lock mounted to the back)		\$240
WL15	<input type="checkbox"/>	Unilateral wheel lock - pull to lock <sup>2</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>		\$335
WL99	<input type="checkbox"/>	Omit wheel locks		NC

1- Only available with Newton push to lock (WL08), pull to lock (WL02) and Newton Grade aid (WL04); Only the extension for the Newton push to lock is available with the vision package (VP01)

2- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:

[https://www.motioncomposites.com/en\\_ca/support-and-education/how-to-documents](https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents)

FOOTREST									
<b>SWING-AWAY FOOTREST</b>									
FR01-01	<input type="checkbox"/>	60° swing away footrest	13.5 - 17.5"					Set at*:	\$270
FR01-02	<input type="checkbox"/>	60° swing away footrest	15.5 - 19.5"					Set at*:	\$270
FR01-03	<input type="checkbox"/>	60° swing away footrest	17 - 21"					Set at*:	\$270
FR02-01	<input type="checkbox"/>	70° swing away footrest	13.5 - 17.5"					Set at*:	\$270
FR02-02	<input type="checkbox"/>	70° swing away footrest	15.5 - 19.5"					Set at*:	\$270
FR02-03	<input type="checkbox"/>	70° swing away footrest	17 - 21"					Set at*:	\$270
<b>FOOTREST OPTIONS</b>									
FR09	<input type="checkbox"/>	High mount <i>(footplate attached to the hanger with clamps)</i>	5" - 13.5"					Set at*:	\$105
FR10	<input type="checkbox"/>	Elevating legrest with calf pad, Length adjustable <sup>1</sup>	Left <input type="checkbox"/> Right <input type="checkbox"/>	14" - 22"				Set at:	\$395
FR99	<input type="checkbox"/>	Omit Footrest							NC
<i>1- The elevating legrest replace the footrest side selected; Not available with high mount</i>									
<i>* +/- 1/2"</i>									
<b>FOOTPLATES</b>									
FP01	<input type="checkbox"/>	Standard - Composite							STD
FP02	<input type="checkbox"/>	Newton adjustable angle - Composite <sup>*VP01</sup>							\$155
FP03	<input type="checkbox"/>	Newton adjustable angle - Aluminum							\$175
FP23	<input type="checkbox"/>	Newton adjustable angle extra deep - Aluminum							\$175
<b>FOOTPLATE ACCESSORIES</b>									
FP08	<input type="checkbox"/>	Length adjustable heel loop							\$26
FP09	<input type="checkbox"/>	Calf strap							\$62
FP30	<input type="checkbox"/>	Bodypoint Aeromesh Padded calf strap							\$145
<b>SWING AWAY RESIDUAL LIMB SUPPORT</b>									
FR27	<input type="checkbox"/>	Residual limb support - flat 6"x8"	Left <input type="checkbox"/> Right <input type="checkbox"/>						\$505
FR28	<input type="checkbox"/>	Residual limb support - flat 8"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>						\$505
FR29	<input type="checkbox"/>	Residual limb support - flat 10"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>						\$505
FR30	<input type="checkbox"/>	Residual limb support - flat 10"x14"	Left <input type="checkbox"/> Right <input type="checkbox"/>						\$505
FR32	<input type="checkbox"/>	Residual limb support - contour 6"x8"	Left <input type="checkbox"/> Right <input type="checkbox"/>						\$505
FR33	<input type="checkbox"/>	Residual limb support - contour 8"x10"	Left <input type="checkbox"/> Right <input type="checkbox"/>						\$505
<b>SEATING AND BACKREST</b>									
<b>BACK UPHOLSTERY</b>									
BH05	<input type="checkbox"/>	Slip-on Nylon back upholstery							NC
BH06	<input type="checkbox"/>	Adjustable tension Nylon back upholstery							\$265
BH99	<input type="checkbox"/>	Omit back upholstery							STD
	<input type="checkbox"/>	NXT back Support						See price list	
<b>SEAT CUSHIONS</b>									
SC01	<input type="checkbox"/>	2" Cushion							\$185
SC02	<input type="checkbox"/>	3" Cushion							\$205
SC99	<input type="checkbox"/>	Omit seat cushion							NC
	<input type="checkbox"/>	NXT cushion						See price list	
<b>FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE</b>									
BC01-16	<input type="checkbox"/>	8° bend aluminum back cane	16"						NC
BC01-18	<input type="checkbox"/>	8° bend aluminum back cane	18"						NC
BC01-20	<input type="checkbox"/>	8° bend aluminum back cane	20"						NC
BC59	<input type="checkbox"/>	8° bend aluminum back cane	21" - 24"	increment of 1"				Set at <sup>2</sup> :	\$100
<b>ADJUSTABLE ANGLE BACK WITH INTEGRATED PUSH HANDLE</b>									
BC02-16	<input type="checkbox"/>	8° bend aluminum back cane	16"						\$330
BC02-20	<input type="checkbox"/>	8° bend aluminum back cane	20"						\$330
		Select back cane angle:	90° <input type="checkbox"/> 95° <b>STD</b> <input type="checkbox"/> 100° <input type="checkbox"/> 105° <input type="checkbox"/> 110° <input type="checkbox"/>						
<b>STROLLER HANDLE AND STABILIZER BAR <sup>1</sup></b>									
PH04	<input type="checkbox"/>	Stroller handle <i>(not available with the HD option)</i>							\$220
PH03	<input type="checkbox"/>	Newton folding stabilizer bar <i>(not available with stroller handle)</i>							\$275
<i>1- Stroller handle and stabilizer bar can't be combined 2- If not indicated, will be at the lowest position</i>									
<b>ARMRESTS</b>									
<b>HEIGHT ADJUSTABLE FLIP BACK ARMREST <i>(Convertible to "T" armrest)</i><sup>2,4</sup> *VP01</b>									
AR01-812	<input type="checkbox"/>	10" Desk length armpad	8" - 12"	Locked <input type="checkbox"/> Unlocked <input type="checkbox"/>					STD \$405
AR01-1014	<input type="checkbox"/>	10" Desk length armpad	10" - 14"	Locked <input type="checkbox"/> Unlocked <input type="checkbox"/>					STD \$405
AR02-812	<input type="checkbox"/>	14" Full length armpad	8" - 12"	Locked <input type="checkbox"/> Unlocked <input type="checkbox"/>					STD \$405
AR02-1014	<input type="checkbox"/>	14" Full length armpad	10" - 14"	Locked <input type="checkbox"/> Unlocked <input type="checkbox"/>					STD \$405
<b>HEIGHT ADJUSTABLE "T" ARMREST <sup>3</sup> *VP01</b>									
AR13-610	<input type="checkbox"/>	10" Desk length armpad <sup>4</sup>	6" - 10"						\$385
AR03-812	<input type="checkbox"/>	10" Desk length armpad	8" - 12"						\$385
AR03-1014	<input type="checkbox"/>	10" Desk length armpad	10" - 14"						\$385
AR14-610	<input type="checkbox"/>	14" Full length armpad <sup>4</sup>	6" - 10"						\$385
AR04-812	<input type="checkbox"/>	14" Full length armpad	8" - 12"						\$385
AR04-1014	<input type="checkbox"/>	14" Full length armpad	10" - 14"						\$385
<b>ARMREST OPTIONS</b>									
AR12	<input type="checkbox"/>	14" Full Length armpad - Gel Ovation (pair) <sup>3</sup>							\$160
AR99	<input type="checkbox"/>	Omit armrest							NC

1- May have 1/2" adjust. restriction on both limits of the range; 2- Not available with attendant lock; 3- Only available with full length armrest; Not available with swing away and desk length armrests

4- May have restriction with seat-to-floor height and rear wheel

ACCESSORIES

POSITIONING BELTS			
SB01-S	<input type="checkbox"/>	Velcro adjustable - <b>Small 35"</b> <sup>1</sup>	\$50
SB01-L	<input type="checkbox"/>	Velcro adjustable - <b>Long 60"</b>	\$50
SB02-S	<input type="checkbox"/>	Auto buckle - <b>Small 35"</b> <sup>1</sup>	\$80
SB02-M	<input type="checkbox"/>	Auto buckle - <b>Medium 48"</b>	\$80
SB02-L	<input type="checkbox"/>	Auto buckle - <b>Long 60"</b>	\$80
SB04-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - <b>Medium 56"</b>	\$105
SB04-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point not padded - <b>Long 61"</b>	\$105
SB05-S	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Small 51"</b>	\$195
SB05-M	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Medium 56"</b>	\$195
SB05-L	<input type="checkbox"/>	BodyPoint Hip Belt - 2 point padded - <b>Long 61"</b>	\$195
SB06-S	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Small 51"</b>	\$220
SB06-M	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Medium 56"</b>	\$220
SB06-L	<input type="checkbox"/>	BodyPoint Hip Belt - 4 point padded - <b>Long 61"</b>	\$220
SB03	<input type="checkbox"/>	Hardware attachment for 2 extra points (1 pair)	\$15
SB08-M	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Medium 18"</b>	\$365
SB08-L	<input type="checkbox"/>	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - <b>Long 24"</b>	\$365
SB99	<input type="checkbox"/>	Omit positioning belt	STD

1- Not recommended for seat width of 18" and more

HEADREST <sup>1</sup>			
AC01-S	<input type="checkbox"/>	Profiled Headrest - <b>Small (3.5"X6")</b>	\$215
AC01-M	<input type="checkbox"/>	Profiled Headrest - <b>Medium (4.5"X7")</b>	\$215
AC01-L	<input type="checkbox"/>	Profiled Headrest - <b>Large (5.5"X8")</b>	\$215
AC28	<input type="checkbox"/>	Headrest mounting hardware with removable non folding stabilizer bar <sup>1</sup>	\$365

1- Not available with a stroller handle or a folding stabilizer bar

ACCESSORIES			
AC03	<input type="checkbox"/>	Anti-tippers	\$125
AC29	<input type="checkbox"/>	Swing-away anti-tipper <sup>1</sup> Left <input type="checkbox"/> Right <input type="checkbox"/>	\$190
AC27	<input type="checkbox"/>	Tip assist <sup>2</sup>	\$70
AC06	<input type="checkbox"/>	Cane and crutch holder      Left <input type="checkbox"/> Right <input type="checkbox"/>	\$115
AC07	<input type="checkbox"/>	Oxygen tank holder <sup>4</sup> - D and E Type - Installed to the left	\$398
AC08	<input type="checkbox"/>	Spoke guards <sup>3</sup> Available: 20", 22", 24", 25", 26"	\$165
AC09	<input type="checkbox"/>	Footrest neoprene impact guard (pair) - Black	\$90
AC10	<input type="checkbox"/>	Seat pouch	\$105
AC30	<input type="checkbox"/>	Newton multi-purpose pouch	\$80
AC26	<input type="checkbox"/>	BodyPoint Mobility Bag	\$150
AC23	<input type="checkbox"/>	Newton Backpack	\$135
AC12	<input type="checkbox"/>	Velcro Closing Strap	\$60
AC13	<input type="checkbox"/>	Tool kit	\$60

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"; Not available with amputee axle plate;

2- Cannot be combined with Anti Tippers; 3- Not available with Mag wheels and Spinerly; the 24" spoke guard will be installed on the 25" wheel

4 - Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down


Colors

Options <sup>1</sup>	Black	Sunkissed Orange	Burgundy	Sapphire Blue	Monster Green	Charcoal	Fuchsia	Steel Blue	Ferrari Red	Acid Green	Electric Purple	Caribbean Blue	White	Pink Pearl	Plum Purple	Sky Blue	Black Cherry	Khaki Green (matte)	Desert Sand (matte)		
AirGrip <input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$350
Frame <sup>2</sup> <input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	\$50
FC99 <input type="checkbox"/>	Custom Color <a href="#">Specify PPG Paintit color:</a> <i>additional 10 days lead-time</i>																	\$422			

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

2- with black crossbrace

### Fork and Stem bolt combinations



### Frame length chart

Frame length	Seat depth		
	16"	18"	20"
STD			
+1			N/A
+2			N/A
-1			
-2			

### Unilateral wheel lock restrictions

Seat depth: shall be 16 in and greater		
Rear wheel: shall be 24 in and less		
"T" or flip-back armrest only compatible with 0° Camber		
One arm drive: Not compatible with "T" or flip-back armrest		
"T" or flip-back armrest may require to move the rear wheels outwards		
Not compatible with wheel lock lever extension		
Rear wheel diameter	Seat depth (in)	COG Compatibility Range (in)
20 in	16+	2
22 in	16	0.75 to 2
	17+	0.75 to 2
24 in	16	0.75
	17	0.75 to 2
	18+	0.75 to 2



### CUSTOMIZATION STARTS HERE!

With EVO Custom, whether it's a subtle tweak for comfort or mobility, or a fully reimagined component tailored to your unique needs, we're here to support you.



Explore the **EVO Custom Catalog**.

You can include your request or reference an existing customization from our catalog in the notes section below.



Have a customization in mind or a specific request? The **EVO Custom Request Form** is the best place to start.

Notes :



# FOLDING ON CHAIR – CDN ORDER FORM

 <p>BioFit™ CUSHION</p> <p>Use wheelchair configuration Or use this size:</p> <p>NXCABF _____ (W" x D")</p>	<b>WIDTH</b>												
	<b>DEPTH</b>	<b>12"</b>	<b>13"</b>	<b>14"</b>	<b>15"</b>	<b>16"</b>	<b>17"</b>	<b>18"</b>	<b>19"</b>	<b>20"</b>	<b>21"</b>	<b>22"</b>	
		<b>12"</b>	\$625	\$625	\$625		N/A	N/A	N/A	N/A	N/A	N/A	
		<b>13"</b>						N/A	N/A	N/A	N/A	N/A	
		<b>14"</b>	\$625	\$625	\$625								
		<b>15"</b>		\$625	\$625	\$625	\$625						
		<b>16"</b>		\$625	\$625	\$625	\$625	\$625	\$625	\$655	\$655	\$655	
		<b>17"</b>	N/A			\$625	\$625	\$625	\$625	\$655	\$655		
		<b>18"</b>	N/A			\$625	\$625	\$625	\$625	\$655	\$655	\$655	\$655
		<b>19"</b>	N/A			\$625	\$625	\$625	\$625	\$655	\$655		
		<b>20"</b>	N/A			\$625	\$625	\$625	\$655	\$655	\$655	\$655	\$655
		Option: Change cover for incontinent/infection control outer cover											N/C

Green boxes are stocked sizes; the red boxes require customization\*; the N/A boxes are not available sizes.  
\*Customization (red boxes) is available with an upcharge of \$160.00 and has a 4 to 6 weeks additional lead time.

## SEATING AND CUSHION ACCESSORIES

	MODEL #	Description	Fits Back Width	Retail	QTY
	RSSI 1600	Solid Seat Insert	Fits cushions 16" & 17" wide (40.6-43.2 cm) and 14"-18" deep (35.6-45.7 cm)	\$80	
RSSI 1800	Solid Seat Insert	Fits cushions 18" & 19" wide (45.7-48.3 cm) and 14"-18" deep (35.6-45.7 cm)	\$80		
RSSI 2000	Solid Seat Insert	Fits cushions 20"-22" wide (50.8-55.9 cm) and 18"-20" deep (45.7-50.8 cm)	\$80		

<b>THORACIC SHAPE</b>	 <p><b>Optima™</b> THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)</p>					
	WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM
	14"	13"	NOCFACB1413	\$1,150	NOACB1413	\$725
	14"	16"	NOCFACB1416	\$1,150	NOACB1416	\$725
	15"	13"	NOCFACB1513	\$1,150	NOACB1513	\$725
	15"	16"	NOCFACB1516	\$1,150	NOACB1516	\$725
	16"	13"	NOCFACB1613	\$1,150	NOACB1613	\$725
	16"	16"	NOCFACB1616	\$1,150	NOACB1616	\$725
	16"	19"			NOACB1619	\$725
	17"	13"	NOCFACB1713	\$1,150	NOACB1713	\$725
	17"	16"	NOCFACB1716	\$1,150	NOACB1716	\$725
	18"	13"	NOCFACB1813	\$1,150	NOACB1813	\$725
	18"	16"	NOCFACB1816	\$1,150	NOACB1816	\$725
	18"	19"			NOACB1819	\$725
	19"	16"			NOACB1916	\$725
	19"	19"			NOACB1919	\$725
	20"	13"	NOCFACB2013	\$1,150	NOACB2013	\$725
	20"	16"	NOCFACB2016	\$1,150	NOACB2016	\$725
	20"	19"			NOACB2019	\$725
	21"	16"			NOACB2116	\$755
	21"	19"			NOACB2119	\$755
	22"	13"			NOACB2213	\$755
	22"	16"			NOACB2216	\$755
	22"	19"			NOACB2219	\$755
	ITEM		WIDTH X DEPTH		RETAIL	QTY
	Replacement Cover & Foam		NOACC		\$130	

 <p><b>Optima™ DEEP</b> THORACIC BACK SUPPORT CARBON FIBER SHELL (WEIGHT SAVINGS – 1lb)</p>					
WIDTH	HEIGHT	MODEL # CARBON FIBER	\$ CARBON	MODEL # ALUMINUM	\$ ALUMINUM
14"	13"	NOCFPCB1413	\$1,250	NOPCB1413	\$825
14"	16"	NOCFPCB1416	\$1,250	NOPCB1416	\$825
15"	13"	NOCFPCB1513	\$1,250	NOPCB1513	\$825
15"	16"	NOCFPCB1516	\$1,250	NOPCB1516	\$825
16"	13"	NOCFPCB1613	\$1,250	NOPCB1613	\$825
16"	16"	NOCFPCB1616	\$1,250	NOPCB1616	\$825
16"	19"			NOPCB1619	\$825
17"	13"	NOCFPCB1713	\$1,250	NOPCB1713	\$825
17"	16"	NOCFPCB1716	\$1,250	NOPCB1716	\$825
18"	13"	NOCFPCB1813	\$1,250	NOPCB1813	\$825
18"	16"	NOCFPCB1816	\$1,250	NOPCB1816	\$825
18"	19"			NOPCB1819	\$825
19"	16"			NOPCB1916	\$825
19"	19"			NOPCB1919	\$825
20"	13"	NOCFPCB2013	\$1,250	NOPCB2013	\$825
20"	16"	NOCFPCB2016	\$1,250	NOPCB2016	\$825
20"	19"			NOPCB2019	\$825
21"	16"			NOPCB2116	\$860
21"	19"			NOPCB2119	\$860
22"	13"			NOPCB2213	\$860
22"	16"			NOPCB2216	\$860
22"	19"			NOPCB2219	\$860
ITEM		WIDTH X DEPTH		RETAIL	QTY
Replacement Cover & Foam		NOPCC		\$130	

PLEASE NOTE: Each back will fit widths up to 2" wider.

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at [www.motioncomposites.com](http://www.motioncomposites.com)

**Xtend™ LOW HEIGHT ADJUSTABLE THORACIC BACK SUPPORT**

THORACIC SHAPE	xtend™				
	MODEL #	WIDTH	HEIGHT	RETAIL	✓
	NHAACB (L)1213/17	12"	13"-17"	\$725	
	NHAACB (L)1413/17	14"	13"-17"	\$725	
	NHAACB (L)1613/17	16"	13"-17"	\$725	
	NHAACB (L)1813/17	18"	13"-17"	\$755	
	NHAACB (L)2013/17	20"	13"-17"	\$755	
NHAACB (L)2213/17	22"	13"-17"	\$755		
ITEM		WIDTH X DEPTH	RETAIL	QTY	
Replacement Cover & Foam		N3LHAACC	\$130		

**Xtend™ DEEP HEIGHT ADJUSTABLE THORACIC BACK SUPPORT**

	xtend™ Deep				
	MODEL #	WIDTH	HEIGHT	RETAIL	✓
	NHAPCB1616/20	16"	16"-20"	\$725	
	NHAPCB1816/20	18"	16"-20"	\$725	
	NHAPCB2016/20	20"	16"-20"	\$725	
NHAPCB2216/20	22"	16"-20"	\$755		
ITEM		WIDTH X DEPTH	RETAIL	QTY	
Replacement Cover & Foam		N3HAPCC	\$130		

PLEASE NOTE: Each back will fit widths up to 2" wider.

**Classic™ PELVIC BACK SUPPORT**

PELVIC SHAPE	Classic™				
	MODEL #	WIDTH	HEIGHT	RETAIL	✓
	NACB1616	16"	16"	\$725	
	NACB1618	16"	18"	\$725	
	NACB1816	18"	16"	\$725	
	NACB1818	18"	18"	\$725	
	NACB1820	18"	20"	\$725	
	NACB2020	20"	20"	\$725	
	NACB2216	22"	16"	\$755	
	NACB2218	22"	18"	\$755	
NACB2220	22"	20"	\$755		
ITEM		WIDTH X DEPTH	RETAIL	QTY	
Replacement Cover & Foam		NACC	\$130		

**Classic™ DEEP PELVIC BACK SUPPORT**

	Classic™ Deep				
	MODEL #	WIDTH	HEIGHT	RETAIL	✓
	NPCB1616	16"	16"	\$830	
	NPCB1618	16"	18"	\$830	
	NPCB1818	18"	18"	\$830	
	NPCB1820	18"	20"	\$830	
	NPCB2020	18"	20"	\$830	
	NPCB2216	20"	16"	\$860	
	NPCB2218	22"	18"	\$860	
	NPCB2220	22"	20"	\$860	
ITEM		WIDTH X DEPTH	RETAIL	QTY	
Replacement Cover & Foam		NPCC	\$130		

PLEASE NOTE: Each back will fit widths up to 2" wider.

**BACK SUPPORT ACCESSORIES**

	MODEL #	DESCRIPTION	FITS BACK WIDTH	RETAIL	QTY
	Custom - Cover	Optional Infection Control Wipe Down Back Support Cover	Specify back model # when ordering	\$205	
	NIB-1938	Extra Small Privacy Flap for 11" - 13" wide back supports	11" - 13" (27.9 - 33.0)	\$70	
	NIB-1078	Small Privacy Flap for 14"-16" wide back supports	14" - 16" (35.6-40.6 cm)	\$70	
	NIB-1079	Medium Flap for 17"-19" wide back supports	17" - 19" (43.2-48.3 cm)	\$70	
	NIB-1080	Large Privacy Flap for 20"-22" wide back supports	20" - 22" (50.8-55.9 cm)	\$70	
	NIB-1081	Extra Large Privacy Flap for 23" - 26" wide back supports	23" - 26" (58.4 - 66.0 cm)	\$70	

**MultiFit™ HEAD SUPPORT – MULTI-ADJUSTABLE PADS**

MODEL #	Description	RETAIL	✓
NIB-0784	10" x 5" Large Rectangular	\$205	
NIB-0818	14" x 5" Extra-Large Rectangular	\$205	

**MultiFit™ HEAD SUPPORT – FIXTURES**

MODEL	Description	RETAIL	✓
NIB-0747	Multi-Axis Hardware	\$150	
NIB-0823	Head Support Receiver Bracket	\$60	

**EXTRA COVERS**

Extra Outer smartx3D Cover



Extra Outer Incontinent Cover



ITEM	WIDTH X DEPTH	RETAIL	QTY
Extra Outer smartx <sup>3D</sup> Cover – For a cover up to 18" x 18"	NXCCBF _____ v	\$130	
Extra Outer smartx <sup>3D</sup> Cover – For a cover 18" x 20" and above	NXCCBF _____ v	\$150	
Extra Outer Incontinent Cover – For a cover up to 18" x 18"	NXCCAF _____ v	\$130	
Extra Outer Incontinent Cover – For a cover 18" x 20" and above	NXCCAF _____ v	\$150	

ALL PRICES QUOTED IN CAD PRICES ARE SUBJECT TO CHANGE WITHOUT NOTICE.

For additional NXT models, sizes and hardware not shown, please refer to the complete NXT Seating Order Form at [www.motioncomposites.com](http://www.motioncomposites.com)