



1- Purchasing information

P.O. number:* _____ Date: _____

Company Name:* _____

Purchasing contact:* _____

Phone:* _____

Email:* _____

Contact preference:* Phone ☐ Email ☐

Marked for: _____

3- Shipping information

Ship to:* **Manitoba Wheelchair Program**

Name: _____

Address: **1857 Notre Dame Ave.**

City: **Wpg** Province: **MB**

Postal Code: **R3E 3E7**

Phone: **204-975-3250**

NC = NO CHARGE (option available through Manitoba Possible)

J = JUSTIFICATION REQUIRED (provide clinical rationale)

OTP = OPTION TO PURCHASE (client responsible for cost)

N/A = NOT AVAILABLE (option not available through Manitoba Possible)

MA080 
 English Manual (STD)

Taxes not included

| MODEL | | | | | | |
|-------------------------|--------------------------|---|---|---|--|----|
| WC01 | <input type="checkbox"/> | Ultralight MOVE wheelchair aluminum. (265 lb weight capacity) | | | | NC |
| WC09 | <input type="checkbox"/> | Heavy duty kit (350 lb weight capacity) | | | | NC |
| Frame | | | | | | |
| SEAT WIDTH ³ | | | | | | |
| SW | <input type="checkbox"/> | 14" <input type="checkbox"/> | 16" <input type="checkbox"/> | 18" ¹ <input type="checkbox"/> | | NC |
| | <input type="checkbox"/> | 20" ¹ <input type="checkbox"/> | 22" ² <input type="checkbox"/> | | | NC |
| SEAT DEPTH | | | | | | |
| SD | <input type="checkbox"/> | 16" <input type="checkbox"/> | 18" <input type="checkbox"/> | | | NC |
| | <input type="checkbox"/> | 20" <input type="checkbox"/> | | | | NC |
| TRANSIT TIE DOWN | | | | | | |
| AC04 | <input type="checkbox"/> | Transit Tie Down (Unoccupied) ⁴ | | | | NC |
| AC05-S | <input type="checkbox"/> | WC-19 Transit Tie Down w/Q'straint belt ⁵ | Short belt 20"-41" | | | J |
| AC05-L | <input type="checkbox"/> | WC-19 Transit Tie Down w/Q'straint belt ⁵ | Long belt 20"-46" | | | J |
| AC35 | <input type="checkbox"/> | WC-19 Transit Tie Down without belt ⁵ | | | | NC |

1- Available with HD option; 2- Only available with HD option [3. Width calculator](#) 4- required when ordering HD; 5- May interfere with swing away armrest and Clamp on height adjustable push handle; Not available with HD option

Move FRONT Seat-to-floor Height - Please select your choice

| | | 4" caster (OTP) | | | 5" caster (J) | | | 6" caster (NC) | | | 7" caster (OTP) | | | 8" caster (NC) | | |
|---------------|----|------------------|---------|---------|------------------|---------|---------|------------------|---------|---------|------------------|---------|---------|------------------|---------|---------|
| | | Stem bolt length | | | Stem bolt length | | | Stem bolt length | | | Stem bolt length | | | Stem bolt length | | |
| Fork/Position | | STD | +1" | +2" | STD | +1" | +2" | STD | +1" | +2" | STD | +1" | +2" | STD | +1" | +2" |
| 4" | P1 | 12 1/2" | 13 1/2" | 14 1/2" | | | | | | | | | | | | |
| | P2 | 13" | 14" | 15" | 13 1/2" | 14 1/2" | 15 1/2" | | | | | | | | | |
| 5" | P1 | | | | | | | | | | | | | | | |
| | P2 | 13" | 14" | 15" | 13 1/2" | 14 1/2" | 15 1/2" | | | | | | | | | |
| | P3 | 13 1/2" | 14 1/2" | 15 1/2" | 14" | 15" | 16" | 14 1/2" | 15 1/2" | 16 1/2" | | | | | | |
| | P4 | 14" | 15" | 16" | 14 1/2" | 15 1/2" | 16 1/2" | 15" | 16" | 17" | 15 1/2" | 16 1/2" | 17 1/2" | | | |
| 7" | P1 | | | | | | | 15 1/2" | 16 1/2" | 17 1/2" | 16" | 17" | 18" | | | |
| | P2 | | | | 15 1/2" | 16 1/2" | 17 1/2" | 16" | 17" | 18" | 16 1/2" | 17 1/2" | 18 1/2" | | | |
| | P3 | | | | 16" | 17" | 18" | 16 1/2" | 17 1/2" | 18 1/2" | 17" | 18" | 19" | 17 1/2" | 18 1/2" | 19 1/2" |
| | P4 | 16" | 17" | 18" | 16 1/2" | 17 1/2" | 18 1/2" | 17" | 18" | 19" | 17 1/2" | 18 1/2" | 19 1/2" | 18" | 19" | 20" |
| | P5 | 16 1/2" | 17 1/2" | 18 1/2" | 17" | 18" | 19" | 17 1/2" | 18 1/2" | 19 1/2" | 18" | 19" | 20" | 18 1/2" | 19 1/2" | 20 1/2" |

Height available + - 1/4" (see stem bolt and fork combination at the end of the order form)

| CASTERS | | | | Available | | |
|--------------------|--------------------------|--------------------------------------|--|--------------------|--|-----|
| CA01-04 to CA01-08 | <input type="checkbox"/> | Composite wheel w/ PU tire (1") | | 4", 5", 6", 7", 8" | | NC |
| CA10 to CA11 | <input type="checkbox"/> | Pneumatic casters (1 ¼") | | 6", 8" | | OTP |
| CA14 to CA16 | <input type="checkbox"/> | Newton UltraCaster Composite (1 ½ ") | | 4", 5", 6" | | OTP |
| CA17 to CA19 | <input type="checkbox"/> | Newton UltraCaster Aluminum (1 ½ ") | | 4", 5", 6" | | OTP |

Move REAR Seat-to-floor Height - Please circle your choice

| 20" (451) (OTP) | 22" (501) (NC) | 24" (540) (NC) | 25" (559) (OTP) | 26" (590) (OTP) |
|-----------------|----------------|----------------|-----------------|-----------------|
| 12" | 13" | 14" | 14 1/2" | 15" |
| 13" | 14" | 15" | 15 1/2" | 16" |
| 14" | 15" | 16" | 16 1/2" | 17" |
| 15" | 16" | 17" | 17 1/2" | 18" |
| 16" | 17" | 18" | 18 1/2" | 19" |
| 17" | 18" | 19" | 19 1/2" | 20" |

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

| Footrest | | | | | | |
|---------------------|--------------------------|--|--|----------|--|-----|
| SWING-AWAY FOOTREST | | | | | | |
| FR01-01 | <input type="checkbox"/> | 60° swing away footrest | 13.5 - 17.5" | Set at:* | | OTP |
| FR01-02 | <input type="checkbox"/> | 60° swing away footrest | 15.5 - 19.5" | Set at:* | | OTP |
| FR01-03 | <input type="checkbox"/> | 60° swing away footrest | 17 - 21" | Set at:* | | OTP |
| FR02-01 | <input type="checkbox"/> | 70° swing away footrest | 13.5 - 17.5" | Set at:* | | NC |
| FR02-02 | <input type="checkbox"/> | 70° swing away footrest | 15.5 - 19.5" | Set at:* | | NC |
| FR02-03 | <input type="checkbox"/> | 70° swing away footrest | 17 - 21" | Set at:* | | NC |
| FOOTREST OPTIONS | | | | | | |
| FR09 | <input type="checkbox"/> | High mount <i>(footplate attached to the hanger with clamps)</i> | 5" - 13.5" | Set at:* | | J |
| FR10 | <input type="checkbox"/> | Elevating legrest with calf pad, Length adjustable ¹ | Left <input type="checkbox"/> Right <input type="checkbox"/> 14" - 22" | Set at: | | J |
| FR99 | <input type="checkbox"/> | Omit Footrest | | | | J |

1- The elevating legrest replace the footrest side selected; Not available with high mount

* +/- 1/2"

| FOOTPLATES | | | | | | |
|-----------------------|-------------------------------------|--------------------------------------|--|--|--|-----|
| FP01 | <input type="checkbox"/> | Standard - Composite | | | | NC |
| FP02 | <input type="checkbox"/> | Newton adjustable angle - Composite | | | | J |
| FP03 | <input type="checkbox"/> | Newton adjustable angle - Aluminum | | | | OTP |
| FOOTPLATE ACCESSORIES | | | | | | |
| FP08 | <input checked="" type="checkbox"/> | Length adjustable heel loop | | | | NC |
| FP09 | <input type="checkbox"/> | Calf strap | | | | OTP |
| FP30 | <input type="checkbox"/> | Bodypoint Aeromesh Padded calf strap | | | | OTP |

| SWING AWAY RESIDUAL LIMB SUPPORT | | | | | | |
|----------------------------------|--------------------------|--|--|--|--|-----|
| FR27 | <input type="checkbox"/> | Residual limb support - flat 6"x8" | Left <input type="checkbox"/> Right <input type="checkbox"/> | | | OTP |
| FR28 | <input type="checkbox"/> | Residual limb support - flat 8"x10" | Left <input type="checkbox"/> Right <input type="checkbox"/> | | | OTP |
| FR29 | <input type="checkbox"/> | Residual limb support - flat 10"x10" | Left <input type="checkbox"/> Right <input type="checkbox"/> | | | OTP |
| FR30 | <input type="checkbox"/> | Residual limb support - flat 10"x14" | Left <input type="checkbox"/> Right <input type="checkbox"/> | | | OTP |
| FR32 | <input type="checkbox"/> | Residual limb support - contour 6"x8" | Left <input type="checkbox"/> Right <input type="checkbox"/> | | | OTP |
| FR33 | <input type="checkbox"/> | Residual limb support - contour 8"x10" | Left <input type="checkbox"/> Right <input type="checkbox"/> | | | OTP |

| Backrest | | | | | | | | | |
|---|--------------------------|--|--|------------------------------|----------------------------------|--|-------------------------------|-------------------------------|-----|
| FIXED ANGLE BACK CANES WITH INTEGRATED PUSH HANDLE | | | | | | | | | |
| BC01-16 | <input type="checkbox"/> | 8° bend aluminum back cane | | 16" | | | | | NC |
| BC01-18 | <input type="checkbox"/> | 8° bend aluminum back cane | | 18" | | | | | NC |
| BC01-20 | <input type="checkbox"/> | 8° bend aluminum back cane | | 20" | | | | | NC |
| BC59 | <input type="checkbox"/> | 8° bend aluminum back cane | | 21" - 24" | | increment of 1" | | Set at ² : | OTP |
| ADJUSTABLE ANGLE BACK WITH INTEGRATED PUSH HANDLE | | | | | | | | | |
| BC02-16 | <input type="checkbox"/> | 8° bend aluminum back cane 16" | | | | | | | J |
| BC02-20 | <input type="checkbox"/> | 8° bend aluminum back cane 20" | | | | | | | J |
| | Select back cane angle: | | | 90° <input type="checkbox"/> | 95° STD <input type="checkbox"/> | 100° <input type="checkbox"/> | 105° <input type="checkbox"/> | 110° <input type="checkbox"/> | |
| STROLLER HANDLE AND STABILIZER BAR ¹ | | | | | | | | | |
| PH04 | <input type="checkbox"/> | Stroller handle (not available with the HD option) | | | | | | | OTP |
| PH03 | <input type="checkbox"/> | Newton folding stabilizer bar (not available with stroller handle) | | | | | | | OTP |
| 1- Stroller handle and stabilizer bar can't be combined 2- If not indicated, will be at the lowest position | | | | | | | | | |
| BACK UPHOLSTERY | | | | | | | | | |
| BH05 | <input type="checkbox"/> | Slip-on Nylon back upholstery | | | | | | | NC |
| BH06 | <input type="checkbox"/> | Adjustable tension Nylon back upholstery | | | | | | | OTP |
| BH99 | <input type="checkbox"/> | Omit back upholstery | | | | | | | NC |
| | <input type="checkbox"/> | NXT back Support | | | | See pricelist | | | OTP |
| SEAT CUSHIONS | | | | | | | | | |
| SC01 | <input type="checkbox"/> | 2" Cushion | | | | | | | OTP |
| SC02 | <input type="checkbox"/> | 3" Cushion | | | | | | | OTP |
| SC99 | <input type="checkbox"/> | Omit seat cushion | | | | | | | NC |
| | <input type="checkbox"/> | NXT cushion | | | | | | | OTP |
| Rear wheels | | | | | | | | | |
| ONE ARM DRIVE ¹ | | | | | | Only available with Newton One rear wheels, permanent axles, 0° camber | | | |
| HR22 | <input type="checkbox"/> | One arm drive - aluminum inner handrim ^{2, 3} | | 24" | Left <input type="checkbox"/> | Right <input type="checkbox"/> | | | J |
| HR60 | <input type="checkbox"/> | One arm drive - plastic coated inner handrim ^{2, 3} | | 24" | Left <input type="checkbox"/> | Right <input type="checkbox"/> | | | OTP |
| HR98 | <input type="checkbox"/> | Omit handrim on opposite side of the One Arm Drive (OAD) | | | | | | | OTP |
| REAR WHEELS | | | | | | Available | | | |
| RW01 | <input type="checkbox"/> | Mag ⁴ | | | | 20", 22" 24" | | | OTP |
| RW02 | <input type="checkbox"/> | Newton One - spoke wheel ⁵ | | | | 20", 22" 24" | | | NC |
| RW03 | <input type="checkbox"/> | Newton Gravity - ultralight wheel | | | | 20", 22", 24", 25", 26" | | | OTP |
| TIRES | | | | | | Available | | | |
| Solid | | | | | | | | | |
| TI01 | <input type="checkbox"/> | Soft Urethane 1 3/8" | | | | 20", 22" 24" | | | NC |
| TI02 | <input type="checkbox"/> | Soft Urethane knobby 1 3/8" ⁶ | | | | 22", 24", 25", 26" | | | OTP |
| TI04 | <input type="checkbox"/> | Soft Urethane 1" - Shox | | | | 22", 24", 25", 26" | | | OTP |
| Pneumatic | | | | | | | | | |
| TI05 | <input type="checkbox"/> | Pneumatic 1 3/8" | | | | 20", 22", 24", 25", 26" | | | NC |
| TI08 | <input type="checkbox"/> | Pneumatic with Airless insert 1 3/8" | | | | 22", 24", 26" | | | OTP |
| TI06 | <input type="checkbox"/> | Pneumatic HP, Low Tread, Puncture resistant, black 1" - SpeedLite | | | | 20", 22", 24", 25", 26" | | | OTP |
| TI07 | <input type="checkbox"/> | Pneumatic HP, Medium Tread, Puncture resistant, black 1" - TrailBlazer | | | | 24", 25", 26" | | | OTP |
| 1- Offered with handrims on each side of the wheelchair; select the outer handrim in the handrims section; 2- Only standard handrim position available; 3- Not available with HD option | | | | | | | | | |
| 4- Not compatible with HP tires ; Only available with aluminum or plastic coated handrim in standard position; 5- Not compatible with Shox tires; 6- The 22" is only available with Mags and Newton One rear wheels | | | | | | | | | |
| Handrims | | | | | | | | | |
| HANDRIM | | | | | | Available | | | |
| HR01 | <input type="checkbox"/> | Aluminum anodized | | | | 20", 22", 24", 25", 26" | | | NC |
| HR61 | <input type="checkbox"/> | Aluminum black hard anodized ² | | | | 20", 22", 24", 25" | | | OTP |
| HR02 | <input type="checkbox"/> | Plastic coated ¹ | | | | 20", 22", 24", 25", 26" | | | OTP |
| HR03 | <input type="checkbox"/> | High friction coated ^{1,2, 3} | | | | 22", 24", 25", 26" | | | OTP |
| HR07-HR18 and HR65-HR71 | <input type="checkbox"/> | Newton Air Grip ² | | | | 20", 22", 24", 25", 26" | | | J |
| | | See page 5 for the color selection and price | | | | | | | |
| HR28 | <input type="checkbox"/> | Surge ^{1,2,3,4} (1 1/2" Oval handrim w/ Gripton high friction strip) | | | | 22", 24", 25", 26" | | | OTP |
| HR05 | <input type="checkbox"/> | Surge LT ^{1,2,3,4} (1 1/8" Oval handrim w/ Gripton high friction strip) | | | | 22", 24", 25", 26" | | | OTP |
| HR04 | <input type="checkbox"/> | Natural Fit ^{1,2,3,4} | | | | 22", 24", 25", 26" | | | OTP |
| HR62 | <input type="checkbox"/> | Simi H ^{1,2,3} (1 1/4" Oval handrim with high friction strip) | | | | 22", 24", 25", 26" | | | OTP |
| HR63 | <input type="checkbox"/> | Nova H ^{1,2,3} (1 1/4" Oval handrim with high friction strip) | | | | 22", 24", 25", 26" | | | OTP |
| HR64 | <input type="checkbox"/> | Optimum H ^{1,2,3} (1 1/8" heart-shaped handrim) | | | | 22", 24", 25", 26" | | | OTP |
| HR99 | | Omit handrim | | | | | | | N/A |

1- Only standard and narrow handrim position. By default, the tab won't be cut with the narrow position; 2- Not available with Mag wheels; 3 - Not available with one arm drive;
4- 25" and 26" are not available in narrow position; 24" is not available with the Spinerigy Spox, LX and CLX in narrow position; CLX rear wheel are not available with thumb piece

Handrims option

ASSEMBLY POSITION

| | | | | | | |
|-------|--------------------------|---|--|--|--|-----|
| HR29 | <input type="checkbox"/> | Standard Handrim position | Full tab | | | NC |
| HR030 | <input type="checkbox"/> | Narrow Handrim position no cut | Full tab | | | J |
| HR30 | <input type="checkbox"/> | Narrow Handrim position with cut ⁵ | Short tab 2 positions | | | OTP |
| HR31 | <input type="checkbox"/> | Super Narrow Handrim ⁵ | Short tab 2 positions STD <input type="checkbox"/> Short tab 1 position <input type="checkbox"/> | | | J |

THUMB PIECE SELECTION FOR NATURAL FIT, SURGE, SURGE LT HANDRIMS

| | | | | | |
|-------|--------------------------|------------------|--|--|-----|
| HR23 | <input type="checkbox"/> | Standard grip | | | OTP |
| HR27 | <input type="checkbox"/> | Super grip | | | OTP |
| HR099 | <input type="checkbox"/> | Omit thumb piece | | | NC |

5- Only available with Aluminum and Air Grip handrims;

AXLES (0° camber) ¹


| | | | | | |
|------|--------------------------|--------------------|--|--|-----|
| AX01 | <input type="checkbox"/> | Quick release axle | | | NC |
| AX02 | <input type="checkbox"/> | Permanent axle | | | J |
| AX03 | <input type="checkbox"/> | Quad release axle | | | OTP |

AXLE PLATE

| | | | | | |
|------|--------------------------|--|--|--|-----|
| AX04 | <input type="checkbox"/> | Amputee adjustable axle plate ² | | | OTP |
|------|--------------------------|--|--|--|-----|

1- Not required with one arm drive 2- May limit seat to floor height available

CENTER OF GRAVITY



Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

Default setting is at 3/4"

Most stable

Most responsive

AMPUTEE ¹

3/4"

2"

☐

☐

☐

1- Only available with amputee axle plate

WHEEL LOCKS

| | | | | | | |
|------|--------------------------|---|--|--|--|-----|
| WL08 | <input type="checkbox"/> | Newton Wheel lock - push to lock | PRODUCT IMPROVED | | | NC |
| WL14 | <input type="checkbox"/> | 6" removable extension handle for push to lock | | | | J |
| WL16 | <input type="checkbox"/> | Aluminum push to lock with extensions | | | | OTP |
| WL02 | <input type="checkbox"/> | Pull to lock | | | | J |
| WL04 | <input type="checkbox"/> | Newton grade aid push to lock (works best with pneumatic or high pressure tire) | | | | OTP |
| WL13 | <input type="checkbox"/> | 6" removable extension handle for pull to lock and grade aid | | | | J |
| WL06 | <input type="checkbox"/> | Attendant Lock (additional set of wheel lock mounted to the back) | | | | OTP |
| WL15 | <input type="checkbox"/> | Unilateral wheel lock - pull to lock ¹ | Left <input type="checkbox"/> Right <input type="checkbox"/> | | | OTP |
| WL99 | <input type="checkbox"/> | Omit wheel locks | | | | N/A |

1- For more information about incompatibilities and restriction, refer to the table at page 7 and the "How-to documents" page of our website:
https://www.motioncomposites.com/en_ca/support-and-education/how-to-documents

Armrests

HEIGHT ADJUSTABLE FLIP BACK ARMREST (Convertible to "T" armrest) ^{1,2}

| | | | | | | | |
|-----------|--------------------------|------------------------|-----------|---|--|--|----|
| AR01-812 | <input type="checkbox"/> | 10" Desk length armpad | 8" - 12" | Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/> | | | NC |
| AR01-1014 | <input type="checkbox"/> | 10" Desk length armpad | 10" - 14" | Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/> | | | NC |
| AR02-812 | <input type="checkbox"/> | 14" Full length armpad | 8" - 12" | Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/> | | | NC |
| AR02-1014 | <input type="checkbox"/> | 14" Full length armpad | 10" - 14" | Locked <input type="checkbox"/> Unlocked STD <input type="checkbox"/> | | | NC |

HEIGHT ADJUSTABLE "T" ARMREST ¹

| | | | | | | | |
|-----------|--------------------------|-------------------------------------|-----------|--|--|--|----|
| AR13-610 | <input type="checkbox"/> | 10" Desk length armpad ⁴ | 6" - 10" | | | | NC |
| AR03-812 | <input type="checkbox"/> | 10" Desk length armpad | 8" - 12" | | | | NC |
| AR03-1014 | <input type="checkbox"/> | 10" Desk length armpad | 10" - 14" | | | | NC |
| AR14-610 | <input type="checkbox"/> | 14" Full length armpad ⁴ | 6" - 10" | | | | NC |
| AR04-812 | <input type="checkbox"/> | 14" Full length armpad | 8" - 12" | | | | NC |
| AR04-1014 | <input type="checkbox"/> | 14" Full length armpad | 10" - 14" | | | | NC |

ARMREST OPTIONS

| | | | | | | |
|------|--------------------------|--|--|--|--|-----|
| AR12 | <input type="checkbox"/> | 14" Full Length armpad - Gel Ovation (pair) ³ | | | | OTP |
| AR99 | <input type="checkbox"/> | Omit armrest | | | | NC |

1- May have 1/2" adjust. restriction on both limits of the range; 2- Not available with attendant lock; 3- Not available with short flip back armrest and short "T" armrest; 4- May have restriction with seat-to-floor height and rear wheel

MOVE Orderform Manitoba 6.3

February 5th, 2024

4 of 6

| POSITIONING BELTS | | | | | | | |
|-------------------|--------------------------|---|--|--|--|--|-----|
| SB01-S | <input type="checkbox"/> | Velcro adjustable - Standard 35" ¹ | | | | | OTP |
| SB01-L | <input type="checkbox"/> | Velcro adjustable - Long 60" | | | | | OTP |
| SB02-S | <input type="checkbox"/> | Auto buckle - Standard 35" ¹ | | | | | NC |
| SB02-M | | Auto buckle - Medium 48" | | | | | NC |
| SB02-L | <input type="checkbox"/> | Auto buckle - Long 60" | | | | | NC |
| SB04-M | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point not padded - Medium 56" | | | | | NC |
| SB04-L | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point not padded - Long 61" | | | | | OTP |
| SB05-S | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point padded - Small 51" | | | | | OTP |
| SB05-M | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point padded - Medium 56" | | | | | OTP |
| SB05-L | <input type="checkbox"/> | BodyPoint Hip Belt - 2 point padded - Long 61" | | | | | OTP |
| SB06-S | <input type="checkbox"/> | BodyPoint Hip Belt - 4 point padded - Small 51" | | | | | OTP |
| SB06-M | <input type="checkbox"/> | BodyPoint Hip Belt - 4 point padded - Medium 56" | | | | | OTP |
| SB06-L | <input type="checkbox"/> | BodyPoint Hip Belt - 4 point padded - Long 61" | | | | | OTP |
| SB03 | <input type="checkbox"/> | Hardware attachment for 2 extra points (1 pair) | | | | | OTP |
| SB08-M | <input type="checkbox"/> | BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Medium 18" | | | | | OTP |
| SB08-L | <input type="checkbox"/> | BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps - Long 24" | | | | | OTP |
| SB99 | <input type="checkbox"/> | Omit positioning belt | | | | | NC |




















1- Not available with wheelchair width of 18" and more

| ACCESSORIES | | | | | | | |
|-------------|-------------------------------------|--|-------------------------------|--------------------------------|------------------------------|--------------------------------|---------------------------------|
| AC03 | <input checked="" type="checkbox"/> | Anti-tippers | | | | | NC |
| AC29 | <input type="checkbox"/> | Swing-away anti-tipper ¹ | Left <input type="checkbox"/> | Right <input type="checkbox"/> | | | OTP |
| AC27 | <input type="checkbox"/> | Tip assist ² | | | | | OTP |
| AC06 | <input type="checkbox"/> | Cane and crutch holder | Left <input type="checkbox"/> | Right <input type="checkbox"/> | | | OTP |
| AC07 | <input type="checkbox"/> | Oxygen tank holder ⁴ - D and E Type - Installed to the left | | | | | OTP |
| AC08 | <input type="checkbox"/> | Spoke guards ³ | Available: 20", 22", 24", 26" | | | | OTP |
| AC09 | <input type="checkbox"/> | Footrest neoprene impact guard (<i>pair</i>) - Black | | | | | OTP |
| AC10 | <input type="checkbox"/> | Seat pouch | | | | | OTP |
| AC30 | <input type="checkbox"/> | Newton multi-purpose pouch | | | | | OTP |
| AC26 | <input type="checkbox"/> | BodyPoint Mobility Bag | | | | | OTP |
| AC23 | <input type="checkbox"/> | Newton Backpack | | | | | OTP |
| AC12 | <input type="checkbox"/> | Velcro Closing Strap | | | | | OTP |
| AC13 | <input type="checkbox"/> | Tool kit | | | | | OTP |
| AC16 | <input type="checkbox"/> | Motion Composites Polo Shirt | W/S <input type="checkbox"/> | W/M <input type="checkbox"/> | W/L <input type="checkbox"/> | | OTP |
| | | | M/S <input type="checkbox"/> | M/M <input type="checkbox"/> | M/L <input type="checkbox"/> | M/XL <input type="checkbox"/> | |
| | | | | | | M/XXL <input type="checkbox"/> | M/XXXL <input type="checkbox"/> |

1- Must select the pair (left and right) with HD option; Only available with rear seat to floor height between 13 1/2" and 19 3/4"


2- Cannot be combined with Anti Tippers

3- Not available with Mag wheels and Spinergy; 4- Not compatible with swing away armrest, amputee axle plate and WC-19 transit tie-down

| Colors | | | | | | | | | | | | | | | | | | | | |
|----------------------|---|---|---|---|---|---|---|---|---|---|---|---|--|---|---|---|---|---|---|-----|
| Options ¹ | Black | Sunkissed Orange | Burgundy | Sapphire Blue | Monster Green | Charcoal | Fuchsia | Steel Blue | Ferrari Red | Acid Green | Electric Purple | Caribbean Blue | White | Pink Pearl | Plum Purple | Sky Blue | Black Cherry | Khaki Green (matte) | Desert Sand (matte) | |
| |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  |  | |
| AirGrip | <input type="checkbox"/> J | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | OTP |
| Frame | <input type="checkbox"/> | | | | | | | | | | | | | | | | | | | NC |

1- Color representation may vary from one monitor to another. Therefore, the actual paint color on your wheelchair might differ slightly from what you see on your screen.; Finish between Airgrip and frame might be different;

Fork and Stem bolt combinations



Unilateral wheel lock restrictions

| | | |
|---|-----------------|------------------------------|
| Seat depth: shall be 16 in and greater | | |
| Rear wheel: shall be 24 in and less | | |
| "T" or flip-back armrest only compatible with 0° Camber | | |
| One arm drive: Not compatible with "T" or flip-back armrest | | |
| "T" or flip-back armrest may require to move the rear wheels outwards | | |
| Not compatible with wheel lock lever extension | | |
| Rear wheel diameter | Seat depth (in) | COG Compatibility Range (in) |
| 20 in | 16+ | 2 |
| 22 in | 16 | 0.75 to 2 |
| | 17+ | 0.75 to 2 |
| 24 in | 16 | 0.75 |
| | 17 | 0.75 to 2 |
| | 18+ | 0.75 to 2 |

Notes :