



160 Armand-Majeau Sud  
 Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0  
 T. +1-866-650-6555 F. +1-888-966-6555  
[quotes@motioncomposites.com](mailto:quotes@motioncomposites.com)  
[orders@motioncomposites.com](mailto:orders@motioncomposites.com)  
[www.motioncomposites.com](http://www.motioncomposites.com)



<b>Move</b>	<b>MVQs</b>
<b>Quick Ship Order Form</b>	
Canada	6.1

P.O. number:\* \_\_\_\_\_ Date: \_\_\_\_\_

Company Name:\* \_\_\_\_\_

Purchasing contact:\* \_\_\_\_\_ Phone:\* \_\_\_\_\_

Clinician/Sales rep: \_\_\_\_\_

Email: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to:\* \_\_\_\_\_

Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ Province \_\_\_\_\_

Postal Code: \_\_\_\_\_

Phone: \_\_\_\_\_

<b>MOVE</b>				<i>Taxes not included</i>	
Ultralight MOVE wheelchair aluminum	Basalt Black	WM0001196		\$1 425	<input type="checkbox"/>
Included: seating upholstery / Excluded back upholstery and cushion	265 lbs. weight capacity				

### Front Seat-to-Floor Height Configuration

13"				14"				16 1/2"							
M6 16x16	NC	<input type="checkbox"/>		M1 18x16	NC	<input type="checkbox"/>		M11 16x18	NC	<input type="checkbox"/>					
M7 16x18	NC	<input type="checkbox"/>		M2 18x18	NC	<input type="checkbox"/>		M12 18x18	NC	<input type="checkbox"/>					
M9 18x18	NC	<input type="checkbox"/>		M3 18x20	WAMB	\$235	<input type="checkbox"/>	M13 18x20	WAMB	\$235	<input type="checkbox"/>				
M10 20x18	WAMA	\$315	<input type="checkbox"/>	M4 20x20	WAMA WAMB	\$550	<input type="checkbox"/>	M14 20x18	WAMA	\$315	<input type="checkbox"/>				
WHEELS Front, Rear & Handrim				WHEELS Front, Rear & Handrim				WHEELS Front, Rear & Handrim							
4"x1" Composites caster w/PU tire	NC	<input type="checkbox"/>		5"x1" Composites caster w/PU tire	NC	<input type="checkbox"/>		6"x1" Composites caster w/PU tire	NC	<input type="checkbox"/>					
5" fork + Std stem bolt	NC	<input type="checkbox"/>		5" fork + Std stem bolt	NC	<input type="checkbox"/>		7" fork + Std stem bolt	NC	<input type="checkbox"/>					
22" Newton One rear wheel	NC	<input type="checkbox"/>		22" Newton One rear wheel	NC	<input type="checkbox"/>		24" Newton One rear wheel	NC	<input type="checkbox"/>					
Soft urethane tire 1 <sup>3/8</sup> "	NC	<input type="checkbox"/>		Soft urethane tire 1 <sup>3/8</sup> "	NC	<input type="checkbox"/>		Soft urethane tire 1 <sup>3/8</sup> "	NC	<input type="checkbox"/>					
Aluminum anodized handrim	NC	<input type="checkbox"/>		Aluminum anodized handrim	NC	<input type="checkbox"/>		Aluminum anodized handrim	NC	<input type="checkbox"/>					
SWING-AWAY FOOTREST				SWING-AWAY FOOTREST				SWING-AWAY FOOTREST							
60° Swing-away footrests	\$270	<input type="checkbox"/>		70° Swing-away footrests	\$270	<input type="checkbox"/>		70° Swing-away footrests	\$270	<input type="checkbox"/>					
Newton Composite Adjustable angle footplate w/ heel loop	WAMJ	\$181	<input type="checkbox"/>	Newton Composite Adjustable angle footplate w/ heel loop	WAMJ	\$181	<input type="checkbox"/>	Newton Composite Adjustable angle footplate w/ heel loop	WAMJ	\$181	<input type="checkbox"/>				
This configuration:				This configuration:				This configuration:							
Front STF Height	Rear STF Height			Front STF Height	Rear STF Height			Front STF Height	Rear STF Height						
13" - 14"	13" - 18"			13 3/4" - 14 3/4"	13" - 18"			15 1/2" - 17 1/2"	14" - 19"						
Possible configurations:															
Front STF Height				12 1/2" - 19 1/2"				Rear STF Height				12 1/4" - 20"			

PARTS & ACCESSORIES on all configurations			
<b>8° BEND BACK</b> (with fixed height push handles)			
Aluminum back cane	20" height	NC	<input type="checkbox"/>
<b>FLIP BACK ARMREST</b> (Convertible to "T" armrest)			
14" Full Length, 8 - 12" height adjustable		\$395	<input type="checkbox"/>
<b>AXLE</b> (with 0° camber)			
Quick release axle		WAMU	NC <input type="checkbox"/>
<b>WHEEL LOCK</b>			
Push to lock with 6" removable extension handle		\$73	<input type="checkbox"/>
<b>POSITIONING BELT</b>			
Additional attachment hardware for 4 point seatbelt		\$15	<input type="checkbox"/>
<b>ACCESSORIES</b>			
Anti tippers		\$120	<input type="checkbox"/>

# Quick Ship Parts and Accessories Program

Motion Composites' Quick Ship Parts and Accessories Program adds flexibility to the convenience of our Quick Ship wheelchairs

These specially discounted parts will be added to your Quick Ship wheelchair order for you to install and customize upon receipt

Invoice code					Invoice code				
QP0004	Aluminum straight back cane 18-21"	PR	\$120	<input type="checkbox"/>	QP0033	Standard stem bolt (anti-flutter)	PR	\$21	<input type="checkbox"/>
QP0005	Aluminum back cane 8° bend 18-21"	PR	\$140	<input type="checkbox"/>	QP0034	+1" stem bolt (anti-flutter)	PR	\$21	<input type="checkbox"/>
QP0006	Angle adjustable - Aluminium straight back cane 18-21"	PR	\$240	<input type="checkbox"/>	QP0035	+2" stem bolt (anti-flutter)	PR	\$21	<input type="checkbox"/>
QP0098	Push handles	PR	\$110	<input type="checkbox"/>	QP0076	Length adjustable heel loop	PR	\$41	<input type="checkbox"/>
QP0134	Slip-On Nylon Back UPH, 16x18-21in	EA	\$67	<input type="checkbox"/>	QP0096	Newton Airgrip 22" black handrim	PR	\$160	<input type="checkbox"/>
QP0135	Slip-On Nylon Back UPH, 18x18-21in	EA	\$67	<input type="checkbox"/>	QP0097	Newton Airgrip 24" black handrim	PR	\$160	<input type="checkbox"/>
QP0136	Slip-On Nylon Back UPH, 20x18-21in	EA	\$67	<input type="checkbox"/>	QP0036	4" fork (2 holes)	PR	\$62	<input type="checkbox"/>
QP0137	Adj Tension Nylon Back UPH, 16x18-21in	EA	\$260	<input type="checkbox"/>	QP0037	5" fork (4 holes)	PR	\$62	<input type="checkbox"/>
QP0138	Adj Tension Nylon Back UPH, 18x18-21in	EA	\$260	<input type="checkbox"/>	QP0038	7" fork (5 holes)	PR	\$62	<input type="checkbox"/>
QP0139	Adj Tension Nylon Back UPH, 20x18-21in	EA	\$260	<input type="checkbox"/>	QP0114	Newton One wheels w/Urethane tire & Alu Pushrim - 22"	PR	\$450	<input type="checkbox"/>
QP0167	Swing Out residual limb support flat cushion - L 6"x8"	E1020 EA	\$505	<input type="checkbox"/>	QP0115	Newton One wheels w/Urethane tire & Alu Pushrim - 24"	PR	\$450	<input type="checkbox"/>
QP0168	Swing Out residual limb support flat cushion - R 6"x8"	E1020 EA	\$505	<input type="checkbox"/>	QP0039	Mag wheels w/Urethane tire & Alu Pushrim - 20"	PR	\$450	<input type="checkbox"/>
QP0169	Swing Out residual limb support flat cushion - L 8"x10"	E1020 EA	\$505	<input type="checkbox"/>	QP0040	Mag wheels w/Urethane tire & Alu Pushrim - 22"	PR	\$450	<input type="checkbox"/>
QP0170	Swing Out residual limb support flat cushion - R 8"x10"	E1020 EA	\$505	<input type="checkbox"/>	QP0041	Mag wheels w/Urethane tire & Alu Pushrim - 24"	PR	\$450	<input type="checkbox"/>
QP0171	Swing Out residual limb support flat cushion - L 10"x10"	E1020 EA	\$505	<input type="checkbox"/>	QP0044	Wheel lock - pull to lock	PR	\$105	<input type="checkbox"/>
QP0172	Swing Out residual limb support flat cushion - R 10"x10"	E1020 EA	\$505	<input type="checkbox"/>	QP0045	Wheel lock - pull to lock with removable extension handle	PR	\$175	<input type="checkbox"/>
QP0173	Swing Out residual limb support flat cushion - L 10"x14"	E1020 EA	\$505	<input type="checkbox"/>	QP0046	Wheel lock - Newton grade aid push to lock	PR	\$305	<input type="checkbox"/>
QP0174	Swing Out residual limb support flat cushion - R 10"x14"	E1020 EA	\$505	<input type="checkbox"/>	QP0048	Positioning belt - velcro adjustable	EA	\$50	<input type="checkbox"/>
QP0175	Swing Out residual limb support contour cushion - L 6"x8"	E1020 EA	\$505	<input type="checkbox"/>	QP0049	Positioning belt - auto buckle	EA	\$60	<input type="checkbox"/>
QP0176	Swing Out residual limb support contour cushion - R 6"x8"	E1020 EA	\$505	<input type="checkbox"/>	QP0050	Bodypoint 2 point Non Padded - medium (56")	EA	\$90	<input type="checkbox"/>
QP0177	Swing Out residual limb support contour cushion - L 8"x10"	E1020 EA	\$505	<input type="checkbox"/>	QP0051	Bodypoint 2 point Non Padded - large (61")	EA	\$90	<input type="checkbox"/>
QP0178	Swing Out residual limb support contour cushion - R 8"x10"	E1020 EA	\$505	<input type="checkbox"/>	QP0052	Bodypoint 2 point Padded - small (51")	EA	\$165	<input type="checkbox"/>
QP0149	Swing Away ELR/calf pad/std FP, Left & heel loop (14-16in)	EA	\$410	<input type="checkbox"/>	QP0053	Bodypoint 2 point Padded - medium (56")	EA	\$165	<input type="checkbox"/>
QP0150	Swing Away ELR/calf pad/std FP, Left & heel loop (17-18in)	EA	\$410	<input type="checkbox"/>	QP0054	Bodypoint 2 point Padded - large (61")	EA	\$165	<input type="checkbox"/>
QP0151	Swing Away ELR/calf pad/std FP, Left & heel loop (19in +)	EA	\$410	<input type="checkbox"/>	QP0055	Bodypoint 4 point Padded - small (51")	EA	\$205	<input type="checkbox"/>
QP0152	Swing Away ELR/calf pad/std FP, Right & heel loop (14-16in)	EA	\$410	<input type="checkbox"/>	QP0056	Bodypoint 4 point Padded - medium (56")	EA	\$205	<input type="checkbox"/>
QP0153	Swing Away ELR/calf pad/std FP, Right & heel loop (17-18in)	EA	\$410	<input type="checkbox"/>	QP0057	Bodypoint 4 point Padded - large (61")	EA	\$205	<input type="checkbox"/>
QP0154	Swing Away ELR/calf pad/std FP, Right & heel loop (19in +)	EA	\$410	<input type="checkbox"/>	QP0058	Cane carrier - Left	EA	\$115	<input type="checkbox"/>
QP0013	Calf strap (wheelchair width: 15 to 17 in)	EA	\$57	<input type="checkbox"/>	QP0059	Cane carrier - Right	EA	\$115	<input type="checkbox"/>
QP0014	Calf strap (wheelchair width: 18 to 20 in)	EA	\$57	<input type="checkbox"/>	QP0060	Spoke guards - 22"	PR	\$155	<input type="checkbox"/>
QP0125	Bodypoint Aeromesh Padded calf strap (wheelchair Width: 12 to 16 in)	EA	\$140	<input type="checkbox"/>	QP0061	Spoke guards - 24"	PR	\$155	<input type="checkbox"/>
QP0126	Bodypoint Aeromesh Padded calf strap (wheelchair width: 17 to 20 in)	EA	\$140	<input type="checkbox"/>	QP0062	Footrest neoprene impact guard (pr) - Black	PR	\$85	<input type="checkbox"/>
QP0127	Bodypoint Aeromesh Padded calf strap (wheelchair width: 21 to 24 in)	EA	\$140	<input type="checkbox"/>	QP0063	Seat pouch	EA	\$105	<input type="checkbox"/>
QP0017	Composite caster 3" x 1"	PR	\$85	<input type="checkbox"/>	QP0103	Newton Backpack	EA	\$130	<input type="checkbox"/>
QP0018	Composite caster 4" x 1"	PR	\$85	<input type="checkbox"/>	QP0106	Motion Composites Polo-Shirt Men small	EA	\$36	<input type="checkbox"/>
QP0019	Composite caster 5" x 1"	PR	\$100	<input type="checkbox"/>	QP0107	Motion Composites Polo-Shirt Men medium	EA	\$36	<input type="checkbox"/>
QP0020	Composite caster 6" x 1"	PR	\$100	<input type="checkbox"/>	QP0108	Motion Composites Polo-Shirt Men large	EA	\$36	<input type="checkbox"/>
QP0021	Composite caster 7" x 1"	PR	\$110	<input type="checkbox"/>	QP0109	Motion Composites Polo-Shirt Men X-large	EA	\$36	<input type="checkbox"/>
QP0022	Composite caster 8" x 1"	PR	\$110	<input type="checkbox"/>	QP0110	Motion Composites Polo-Shirt Men XX-large	EA	\$36	<input type="checkbox"/>
QP0077	Newton Ultracaster Composite caster 4" x 1½" Soft-Roll	PR	\$105	<input type="checkbox"/>	QP0118	Motion Composites Polo-Shirt Men XXX-large	EA	\$36	<input type="checkbox"/>
QP0078	Newton Ultracaster Composite caster 5" x 1½" Soft-Roll	PR	\$105	<input type="checkbox"/>	QP0111	Motion Composites Polo-Shirt Women small	EA	\$36	<input type="checkbox"/>
QP0079	Newton Ultracaster Composite caster 6" x 1½" Soft-Roll	PR	\$105	<input type="checkbox"/>	QP0112	Motion Composites Polo-Shirt Women medium	EA	\$36	<input type="checkbox"/>
QP0080	Newton Ultracaster Aluminum caster 4" x 1½" Soft-Roll	PR	\$185	<input type="checkbox"/>	QP0113	Motion Composites Polo-Shirt Women large	EA	\$36	<input type="checkbox"/>
QP0081	Newton Ultracaster Aluminum caster 5" x 1½" Soft-Roll	PR	\$195	<input type="checkbox"/>	QP0071	Tool Kit	EA	\$85	<input type="checkbox"/>
QP0082	Newton Ultracaster Aluminum caster 6" x 1½" Soft-Roll	PR	\$205	<input type="checkbox"/>					
QP0031	Pneumatic caster 6" x 1½"	PR	\$170	<input type="checkbox"/>					
QP0032	Pneumatic caster 8" x 1½"	PR	\$170	<input type="checkbox"/>					

### Move REAR Seat-to-floor Height

20"	22"	24"	25"	26"
12 1/4"	13"	14"	14 1/4"	15"
13 1/4"	14"	15"	15 1/4"	16"
14 1/4"	15"	16"	16 1/4"	17"
15 1/4"	16"	17"	17 1/4"	18"
16 1/4"	17"	18"	18 1/4"	19"
17 1/4"	18"	19"	19 1/4"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- 1/2" of selected seat to floor height"

Table 1

SLIP ON BACK UPHOLSTERY		LENGTH inches (cm)
		18-21 (45.7-53.3)
Width Inches (cm)	16 (40.6)	9001832
	18 (45.7)	9001842
	20 (50.8)	9001852

### Move FRONT Seat-to-floor Height

		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
Fork/Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2	13 1/2	14 1/2												
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1															
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1							15 1/2	16 1/2	17 1/2	16	17	18			
	P2				15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3				16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +/- 1/4"

Table 2

ADJUSTABLE TENSION BACK UPHOLSTERY		LENGTH inches (cm)
		18-21 (45.7-53.3)
Width Inches (cm)	16 (40.6)	9001777
	18 (45.7)	9001787
	20 (50.8)	9001797

### Notes :