



160 Armand-Majeau Sud
Saint-Roch-de-l'Achigan, QC, Canada J0K 3H0
T. +1-866-650-6555 F. +1-888-966-6555
quotes@motioncomposites.com
orders@motioncomposites.com
www.motioncomposites.com



Move

MVQs

Quick Ship Order Form

Canada

6.4

P.O. number: *

Company Name: *

Purchasing contact: *

Phone: *

Clinician/Sales rep:

Email:

Marked for:

* Required

Ship to: *

Name:

Address:

City:

Province

Postal Code:

Phone:

Taxes not included

MOVE						
Ultralight MOVE wheelchair aluminum	Basalt Black	WM0001196		\$1,425	☐	
Included: seating upholstery / Excluded back upholstery and cushion	265 lbs. weight capacity					

Front Seat-to-Floor Height Configuration

13"				14"				16 1/2"			
M6	16x16	NC	☐	M1	18x16	NC	☐	M11	16x18	NC	☐
M7	16x18	NC	☐	M2	18x18	NC	☐	M12	18x18	NC	☐
M9	18x18	NC	☐	M3	18x20	WAMB \$235	☐	M13	18x20	WAMB \$235	☐
M10	20x18	WAMA \$315	☐	M4	20x20	WAMA WAMB \$550	☐	M14	20x18	WAMA \$315	☐
WHEELS Front, Rear & Handrim				WHEELS Front, Rear & Handrim				WHEELS Front, Rear & Handrim			
4"x1" Composites caster w/PU tire	NC	☐		5"x1" Composites caster w/PU tire	NC	☐		6"x1" Composites caster w/PU tire	NC	☐	
5" fork + Std stem bolt	NC	☐		5" fork + Std stem bolt	NC	☐		7" fork + Std stem bolt	NC	☐	
22" Newton One rear wheel	NC	☐		22" Newton One rear wheel	NC	☐		24" Newton One rear wheel	NC	☐	
Soft urethane tire 1 3/8"	NC	☐		Soft urethane tire 1 3/8"	NC	☐		Soft urethane tire 1 3/8"	NC	☐	
Aluminum anodized handrim	NC	☐		Aluminum anodized handrim	NC	☐		Aluminum anodized handrim	NC	☐	
SWING-AWAY FOOTREST				SWING-AWAY FOOTREST				SWING-AWAY FOOTREST			
60° Swing-away footrests	\$270	☐		70° Swing-away footrests	\$270	☐		70° Swing-away footrests	\$270	☐	
Newton Composite Adjustable angle footplate w/ heel loop	WAMJ \$181	☐		Newton Composite Adjustable angle footplate w/ heel loop	WAMJ \$181	☐		Newton Composite Adjustable angle footplate w/ heel loop	WAMJ \$181	☐	
This configuration:				This configuration:				This configuration:			
Front STF Height	Rear STF Height			Front STF Height	Rear STF Height			Front STF Height	Rear STF Height		
13" - 14"	13" - 18"			13 3/4" - 14 3/4"	13" - 18"			15 1/2" - 17 1/2"	14" - 19"		
Possible configurations:											
Front STF Height		12 1/2" - 19 1/2"				Rear STF Height		12 1/4" - 20"			

PARTS & ACCESSORIES on all configurations			
8° BEND BACK (with fixed height push handles)			
Aluminum back cane	20" height	NC	☐
FLIP BACK ARMREST (Convertible to "T" armrest)			
14" Full Length, 8 - 12" height adjustable		\$395	☐
AXLE (with 0° camber)			
Quick release axle		WAMU NC	☐
WHEEL LOCK			
Push to lock with 6" removable extension handle		\$73	☐
POSITIONING BELT			
Additional attachment hardware for 4 point seatbelt		\$15	☐
ACCESSORIES			
Anti tippers		\$120	☐

Quick Ship Parts and Accessories Program

Motion Composites' Quick Ship Parts and Accessories Program adds flexibility to the convenience of our Quick Ship wheelchairs

These specially discounted parts will be added to your Quick Ship wheelchair order for you to install and customize upon receipt

Invoice code					Invoice code				
QP0004	Aluminum straight back cane 18-21"	PR	\$120	☐	QP0033	Standard stem bolt (anti-flutter)	PR	\$21	☐
QP0005	Aluminum back cane 8" bend 18-21"	PR	\$140	☐	QP0034	+1" stem bolt (anti-flutter)	PR	\$21	☐
QP0006	Angle adjustable - Aluminium straight back cane 18-21"	PR	\$240	☐	QP0035	+2" stem bolt (anti-flutter)	PR	\$21	☐
QP0098	Push handles	PR	\$110	☐	QP0076	Length adjustable heel loop	PR	\$41	☐
QP0134	Slip-On Nylon Back UPH, 16x18-21in	EA	\$67	☐	QP0096	Newton Airgrip 22" black handrim	PR	\$160	☐
QP0135	Slip-On Nylon Back UPH, 18x18-21in	EA	\$67	☐	QP0097	Newton Airgrip 24" black handrim	PR	\$160	☐
QP0136	Slip-On Nylon Back UPH, 20x18-21in	EA	\$67	☐	QP0036	4" fork (2 holes)	PR	\$62	☐
QP0137	Adj Tension Nylon Back UPH, 16x18-21in	EA	\$260	☐	QP0037	5" fork (4 holes)	PR	\$62	☐
QP0138	Adj Tension Nylon Back UPH, 18x18-21in	EA	\$260	☐	QP0038	7" fork (5 holes)	PR	\$62	☐
QP0139	Adj Tension Nylon Back UPH, 20x18-21in	EA	\$260	☐	QP0114	Newton One wheels w/Urethane tire & Alu Pushrim - 22"	PR	\$450	☐
QP0167	Swing Out residual limb support flat cushion - L 6"x8" E1020	EA	\$505	☐	QP0115	Newton One wheels w/Urethane tire & Alu Pushrim - 24"	PR	\$450	☐
QP0168	Swing Out residual limb support flat cushion - R 6"x8" E1020	EA	\$505	☐	QP0039	Mag wheels w/Urethane tire & Alu Pushrim - 20"	PR	\$450	☐
QP0169	Swing Out residual limb support flat cushion - L 8"x10" E1020	EA	\$505	☐	QP0040	Mag wheels w/Urethane tire & Alu Pushrim - 22"	PR	\$450	☐
QP0170	Swing Out residual limb support flat cushion - R 8"x10" E1020	EA	\$505	☐	QP0041	Mag wheels w/Urethane tire & Alu Pushrim - 24"	PR	\$450	☐
QP0171	Swing Out residual limb support flat cushion - L 10"x10 E1020	EA	\$505	☐	QP0044	Wheel lock - pull to lock	PR	\$105	☐
QP0172	Swing Out residual limb support flat cushion - R 10"x10 E1020	EA	\$505	☐	QP0045	Wheel lock - pull to lock with removable extension handle	PR	\$175	☐
QP0173	Swing Out residual limb support flat cushion - L 10"x14 E1020	EA	\$505	☐	QP0046	Wheel lock - Newton grade aid push to lock	PR	\$305	☐
QP0174	Swing Out residual limb support flat cushion - R 10"x14 E1020	EA	\$505	☐	QP0048	Positioning belt - velcro adjustable - 35"	EA	\$50	☐
QP0175	Swing Out residual limb support contour cushion - L 6"x8" E1020	EA	\$505	☐	QP0049	Positioning belt - auto buckle - 35"	EA	\$60	☐
QP0176	Swing Out residual limb support contour cushion - R 6"x8" E1020	EA	\$505	☐	QP0195	Positioning belt - auto buckle - 48"	EA	\$60	☐
QP0177	Swing Out residual limb support contour cushion - L 8"x10" E1020	EA	\$505	☐	QP0050	Bodypoint 2 point Non Padded - medium (56")	EA	\$95	☐
QP0178	Swing Out residual limb support contour cushion - R 8"x10" E1020	EA	\$505	☐	QP0051	Bodypoint 2 point Non Padded - large (61")	EA	\$95	☐
QP0149	Swing Away ELR/calf pad/std FP, Left & heel loop (14-16in)	EA	\$410	☐	QP0052	Bodypoint 2 point Padded - small (51")	EA	\$185	☐
QP0150	Swing AwayELR/calf pad/std FP, Left & heel loop (17-18in)	EA	\$410	☐	QP0053	Bodypoint 2 point Padded - medium (56")	EA	\$185	☐
QP0151	SwingAway ELR/calf pad/std FP, Left & heel loop (19in +)	EA	\$410	☐	QP0054	Bodypoint 2 point Padded - large (61")	EA	\$185	☐
QP0152	Swing Away ELR/calf pad/std FP, Right & heel loop (14-16in)	EA	\$410	☐	QP0055	Bodypoint 4 point Padded - small (51")	EA	\$210	☐
QP0153	Swing Away ELR/calf pad/std FP, Right & heel loop (17-18in)	EA	\$410	☐	QP0056	Bodypoint 4 point Padded - medium (56")	EA	\$210	☐
QP0154	Swing Away ELR/calf pad/std FP, Right & heel loop (19in +)	EA	\$410	☐	QP0057	Bodypoint 4 point Padded - large (61")	EA	\$210	☐
QP0013	Calf strap (wheelchair width: 15 to 17 in)	EA	\$57	☐	QP0058	Cane carrier - Left	EA	\$115	☐
QP0014	Calf strap (wheelchar width: 18 to 20 in)	EA	\$57	☐	QP0059	Cane carrier - Right	EA	\$115	☐
QP0125	Bodypoint Aeromesh Padded calf strap (wheelchair Width: 12 to 16 in)	EA	\$140	☐	QP0060	Spoke guards - 22"	PR	\$155	☐
QP0126	Bodypoint Aeromesh Padded calf strap (wheelchar width: 17 to 20 in)	EA	\$140	☐	QP0061	Spoke guards - 24"	PR	\$155	☐
QP0127	Bodypoint Aeromesh Padded calf strap (wheelchar width: 21 to 24 in)	EA	\$140	☐	QP0062	Footrest neoprene impact guard (pr) - Black	PR	\$85	☐
QP0017	Composite caster 3" x 1"	PR	\$85	☐	QP0063	Seat pouch	EA	\$105	☐
QP0018	Composite caster 4" x 1"	PR	\$85	☐	QP0103	Newton Backpack	EA	\$130	☐
QP0019	Composite caster 5" x 1"	PR	\$100	☐	QP0106	Motion Composites Polo-Shirt Men small	EA	\$36	☐
QP0020	Composite caster 6" x 1"	PR	\$100	☐	QP0107	Motion Composites Polo-Shirt Men medium	EA	\$36	☐
QP0021	Composite caster 7" x 1"	PR	\$110	☐	QP0108	Motion Composites Polo-Shirt Men large	EA	\$36	☐
QP0022	Composite caster 8" x 1"	PR	\$110	☐	QP0109	Motion Composites Polo-Shirt Men X-large	EA	\$36	☐
QP0077	Newton Ultracaster Composite caster 4" x 1½" Soft-Roll	PR	\$105	☐	QP0110	Motion Composites Polo-Shirt Men XX-large	EA	\$36	☐
QP0078	Newton Ultracaster Composite caster 5" x 1½" Soft-Roll	PR	\$105	☐	QP0118	Motion Composites Polo-Shirt Men XXX-large	EA	\$36	☐
QP0079	Newton Ultracaster Composite caster 6" x 1½" Soft-Roll	PR	\$105	☐	QP0111	Motion Composites Polo-Shirt Women small	EA	\$36	☐
QP0080	Newton Ultracaster Aluminum caster 4" x 1½" Soft-Roll	PR	\$198	☐	QP0112	Motion Composites Polo-Shirt Women medium	EA	\$36	☐
QP0081	Newton Ultracaster Aluminum caster 5" x 1½" Soft-Roll	PR	\$198	☐	QP0113	Motion Composites Polo-Shirt Women large	EA	\$36	☐
QP0082	Newton Ultracaster Aluminum caster 6" x 1½" Soft-Roll	PR	\$198	☐	QP0071	Tool Kit	EA	\$85	☐
QP0031	Pneumatic caster 6" x 1¼"	PR	\$170	☐	QP0192	Unoccupied transit tie-down	pr	\$120	☐
QP0032	Pneumatic caster 8" x 1¼"	PR	\$170	☐					

Table 1

SLIP ON BACK UPHOLSTERY		LENGTH inches (cm)
		18-21
		(45.7-53.3)
Width Inches (cm)	16 (40.6)	9001832
	18 (45.7)	9001842
	20 (50.8)	9001852

Move REAR Seat-to-floor Height

20"	22"	24"	25"	26"
12 1/4"	13"	14"	14 1/4"	15"
13 1/4"	14"	15"	15 1/4"	16"
14 1/4"	15"	16"	16 1/4"	17"
15 1/4"	16"	17"	17 1/4"	18"
16 1/4"	17"	18"	18 1/4"	19"
17 1/4"	18"	19"	19 1/4"	20"

Note: Tire size, center of gravity setting and seat slope may result in a variation of +/- ½" of selected seat to floor height"

Table 2

ADJUSTABLE TENSION BACK UPHOLSTERY		LENGTH inches (cm)
		18-21
		(45.7-53.3)
Width Inches (cm)	16 (40.6)	9001777
	18 (45.7)	9001787
	20 (50.8)	9001797

Move FRONT Seat-to-floor Height

		4" caster			5" caster			6" caster			7" caster			8" caster		
		Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length			Stem bolt length		
Fork/Position		STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"	STD	+1"	+2"
4"	P1	12 1/2	13 1/2	14 1/2												
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
5"	P1															
	P2	13	14	15	13 1/2	14 1/2	15 1/2									
	P3	13 1/2	14 1/2	15 1/2	14	15	16	14 1/2	15 1/2	16 1/2						
	P4	14	15	16	14 1/2	15 1/2	16 1/2	15	16	17	15 1/2	16 1/2	17 1/2			
7"	P1							15 1/2	16 1/2	17 1/2	16	17	18			
	P2				15 1/2	16 1/2	17 1/2	16	17	18	16 1/2	17 1/2	18 1/2			
	P3				16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2
	P4	16	17	18	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20
	P5	16 1/2	17 1/2	18 1/2	17	18	19	17 1/2	18 1/2	19 1/2	18	19	20	18 1/2	19 1/2	20 1/2

Height available +/- 1/4"

Notes :