



# RAZ-AP600 ORDER FORM

EFFECTIVE MAY 2026 (CDN MSRP)

Date: \_\_\_\_\_

Dealer / Location: \_\_\_\_\_ T: \_\_\_\_\_ F: \_\_\_\_\_

PO# / Client: \_\_\_\_\_ Contact: \_\_\_\_\_ Email: \_\_\_\_\_

**Raz-AP600** Heavy Duty Mobile Shower Commode Chair with stainless steel frame, flip-up padded arm supports, heavy duty foot supports, four dual-locking 5" (125mm) casters and commode pan. Also included is a 20"H Multi-Position Back with 2" fore/aft adjustment range, two back angles (90°/97°) and tension-adjustable vinyl upholstery. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. User weight over 450 lbs. requires 6" (150mm) casters (ZCUAP6). **User weight capacity 275kgs.**  
For help with justification, please see our [Evaluation and Justification Form](#).

**Need help with this order form? Call customer service at (866) 875-2877 or (519) 942 8441**

## 1. Raz-AP600 codes by size – Must Choose One

Chair Code	Seat Frame Width	MSRP
Z16022	22" W Frame	\$3,615
Z16024	24" W Frame	\$3,925
Z16026	26" W Frame	\$4,260
Z16028	28" W Frame	\$4,600
Z16030	30" W Frame – Requires custom seat	\$5,175

## 2. Seats – Must Choose One

Molded and Visco Foam Interface Seats			
ZMS23FRO	Molded 23"W x 18"D – Front Access	Fits 22" and 24"W <sup>1</sup> frames only	N/C
ZMS26FRO	Molded 26"W x 18"D – Front Access	Fits 26" and 28"W <sup>1</sup> frames only	N/C
ZVFI23FRO	Visco Foam Interface 23"W x 18"D – Front Access	Fits 22" and 24"W <sup>1</sup> frames only	\$835
ZVFI26FRO	Visco Foam Interface 26"W x 18"D – Front Access	Fits 26" and 28"W <sup>1</sup> frames only	\$860
ZSSKASHD	Seat Support Kit for Raz-AP600 – for seats that are narrower than seat frames		\$225

<sup>1</sup> Requires a Seat Support Kit (ZSSKASHD)

### Custom Seats

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. [See Custom Seat Assessment Form – Resources tab at razdesigninc.com](#)

ZCSEATHD	Heavy Duty Custom Seat (enter specs below) <b>Select if both width and depth are ≤ 28"</b>	\$1,545
ZCSEATXHD	Heavy Duty Custom Seat for widths 28" and 2.5" (enter specs below) <b>Select if width and/or depth is ≥ 28"</b>	\$1,745

Seat Frame Width: \_\_\_ 22" / 24" / 26" / 28" / 30"  
 Custom Seat Width: \_\_\_ 22" / 24" / 26" / 28" / 30" / 32" (cannot be more than 2" wider than frame width)  
 Seat Depth: \_\_\_ 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" [≥ 21" requires custom longer frame]  
 Access: \_\_\_ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged  
 Aperture Width: \_\_\_ 6.5" recommended

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 6.5" for custom seat.
- Custom seat widths can be up to 2" wider than the frame width to provide additional thigh support without increasing overall width of the chair.
- For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 0,75" fore / aft adjustment from specified custom seat depth.
- Access for Rear-Open seats is limited due to rear seat support tube.
- Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc.



Motion Composites is the exclusive Canadian Distributor of RAZ DESIGN INC.

[motioncomposites.com](http://motioncomposites.com)

**3. Symphony Back** (Recessed solid back that replaces the standard Tension-Adjustable Vinyl Upholstery [recommended for infection-control environments]). All Symphony Back panels are 13" H.

ZSB22	Symphony Back – Recessed Solid Back Support – 22" W	\$525
ZSB24	Symphony Back – Recessed Solid Back Support – 24" W	\$555
ZSB26	Symphony Back – Recessed Solid Back Support – 26" W	\$580
ZSB28	Symphony Back – Recessed Solid Back Support – 28" W	\$745
ZSB30	Symphony Back – Recessed Solid Back Support – 30" W	\$770

**4. Arm Supports Accessories**

ZALKHD	Arm Support Locks – HD (pr)	\$110
ZFAPLL	Long Molded Flat Armpad <sup>1</sup> – Left (14"L x 3.25"W)	\$95
ZFAPLR	Long Molded Flat Armpad <sup>1</sup> – Right (14"L x 3.25"W)	\$95
ZTRML	Arm Trough <sup>2</sup> - Medium Left (18.5" L x 5" W)	\$130
ZTRMR	Arm Trough <sup>2</sup> - Medium Right (18.5" L x 5" W)	\$130
ZTRSL	Arm Trough <sup>2</sup> - Large Left (21" L x 5.63.5cm W) <sup>NEW</sup>	\$130
ZTRSR	Arm Trough <sup>2</sup> - Large Right (21" L x 5.63.5cm W) <sup>NEW</sup>	\$130
ZPAML	Pivoting Arm Mount <sup>3</sup> – Left (for ZFAPLL / ZFAPLR / ZTRSL / ZTRML / ZTRLL)	\$195
ZPAMR	Pivoting Arm Mount <sup>3</sup> – Right (for ZFAPLR / ZFAPLR / ZTRSR / ZTRMR / ZTRLR)	\$195
ZAPSB_*	Anterior Posture Support Bar <sup>4</sup> * Enter frame width to complete code (22" / 24")	\$260
ZAPSB_*	Anterior Posture Support Bar <sup>4</sup> * Enter frame width to complete code (26" / 28" / 30")	\$305

<sup>1</sup>May require ZPAM for Arm Support to flip up with Flat Arm Pad  
<sup>2</sup> ZPAM is required for Arm Support to flip up with Arm Trough  
<sup>3</sup> Recommended when Lateral Thoracic Support (ZTLS) and Molded Arm Throughs (ZTR) are both ordered  
<sup>4</sup> Not compatible with ZFAP / ZTR

**5. Foot and Leg Supports**

ZFSHDL	Foot Support Left – 15,75" - 20,5"	Standard
ZFSHDR	Foot Support Right – 15,75" - 20,5"	Standard
ZFET2L	Footrest Extension Tube <sup>1</sup> Left – 2" longer	\$30
ZFET2R	Footrest Extension Tube <sup>1</sup> Right – 2" longer	\$30
ZADML	Adj Angle/Depth Molded Footplate - Left	\$110
ZADMR	Adj Angle/Depth Molded Footplate - Right	\$110
ZFET2AL	Footrest Extension Tube Left for ZADML – 2" longer	\$30
ZFET2AR	Footrest Extension Tube Right for ZADMR – 2" longer	\$30
ZFBT	Flip-Back Foot Support (1-piece) 16,5" - 20" with calf pads	\$525
ZELRL	Elevating Leg Support – Left (includes ZADML)	\$570
ZELRR	Elevating Leg Support – Right (includes ZADMR)	\$570
ZRLSL	Universal Residual Limb Support – Left	\$225
ZRLSR	Universal Residual Limb Support – Right	\$225

<sup>1</sup> Not compatible with ZADM / ZVS / ZFBT

## 6. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$110
ZCS_ _*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (22"/ 24")	\$155
ZCS_ _*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (26"/ 28" / 30")	\$310
ZCSIC_ _*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (22" / 24" / 26" / 28" / 30")	\$165
ZHS2224	H-Strap – for 22"/ 24"W frames – hook and loop attachment	\$130
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBT)	\$70
ZNC	Neoprene Cover for Standard Footplate (Each per Side)	\$55
ZNCA	Neoprene Cover for Angle/Depth Adjustable (Each Per Side)	\$55
ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$115
ZLOR	Lateral Offset Foot Support Receiver (ea) – moves foot support 5cm laterally or medially. Specify 1 or 2 ___	\$135
ZMOR	Medial Offset Foot Support Receiver (ea) – moves foot support 5cm laterally or medially. Specify 1 or 2 ___	\$135
ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 6.5cm to accomodate deeper seats and reduced knee flexion range. Specify 1 or 2 ___	\$140

<sup>1</sup>Not compatible with ZFBT / ZVS

## 7. Positioning Pads – Pad Sizes: Small 5,5" x 3,5" / Large – 17cm x 11cm.

<b>Head Supports</b>		
ZHSAS1	Head Support with small pad	\$570
ZHSAS2	Head Support with large pad	\$630
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	\$90
<b>Lateral Thoracic Supports</b> – Pads on Lateral Thoracic Supports can be adjusted medially up to 2.75" per side. For up to 4.3" medial offset, specify ZLAEP extension plate.		
ZLTS2LHD	Lateral Thoracic Support w/ Large Pad for HD chairs – Left	\$380
ZLTS2RHD	Lateral Thoracic Support w/ Large Pad for HD chairs – Right	\$380
<b>Leg Adductors</b>		
ZLGAS2L	Leg Adductor w/ Large Pad – Left	\$380
ZLGAS2R	Leg Adductor w/ Large Pad – Right	\$380
<b>Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors</b>		
These plates allow for an additional 4cm medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors		
ZLAEPL	Lateral/Adductor Extension Plate - Left (ea)	\$30
ZLAEPR	Lateral/Adductor Extension Plate - Right (ea)	\$30
Lateral Thoracic Support	Left      Right	
Leg Adductor	Left      Right	
Hip Adductor	Left      Right	
Total Quantity ___		
<b>Medial Thigh Supports (Pommel) Large Pad Size – 17cm x 11cm</b>		
ZMTSAS2	Medial Thigh Support w/ Large Pads	\$690

## 8. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
ZPBXLHDAS	Pelvic Belt – XL (seat frame widths 22"-30", 2-piece w/ side release buckle)	\$90	ZCBXL	Chest Belt – XL <sup>2</sup> (back frame widths 22"-30", 1-piece w/ hook and loop closure)	\$70
ZPBXXLHDAS	Pelvic Belt – XXL (seat frame widths 28" -2.5", 2-piece w/ side-release buckle)	\$110	ZCBXXL	Chest Belt – XXL <sup>2</sup> (back frame widths 26" - 2.5", 1-piece w/ hook and loop closure)	\$75
ZPBICXLHDAS	Infection Control Pelvic Belt – XL <sup>1</sup> (seat frame widths 24" -26" , 2-piece w/ side-release buckle)	\$145	ZCBICXL	Infection Control Chest Belt – XL <sup>1,2</sup> (back frame widths 22"-30", 1-piece w/ side-release buckle)	\$125
ZPBICXXLHDAS	Infection Control Pelvic Belt – XXL <sup>1</sup> (seat frame widths 26" - 30", 2-piece w/ side-release buckle)	\$145	ZCBICXXL	Infection Control Chest Belt – XXL <sup>1,2</sup> (back frame widths 26" - 30", 1-piece w/ side-release buckle)	\$120
ZPBBPXLHDAS	Pelvic Belt – Bodypoint Aeromesh – XL	\$275	ZCBBPXL	Chest Belt – Bodypoint Aeromesh – XLarge	\$310
<b>Recommended sizes for Bodypoint belts</b>	Chest belt	ZCBBPXL	22"- 30" W back frame		
	Pelvic belt	ZPBBPXLHDAS	22"- 30" W base frame		
<sup>1</sup> Non-absorbent belts – recommended for infection-control environments. <sup>2</sup> Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCB). When the chair is supplied with a Symphony Back, Chest Belt loops are provided.					

## 9. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$75
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$170
ZSC23HD	Splash Curtain <sup>1</sup> for 23" W HD Molded and VFI Seats (surrounds entire aperture, raises seat 0.4").	\$410
ZSC26HD	Splash Curtain <sup>1</sup> for 26" W HD Molded and VFI Seats (surrounds entire aperture, raises seat 0.4").	\$440
ZSCCHD	Custom Splash Curtain for Custom Seats – Raises seat 0.4"	\$420
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$390
ZBUDBH	Bard urine drainage bag holder for 22" / 24" / 26" / 28" / 30"W frames	\$185
<sup>1</sup> Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.		

**10. Other**

ZDCLA	Directional Caster Lock – fits on left rear caster (facilitates steering)	\$165	ZPSK1	Commode Pan Spacer Kit <sup>1</sup> – 1" lower	\$225
ZATAS5	Anti-Tippers <sup>1</sup> – Rear (pr)	\$395	ZPSK2	Commode Pan Spacer Kit <sup>1</sup> – 2" lower	\$245
ZCUAP6	6" Caster Upgrade (4x)	\$395	ZPSK3	Commode Pan Spacer Kit <sup>1</sup> – 3" lower	\$285
			ZPSK4	Commode Pan Spacer Kit <sup>1</sup> – 4" lower	\$285

<sup>1</sup> Commode Pan Spacer Kits provide extra clearance for perianal care with pan in place. Increases height required to clear toilet.

**11. Credits and Omissions**

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z907HD	Omit HD Foot Supports – if no other Foot / Leg Support Option is selected	-\$45
------	--	-------	--------	---	-------

Some assembly of Raz products is required.

**12. Notes:**

--

Specifications	Raz-AP600 (Attendant Propel)	
Weight Capacity	600 lbs	
Overall Width	Seat frame width plus 4.5" (with standard arm pad)	
Overall Length (with foot supports)	40.5" (with rear casters trailing), 37.25" (to rear of back canes)	
Seat Depth	18"-20"	
Toilet Clearance Width	seat frame width minus 6"	
Toilet Clearance Height (under seat)	with 5" casters	with 6" casters
- with pan brackets	14.3"-18.3"	15.3"-19.3"
- without pan brackets	15.1"-19.1"	16.1"-20.1"
Arm Support Height - above molded seat	10"	
Seat-to-Floor Height (with molded seat)	19"-23" (with 5" casters)	
Back Height (Molded Seat to top of upholstery)	20"	
Product Weight (with 22"W seat frame)	54.5 lbs	
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)	