



RAZ-AT600 ORDER FORM

EFFECTIVE MAY 2026 (CDN MSRP)

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-AT600 Heavy Duty Mobile Shower Commode Chair with stainless steel frame and components, 25° tilt range, Fixed Back, four dual-locking 6" (150mm) casters, flip-up arm supports, heavy duty foot supports, commode pan, head support with small pad and the proprietary Weight Assist Spring System (WASS). Also included is a 20"H Multi-Position Back with 2" fore/aft adjustment range, two back angles (90°/97°) and tension-adjustable vinyl upholstery. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2" of fore/aft seat adjustment. **User weight capacity 275 kgs.** For help with justification, please see our [Evaluation and Medical Justification Form](#).

Need help with this order form? Call customer service at (866) 875-2877 or (519) 942 8441

1. Raz-AT600 codes by size – Must Choose One

Z36022	22"W Seat and Back Frames	\$7,275
Z36024	24"W Seat and Back Frames	\$7,590
Z36026	26"W Seat and Back Frames	\$8,165
Z36028	28"W Seat and Back Frames	\$8,730
Z36030	30"W Seat and Back Frames – Requires custom seat	\$9,310

2. User Weight Range

Z36WR450	WASS User Weight Range 159 - 204 kgs.	Std	Z36WR600	WASS User Weight Range 181 - 275 kgs.	\$290
----------	---------------------------------------	-----	----------	---------------------------------------	-------

3. Seats – Must Choose One

Molded and Visco Foam Interface Seats – Seat Support Kit (ZSSKT) required if seat is narrower than seat frame.			
ZMS23FRO	Molded 58.5cm W x 18" D – Front Access	Fits 22" and 24"W ¹ frames only	N/C
ZMS26FRO	Molded 26" W x 18" D – Front Access	Fits 26" and 28"W ¹ frames only	N/C
ZVFI23FRO	Visco Foam Interface 58.5cm W x 18" D – Front Access	Fits 22" and 24"W ¹ frames only	\$835
ZVFI26FRO	Visco Foam Interface 26" W x 18" D – Front Access	Fits 26" and 28"W ¹ frames only	\$860
ZSSKTHD	Seat Support Kit for Raz-AT600 – for seats that are narrower than seat frames		\$165

¹ Requires a Seat Support Kit (ZSSKTHD)

Custom Seats

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. See [Custom Seat Assessment Form – Resources tab at razdesigninc.com](#)

ZCSEATHD	Custom Dual-Visco Foam Seat – HD (enter specs below) Select if both width and depth are ≤ 28"	\$1,545
ZCSEATXHD	Custom Dual-Visco Foam Seat – XHD" (enter specs below) Select if width and/or depth is ≥ 28"	\$1,745
Base Frame Width:	__ 22" / 24" / 26" / 28" / 30"	
Custom Seat Width:	__ 22" / 24" / 26" / 28" / 30" / 32" (cannot be more than 2" wider than seat frame width)	
Seat Depth:	__ 16" / 17" / 18" / 19" / 20" / 21" / 22" / 23" / 24" (≥ 21" requires custom longer frame)	
Access:	__ FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged	
Aperture Width:	__ 6.5" recommended	

- Unless otherwise requested, the width of the aperture (inside width dimension of the commode opening) will be 16.5cm for custom seats.
- Custom seat widths can be up to 5cm wider than seat the frame width to provide additional thigh support without increasing overall width of the chair.
- The Ischial Pelvic Alignment System (IPAS) will allow 5cm fore/aft adjustment of the seat position (and therefore seat depth and aperture position) on all seats.
- Access for Rear-Open seats is limited due to rear seat support tube.



Motion Composites is the exclusive Canadian Distributor of RAZ DESIGN INC.

motioncomposites.com

4. Symphony Back (Recessed solid back that replaces the standard Tension-Adjustable Vinyl Upholstery [recommended for infection-control environments]). All Symphony Back panels are 13" H.

ZSB22	Symphony Back – 22"W	\$525
ZSB24	Symphony Back – 24"W	\$555
ZSB26	Symphony Back – 26"W	\$580
ZSB28	Symphony Back – 28"W	\$745
ZSB30	Symphony Back – 30"W	\$770

5. Arm Supports Accessories

ZALKHD	Arm Support Locks – HD (pr)	\$110
ZFAPLL	Long Molded Flat Armpad ¹ – Left (14"L x 3.25"W))	\$95
ZFAPLR	Long Molded Flat Armpad ¹ – Right (14"L x 3.25"W)	\$95
ZTRML	Arm Trough ² - Medium Left (18.5" L x 5" W)	\$130
ZTRMR	Arm Trough ² - Medium Left (18.5" L x 5" W)	\$130
ZTRSL	Arm Trough ² - Large Left (21" L x 5" W) NEW	\$130
ZTRSR	Arm Trough ² - Large Right (21" L x 5" W) NEW	\$130
ZPAML	Pivoting Arm Mount ³ – Left (for ZFAPLL / ZFAPLR / ZTRSL / ZTRML / ZTRLL)	\$195
ZPAMR	Pivoting Arm Mount ³ – Right (for ZFAPLR / ZFAPLR / ZTRSR / ZTRMR / ZTRLR)	\$195
ZAPSB_*	Anterior Posture Support Bar ⁴ * Enter frame width to complete code (22" / 24")	\$260
ZAPSB_*	Anterior Posture Support Bar ⁴ * Enter frame width to complete code (26" / 28" / 30")	\$305

¹May require ZPAM for Arm Support to flip up with Flat Arm Pad
² ZPAM is required for Arm Support to flip up with Arm Trough
³ Recommended when Lateral Thoracic Support (ZTLS) and Molded Arm Throughs (ZTR) are both ordered
⁴ Not compatible with ZFAP / ZTR

6. Leg Supports

ZFSHDL	Foot Support Left – 15.75" - 20.5"	Standard
ZFSHDR	Foot Support Right – 15.75" - 20.5"	Standard
ZFET2L	Footrest Extension Tube ¹ Left – 5cm longer	\$30
ZFET2R	Footrest Extension Tube ¹ Right – 5cm longer	\$30
ZADML	Adj Angle/Depth Molded Footplate - Left	\$110
ZADMR	Adj Angle/Depth Molded Footplate - Right	\$110
ZFET2AL	Footrest Extension Tube Left for ZADML – 5cm longer	\$30
ZFET2AR	Footrest Extension Tube Right for ZADMR – 5cm longer	\$30
ZFBT	Flip-Back Foot Support (1-piece) 16.5" - 20" with calf pads	\$525
ZELRL	Elevating Leg Support – Left (includes ZADML)	\$570
ZELRR	Elevating Leg Support – Right (includes ZADMR)	\$570
ZRLSL	Universal Residual Limb Support – Left	\$225
ZRLSR	Universal Residual Limb Support – Right	\$225

¹ Not compatible with ZADM / ZVS / ZFBT

7. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$110
ZCS_*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (22"/ 24")	\$155
ZCS_*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (26"/ 28" / 30")	\$310
ZCSIC_*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (22" / 24" / 26" / 28" / 30")	\$165
ZHS2224	H-Strap – for 22" / 24" W frames – hook and loop attachment	\$130
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBT)	\$70
ZNC	Neoprene Cover for Standard Footplate (Each)	\$55
ZNCA	Neoprene Cover for Angle/Depth Adjustable (Each)	\$55
ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$115
ZLOR	Lateral Offset Foot Support Receiver (ea) – moves foot support 5cm laterally or medially. Specify 1 or 2 ___	\$135
ZMOR	Medial Offset Foot Support Receiver (ea) – moves foot support 5cm laterally or medially. Specify 1 or 2 ___	\$135
ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 6.5cm to accomodate deeper seats and reduced knee flexion range. Specify 1 or 2 ___	\$140

¹ Not compatible with ZFBT / ZVS

8. Positioning Pads – Pad Sizes: Small 5.5" x 3.5" / Large – 6.75" x 4.3".

Head Supports		
ZHSP2	Headrest Pad – Large (12" x 6")	\$60
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	\$90
Lateral Thoracic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 2.75" per side. For up to 4.3" medial offset, specify ZLAEPL extension plate.		
ZLTS2LHD	Lateral Thoracic Support w/ Large Pad for HD chairs – Left	\$380
ZLTS2RHD	Lateral Thoracic Support w/ Large Pad for HD chairs – Right	\$380
Leg Adductors		
ZLGT2L	Leg Adductor w/ Large Pad – Left	\$380
ZLGT2R	Leg Adductor w/ Large Pad – Right	\$380
Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors		
These plates allow for an additional 4cm medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors		
ZLAEPL	Lateral/Adductor Extension Plate - Left (ea)	\$30
ZLAEPR	Lateral/Adductor Extension Plate - Right (ea)	\$30
Lateral Thoracic Support	Left Right	
Leg Adductor	Left Right	
Hip Adductor	Left Right	
Total Quantity ___		
Medial Thigh Supports (Pommel) Large Pad Size – 6.7" x 4.3"		
ZMTST2	Medial Thigh Support w/ Large Pads	\$515

9. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts			Chest Belts		
ZPBXLHDT	Pelvic Belt – XLarge (seat frame widths 22"-30", 2-piece w/ side-release buckle)	\$90	ZCBXL	Chest Belt – XLarge ² (back frame widths 22"-30", 1-piece w/ hook and loop closure)	\$70
ZPBXXLHDT	Pelvic Belt – XXL Large (seat frame widths 28" -30" , 2-piece w/ side-release buckle)	\$110	ZCBXXL	Chest Belt – XXL Large ² (back frame widths 26" -30" , 1-piece w/ hook and loop closure)	\$75
ZPBICXLHDT	Infection Control Pelvic Belt – XLarge ¹ (seat frame widths 22"-30", 2-piece w/ side-release buckle)	\$155	ZCBICXL	Infection Control Chest Belt – XLarge ^{1,2} (back frame widths 24" -26" , 1-piece w/ side-release buckle)	\$125
ZPBICXXLHDT	Infection Control Pelvic Belt – XX-Large ¹ (seat frame widths 26" -30" , 2-piece w/ side-release buckle)	\$170	ZCBICXXL	Infection Control Chest Belt – XLarge ^{1,2} (back frame widths 26" -30" , 1-piece w/ side-release buckle)	\$120
ZPBBPXLHDT	Pelvic Belt – Bodypoint Aeromesh – XLarge	\$235	ZCBBPXL	Chest Belt – Bodypoint Aeromesh – XLarge	\$310
Recommended sizes for Bodypoint belts	Pelvic belt	ZPBBPXLHDT	22"-30"W base frame		
	Chest belt	ZCBBPXL	22"-30"W back frame		
¹ Non-absorbent belts – recommended for infection-control environments. ² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCB). When the chair is supplied with a Symphony Back, Chest Belt loops are provided.					

10. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.	\$80
ZSS	SplaShield – Molded urine deflector (must be held in place during use)	\$170
ZSC23HD	Splash Curtain ¹ for 23" W HD Molded and VFI Seats (surrounds entire aperture, raises seat 1.25cm).	\$410
ZSC26HD	Splash Curtain ¹ for 26" W HD Molded and VFI Seats (surrounds entire aperture, raises seat 0.4").	\$440
ZSCCHD	Custom Splash Curtain for Custom Seats – Raises seat 0.4"	\$420
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets	\$390
ZBUDBH	Bard urine drainage bag holder for 22" / 24" / 26" / 28" / 30"W frames	\$185
¹ Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for pricing.		

11. Other

ZDCLT	Directional Caster Lock – fits on left rear caster (facilitates steering)	\$195	ZPSK1	Commode Pan Spacer Kit ² – 1" lower	\$225
ZATTRR	Anti-Tippers – Rear (pr)	\$395	ZPSK2	Commode Pan Spacer Kit ² – 2" lower	\$245
ZATFFR	Anti-Tippers – Front (pr)	\$395	ZPSK3	Commode Pan Spacer Kit ² – 3" lower	\$285
ZTLS	Tilt Limit Spacers Pair – 1/2" spacers snap onto gas struts to limit tilt (1/2" reduces tilt range approx. 5°. Can be stacked for additional tilt limiting. Quantity ____	N/C	ZPSK4	Commode Pan Spacer Kit ² – 4" lower	\$285
ZSPKTHD24	AT600 Self-Propel Kit – 24" wheels – includes rear anti-tippers	\$1,440	¹ Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.		

12. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z911	Omit Molded Head Support Pad	-\$20
Z907HD	Omit HD Foot Supports – if no other Foot / Leg Support Option is selected	-\$45	Z912	Omit Complete Head Support Assembly	-\$50

Some assembly of Raz products is required.

13. Notes:

--

Specifications	Raz-AT600 (Attendant Tilt)
Weight Capacity	600 lbs
Seat Depth (with Standard Molded Seat)	18"-20"
Overall Width	Seat frame width plus 4.5" (with standard arm pad)
Toilet Clearance Width	17"
Overall Length (with foot supports and headrest adjusted to middle positions)	
- at 0° tilt and standard foot supports	42"
- at 25°	52"
Clearance Height (under seat)	with 6" casters
- with pan brackets	17.75"-18.75"
- without pan brackets	18.5"-19.5"
Arm Support Height - above molded seat	9.25"
Seat-to-Floor Height (with molded seat)	22.25"
Back Height (Molded Seat to top of upholstery)	20"
Product Weight (with 22"W seat frame)	84 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)