



RAZ-CAT ORDER FORM

EFFECTIVE MAY 2026 (CDN MSRP)

Date: _____

Dealer / Location: _____ T: _____ F: _____

PO# / Client: _____ Contact: _____ Email: _____

Raz-CAT –Compact Attendant Tilt Mobile Shower Commode Chair with 40° tilt range, head support fixture with small pad, height adjustable stainless steel frame, Fixed Back, flip-up padded arm supports, foot supports, four dual-locking 5” (125mm) casters and commode pan. Exclusive to Raz, IPAS helps reduce the possibility of pressure injury at key pressure points by providing up to 2” of fore/aft seat adjustment. **User weight capacity 136 kgs.**

For help with justification, please see our [Evaluation and Medical Justification Form](#).

Need help with this order form? Call customer service at (866) 875-2877 or (519) 942 8441

1. Standard Raz-CAT Frame ((18”W) – If you require a special width, depth and / or height see [Sec. 2]


Z30318	Raz-CAT w/ 18”W base frame – Standard base frame width, will fit over most toilets	\$3,925
--------	--	---------

2. Custom / Special Frame Widths, Depths and Heights – Adjusta-Back must be selected for all custom widths see [Sec.7]

Z30316	Raz-CAT w/ 16”W frame width ¹ (base chair price – custom upcharge below will be added)	\$3,925
Z3CSW16	Upcharge for custom 16”W seat and base frame – will not fit over most toilets	\$690
Z30320	Raz-CAT w/ 20”W frame width ² (base chair price – custom upcharge below will be added)	\$3,925
Z3CSW20	Upcharge for custom 20”W seat and base frame	\$690
Custom Frame Configuration for 18” W frame. (Call Raz for codes and pricing for 16” / 20” widths)		
Z3CLOWF5	Upcharge for custom Low Frame (16”-20” STF)	\$690

¹ 18”W Visco-Foam Interface seats can be modified to fit.
² 18”W seats can be selected but must also select Seat Support Kit (Z4SSKCT) [Sec.5]

3. Seats – Standard – Please note that the Molded and VFI seats, formerly designated 15.75” D, are now designated 17” D, which is the center of the IPAS range and the position in which they are shipped from the factory.

Standard CAT Seat – Visco-Foam w/ Tapered Oblong Aperture (thickness ~2.5”)			
ZCSC1816	18”W x 11”-16”D (adjustable)	Standard w/ 18”W frame	
ZCSC1616	16”W x 11”-16”D (adjustable)	Standard w/ 16”W frame	
ZCSC2016	20”W x 11”-16”D (adjustable)	Standard w/ 20”W frame	
Set seat depth at: 11” 12” 13” 14” 15” 16”			
Standard aperture width: 3.5” at front, 4.5” at rear			
Other Seats Available – for seat depth adjustment up to 17” – Modified to fit CAT seat frame			
Molded Seats (thickness ~6.2”-2.75”)		Visco Foam Interface Seats (thickness ~9cm)	
ZMSC17FRO	18”W x 16”D – Front Access N/C	ZVFIC17FRO	18”W x 17”D – Front Access \$495
ZMSC17BRI	18”W x 16”D – Front Bridged N/C	ZVFIC17BRI	18”W x 17”D – Front Bridged \$495
		ZVFIC17LTO	18”W x 17”D – Left Access \$495
		ZVFIC17RTO	18”W x 17”D – Right Access \$495
		ZVFIC17RRO	18”W x 17”D – Rear Access \$495

4. Seats – Custom

For custom seats, please specify seat width and depth. Indicate if the access opening is to be in the Front, Left, Right or Rear. Select Bridged if an access opening is not required. 18" is the standard width for CAT seat frames.

See [Custom Seat Assessment Form](#) – Resources tab at [razdesigninc.com](#)

ZCSEAT	Custom Dual-Visco Foam Seat (enter specs below)	\$1,190
Seat Frame Width:	-- 16" / 18" / 20"	
Custom Seat Width:	-- 16" / 17" / 18" / 19" / 20"	
Custom Seat Depth:	-- 11" / 12" / 13" / 14" / 15" / 16" / 17"	
Access Opening:	-- FR – Front / LT – Left / RT – Right / RR – Rear / BR – Bridged	
Aperture:	Tapered Oblong: __"W (at front) __"W (at rear) [standard is 3.5"W at front and 4.5"W at rear] Keyhole: __"W (3.0" - 8.0"W – enter "NA" if aperture is not required)	
<ul style="list-style-type: none"> ● Oblong-shaped apertures can be specified in conjunction with any of the access openings listed. ● For Custom Seats, the Ischial Pelvic Alignment System (IPAS) allows +/- 2.2" fore / aft adjustment from specified custom seat depth. ● Special modifications on Custom Seats include: leg-length discrepancy, soft spots, pressure-relief areas, contouring, pelvic obliquity accommodation, etc. 		

5. Back Frames – for standard 18"W seat frames (back heights measured from top of seat to top of upholstery)

Multi-Position Adjusta-Back: The "Multi-Position" feature allows for fore/aft adjustment of the backrest to provide seat depth adjustability from 23.15" - 15.75". Includes tension-adjustable vinyl back support, flip-up padded arm supports, back angle adjustment (90° / 97°) with cane-style back posts. Back height is measured from top of molded seat to top of back support. Hip Adductors [ZHPG] are not compatible with Fixed Back.

Back frames for standard 18"W seat frame (Z30318)			Reclining Backs for 18"W seat frame		
ZCAB845	14" Back Width / 15" Back Height	N/C	ZRBC848	Recliner Back for CAT – 14"W x 18" H ^{2,3}	\$1,145
ZCAB848	14" Back Width / 18" Back Height	\$110	ZRBC868	Recliner Back for CAT – 16"W x 18" H ^{2,3}	\$1,145
ZCAB865	16" Back Width / 15" Back Height	N/C	ZRBCR888	Recliner Back for CAT – 18"W x 18" H ^{2,3}	\$1,145
ZCAB868	16" Back Width / 18" Back Height	\$110			
ZCAB885	18" Back Width / 15" Back Height	N/C			
ZCAB888	18" Back Width / 18" Back Height	\$110			
ZDTL	Dual Triggers – left side ¹	N/C			
ZDTR	Dual Triggers – right side ¹	N/C			

¹ Standard trigger configuration for tilt-only is one on each side. Dual Triggers option positions both triggers on one side, allowing for single-handed operation (tilt only).

² Standard trigger locations for Raz-CAT with recline are: tilt triggers on the right and recline triggers on the left.

³ Tilt is limited to 25°

6. Back Frames for 16"W and 20"W Seat Frames (back heights measured from top of seat to top of upholstery)

Multi-Position Adjusta-Back: The "Multi-Position" feature allows for fore/aft adjustment of the backrest to provide seat depth adjustability from 11"-16". Includes tension-adjustable vinyl back support, back angle adjustment (90° / 97°) and flip-up padded arm supports. Back upholstery can be positioned 2.5" higher than factory settings stated below.

16"W Seat Frames			Reclining Backs for 16"W seat frame		
ZCAB645	14" Back Width / 15" Back Height	N/C	2-Point Arms (Z2PLT / Z2PST) are included with Reclining Backs.		
ZCAB648	14" Back Width / 18" Back Height	\$110	ZRBC648	Recliner Back for CAT –14"W x 18"H ^{2,3}	\$1,145
ZCAB665	16" Back Width / 15" Back Height	N/C	ZRBC668	Recliner Back for CAT – 16"W x 18"H ^{2,3}	\$1,145
ZCAB668	16" Back Width / 18" Back Height	\$110			
20"W Seat Frame			Reclining Backs for 20"W seat frame		
ZCAB065	16" Back Width / 15" Back Height	N/C	2-Point Arms (Z2PLT / Z2PST) are included with Reclining Backs.		
ZCAB068	16" Back Width / 18" Back Height	\$110	ZRBC088	Recliner Back for CAT – 18"W x 18"H ^{2,3}	\$1,145
ZCAB085	18" Back Width / 15" Back Height	N/C			
ZCAB088	18" Back Width / 18" Back Height	\$110			
ZCAB005	20" Back Width / 15" Back Height	N/C			
ZCAB008	20" Back Width / 18" Back Height	\$110			
ZDTL	Dual Triggers – left side ¹	N/C			
ZDTR	Dual Triggers – right side ¹	N/C			
¹ Standard trigger configuration for tilt-only is one on each side. Dual Triggers option positions both triggers on one side, allowing for single-handed operation (tilt only).			² Standard trigger locations for Raz-CAT with recline are tilt triggers on the right and recline triggers on left. ³ Tilt is limited to 25°		

7. Back Support Options – These back supports replace tension-adjustable upholstery on an Adjusta-Back. Symphony Backs are recommended for infection-control environments (for 18"H backs only).

ZSB14	Symphony Back – recessed solid back support – 14" back frame width	\$340
ZSB16	Symphony Back – recessed solid back support – 16" back frame width	\$340
ZSB18	Symphony Back – recessed solid back support – 18" back frame width	\$340
ZSB20	Symphony Back – recessed solid back support – 20" back frame width	\$460

8. Arm Supports and Accessories

ZARMTS	Short Arm Support Tubes (pr) ¹	Standard
ZARMTL	Long Arm Support Tubes (pr)	N/C
ZARMSL	Short Arm Support Pad 10" L – Left	Standard
ZARMSR	Short Arm Support Pad 10" L – Right	Standard
ZARMLL	Long Arm Support Pad 13" L – Left	N/C
ZARMLR	Long Arm Support Pad 13" L – Right	N/C
ZALK	Arm Support Locks (pr) for Fixed Back (ZFB1818)	\$90
Z2PLT	2-Point Arm Supports ² (pr) with 13" L pads (std with reclining backs)	\$525
Z2PST	2-Point Arm Supports ² (pr) with 10" L pads (std with reclining backs)	\$525
ZFAPLL	Long Molded Flat Armpad ³ – Left (13.75" L x 3.15" W – fits ZARM13)	\$95
ZFAPLR	Long Molded Flat Armpad ³ – Right (13.75" L x 3.15" W – fits ZARM13)	\$95
ZFAPSL	Short Molded Flat Armpad ³ – Left (10" L x 2.75" W – fits ZARM10)	\$95
ZFAPSR	Short Molded Flat Armpad ³ – Right (10" L x 2.75" W – fits ZARM10)	\$95
ZTRSL	Arm Trough ⁴ - Small Left (15.75" L x 4.3" W) NEW	\$130
ZTRSR	Arm Trough ⁴ - Small Right (15.75" L x 4.3" W) NEW	\$130
ZTRML	Arm Trough ⁴ - Medium Left (18.5" L x 4.75" W)	\$130
ZTRMR	Arm Trough ⁴ - Medium Right (18.5" L x 4.75" W)	\$130
ZTRLL	Arm Trough ⁴ - Large Left (21" L x 5.1" W) NEW	\$130
ZTRLR	Arm Trough ⁴ - Large Right (21" L x 5.1" W) NEW	\$130
ZPAML	Pivoting Arm Mount ⁵ – Left (for ZFAPLL / ZFAPLR / ZTRSL / ZTRML / ZTRLL)	\$190
ZPAMR	Pivoting Arm Mount ⁵ – Right (for ZFAPLR / ZFAPLL / ZTRSR / ZTRMR / ZTRLR)	\$190
ZASK05	Arm Support Spacer Kit ⁶ – 1.22" per side for all Adjusta-Backs (pr)	\$70
ZASK10	Arm Support Spacer Kit ⁶ – 2.2" per side for all Adjusta-Backs (pr)	\$135
ZPHGL	Pivoting Hand Grips (pr) for 13" L arms ⁷	\$460
ZPHGS	Pivoting Hand Grips (pr) for 10" L arms ⁷	\$460
ZAPSB_*	Anterior Posture Support Bar ^{8,9} * Enter frame width to complete code (14" / 16" / 18" / 20")	\$225
ZAPSB_*	Anterior Posture Support Bar ^{8,9} * Enter frame width to complete code (22" / 24")	\$260
Reclining Back Arm Support Accessories		
Z2PLTR	2-Point Flip-up Arm Supports (Long) (pr)	N/C
Z2PSTR	2-Point Flip-up Arm Supports (Short) (pr)	N/C
ZASKR05	Arm Support Spacers 1.22" per side (pr) for Reclining Back Support	\$70
ZASKR10	Arm Support Spacers 2.2" per side (pr) for Reclining Back Support	\$70
<p>¹ Not compatible with ZPAM / ZFAPL / ZFAPS / ZTRM / ZTRL</p> <p>² For chairs with same back width and frame width only</p> <p>³ May require ZPAM for Arm Support to flip up with Flat Arm Pad</p> <p>⁴ ZPAM is required for Arm Support to flip up with Arm Trough</p> <p>⁵ Recommended when Lateral Thoracic Support (ZTLS) and Molded Arm Throughs (ZTR) are both ordered.</p> <p>⁶ ZASK05 or ZASK10 is required for Standard Arm Support to flip up with lateral pad in forward position</p> <p>⁷ ZPHGL ZPHGS combined with 2 Point Arms requires a custom up charge</p> <p>⁸ Not compatible with ZFAP / ZTR</p> <p>⁹ Special length required with ZASK05 and ZASK10 (contact customer service)</p>		

9. Foot and Leg Supports

ZFSL	Foot Support Left – 15.75” - 20.5”	N/C
ZFSR	Foot Support Right – 15.75” - 20.5”	N/C
ZFET2L	Footrest Extension Tube ¹ Left – 2” longer	\$30
ZFET2R	Footrest Extension Tube ¹ Right – 2” longer	\$30
ZMFX8CL	MFX8 Foot Support Left – 8.25” - 12”	N/C
ZMFX8CR	MFX8 Foot Support Right – 8.25” - 12”	N/C
ZMFX12CL	MFX12 Foot Support Left – 12” - 16.5”	Standard
ZMFX12CR	MFX12 Foot Support Right – 12” - 16.5”	Standard
ZADABS18L	Angle/Depth Adj ABS Footplate Left For 18” W frame	\$110
ZADABS18R	Angle/Depth Adj ABS Footplate Right For 18” W frame	\$110
ZADML	Adj Angle/Depth Molded Footplate ² - Left	\$110
ZADMR	Adj Angle/Depth Molded Footplate ² - Right	\$110
ZFET2AL	Footrest Extension Tube Left for ZADML – 2” longer	\$30
ZFET2AR	Footrest Extension Tube Right for ZADMR – 2” longer	\$30
ZVS__*	V-Style Foot Support (1-piece) 15.75” - 20.5” * Enter frame width to complete code (16” / 18”)	\$340
ZVS20	V-Style Foot Support (1-piece) 15.75” - 20.5” for 20”W frame	\$375
ZFBT12	MFX12 Flip-Back Foot Support (1-piece) 12” - 16.5” with calf pads	\$525
ZFBT	Flip-Back Foot Support (1-piece) 16.5” - 51cm with calf pads	\$525
ZELRL	Elevating Leg Support ³ – Left (includes ZADML)	\$570
ZELRR	Elevating Leg Support ³ – Right (includes ZADMR)	\$570
ZRLSL	Universal Residual Limb Support – Left	\$225
ZRLSR	Universal Residual Limb Support – Right	\$225
Foot Support Options for 16”W Frames		
ZFS16L	Foot Support for 16”W frame Left – (15.75” - 20.5”)	N/C
ZFS16R	Foot Support for 16”W frame Right – (15.75” - 20.5”)	N/C
ZADABS16L	Angle/Depth Adj ABS Footplate Left for 16”W frame	\$110
ZADABS16R	Angle/Depth Adj ABS Footplate Right For 16”W frame	\$110
ZMFX816L	MFX8 Foot Support for 15.75” W Frame Left – 8.25” - 12”	\$45
ZMFX816R	MFX8 Foot Supports for 15.75” W Frame Right – 8.25” - 12”	\$45
ZMFX1216L	MFX12 Foot Support for 15.75” W Frame Left – 12” - 16.5”	\$55
ZMFX1216R	MFX12 Foot Support for 15.75” W Frame Right – 12” - 16.5”	\$55
¹ Not compatible with ZMFX12/ZADM/ZVS/ZFBT		
² Not compatible with Z30016		
³ ZELRL/R compatible with frame widths ≥ 18” and adjusted seat depths ≤ 13” or ≤ 17” w/ forward offset footrest receivers (ZFOR)		

10. Foot and Leg Support Accessories

ZFSLK	Foot Support Lock (pr)	\$110
ZCS_*	Calf Strap, Snap-On – Bodypoint Aeromesh * Enter frame width to complete code (16" / 18" / 20")	\$155
ZCSIC_*	Calf Strap, Snap-On – Infection Control * Enter frame width to complete code (16" / 18" / 20")	\$145
ZCPL	Calf Panel – Bodypoint Aeromesh ¹ Fits frame widths (16" / 18" / 20")	\$210
ZHS1820	H-Strap – for 18" / 20"W frames – hook and loop attachment	\$105
ZAHL	Adjustable Heel Loops (pr) (not available for ZVS / ZFBT)	\$70
ZNC	Neoprene Cover for Standard Footplate (Each)	\$55
ZNCA	Neoprene Cover for Angle/Depth Adjustable (Each)	\$55
ZNC16	Neoprene Covers for Narrow Molded Footplates –16"W frame (Each)	\$55
ZNCABS16	Neoprene Cover for Narrow Molded ABS Angle/Depth Adjustable (Each)	\$55
ZNCVS	Neoprene Cover for V-Style Foot Support	\$55
ZNCFB	Neoprene Cover for Flip-Back Foot Support	\$55
ZLOR	Lateral Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ___	\$135
ZMOR	Medial Offset Foot Support Receiver (ea) – moves foot support 2" laterally or medially. Specify 1 or 2 ___	\$135
ZFOR	Forward Offset Foot Support Receiver – (ea) Moves foot support forward 6.2" to accommodate deeper seats and reduced knee flexion range. Specify 1 or 2 ___	\$140

¹ Not compatible with ZMFX12 / ZFBT / ZVS

11. Positioning Pads – Hip Adductors are not compatible with Reclining Backrest unless ordered on-chair.

Head Supports		
ZHSP1	Head Support Pad – Small (10" x 10")	Standard
ZHSP2	Head Support Pad – Large (12" x 12")	\$60
ZOHIP	Lateral Offset Head Support Interface Plate – up to 22° of adjustment	\$90
Lateral Thoracic Supports – Pads on Lateral Thoracic Supports can be adjusted medially up to 2.75" per side. For up to 4.3" medial offset, specify ZLAEP extension plate.		
ZLTS1L	Lateral Thoracic Support w/ Small Pad – Left	\$380
ZLTS1R	Lateral Thoracic Support w/ Small Pad – Right	\$380
ZLTS2L	Lateral Thoracic Support w/ Large Pad – Left	\$380
ZLTS2R	Lateral Thoracic Support w/ Large Pad – Right	\$380
Leg Adductors		
ZLGT1L	Leg Adductor w/ Small Pad – Left	\$380
ZLGT1R	Leg Adductor w/ Small Pad – Right	\$380
ZLGT2L	Leg Adductor w/ Large Pad – Left	\$380
ZLGT2R	Leg Adductor w/ Large Pad – Right	\$380
Hip Adductors for 4" narrower backs		
ZHP4C1L	Hip Adductor for 4" narrower back w/ Small Pad– Left	\$380
ZHP4C1R	Hip Adductor for 4" narrower back w/ Small Pad– Right	\$380
ZHP4C2L	Hip Adductor for 4" narrower back Large Pad– Left	\$380
ZHP4C2R	Hip Adductor for 4" narrower back w/ Large Pad– Right	\$380
Hip Adductors for backs that are same width as seat frame or 2" narrower		
ZHP02C1L	Hip Adductor for 0 - 2" narrower back w/ Small Pad – Left	\$380
ZHP02C1R	Hip Adductor for 0 - 2" narrower back w/ Small Pad – Right	\$380
ZHP02C2L	Hip Adductor for 0 - 2" narrower back w/ Large Pad – Left	\$380
ZHP02C2R	Hip Adductor for 0 - 2" narrower back w/ Large Pad – Right	\$380
Hip Adductors For Reclining Backs		
ZHPRC1L	Hip Adductor for Reclining Back w/ Small Pad – Left	\$380
ZHPRC1R	Hip Adductor for Reclining Back w/ Small Pad – Right	\$380
ZHPRC2L	Hip Adductor for Reclining Back w/ Large Pad – Left	\$380
ZHPRC2R	Hip Adductor for Reclining Back w/ Large Pad – Right	\$380
Extension Plates for Lateral Thoracic Supports and Hip / Leg Adductors		
These plates allow for an additional 2.2" medially and can be used on Lateral Thoracic Support, Hip and Leg Adductors		
ZLAEPL	Lateral/Adductor Extension Plate – Left	\$30
ZLAEPR	Lateral/Adductor Extension Plate – Right	\$30
Lateral Thoracic Support	Left	Right
Leg Adductor	Left	Right
Hip Adductor	Left	Right
Total Quantity ____		
Medial Thigh Supports (Pommel)		
ZMTST1	Medial Thigh Support w/ Small Pads	\$515
ZMTST2	Medial Thigh Support w/ Large Pads	\$515

12. Pelvic and Chest Belts – A CHEST BELT MUST NEVER BE USED WITHOUT A PELVIC BELT!

Pelvic Belts				Chest Belts			
ZPBS	Pelvic Belt – Small (frame widths 13.75" - 15.75", 2-piece w/ side-release buckle)		\$50	ZCBS	Chest Belt ² – Small (back widths 13.75" - 15.75", 1-piece w/ hook and loop closure)		\$40
ZPBM	Pelvic Belt – Medium (frame widths 15.75" - 18", 2-piece w/ side-release buckle)		\$60	ZCBM	Chest Belt ² – Medium (back widths 15.75" - 18", 1-piece w/ hook and loop closure)		\$45
ZPBL	Pelvic Belt – Large (frame widths 18" - 51cm, 2-piece w/ side-release buckle)		\$55	ZCBL	Chest Belt ² – Large (back widths 18" - 51cm, 1-piece w/ hook and loop closure)		\$55
ZPBICS	Pelvic Belt ¹ – Infection Control, Small 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 13.75" - 15.75")		\$110	ZCBICS	Chest Belt ^{1,2} – Infection Control, Small 1-piece w/ side-release buckle (back widths 13.75" - 15.75")		\$90
ZPBICM	Pelvic Belt ¹ – Infection Control, Medium 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 15.75" - 18")		\$110	ZCBICM	Chest Belt ^{1,2} – Infection Control, Medium 1-piece w/ side-release buckle (back widths 15.75" - 18")		\$90
ZPBICL	Pelvic Belt ¹ – Infection Control, Large 2-piece w/ frame clamps, mounting plates, side-release buckle (frame widths 18" - 51cm)		\$110	ZCBICL	Chest Belt ^{1,2} – Infection Control, Large 1-piece w/ side-release buckle (back widths 18" - 51cm)		\$100
ZPBBPM	Pelvic Belt – Bodypoint Aeromesh, Med		\$295	ZCBBPM	Chest Belt – Bodypoint Aeromesh, Med		\$310
ZPBBPL	Pelvic Belt – Bodypoint Aeromesh, Large		\$295	ZCBBPL	Chest Belt – Bodypoint Aeromesh, Large		\$310
Recommended sizes for Bodypoint belts	Pelvic belt	ZPBBPS	N/A	ZPBBPM	15.75" W base frame	ZPBBPL	18"-22" W base frame
	Chest belt	ZCBBPS	13.75" W back frame	ZCBBPM	15.75" W back frame	ZCBBPL	18"-22" W back frame
¹ Non-absorbent belts – recommended for infection-control environments ² Tension-Adjustable Upholstery includes loops to retain a chest belt (ZCBM/ZCBL/ZCBLXL). When the chair is supplied with a Harmony or Symphony Back, Chest Belt loops are provided.							

13. Urine Management

ZIHW	Whizard – winged urine deflector – Has patent-pending lateral wings that are positioned under the medial aspect of the thighs to keep the Whizard in place.		\$75
ZSS	SplaShield – Molded urine deflector (must be held in place while using)		\$170
ZSC17	Splash Curtain ¹ for 17" D Bridged Front Seats, Molded and VFI (surrounds entire aperture, raises seat 1.22").		\$345
ZSCC	Custom Splash Curtain for Custom Seats – Raises seat 1.22"		\$420
ZSCSI	Slide-In Splash Curtain – Vinyl curtain, slides into commode pan brackets		\$390
ZBUDBH	Bard Urine Drainage Bag Holder 18"		\$190
¹ Bridged Front Seats only. For Left, Right, Front, and Rear Open, call for order codes and pricing.			

14. Other

ZDCLT	Directional Caster Lock – fits on one rear caster (facilitates steering)	\$195	ZPSK1	Commode Pan Spacer Kit ¹ – 2.2" lower	\$225
ZATTRR	Anti-Tippers – Rear (pr)	\$395	ZPSK2	Commode Pan Spacer Kit ¹ – 2" lower	\$245
ZATFR	Anti-Tippers – Front (pr)	\$395	ZPSK3	Commode Pan Spacer Kit ¹ – 2.75" lower	\$285
ZCUT4	4" Caster Upgrade (4x)	\$140	ZPSK4	Commode Pan Spacer Kit ¹ – 4" lower	\$285
ZCUT6	6" Caster Upgrade (4x)	\$400	ZSPKC22	CAT Self-Propel Kit with 56cm wheels – incl anti-tippers (tilt limited to 25°)	\$1,765
ZTLS	Tilt Limit Spacers – 1.22" spacers snap onto gas struts to limit tilt (each reduces tilt range approx. 5°. Can be stacked for additional tilt limiting. Quantity ____	N/C	ZSPKC24	CAT Self-Propel Kit with 61cm wheels – incl anti-tippers (tilt limited to 25°)	\$1,765

¹ Commode Pan Spacer Kits provide extra clearance for self-care with pan in place. Increases height required to clear toilet.

15. Credits and Omissions

Z905	Omit Commode Pan and Mounting Brackets	-\$20	Z911	Omit Molded Head Support Pad	-\$20
Z907	Omit Foot Supports	-\$45	Z912	Omit Complete Head Support Assembly	-\$50

Some assembly of Raz products is required.

16. Notes

Specifications	Raz-CAT (Attendant Tilt)
Weight Capacity	300 lbs
Overall Width	Seat width plus 4.5"
Overall Length (with MFX12 foot supports and headrest adjusted to middle positions)	
• at 0° tilt	36.75" (rear casters trailing), 33.75" (rear casters leading)
• at 40°	44"
Seat Depth	11" or 16" depending on seat/cushion length
Clearance Height (under seat)	
• with pan brackets	17"-21" (1" increments)
• without pan brackets	17.75"-21.75 (1" increments)
Clearance Width (to fit around toilet)	15.75" for standard 18"W seat frame. For other widths, seat frame width minus 2.25"
Arm Support Height - above molded seat (1/2" increments)	
• 15"H Adjusta-Back	5.5"-8"
• 18"H Adjusta-Back	6.5"-9"
• 18"H Reclining Adjusta-Back	6.5"-9"
Seat-to-Floor Height (with molded seat)	21"
Product Weight	55 lbs
Frame and Hardware Material	Stainless Steel (lifetime warranty against corrosion)